1

SONY

LMD-X2710MD

4K 27-inch Surgical LCD Monitor

Eil. Nr. 1.1.2.
27" jstrižainēs

High quality 4K Medical Monitor with extensive inputs to support the evolving surgical system of 4K/HD and HDR

- Wide Color Gamut complying with BT.2020
- HDR Gamma support (Hybrid Log-Gamma)
- 12G/3G/HD/SD-SDI inputs
- Variety of display modes (side by side, PinP, POP)
- Easy to clean ergonomic design

SUPER PICTURE QUALITY

4K Ultra HD Resolution

4K resolution on a 27-inch (3840 x 2160)* screen offers high picture quality.

* The 27-inch (684-mm) screen is measured diagonally.


HD/SD to 4K Upscaling

Thanks to Sony's unique upscaling technology, LMD-X2710MD provides a natural, sharp 4K view when upscaling HD/SD images to 4K.

Wide Color Gamut

The LCD panel and signal processing technology provide a wide color gamut which complies with ITU-R Recommendation BT.2020.

HDR Gamma Support

HDR (high dynamic range) enables the display of images with a wider range of brightness levels, greater contrast and also rich colors. The LMD-X2710MD offers HLG (Hybrid Log-Gamma) for supporting HDR.

SDR


HDR

*Simulated image

INTUITIVE OPERATION

Variety of Display Modes

The LMD-X2710MD can display images from multiple sources quickly from the menu in a variety of display modes; side-by-side, picture-in-picture (PIP) and picture-out-picture (POP). Also flip patterns (mirror and flip rotation) can be selected easily by pressing buttons.

Easy-to-use Control Panel with LED Lighting Navigation

Operation is simple with intuitive, easy-to-use control panel. The LED backlighting navigation is useful, especially in dark environments. In addition, three custom buttons are available to assign user's preferred functions.


Auto Signal-Detect and Fail-Safe Modes

Two modes for auto input selection are provided. Mode 1 is to accept input automatically, without any requirement to select the input. Mode 2 is for failsafe operation; the monitor can switch to a backup input automatically by detecting no signal input. This is particularly useful if input signals are accidentally interrupted.

EASY INSTALLATION

Installation-Friendly Cabling

All the connectors face downwards, allowing for easy and organized cable connection. The cables can also be pulled in both upward and downward directions, to suit your various installation needs. Single cable of 12G-SDI makes it far easier to handle.

Installation Flexibility via Direct AC Input or AC Adaptor

To suit different installation environments, two types of power supply are available: AC input and DC input using an AC adaptor. * The AC-300MD AC adaptor is sold separately

VESA Mounting

The VESA-mounting standard (100 x 100 mm) enables these devices to be used in a variety of medical installations including cart or boom-mounted articulating arms.



Compact Design and Easy Cleaning

The narrow bezel achieves a wider display area, even within the compact body. The monitor's slim and easy-to-hold design facilitates simple user adjustment of the monitor position. The edge-to-edge flat surface is helpful to easily wipe off liquids and gels from the LCD panel and control buttons, facilitating cleanliness and disinfection.

Various inputs and outputs

Various inputs and outputs are available to meet a wide variety of user needs; DP, HDMI, DVI, 12G/3G/HD/SD-SDI

- 1 Cable clamp
- 2 SERIAL REMOTE RS-232C connector (D-sub 9-pin, female)
- 3 REMOTE connector (Stereo mini jack)
- 4 SERIAL REMOTE connector (RJ-45)
- 5 DVI-D input/output connector (DVI-D)
- 6 HDMI input connector
- 7 DP1 (Display Port 1) input connector
- 8 DP2 (Display Port 2) input/output connector
- 9 12 V 2.5 A (DC output) connector
- 10 5 V 2.0 A (DC output) connector
- 11 12G/3G/HD/SD-SDI input/output connector (BNC-type)
- 12 Equipotential terminal
- 13 AC input connector
- 14 DC input connector

Specifications

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
LCD Panel	
Picture Size (Diagonal)	684 mm (26.93 inches)
Effective Picture Size	596.2 x 335.3 mm (23 1/2 x 13 1/4 inches)
Pixel pitch	0.15525 x 0.15525 mm
Resolution	3840 x 2160 pixels
Aspect Ratio	16:9
Pixel Efficiency	99.99%
Backlight	LED
Luminance (Panel Specification)	800 cd/m ² (typical)
Contrast Ratio	1000:1
Colors	Approx. 1.07 billion colors
Viewing Angle (Panel Specification)	89°/89°/89°/89° (typical) (up/down/left/right, contrast > 10%)
Gamma	1.8, 2.0, 2.2, 2.4, 2.6, DICOM, HLG

Eil.Nr. 1.1.1
Ekranu raiška 3840 x 2160 taškų

Eil.Nr. 1.1.3
Ekranu šviesumas 800 cd/m²

Input	
HDMI Input	HDMI connector (x1), HDCP1.4 correspondence
DVI-D Input	DVI-D connector (x1), TMDS single link
SDI Input	12G/3G/HD/SD-SDI connector, BNC type (x1)
Display Port	Display Port connector (x2), 1.2, 1.3, 1.4, HDCP1.3 correspondence
Serial Remote (LAN)	D-sub 9-pin (RS-232C) (x1), RJ-45 modular connector (ETHERNET) (x1)
Remote	Stereo mini jack (x1)
AC input	AC input connector (x1), 100 V to 240 V, 50/60 Hz
DC Input	DC input connector (x1), DC 26 V
Output	
DVI-D Output	DVI-D connector (x1)
SDI Output	12G/3G/HD/SD-SDI connector, BNC type (x1)
Display Port Output	Display Port connector (x1)
DC 5 V / 12 V Output	5 V Output (x1) up to 2 A, Round type pin (female) (x1) 12 V Output (x1) up to 2.5 A, Round type pin (female) (x1)

Eil.Nr. 1.1.4
Vaizdo signalo įvestys - DisplayPort, HDMI, DVI, SDI


Unit: mm (inches)

General	
Power Requirements	AC IN: 100 V - 240 V, 50/60 Hz, 1.5 A - 0.7 A DC IN: 26 V, 5.2 A (Supplied from AC adaptor)
Power Consumption	Approx. 149 W (max.)
Operating Temperature	0°C to 35°C (32°F to 95°F)
Operating Humidity	30% to 85% (no condensation allowed)
Storage/Transport Temperature	-20°C to +60°C (-4°F to +140°F)
Storage/Transport Humidity	20% to 90%
Operating/Storage/Transport Pressure	700 hPa to 1060 hPa
Dimensions (W x H x D)	658.8 x 426.8 x 80 mm (26 x 16 7/8 x 3 1/4 inches) 658.8 x 508.6 x 319 mm (26 x 20 1/8 x 12 5/8 inches) (with SU-600 optional stand)
Mass	Approx. 9.2 Kg (Approx. 20 lb 4.5 oz) (when the optional stand is not installed)
Mounting	100 x 100 mm VESA AC plug holder (x2)
Supplied Accessories	Before Using This Unit (x1) CD-ROM (including the Instructions for Use) (x1) Service Contact List (x1)
Optional Accessories	AC adaptor AC-300MD Monitor stand SU-600MD Foot Switch FS-24

These products are distributed to US and EU as medical devices. They satisfy product safety standards (e.g. IEC 60601-1). For more details, please contact your nearest Sony sales office or an authorized dealer.

Side
 When an optional stand SU-600MD is attached

Optional Accessories

AC-300MD
AC Adaptor

SU-600MD
Monitor Stand

FS-24
Foot Switch

Distributed by

MK20476V20HB22APR

©2022 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features, design, and specifications are subject to change without notice. The values for mass and dimension are approximate. Some images in this document are simulated. "SONY" is a registered trademark of Sony Corporation. "NUCLEUS", "A.I.M.E." and "BRAVIA" are trademarks of Sony Corporation. All other trademarks are the property of their respective owners. Please visit Sony's professional website or contact your Sony representative for specific models available in your region.

EVIS X1 Video System Center

CV-1500

A Unified Platform with 5 LED Spectrum Technology

1,2,3

A Unified Platform with 5 LED Spectrum Technology

By integrating the LED light source with the video processor, Olympus has developed a powerful system that is much more compact and lightweight than the predecessors*1.

Broad Compatibility

The CV-1500 can be connected to many different types of endoscopes, providing access to a wide variety of endoscopy-supporting functions.

Enhanced Observations

1,2,4,1

1,2,4,2

In addition to conventional white light and NBI (Narrow Band Imaging) and AFI (Auto Fluorescence Imaging) observation, the CV-1500 offers three other powerful enhanced observations to improve diagnostic and therapeutic capability:

1,2,4,3

- TXI (Texture & Color Enhancement Imaging) optimizes the structure, color tone and brightness of the mucosal surface.
- RDI (Red Dichromatic Imaging) improves visibility of deep blood vessels and bleeding points.
- BAI-MAC (Brightness Adjustment Imaging with MAintenance of Contrast) improves brightness in darker portions.

Intuitive, User-friendly Functions

With One-Touch Connector for quick, easy connection and no need for white balance adjustment**2, setup is simplified, with the aim of streamlining workflow and accelerating procedure time. Touch-sensitive panel facilitates intuitive operation, while convenient functions like Pre-freeze and MyCV mode ensure user-friendly working environment. Downtime is reduced thanks to the use of LED bulbs that last years without needing replacement.

*1 Combination of EVIS EXERA III/EVIS LUCERA ELITE series light source and processor *2 Olympus 1100/1200/1500 series endoscopes only

Specifications		
Power Supply	Rated voltage	100-240 V AC; Within ±10%
	Frequency	50/60 Hz; within ±3 Hz
	Rated input	600 VA
Size	Dimensions (W x H x D)	370 x 198 x 488 mm; 398 x 218 x 580 mm (maximum)
	Weight	19.4 kg
Classification (Medical Electrical Equipment)	Type of protection against electric shock	Class I
	Degree of protection against electric shock of applied part	Depend on applied part. (The degree of protection against electric shock of this product is BF type if the mounting part to be connected to this product is BF type. However CF type is not subject to combination in this product.)
	Degree or protection against explosion	The video system center should be kept away from flammable gases.
	Analog signal output	VBS composite and Y/C; simultaneous outputs possible.
Observation	Digital signal output	12G-SDI (SMPTE ST 2082), 3G-SDI (SMPTE424M), HD-SDI (SMPTE292M), SD-SDI (SMPTE259M)
	User settings	The function settings for up to 20 users can be stored.
	Color tone adjustment	Adjust the color tone of each endoscopic image for White light observation mode, NBI observation mode, and RDI observation mode. · Red adjustment : ±8 steps · Blue adjustment : ±8 steps · Chroma adjustment : ±8 steps
	Automatic gain control (AGC)	The image can be electronically amplified when the light is inadequate due to the distal end of the endoscope being too far from the object.
	Contrast	· H (High) : Darkens the dark part and brightens the bright part. · L (Low) : Brightens the dark part and darkens the bright part.
	BAI-MAC	Brightness adjustment with maintenance of contrast
	Iris	The iris modes can be switched. · Auto : The brightness is adjusted based on the brightest part of the central part and the average brightness of the periphery part. · Peak : The brightness is adjusted based on the brightest part of the endoscopic image. · Average : The brightness is adjusted based on the average brightness of the endoscopic image.
	Image enhancement settings	Fine patterns or edges in the endoscopic images can be enhanced electrically to increase the image sharpness. · Enhancement type A : Emphasizes the pattern and contour of the endoscopic image. · Enhancement type B : Emphasizes the finer parts than structure emphasis type A.
	Switching the enhancement modes	The enhancement level can be selected from 3 levels (OFF, 1, 2, and 3)
	Image size selection	The size of the endoscopic image can be selected from 2 modes. (Except SDTV)
Electric zoom	Switch between mode 1, mode 2, and mode 3.	
PIP/POP	Switch between PIP and POP.	
Aspect ratio	Switch between 16:9 and 4:3. (Except SDTV)	
Freeze	Freeze the endoscopic image.	
Pre-freeze	The image with the least blur is selected from the images captured in the set time period before freeze operation and displayed.	
Optical-digital observation	The optical-digital observation can be performed. The endoscope compatible with the optical-digital observation is required. · NBI observation : This observation mode uses the narrow band light. · RDI observation : This observation mode uses the red dichromatic lights. · AFI observation : This observation mode uses the blue light. · TXI observation : This observation mode enhances color, texture and brightness.	
Beginning and ending examination	Beginning and ending examination timing can be set interlock with the particular operation.	
Custom switch	Assign specific functions to the following buttons. · Remote switches (Up to 5) · Foot switches (Up to 2) · Keyboard custom key (Up to 4) · Touch panel custom button of basic functions screen (Up to 3) · Touch panel custom button of custom functions screen (Up to 10)	
MyCV mode	Switch setting values of multiple functions at once.	
Remote control	The following peripheral device can be controlled (specified models only). · Portable memory · Video recorder · Color video printer · Image filing system · Server	
Patient information	The following data can be displayed on the monitor. · Patient ID · Patient name · Gender · Age · Date of birth · Comment	
Documentation	Displaying the record state	The recording state of the following peripheral device can be displayed on the monitor. · Portable memory : Remaining capacity · Video recorder : Number of shots / Recording status · Color video printer : Number of shots · Image filing system : Number of shots
	Displaying the image information	The following data can be displayed on the monitor. · Image enhancement · Electric zoom ratio · Color mode · Focus · Observation mode
	Advanced registration of patient information	Up to 50 patient information can be registered. · Patient ID · Patient name · Gender · Age · Date of birth
	Recording format	Standard image quality: TIFF; Low image quality: JPEG
Memory Backup	Memorization of user settings	The settings are held in memory even after the video system center is turned OFF.
	White balance	The white balance that is once set is held in memory (only when using the compatible endoscope).

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
Endoscopy CAD System OIP-1

ENDO-AID


Welcome to the AI Future in Endoscopy

ENDO-AID

ENDO-AID is Olympus' Intelligent Platform, which can be seamlessly integrated into our latest EVIS X1 endoscopy system. It comes with ENDO-AID CADe, a high-performing computer-aided detection application for colonoscopy that uses artificial intelligence (AI) to suggest the potential presence of lesions, such as colonic polyps, malignant neoplasms and adenomas.

ENDO-AID CADe can alert the endoscopist when a potential lesion appears, surrounding it with an intuitive marker.

Get into AI and Enjoy Its Many Benefits

- By supporting the detection of potential lesions, the system aims to increase the adenoma detection rate (ADR).
- Real-time display: The system can alert abnormalities virtually in real time (no video delays).
- Physicians may be more confident during colonoscopies, thanks to additional support by AI.
- The system may contribute to the overall improvement of clinical outcomes, irrespective of operator experience.

Product Specifications

OIP-1

Compatibility	Compatible Processor	
	Compatible Scopes	Colonovideoscope (1500/1200/1100/290/190/185 series)
Size	Dimensions	W370xH85xD444 / W390xH87xD467 (Maximum)
	Weight	Approx.10 kg
Power Supply	Power Rating	Less than 275 VA, 100-240 VAC ±10%, 50/60 Hz
Connection	Memory	Front: 1x portable memory
	Function	Colon CADe
	Input/Output	In: SDI (1080p/4k) x1 Out: SDI (1080i/1080p/4k) x2
	Frame Rate	50fps, 59.94fps *2
	Compatible Recording Device	IMH-200, IMH-20, IMH-10
	Compatible Observation Mode	WLI, TXI
	Sub Image	On (large/small), off
Observation	Display Mode	Normal mode, target mode
	Remote Control	"CUSTOM" key on the keyboard, remote switch on the endoscope, remote switch on the foot switch and custom button on the touch panel of the video system center
	Classification (Medical Electrical Equipment)	Type of Protection against Electric Shock Class I

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*1 EVIS X1 VIDEO SYSTEM CENTER OLYMPUS CV-1500
 *2 Frame rate depends on region

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.

OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 23765
 www.olympus-europa.com

GIF-H185

Gastrointestinal videoscope – outstanding HDTV imaging for routine endoscopy.


GIF-H185

Main features


HDTV image quality

With the new EVIS EXERA III system, HDTV image quality delivers high-definition observation capabilities, even with the slim scope design.

NBI (Narrow Band Imaging)

NBI in EVIS EXERA III 185 series scopes provides twice the viewable distance of EVIS EXERA II 180 series scopes and offers much greater contrast between blood vessels and mucosa. The greatly improved performance of NBI opens up exciting new clinical applications and reinforces NBI's position as the standard of care for GI endoscopy.

Close focus

Close focus enables you to obtain an enlarged, close-up image simply by moving the scope tip as close as 2 mm from the mucosa.

Slim design

This scope offers an excellent balance between size and performance, with HDTV image quality in a slim 9.2 mm diameter size.

Waterproof One-touch Connector


1.3.1

A new connector design minimises the effort required for set-up prior to and in between cases. In addition, it is fully submersible and eliminates the need for a water-resistant cap and the associated risk of an expensive repair due to accidental immersion.

Specifications

Optical system	Field of view	140°	- 1.3.2
	Direction of view	Forward viewing	
	Depth of field	2-100 mm	- 1.3.3
Insertion section	Distal end outer diameter	9.2 mm	- 1.3.6
	Distal end enlarged		
Instrument channel	Insertion tube outer diameter	9.2 mm	- 1.3.5
	Working length	1030 mm	- 1.3.8
	Channel inner diameter	2.8 mm	- 1.3.7
	Minimum visible distance	3.0 mm from the distal end	
	Direction from which endotherapy accessories enter and exit the endoscopic image		

Bending section	Angulation range	Up 210°	} 1.3.4
		Down 90°	
		Right 100°	
		Left 100°	
Total length	1350 mm		
Compatible EVIS EXERA system	Video system center OLYMPUS CV-190		
	Xenon light source OLYMPUS CLV-190		

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS EUROPA HOLDING GMBH
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 (0)40 237 730, Fax: +49 (0)40 230 761
 www.olympus-europa.com

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EVIS EXERA III

CF-H185L/I

Routine colonoscopy at its best – featuring HDTV and variable stiffness.


CF-H185L/I

Main features

HDTV image quality

With the new EVIS EXERA III system, HDTV image quality enables observation of the mucosa and capillaries in much more detail.

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Ty

Variable stiffness

Variable stiffness helps to prevent the endoscope from re-looping, for example at the sigmoid colon, and also allows the stiffness of the scope to be adjusted on a case-by-case basis in order to meet the unique anatomical needs of each patient or the handling preferences of the physician.

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
Specifications

Optical system	Field of view	140°	- 1.4.2
	Direction of view	Forward viewing	
	Depth of field	2-100 mm	- 1.4.3
Insertion section	Distal end outer diameter	12.8 mm	- 1.4.4
	Distal end enlarged		
	Objective lens		1.4.8
	Insertion tube outer diameter	12.8 mm	- 1.4.5
	Working length	L: 1680 mm I: 1330 mm	1.4.9
Instrument channel	Channel inner diameter	3.7 mm	- 1.4.7
	Minimum visible distance	3.0 mm from the distal end	
	Direction from which endotherapy accessories enter and exit the endoscopic image		

Bending section	Angulation range	Up 180°
		Down 180°
		Right 160°
		Left 160°
Total length	L: 2005 mm I: 1655 mm	1.4.4
Compatible EVIS EXERA system	Video system center	OLYMPUS CV-190
	Xenon light source	OLYMPUS CLV-190

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
WM-NP3 Workstation

Features

The WM-NP3 has been designed and manufactured to enhance the user experience and add value to the Olympus Imaging Platform systems. Using a slim line profile and offering increased loading capacity while taking up less floor space and lowering running costs, the WM-NP3 will support Olympus systems without dominating the clinical environment, making it a valuable successor to WM-NP2.

Customization

Dedicated accessories are available to optimize system and procedural efficiency.

1,5,4

Electrical Safety

The WM-NP3 is supplied with a separation transformer as standard. Output from the transformer is controlled using the central on/off switch, which allows all equipment to be powered up simultaneously.

Inrush Current Control and Low-Power Standby

The WM-NP3 uses a level detector to detect mains distortion and protect the transformer from inrush current. The inrush current protection system is powered from a medically approved AC-DC power supply. To reduce the environmental impact of the WM-NP3, the low-power standby mode reduces standby power consumption giving greater efficiency and energy savings.

Easy to Move

Ergonomically designed handles and twin-wheeled castors allow the workstation to be moved into the ideal position in the clinical environment.

Cable Management

Optimized cable management has been achieved by locating the cable management modules within the workstation hoop, increasing the cable capacity and providing access to route and remove cables easily.

1,5,7


Imaging

The WM-NP3 workstation is offered with a choice of monitor arms that optimize the range of compatible monitors and enhance the procedural view with both swivel and tilt functionality. The MAJ-2216 has a weight range between 6.5 and 12 kg, and the MAJ-2217 has a weight range between 12 and 14 kg. Using a VESA fitting, both arms will accommodate a monitor up to 32".

WM-NP3 Workstation

Accessories

MAJ-2154 Consumable Storage Holder

Article number — K10030423

MAJ-2158 Suction Jar Holder

Article number — K10028139

MAJ-2159 Scope Pole Kit

Article number — K10028140

MAJ-2160 Irrigation Tube Holder

Article number — K10028141

MAJ-2165 Drawer Unit

Article number — K10030178

MAJ-2166 Sterile Water Holder

Article number — K10030036

MAJ-2167 Nurses Control Arm

Article number — K10030210

MAJ-2173 Side Shelf

Article number — K10030389

MAJ-2211 Side Handles

Article number — K10035108

MAJ-2216 LCD Monitor Arm 6.5 kg - 12 kg

Article number — K10035789

WM-NP3

Standard Sets Include:

WM-NP3 Workstation plus items marked (*)

Article Number	Part Description	Region	Basic Set	GI Standard	URO Set
		EUROPE ROW	K10035360 K10035363	K10035365 K10035368	K10035370 K10035373
	Workstation Fitting Kit		•	•	•
	Transformer (market-specific)		•	•	•
K10028141	MAJ-2160 Irrigation Tube Holder		•	•	•
K10027573	MAJ-2146 Sliding Keyboard Tray			•	•
K10028140	MAJ-2159 Scope Pole Kit			•	
K10035789	MAJ-2216 LCD Monitor Arm 6.5-12 kg			•	
K10035108	MAJ-2211 Side Handles			•	•
K10009210	IEC Lead Set		•	•	•

Optional Compatible Accessories

Article Number	Part Description
K10021041	MAJ-1639 CO ₂ Cylinder Holder-140
K10016952	MAJ-1642 IV Pole
K10021042	MAJ-1650 CO ₂ Cylinder Holder-205
K10021043	MAJ-1653 CO ₂ Cylinder Holder — Double
K10021352	MAJ-1654 Equipotential Terminal Strip
K10021791	MAJ-1657 Keyboard Arm Side-Mounted
K10021795	MAJ-1661 Side-Mounted LCD Monitor Arm
K10021797	MAJ-1663 EUS Arm-Mount Kit
K10022056	MAJ-1665 Camera Head Holder
K10027573	MAJ-2146 Sliding Keyboard Tray
K10027575	MAJ-2149 Dual Monitor Arm

Specifications

WM-NP3	
Dimensions	1400 mm (H) x 675 mm (D) x 665 mm (W)
Weight	84 kg (unladen) including fitted transformer
Maximum Load Capacity	Maximum 20 kg on top tray
	Maximum 31 kg per shelf
	Maximum 35 kg on base panel
Castors	125 mm diameter
	2x braked conductive twin-wheeled castors
	2x braked twin-wheeled castors

Voltage

The WM-NP3 is supplied with a separation transformer based on one of the following voltage options:

MAJ-2155	WM-T3 Transformer 100 V
MAJ-2156	WM-T3 Transformer 110-120 V
MAJ-2157	WM-T3 Transformer 220-240 V

Construction

Workstation	
Hoop	Extruded profiled aluminum
Shelves	Mild steel
Top Tray	Injection-molded ABS

Olympus reserves the right of errors, modification and changes of the service and/or product offerings.

OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 233765
 www.olympus-europa.com

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
SUCTION PUMP

KV-6

The KV-6 suction pump is the successor to the KV-5 pump. It is designed for endoscopic aspiration and for non-endoscopy applications. The unique flow change buttons allow a fast vacuum build-up with a single touch. It has a “whisper mode” for quiet operation and “turbo mode” for faster suction. The pump connects to an extensive fluid collection system.

Key Benefits

- **High Suction Performance**

The KV-6 suction pump has vacuum capabilities of up to **-95 kPa**. There are three pre-set speed settings (40 l/min, 50 l/min, and **60 l/min**), which enables the operator the flexibility to control the rate of vacuum build-up at the touch of a button.

- **Cleaning and Hygiene**

- Single-use collection system allows easy and safe disposal of clinical waste.
- Touch-sensitive buttons are easy to clean.

- **Compatibility**

Compatible with existing KV-5 accessories and Olympus workstations.

The KV-6 suction pump is ideal for use either mounted on a procedure trolley or as a stand-alone unit. The variable suction control allows adjustment to meet individual requirements. A high-performance pump and motor combination ensure that the unit is “whisper quiet” for any working environment.

The KV-6 may be adapted to suit individual requirements by the use of optional components, to allow a single-use or reusable suction system. The single-use liners are conveniently supplied with solidifier for safe and simple use and disposal.

Product Details

Compact Design for Convenient Use

The KV-6's compact and easy-to-use design means that it is suitable for use in office space or the operating theater.

Simple and Easy to Operate with Clearly Labelled Controls

A clear control panel ensures that the KV-6 suction pump can be used safely and efficiently by all staff.

Trolley-Mountable for Convenience

The KV-6 is ideal for use either mounted on a procedure trolley or as a stand-alone unit for maximum flexibility of use.

Use with 1.5- or 2.5-Liter Jars and Liners

The KV-6 can be used with either 1.5- or 2.5-liter jars, enabling its use in a variety of procedures.

Suction Jar Mountable on Olympus Workstations

As the suction jar can be mounted on Olympus workstations, it is compatible with the Olympus range of devices and equipment for the office or operating theater.

Overview of Accessories

Suction Tubing

14921  **K10026130**
with 2 coupling pieces (1 x 90°), 7 x 12 mm, 60 cm

14903  **K10026804**
with 2 coupling pieces, 7 x 12 mm, 60 cm

Patient Connecting Tubing

14919  **K7503432**
50 pack

Disposable Liners

14908  **K10021894**
K10021895
1.5 liter

14909  **K10021897**
K10021898
2.5 liter

Reusable Lids

14912  **K10026127**
0.25 liter

14911  **K7503760**
1 & 2 liter

Jars for Liners

14906  **K10021893**
1.5 liter

14907  **K10021896**
2.5 liter

Reusable Jars

14943  **K10026128** 0.25 liter **K7503762** 1 liter **K7503763** 2 liter

1.6.3

Filters

14945  **K10027321**
Disposable bacteria/odor filter 10 pack


14946  **K10027322**
Disposable bacteria filter 10 pack


14913  **K10026130**
K10026127

14915  **K10026128**

Clamp Holder

14901  **K10002819**
MAJ-191

14922  **K10021793**
MAJ-1659
K10001115
MAJ-190

14944  **K7503556**
WM-60 / D-60

Two or More Jars

Disposable collection system

14900  **K10021885**
MAJ-1673

14889  **K10021188**
10 Pack

14902  **K10021878**
MAJ-1672

Disposable/reusable collection system

14920  **K10021882**
MAJ-1674

Specifications

Classification (Electromedical Equipment)	Standards Compliance	Complies with EN ISO 10079-1 and EN IEC 60601-1 and UL 60601-1.		
	Electromagnetic Compatibility	This product complies with the requirements of EN IEC 60601-1-2 for emissions and immunity, and as such its operation is unlikely to be affected by, or cause interference with, equipment meeting appropriate EMC standards. As a precaution, equipment which may be sensitive to interference outside the limits specified by EN IEC 60601-1-2 should not be placed in close proximity to the KV-6.		
	Type of Protection against Electrical Shock	Classification according to EN IEC 60601-1 and UL 60601-1: Class 1 with Type CF applied part.		
	Degree of Protection against Electrical Shock	In accordance with EN IEC 60601-1 and UL 60601-1, the KV-6 is marked with the symbol to indicate the provision of an adequate degree of protection against electrical shock and that it has an applied part isolated from all other parts of the equipment.		
	Degree of Protection against Explosion	None: the Olympus KV-6 must NOT be used in the zone of risk of flammable anesthetic gases.		
	Mode of Operation	Continuous.		
Regulatory Status	European Economic Area (EEA)	CE: This mark on the product indicates compliance with Directive 93/42/EEC relating to medical devices, Class IIa.		
Dimensions and Weight	Dimensions	Height: 210 mm	Width: 305 mm	Depth: 375 mm
	Weight	9.3 kg		
Power Switch	Marking	The power switch is marked. A white LED in the center signalizes it is turned on.		
Power Requirements	Power Supply	100–240 V~	Fusing	2 x T1.6 AH
	Frequency	50/60 Hz	Power Rating	120 W
Environmental Conditions	Ambient Temperature	Operational:	+5°C to +40°C (+41°F to +104°F)	
		Storage:	–20°C to +50°C (–4°F to +122°F)	
	Relative Humidity	Operational:	30–75%	
		Storage:	20–95%	
Atmospheric Pressure	Operational/ Storage:	70.0 to 106 kPa		
Fluid Ingress	In accordance with EN IEC 60601-1 and UL 60601-1, the product is marked IP21 to indicate that it is provided with an enclosure that prevents entry of such an amount of falling liquid as might interfere with its safe and satisfactory operation when correctly positioned.			
Resistance to Chemicals	The external surfaces of the Olympus KV-6 are resistant to: 2% aqueous neutral detergent, 70% ethyl alcohol, isopropyl alcohol, water.			
Pump Specification	Nominal Vacuum:	95 kPa	Nominal Free Air Flow Rate: 50 l/min, adjustable to 40 or 60 l/min.	
	In accordance with ISO 10079-1, the unit is marked with the words "HIGH VACUUM, HIGH FLOW" to indicate it has a vacuum of at least 60 kPa in 10 seconds and a free air flow rate of 20 l/min or greater.			
Thermal Cut-Out	The temperature is observed by a sensor in the housing. If the temperature in the housing exceeds 80°C (176°F), the software stops the motor and turns on the yellow LED. It is possible to restart the motor as soon as the temperature in the housing drops below 70°C (158°F).			
Single-Use Suction Liner	To prevent fluid from being drawn up to the pump the single-use suction liner contains a mechanism to shut off the flow when the liner is overfilled. This mechanism is activated when contact is made with aspirated fluid.			
Suction Jar	The suction jars can be autoclaved up to 135°C or disinfected. The jar is impact-resistant.			
Suction Jar Lid	The suction jar lid clips and where applicable the float mechanism are autoclavable up to 135°C. The inlet connection (from the endoscope) is marked with the word PATIENT, and the outlet connection (to the pump) is marked with the word VACUUM.			
Tubing	Use only tubing specified in the KV-6 Instructions for Use, otherwise equipment malfunction may result.			
Vacuum Gauge	Range:	0 to 100 kPa	Accuracy:	±2.5%

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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080.0023
1 ltr

080.0024
2 ltr

Suction jar made of polysulfone (PSU)

- The suction jars are made of polysulfone and are suitable for all normal disinfection and sterilization processes, e.g. including autoclaving at up to 134°C for 18 minutes or 137°C for 3 minutes. Prior to autoclaving, any chemical disinfectant residues must be removed by thorough rinsing. The suction jars must not be stacked during autoclaving.
- Important: Do not use any cleaning/disinfecting agents based on phenol.
- When preparing chemical disinfectants for use, do not exceed the concentrations, soaking times or temperatures recommended by the manufacturer.
- Prior to use, always check the suction jars for cracks, brittle and flawed spots. As an additional safety test, evacuate the suction jars to maximum vacuum level several times before actual use.

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KEYMED LTD.

KeyMed House, Stock Road, Southend-on-Sea,
Essex SS2 5QH, United Kingdom
Fax: (01702) 485677, Telephone: (01702) 616333

OLYMPUS SINGAPORE PTE LTD.

491B, River Valley Road #12-01/04,
Valley Point Office Tower, Singapore 248373
Fax: 6834-2438, Telephone: 6834-0010

OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.

A6F, Ping An International Financial Center, No. 1-3,
Xinyuan South Road, Chaoyang District, Beijing, 100027 P.R.C.
Fax: (86)10-5976-1299 Telephone: (86)10-5819-9000

OLYMPUS MOSCOW LIMITED LIABILITY COMPANY

117071, Moscow, Malaya Kaluzhskaya 19,
bld. 1, fl.2, Russia
Fax: (095) 958-2277 Telephone: (095) 958-2245

OLYMPUS AUSTRALIA PTY. LTD.

3 Acacia Place, Notting Hill, VIC 3168, Australia
Fax: (03)9543-1350 Telephone: (03)9285-5400

OLYMPUS LATIN AMERICA INC.

5301 Blue Lagoon Drive, Suite 290
Miami, FL 33126-2097, U.S.A.
Fax: (305)261-4421 Telephone: (305)266-2332


OLYMPUS KOREA CO., LTD.

Olympus Tower 9F, 446, Bongeunsa-ro,
Gangnam-gu, Seoul, Korea 135-509
Fax: (02)6255-3494 Telephone: (02)6255-3210

OLYMPUS[®]

Your Vision, Our Future

OLYMPUS FLUSHING PUMP OFF-2

With its powerful but controllable **peristaltic pump** the OFP-2 will irrigate fluid down either the instrumentation or ~~auxiliary~~ water channel of Olympus Gastrointestinal EVIS ACTERA, EXERA, LUCERA, SMARTAGE and EUS endoscopes.



Eil.Nr. 1.7.1.
Pompos tipas - perstaltinis plovimo siurblys

The OFP-2 facilitates two main functions:

- The washing of gastric and colonic mucosa resulting in improved visualisation, diagnosis and therapy during endoscopic procedures
- To aid the use of trans-endoscopic ultrasound probes by providing rapid filling of the organ to be examined

Main Unit:

- Designed to prevent over pressurisation of Olympus endoscopes
- Microprocessor control ensures accurate and precise, flow rate delivery
- Ergonomically designed touch buttons for simple use
- Bright LED display for easy verification of selected settings

Practicality:

- Use of the MAJ-1606 irrigation adaptor allows either irrigation or a combination of irrigation and the passage of an EndoTherapy device through the endoscopes instrumentation channel

Safety Features:

- Pump Head automatically cuts out if opened during operation preventing accidental user injury
- Standby mode allows quick and safe exchange of water channel tubes without the need to operate the mains switch
- A 20 second automatic cut-off timer allows controlled irrigation

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

Water Container:

- Holds two litres with markings to indicate fluid levels
- A specially designed lid guides the water tubing to optimise fluid usage
- The water container and lid are autoclavable for effective sterilisation
- Washer disinfectant compatible

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Size and functionality:

- Improved functionality allows use with either the instrumentation or auxiliary water channel, including EUS ultrasound endoscopes for the first time
- Compact design for essential space saving on workstations
- Operates via either remote controlled* scope operation or foot switch

Eil.Nr.1.7.3

Pompos valdymas - valdoma pedalu

* requires the use of MAJ-920 remote control cable

Consumables:

- The MAJ-1607 and MAJ-1608 channel water tubes and the MAJ-1606 adaptor are supplied sterile for immediate use
- The MAJ-1608 auxiliary channel water tube is designed for one single days use
- Instrument channel adaptor and MAJ-1607 water tube are single patient use items

MAJ-1606
Single patient use

Supplied sterile
Instrument channel adaptor
Pack 100: Part number - K10007072
Pack 10: Part number - K10016091

MAJ-1607
Single patient use

Supplied sterile
Instrument channel water tube
Pack 50: Part number - K10001146
Pack 10: Part number - K10016136


MAJ-1608
Single days use

Supplied sterile
Auxiliary channel water tube
Pack 50: Part number - K10001147
Pack 10: Part number - K10016135

Accessories:

MAJ-1603
Pack of three

Water container two litre
Part number - K10007071

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Olympus reserves the right to alter the specification without notice.


OLYMPUS KEYMED GROUP
KeyMed House, Stock Road, Southend-on-Sea,
Essex, SS2 5DH, United Kingdom
Telephone: +44 (0)1702 616393 Fax: +44 (0)1702 465077
e-mail: info@ayntrus.co.uk www.olympus.co.uk

OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914,
OLYMPUS EUROPA HOLDING GMBH
Premises/Goods delivery Windenerstraße 14-16, 22691 Hamburg, Germany
(Leiters) Postfach 10 49 06, D-20094 Hamburg, Germany
OLYMPUS AMERICA INC.
3500 Corporate Parkway, PO Box 610, Carlsbad, PA 18034-0610, USA
OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 246373, Singapore.

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