
SONY

LMD-X2710MD


4K 27-inch Surgical LCD Monitor

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High quality 4K Medical Monitor with extensive inputs to support the evolving surgical system of 4K/HD and HDR

- Wide Color Gamut complying with BT.2020
- HDR Gamma support (Hybrid Log-Gamma)
- 12G/3G/HD/SD-SDI inputs
- Variety of display modes (side by side, PinP, POP)
- Easy to clean ergonomic design

SUPERB PICTURE QUALITY

4K Ultra HD Resolution

4K resolution on a 27-inch (3840 x 2160)* screen offers high picture quality.

* The 27-inch (684-mm) screen is measured diagonally.


HD/SD to 4K Upscaling

Thanks to Sony's unique upscaling technology, LMD-X2710MD provides a natural, sharp 4K view when upscaling HD/SD images to 4K.

Wide Color Gamut

The LCD panel and signal processing technology provide a wide color gamut which complies with ITU-R Recommendation BT.2020.

HDR Gamma Support

HDR (high dynamic range) enables the display of images with a wider range of brightness levels, greater contrast and also rich colors. The LMD-X2710MD offers HLG (Hybrid Log-Gamma) for supporting HDR.

SDR


HDR

*Simulated image

INTUITIVE OPERATION

Variety of Display Modes

The LMD-X2710MD can display images from multiple sources quickly from the menu in a variety of display modes; side-by-side, picture-in-picture (PIP) and picture-out-picture (POP). Also flip patterns (mirror and flip rotation) can be selected easily by pressing buttons.

Easy-to-use Control Panel with LED Lighting Navigation

Operation is simple with intuitive, easy-to-use control panel. The LED backlighting navigation is useful, especially in dark environments. In addition, three custom buttons are available to assign user's preferred functions.


Auto Signal-Detect and Fail-Safe Modes

Two modes for auto input selection are provided. Mode 1 is to accept input automatically, without any requirement to select the input. Mode 2 is for failsafe operation; the monitor can switch to a backup input automatically by detecting no signal input. This is particularly useful if input signals are accidentally interrupted.

EASY INSTALLATION

Installation-Friendly Cabling

All the connectors face downwards, allowing for easy and organized cable connection. The cables can also be pulled in both upward and downward directions, to suit your various installation needs. Single cable of 12G-SDI makes it far easier to handle.

Installation Flexibility via Direct AC Input or AC Adaptor

To suit different installation environments, two types of power supply are available: AC input and DC input using an AC adaptor. * The AC-300MD AC adaptor is sold separately

VESA Mounting

The VESA-mounting standard (100 x 100 mm) enables these devices to be used in a variety of medical installations including cart or boom-mounted articulating arms.



Compact Design and Easy Cleaning

The narrow bezel achieves a wider display area, even within the compact body. The monitor's slim and easy-to-hold design facilitates simple user adjustment of the monitor position. The edge-to-edge flat surface is helpful to easily wipe off liquids and gels from the LCD panel and control buttons, facilitating cleanliness and disinfection.

Various inputs and outputs

Various inputs and outputs are available to meet a wide variety of user needs; DP, HDMI, DVI, 12G/3G/HD/SD-SDI

- 1 Cable clamp
- 2 SERIAL REMOTE RS-232C connector (D-sub 9-pin, female)
- 3 REMOTE connector (Stereo mini jack)
- 4 SERIAL REMOTE connector (RJ-45)
- 5 DVI-D input/output connector (DVI-D)
- 6 HDMI input connector
- 7 DP1 (Display Port 1) input connector
- 8 DP2 (Display Port 2) input/output connector
- 9 12 V 2.5 A (DC output) connector
- 10 5 V 2.0 A (DC output) connector
- 11 12G/3G/HD/SD-SDI input/output connector (BNC-type)
- 12 Equipotential terminal
- 13 AC input connector
- 14 DC input connector

Specifications

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| LMD-X2710MD | |
|--------------------------------------|---|
| Picture Performance | |
| LCD Panel | a-Si TFT Active Matrix LCD |
| Picture Size (Diagonal) | 684 mm (26.93 inches) |
| Effective Picture Size | 596.2 x 335.3 mm (23 1/2 x 13 1/4 inches) |
| Pixel pitch | 0.15525 x 0.15525 mm |
| Resolution | 3840 x 2160 pixels |
| Aspect Ratio | 16:9 |
| Pixel Efficiency | 99.99% |
| Backlight | LED |
| Luminance (Panel Specification) | 800 cd/m ² (typical) |
| Contrast Ratio | 1000:1 |
| Colors | Approx.1.07 billion colors |
| Viewing Angle (Panel Specification) | 89°/89°/89°/89° (typical) (up/down/left/right, contrast > 10%) |
| Gamma | 1.8, 2.0, 2.2, 2.4, 2.6, DICOM, HLG |
| Input | |
| HDMI Input | HDMI connector (x1), HDCP1.4 correspondence |
| DVI-D Input | DVI-D connector (x1), TMDS single link |
| SDI Input | 12G/3G/HD/SD-SDI connector, BNC type (x1) |
| Display Port | Display Port connector (x2), 1.2, 50T, HDCP1.3 correspondence |
| Serial Remote (LAN) | D-sub 9-pin (RS-232C) (x1), RJ-45 modular connector (ETHERNET) (x1) |
| Remote | Stereo mini jack (x1) |
| AC input | AC input connector (x1), 100 V to 240 V, 50/60 Hz |
| DC Input | DC input connector (x1), DC 26 V |
| Output | |
| DVI-D Output | DVI-D connector (x1) |
| SDI Output | 12G/3G/HD/SD-SDI connector, BNC type (x1) |
| Display Port Output | Display Port connector (x1) |
| DC 5 V / 12 V Output | 5 V Output (x1) up to 2 A, Round type pin (female) (x1) 12 V Output (x1) up to 2.5 A, Round type pin (female) (x1) |
| General | |
| Power Requirements | AC IN: 100 V - 240 V, 50/60 Hz, 1.5 A - 0.7 A DC IN: 26 V, 5.2 A (Supplied from AC adaptor) |
| Power Consumption | Approx. 149 W (max.) |
| Operating Temperature | 0°C to 35°C (32°F to 95°F) |
| Operating Humidity | 30% to 85% (no condensation allowed) |
| Storage/Transport Temperature | -20°C to +60°C (-4°F to +140°F) |
| Storage/Transport Humidity | 20% to 90% |
| Operating/Storage/Transport Pressure | 700 hPa to 1060 hPa |
| Dimensions (W x H x D) | 658.8 x 426.8 x 80 mm (26 x 16 7/8 x 3 1/4 inches) 658.8 x 508.6 x 319 mm (26 x 20 1/8 x 12 5/8 inches) (with SU-600 optional stand) |
| Mass | Approx. 9.2 Kg (Approx. 20 lb 4.5 oz) (when the optional stand is not installed) |
| Mounting | 100 x 100 mm VESA AC plug holder (x2) |
| Supplied Accessories | Before Using This Unit (x1) CD-ROM (including the Instructions for Use) (x1) Service Contact List (x1) |
| Optional Accessories | AC adaptor AC-300MD Monitor stand SU-600MD Foot Switch FS-24 |



Front
When an optional stand SU-600MD is attached

Side
When an optional stand SU-600MD is attached

Eil.Nr.2.1.1
Ekranas raiška 3840 x 2160 taškų


Eil.Nr.2.1.3
Ekranas šviesumas 800 cd/m²

Eil.Nr.2.1.4
Vaizdo signalo įvestys - DisplayPort, HDMI, DVI, SDI

Unit: mm (inches)

Optional Accessories

AC-300MD
AC Adaptor

SU-600MD
Monitor Stand

FS-24
Foot Switch

These products are distributed to US and EU as medical devices. They satisfy product safety standards (e.g. IEC 60601-1). For more details, please contact your nearest Sony sales office or an authorized dealer.

Distributed by

MK20476V20HB22APR

©2022 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features, design, and specifications are subject to change without notice. The values for mass and dimension are approximate. Some images in this document are simulated. "SONY" is a registered trademark of Sony Corporation. "NUCLEUS", "A.I.M.E." and "BRAVIA" are trademarks of Sony Corporation. All other trademarks are the property of their respective owners. Please visit Sony's professional website or contact your Sony representative for specific models available in your region.

EVIS X1 Video System Center CV-1500

A Unified Platform with 5 LED Spectrum Technology

2.2.3

A Unified Platform with 5 LED Spectrum Technology

By integrating the LED light source with the video processor, Olympus has developed a powerful system that is much more compact and lightweight than the predecessors*1.

Broad Compatibility

The CV-1500 can be connected to many different types of endoscopes, providing access to a wide variety of endoscopy-supporting functions.

Enhanced Observations

2.2.4.1

2.2.4.2

In addition to conventional white light and NBI (Narrow Band Imaging) and AFI (Auto Fluorescence Imaging) observation, the CV-1500 offers three other powerful enhanced observations to improve diagnostic and therapeutic capability:

- TXI (Texture & Color Enhancement Imaging) optimizes the structure, color tone and brightness of the mucosal surface.
- RDI (Red Dichromatic Imaging) improves visibility of deep blood vessels and bleeding points.
- BAI-MAC (Brightness Adjustment Imaging with MAintenance of Contrast) improves brightness in darker portions.

2.2.4.3

T1
T2

Intuitive, User-friendly Functions

With One-Touch Connector for quick, easy connection and no need for white balance adjustment*2, setup is simplified, with the aim of streamlining workflow and accelerating procedure time. Touch-sensitive panel facilitates intuitive operation, while convenient functions like Pre-freeze and MyCV mode ensure user-friendly working environment. Downtime is reduced thanks to the use of LED bulbs that last years without needing replacement.

*1 Combination of EVIS EXERA III/EVIS LUCERA ELITE series light source and processor *2 Olympus 1100/1200/1500 series endoscopes only

| Specifications | | |
|---|---|--|
| Power Supply | Rated voltage | 100-240 V AC; Within ±10% |
| | Frequency | 50/60 Hz; within ±3 Hz |
| | Rated input | 600 VA |
| Size | Dimensions (W x H x D) | 370 x 198 x 488 mm; 398 x 218 x 580 mm (maximum) |
| | Weight | 19.4 kg |
| Classification (Medical Electrical Equipment) | Type of protection against electric shock | Class I |
| | Degree of protection against electric shock of applied part | Depend on applied part. (The degree of protection against electric shock of this product is BF type if the mounting part to be connected to this product is BF type. However CF type is not subject to combination in this product.) |
| | Degree or protection against explosion | The video system center should be kept away from flammable gases. |
| | Analog signal output | VBS composite and Y/C; simultaneous outputs possible. |
| Digital signal output | | 12G-SDI (SMPT E ST 2082), 3G-SDI (SMPT E 424M), HD-SDI (SMPT E 292M), SD-SDI (SMPT E 259M) |
| | User settings | The function settings for up to 20 users can be stored. |
| Observation | Color tone adjustment | Adjust the color tone of each endoscopic image for White light observation mode, NBI observation mode, and RDI observation mode. · Red adjustment : ±8 steps · Blue adjustment : ±8 steps · Chroma adjustment : ±8 steps |
| | Automatic gain control (AGC) | The image can be electronically amplified when the light is inadequate due to the distal end of the endoscope being too far from the object. |
| | Contrast | · H (High) : Darkens the dark part and brightens the bright part. · L (Low) : Brightens the dark part and darkens the bright part. |
| | BAI-MAC | Brightness adjustment with maintenance of contrast |
| | Iris | The iris modes can be switched. · Auto : The brightness is adjusted based on the brightest part of the central part and the average brightness of the periphery part. · Peak : The brightness is adjusted based on the brightest part of the endoscopic image. · Average : The brightness is adjusted based on the average brightness of the endoscopic image. |
| | Image enhancement settings | Fine patterns or edges in the endoscopic images can be enhanced electrically to increase the image sharpness. · Enhancement type A : Emphasizes the pattern and contour of the endoscopic image. · Enhancement type B : Emphasizes the finer parts than structure emphasis type A. |
| | Switching the enhancement modes | The enhancement level can be selected from 3 levels (OFF, 1, 2, and 3) |
| | Image size selection | The size of the endoscopic image can be selected from 2 modes. (Except SDTV) |
| | Electric zoom | Switch between mode 1, mode 2, and mode 3. |
| | PIP/POP | Switch between PIP and POP. |
| Documentation | Aspect ratio | Switch between 16:9 and 4:3. (Except SDTV) |
| | Freeze | Freeze the endoscopic image. |
| | Pre-freeze | The image with the least blur is selected from the images captured in the set time period before freeze operation and displayed. |
| | Optical-digital observation | The optical-digital observation can be performed. The endoscope compatible with the optical-digital observation is required. · NBI observation : This observation mode uses the narrow band light. · RDI observation : This observation mode uses the red dichromatic lights. · AFI observation : This observation mode uses the blue light. · TXI observation : This observation mode enhances color, texture and brightness. |
| | Beginning and ending examination | Beginning and ending examination timing can be set interlock with the particular operation. |
| | Custom switch | Assign specific functions to the following buttons. · Remote switches (Up to 5) · Foot switches (Up to 2) · Keyboard custom key (Up to 4) · Touch panel custom button of basic functions screen (Up to 3) · Touch panel custom button of custom functions screen (Up to 10) |
| | MyCV mode | Switch setting values of multiple functions at once. |
| | Remote control | The following peripheral device can be controlled (specified models only). · Portable memory · Video recorder · Color video printer · Image filling system · Server |
| | Patient information | The following data can be displayed on the monitor. · Patient ID · Patient name · Gender · Age · Date of birth · Comment |
| | Displaying the record state | The recording state of the following peripheral device can be displayed on the monitor. · Portable memory : Remaining capacity · Video recorder : Number of shots / Recording status · Color video printer : Number of shots · Image filling system : Number of shots |
| Memory Backup | Displaying the image information | The following data can be displayed on the monitor. · Image enhancement · Electric zoom ratio · Color mode · Focus · Observation mode |
| | Advanced registration of patient information | Up to 50 patient information can be registered. · Patient ID · Patient name · Gender · Age · Date of birth |
| | Recording format | Standard image quality: TIFF; Low image quality: JPEG |
| Memory Backup | Memorization of user settings | The settings are held in memory even after the video system center is turned OFF. |
| | White balance | The white balance that is once set is held in memory (only when using the compatible endoscope). |

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


Endoscopy CAD System OIP-1

ENDO-AID


Welcome to the AI Future in Endoscopy

ENDO-AID

ENDO-AID is Olympus' Intelligent Platform, which can be seamlessly integrated into our latest EVIS X1 endoscopy system. It comes with ENDO-AID CADe, a high-performing computer-aided detection application for colonoscopy that uses artificial intelligence (AI) to suggest the potential presence of lesions, such as colonic polyps, malignant neoplasms and adenomas.

ENDO-AID CADe can alert the endoscopist when a lesion appears, surrounding it with an intuitive marker.

Get into AI and Enjoy Its Many Benefits

- By supporting the detection of potential lesions, the system aims to increase the adenoma detection rate (ADR).
- Real-time display: The system can alert abnormalities virtually in real time (no video delays).
- Physicians may be more confident during colonoscopies, thanks to additional support by AI.
- The system may contribute to the overall improvement of clinical outcomes, irrespective of operator experience.

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Product Specifications

OIP-1

| | | |
|--|--|--|
| Compatibility | Compatible Processor | CV-1500 |
| | Compatible Scopes | Colonovideoscope (1500/1200/1100/290/190/185 series) |
| Size | Dimensions | W370×H85×D444 / W390×H87×D467 (Maximum) |
| | Weight | Approx.10 kg |
| Power Supply | Power Rating | Less than 275 VA, 100-240 VAC ±10%, 50/60 Hz |
| Connection | Memory | Front: 1× portable memory |
| | Function | Colon CADe |
| Observation | Input/Output | In: SDI (1080p/4k) ×1 Out: SDI (1080i/1080p/4k) ×2 |
| | Frame Rate | 50fps, 59.94fps *2 |
| | Compatible Recording Device | IMH-200, IMH-20, IMH-10 |
| | Compatible Observation Mode | WLI, TXI |
| | Sub Image | On (large/small), off |
| | Display Mode | Normal mode, target mode |
| | Remote Control | "CUSTOM" key on the keyboard, remote switch on the endoscope, remote switch on the foot switch and custom button on the touch panel of the video system center |
| Classification (Medical Electrical Equipment) | Type of Protection against Electric Shock | Class I |

*1 EVIS X1 VIDEO SYSTEM CENTER OLYMPUS CV-1500
*2 Frame rate depends on region

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com

GIF-H185

Gastrointestinal videoscope – outstanding HDTV imaging for routine endoscopy.

GIF-H185

Main features

HDTV image quality

With the new EVIS EXERA III system, HDTV image quality delivers high-definition observation capabilities, even with the slim scope design.

NBI (Narrow Band Imaging)

NBI in EVIS EXERA III 185 series scopes provides twice the viewable distance of EVIS EXERA II 180 series scopes and offers much greater contrast between blood vessels and mucosa. The greatly improved performance of NBI opens up exciting new clinical applications and reinforces NBI's position as the standard of care for GI endoscopy.

Close focus

Close focus enables you to obtain an enlarged, close-up image simply by moving the scope tip as close as 2 mm from the mucosa.



Slim design

This scope offers an excellent balance between size and performance, with HDTV image quality in a slim 9.2 mm diameter size.

Waterproof One-touch Connector



2.3.1

A new connector design minimises the effort required for set-up prior to and in between cases. In addition, it is fully submersible and eliminates the need for a water-resistant cap and the associated risk of an expensive repair due to accidental immersion.

Specifications

| | | |
|--------------------------|---------------------------|-----------------|
| Optical system | Field of view | 140° |
| | Direction of view | Forward viewing |
| | Depth of field | 2-100 mm |
| Insertion section | Distal end outer diameter | 9.2 mm |
| | Distal end enlarged | |

| | | |
|---------------------------|--|----------------------------|
| Insertion tube | Insertion tube outer diameter | 9.2 mm |
| | Working length | 1030 mm |
| Instrument channel | Channel inner diameter | 2.8 mm |
| | Minimum visible distance | 3.0 mm from the distal end |
| | Direction from which endotherapy accessories enter and exit the endoscopic image | |

| | | |
|--------------------------|------------------------------------|------------|
| Bending section | Angulation range | Up 210° |
| | | Down 90° |
| | | Right 100° |
| | | Left 100° |
| Total length | 1350 mm | |
| Compatible | Video system center OLYMPUS CV-190 | |
| EVIS EXERA system | Xenon light source OLYMPUS CLV-190 | |

2.3.4

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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EVIS EXERA III

CF-H185L/I

Routine colonoscopy at its best – featuring HDTV and variable stiffness.


CF-H185L/I

Main features

HDTV image quality

With the new EVIS EXERA III system, HDTV image quality enables observation of the mucosa and capillaries in much more detail.

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Variable stiffness

T₅ Variable stiffness helps to prevent the endoscope from re-looping, for example at the sigmoid colon, and also allows the stiffness of the scope to be adjusted on a case-by-case basis in order to meet the unique anatomical needs of each patient or the handling preferences of the physician.

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
Specifications

| | | | |
|---------------------------|--|----------------------------|---------|
| Optical system | Field of view | 140° | - 2.4.2 |
| | Direction of view | Forward viewing | |
| | Depth of field | 2-100 mm | - 2.4.3 |
| Insertion section | Distal end outer diameter | 12.8 mm | - 2.4.6 |
| | Distal end enlarged | | |
| | | | |
| Instrument channel | Channel inner diameter | 3.7 mm | - 2.4.7 |
| | Minimum visible distance | 3.0 mm from the distal end | |
| | Direction from which endotherapy accessories enter and exit the endoscopic image | | |

| | | | |
|--------------------------|------------------------------------|------------|-------|
| Bending section | Angulation range | Up 180° | 2.4.4 |
| | | Down 180° | |
| | | Right 160° | |
| | | Left 160° | |
| Total length | L: 2005 mm I: 1655 mm | | |
| Compatible | Video system center OLYMPUS CV-190 | | |
| EVIS EXERA system | Xenon light source OLYMPUS CLV-190 | | |

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

WM-NP3 Workstation

Features

The WM-NP3 has been designed and manufactured to enhance the user experience and add value to the Olympus Imaging Platform systems. Using a slim line profile and offering increased loading capacity while taking up less floor space and lowering running costs, the WM-NP3 will support Olympus systems without dominating the clinical environment, making it a valuable successor to WM-NP2.

Customization

Dedicated accessories are available to optimize system and procedural efficiency.

Electrical Safety

The WM-NP3 is supplied with a separation transformer as standard. Output from the transformer is controlled using the central on/off switch, which allows all equipment to be powered up simultaneously.

Inrush Current Control and Low-Power Standby

The WM-NP3 uses a level detector to detect mains distortion and protect the transformer from inrush current. The inrush current protection system is powered from a medically approved AC-DC power supply. To reduce the environmental impact of the WM-NP3, the low-power standby mode reduces standby power consumption giving greater efficiency and energy savings.

Easy to Move


Ergonomically designed handles and twin-wheeled castors allow the workstation to be moved into the ideal position in the clinical environment.

Cable Management

Optimized cable management has been achieved by locating the cable management modules within the workstation hoop, increasing the cable capacity and providing access to route and remove cables easily.

Imaging

The WM-NP3 workstation is offered with a choice of monitor arms that optimize the range of compatible monitors and enhance the procedural view with both swivel and tilt functionality. The MAJ-2216 has a weight range between 6.5 and 12 kg, and the MAJ-2217 has a weight range between 12 and 14 kg. Using a VESA fitting, both arms will accommodate a monitor up to 32".

WM-NP3 Workstation

Accessories

MAJ-2154

Consumable Storage Holder

Article number — K10030423

MAJ-2160

Irrigation Tube Holder

Article number — K10028141

MAJ-2158

Suction Jar Holder

Article number — K10028139

MAJ-2165

Drawer Unit

Article number — K10030178

MAJ-2167

Nurses Control Arm

Article number — K10030210

MAJ-2159

Scope Pole Kit

Article number — K10028140

MAJ-2166

Sterile Water Holder

Article number — K10030036

MAJ-2211

Side Handles

Article number — K10035108

MAJ-2173

Side Shelf

Article number — K10030389

MAJ-2216

LCD Monitor Arm 6.5 kg - 12 kg

Article number — K10035789

WM-NP3

Standard Sets Include:

| WM-NP3 Workstation plus items marked (•) | | Region | Basic Set | GI Standard | URO Set |
|--|------------------------------------|---------------|------------------------|------------------------|------------------------|
| Article Number | Part Description | EUROPE ROW | K10035360 K10035363 | K10035365 K10035368 | K10035370 K10035373 |
| | Workstation Fitting Kit | | • | • | • |
| | Transformer (market-specific) | | • | • | • |
| K10028141 | MAJ-2160 Irrigation Tube Holder | | • | • | • |
| K10027573 | MAJ-2146 Sliding Keyboard Tray | | | • | • |
| K10028140 | MAJ-2159 Scope Pole Kit | | | • | |
| K10035789 | MAJ-2216 LCD Monitor Arm 6.5-12 kg | | | • | |
| K10035108 | MAJ-2211 Side Handles | | | • | • |
| K10009210 | IEC Lead Set | | • | • | • |

Optional Compatible Accessories

| Article Number | Part Description |
|----------------|--|
| K10021041 | MAJ-1639 CO ₂ Cylinder Holder-140 |
| K10016952 | MAJ-1642 IV Pole |
| K10021042 | MAJ-1650 CO ₂ Cylinder Holder-205 |
| K10021043 | MAJ-1653 CO ₂ Cylinder Holder -- Double |
| K10021352 | MAJ-1654 Equipotential Terminal Strip |
| K10021791 | MAJ-1657 Keyboard Arm Side-Mounted |
| K10021795 | MAJ-1661 Side-Mounted LCD Monitor Arm |
| K10021797 | MAJ-1663 EUS Arm-Mount Kit |
| K10022056 | MAJ-1665 Camera Head Holder |
| K10027573 | MAJ-2146 Sliding Keyboard Tray |
| K10027575 | MAJ-2149 Dual Monitor Arm |

Specifications

| WM-NP3 | |
|-----------------------|--|
| Dimensions | 1400 mm (H) × 675 mm (D) × 665 mm (W) |
| Weight | 84 kg (unladen) including fitted transformer |
| Maximum Load Capacity | Maximum 20 kg on top tray |
| | Maximum 31 kg per shelf |
| | Maximum 35 kg on base panel |
| Castors | 125 mm diameter |
| | 2× braked conductive twin-wheeled castors |
| | 2× braked twin-wheeled castors |

Voltage

The WM-NP3 is supplied with a separation transformer based on one of the following voltage options:

| | |
|----------|-----------------------------|
| MAJ-2155 | WM-T3 Transformer 100 V |
| MAJ-2156 | WM-T3 Transformer 110-120 V |
| MAJ-2157 | WM-T3 Transformer 220-240 V |

Construction

| Workstation | |
|-------------|----------------------------|
| Hoop | Extruded profiled aluminum |
| Shelves | Mild steel |
| Top Tray | Injection-molded ABS |

Olympus reserves the right of errors, modification and changes of the service and/or product offerings.

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
OLYMPUS FLUSHING PUMP OFF-2

With its powerful but controllable **peristaltic pump** the OFP-2 will irrigate fluid down either the instrumentation or a ~~an~~ auxiliary water channel of Olympus Gastrointestinal EVIS ACTERA, EXERA, LUCERA, SMARTAGE and EUS endoscopes.



Eil.Nr. 2.6.1.
Pompos tipas - peristaltinis plovimo siurblys

The OFP-2 facilitates two main fu

- The washing of gastric and colonic mucosa resulting in improved visualisation, diagnosis and therapy during endoscopic procedures
- To aid the use of trans-endoscopic ultrasound probes by providing rapid filling of the organ to be examined

Main Unit:

- Designed to prevent over pressurisation of Olympus endoscopes
- Microprocessor control ensures accurate and precise, flow rate delivery
- Ergonomically designed touch buttons for simple use
- Bright LED display for easy verification of selected settings

Practicality:

- Use of the MAJ-1606 irrigation adaptor allows either irrigation or a combination of irrigation and the passage of an EndoTherapy device through the endoscopes instrumentation channel

Safety Features:

- Pump Head automatically cuts out if opened during operation preventing accidental user injury
- Standby mode allows quick and safe exchange of water channel tubes without the need to operate the mains switch
- A 20 second automatic cut-off timer allows controlled irrigation

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

Water Container:

- Holds two litres with markings to indicate fluid levels
- A specially designed lid guides the water tubing to optimise fluid usage
- The water container and lid are autoclavable for effective sterilisation
- Washer disinfectant compatible

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Size and functionality:

- Improved functionality allows use with either the instrumentation or auxiliary water channel, including EUS ultrasound endoscopes for the first time
- Compact design for essential space saving on workstations
- Operates via either remote controlled* scope operation or foot switch

Eil.Nr.2.6.3

Pompos valdymas - valdoma pedalu

* requires the use of MAJ-920 remote control cable

Consumables:

- The MAJ-1607 and MAJ-1608 channel water tubes and the MAJ-1606 adaptor are supplied sterile for immediate use
- The MAJ-1608 auxiliary channel water tube is designed for one single days use
- Instrument channel adaptor and MAJ-1607 water tube are single patient use items

MAJ-1606
Single patient use

Supplied sterile
Instrument channel adaptor
Pack 100: Part number - K10007072
Pack 10: Part number - K10016091

MAJ-1607
Single patient use

Supplied sterile
Instrument channel water tube
Pack 50: Part number - K10001146
Pack 10: Part number - K10016136

MAJ-1608
Single days use

Supplied sterile
Auxiliary channel water tube
Pack 50: Part number - K10001147
Pack 10: Part number - K10016135

Accessories:

MAJ-1603
Pack of three

Water container two litre
Part number - K10007071

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Olympus reserves the right to alter the specification without notice.

OLYMPUS

OLYMPUS KEYMED GROUP
KeyMed House, Stock Road, Southend-on-Sea,
Essex, SS2 5GH, United Kingdom
Telephone: +44 (0)1702 616933 Fax: +44 (0)1702 405577
e-mail: info@olympus.co.uk www.olympus.co.uk


OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Municipal 2-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0294
OLYMPUS EUROPA HOLDING GMBH
(Prämiestraße/Goods delivery) Wendenstraße 14-16, 20087 Hamburg, Germany
(Leibniz) Postfach 10 49 98, D-20024 Hamburg, Germany
OLYMPUS AMERICA INC.
3500 Corporate Parkway, PO Box 610, Canal Valley, PA 19034-0610, USA
OLYMPUS SINGAPORE PTE LTD.
401B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 244573, Singapore.

ISO 9001 Registered
Certificate No. 78 0000
Quality Management System
ISO 14001 Registered
Certificate No. 78 0000
Environmental Management System

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