Helix Locking Tool

Model A001HELX

User's Manual

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2023 Abbott. All Rights Reserved.

Helix Locking Tool


User's Manual

Intended Use

The Helix Locking Tool (Model A001HELX) is designed to aid with the Tendril™ STS (Model 2088TC) and is a non-implantable accessory tool designed to aid with lead fixation by providing control over the extension and retraction of the helix. This can be an alternative to the clip-on tool. The Helix Locking Tool is placed on the lead by inserting the tool on the connector pin. The longitudinal slot feature on the Helix Locking Tool allows for placement on the lead with a stylet already in place.

Refer to the appropriate lead user's manual for a complete list of lead-specific indications, contraindications, warnings, precautions, potential adverse effects and instructions for use.

Figure 1. Helix Locking Tool

Clinical Benefits

The Helix Locking Tool is primarily used to facilitate the implantation of implantable transvenous pacing leads. Thus, this accessory is contributing to the same clinical benefits associated with the pacing leads. When this device performs as intended, it provides the associated clinical benefits to the target patients. When the leads and the device system perform as intended, they provide the appropriate pacing therapy to manage their heart conditions.

Securing the Tip


Securing the Tip with the Helix Locking Tool

As an alternative to the clip-on tool, the Helix Locking Tool may be used to extend or retract the helix.

1. Check the Use-By date on the packaging. Do not use tool if the Use-By date has passed.
2. Aseptically remove tool from packaging.
Do not use tool if:
 - Dropped from packaging
 - Sterility is breached
3. Slide Helix Locking Tool onto the connector end of lead. The tool should be fully inserted as shown in illustration.

Figure 2. Helix Locking Tool Fully Inserted

1. Fully inserted to connector pin


4. After the fixation site has been selected, hold the lead body stationary in one hand and turn the Helix Locking Tool in the direction of the arrow marked "EXTEND HELIX" (clockwise).

NOTE: See Technical Specifications for the approximate number of turns required. On the fluoroscopic image, the helix will be extended beyond the marker ring.

5. PSA cables may be attached in the zones depicted below.


Figure 3. PSA Cable Attachment Zones

1. Attachment Zones

6. If necessary, retract the helix by holding the lead body stationary in one hand and turning the Helix Locking Tool in the direction marked "RETRACT HELIX" (counterclockwise).

Figure 4. Helix Locking Tool retraction

7. After completion, slide tool backward to remove.

Atrial Lead Placement

1. Using a straight stylet, introduce the lead into the atrium so that it rests on the floor of the atrial chamber.
2. Replace the straight stylet with a J-shaped stylet, or withdraw the existing stylet, bend it into a soft J-shape, and reinsert the curved stylet into the lead.
3. As the stylet approaches the electrode tip, introduce more lead to ensure that the tip remains in the atrium as the lead takes its "J" shape.
4. Retract the lead as necessary to ensure that the electrode tip slides into the atrial appendage. Observe the fluoroscopy monitor to verify that the "J" is straightening.
5. When the lead is past the appendage and in the chamber, feed more lead into the heart so that it regains its "J" shape.
6. Take a firm grip on the stylet, then introduce more of the lead so that the electrode tip goes as far as possible into the atrium. On fluoroscopy, the electrode tip will "tilt over" as proof that it can go no further.
7. With the Helix Locking Tool, extend the helix so that the lead is fixed to the atrial wall.
8. Retract the entire stylet from the lead with a smooth and steady motion.
9. Check that the lead is properly anchored by introducing more of it into the heart until the loop that forms either lies on the bottom of the atrium, or is about to enter the inferior vena cava or the right ventricle. Retract any excess lead until it acquires the correct "J" shape.

Figure 5. Atrial lead placement


10. Ask the patient to breathe deeply and check that the lead keeps its “J” shape.
11. Ask the patient to cough to ensure that the electrode is securely anchored.

Ventricular Lead Placement

1. Advance the lead into the atrium.
2. Pull the stylet back a few centimeters to reduce the risk of the lead damaging the valves or penetrating the heart muscle when it continues down into the ventricle.
3. Continue to advance the lead. When the tip reaches the apex, retract the stylet an additional ten centimeters or more.
4. With the Helix Locking Tool, extend the helix to fix the lead tip to the ventricular wall. If the tip is correctly secured, the lead will be felt to jerk slightly.
5. Remove the stylet completely. Adjust the lead so that it lies in the desired position in the ventricle.

Figure 6. Ventricular lead placement

Device Storage

- Store the accessory at room temperature of 15°C and 30°C (59°F and 86°F). During transportation and handling, the accessory can be exposed to temperature excursions between -20°C and 55°C (-4°F and 131°F).
- Before the package is opened, inspect it visually for any damage that may have compromised sterility.

Sterilization Instructions

- The package content has been sterilized with ethylene oxide before shipment. This accessory is for single use only and is not intended to be resterilized.
- If the sterile package has been compromised, contact Abbott Medical.

Technical Support

Abbott Medical maintains 24-hour phone lines for technical questions and support:

- 1 818 362 6822
- 1 800 722 3774 (toll-free within North America)
- + 46 8 474 4147 (Sweden)
- + 61 2 9936 1200 (Australia)
- medical.abbott/manuals

For additional assistance, call your local Abbott Medical representative.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical.abbott/manuals.

Symbol	Description
 medical.abbott/manuals	Follow instructions for use on this website
	Medical Device
	Unique Device Identifier
	Storage temperature. Temperature value is indicated adjacent to the symbol.
	Double sterile barrier

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