Thermo Scientific B·R·A·H·M·S KRYPTOR compact PLUS

Vartotojo vadovas | Version 6

Data: 20.11.2013

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thermoscientific.com/brahms

© 2014 Thermo Fisher Scientific Inc. Všechna práva vyhrazena. Windows is a registered trademark of Microsoft Corporation in the United States and other countries. KRYPTOR and TRACE are registered trademarks of CIS bio international, licensed for use by B-R-A-H-M-S, a part of Thermo Fisher Scientific. CEZANNE SAS hold exclusive rights in the B-R-A-H-M-S KRYPTOR compact software. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. All data regarding specifications, terms and pricing correspond to the existing knowledge at the time of the printing. We are not responsible for any errors, misprints or changes. Reprint, also in parts, solely with prior written consent of B-R-A-H-M-S GmbH.

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Contents

Labelling	4
Safety Instructions	5
Introduction	7
Instrument Components	8
Description of the System Status Window	10
Main Menus	11
The Tool Palette	12
Work Surface Colour Codes	13
Sample Carousel	13
Reagent Area	14
Reaction Area	15
Dilution Plate	15
B·R·A·H·M·S KRYPTOR compact SOLUTIONS 1 to 4	16
Fluidic System	16
Routine Work	17
Start-up	17
Maintenance	17
Reagent registration	18
Reconstitution of reagents	18
Perform a calibration	18
Run controls	19
Run patient samples with barcode	19
End of day	19
Analytes with Pre-incubation	20
Important User Information	22
Barcodes Specifications	24
Troubleshooting Guide	25
Result Window Messages	25
Degraded Mode	27

Labelling

Symbol	Explanation
	Laser, danger to eyes Location: on fluidic hood, on barcode readers
	Electrical risk Location: high voltage power supply, inside electronic boxes
	Biohazard Location: waste bottle plug, wash station, waste pump, clot detection board
	Thermal risk Location: reaction area entrance
	Electrostaticsensitive part Location: inside electronic boxes
	User hazard Location: on sample and cooling trays
	Caution - Class 3B invisible laser radiation when open. Avoid exposure to the beam Location: on reader head shielding
	Caution - Class 3B invisible laser radiation when open and interlocks defeated. Avoid exposure to the beam Location: on side panel, access to the laser
	Laser light Do not stare into the beam. Class 2 laser product Location: on fluidic hood
	Serial number of Reading Module
	Serial number of Pipetting Module
	Identification of the machine
	Alternative current (AC)

Symbol	Explanation
	Humidity range
	Temperature range
	Medical device conforms to IVD Directive 98/79/EC
	Consult instructions for use
	Medical device is CE certified to conform with the IVD Directive
	Date of manufacturing
	Name and address of the Manufacturer
	Protective nominal value (Rating of electrical fuse)
	Electric and electronic equip- ments have to be selectively collected under the manufacturer responsibility B·R·A·H·M·S GmbH
106172	Product reference number of the assembly

Labelling of B·R·A·H·M·S KRYPTOR compact SOLUTIONS

	Xn		SGH 06 harmful
	Xi		SGH 07 irritant

Safety Instructions

The safety instructions of the User Manual allow the user to avoid injury to persons, material damages and environmental contamination.

It is mandatory for users of Thermo Scientific™ B·R·A·H·M·S™ KRYPTOR™ compact PLUS to pay particular attention to SAFETY INSTRUCTIONS written in the B·R·A·H·M·S KRYPTOR compact PLUS Manual. Installation of the B·R·A·H·M·S KRYPTOR compact PLUS can only be performed by a properly trained service engineer. At the time of installation all performance specifications will be verified. Any attempt to install, repair or modify the instrument by unauthorized personnel is not allowed.

This *in vitro* device must only be used by qualified and trained personnel according to GMP guide lines. Local health and safety regulations must also be taken into account.

Reading and interpretation of results must be done by a qualified user.

B·R·A·H·M·S KRYPTOR compact PLUS must only be used with materials, equipment and accessories specified in the B·R·A·H·M·S KRYPTOR compact PLUS Manual.

B·R·A·H·M·S KRYPTOR compact PLUS analyser is a Class 2 laser product.

Caution: Use of controls or adjustments or performing any procedures other than those specified herein may result in hazardous radiation exposure.

Precaution for installation

Indoor use only

Altitude	Up to 2000 m
Temperature range	18...30 °C
Humidity	20-85% (non condensing)
Corrosion	Protection against 1N HCl, 1N NaOH, decontamination solution, bleached reagent
Electrical mains	100-240 V~, mains supply voltage fluctuations up to ± 10% of the nominal voltage 50-60 Hz Power 465 VA Transient overvoltage: category II
Placement	Clearance distance: Behind: 5 cm Side-left: 10 cm Side-right: 20 cm Disconnecting device should be accessible.
Transport conditions	-20...70 °C
Long term storage	0...50 °C

Input and output connections

Caution:

On the serial communication ports, connect only a RS-232 link (Very low safety voltage).

On the USB port, connect only the specific USB cable provided (Very low safety voltage).

On Waste bottle cap (liquid level and mechanical sensor), connect only the specific cable provided (Very low safety voltage).

Electromagnetic compatibility

Changes or modifications not expressly approved by B·R·A·H·M·S GmbH could void the user's authority to operate the equipment.

B·R·A·H·M·S KRYPTOR compact PLUS is compliant with class B product requirements as defined in IEC 61326-2-6 standard.

B·R·A·H·M·S KRYPTOR compact PLUS complies with the emission and immunity requirements described in EN 61326-2-6. The electromagnetic environment should be evaluated prior to operation of the device. Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g. unshielded intentional RF sources), as these may interfere with the proper operation.

Electrical security

Do not connect B·R·A·H·M·S KRYPTOR compact PLUS to a power supply before ensuring that the voltage setting is correct.

The analyzer can be used with a power supply (mains) voltage of 100-240 V~ (50-60 Hz). Verify the voltage of the local power supply (mains) to be used. Always plug the analyzer into a grounded outlet. Operating technicians and maintenance personnel are urged to follow sound electrical safety practices at all times.

Although all metal parts of the analyzer are at ground potential (zero voltage), they should never be touched with one hand while touching a plumbing fixture, radiator, AC-operated device or other grounded object with the other.

Before opening the analyzer, remove the power cable from the power outlet. Do not replace components or attempt any repairs with the analyzer switched on.

Do not operate the analyzer in an atmosphere containing explosive gases since components of the analyzer could generate sparks.

Avoid spilling fluid on or into the analyzer at any time. All spills should be wiped up promptly.

If the equipment is used in a manner not specified by B·R·A·H·M·S GmbH, the protection provided by the equipment may be impaired.

Biohazardous

When working with human serum, controls and calibrators all accessible parts of the analyzer must be considered as biohazardous. The pipette tip, sample cassettes, reagent cassettes, carousel drip pan, and analyzer deck should be routinely disinfected.

Concerning all areas exposed to patient potential infectious material or to user contact follow the cleaning requirements and use 5% Hypochloride solution for decontamination. (See B·R·A·H·M·S KRYPTOR compact PLUS Manual, chapter maintenance.)

It is strongly recommended to wear gloves and a special coat.

Waste

The reagents and the waste must be disposed of as potential infectious laboratory waste according to local regulations.

Dilution plate, reaction plate

Before disposing of the dilution plate as well as the reaction plate stick the adhesive cover film with biohazard sign on the plate. The adhesive covers are supplied with the plates.

B·R·A·H·M·S KRYPTOR compact SOLUTIONS 1 to 4

For detailed information read the safety instructions for use inserts.

Safety data sheets are available at International Product Support on request.

Laser

Concerning the sample carousel barcode reader (Class 2 laser) and the manual barcode reader do not look at the laser. Avoid direct eye exposure.

Class 2 lasers are limited to a maximum output power of 1 milliwatt and the beam must have a wavelength between 400 and 700 nm. A person receiving an eye exposure from a Class 2 laser beam, either accidentally or as a result of someone else's deliberate action (misuse) will be protected from injury by their own natural aversion response. This is a natural involuntary response that causes the individual to blink and avert their head thereby terminating the eye exposure. Repeated, deliberate exposure to the laser beam may not be safe. The laser inside the instrument (Class 3B laser) emits an invisible radiation whose characteristics are the following:

Laser Type	LTB	SRS
Beam Deviation (mrad)	3*3	5*8
Pulse Length (ns)	2.5	<3.5
Max Peak Power (kW)	100	45
Repetition Rate (Hz)	20	20

Class 3B lasers may have sufficient power to cause an eye injury, both from the direct beam and from reflections. The higher the output power of the device the greater the risk of injury. Class 3B lasers are therefore considered hazardous to the eye.

The laser is only accessible by opening the instrument cover. The instrument cover (side access panel) is equipped with a safety interlock, which stops laser beam immediately.

For further information see latest version of B·R·A·H·M·S KRYPTOR compact PLUS Manual and latest version of B·R·A·H·M·S KRYPTOR instructions for use for reagents and consumables.

Introduction

Intended use

B·R·A·H·M·S KRYPTOR compact PLUS is a fully automated system for *in vitro* diagnostic.

It is used for measurement of patient samples in random access mode routinely.

B·R·A·H·M·S KRYPTOR compact PLUS is a closed system and can only operate utilising special reagents offered by B·R·A·H·M·S GmbH. The system is based on TRACE technology (Time-Resolved Amplified Cryptate Emission). The system and the reagents are under continuous development by B·R·A·H·M·S GmbH.

TRACE – the unique measuring principle of B·R·A·H·M·S KRYPTOR compact PLUS

The measurement principle of B·R·A·H·M·S KRYPTOR/ KRYPTOR compact is based on TRACE technology, which measures the signal that is emitted from an immuno complex with time delay.

The basis of the TRACE technology is non-radiative energy transfer from a donor to an acceptor.

The proximity of donor and acceptor when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the cryptate and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.



Precise measurement of analyte concentration

When the sample is excited with a nitrogen laser at 337 nm, the donor emits a long-life fluorescent signal in the milli-second range at 620 nm, while the acceptor generates a short-life signal in the nanosecond-range at 665 nm or 707 nm depending on the type of acceptor. When the two components are bound in an immuno-complex, both the signal amplification and the prolonged life span of the acceptor signal occur at 665 nm or 707 nm, so that it can be measured over μ -seconds. This long-life signal is proportional to the concentration of the analyte to be measured.

Reliable prevention of interference

Non-specific signals are eliminated by the internally calculated ratio of the intensities at these wavelengths (665/620 or 707/620).

The signal generated by the cryptate at 620 nm serves as an internal reference and is measured simultaneously with the long-life acceptor signal at 665 nm or 707 nm. Interfering influences, *e.g.* from turbid sera, are automatically corrected.

Instrument Components

This instrument has two parts: pipetting module and reading module.

Pipetting Module

This part comprises a pipetting unit **1**, fluidic system **2** and carousel **3**.

1 Pipetting unit

The pipetting unit consists of a heated tip sample probe and a wash station. The heated tip uses a capacitive liquid level detection system to aspirate all components required to perform the biological analysis. This includes the sample, diluent and reagents. After aspirating the sample and reagent components into the tip, the tip is heated to reaction temperature and acts as a pre-reaction incubator. This ensures that the fluids are at reaction temperature before dispensing in to the reaction area. At the end of dispense, the tip is washed to avoid carry-over contamination between samples.

Caution, any fluidic hood opening stops tip motion.


2 Fluidic system

The fluidic system includes 3 bottles of 5 liters. One for the PBS wash liquid (BUFFER; order code for the buffer: B·R·A·H·M·S KRYPTOR BUFFER) used to wash the tubing system after pipetting. A second bottle with distilled or demineralized water used as liquid system (DISTILLED WATER) by priming before pipetting and to keep the tubing line in water during instrument stand by or switch off. A third bottle is used to collect liquid waste (WASTE). Only the waste bottle contains a float sensor on the cap to manage the liquid level and a mechanical sensor for bottle

presence. Levels of buffer and distilled water are managed by float switches within the main body of the instrument. When the bottle becomes red on the VDU screen, it is possible to refill the bottles without interruption to the current run. There is no fluid sensor on these two bottles.

Distilled or demineralized water tank with floating switch

PBS Tank with floating switch

3 Carousel

The carousel is divided into 5 positions. Positions 1, 2 and 3 are hybrid locations which allow to put reagent or sample trays. Positions 4 and 5 are dedicated

only for sample trays. When the hood is open, it is possible to turn the carousel by hand to simplify access to each position.

Reagents tray

This accessory contains up to 4 reagent kits. The cooling system (2...8 °C) is operational only when the tray is installed on the carousel and the instrument is turned on. An infrared system is used to detect reagent segments during movement of the carousel. Barcoded identification of reagent kits is automatically scanned when the carousel rotates. Caution: Put the kit in right order, barcode must be visible within the perspex window (refer to photo). Inventory for all reagent kits and reagent segment temperature is fully managed by the software with use of colour codes.

Tray for samples, dilution plates and SOLUTIONS 1-4

This tray has space for 16 sample tubes which can be any combination of primary or secondary tubes at a height between 60 mm and 120 mm and a diameter between 11 mm and 17 mm. Calibrators, controls and micro-cups in a specific metallic adaptor can also be installed. This tray also has space to install SOLUTIONS 1-4 or dilution plates. Sample trays are detected on board by barcode identification placed on the top and the middle of the segment.

SOLUTIONS 1, 2, 3 and 4 (order code B·R·A·H·M·S KRYPTOR compact SOLUTION 1-4)

The SOLUTIONS 1 and 2 are used for reconstitution of freeze dried reagents. The SOLUTIONS 3 and 4 are needed for different washing steps between pipetting reagents and samples. Their use depends on the analytes. Do not top-up the bottles. Each solution bottle is managed on the instrument by the software through a barcode placed on the top of the bottle. SOLUTIONS 1 and 2 are mandatory to launch reconstitution of reagent kits. SOLUTIONS 3 and 4 are mandatory to launch pipetting sequence for samples.

DILUTION PLATES (order code B·R·A·H·M·S KRYPTOR compact DILCUP)

The dilution plate, which consists of 24 wells, is used by the system to reduce a high concentration or make a recovery test by using a diluent or a specific reagent. A barcode placed on the convex side of the plate, allows the software to monitor the availability of wells and lifetime on board the instrument.

4 Reading Module


Power on or off of the complete instrument is done with the button O/I placed on right side of reading module. This module consists of the reaction area and optical reading system to measure the signal emitted by the immuno complex. The reaction plate contains 96 wells (order code B·R·A·H·M·S KRYPTOR compact REACT) and its life onboard the instrument is monitored by the software through a barcode placed on each reaction plate and transparent damper closing. The reaction plate is heated at $37\text{ °C} \pm 0.5\text{ °C}$.

5 Handheld barcode reader

All relevant information concerning calibrators, controls, reagent lots and reaction plate barcodes are registered on B·R·A·H·M·S KRYPTOR compact PLUS by the barcode reader as well as reaction plate barcode.

Description of the System Status Window

① ACCESS TO MAIN MENUS

② B·R·A·H·M·S KRYPTOR compact PLUS TOOL BAR

③ WINDOWS COMMAND BAR

④ WINDOWS START BUTTON

⑤ WINDOWS TASK BAR

⑥ STATUS BAR

⑦ B·R·A·H·M·S KRYPTOR compact PLUS WORK SURFACE

① Waste

② Distilled or demineralized water: Liquid system

③ Buffer PBS: Phosphate Buffer Saline: Wash liquid

④ B·R·A·H·M·S KRYPTOR compact Wash Solutions

⑤ Dilution plates

⑥ Sample tray

⑦ Reagents tray

⑧ Reaction plate

Main Menus

These menus are available by clicking on: 

The main menus bar appears

To close it, click again on: 

SYSTEM	DATA	INSTRUMENT	REGISTRATION	ADMINISTRATION	HELP
* Logon	* Worklist	* Rescan Carousel	* Calibrator/Standard	* Panels	* User Manual
* Shutdown	* Results	* Start Processing	* Control	* Analytes	* About
* Maintenance	* QC Functions	* Pause Processing	* Reagent Lot	* Combined Analytes	
* Maintenance Log	* Work Analysis	* Query All		* Reflex Testing	
* Session Log		* Prime		* User Accounts	
* Service Diagnostics		* Reconstitute Kits		* Instrument Params	
* Printer Setup		* Reaction Plate		* LIS Interface	
* Close		* Sample Carousel		* Preferences	
		* Reagents			

The Tool Palette

ICON	COMMAND	FUNCTION	MENU OPTION	DESCRIPTION
	Scan Carousel	F1	Instrument Rescan Carousel	Identification of samples without starting a run
	Start Processing	F2	Instrument Start Processing	This icon disappears after click on it to give access to Pause processing.
	Pause Processing	F3	Instrument Pause Processing	This icon appears after Processing has started. Addition of more samples to the carousel. If you push the blue button on the front of the instrument, at the end of pause, carousel hood is open automatically. If you use this icon, at the end of pause it is needed to click on the blue button to open the hood.
	Manage Worklist	F5	Data Worklist	
	Validate Results	F6	Data Results	
	QC Functions	F7	Data QC Functions	
	Reagent Lot Manager	F8	Registration Reagent Lot	
	Calibrator Lot Manager		Registration Calibrator/Standard	
	Prime Instrument	F9	Instrument Prime	
	Reconstitute Reagents	F10	Instrument Reconstitute Kits	
	Advanced Menu			To close it, click again on the icon.
	User Manual		Help	
	Load/unload the Reaction plate	F11	Instrument Reaction Plate	Allow to load reaction plate by plate ID manual entry. Allow to unload reaction plate.
	System Logon procedures		System Logon	
	System Shutdown procedures		System Shutdown	

Work Surface Colour Codes


Sample Carousel

At the end of carousel scan, all samples with a barcode and tests in the worklist will be shown in green, blue, yellow or red. If the system is connected to LIS, there are different colours for the sample status (see Administration Menu/ LIS interface in B·R·A·H·M·S KRYPTOR compact PLUS Manual).

An information window is available by clicking with the left button of the mouse on a selected sample.

In case of sample tubes without barcodes, it is possible to identify a location in the sample tray by manual entry through the **Carousel Status** window, or the worklist.

Double click on sample tube to open **Carousel Status** window.

COLOUR	EXPLANATION
	Black: no sample present
	Dark green: the sample is ready to be tested
	Dark green with a white point: sample during pipetting process
	Dark green with a red circle: STAT sample
	Blue: sample to be removed (all the tests of this sample have processed successfully including the out of range detection)
	Yellow: sample with barcode without test in worklist or presence on carousel of same ID on different sample tube
	Red: problem with sample processing (e.g. reagent related, consumables related; clot, insufficient sample volume)
	Red with a white point: sample during pipetting process with at least one test with problem

Reagent Area

Each reagent cassette has its own identification. An information area appears by clicking (left button of the mouse) on reagent tray screen image, this shows cassette ID and temperature in degrees of the kit.



A **colour code** on the circle portion exists as a quick visual check of the right temperature in the reagent area or a reagent cassette's fan fault.

COLOUR	TEMPERATURE BAR ON THE REAGENT AREA
	Green: 2.0...8.0 °C
	Blue: < 2.0 °C
	Red: > 8.0 °C or communication failed
	Yellow: one or more reagent cassette's fan is out of order

All reagent kits detected by barcode are shown on screen in the appropriate reagent cassette. A more detailed information Window is available by clicking (left button of the mouse) on a reagent kit.

The colour of the background indicates the availability of the kits and the colour of the rectangle indicates the number of remaining tests.

COLOUR	EXPLANATION
	White background: Kit OK (calibration and reconstitution OK) Red rectangle: Number of tests remaining: ≤ 5
	White background: Kit OK (calibration and reconstitution OK) Yellow rectangle: Number of tests remaining: > 5 – ≤ 10
	White background: Kit OK (calibration and reconstitution OK) Green rectangle: Number of tests remaining: > 10
	Yellow background: Kit or calibration expires tomorrow Green rectangle: Number of tests remaining: > 10
	Red background: Kit to reconstitute, to calibrate or expired kit or no test remaining

Double clicking anywhere in the reagent area opens the **Reagent Status** window. This window enables you to display information about each reagent kit installed in the reagent area of the instrument and to request calibration for each reagent kit.

Reaction Area

Each reaction plate is barcode labelled to allow the loading of the plate. The reaction plate has to be loaded to activate its preheating and to allow the dispensing of samples.

Click on the plate to view the exact number of tests dispensed and completed. When the reaction plate is full, there is no automatic unloading, it is necessary to request by clicking on reaction plate. A reaction plate not used but unloaded could be loaded again during 2 hours after its first loading.

There is no automatic shutdown (machine remains in ready stage). As long as the reaction plate is on board, it is ready to be used, until 7 days are passed or plate is full. After 7 days, the plate is automatically unloaded at the next end of day (first ready or initialized stage after midnight).

Time validity status of reaction plate is available when clicking with the left mouse button on the reaction plate on the screen. If there is an "xy position error", the plate is automatically unloaded, and cannot be used any more.

COLOUR	EXPLANATION
	Green: Proportion of tests running and not completed
	Blue: Proportion of tests completed

	Green: 36.2...38.0 °C
	Blue: < 36.2 °C
	Red: > 38.2 °C

Pipetting sequence does not start if the reaction area temperature is below to 35.7 °C or over 38.5 °C.

Dilution Plate

Each dilution plate is barcode labelled to manage the wells used and the shelf life on board (365 days). Check dilution plate cleanness during daily maintenance.

COLOUR	EXPLANATION
	White well = available well / coloured well = used well
	Green: More than 10 wells available
	Yellow: From 5 to 10 wells available
	Red: Less than 5 wells available or expired

B·R·A·H·M·S KRYPTOR compact SOLUTIONS 1 to 4

For SOLUTIONS 1 and 2, the dead volume has been defined to finish a reconstitution started even if red level is detected; it is normal to have around 12 mL in SOLUTION 1 and approximately 8 mL in SOLUTION 2 when red status appears on the interface screen. For SOLUTIONS 3 and 4 the dead volume is around 4 mL, defined to finish wash step.

COLOUR	EXPLANATION
	Green: OK
	Yellow: Take care to change the bottle soon / solution expires tomorrow
	Red: Bottle to be replaced = dead volume reached or solution expired

Fluidic System

An intermediate tank system enables you to refill B·R·A·H·M·S KRYPTOR BUFFER or distilled/demineralized water bottle when red status is detected without stopping the current run by use of dead volume in intermediate tank system. If the dead volume is used then waste red level is detected and the instrument stops at the end of current run or prime. A prime is mandatory after refilling a bottle.

COLOUR	EXPLANATION
	Green: OK
	Red: Waste bottle full (void), buffer or water bottle empty (filling), intermediate waste collector full (Call Hot Line)

Routine Work

Start-up

Fill the appropriate bottles with distilled or demineralized water and buffer PBS, empty the waste bottle. Ensure B·R·A·H·M·S KRYPTOR compact SOLUTIONS 1 to 4 are on board. If B·R·A·H·M·S KRYPTOR compact SOLUTIONS 1 to 4 are closed, open the caps. **Put caps into specific locations and inscribe the caps with number of solution.**

Check the fluidic and carousel hoods are closed.

After a power off


- Check if carousel hood is closed.
- Switch on the instrument with O/I (ON/OFF) button found on right side of the instrument.
- Switch on the XPC, monitor, printer.
- Click on Windows Start button and choose

icon to launch the B·R·A·H·M·S KRYPTOR compact software.

2 programs will be started:

- 1) XPC program is the user interface.
- 2) Check if XIPC program is active by icon presence on clock bar.

After a Shutdown start here

- Check if carousel hood is closed.
- Click on  icon.
- Select **Name** and **Password** (e.g. user: Admin and password: Admin) in the **System LogOn** window.
- An automatic scan is launched after motor initialization. System is ready to use at the end of scan.
- An automatic scan is done every two minutes when carousel hood is closed to manage reagent tray temperature.
- Do **maintenances** as requested.
- Place a new reaction plate, close the transparent damper and load it by scanning the barcode.
- Request carousel hood opening by pressing on the blue button in front of carousel.
- If there is no reagent cassette, load one or more on specific positions.
- Check status of dilution plate and B·R·A·H·M·S KRYPTOR compact SOLUTIONS 1 to 4 on interface, according to necessity, change or load new consumables.
- Load kits on reagent tray when temperature status is green.
- Ensure the reagent lid is well closed.

After an automatic change of day

The automatic change of day procedure is launched every day, on the first ready or initialized stage after midnight. This will perform automatically all the initializations (init pipettor + init reader).

Nevertheless, the user will have a logon to perform manual maintenances and database maintenance.

User is warned that he has to logon with the message "KRYPTOR compact – Change of day made – Maintenance needed" in the title of the window.

User can go on using the analyser, without logon, in case of use by night team, who wants that maintenance is performed by "morning team".

Maintenance

Items on the maintenance screens only need to be performed if they are marked ***Expired***.

(!) These procedures will be performed automatically when the system has been started.

Daily maintenance

- **Initialize pipettor (!)**. To be done only if problems occur during run on pipettor module.
- **Initialize reader (!)**. To be done only if problems occur during run on reader module.
- **Prime Liquid Handling System (!)**. Check visually the absence of leaks and bubbles during the priming procedure.
- **Check and clean reader head window** by using a cotton-swab first with water then with alcohol.
- **Check and clean condensation in Reagent cassettes** with absorbent paper.
- **LIS end of day**. This action enables to purge LIS files and avoid the occurrence of data transmission problems.
- **Check dilution plate cleanness**. Look at dilution plate already on board if there is no dirty particles in dilution well.

Weekly maintenance

- **Check for liquid leaks**. Open fluidic hood to check the tubing system (tubing connectors, syringes). When fluidic hood is closed, pipettor initialization is automatically done.
- **Tip path cleaning with decontamination solution** to clean up serum splashes with a paper towel slightly wet with water, then repeat with alcohol or specific product to decontaminate.
- **Backup DataBases**. Backup of all databases older than 3 days.

- **Backup Log Files.** Backup of log of run and maintenance log files.
- **Remove messages from box office.** This action enables the purge of the box office and avoid the occurrence of data transmission problems.
- **Cleaning of water bottle.** Empty and replace distilled/demineralized water if it is more than seven days old.
- **Clean carousel pan.** Clean the tray under the sample carousel with a paper towel slightly wet with water, then repeat with alcohol or specific product to decontaminate.

Monthly maintenance

- **Bottles decontamination.** Empty the bottles and pour in 1 liter of 5% sodium hypochlorite solution. Fit a bottle cap and swirl the liquid inside the bottle so that it comes into contact with all the internal surfaces. Allow to stand for 15 minutes, then empty and rinse the bottle with water.
- **Secure tip cleaning.** Tip parking over night or if instrument is not in use for a longer time. B·R·A·H·M·S KRYPTOR compact SOLUTION 3 and 4 must be available.
- **Automatic check dot.** Tip coming on dot point. If adjustment is correct, tip coming back to wash cup. In case of adjustment problem, user should open the fluidic hood and check if tip can be slightly bent to enter the hole.

Reagent registration

For the following working steps installation of valid version of K-DISK ANA is needed.

Necessary when using new lots only:

Reagent

- Click on **Reagent Lot Manager:** 
- Click on **Register.**
- Scan the barcode sheet.
- Click on **OK** and confirm the registration with **Yes.**

Calibrator

- Click on **Calibrator/Standard Manager:** 
- Click on **Register.**
- Scan the barcode sheet.
- Click on **OK** and confirm the registration with **Yes.**

Control

- Click on: 
- Click on: 
- Click on **Register new Controls.**
- Scan the barcode sheet.
- Click on **OK.**

Reconstitution of reagents

- Click on  to start reconstitution process.
- All freeze dried reagent kits placed in the instrument will be reconstituted.
- In the status bar the progression of the reconstitution in % is visible.

Perform a calibration

Calibrator

- Double click on the kit you want to calibrate.
- Select a calibrator lot from the drop-down arrow list.
- Click on **Add to worklist, OK, Close.**
- Request a carousel hood opening by pressing on blue button in front of instrument.
- Place the calibrator tube on the sample carousel making sure the barcode fills the slot.
- Close the carousel hood.
- Start the instrument: 

Results

- Open results list: 
- Select the calibrator result and click on **Validate curve.**
- The software indicates if the new calibration curve is acceptable or not:
if Yes click on: 
if No click on: 

Run controls

Control

Adding a Control through Worklist

- Select **Add control** in the **Data, Worklist:** 
- Select a control from the list, press **Add to worklist, Close, Close.**

Adding a Control to be tested on specific reagent kit

- Double click on the kit which should be used to measure the controls.
- Select **Control lot** from the drop-down arrow list.
- Select **Low volume** and/or **Stat status** for this control.
- Click on **Add to worklist, OK, Close.**

- Request a carousel hood opening by pressing on blue button in front of instrument.
- Place the control tube on the sample carousel making sure the barcode fills the slot.
- Close the carousel hood.
- Start the instrument: 

Results

- Open results list: 
- If result is in "needs resolution" status related to concentration out of 2SD range (shown in red) or with flag, check the information and validate, reject or relaunch the test. Accepting the control automatically sends the result to the QC program.
- Select the B·R·A·H·M·S KRYPTOR compact Quality control program  to follow control life.

Run patient samples with barcode

Patient Samples with barcode

- Select **Add sample** from **Data, Worklist:** 
- Enter the Sample ID. Cassette ID and Tube position are mandatory in case of sample without barcode.
- Select tests and dilutions from the displayed list, **Add.**
- Click on **Stat** to pipette a sample in priority.
- Click on **Low volume sample vial** if sample is in microcup.
- Select **Save sample** and continue with required samples, **Close.**
- Request a carousel hood opening by pressing on blue button in front of carousel part.
- Place sample tubes on the carousel making sure the barcode fills the slot.
- Close the carousel hood.
- Start the instrument: 

Results

- Open results list: 
- Select a test.
- It is possible to **Cancel** tests in detecting, or counting step or in pending, only in "Ready/..." status.
- Results needing resolution can be: **Accepted, Re-run** or **Rejected.**
- Select the results then **Print Status** or **Print Report** to print results.

End of day

Shutdown (stand by modus)

- Check inventory of reagents on board for tomorrow's workload.
- Check the level of waste, distilled or demineralized water and B·R·A·H·M·S KRYPTOR BUFFER.
- Check if the B·R·A·H·M·S KRYPTOR compact SOLUTION 3 and 4 is on board.
- Select  and confirm with **Yes.**
- Unload the reaction plate and confirm the popup window with **OK.**
- Close the transparent damper and confirm the popup window with **OK.**

To power OFF the instrument

- Remove dilution plate(s), samples, and reagents.
- Close carousel hood.
- Select  and confirm with **Yes.**
- Unload the reaction plate and confirm the popup window with **OK.**
- Close the transparent damper and confirm the popup window with **OK.**
- Instrument is washing, wait for the completion of steps (SOLUTION 3 and water washing) in case of instrument power off.
- Close B·R·A·H·M·S KRYPTOR compact software: 
- Close Windows by clicking on **Start** and **Turn off.**
- Switch OFF the PC, screen, printer and instrument.

Analytes with Pre-incubation

The assay of an analyte with a pre-incubation is composed of two phases:

Pre-incubation or Phase I: the system aspirates one of the 2 conjugates and the antigen (sample). The system dispenses the mixture into the reaction well and the pre-incubation starts.

Phase II: the system aspirates the second conjugate and dispenses it into the reaction well. At this moment the Out of Range detection starts.


This leads to a specific operation of the interface to avoid losing tests in Phase I distributed into the reaction area.

The final test result will be achieved through Phase I and Phase II using the same reagent kit.


To alert the user that a kit is in use for Phase I and is required for Phase II, a padlock will appear on the user interface at the edge of the reagent kit as shown below. The kit should not be removed while the padlock is displayed.

At the same time, a flashing red square will appear in the top left corner of the screen.

The message behind the red square indicates the estimated time in minutes and seconds remaining of Phase I. There is also a number in parentheses to show the number of tests in Phase I.


By reading this message the user can estimate the time available to open and close the hood.


The padlock is removed when the reagent kit has no tests remaining in Phase I, however the RED square can remain active if other tests in Phase I are underway.

If the buffer bottle, water bottle, one of the solution bottles or the waste bottle is detected in red status during Phase I or Phase II, pipetting will stop and a popup window on interface informs to replace the solution. In Phase II the carousel hood will be opened by the system software and a beeping sound will alert the user.

A user request to open the carousel hood during Phase I is subject to confirmation by the user via a message appearing on the screen.


Clicking on the "Yes" button allows the hood to be opened with a beeping sound to alert the user and a new text message next to the red flashing square indicating the time available for the user to close the hood. Clicking on "No" led to the abandonment of the application to open the carousel hood.


When closing the carousel hood, the text indicates the estimated time in minutes and seconds before the distribution of Phase II.

When closing the hood, a carousel scan is performed. If a reagent kit containing a pre-incubation test was not found on board of the instrument, the user is notified via a message displayed on the interface:


If the kit is not returned before the start of Phase II, any tests in progress with this kit will be declared failed. At the results list, the user can view the status of a test.

During Phase I, tests appear in blue with the message pre-incubating:

The start time of a test is given after the start of Phase II. If the test is failed during Phase I, the test is indicated in red with an alarm: "System error: Pre-incubation has failed."


Important User Information


Detailed information about consumables and patient samples

By clicking and holding the left mouse button an information window appears.


1. Buffer


2. SOLUTIONS 1-4


3. Dilution plates


4. Reagent cassette


5. Patient samples



6. Reaction plate

7. Reagent


When a consumable is loaded by manual entry, it is marked by a symbol  (not reaction plate) and on session log.

If during a scan, a barcode is detected in place of manual entry ID, scan reading has priority, ID consumable detected by scan erases manual entry ID. This action is indicated in session log.

Interface drawing of tray positioning may differ from actual tray position. For manual entry make sure that placements of reagent kits are done correctly as illustrated below.

Reagent kit placement

Check the kit's correct position in reagent tray with kit notch towards to the transparent window. A kit mispositioning could impact data results.

Sample ID barcode

It is forbidden to use a sample barcode starting with "R" or "_". In these cases, problems occur on interface.

Manual entry of consumables and reagents

In case of barcode reading problem, it is possible to load barcode information by manual entry with the keyboard. Click with the right mouse button on the corresponding area. A specific window for manual entry of the barcode data will open.

Caution

Manual identification of consumables or sample tubes is under user responsibility. Entering incorrect data may lead to wrong result. This option is to be used only in case of barcode reading fault while waiting for repair of problem cause.

Pipetting / Reading sequence

During pipetting/reading sequence it is forbidden to close interface to avoid to lose samples results or to have bad behavior of results list management.

Clot management

A clot management is available on the system during aspiration step on sample and dilution. The user is warned of clot presence by two kinds of information.


Case 1:

Clot has been released into the tube, tip has been washed before continuing the pipetting sequence. Sample tube with clot is in red status with tests not pipetted in cancel status, and test with clot problem with flag "System error: clot detected".

User intervention should be done on tube (clot removal or sample tube change) before to launching new tests.

Case 2:

Clot could not be released, tip is blocked. This popup appears to inform that pipetting sequence is stopped with carousel hood open.

Tip cleaning after fatal clot detection

- Open fluidic hood, place tip on Zhome position.
- Clean point of tip.
- Unscrew tubing on tip, liquid must flow out into wash bowl.
- Disconnect tip electronic board, Caution to not bend the tip point.
- Pass through a fish line (dia 0.45 mm) into the tip from the small end towards the tip board.
- Clean by blowing with a neutral gas through the tube from the small end.
- Connect tip electronic board on specific location.
- Close fluidic hood.
- Close carousel hood.
- Check dot point through monthly maintenance "Automatic check dot".
- Launch a prime.

Help

Some help is available, in three languages (English, French and German), through  button or Help/User Manual menu.

Barcodes Specifications

Available symbologies on KRYPTOR instruments:


- Code 128
- Code 39
- Codabar
- Code 2/5 interleaved
- Code UPC/EAN

Other specifications:

- Minimum resolution: 0.21 mm
(size of narrowest bar (black or white bar))
- Minimum ratio = 2.5
(ratio of narrow bars to wide bars)
- Maximum ratio = 3.0


Silence section (blank section before and after the code): at least 10 times the resolution (a silent section before and after the code is mandatory for a good reading otherwise the code cannot be recognized).

Barcode positioning: We recommend to stick the label 5 mm above the bottom of the tube (if the code is long and the label is stuck too high the barcode reader may have some difficulties to read the bars located at the top, specially if the resolution is low (thin bars)). At least a few bars from the barcode must be in the window formed by the tube holder: the window will determine the barcode position within the cassette. The code must not be stuck entirely above the tube holder because there is a risk of misidentification in this case (a tube in position 2 may be seen in position 3 for example).

- A:** Barcode position is acceptable but it's not the best position.
- B:** Best barcode position.
- C:** This barcode position is forbidden. At least a few bars must be in the tube holder window.

Tube positioning: The tube must be straight and the bars must be visible throughout the window width.


- A:** Good position
- B-C:** Bad position, the barcode reader may not be able to scan the code

Paper quality: A paper having a too high ink absorption will give a blurred result and will decrease the amount of good readings (this problem will have a higher impact on barcodes having a low resolution).

Paper brightness: The barcode readers are very sensitive to the light reflection. Even though the barcode reader is mechanically adjusted in order to limit the light reflection a too glossy paper may affect the reading performances of the codes.

Printing quality: The printer quality has an influence on the reading performances. The best printers for barcodes are thermal printers (recommended for low resolution barcodes) but most of the time a laser printer is sufficient.

- A: Printing quality is good.**
Good contrast. Good resolution (accurate printing).
- B: Printing quality is low.**
The printer resolution is not good enough. The contrast is not sufficient.

Troubleshooting Guide

Result Window Messages


FLAG NO.	RESULT WINDOW MESSAGE	WHAT TO DO
0	"System Warning: Low Laser Power"	<ul style="list-style-type: none"> • Call the Hot Line
1	"System Error: Reaction Plate position"	<ul style="list-style-type: none"> • Request a reader initialization in Daily Maintenance or restart the system completely. • Unload the reaction plate, check that the reaction plate was correctly loaded and load a new reaction plate. • If the problem persists, call the Hot Line.
2	"System Warning: Reagent Cooler out of Range"	<ul style="list-style-type: none"> • Adjust the room temperature within the specification (between 18 - 30°C). • Check if a reagent cassette fan is out of order. • If the problem persists, call the Hot Line.
3	"System Warning: Incubator Out of Range"	<ul style="list-style-type: none"> • Stop the B·R·A·H·M·S KRYPTOR compact PLUS and wait a few minutes before switching on the instrument. • Adjust the room temperature within the specification (between 18 - 30°C). • If the problem persists, call the Hot Line.
4	"System Error: Missed Reading"	<ul style="list-style-type: none"> • Call the Hot Line.
6	"Error: Insufficient Sample Volume"	<p>Check if the trouble occurs on one or several samples tubes, cup or calibrator vials.</p> <p>If it concerns one tube:</p> <ul style="list-style-type: none"> • check the liquid level in the sample cup, the sample tube or the calibrator vial. • check the tube positioning in the carousel. <p>If it concerns several tubes:</p> <ul style="list-style-type: none"> • Call the Hot Line.
7	"System Error: Insufficient Reagent Volume"	<ul style="list-style-type: none"> • Check the liquid level in reagent bottles. • Call the Hot Line.
9	"System Warning: Late T0/TM"	<ul style="list-style-type: none"> • If the problem occurs more than once in a series, shut down, and then restart the instrument after 1 or 2 minutes. If the problem persists, call the Hot Line. • If the problem occurs sporadically while a significant work list is loaded on B·R·A·H·M·S KRYPTOR compact PLUS, the problem is not relevant. It is then advisable to relaunch the sample displaying this message.
11	"System Warning: Late TE"	<ul style="list-style-type: none"> • If the problem occurs more than once in a series, shut down, and then restart the instrument after 1 or 2 minutes. If the problem persists, launch an intervention. If the problem persists, call the Hot Line. • If the problem occurs sporadically while a significant work list is loaded on B·R·A·H·M·S KRYPTOR compact PLUS, the problem is not relevant. It is then advisable to relaunch the sample displaying this message.
12	"System Error: Pipetting"	<ul style="list-style-type: none"> • Check obstacle in the path (Cap on sample tube...). <p>Remark: In that case, the carousel hood is open and the run is interrupted.</p> <ul style="list-style-type: none"> • Check the liquid level in vial. <p>Remark: In that case, test is in pending status. If following popup appears <Error: tip level sense baseline. Contact technical support>, the carousel hood is open and the run is interrupted.</p>
13	"System Error: Clot Detected"	<ul style="list-style-type: none"> • Clean the tip (fish line, Secure tip cleaning from the monthly maintenance). • Check the sample and relaunch the test. • If the trouble persists, call the Hot Line.
14	"System Error: Pipetting in Sample area"	<ul style="list-style-type: none"> • Check obstacle in the path (Cap...). • Check the sample tube positioning.
15	"System Error: Pipetting in Dilution area"	<ul style="list-style-type: none"> • Check the dilution plate positioning.
16	"System Error: Pipetting in Reagent area"	<ul style="list-style-type: none"> • Check obstacle in the path (reagent area lid, aluminium foil...). • Check the reagent kit positioning.
17	"System Error: Pipetting in Wash area"	<ul style="list-style-type: none"> • Check if the SOLUTION bottles are present, open and if sufficient liquid remains. • Check if there is any obstruction on the path. • If the problem persists, call the Hot Line.
18	"System Error: Pipetting in Reaction area"	<ul style="list-style-type: none"> • Check if there is any obstruction in the path. • If the problem persists, call the Hot Line.

FLAG NO.	RESULT WINDOW MESSAGE	WHAT TO DO
19	"System Warning: Heated Tip Out of Range"	<ul style="list-style-type: none"> • Prime the system to remove bubbles. • If the problem persists, call the Hot Line.
20	"System Error: Interrupted"	<ul style="list-style-type: none"> • Reboot the system. • If the problem persists, call the Hot Line.
21	"System Warning: Not Performed"	<ul style="list-style-type: none"> • Start testing or delete it from the work list.
24	"Check Results"	<ul style="list-style-type: none"> • The software is working properly. If the "Flag report" option in the combined analytes is checked, any flag on the individual tests A1 and A2 generates a "Check results" message on the combined result which makes it only acceptable manually.
25	"Math Error"	<ul style="list-style-type: none"> • Call the Hot Line.
26	"Reflex Test Launched"	<ul style="list-style-type: none"> • It is not a problem, this is an information only.
29	"Clot detected or insufficient sample volume"	<p>Message occurs when a clot is detected too close to sample's zmax. The current test will be flagged in red.</p> <ul style="list-style-type: none"> • Clean the tip (fish line, Secure tip cleaning from the monthly maintenance). • Check the sample and relaunch the test. • If the trouble persists, call the Hot Line.
30	"Data Error: Unknown Error"	<ul style="list-style-type: none"> • Call the Hot Line.
31	"Data Error: Ratio"	<ul style="list-style-type: none"> • Call the Hot Line.
32	"Data Error: Response"	<ul style="list-style-type: none"> • Call the Hot Line.
33	"Data Warning: Abnormal"	<ul style="list-style-type: none"> • Prime the system to remove bubbles. • Check the samples for bubbles. • Check sample and reaction plate for dust. • Check if there are any bubbles in the reaction plate. • If the problem persists, call the Hot Line.
34	"Data Warning: Out of Range"	<ul style="list-style-type: none"> • Automatic dilution is started depending on the Preferences configuration in the Administration menu.
35	"Data Warning: Detection Limit"	<ul style="list-style-type: none"> • Result lower than the lowest detectable concentration.
36	"Data Error: Above max. Range"	<ul style="list-style-type: none"> • Relaunch the test manually with a small dilution factor (1/2 or 1/5).
38	"Data Warning: Calibrator Warning"	<ul style="list-style-type: none"> • Start a calibration not later than the following day.
39	"Data Warning: Below Normal"	<ul style="list-style-type: none"> • Information that sample is below the cut off.
40	"Data Warning: Above Normal"	<ul style="list-style-type: none"> • Information that sample is above the cut off.
42	"Concentration not Consistent with Dilution Used"	<ul style="list-style-type: none"> • Relaunch the test with another dilution. Use a lower dilution factor.
43	"Sample Pipetting Problem - Check Sample"	<ul style="list-style-type: none"> • Check the sample for bubbles or foam and relaunch the test.
44	"Data Warning: Abnormal Kinetics"	<ul style="list-style-type: none"> • Prime the system and relaunch the test. • Check the samples for bubbles. • Check sample and reaction plate for dust. • Check if there are any bubbles in the reaction plate. • If the problem persists, call the Hot Line.
45	"_Data Warning: Second Response Used"	<ul style="list-style-type: none"> • Information only when logon as Service.
46	"_Data Warning: Short Lived Fluorescence"	<ul style="list-style-type: none"> • Information only when logon as Service.
48	"Data Warning: Control out of 2 SD range"	<ul style="list-style-type: none"> • Repeat measurement with a new aliquot of the control. • Recalibrate the reagent. • Use a fresh reagent box and reconstitute a fresh control vial.
51	"Inconsistent incubation time"	<ul style="list-style-type: none"> • Check in the result list the time of "end of test" and "start of test". If the real incubation time = theoretical incubation time + <4 minutes => the result can be accepted. If not check if system is disconnected. • Close the KRYPTOR compact PLUS program or restart the computer.
52	"System Error: Missed flashes"	<ul style="list-style-type: none"> • Call the Hot Line.
69	"Pre-incubating"	<ul style="list-style-type: none"> • Preincubation test, awaiting Phase II incubation.
70	"System Error: Pre-incubation has failed"	<ul style="list-style-type: none"> • Pipetting problem during the Phase II. • Reconstitution is required during the pre-incubation. • No reagent unit on board at the end of the pre-incubation. • The carousel hood stills open at the end of the pre-incubation. • Kit discrepancy because the diluent is empty (Insufficient reagent volume).

Degraded Mode

When one of the following sensors has a malfunction the local hotline can disable it to finalize the ongoing measurement protocol. The sensors listed are installed to warrant safe routine use. As the sensor malfunctions do not affect the patient result they can be disabled. In these exceptional cases the user has to follow the precautions described below. This allows use of the instrument until service engineer intervention.

After the local hotline has disabled a sensor, the window "Securities disabled" appears on the interface.

When you click on the window, detail of sensors disabled is shown.

Specific conditions to use the instrument

Caution

The user must be careful when using the instrument because the securities are off and some specific conditions are mandatory. The user must follow the advice given.

Fluidic Hood: The tip and the carousel may move while the fluidic hood is open. The user must be vigilant at each requested action on the interface (starting a pipetting sequence, a carousel scan, maintenances, etc ...).

Lower Barcode Reader: Sample tubes and reagent kits must be identified by manual entry. Refer to „Important User informations“ for precautions.

Upper Barcode Reader: Sample cassettes are not detected. Their positions are imposed on the carousel in positions 4 and 5. The dilution plates and SOLUTIONS 1-4 must be identified manually. Refer to the "Important User Information" for precautions. Consumables status (liquid level, dilution wells available,...) will be managed as usual.

Plexi Door: The Plexi door is the clear flap which is lowered when loading a new reaction plate. The software no longer identifies if it is open or closed. Check that the plexi door is closed before loading the reaction plate to prevent damage to the reaction plate carriage.

Reader Front Door: This door is used to load the reaction plate and the software no longer identifies if it is open or closed. Check that the door is closed after loading the reaction plate to avoid allowing light into the reaction area and maintaining temperature control of the incubator.

Silica Window: This is the Reader Head Window and the heating is disabled. There are no special precautions to be applied in addition to daily maintenance.

Clot Detection: The clot detection for samples is disabled, check samples before placing them on the carousel.

Waste Status: The liquid level in the Waste bottle is no longer managed as well as in the waste collector placed below the washbowl. The status of the waste is in green colour on the screen. Check waste level before starting a pipetting sequence. Check that no liquid comes out underneath the instrument.

IRDA Communication: IRDA is used to inform the system about temperature of the reagent cassettes. The temperature is still controlled accurately but it is not reported to the system which leads to display of the alarm "System Warning: reagent cooler out of range" on all test results. The temperature status of the reagent cassette appears in red on the screen. When IRDA communication is OFF, the position of the reagent cassette is imposed at location 1 of the carousel. The identification of the reagent cassette is R1. Reagent kits are always detected through the lower bar code reader.

B-R-A-H-M-S GmbH
Neuendorfstr. 25
16761 Hennigsdorf
Germany

+49-3302-883-0
+49-3302-883-100 fax
info.brahms@thermofisher.com
www.thermoscientific.com/brahms

www.thermoscientific.com/copeptin
www.thermoscientific.com/proadrenomedullin
www.thermoscientific.com/procalcitonin
www.thermoscientific.com/kryptor

International Product Support

CEZANNE SAS
Allée Graham Bell
Parc Scientifique Georges Besse
30035 Nîmes Cedex 01
France

+33-466-365-246
+33-466-365-261 fax
productsupport.brahms.frnim@thermofisher.com

Date: 20.11.2013
This version supersedes all earlier User Manuals.
Software 7.05.05

thermoscientific.com/brahms

© 2013 Thermo Fisher Scientific Inc. All rights reserved. Windows is a registered trademark of Microsoft Corporation in the United States and other countries. KRYPTOR and TRACE are registered trademarks of CIS bio international, licensed for use by B-R-A-H-M-S, a part of Thermo Fisher Scientific. CEZANNE SAS hold exclusive rights in the B-R-A-H-M-S KRYPTOR compact PLUS software. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. All data regarding specifications, terms and pricing correspond to the existing knowledge at the time of the printing. We are not responsible for any errors, misprints or changes. Reprint, also in parts, solely with prior written consent of B-R-A-H-M-S GmbH.

Thermo Fisher Scientific products are distributed worldwide; not all intended uses and applications mentioned in this printing are registered in every country.

Clinical Diagnostics

Thermo Fisher Scientific
B-R-A-H-M-S GmbH
Neuendorfstr. 25
16761 Hennigsdorf
Germany

+49 (0)3302 883 0
+49 (0)3302 883 100 fax
info.brahms@thermofisher.com

www.thermoscientific.com/brahms
www.thermoscientific.com/copeptin
www.thermoscientific.com/proadrenomedullin
www.thermoscientific.com/procalcitonin
www.thermoscientific.com/kryptor

Thermo
S C I E N T I F I C
Part of Thermo Fisher Scientific