Naudojimo instrukcija

B·R·A·H·M·S MR-proADM KRYPTOR

IVD
tik profesionaliam
naudojimui
TRACE™

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Bibliografija

Analitė

- [1] T. Eto. A review of the biological properties and clinical implications of adrenomedullin and proadrenomedullin N-terminal 20 peptide (PAMP), hypotensive and vasodilating peptides. *Peptides* 2001;22 (11): 1693-711.
- [2] J. Struck, C. Tao, N. G. Morgenthaler and A. Bergmann. Identification of an Adrenomedullin precursor fragment in plasma of sepsis patients. *Peptides* 2004;25 (8): 1369-72.
- [3] Caruhel P, et al. (2009): Homogeneous time-resolved fluoroimmunoassay for the measurement of midregional proadrenomedullin in plasma on the fully automated system B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD *Clin Biochem* 42: (7-8); 725-8.
- [4] Kruger S, Ewig S, Giersdorf S, Hartmann O, Suttrop N, Welte T. Cardiovascular and inflammatory biomarkers to predict short- and long-term survival in community-acquired pneumonia: Results from the German Competence Network, CAPNETZ. *Am J Respir Crit Care Med*. Dec 1 2010;182(11):1426-1434.
- [5] Renaud B, Schuetz P, Claessens YE, Labarere J, Albrich W, Mueller B. Proadrenomedullin improves Risk of Early Admission to ICU score for predicting early severe community-acquired pneumonia. *Chest*. Dec 2012;142(6):1447-1454.
- [6] Gouya G, Sturm G, Lamina C, et al. The association of mid-regional pro-adrenomedullin and mid-regional pro-atrial natriuretic peptide with mortality in an incident dialysis cohort. *PLoS One*. 2011;6(3):e17803.
- [7] Litke A, Bossart R, Regez K, et al. The potential impact of biomarkerguided triage decisions for patients with urinary tract infections. *Infection*. Aug 2013;41(4):799-809.
- [8] Abraham E, Laterre PF, Garg R, et al. Drotrecogin alfa (activated) for adults with severe sepsis and a low risk of death. *N Engl J Med*. Sep 29 2005;353(13):1332-1341.
- [9] Christ-Crain M, Morgenthaler NG, Struck J, Harbarth S, Bergmann A, Muller B. Mid-regional pro-adrenomedullin as a prognostic marker in sepsis: an observational study. *Crit Care*. 2005;9(6):R816-824.
- [10] Guignat C, Voirin N, Venet F, et al. Assessment of pro-vasopressin and pro-adrenomedullin as predictors of 28-day mortality in septic shock patients. *Intensive Care Med*. Nov 2009;35(11):1859-1867.
- [11] Hagag AA, Elmahdy HS, Ezzat AA. Prognostic value of plasma proadrenomedullin and antithrombin levels in neonatal sepsis. *Indian Pediatr*. Jun 2011;48(6):471-473.
- [12] Klip IT, Voors AA, Anker SD, et al. Prognostic value of mid-regional pro-adrenomedullin in patients with heart failure after an acute myocardial infarction. *Heart*. Jun 2011;97(11):892-898.
- [13] Adlbrecht C, Hulsmann M, Strunk G, et al. Prognostic value of plasma midregional pro-adrenomedullin and C-terminal-pro-endothelin-1 in chronic heart failure outpatients. *Eur J Heart Fail*. Apr 2009;11(4): 361-366.
- [14] Nickel CH, Messmer AS, Geigy N, et al. Stress markers predict mortality in patients with nonspecific complaints presenting to the emergency department and may be a useful risk stratification tool to support disposition planning. *Acad Emerg Med*. Jul 2013;20(7): 670-679.

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Alexa Fluor® is a registered trademark of Molecular Probes, licensed for use by B·R·A·H·M·S.

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- New version

Intended Use

B·R·A·H·M·S CgA II KRYPTOR is an automated immunofluorescent assay for the quantitative determination of the concentration of Chromogranin A in human serum.

Instruments

B·R·A·H·M·S CgA II KRYPTOR may be used on:

- ☒ B·R·A·H·M·S KRYPTOR
- ☒ B·R·A·H·M·S KRYPTOR compact
- ☒ B·R·A·H·M·S KRYPTOR compact PLUS

Introduction

The B·R·A·H·M·S CgA II KRYPTOR test is indicated as an aid to be used in conjunction with clinical evaluation in the follow up / monitoring of patients with Neuroendocrine Tumors (NET).

Chromogranin A (CgA) is an acidic, hydrophilic protein of 49 kDa present in chromaffin granules of the neuroendocrine cells and is a member of the granin family.

High levels of CgA could also be found in cases of benign diseases (gastro-intestinal disorders, kidney failure, cardiovascular disorders).

Chromogranin A values may rise during treatment with proton pump inhibitors (PPI). It is recommended to stop PPI treatment for at least two weeks before determination of CgA.

Contents

Kit

B·R·A·H·M·S CgA II KRYPTOR

CONT

REF 839.050

50

2...8 °C

see label for expiry date

Name	Quantity	Quality	Description
K-CONJUGATE	VIAL 1	LYOPH	Anti-CgA monoclonal antibody conjugated with cryptate, buffer, bovine albumin, bovine immunoglobulins, mouse immunoglobulins, trehalose, mannitol

Name	Quantity	Quality	Description
DILUENT	VIAL 1	ready for use	Human serum, Kathon, EDTA
XL-CONJUGATE	VIAL 1	ready for use	Anti-CgA monoclonal antibody conjugated with Alexa Fluor® 647, buffer, bovine albumin, bovine immunoglobulins, mouse immunoglobulins, potassium fluoride
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS User Manual. The bar code card contains all information needed for the registration of a new reagent lot.

Accessories

B·R·A·H·M·S CgA II KRYPTOR CAL

CAL

Not supplied with the kit.

REF 83991

2...8 °C


see label for expiry date

Intended Use: To readjust the standard curve memorized in B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS.

Name	Quantity	Quality	Description
B·R·A·H·M·S CgA II KRYPTOR CAL	VIAL 1	LYOPH	Recombinant Chromogranin A diluted in horse serum
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS User Manual. The bar code card contains information related to the calibrator lot including its concentration

Not supplied with the kit.

REF 83992

see label for expiry date

Intended Use: Quality control on board of B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS for the assay.

Name	Quantity	Quality	Description
B-R-A-H-M-S CgA II KRYPTOR QC Control 1	VIALS 3	LYOPH	Recombinant Chromogranin A diluted in horse serum
B-R-A-H-M-S CgA II KRYPTOR QC Control 2	VIALS 3	LYOPH	Recombinant Chromogranin A diluted in horse serum
Bar code card	1	ready for use	See the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS User Manual. The bar code card contains information related to the control lot, particularly the target concentrations, the standard deviations obtained and the concentration acceptance ranges. This information is visible on the B-R-A-H-M-S KRYPTOR/ KRYPTOR compact/KRYPTOR compact PLUS monitor screen in the quality control section.
Bar code stick-on labels	20 for each control	ready for use	See the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS User Manual. The bar code stick-on labels are used for identifying the controls when assayed on B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS.

Other requisites

Not supplied with the kit.

B-R-A-H-M-S KRYPTOR Consumables

Name	REF
B-R-A-H-M-S KRYPTOR BUFFER	89970
B-R-A-H-M-S KRYPTOR SOLUTION 1	89971
B-R-A-H-M-S KRYPTOR SOLUTION 2	89972
B-R-A-H-M-S KRYPTOR SOLUTION 3	89973
B-R-A-H-M-S KRYPTOR SOLUTION 4	89974
B-R-A-H-M-S KRYPTOR DILCUP	89975
B-R-A-H-M-S KRYPTOR REACT	89976

B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS Consumables

Name	REF
B-R-A-H-M-S KRYPTOR BUFFER	89970
B-R-A-H-M-S KRYPTOR compact SOLUTION 1	89981
B-R-A-H-M-S KRYPTOR compact SOLUTION 2	89982
B-R-A-H-M-S KRYPTOR compact SOLUTION 3	89983
B-R-A-H-M-S KRYPTOR compact SOLUTION 4	89984
B-R-A-H-M-S KRYPTOR compact DILCUP	89985
B-R-A-H-M-S KRYPTOR compact REACT	89986

■ Precautions

Respect the user information given in the User Manual of the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS and in this document.

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■ Principle

The measurement principle of B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS is based on TRACE™ Technology (Time-Resolved Amplified Cryptate Emission), which measures the signal that is emitted from an immunocomplex with time delay. The basis of the TRACE™ Technology is non-radiative energy transfer from a donor (a cage-like structure with an europium ion in the center [Cryptate]) to an acceptor (Alexa Fluor® 647). The proximity of donor (Cryptate) and acceptor (Alexa Fluor® 647) when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the cryptate and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.

Precise measuring of analyte concentration: When the sample is excited with a nitrogen laser at 337 nm, the donor (Cryptate) emits a long-life fluorescent signal in the milli-second range at 620 nm, while the acceptor (Alexa Fluor® 647) generates a short-life signal in the nanosecond-range at 647 nm.

When the two components are bound in an immunocomplex, both the signal amplification and the prolongation of the life span of the acceptor signal occur at 647 nm, so that it can be measured over μ-seconds.

Reliable prevention of interference: Non-specific signals, e.g. the signals of the short-life and unbound acceptor Alexa Fluor® 647 and the medium-specific interference signals conditional upon the natural fluorescence of the sample, are eliminated by temporal delay of the fluorescence measurement. The signal generated by the cryptate at 620 nm serves as an internal reference and is measured simultaneously with the long-life acceptor signal at 647 nm which is the specific signal. Interfering influences, e.g. from turbid sera, are automatically corrected by means of the internally calculated ratio of the intensities at these wavelengths.

■ Instructions

Sample volume.....	14 µL
Incubation time.....	29 min
Results are given in.....	ng/mL
Conversion factor.....	1 ng/mL = 0.0204 nmol/L
Direct measuring range.....	11.8...3 000 ng/mL
Measuring range with automatic dilution.....	11.8...1 000 000 ng/mL
Sample type.....	serum, EDTA plasma
Kit stability on board.....	29 days
Calibrator.....	1 point
Calibration stability.....	15 days
Assay principle.....	sandwich

- If the assay is carried out within 48 hours of blood collection, store samples at room temperature or at 2...8 °C. If not, samples shall be aliquoted and stored frozen (-20 °C). Up to four freezing/thawing cycles are possible.
- Citrated plasma should not be used.
- If the samples are taken on EDTA plasma, the values are systematically lower compared to serum.
- Place the sample in a tube suited for use on B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS (11–17 mm diameter). This might be the primary tube.
- The sample tube must contain a minimum volume which will vary depending on the diameter of the sample tube. A 13 mm diameter tube will require an additional 150 µL of sample.
- Should a dilution be requested either automatically or by the user, the volume of sample necessary will be 75 µL maximum.
- Icteric, hemolytic or hyperlipemic samples, or samples which are turbid or contain fibrin may yield imprecise results. Such samples are signaled by B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS.


The operation and maintenance of B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS are described in the related User Manual.

After it has been opened the reagent unit may be stored on the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS in the space provided. Each reagent unit is individually identified (bar code) and its maximum period of use after opening is controlled by the instrument.

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Take care to remove all the foil.

Opening the Kit

Note: Always handle the reagent unit carefully to avoid generation of foam or bubbles. Such foam/bubbles may interfere with proper detection of reagent

but also correct dispensing of reagents. To avoid the presence of foam/ bubbles allow the kit to stand 5...10 min within B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS prior to use.

A standard curve does not need to be constructed on B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS.

The instrument stores the required information after reading the bar code from the reagent card. A calibration must be carried out for every new reagent lot, then repeated on a regular basis. B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS automatically indicate when a calibration is required. Both the stored and recalibrated standard curve may be displayed on the screen.

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- Conjugates and sample are dispensed into the reaction plate and the signal emitted is measured periodically.
- After measurement of the fluorescent signal, the data obtained from the software are compared to the memorized standard curve.

Calibration CAL

- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) indicated on the vial label.
- Mix gently after reconstitution. Allow 15 min for the complete dissolution of the lyophilisate.
- Calibration must be carried out with every new reagent kit lot, it is then repeated on a regular basis automatically managed by B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS in order to re-adjust the standard curve.
- It is recommended that controls be run after each calibration.
- Use the calibrator/calibrators only once.
- Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5 hours.
- The calibrator bar code card must be read in for each new lot of calibrator.
- For further information see the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS User Manual.

Control CONTROL

- It is recommended that controls be run at least once a day and after each calibration.
- A control tube is processed like a sample tube.
- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) given on the vial label.
- Allow 15 min. for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- After reconstitution, do not keep a vial more than 2 hours at 18...25 °C or 24 hours at 2...8 °C.
- It is recommended that the contents of a reconstituted vial be divided into aliquots which may then be stored frozen at -20 °C for a maximum period of 1 month.
- Use one of the tubes immediately for measurement.
- After thawing an aliquot, mix gently and use immediately for measurement.
- Once thawed, a control aliquot must not be refrozen.
- The bar code stick-on labels are used for identifying the controls when assayed on B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS.
- The control kit bar code card must be scanned in for each new lot of control.
- For further information see the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS User Manual.

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If desired, B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS can automatically check the quality of assays at intervals, by statistical analysis on the basis of Levey Jennings graphs.

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■ Reference Range

Serum samples: A study of apparently healthy population (n=176, both sexes) showed that 95 % of the samples were below 101.9 ng/mL (95th percentile; median: 45.8 ng/mL).

EDTA plasma: A study of apparently healthy population (n=116, both sexes) showed that 95 % of the samples were below 76.3 ng/mL (95th percentile; median: 34.5 ng/mL).

Note: It is recommended that each laboratory establishes its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values have been established on serum and are given only as an example.

■ Assay Performance

Detection Limit

The **Detection Limit (LOD)**, calculated using the imprecision profile, is estimated in accordance with CLSI guidelines EP17-AE (Protocols for Determination of Limits of Detection and Limits of quantitation). It has been assessed as being 11.8 ng/mL.

Sensitivity

The **Functional Assay Sensitivity (FAS)**, detected by inter-assay precision of 20 %, is estimated in accordance with CLSI guidelines EP5-A2 (Evaluation of Precision Performance of Quantitative Measurement Methods). It has been assessed as being 13.7 ng/mL.

The **Limit of Quantitation (LOQ)** is the lowest reported concentration level with total error (imprecision and bias) at 40 %. The LOQ Is estimated in accordance with CLSI guidelines EP17-AE (Protocols for Determination of Limits of Detection and Limits of quantitation). It has been assessed as being 13.71 ng/mL.

Specificity

The following cross-reactivities were found:

cross-reactant	Amino acids	concentration tested	max. cross reactivity
Parastatin (porcine)	347 to 365	100 nmol/L	0.08 %
Catestatin (human)	352 to 372	452 nmol/L	0.15 %
Pancreastatin (human)	250 to 301	182 nmol/L	2.1 %
Vasostatin I (human)	17 to 76	9 nmol/L	6.7 %
Vasostatin II C-term (human)	97 to 131	15 nmol/L	1.8 %
Vasostatin II (human)	19 to 131	5 nmol/L	21.6 %
Chromostatin (bovin)	124 to 143	10 nmol/L	11.0 %
Chromogranin A protein fragment (human)	260 to 454	217 nmol/L	0.19 %
Chromogranin B (Secretogranin 1) (human)	21 to 677	72 nmol/L	1.3 %
Chromogranin C (Secretogranin 2) (human)	1 to 617	148 nmol/L	0.72 %
WE14 (human)	315 to 329	606 nmol/L	0.18 %

Accuracy/ Linearity

Dilution of highly concentrated samples revealed recovery rates between 80 % and 120 %. 0.02 % of the measured samples did not dilute linearly.

Accuracy/ "High Dose Hook" Effect

– not applicable –

Precision/ Reproducibility for Intra Assay CV

This was evaluated using samples of different concentration.

Concentration range	Intra-Assay-CV
30-80 ng/mL	≤ 7 %
80–120 ng/mL	≤ 3 %
120–1000 ng/mL	≤ 3 %
> 1000 ng/mL	≤ 3 %
Out of range samples (> 3000 ng/mL)	≤ 5 %

Precision/ Reproducibility for Inter Assay CV

This was evaluated using samples of different concentration.

Concentration range	Inter-Assay-CV
30-80 ng/mL	≤ 12 %
80–120 ng/mL	≤ 8 %
120–1000 ng/mL	≤ 8 %
> 1000 ng/mL	≤ 6 %
Out of range samples (> 3000 ng/mL)	≤ 9 %

Disturbing Factors

Factor	Description
Hemoglobin	No significant effect up to 5 mg/mL
Bilirubin	No significant effect up to 0.5 mg/mL
Triglycerides	No significant effect up to 5 mg/mL
Human albumin	No significant effect up to 50 mg/mL

Traceability

A method comparison of B-R-A-H-M-S CgA II KRYPTOR with a commercial available assay (B-R-A-H-M-S Chromogranin A KRYPTOR; REF: 823.050) using fresh routine samples gave the following correlation:

number of samples measured: 98

Passing/ Bablok: $y=1.03x + 6.93$

r spearman: $r=0.99$

The samples concentrations were between 17.9 to 2 197 ng/mL.

■ Bibliography

Analyte

[1] Hsiao RJ, Seeger RC, Yu AL, O'Connor DT. Chromogranin A in children with neuroblastoma. J Clin Invest. 1990;85:1555-9.

[2] Nobels FRE, Kwekkeboom DJ, Coopmans W, et al. Chromogranin A as serum marker for neuroendocrine neoplasia: comparison with neuron-specific enolase and the a-subunit of glycoprotein hormones. J Clin Endocrinol Metab. 1997;82:2622-8.

[3] D'Herbomez M, Forzy G, Bauters C et al. An analysis of the biochemical diagnosis of 66 pheochromocytomas. Eur J Endocrinol. 2007.156(5):569-75.

[4] Cazaentre T, d'Herbomez M, Rouaix N, Bauters C, Wémeau JL, Marchandise X . Chromogranin A in the follow-up of pheochromocytomas. Médecine Nucléaire. 2008.32:354-60.

[5] Di Comite G, Morganti A. Chromogranin A: a novel factor acting at the cross road between the neuroendocrine and the cardiovascular systems. J Hypertens. 2011 Mar;29(3):409-14.

- [6] Stokkel MP, Rietbergen DD, Korse CM, Taal BG. Somatostatin receptor scintigraphy and chromogranin A assay in staging and follow-up of patients with well-differentiated neuroendocrine tumors. Nucl Med Commun. 2011 Aug;32(8):731-7.
- [7] Arnold R. Endocrine tumours of the gastrointestinal tract. Introduction: definition, historical aspects, classification, staging, prognosis and therapeutic options. Best Pract Res Clin Gastroenterol. 2005 Aug; 19(4):491-505.
- [8] Korse CM, Bonfrer J, Aaronson NK, Hart A, Taal B. Chromogranin A as an Alternative to 5-Hydroxyindoleacetic Acid in the Evaluation of Symptoms during Treatment of Patients with Neuroendocrine Tumors. Neuroendocrinology 2009;89:296–301.
- [9] Irvin M. Modlin. Chromogranin A—Biological Function and Clinical Utility in Neuro Endocrine Tumor Disease. Ann Surg Oncol (2010) 17:2427–2443
- [10] Davide Campana, Francesca Nori, Lidya Piscitelli. Chromogranin A: Is It a Useful Marker of Neuroendocrine Tumors? J Clin Oncol 25, 2007 : 1967-1973
- [11] Aaron I. Vinik, MD, PhD, Maria P. Silva, MD. Biochemical Testing for Neuroendocrine Tumors. Pancreas 2009;38: 876-889

Technology

- [A] Mathis G, Lehn JM.: Trace – Another Story of Time. Isotopics 1995; Vol. 9.
- [B] Mathis G. , Clin. Chem. 1993; 39:1953-9.
- [C] B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS User Manual.
- [D] B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS Manual.

■ Symbols

Symbols used in Instruction for Use and Product Labelling of B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
Intended Use	Reference to the Intended use of the Medical Device	IVD	In Vitro Diagnostic Medical Device	LOT	Batch Code
CONT	Contents	CAL	Calibrator	CONTROL	Control
BUF	Buffer	SOLN 1	B-R-A-H-M-S KRYPTOR SOLUTION 1/B-R-A-H-M-S KRYPTOR compact SOLUTION 1	SOLN 2	B-R-A-H-M-S KRYPTOR SOLUTION 2/B-R-A-H-M-S KRYPTOR compact SOLUTION 2
SOLN 3	B-R-A-H-M-S KRYPTOR SOLUTION 3/B-R-A-H-M-S KRYPTOR compact SOLUTION 3	SOLN 4	B-R-A-H-M-S KRYPTOR SOLUTION 4/B-R-A-H-M-S KRYPTOR compact SOLUTION 4	CONT BAGS	Bags contained
BAGS	Bags	CONT PLATES	Plates contained	PLATES	Plates
CONT VIALS	Vials contained	VIALS	Vials	VIAL	Vial
H₂O	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL	LYOPH	Lyophilized, Freeze Dried	RCNS	Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
®	Registered Trade Mark	TM	Trade Mark	REF	Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		Harmful
	Irritant	TRACE	Trade Mark for TRACE™-technology	CE	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices

Symbol	Usage	Symbol	Usage	Symbol	Usage
	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg.Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste

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[6] Stokkel MP, Rietbergen DD, Korse CM, Taal BG. Somatostatin receptor scintigraphy and chromogranin A assay in staging and follow-up of patients with well-differentiated neuroendocrine tumors. Nucl Med Commun. 2011 Aug;32(8):731-7.

[7] Arnold R. Endocrine tumours of the gastrointestinal tract. Introduction: definition, historical aspects, classification, staging, prognosis and therapeutic options. Best Pract Res Clin Gastroenterol. 2005 Aug; 19(4):491-505.

[8] Korse CM, Bonfrer J, Aaronson NK, Hart A,Taal B. Chromogranin A as an Alternative to 5-Hydroxyindoleacetic Acid in the Evaluation of Symptoms during Treatment of Patients with Neuroendocrine Tumors. Neuroendocrinology 2009;89:296–301.

[9] Irvin M. Modlin.Chromogranin A—Biological Function and Clinical Utility in Neuro Endocrine Tumor Disease.Ann Surg Oncol (2010) 17:2427–2443

[10] Davide Campana, Francesca Nori, Lidya Piscitelli.Chromogranin A: Is It a Useful Marker of Neuroendocrine Tumors?.J Clin Oncol 25, 2007 : 1967-1973

[11] Aaron I. Vinik, MD, PhD, Maria P. Silva, MD. Biochemical Testing for Neuroendocrine Tumors. Pancreas 2009;38: 876-889

Technologija

[A] Mathis G, Lehn JM.: Trace – Another Story of Time. Isotopics 1995; Vol. 9.

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TRACE™

The manufacture and/or use of this product is covered by one or more of the following patents:
AU2003265294; EP1575416, EP2172220; JP4711972; KR1215701; MX297152, MX275242;
PH12005500130; SG109728; US7407659, US7435419, US7947449, US5627027; ZL03822180.2

Date This version supersedes all earlier versions.

30.05.2013 Changes versus previous version:

- New version

Intended Use

B·R·A·H·M·S PIGF KRYPTOR is an automated immunofluorescent assay for the quantitative determination of the concentration of human PIGF (Placental Growth Factor) in human serum.

Instruments

B·R·A·H·M·S PIGF KRYPTOR may be used on:

- ☐ B·R·A·H·M·S KRYPTOR
- ☐ B·R·A·H·M·S KRYPTOR compact
- ☒ B·R·A·H·M·S KRYPTOR compact PLUS

Introduction

The B·R·A·H·M·S PIGF KRYPTOR test is indicated as an aid to be used in conjunction with other biological and clinical data in the screening of pregnant women to assess the probability of developing pre-eclampsia.

Reduced levels of free PIGF imply an increased probability to subsequently develop pre-eclampsia.

PIGF is an angiogenic factor and belongs to the vascular endothelial growth factors (VEGF) family. PIGF acts as a vasodilator that increases the diameter of maternal spiral arteries during placentation. In the blood circulation, PIGF is mainly bound to VEGF receptor-1 (known as fms-like tyrosine kinase-1 or Flt-1) but exists also as a free form.

Pre-eclampsia, which is characterized by gestational hypertension and proteinuria, is a major cause of perinatal and maternal morbidity and mortality and complicates about 2-8 % of pregnancies worldwide.

Clinically pre-eclampsia may vary from a mild or late form to a severe or early form requiring delivery before 34 weeks of gestation.

Contents

Kit

B·R·A·H·M·S PIGF KRYPTOR

CONT

REF 844.075

see label for expiry date

Name	Quantity	Quality	Description
K-CONJUGATE	VIAL 1	LYOPH	Anti human PIGF polyclonal antibody conjugated with Lumi4®-Tb, buffer, bovine albumin, rat immunoglobulins, bovine immunoglobulins goat immunoglobulins, trehalose, mannitol
XL-CONJUGATE	VIAL 1	ready for use	Anti human PIGF monoclonal antibody conjugated with Cy5.5, buffer, bovine albumin, rat immunoglobulins, bovine immunoglobulins, goat immunoglobulins, kathon
Bar code card	1	ready for use	see the B·R·A·H·M·S KRYPTOR compact PLUS User Manual. The bar code card contains all information needed for the registration of a new reagent lot.


Accessories

B-R-A-H-M-S PIGF KRYPTOR CAL

CAL

Not supplied with the kit.

REF 84491

2...8 °C

see label for expiry date

Intended Use: To readjust the standard curve memorized in B-R-A-H-M-S KRYPTOR compact PLUS.

Name	Quantity	Quality	Description
B-R-A-H-M-S PIGF KRYPTOR CAL	VIALS 6 (0.85 mL)	LYOPH	Defibrinated normal human plasma heat treated, recombinant human PIGF, kathon
Bar code card	1	ready for use	See the B-R-A-H-M-S KRYPTOR compact PLUS User Manual. The bar code card contains information related to the calibrator lot including its concentration.

B-R-A-H-M-S PIGF KRYPTOR QC

CONTROL

Not supplied with the kit.

REF 84492

2...8 °C

see label for expiry date

Intended Use: Quality control on board of B-R-A-H-M-S KRYPTOR compact PLUS for the assay.

Name	Quantity	Quality	Description
B-R-A-H-M-S PIGF KRYPTOR Control 1	VIALS 2 (2 mL)	LYOPH	Defibrinated normal human plasma heat treated, recombinant human PIGF, kathon
B-R-A-H-M-S PIGF KRYPTOR Control 2	VIALS 2 (2 mL)	LYOPH	Defibrinated normal human plasma heat treated, recombinant human PIGF, kathon
B-R-A-H-M-S PIGF KRYPTOR Control 3	VIALS 2 (2 mL)	LYOPH	Defibrinated normal human plasma heat treated, recombinant human PIGF, kathon
Bar code card	1	ready for use	See the B-R-A-H-M-S KRYPTOR compact PLUS user's manual. The bar code card contains information related to the control lot, particularly the target concentrations, the standard deviations obtained and the concentration acceptance ranges. This information is visible on the B-R-A-H-M-S KRYPTOR compact PLUS monitor screen in the quality control section.
Bar code stick-on labels	20 for each control	ready for use	See the B-R-A-H-M-S KRYPTOR compact PLUS User Manual. The bar code stick-on labels are used for identifying the controls when assayed on B-R-A-H-M-S KRYPTOR compact PLUS

Other requisites

Not supplied with the kit.

B-R-A-H-M-S KRYPTOR compact PLUS Consumables

Name	REF
B-R-A-H-M-S KRYPTOR BUFFER	89970
B-R-A-H-M-S KRYPTOR compact SOLUTION 1	89981
B-R-A-H-M-S KRYPTOR compact SOLUTION 2	89982
B-R-A-H-M-S KRYPTOR compact SOLUTION 3	89983
B-R-A-H-M-S KRYPTOR compact SOLUTION 4	89984
B-R-A-H-M-S KRYPTOR compact DILCUP	89985
B-R-A-H-M-S KRYPTOR compact REACT	89986

■ Precautions

Respect the user information given in the User Manual of the B-R-A-H-M-S KRYPTOR compact PLUS and in this document.

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■ Principle

The measurement principle of B·R·A·H·M·S KRYPTOR compact PLUS is based on TRACE™ Technology (Time-Resolved Amplified Cryptate Emission), which measures the signal that is emitted from an immunocomplex with time delay. The basis of the TRACE™ Technology is non-radiative energy transfer from a donor (Cryptate) to an acceptor (XL). The proximity of donor (Cryptate) and acceptor (XL) when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the cryptate and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.

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■ Instructions

Sample volume.....	70 μ L
Incubation time	29 min
Results are given in	pg/mL
Conversion factor	not applicable
Direct measuring range	3.6...7 000 pg/mL
Measuring range with automatic dilution	not applicable
Sample type	serum
Kit stability on board	15 days
Calibrator	1 point
Calibration stability	7 days
Assay principle	sandwich

- If the assay is carried out within 24 hours of blood collection, store samples at 2...8 °C. If not, samples shall be aliquoted and stored frozen (-20 °C).
- Citrated plasma, heparin and EDTA should not be used.
- Place the sample in a tube suited for use on B·R·A·H·M·S KRYPTOR compact PLUS (11–17 mm diameter). This might be the primary tube.
- The sample tube must contain a minimum volume which will vary depending on the diameter of the sample tube. A 13 mm diameter tube will require an additional 150 μ L of sample.
- PIGF concentrations in patient samples are not expected to exceed the direct measuring range. Samples cannot be diluted.
- Icteric, hemolytic or hyperlipemic samples, or samples which are turbid or contain fibrin may yield imprecise results. Such samples are signaled by B·R·A·H·M·S KRYPTOR compact PLUS.

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
The instrument stores the required information after reading the bar code from the reagent card. A calibration must be carried out for every new reagent lot, then repeated on a regular basis. B·R·A·H·M·S KRYPTOR compact PLUS automatically indicates when a calibration is required. Both the stored and recalibrated standard curve may be displayed on the screen.

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Calibration CAL

- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 μ S/cm is recommended) indicated on the vial label.
- Allow 15 min for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- Calibration must be carried out with every new reagent kit lot, it is then repeated on a regular basis automatically managed by B·R·A·H·M·S KRYPTOR compact PLUS in order to readjust the standard curve.
- Use the calibrator/calibrators only once.
- Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5 hours.
- The calibrator bar code card must be read in for each new lot of calibrator.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS User Manual.

Control CONTROL

- It is recommended that controls be run at least once a day and after each calibration.
- A control tube is processed like a sample tube.
- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 μ S/cm is recommended) given on the vial label.
- Allow 15 min. for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- After reconstitution, do not keep a vial more than 5 hours at 18...25 °C or 24 hours at 2...8 °C.
- It is recommended that the contents of a reconstituted vial be divided into aliquots which may then be stored frozen at -20 °C for a maximum period of 1 month.
- Use one of the tubes immediately for measurement.
- After thawing an aliquot, mix gently and use immediately for measurement.
- Once thawed, a control aliquot must not be refrozen.
- The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR compact PLUS.
- The control kit bar code card must be scanned in for each new lot of control.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS User Manual.

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■ Reference Range

Normal non pregnant subjects

A study of 100 presumably healthy male and non pregnant female subjects showed that 95 % of the samples were found below 16.5 pg/mL. The median value was 10.3 pg/mL.

Normal values during pregnancy

A continuous study of the Fetal Medicine Foundation in U.K., based on 798 pregnant women, reveals the following results using the B·R·A·H·M·S PIGF KRYPTOR assay:

Gestational age	n	5th percentile (pg/mL)	Median (pg/mL)	95th percentile (pg/mL)
8	99	9.3	17.1	28.9
9	100	11.7	22.7	37.8
10	150	16.2	28.2	55.0
11	83	19.0	32.7	61.0
12	231	21.3	42.3	96.7
13	122	25.2	54.7	100.4
14	13	34.3	61.6	132.0

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■ Assay Performance

Detection Limit

The **Limit of Detection (LOD)**, and **Limit of Blank (LOB)** have been determined following CLSI EP17-A2 guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). The LOB has been assessed as being 0.0 pg/mL. The LOD has been assessed as being 3.6 pg/mL.

Sensitivity

The **Functional Assay Sensitivity (FAS)**, detected by inter-assay precision of 20 %, is estimated in accordance with CLSI EP05-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods). It has been assessed as being 6.7 pg/mL.

The **Limit of Quantitation (LOQ)** is the lowest reported concentration level with total error (imprecision and bias) at 40 %. The LOQ is estimated in accordance with CLSI EP17-A guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). It has been assessed as being 6.9 pg/mL.

Specificity

The antibodies used in this assay are specific for the measurement of free PIGF-1.

The **cross reactivity** is determined in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with PIGF concentrations of 12 pg/mL, 35 pg/mL or 480 pg/mL and found to cross-react at the concentrations indicated below:

Tested substance	Tested concentration	Observed cross reactivity
PIGF-2 (glycosylated)	10 000 pg/mL	13 %
PIGF-3	10 000 pg/mL	4 %
VEGF 165	10 000 pg/mL	No cross reactivity
VEGF/PIGF heterodimer	10 000 pg/mL	< 1 %

The **interference testing** is performed in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with PIGF concentrations of 12 pg/mL, 35 pg/mL or 480 pg/mL and found not to interfere at the concentrations indicated below:

Tested substance	Tested concentration
Acetaminophen	200 mg/L
Acetylsalicylic Acid	652 mg/L
Ascorbic Acid	30 mg/mL
Caffeine	59 mg/L
Calcium	200 mg/L
Ethanol (99.5 %)	5 %
Folic acid	18 µg/mL
Gentamicin	10 mg/L
Heparin	3 000 U/L
Ibuprofen	500 mg/L

Accuracy/ Linearity

The linearity is determined in accordance with the CLSI EP6-A guideline (Evaluation of the Linearity of Quantitative Measurement Procedures); based on dilution with recombinant PIGF the linearity of the assay is demonstrated from the LOQ up to 7 000 pg/mL. Patient samples cannot be diluted.

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Concentration range	Within-assay CV %
7 – 15 pg/mL	≤ 20
15 – 40 pg/mL	≤ 8
40 – 150 pg/mL	≤ 5
> 150 pg/mL	≤ 3

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Concentration range	Within-assay CV %
7 – 15 pg/mL	≤ 22
15 – 40 pg/mL	≤ 10
40 – 150 pg/mL	≤ 7
> 150 pg/mL	≤ 5

Disturbing Factors

This is determined in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with PIGF concentrations of 12 pg/mL, 35 pg/mL or 480 pg/mL and found not to interfere at the concentrations indicated below:

Factor	Description
Bilirubin	No significant effect to 250 mg/L
Hemoglobin	No significant effect to 5 g/L
Human serum albumin	No significant effect to 60 g/L
Triglycerids	No significant effect to 30 g/L

Traceability

An International PIGF- Reference Preparation is not available. Therefore the assay is calibrated using a reference preparation of recombinant human PIGF as highest available standard. Dilutions of this preparation are checked by regression analysis and compared to theoretical expected results. Quality Control of new calibrator manufacturing ensures reproducibility and stability of calibration along time (according to internal procedure).

Bibliography

Analyte

[1] Akolekar, R., A. Syngelaki, et al. Prediction of early, intermediate and late pre-eclampsia from maternal factors, biophysical and biochemical markers at 11-13 weeks. Prenat Diagn 2011, 31(1): 66-74.

[2] Akolekar, R., E. Zaragoza, et al. Maternal serum placental growth factor at 11 + 0 to 13 + 6 weeks of gestation in the prediction of pre-eclampsia. Ultrasound Obstet Gynecol 2008, 32(6): 732-9.

Symbols

Symbols used in Instruction for Use and Product Labelling of B-R-A-H-M-S KRYPTOR compact PLUS products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
Intended Use	Reference to the Intended use of the Medical Device	IVD	In Vitro Diagnostic Medical Device	LOT	Batch Code
CONT	Contents	CAL	Calibrator	CONTROL	Control
BUF	Buffer	SOLN 1	B-R-A-H-M-S KRYPTOR SOLUTION 1/B-R-A-H-M-S KRYPTOR compact SOLUTION 1	SOLN 2	B-R-A-H-M-S KRYPTOR SOLUTION 2/B-R-A-H-M-S KRYPTOR compact SOLUTION 2
SOLN 3	B-R-A-H-M-S KRYPTOR SOLUTION 3/B-R-A-H-M-S KRYPTOR compact SOLUTION 3	SOLN 4	B-R-A-H-M-S KRYPTOR SOLUTION 4/B-R-A-H-M-S KRYPTOR compact SOLUTION 4	CONT BAGS	Bags contained
BAGS	Bags	CONT PLATES	Plates contained	PLATES	Plates
CONT VIALS	Vials contained	VIALS	Vials	VIAL	Vial
H₂O	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL	LYOPH	Lyophilized, Freeze Dried	RCNS	Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
	Registered Trade Mark	TM	Trade Mark	REF	Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses

[3] Andraweera, P. H., G. A. Dekker, et al. The vascular endothelial growth factor family in adverse pregnancy outcomes. Hum Reprod Update 2012, 18(4): 436-57.

[4] Cowans, N. J., H. Alfthan, et al. Stability of first trimester placental growth factor in serum and whole blood. Prenat Diagn 2011, 31(12): 1193-7

[5] Crispi, F., E. Llurba, et al. Predictive value of angiogenic factors and uterine artery Doppler for early- versus late-onset pre-eclampsia and intrauterine growth restriction. Ultrasound Obstet Gynecol 2008, 31(3): 303-9.

[6] Lai, J., A. Pinas, et al. Maternal Serum Placental Growth Factor, Pregnancy-Associated Plasma Protein-A and Free beta-Human Chorionic Gonadotrophin at 30-33 Weeks in the Prediction of Pre-Eclampsia. Fetal Diagn Ther 2013, 33(3): 164-72.

Technology

[A] Mathis G, Lehn JM.: Trace – Another Story of Time. Isotopics 1995; Vol. 9.

[B] Mathis G. , Clin. Chem. 1993; 39:1953-9.

[C] B-R-A-H-M-S KRYPTOR compact PLUS User Manual.

[D] B-R-A-H-M-S KRYPTOR compact PLUS Manual.

Symbol	Usage	Symbol	Usage	Symbol	Usage
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		Harmful
	Irritant	TRACE	Trade Mark for TRACE™-technology		CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices
	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg.Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste

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Bibliografija

Analitė

- [1] Akolekar, R., A. Syngelaki, et al. Prediction of early, intermediate and late pre-eclampsia from maternal factors, biophysical and biochemical markers at 11-13 weeks. Prenat Diagn 2011, 31(1): 66-74.
- [2] Akolekar, R., E. Zaragoza, et al. Maternal serum placental growth factor at 11 + 0 to 13 + 6 weeks of gestation in the prediction of pre-eclampsia. Ultrasound Obstet Gynecol 2008, 32(6): 732-9.

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- [3] Andraweera, P. H., G. A. Dekker, et al. The vascular endothelial growth factor family in adverse pregnancy outcomes. Hum Reprod Update 2012, 18(4): 436-57.
- [4] Cowans, N. J., H. Alfthan, et al. Stability of first trimester placental growth factor in serum and whole blood. Prenat Diagn 2011, 31(12): 1193-7
- [5] Crispi, F., E. Llurba, et al. Predictive value of angiogenic factors and uterine artery Doppler for early- versus late-onset pre-eclampsia and intrauterine growth restriction. Ultrasound Obstet Gynecol 2008, 31(3): 303-9.
- [6] Lai, J., A. Pinas, et al. Maternal Serum Placental Growth Factor, Pregnancy-Associated Plasma Protein-A and Free beta-Human Chorionic Gonadotrophin at 30-33 Weeks in the Prediction of Pre-Eclampsia. Fe-tal Diagn Ther 2013, 33(3): 164-72.

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Lumi4®-Tb is a registered trademark of Lumiphore, Inc., licensed for use by B·R·A·H·M·S.

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- Changes versus previous version
- KRYPTOR GOLD implemented
 - Adjustment of Limit of Blank (LOB)
 - Completion of existing precision profile
 - Addition of Reference (Salahuddin et al)

■ Intended Purpose/Intended Use

B·R·A·H·M·S sFlt-1 KRYPTOR is an automated immunofluorescent assay for the quantitative determination of the concentration of sFlt-1 (soluble fms-like tyrosine kinase 1, also known as VEGF receptor-1) in human serum and EDTA plasma.

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■ Introduction

The B·R·A·H·M·S sFlt-1 KRYPTOR assay is intended to be used - together with the B·R·A·H·M·S PIGF plus KRYPTOR assay and additional clinical and diagnostic findings - as an aid in the diagnosis of pre-eclampsia and/or short-term prognosis of pre-eclampsia in pregnant women with suspicion of pre-eclampsia.

Pre-eclampsia, which affects about 2-8 % of pregnancies, is a major cause of perinatal and maternal morbidity and mortality. First clinical symptoms of pre-eclampsia such as new onset of hypertension and proteinuria are observed after 20 weeks of gestation. Clinically, pre-eclampsia may vary from mild to severe forms, requiring delivery before 34 weeks of gestation. The severe form of pre-eclampsia, HELLP syndrome (Hemolysis, elevated liver enzymes, low platelets), occurs in about 20 % of women affected by pre-eclampsia.^{[1] [2]}

Although the cause of pre-eclampsia remains unclear, the syndrome may be initiated by an imbalance of placental factors that induce endothelial dysfunction. The antiangiogenic factor sFlt-1 (soluble fms-like tyrosine kinase) acts as potent antagonist of proangiogenic factors such as Placental Growth Factor (PIGF) by adhering to the receptor-binding domains, preventing the interaction with endothelial receptors and thereby inducing endothelial dysfunction.^[3]

During pregnancy, PIGF levels increase progressively in first and second trimester and decrease towards term. In women with clinical pre-eclampsia,

sFlt-1 levels are significantly increased while concentrations of circulating free PIGF are significantly decreased. Significantly reduced levels of free PIGF have been shown in pregnant women who subsequently develop pre-eclampsia.^{[4] [5]}

Measurement of sFlt-1 levels together with PIGF levels in maternal serum starting in mid pregnancy can improve the current diagnostic workup of patients with suspected pre-eclampsia, which includes clinical symptoms, proteinuria and Uterine artery Doppler velocimetry. The sFlt-1/PIGF ratio is a tool for the diagnosis of pre-eclampsia by discriminating between pre-eclamptic pregnancies and other forms of pregnancy-related hypertensive disorders. In addition, the sFlt-1/PIGF ratio may be a valuable prognostic tool of subsequent maternal and fetal adverse outcome in women diagnosed with pre-eclampsia.^{[6] [7] [8] [9]}

■ Contents

Kit

B·R·A·H·M·S sFit-1 KRYPTOR CONT


REF 845.075    see label for expiry date

Name	Quantity	Quality	Description
K-CONJUGATE	VIAL 1	ready for use	1 Anti-sFit-1 monoclonal antibody conjugate with Lumi4®-Tb, buffer, bovine albumin, bovine immunoglobulins, mouse immunoglobulins, dextran, kathon
XL-CONJUGATE	VIAL 1	ready for use	1 Anti-sFit-1 monoclonal antibody conjugate with Cyanin5.5, buffer, bovine albumin, bovine immunoglobulins, mouse immunoglobulins, dextran, kathon
Bar code card	1	ready for use	see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual. The bar code card contains all information needed for the registration of a new reagent lot.

Accessories

B·R·A·H·M·S sFit-1 KRYPTOR CAL CAL

Not supplied with the kit.


REF 84591   see label for expiry date

Intended Purpose/Intended Use: To readjust the standard curve memorized in B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD.

Name	Quantity	Quality	Description
B·R·A·H·M·S sFit-1 KRYPTOR CAL 1	VIALS 6	LYOPH	Recombinant sFit-1, human defibrinated plasma, kathon, acid citric
B·R·A·H·M·S sFit-1 KRYPTOR CAL 2	VIALS 6	LYOPH	Recombinant sFit-1, human defibrinated plasma, kathon, acid citric
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual. The bar code card contains information related to the calibrator lot including its concentration

B·R·A·H·M·S sFit-1 KRYPTOR QC CONTROL

Not supplied with the kit.

REF 84592   see label for expiry date

Intended Purpose/Intended Use: Quality control on board of B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD for the assay.

Name	Quantity	Quality	Description
B·R·A·H·M·S sFit-1 KRYPTOR Control 1	VIALS 2	LYOPH	Recombinant sFit-1, human defibrinated plasma, kathon, acid citric
B·R·A·H·M·S sFit-1 KRYPTOR Control 2	VIALS 2	LYOPH	Recombinant sFit-1, human defibrinated plasma, kathon, acid citric
B·R·A·H·M·S sFit-1 KRYPTOR Control 3	VIALS 2	LYOPH	Recombinant sFit-1, human defibrinated plasma, kathon, acid citric
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual. The bar code card contains information related to the control lot, particularly the target concentrations, the standard deviations obtained and the concentration acceptance ranges. This information is visible on the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD monitor screen in the quality control section.
Bar code stick-on labels	32 for each control	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual. The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD

Other requisites

Not supplied with the kit.

B·R·A·H·M·S KRYPTOR compact PLUS Consumables

Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89981
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89982
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89983
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89984
B·R·A·H·M·S KRYPTOR compact DILCUP	89985
B·R·A·H·M·S KRYPTOR compact REACT	89986

B-R-A-H-M-S Consumables KRYPTOR GOLD

Name	REF
B-R-A-H-M-S KRYPTOR BUFFER	89970
B-R-A-H-M-S Solution 1 KRYPTOR GOLD	89991
B-R-A-H-M-S Solution 2 KRYPTOR GOLD	89992
B-R-A-H-M-S Solution 3 KRYPTOR GOLD	89993
B-R-A-H-M-S Solution 4 KRYPTOR GOLD	89994
B-R-A-H-M-S Dilution Plates KRYPTOR GOLD	89995
B-R-A-H-M-S Reaction Plates KRYPTOR GOLD	89996

■ Precautions

	Follow the information given in the user manual of the B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD and in this document.
	<ul style="list-style-type: none"> Check the expiration dates specified on the label. Reagents from different reagent kits should not be mixed. Do not separate the elements of a unit. Avoid any microbial contamination of the reagents.
	Raw materials of animal and/or human origin contained in the reagents have been tested and/or processed through virus inactivation procedures as far as possible. However, as it is impossible to strictly warrant that such products will not transmit any infectious agents, all raw materials of human or animal origin including the samples to be assayed must be treated as potentially infectious.
 	<p>The generally acknowledged safety precautions and laboratory techniques must be observed when handling reagents and patient samples.</p> <ul style="list-style-type: none"> Do not pipet by mouth. Wash hands after work.
 	<ul style="list-style-type: none"> Wear protective clothing, protective gloves and safety glasses at work.
 	<ul style="list-style-type: none"> Do not eat, drink or smoke in areas where samples or kit-reagents are handled.
	<ul style="list-style-type: none"> Remove spills with absorbing paper. All the material used for cleaning up must be disposed of as infectious laboratory waste. Prevent from getting into sewage, water, ground.
	<ul style="list-style-type: none"> Used reagent plates and reagent kits are disposed of as potential infectious laboratory waste according to local regulations. Empty containers could be recycled.

■ Principle

The measurement principle of B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD is based on TRACE™ Technology (Time-Resolved Amplified Cryptate Emission), which measures the signal that is emitted from an

immunocomplex with time delay. The basis of the TRACE™ Technology is non-radiative energy transfer from a donor (cryptate) to an acceptor (XL).

The proximity of donor (cryptate) and acceptor (XL) when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the cryptate and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.

Precise measuring of analyte concentration: When the sample is excited with a nitrogen laser at 337 nm, the donor (cryptate) emits a long-life fluorescent signal in the millisecond range at 620 nm, while the acceptor (XL) generates a short-life signal in the nanosecond-range at 707 nm.

When the two components are bound in an immunocomplex, both the signal amplification and the prolongation of the life span of the acceptor signal occur at 707 nm, so that it can be measured over microseconds. This long-life signal is proportional to the concentration of the analyte to be measured.

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■ Instructions

Sample volume.....	8 µL
Incubation time.....	9 min
Results are given in.....	pg/mL
Conversion factor.....	not applicable
Direct measuring range.....	22...90 000 pg/mL
Measuring range with automatic dilution.....	not applicable
Sample type.....	serum, EDTA plasma
Kit stability on board.....	29 days
Calibrator.....	2 points
Calibration stability.....	15 days
Assay principle.....	sandwich

	<ul style="list-style-type: none"> If the assay is carried out within 72 hours of blood collection, store samples at 2...8 °C. If not, samples shall be aliquoted and stored frozen (-20 °C). Citrated and heparin plasma should not be used. Place the sample in a tube suited for use on B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD (11–17 mm diameter). This might be the primary tube. The sample tube must contain a minimum volume which will vary depending on the diameter of the sample tube. A 13 mm diameter tube will require an additional 150 µL of sample. sFit-1 concentrations in patient samples are not expected to exceed the direct measuring range. Samples cannot be diluted. Icteric, hemolytic or hyperlipemic samples, or samples which are turbid or contain fibrin may yield imprecise results. Such samples are signaled by B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD.
---	--


The operation and maintenance of B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD are described in the related user manual.

After it has been opened the reagent unit may be stored on the B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD in the space provided. Each reagent unit is individually identified (bar code) and its maximum period of use after opening is controlled by the instrument.

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Opening the Kit

Note: Always handle the reagent unit carefully to avoid foaming or bubbles. Foam/bubbles may interfere with proper detection of reagent and/or correct dispensing of reagents. Allow the kit to stand 5...10 min within B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD prior to use.

A standard curve does not need to be established on B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD.

The instrument stores the required information after reading the bar code from the reagent card. A recalibration must be carried out for every new reagent lot, and subsequently on a regular basis. B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD automatically indicates when a calibration is required. Both, the stored and recalibrated standard curve may be displayed on the screen.

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- Conjugates and sample are dispensed into the reaction plate and the fluorescent signal emitted is measured periodically.
- After measurement of the fluorescent signal, the data obtained from the software are compared to the memorized standard curve.

Calibration CAL

- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) indicated on the vial label.
- Allow 15 min for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- Calibration must be carried out with every new reagent kit lot. It is then repeated on a regular basis automatically managed by B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD in order to readjust the standard curve.
- Use the calibrator/calibrators only once.
- Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5 hours.
- The calibrator bar code card must be read in for each new lot of calibrator.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual.

Control CONTROL

- It is recommended that controls be run at least once a day and after each calibration.
- A control tube is processed like a sample tube.
- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) as stated on the vial label.
- Allow 15 min. for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- After reconstitution, do not keep a vial more than 5 hours at 18...25 °C or 24 hours at 2...8 °C.
- It is recommended that the contents of a reconstituted vial be divided into aliquots which may then be stored frozen at -20 °C for a maximum period of 1 month.
- Use one of the tubes immediately for measurement.
- After thawing an aliquot, mix gently and use immediately for measurement.
- Once thawed, a control aliquot must not be refrozen.
- The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD.
- The control kit bar code card must be scanned in for each new lot of control.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual.

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■ Reference Range

A) Expected values

From EDTA plasma samples collected from 347 singleton pregnancies with a normal pregnancy outcome, levels of sFit-1 and PIGF were determined in parallel and the sFit-1/PIGF ratio was calculated.^[10]

sFit-1 KRYPTOR (pg/mL)

Gestational age (completed weeks)	n	5 th percentile (pg/mL)	Median (pg/mL)	95 th percentile (pg/mL)
14–19	21	764	1 515	2 545
20–24	60	626	1 319	4 443
25–28	66	612	1 231	3 958
29–33	44	697	1 395	3 722
34–36	39	1 023	2 971	8 492
37-delivery	117	961	3 079	7 524

PIGF plus KRYPTOR (pg/mL)

Gestational age (completed weeks)	n	5 th percentile (pg/mL)	Median (pg/mL)	95 th percentile (pg/mL)
14–19	21	43	101	348
20–24	60	183	391	868
25–28	66	145	459	1 000
29–33	44	71	507	1 157
34–36	39	41	170	752
37-delivery	117	37	131	750

sFit-1/PIGF ratio

Gestational age (completed weeks)	n	5 th percentile	Median	95 th percentile
14–19	21	5	14	26
20–24	60	1	4	17
25–28	66	1	3	9
29–33	44	1	3	35
34–36	39	2	17	158
37-delivery	117	2	29	131

Note: It is recommended that each laboratory establishes its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values have been established on serum and serve only as an example.

B) Clinical sensitivity and specificity

In an external study, using the B-R-A-H-M-S PIGF plus KRYPTOR and B-R-A-H-M-S sFlt-1 KRYPTOR assays in parallel on samples from 109 pregnant women with normal pregnancy outcome and 20 patients with pre-eclampsia, the published cut-off^[7] for the sFlt-1/PIGF ratio of 85 was applied to calculate sensitivity and specificity for this cohort. At a cut-off of 85 the sensitivity was calculated at 95 % and the specificity at 84 %.

All pregnancies were singleton pregnancies. Pre-eclampsia was defined as new onset of hypertension and proteinuria after 20 weeks of gestation. A pre-eclamptic pregnancy was defined as early-onset if clinical signs appeared before 34 weeks of gestation.

C) Method comparison

A comparison of the B-R-A-H-M-S sFlt-1 KRYPTOR assay (y) with a commercially available automated sFlt-1 assay (x) using clinical samples gave the following correlation (pg/mL):

Number of samples measured: 150

Passing/Bablok

$y = 0.93x - 267$

$r = 0.996$

The sample concentrations were between 686 and 43 598 pg/mL.

■ Assay Performance

Detection Limit

The **Limit of Detection (LOD)**, and **Limit of Blank (LOB)** have been determined following CLSI EP17-A2 guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). The LOB has been assessed as being 16.7 pg/mL. The LOD has been assessed as being 22 pg/mL.

Sensitivity

The **Functional Assay Sensitivity (FAS)**, detected by inter-assay precision of 20 %, is estimated in accordance with CLSI EP05-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods). It has been assessed as being 29 pg/mL.

The **Limit of Quantitation (LOQ)** is the lowest reported concentration level with total error (imprecision and bias) at 40 %The LOQ is estimated in accordance with CLSI EP17-A guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). It has been assessed as being 34 pg/mL.

Specificity

The antibodies used in this assay are specific for the measurement of sFlt-1.

The **cross reactivity** is determined in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with sFlt-1 concentrations of 1 586 pg/mL or 10 248 pg/mL and found to cross-react at the concentrations indicated below:

Tested substance	Tested concentration	Observed cross reactivity
PIGF-1	10 000 pg/mL	No cross-reactivity
PIGF-2	10 000 pg/mL	No cross-reactivity
PIGF-3	10 000 pg/mL	No cross-reactivity
VEGF-A	10 000 pg/mL	No cross-reactivity
VEGF-B	10 000 pg/mL	No cross-reactivity
VEGF/PIGF heterodimer	10 000 pg/mL	No cross-reactivity

The **interference testing** is performed in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with sFlt-1 concentrations of 1 586 pg/mL or 10 248 pg/mL and found not to interfere at the concentrations indicated below:

Tested substance	Tested concentration
Acetaminophen	200 mg/L
Acetylsalicylic Acid	652 mg/L
L-Ascorbic acid	30 mg/L
Caffeine	59 mg/L
Calcium	200 mg/L
Ethanol	5 % (v/v)
Folic Acid	18 µg/L
Gentamicin	10 mg/L
Heparin	3000 U/L
Ibuprofen	500 mg/L
A-Methyldopa	15 mg/L
Metoprolol	5 mg/L
Iron	150 µg/dL
Dihydralazin	0.2 g/L

Accuracy/ Linearity

The linearity is determined in accordance with the CLSI EP6-A guideline (Evaluation of the Linearity of Quantitative Measurement Procedures); based on dilution with recombinant sFlt-1 the linearity of the assay is demonstrated from the LOQ up to 90 000 pg/mL. Patient samples cannot be diluted.

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Disturbing Factors

This is determined in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum concentrations of 1 586 pg/mL or 10 248 pg/mL and found not to interfere at the concentrations indicated below:

Factor	Description
Bilirubin	No significant effect up to 250 mg/L
Hemoglobin	No significant effect up to 5 g/L
Human serum albumin	No significant effect up to 60 g/L
Triglycerids	No significant effect up to 30 g/L

Traceability

An International sFit-1- Reference Preparation is not available. Therefore the assay is calibrated using a reference preparation of recombinant human sFit-1 as highest available standard. Dilutions of this preparation are checked by regression analysis and compared to theoretical expected results. Quality Control of new calibrator manufacturing ensures reproducibility and stability of calibration along time (according to internal procedure).

■ Bibliography

Analyte

- [1] Ghulmiyyah L, et al (2012), Seminars in Perinatology; 36:56-59
- [2] Powe C et al (2011), Circulation; 123: 2856-69
- [3] Levine et al. (2004). N Engl J Med 2004;350:672-83.
- [4] Lam C et al (2005). Hypertension; 46: 1077-85
- [5] De Vivo, et al. (2008). Acta Obstet et Gynecol 2008; 87: 837-842
- [6] Thadhani, et al. (2004). J Clin Endo; 89(2):770-775.
- [7] Verlohren, et al (2012). Am J Obstet Gynecol. 206(1):58.e1-8.
- [8] Rana S, et al (2012). Circulation; 125 (7):911-919
- [9] Andersen L et al (2015). J Am Soc Hypertens; 9(2):86-96.
- [10] Salahuddin S et al (2016). Hypertens Pregnancy, 35(3):330-45

Technology

- [A] Mathis G, Lehn JM.: Trace – Another Story of Time. Isotopics 1995; Vol. 9.
- [B] Mathis G. , Clin. Chem. 1993; 39:1953-9.
- [C] B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD User Manual.
- [D] B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD Manual.

■ Symbols

Symbols used in Instruction for Use and Product Labelling of B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD products.

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B·R·A·H·M·S Copeptin proAVP KRYPTOR

IVD
for professional use
only

B·R·A·H·M·S is a registered trademark of B·R·A·H·M·S GmbH.
Lumi4®-Tb is a registered trademark of Lumiphore, Inc., licensed for use by B·R·A·H·M·S.

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Changes versus previous version

- Introduction: standardization of product name
- Removal of duplicated literature reference

Intended Use

B·R·A·H·M·S Copeptin proAVP KRYPTOR is an automated immunofluorescent assay to determine the concentration of Copeptin (C-terminal vasopressin = CT-proAVP) in human serum and human plasma (EDTA, heparin).

Instruments

B·R·A·H·M·S Copeptin proAVP KRYPTOR may be used on:

- ☐ B·R·A·H·M·S KRYPTOR
- ☐ B·R·A·H·M·S KRYPTOR compact
- ☒ B·R·A·H·M·S KRYPTOR compact PLUS

Introduction

The B·R·A·H·M·S Copeptin proAVP KRYPTOR assay measures the C-Terminal precursor fragment of Vasopressin (AVP or ADH). Vasopressin is a key regulator of water balance, has haemodynamic effects and plays a critical role in the regulation of hypothalamo-pituitary-adrenal axis reflecting the individual stress response^[1].

Copeptin and vasopressin (AVP, ADH) are released from the same precursor hormone in equimolar amounts into the bloodstream, hence Copeptin can be used as a surrogate marker for vasopressin^{[2][3][4]}.

B·R·A·H·M·S Copeptin proAVP KRYPTOR is indicated in conjunction with Troponin as an aid in ruling out acute myocardial infarction in patients presenting with chest pain or equivalent symptoms of cardiac origin^{[5][6][7][8]}.

B·R·A·H·M·S Copeptin proAVP KRYPTOR is indicated in conjunction with other clinical data (e.g. osmolality and hydration status) as an aid in the differential diagnosis of patients with water balance disorders, e.g. Polyuria-Polydipsia Syndrome and Syndrome of Inappropriate ADH secretion (SIADH)^{[9][10][11][12]}.

Contents

Kit

B·R·A·H·M·S Copeptin proAVP KRYPTOR

CONT

REF 857.050 50 2...8 °C see label for expiry date

Name	Quantity	Quality	Description
K- conjugate	VIAL 1	LYOPH	Anti CT-proAVP polyclonal sheep antibody conjugated with Lumi4®-Tb, buffer, BSA, bovine IgG, mouse IgG, sheep IgG, Trehalose, D-Mannitol
XL-conjugate	VIAL 1	LYOPH	Anti CT-proAVP monoclonal mouse antibody conjugated with Cy 5.5, buffer, BSA, bovine IgG, mouse IgG, sheep IgG, Trehalose, D-Mannitol
DILUENT	VIAL 1	ready for use	Horse serum, EDTA, Kathon
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS User Manual. The bar code card contains all information needed for the registration of a new reagent lot.


Accessories

B·R·A·H·M·S Copeptin proAVP KRYPTOR CAL

CAL

Not supplied with the kit.

REF 85791

see label for expiry date

Intended Use: To readjust the standard curve memorized in B·R·A·H·M·S KRYPTOR compact PLUS.


Name	Quantity	Quality	Description
B·R·A·H·M·S Copeptin proAVP KRYPTOR CAL	VIALS 6	LYOPH	Horse serum, EDTA, complete, CT-proAVP
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS User Manual. The bar code card contains information related to the calibrator lot including its concentration.

B·R·A·H·M·S Copeptin proAVP KRYPTOR QC

CONTROL

Not supplied with the kit.

REF 85792

see label for expiry date

Intended Use: Quality control on board of B·R·A·H·M·S KRYPTOR compact PLUS for the assay.

Name	Quantity	Quality	Description
B·R·A·H·M·S Copeptin proAVP KRYPTOR Control 1	VIALS 3	LYOPH	Horse serum, EDTA, complete, CT-proAVP
B·R·A·H·M·S Copeptin proAVP KRYPTOR Control 2	VIALS 3	LYOPH	Horse serum, EDTA, complete, CT-proAVP
Bar code card	1	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS user's manual. The bar code card contains information related to the control lot, particularly the target concentrations, the standard deviations obtained and the concentration acceptance ranges. This information is visible on the B·R·A·H·M·S KRYPTOR compact PLUS monitor screen in the quality control section.
Bar code stick-on labels	20 for each control	ready for use	See the B·R·A·H·M·S KRYPTOR compact PLUS User Manual. The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR compact PLUS

Other requisites

Not supplied with the kit.

B·R·A·H·M·S KRYPTOR compact PLUS Consumables

Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89981
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89982
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89983
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89984
B·R·A·H·M·S KRYPTOR compact DILCUP	89985
B·R·A·H·M·S KRYPTOR compact REACT	89986

■ Precautions

Respect the user information given in the User Manual of the B·R·A·H·M·S KRYPTOR compact PLUS and in this document.

- Observe the expiration dates specified on the label.
- Reagents from different reagent kits should not be mixed.
- Do not dissociate the elements of a unit.
- Avoid any microbic contamination of the reagents.

Raw materials of human or animal origin contained in the reagents have been tested with approved kits and found negative for the anti-HIV 1, anti-HIV 2, anti-HCV antibodies and the HBs antigen. However as it is impossible to strictly guarantee that such products will not transmit hepatitis, the HIV virus, or any other viral infection, all raw materials of human or animal origin including the samples to be assayed must be treated as potentially infectious.

The generally acknowledged safety precautions and laboratory techniques must be observed when handling reagents and patient samples.

- Do not pipet by mouth.
- Wash hands after work.

- Wear protective clothing, protective gloves and safety glasses at work.

- Do not eat, drink or smoke in areas where samples or kit-reagents are handled.

- Remove spills with absorbing paper.
- All the material used for cleaning up must be disposed of as infectious laboratory waste.
- Prevent from getting into sewage, water, ground.

- Used reagent plates and reagent kits are disposed of as potential infectious laboratory waste according to local regulations.
- Empty containers should be returned to local recyclers.

■ Principle

The measurement principle of B·R·A·H·M·S KRYPTOR compact PLUS is based on TRACE™ Technology (Time-Resolved Amplified Cryptate Emission), which measures the signal that is emitted from an immunocomplex with time delay. The basis of the TRACE™ Technology is non-radiative energy transfer from a donor (a cage-like structure with a terbium ion in the center) to an acceptor, which is part of a chemical protein (XL).

The proximity of donor (Cryptate) and acceptor (XL) when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the cryptate and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.

Precise measuring of analyte concentration: When the sample is excited with a nitrogen laser at 337 nm, the donor (Cryptate) emits a long-life fluorescent signal in the milli-second range at 620 nm, while the acceptor (XL) generates a short-life signal in the nanosecond-range at 707 nm.

When the two components are bound in an immunocomplex, both the signal amplification and the prolongation of the life span of the acceptor signal occur at 707 nm, so that it can be measured over μ -seconds. This long-life signal is proportional to the concentration of the analyte to be measured.

Reliable prevention of interference: Non-specific signals, e.g. the signals of the short-life and unbound acceptor XL and the medium-specific interference signals conditional upon the natural fluorescence of the sample, are eliminated by temporal delay of the fluorescence measurement. The signal generated by the cryptate at 620 nm serves as an internal reference and is measured simultaneously with the long-life acceptor signal at 707 nm which is the specific signal. Interfering influences, e.g. from turbid sera, are automatically corrected by means of the internally calculated ratio of the intensities at these wavelengths.

■ Instructions

Sample volume.....	50 μ L
Incubation time.....	14 min
Results are given in.....	pmol/L
Conversion factor.....	1 pmol/L = 4.02 pg/mL
Direct measuring range.....	0.7...500 pmol/L
Measuring range with automatic dilution.....	0.7...2 000 pmol/L
Sample type.....	serum, plasma (EDTA, heparin)
Kit stability on board.....	29 days
Calibrator.....	1 point
Calibration stability.....	15 days
Assay principle.....	sandwich

- If the assay is carried out within 24 hours of blood collection, store samples at 2...8 °C. If not, samples shall be aliquoted and stored frozen (-20 °C).
- Samples may be frozen and thawed three times.
- Citrated plasma should not be used.
- Place the sample in a tube suited for use on B·R·A·H·M·S KRYPTOR compact PLUS (11–17 mm diameter). This might be the primary tube.
- The sample tube must contain a minimum volume which will vary depending on the diameter of the sample tube. A 13 mm diameter tube will require an additional 150 μ L of sample.
- Should a dilution be requested either automatically or by the user, the volume of sample necessary will be 50 μ L maximum.
- Icteric, hemolytic or hyperlipemic samples, or samples which are turbid or contain fibrin may yield imprecise results. Such samples are signaled by B·R·A·H·M·S KRYPTOR compact PLUS.


The operation and maintenance of B·R·A·H·M·S KRYPTOR compact PLUS are described in the related User Manual.

After it has been opened the reagent unit may be stored on the B·R·A·H·M·S KRYPTOR compact PLUS in the space provided. Each reagent unit is individually identified (bar code) and its maximum period of use after opening is controlled by the instrument.

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Take care to remove all the foil.

Opening the Kit

Note: Always handle the reagent unit carefully to avoid generation of foam or bubbles. Such foam/bubbles may interfere with proper detection of reagent but also correct dispensing of reagents. To avoid the presence of foam/ bubbles allow the kit to stand 5...10 min within B·R·A·H·M·S KRYPTOR compact PLUS prior to use.

A standard curve does not need to be constructed on B·R·A·H·M·S KRYPTOR compact PLUS.

The instrument stores the required information after reading the bar code from the reagent card. A calibration must be carried out for every new reagent lot, and then repeated on a regular basis. B·R·A·H·M·S KRYPTOR compact PLUS automatically indicates when a calibration is required. Both the stored and recalibrated standard curve may be displayed on the screen.

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- Conjugates and sample are dispensed into the reaction plate and the signal emitted is measured periodically.
- Samples with concentrations higher than the direct measurement range, are identified in the first few minutes of incubation, then diluted automatically and reassayed.
- After measurement of the fluorescent signal, the data obtained from the software are compared to the memorized standard curve.

Calibration CAL

- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 μ S/cm is recommended) indicated on the vial label.
- Allow 15 min for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- Calibration must be carried out with every new reagent kit lot, it is then repeated on a regular basis automatically managed by B·R·A·H·M·S KRYPTOR compact PLUS in order to readjust the standard curve.
- It is recommended that controls be run after each calibration.
- Use the calibrator/calibrators only once.
- Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5 hours.
- The calibrator bar code card must be read in for each new lot of calibrator.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS User Manual.

Control **CONTROL**

- It is recommended that controls be run at least once a day and after each calibration.
- A control tube is processed like a sample tube.
- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) given on the vial label.
- Allow 15 min. for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- After reconstitution, do not keep a vial more than 5 hours at 18...25 °C or 24 hours at 2...8 °C.
- It is recommended that the contents of a reconstituted vial be divided into aliquots which may then be stored frozen at -20 °C for a maximum period of 1 month.
- Use one of the tubes immediately for measurement.
- After thawing an aliquot, mix gently and use immediately for measurement.
- Once thawed, a control aliquot must not be refrozen.
- The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR compact PLUS.
- The control kit bar code card must be scanned in for each new lot of control.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS User Manual.

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If desired, B·R·A·H·M·S KRYPTOR compact PLUS can automatically check the quality of assays at intervals, by statistical analysis on the basis of Levey Jennings graphs.

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
■ Reference Range

Note: It is recommended that each laboratory establishes its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values have been established on serum and EDTA plasma and are given only as an indications.

The results, obtained from this assay should always be assessed in combination with the clinical examination, patients medical history, plasma osmolality, plasma sodium and other findings before momentous actions will be prefaced.


A study of an apparently healthy population (N = 297, both sexes) showed a median of 3.9 pmol/L and the following distribution as listed in the histogram graph below.

I. B·R·A·H·M·S Copeptin proAVP KRYPTOR in early rule out AMI:

Following several observational study that showed the highest Negative Predictive Value of the combination of Copeptin and Troponin for ruling-out AMI,

in a randomized controlled clinical process study for early discharge using single Troponin and Copeptin in patients with suspected ACS (Acute Coronary Syndrome), a Copeptin proAVP cut-off at 10 pmol/L is recommended to rule-out AMI (Acute Myocardial Infarction) [5][6][7]. A value of 10 pmol/L or above is considered as positive result.

Modified from Möckel and Searle (2014) [8]


Note: It is recommended that each laboratory establish its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values above have been established on serum and EDTA plasma and are given only as indications.

II. B·R·A·H·M·S Copeptin proAVP KRYPTOR in the differential diagnosis of patients with water balance disorders

In a cohort of 72 healthy individuals who underwent Copeptin stimulation via water deprivation or hypo- and hypertonic saline infusion, the following Copeptin values were observed [9][12]. Copeptin values according to normal plasma osmolality:

Osmolality (mosmol/kg)	Copeptin (pmol/L)
270-280	0.81-11.6
281-285	1.0-13.7
286-290	1.5-15.3
291-295	2.3-24.5
296-300	2.4-28.2

In a prospective international multicenter clinical study (COSIP) [11] an algorithm was defined based on patients with Polyuria-Polydipsia Syndrome (n=55). Baseline Copeptin value of ≥ 21.4 pmol/L identify patients with Nephrogenic Diabetes Insipidus. A 2nd Copeptin value after fluid deprivation with a cut-off of 4.9 pmol/L can distinguish patients between Complete & Partial Diabetes Insipidus (< 4.9 pmol/L) and patients with Primary Polydipsia (≥ 4.9 pmol/L).

Modified from Timper et al. (2015)^[11]

Note: It is recommended that each laboratory establish its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values above have been established on serum and EDTA plasma and are given only as indications.

Limitations

Following diseases / therapies have an influence on the level on Copeptin proAVP:

- SIRS, Sepsis, Severe Sepsis, Septic Shock
- Cardiovascular Diseases, i.e. Chronic Heart Failure
- Low Respiratory Tract Infections, i.e. Community Acquired Pneumonia (CAP) and COPD
- AVP Receptor Antagonist therapies and other diseases in which AVP has been shown to play an important pathophysiologic role

Assay Performance

Detection Limit

The **Limit of Detection (LOD)** and **Limit of Blank (LOB)** have been determined following CLSI EP17-A2 guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). The LOB has been assessed as being 0.41 pmol/L. The LOD has been assessed as being 0.69 pmol/L.

Sensitivity

The **Functional Assay Sensitivity (FAS)**, detected by inter-assay precision of 20 %, is estimated in accordance with CLSI EP05-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods). It has been assessed as being 1.08 pmol/L.

The **Limit of Quantitation (LOQ)** is the lowest reported concentration level with total error (imprecision and bias) at 40 %. The LOQ is estimated in accordance with CLSI EP17-A guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). It has been assessed as being 1.23 pmol/L.

Specificity

The antibodies used in this assay enable a specific measurement of human C-terminal proVasopressin (CT-proAVP).

Accuracy/ Linearity

The linearity is determined in accordance with the CLSI EP6-A guideline (Evaluation of the Linearity of Quantitative Measurement Procedures). Dilution of highly concentrated samples revealed recovery rates between 80 % and 120 %.

Accuracy/ "High Dose Hook" Effect

No "High Dose Hook" Effect can be observed up to 2 000 pmol/L.

Precision/ Reproducibility for Intra Assay CV

Precision is evaluated using samples of different concentrations in accordance with the CLSI EP5-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods).

Concentration range	Within-assay CV%
2.0 – 4.0 pmol/L	< 15.0 %
4.0 – 15 pmol/L	< 8.0 %
15 – 50 pmol/L	< 4.0 %
> 50 pmol/L	< 3.0 %
Out of range samples (> 500 pmol/L)	< 4.0 %

Precision/ Reproducibility for Inter Assay CV

This was evaluated using samples of different concentrations in accordance with the CLSI EP5-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods).

Concentration range	Inter-assay CV%
2.0 – 4.0 pmol/L	< 18.0 %
4.0 – 15 pmol/L	< 10.0 %
15 – 50 pmol/L	< 5.0 %
> 50 pmol/L	< 5.0 %
Out of range samples (> 500 pmol/L)	< 6.0 %

Disturbing Factors

Interferences are estimated in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry).

Factor	Description
Bilirubin	No significant effect to 5 mg/dL
Hemoglobin	No significant effect to 500 mg/dL
Human serum albumin	No significant effect to 5 g/dL
Triglycerids	No significant effect to 5 g/L

Interferences

Interferences are estimated in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry).

The following drugs have no significant effect up to the following concentrations:

Drug	Concentration up to
Abciximab	20 µg/mL
Acetaminophen	250 µg/mL
Acetylsalicylic Acid	600 µg/mL
Allopurinol	40 µg/mL
Ambroxol	18 µg/mL
Ampicillin	50 µg/mL
Ascorbic Acid	40 µg/mL
Atenolol	10 µg/mL
Atorvastatin	24 µg/mL
Caffeine	100 µg/mL
Captopril	5 µg/mL
Chloramphenicol	25 mg/dL

Drug	Concentration up to
Cinnarizine	2.7 µg/mL
Clopidogrel Bisulfate	2.5 µg/mL
Cocaine	1 µg/mL
Diclofenac	50 µg/mL
Digoxin	7.5 ng/mL
Dopamine	900 ng/mL
Eptifibatide	7 µg/mL
Erythromycin	60 µg/mL
Fondaparinux	12 µg/mL
Furosemide	400 µg/mL
Heparin	8 U/mL
Ibuprofen	500 µg/mL
Lovastatin	20 µg/mL
Low MW Heparin	5 U/mL
Methyl dopa	15 µg/mL
Metoprolol	5 µg/mL
Nicotine	2 mg/dL
Nifedipine	600 ng/mL
Nitrofurantoin	4 µg/mL
Nitroglycerin	0.5 µg/mL
Nystatin	7.5 µg/mL
Oxytetracycline	5 µg/mL
Phenytoin	60 µg/mL
Propranolol	5 µg/mL
Quinidine	20 µg/mL
Ramipril	12 µg/mL
Salicylic Acid	600 µg/mL
Simvastatin	16 µg/mL
Theophylline	60 µg/mL
Tirofiban	1.7 µg/mL
Trimethoprim	40 µg/mL
Verapamil	2 µg/mL
Warfarin	30 µg/mL

Traceability

An International Copeptin proAVP Reference Preparation is not available. Therefore the assay is calibrated using a reference preparation of Copeptin proAVP as highest available standard. Dilutions of this preparation are checked by regression analysis and compared to theoretical expected results. Quality Control of new calibrator manufacturing ensures reproducibility and stability of calibration along time (according to internal procedure).

Bibliography

Analyte

- [1] Morgenthaler.NG, Struck J, Alonso C, Bergmann A. Assay for the measurement of copeptin, a stable peptide derived from the precursor of vasopressin. Clin Chem 2006; 52(1): 112-9
- [2] Holt NF and Haspel KL. Vasopressin: a review of therapeutic applications. J Cardiothorac Vasc Anesth 2010; 24(2): 330-47
- [3] Jochberger S, Zitt M, Luckner G, Mayr VD, Wenzel V, Ulmer H, Morgenthaler NG, Hasibeder WR, Dünser MW. Postoperative vasopressin and copeptin levels in noncardiac surgery patients: a prospective controlled trial. Shock 2009; 31(2): 132-8
- [4] Balanescu S, Kopp P, Gaskill MB, Morgenthaler NG, Schindler C, Rutishauser J. Correlation of plasma copeptin and vasopressin concentration in hypo-, iso-, and hyperosmolar states. J Clin Endocrinol Metab. 2011 Apr;96(4):1046-52.
- [5] Möckel M, Searle J, Hamm C, Slagman A, Blankenberg S, Huber K, Katus H, Liebetrau C, Müller C, Muller R, Peitsmeyer P, von Recum J, Tajsic M, Vollert JO, Giannitsis E. Early discharge using single cardiac troponin and copeptin testing in patients with suspected acute coronary syndrome (ACS): a randomized, controlled clinical process study. Eur Heart J. 2015 Feb 7; 36(6): 369–376
- [6] Post F, Gori T, Giannitsis E, Darius H, Baldus S, Hamm C, Hambrecht R, Hofmeister HM, Katus H, Perings S, Senges J, Münzel T. Criteria of the German Society of Cardiology for the establishment of chest pain units: update 2014. Clin Res Cardiol, 2015: DOI 10.1007/s00392-015-0888-2
- [7] Lipinski MJ, Escárcega RO, D'Ascenzo F, Magalhães MA, Baker NC, Torguson R, Chen F, Epstein SE, Miró O, Llorens P, Giannitsis E, Lotze U, Lefebvre S, Sebbane M, Cristol JP, Chenevier-Gobeaux C, Meune C, Eggers KM, Charpentier S, Twerenbold R, Mueller C, Biondi-Zoccai G, Waksman R. A systematic review and collaborative meta-analysis to determine the incremental value of copeptin for rapid rule-out of acute myocardial infarction. Am J Cardiol 2014;113:1581e1591
- [8] Möckel M, Searle J. Copeptin-marker of acute myocardial infarction. Curr Atheroscler Rep 2014: 16:421
- [9] Fenske W, Quinkler M, Lorenz D, Zopf K, Haagen U, Papassotiriou J, Pfeiffer AFH, Fassnacht M, Stoerk S, Allolio B. Copeptin in the differential diagnosis of the polydipsia-polyuria-syndrome — Revisiting the direct and indirect water deprivation tests. J Clin Endocrinol Metab 2011, 96(5): 1506–1515
- [10] Winzeler B, Zweifel C, Nigro N, Arici B, Bally M, Schuetz P, Blum CA, Kelly C, Berkmann S, Huber A, Gentili F, Zadeh G, Landolt H, Mariani L, Müller B, Christ-Crain M. Postoperative Copeptin Concentration Predicts Diabetes Insipidus After Pituitary Surgery. J Clin Endocrinol Metab 2015; <http://dx.doi.org/10.1210/jc.2014-4527>
- [11] Timper K, Fenske W, Kühn F, Frech N, Arici B, Rutishauser J, Kopp P, Allolio B, Stettler C, Müller B, Katan M, Christ-Crain M. Diagnostic Accuracy of Copeptin in the Differential Diagnosis of the Polyuria-polydipsia Syndrome: A Prospective Multicenter Study. J Clin Endocrinol Metab. 2015 Jun; 100(6):2268-74
- [12] Szinnai G, Morgenthaler NG, Berneis K, Struck J, Müller B, Keller U, Christ-Crain M. Changes in plasma copeptin, the C-terminal portion of arginine vasopressin during water deprivation and excess in healthy subjects. J Clin Endocrinol Metab 2007; 92(10): 3973-8

Technology

- [A] Mathis G, Lehn JM.: Trace – Another Story of Time. Isotopics 1995; Vol. 9.
- [B] Mathis G. , Clin. Chem. 1993; 39:1953-9.
- [C] B·R·A·H·M·S KRYPTOR compact PLUS User Manual.
- [D] B·R·A·H·M·S KRYPTOR compact PLUS Manual.

■ Symbols

Symbols used in Instruction for Use and Product Labelling of B·R·A·H·M·S KRYPTOR compact PLUS products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
Intended Use	Reference to the Intended use of the Medical Device	IVD	In Vitro Diagnostic Medical Device	LOT	Batch Code
CONT	Contents	CAL	Calibrator	CONTROL	Control
BUF	Buffer	SOLN 1	B·R·A·H·M·S KRYPTOR SOLUTION 1/B·R·A·H·M·S KRYPTOR compact SOLUTION 1	SOLN 2	B·R·A·H·M·S KRYPTOR SOLUTION 2/B·R·A·H·M·S KRYPTOR compact SOLUTION 2
SOLN 3	B·R·A·H·M·S KRYPTOR SOLUTION 3/B·R·A·H·M·S KRYPTOR compact SOLUTION 3	SOLN 4	B·R·A·H·M·S KRYPTOR SOLUTION 4/B·R·A·H·M·S KRYPTOR compact SOLUTION 4	CONT BAGS	Bags contained
BAGS	Bags	CONT PLATES	Plates contained	PLATES	Plates
CONT VIALS	Vials contained	VIALS	Vials	VIAL	Vial
H2O	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL	LYOPH	Lyophilized, Freeze Dried	RCNS	Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
®	Registered Trade Mark	TM	Trade Mark	REF	Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		GHS05, Corrosion
	GHS07, Exclamation Mark	TRACE	Trade Mark for TRACE™-technology	CE	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices
CE 0483	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg.Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste

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Bibliografija

Analitė

- [1] Morgenthaler NG, Struck J, Alonso C, Bergmann A. Assay for the measurement of copeptin, a stable peptide derived from the precursor of vasopressin. Clin Chem 2006; 52(1): 112-9
- [2] Holt NF and Haspel KL. Vasopressin: a review of therapeutic applications. J Cardiothorac Vasc Anesth 2010; 24(2): 330-47
- [3] Jochberger S, Zitt M, Luckner G, Mayr VD, Wenzel V, Ulmer H, Morgenthaler NG, Hasibeder WR, Dünser MW. Postoperative vasopressin and copeptin levels in noncardiac surgery patients: a prospective controlled trial. Shock 2009; 31(2): 132-8
- [4] Balanescu S, Kopp P, Gaskill MB, Morgenthaler NG, Schindler C, Rutishauser J. Correlation of plasma copeptin and vasopressin concentration in hypo-, iso-, and hyperosmolar states. J Clin Endocrinol Metab. 2011 Apr;96(4):1046-52.
- [5] Möckel M, Searle J, Hamm C, Slagman A, Blankenberg S, Huber K, Katus H, Liebetrau C, Müller C, Muller R, Peitsmeyer P, von Recum J, Tajsic M, Vollert JO, Giannitsis E. Early discharge using single cardiac troponin and copeptin testing in patients with suspected acute coronary syndrome (ACS): a randomized, controlled clinical process study. Eur Heart J. 2015 Feb 7; 36(6): 369–376
- [6] Post F, Gori T, Giannitsis E, Darius H, Baldus S, Hamm C, Hambrecht R, Hofmeister HM, Katus H, Perings S, Senges J, Münzel T. Criteria of the German Society of Cardiology for the establishment of chest pain units: update 2014. Clin Res Cardiol, 2015: DOI 10.1007/s00392-015-0888-2
- [7] Lipinski MJ, Escárcega RO, D'Ascenzo F, Magalhães MA, Baker NC, Torguson R, Chen F, Epstein SE, Miró O, Llorens P, Giannitsis E, Lotze U, Lefebvre S, Sebbane M, Cristol JP, Chenevier-Gobeaux C, Meune C, Eggers KM, Charpentier S, Twerenbold R, Mueller C, Biondi-Zoccai G, Waksman R. A systematic review and collaborative meta-analysis to determine the incremental value of copeptin for rapid rule-out of acute myocardial infarction. Am J Cardiol 2014;113:1581e1591
- [8] Möckel M, Searle J. Copeptin-marker of acute myocardial infarction. Curr Atheroscler Rep 2014: 16:421
- [9] Fenske W, Quinkler M, Lorenz D, Zopf K, Haagen U, Papassotiriou J, Pfeiffer AFH, Fassnacht M, Stoerk S, Allolio B. Copeptin in the differential diagnosis of the polydipsia-polyuria-syndrome — Revisiting the direct and indirect water deprivation tests. J Clin Endocrinol Metab 2011, 96(5): 1506–1515
- [10] Winzeler B, Zweifel C, Nigro N, Arici B, Bally M, Schuetz P, Blum CA, Kelly C, Berkman S, Huber A, Gentili F, Zadeh G, Landolt H, Mariani L, Müller B, Christ-Crain M. Postoperative Copeptin Concentration Predicts Diabetes Insipidus After Pituitary Surgery. J Clin Endocrinol Metab 2015; <http://dx.doi.org/10.1210/jc.2014-4527>
- [11] Timper K, Fenske W, Kühn F, Frech N, Arici B, Rutishauser J, Kopp P, Allolio B, Stettler C, Müller B, Katan M, Christ-Crain M. Diagnostic Accuracy of Copeptin in the Differential Diagnosis of the Polyuria-polydipsia Syndrome: A Prospective Multicenter Study. J Clin Endocrinol Metab. 2015 Jun; 100(6):2268-74
- [12] Szinnai G, Morgenthaler NG, Berneis K, Struck J, Müller B, Keller U, Christ-Crain M. Changes in plasma copeptin, the C-terminal portion of arginine vasopressin during water deprivation and excess in healthy subjects. J Clin Endocrinol Metab 2007; 92(10): 3973-8

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TRACE™

Date This version supersedes all earlier versions.
13.09.2011 Changes versus previous version:
• Information of TRACE patents deleted

1 Intended Use

B·R·A·H·M·S KRYPTOR compact Consumables is a set of accessories, necessary especially for the proper operation of the B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS. The intended use of each individual KRYPTOR compact Consumable is described with the specific component (chapter 4).

2 Instruments

B·R·A·H·M·S KRYPTOR compact Consumables may be used on:

- ☐ B·R·A·H·M·S KRYPTOR
- ☒ B·R·A·H·M·S KRYPTOR compact
- ☒ B·R·A·H·M·S KRYPTOR compact PLUS

3 Introduction

B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS is a fully automated system for B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS immunometric parameters. It is a closed system and can only operate utilizing special reagents and accessories offered by B·R·A·H·M·S GmbH.

4 Contents

B·R·A·H·M·S KRYPTOR BUFFER

BUF

Not supplied with the kit.

REF 89970

18...30 °C

see label for expiry date

Intended Use: B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS BUFFER is designed to prepare the B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS system liquid

Name	Quantity	Quality	Description
KRYPTOR BUFFER	BAGS 5 (65 g)	powder	potassium phosphate, sodium phosphate and sodium chloride

B·R·A·H·M·S KRYPTOR compact SOLUTION 1

SOLN 1

Not supplied with the kit.

REF 89981

2...8 °C

see label for expiry date

Intended Use: B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS SOLUTION 1 is designed to reconstitute the lyophilized conjugates from

B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS kits and to rinse the tip.

Name	Quantity	Quality	Description
KRYPTOR compact SOLUTION 1	VIALS 4 (55 mL)	ready for use	solution containing a preservative

B·R·A·H·M·S KRYPTOR compact SOLUTION 2

SOLN 2

Not supplied with the kit.

REF 89982

2...8 °C

see label for expiry date

Intended Use: B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS SOLUTION 2 is designed to reconstitute the lyophilized conjugates from B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS kits.

Name	Quantity	Quality	Description
KRYPTOR compact SOLUTION 2	VIALS 4 (55 mL)	ready for use	solution containing potassium fluoride

B·R·A·H·M·S KRYPTOR compact SOLUTION 3

SOLN 3

Not supplied with the kit.

REF 89983

2...8 °C

see label for expiry date


Intended Use: B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS SOLUTION 3 is designed for the washing of the B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS fluidic system.

Name	Quantity	Quality	Description
KRYPTOR compact SOLUTION 3	VIALS 4 (30 mL)	ready for use	solution containing sodium hydroxide and sodium hypochloride

B-R-A-H-M-S KRYPTOR compact SOLUTION 4**SOLN 4**

Not supplied with the kit.

REF 89984

2...8 °C

see label for expiry date


Intended Use: B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS SOLUTION 4 is designed for the washing of the B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS fluidic system.

Name	Quantity	Quality	Description
KRYPTOR compact SOLUTION 4	VIALS 4 (55 mL)	ready for use	solution containing sodium hydroxide.

B-R-A-H-M-S KRYPTOR compact DILCUP**PLATES**

Not supplied with the kit.

REF 89985

18...30 °C

not applicable


Intended Use: B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS DILCUP is designed for automatic dilution of samples.

Name	Quantity	Quality	Description
KRYPTOR compact DILCUP	PLATES 30	ready for use	dilution plates, containing 24 wells each

B-R-A-H-M-S KRYPTOR compact REACT**PLATES**

Not supplied with the kit.

REF 89986

18...30 °C

not applicable

Intended Use: B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS REACT is designed as a reaction plate for incubation of samples and reagents and for the fluorescent measurement.

Name	Quantity	Quality	Description
KRYPTOR compact REACT	PLATES 60	ready for use	reaction plates, containing 96 wells each. Additionally 96 adhesive cover films are included.

5 Precautions

	Respect the user information given in the User Manual of the B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS and in this document.
--	---

	<ul style="list-style-type: none"> Observe the expiration dates specified on the label. Avoid any microbial contamination of the reagents. The results, obtained from this assay should always be assessed in combination with the clinical examination, patients medical history and other findings before momentous actions will be prefaced.
--	--

	<p>For reagents containing Potassium Fluoride KF:</p> <p>Xn Harmful</p> <p>R20/21/22: Harmful by inhalation, in contact with skin and if swallowed.</p> <p>S22: Do not breathe dust.</p> <p>S24: Avoid contact with skin.</p>
--	--

	<p>For reagents containing Sodium Hypochlorite and Sodium Hydroxide:</p> <p>Xi Irritant</p> <p>R31: Contact with acids liberates toxic gas.</p> <p>R36/38: Irritating to eyes and skin.</p> <p>S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.</p>
--	--

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6 Instructions

The operation and maintenance of B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS is described in the B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS User Manual.

B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS BUFFER **BUF**

Paration for 5 L of liquid system:

- Completely empty the liquid system bottle,
- Pour in approximately 2 L of distilled water (conductivity of less than 50 µS/cm is recommended),
- Pour the contents of one B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS BUFFER bag into the liquid system bottle,
- Shake the bottle vigorously to ensure the salt dissolution,
- Top up the volume to 5 L with distilled water and mix,
- Replace the bottle on its support,
- Start up a prime of the fluidic system (see User Manual),
- After dissolution the buffer is stable for 15 days at 18...30 °C,
- For further information see the B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS User Manual.

B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS SOLUTIONS

- Before opening, store at 2...8 °C until the expiry date marked on the label.
- After opening all B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS SOLUTIONS are stable for 4 weeks at 18...30 °C.
- For further information see the B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS User Manual.

B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS REACT

- After loading to the system B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS REACT is stable for 24 hours at 18...30 °C.

B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS DILCUP

- After loading to the system B-R-A-H-M-S KRYPTOR compact/KRYPTOR compact PLUS DILCUP is stable for 3 days at 18...30 °C.

7 Quality Control

Not applicable

8 Reference Range

Not applicable

9 Assay Performance

Complementary information

For data on the assay performance of B·R·A·H·M·S KRYPTOR in-vitro diagnostic assays see the instruction for use of the individual product.

11 Symbols

Symbols used in Instruction for Use and Product Labelling of B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
Intended Use	Reference to the Intended use of the Medical Device	IVD	In Vitro Diagnostic Medical Device	LOT	Batch Code
CONT	Contents	CAL	Calibrator	CONTROL	Control
BUF	Buffer	SOLN 1	B·R·A·H·M·S KRYPTOR SOLUTION 1/B·R·A·H·M·S KRYPTOR compact SOLUTION 1	SOLN 2	B·R·A·H·M·S KRYPTOR SOLUTION 2/B·R·A·H·M·S KRYPTOR compact SOLUTION 2
SOLN 3	B·R·A·H·M·S KRYPTOR SOLUTION 3/B·R·A·H·M·S KRYPTOR compact SOLUTION 3	SOLN 4	B·R·A·H·M·S KRYPTOR SOLUTION 4/B·R·A·H·M·S KRYPTOR compact SOLUTION 4	CONT BAGS	Bags contained
BAGS	Bags	CONT PLATES	Plates contained	PLATES	Plates
CONT VIALS	Vials contained	VIALS	Vials	VIAL	Vial
H₂O	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL	LYOPH	Lyophilized, Freeze Dried	RCNS	Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
®	Registered Trade Mark	TM	Trade Mark	REF	Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		Harmful
	Irritant	TRACE	Trade Mark for TRACE™-technology	CE	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices
CE 0483	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg.Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste

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Instruction for Use

B·R·A·H·M·S MR-proADM KRYPTOR

IVD
for professional use
only
TRACE™

B·R·A·H·M·S is a registered trademark of B·R·A·H·M·S GmbH.

Other product names in this document are used for identification purposes; they may be trademarks and/or registered trademarks of their respective companies.

The manufacture and/or use of this product is covered by one or more of the following patents:

EP1488209, MX255044, JP4602321, CN1759319, HK1089232, JP5021791, US9541549, EP2347266, CN102317790, HK1163813, IN283666, JP5584695, US9012151, JP5711131, EP2097748.

This version supersedes all earlier versions.

Date

11.07.2018

Changes versus previous version

- KRYPTOR GOLD implemented
- Medical indications in section "Intended Use" rephrased for more clarity
- Section "Introduction" expanded for more clarity
- Statement on safety of raw materials of animal and human origin added under section "Precautions"
- Statement on method comparability added under section "Instructions"
- Freeze/Thawing cycles updated in the section "Instructions" to four times
- Existing reference range profile completed: Median of normal subjects updated to 0.38 nmol/L
- Value for Accuracy/ "High Dose Hook" Effect updated to 100 nmol/L
- Existing precision profile completed
- Re-ordering of Literature References

■ Intended Purpose/Intended Use

B·R·A·H·M·S MR-proADM KRYPTOR is an automated immunofluorescent assay for the quantitative determination of midregional proadrenomedullin (MR-proADM) in human EDTA plasma.

The B·R·A·H·M·S MR-proADM KRYPTOR assay is indicated for use as an aid, together with other clinical evaluation and laboratory diagnostic findings, to assess a patient's risk for an adverse outcome in the following indications or settings:

- lower respiratory tract infections (LRTI, comprising of community acquired pneumonia (CAP), chronic, obstructive pulmonary disease (COPD) and bronchitis),
- urinary tract infections (UTI) and kidney disease,
- sepsis, severe sepsis or septic shock,
- acute heart failure (AHF),
- chronic heart failure (CHF),
- myocardial infarction (MI),
- Emergency Department (ED) or Emergency Room (ER)

half-life time)^[2]. The mid-regional fragment of pro-Adrenomedullin (MR-proADM), consisting of the amino acid sequence 45-92, was identified as a stable and reliable surrogate marker for ADM release^[2]. The B·R·A·H·M·S MR-proADM KRYPTOR assay measures the mid-regional amino acid sequence 68-94 of the prohormone of active adrenomedullin^[3].

MR-proADM is a diagnostic parameter which can more accurately determine outcome and risk in patients with conditions such as lower respiratory tract infections^{[4][5]}, urinary tract infections and kidney disease^{[6][7]}, as well as in patients suffering from sepsis, severe sepsis and septic shock^{[8][9][10][11]}. Furthermore, it has been shown to be of value in predicting outcome and risk in forms of heart failure, such as acute^[12] and chronic heart failure^[13], as well as in myocardial infarction^[12]. It is also of use in assessing risk in the Emergency Department setting^[14].

■ Instruments

B·R·A·H·M·S MR-proADM KRYPTOR may be used on:

- ☒ B·R·A·H·M·S KRYPTOR
- ☒ B·R·A·H·M·S KRYPTOR compact
- ☒ B·R·A·H·M·S KRYPTOR compact PLUS
- ☒ B·R·A·H·M·S KRYPTOR GOLD

Availability of instrument depends on registration in the country.

■ Introduction

Adrenomedullin (ADM) is one of the most potent vasodilators, playing a significant role in microcirculation and endothelial (dys)function^[1]. The reliable measurement of ADM is not suitable for clinical routine diagnosis assessment due to its ex vivo instability (immediate binding to receptors, 22 min

■ Contents

Kit



B-R-A-H-M-S MR-proADM KRYPTOR CONT

REF 829.050  2 °C  8 °C  see label for expiry date

Name	Quantity	Quality	Description
CRYPTATE-CONJUGATE	VIAL 1 LYOPH		Anti-pro-ADM sheep polyclonal antibody conjugated with europium cryptate, buffer, bovine albumin, potassium fluoride
DILUENT	VIAL 1 ready for use		Newborn calf serum, EDTA, NaN ₃
XL665-CONJUGATE	VIAL 1 ready for use		Anti-pro-ADM sheep polyclonal antibody conjugated with XL665, buffer, bovine albumin, non-immunized mouse, sheep and bovine immunoglobulins, potassium fluoride.
Bar code card	1 ready for use		see the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS/ KRYPTOR GOLD User Manual. The bar code card contains all information needed for the registration of a new reagent lot.

Accessories



B-R-A-H-M-S MR-proADM KRYPTOR CAL CAL

Not supplied with the kit.
REF 82991  2 °C  8 °C  see label for expiry date

Intended Purpose/Intended Use: To readjust the standard curve memorized in B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/ KRYPTOR GOLD.

Name	Quantity	Quality	Description
B-R-A-H-M-S MR-proADM KRYPTOR CAL	VIALS 6 LYOPH		synthetic pro-ADM in human plasma
Bar code card	1 ready for use		see the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS/ KRYPTOR GOLD User Manual. The bar code card contains information related to the calibrator lot including its concentration.

B-R-A-H-M-S MR-proADM KRYPTOR QC CONTROL

Not supplied with the kit.
REF 82992  2 °C  8 °C  see label for expiry date

Intended Purpose/Intended Use: Quality control on board of B-R-A-H-M-S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD for the assay.

Name	Quantity	Quality	Description
B-R-A-H-M-S MR-proADM KRYPTOR QC Control 1	VIALS 3 LYOPH		synthetic pro-ADM in human plasma
B-R-A-H-M-S MR-proADM KRYPTOR QC Control 2	VIALS 3 LYOPH		synthetic pro-ADM in human plasma
Bar code card	1 ready for use		see the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS/ KRYPTOR GOLD User Manual. The bar code card contains information related to the control lot, particularly the target concentrations, the standard deviations obtained and the concentration acceptance ranges. This information is visible on the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS/ KRYPTOR GOLD monitor screen in the quality control section.
Bar code stick-on labels	32 for each control ready for use		see the B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS/ KRYPTOR GOLD User Manual. The bar code stick-on labels are used for identifying the controls when they are assayed on B-R-A-H-M-S KRYPTOR/KRYPTOR compact/ KRYPTOR compact PLUS/ KRYPTOR GOLD.

Other requisites

Not supplied with the kit.

B-R-A-H-M-S KRYPTOR Consumables

Name	REF
B-R-A-H-M-S KRYPTOR BUFFER	89970
B-R-A-H-M-S KRYPTOR SOLUTION 1	89971
B-R-A-H-M-S KRYPTOR SOLUTION 2	89972
B-R-A-H-M-S KRYPTOR SOLUTION 3	89973
B-R-A-H-M-S KRYPTOR SOLUTION 4	89974
B-R-A-H-M-S KRYPTOR DILCUP	89975
B-R-A-H-M-S KRYPTOR REACT	89976

B·R·A·H·M·S KRYPTOR compact/KRYPTOR compact PLUS Consumables

Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89981
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89982
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89983
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89984
B·R·A·H·M·S KRYPTOR compact DILCUP	89985
B·R·A·H·M·S KRYPTOR compact REACT	89986

B·R·A·H·M·S Consumables KRYPTOR GOLD

Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S Solution 1 KRYPTOR GOLD	89991
B·R·A·H·M·S Solution 2 KRYPTOR GOLD	89992
B·R·A·H·M·S Solution 3 KRYPTOR GOLD	89993
B·R·A·H·M·S Solution 4 KRYPTOR GOLD	89994
B·R·A·H·M·S Dilution Plates KRYPTOR GOLD	89995
B·R·A·H·M·S Reaction Plates KRYPTOR GOLD	89996

■ Precautions

	Follow the information given in the user manual of the B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD and in this document.
---	--

	<ul style="list-style-type: none"> – Check the expiration dates specified on the label. – Reagents from different reagent kits should not be mixed. – Do not separate the elements of a unit. – Avoid any microbial contamination of the reagents. – The results obtained from this assay should always be assessed in combination with the clinical examination, patients medical history and other findings before momentous actions will be prefaced.
--	---

	Raw materials of animal and/or human origin contained in the reagents have been tested and/or processed through virus inactivation procedures as far as possible. However, as it is impossible to strictly warrant that such products will not transmit any infectious agents, all raw materials of human or animal origin including the samples to be assayed must be treated as potentially infectious.
--	---

	The generally acknowledged safety precautions and laboratory techniques must be observed when handling reagents and patient samples.
	<ul style="list-style-type: none"> – Do not pipet by mouth. – Wash hands after work.

	– Wear protective clothing, protective gloves and safety glasses at work.
	

	– Do not eat, drink or smoke in areas where samples or kit-reagents are handled.
	

	<ul style="list-style-type: none"> – Remove spills with absorbing paper. – All the material used for cleaning up must be disposed of as infectious laboratory waste. – Prevent from getting into sewage, water, ground.
---	--

	<ul style="list-style-type: none"> – Used reagent plates and reagent kits are disposed of as potential infectious laboratory waste according to local regulations. – Empty containers could be recycled.
---	--

■ Principle

The measurement principle of B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD is based on TRACE™ Technology (Time-Resolved Amplified Cryptate Emission), which measures the signal that is emitted from an immunocomplex with time delay. The basis of the TRACE™ Technology is non-radiative energy transfer from a donor (a cage-like structure with an europium ion in the center [cryptate]) to an acceptor, which is part of a chemically modified, light-collecting algal protein (XL 665).

The proximity of donor (cryptate) and acceptor (XL 665) when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the cryptate and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.

Precise measuring of analyte concentration: When the sample is excited with a nitrogen laser at 337 nm, the donor (cryptate) emits a long-life fluorescent signal in the millisecond range at 620 nm, while the acceptor (XL 665) generates a short-life signal in the nanosecond-range at 665 nm.

When the two components are bound in an immunocomplex, both the signal amplification and the prolongation of the life span of the acceptor signal occur at 665 nm, so that it can be measured over microseconds.

Reliable prevention of interference: Non-specific signals, e.g. the signals of the short-life and unbound acceptor XL 665 and the medium-specific interference signals conditional upon the natural fluorescence of the sample, are eliminated by temporal delay of the fluorescence measurement. The signal generated by the cryptate at 620 nm serves as an internal reference and is measured simultaneously with the long-life acceptor signal at 665 nm which is the specific signal. Interfering influences, e.g. from turbid sera, are automatically corrected by means of the internally calculated ratio of the intensities at these wavelengths.

■ Instructions

Sample volume.....	26 µL
Incubation time.....	29 min
Results are given in.....	nmol/L
Conversion factor.....	not applicable
Direct measuring range.....	0.05...10 nmol/L
Measuring range with automatic dilution.....	0.05...100 nmol/L
Sample type.....	EDTA plasma
Kit stability on board.....	2 weeks
Calibrator.....	1 point
Calibration stability.....	8 days
Assay principle.....	sandwich

- If the assay is carried out within 24 hours of blood collection, store samples at room temperature or alternatively at 2...8 °C. If not, samples shall be aliquoted and stored frozen (-20 °C).
- Samples may be frozen and thawed four times.
- Other sample types than recommended should not be used.
- Place the sample in a tube suited for use on B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD (11–17 mm diameter). This might be the primary tube.
- The sample tube must contain a minimum volume which will vary depending on the diameter of the sample tube. A 13 mm diameter tube will require an additional 150 µL of sample.
- Should a dilution be requested either automatically or by the user, the volume of sample necessary will be 26 µL maximum.
- Icteric, hemolytic or hyperlipemic samples, or samples which are turbid or contain fibrin may yield imprecise results. Such samples are signaled by B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD.
- The concentration of MR-proADM in a given patient sample, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity.


The operation and maintenance of B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD are described in the related user manual.

After it has been opened the reagent unit may be stored on the B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD in the space provided. Each reagent unit is individually identified (bar code) and its maximum period of use after opening is controlled by the instrument.

To prepare a reagent unit, proceed as follows:

- Remove the security band from the reagent pack.
- Push in the lid by pressing it firmly (see diagram below).

Make sure to remove all the aluminium foil.

Opening the Kit

Note: Always handle the reagent unit carefully to avoid foaming or bubbles. Foam/bubbles may interfere with proper detection of reagent and/or correct dispensing of reagents. Allow the kit to stand 5...10 min within B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD prior to use.

A standard curve does not need to be established on B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD.

The instrument stores the required information after reading the bar code from the reagent card. A recalibration must be carried out for every new reagent lot, and subsequently on a regular basis. B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD automatically indicate when a calibration is required. Both, the stored and recalibrated standard curve may be displayed on the screen.

The following steps are carried out:

- Conjugates and sample are dispensed into the reaction plate and the fluorescent signal emitted is measured periodically.
- Samples with concentrations higher than the direct measurement range, are identified in the first few minutes of incubation, then diluted automatically and reassayed.
- After measurement of the fluorescent signal, the data obtained from the software are compared to the memorized standard curve.

Calibration **CAL**

- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) indicated on the vial label.
- Mix gently after reconstitution. Allow 15 min for the complete dissolution of the lyophilisate.
- Calibration must be carried out with every new reagent kit lot. It is then repeated on a regular basis automatically managed by B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD in order to readjust the standard curve.
- It is recommended that controls be run after each calibration.
- Use the calibrator/calibrators only once.
- Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5 hours.
- The calibrator bar code card must be read in for each new lot of calibrator.
- For further information see the B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD user manual.

Control **CONTROL**

- It is recommended that controls be run at least once a day and after each calibration.
- A control tube is processed like a sample tube.
- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 µS/cm is recommended) as stated on the vial label.
- Allow 15 min. for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- After reconstitution, do not keep a vial more than 5 hours at 18...25 °C or 24 hours at 2...8 °C.
- It is recommended that the contents of a reconstituted vial be divided into aliquots which may then be stored frozen at -20 °C for a maximum period of 1 month.
- Use one of the tubes immediately for measurement.
- After thawing an aliquot, mix gently and use immediately for measurement.
- Once thawed, a control aliquot must not be refrozen.
- The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD.
- The control kit bar code card must be scanned in for each new lot of control.
- For further information see the B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD user manual.

■ Quality Control

Good laboratory practice requires that control samples are measured regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that the results be analysed using appropriate statistical methods.

If desired, B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS can check assay quality at intervals, by statistical analysis on the basis of Levey Jennings graphs.

It is necessary to comply with national quality assurance guidelines for quantitative tests in the medical laboratory (**current version**). For instance, test accuracy and precision should be monitored by means of laboratory in house and/or commercially available control materials. If unacceptable control values are obtained, proceed as outlined in standard laboratory diagnostic procedures to determine the cause and implement corrective measures.

■ Reference Range

Normal subjects

From the data available, the following reference table can be drawn up for the assessment of mid-regional pro-ADM concentrations in healthy subjects (n=144, 87 male, 57 female).

Normal subjects	mid-regional pro-ADM [nmol/L]
97.5 % percentile	0.55 nmol/L
Median	0.38 nmol/L

MR-proADM is a strong prognosticator in patients with acute congestive heart failure and in patients presenting with dyspnea in general; furthermore it is particularly strong in predicting short-term prognosis within 30 days after assessment. Based on recent publications a calculated cut off >1.5 nmol/L indicates an individual high risk concerning severity of disease and prognosis.^{[5][11][13]}

Note: It is recommended that each laboratory establish its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values above have been established on EDTA plasma and are given only as indications.

The result obtained from this assay should always be assessed in combination with the clinical examination, patients medical history and other findings before momentous actions will be prefaced.

■ Assay Performance

Detection Limit

The detection limit, calculated using the imprecision profile, has been assessed as being 0.05 nmol/L.

CLSI EP 17-A

Sensitivity

The **functional assay sensitivity (FAS)**, detected by inter-assay precision of 20 % CV, has been assessed as being 0.25 nmol/L.

The **limit of quantitation (LOQ)**, is estimated in accordance with CLSI guideline EP 17-A (Evaluation of Precision Performance of Clinical Chemistry Devices). It has been assessed as being 0.23 nmol/L.

Specificity

- No cross reactivity of
- cTNI up to 400 ng/mL
 - Myoglobin up to 852 ng/mL
 - CKMB up to 204 ng/mL
 - PCT up to 50.2 ng/mL
 - MR-proANP up to 905 pmol/L
 - CT-proET1 up to 409 pmol/L
 - CRP up to 110 µg/mL

Accuracy/ Linearity

Dilution of highly concentrated samples revealed recovery rates between 80 % and 120 %.

CLSI EP 6-A

Accuracy/ "High Dose Hook" Effect

No "High Dose Hook" Effect can be observed up to 100 nmol/L.

Precision/ Reproducibility for Intra Assay CV

This was evaluated using samples of different concentrations in accordance with the CLSI EP05-A3 guideline (Evaluation of Precision of Quantitative Measurement Procedures).

	Mean [nmol/L]	Within-assay CV %
Sample 1	0.282	10.8
Sample 2	1.63	3.1
Sample 3	2.20	1.2
Sample 4	7.74	3.0
Sample 5	120	6.3

Precision/ Reproducibility for Inter Assay CV

This was evaluated using samples of different concentrations in accordance with the CLSI EP5-A3 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods).

	Mean [nmol/L]	Within-assay CV %
Sample 6	0.282	17.5
Sample 7	0.85	10.4
Sample 8	3.07	7.3
Sample 9	7.35	5.6
Sample 10	120	6.8

Disturbing Factors

Factor	Description
Hemoglobin	No significant effect up to 5 g/L
Bilirubin	No significant effect up to 400 mg/L
Triglycerides	No significant effect up to 10 g/L
Rheumatoid factor	No significant effect up to 1 100 IU/mL
Albumin	No significant effect up to 6 g/dL

No interference of

Factor	Description
Imipenem	No significant effect up to 1.18 mg/mL
Heparin	No significant effect up to 8000 U/L
Furosemide	No significant effect up to 2 mg/dL
Tiotropium	No significant effect up to 21.6 ng/mL
Albuterol	No significant effect up to 1.67 µmol/L
Salmeterol	No significant effect up to 60 ng/mL
Theophylline	No significant effect up to 228 µmol/L
Montelukast	No significant effect up to 6 µg/mL
Epinephrin	No significant effect up to 1.8 µg/mL
Terbutaline	No significant effect up to 0.9 µg/mL
Ipratropium bromide	No significant effect up to 0.9 µg/mL
Methylprednisolone	No significant effect up to 72 µg/mL
Prednisolone	No significant effect up to 8.31 µmol/L

Factor	Description
Prednisone	No significant effect up to 0.84 µmol/L
Cromolyn	No significant effect up to 480 µg/mL
Nedocromil	No significant effect up to 8.4 µg/mL
Beclomethasone dipropionate	No significant effect up to 1 µg/mL
Budesonide	No significant effect up to 0.72 µg/mL
Flunisolide	No significant effect up to 2.4 µg/mL
Fluticasone	No significant effect up to 0.3 µg/mL
Formoterol	No significant effect up to 28.8 ng/mL
Cefotaxim	No significant effect up to 90 mg/dL
Vancomycin	No significant effect up to 3.5 mg/mL
Dopamine	No significant effect up to 13 mg/dL
Dobutamine	No significant effect up to 11.2 µg/mL
Triamcinolone	No significant effect up to 2.4 µg/mL

CLSI EP 7-A2

Traceability

– not applicable –

■ Bibliography

Analyte

- [1] T. Eto. A review of the biological properties and clinical implications of adrenomedullin and proadrenomedullin N-terminal 20 peptide (PAMP), hypotensive and vasodilating peptides. *Peptides* 2001;22 (11): 1693-711.
- [2] J. Struck, C. Tao, N. G. Morgenthaler and A. Bergmann. Identification of an Adrenomedullin precursor fragment in plasma of sepsis patients. *Peptides* 2004;25 (8): 1369-72.
- [3] Caruhel P, et al. (2009): Homogeneous time-resolved fluoroimmunoassay for the measurement of midregional proadrenomedullin in plasma on the fully automated system B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD Clin Biochem 42: (7-8); 725-8.
- [4] Kruger S, Ewig S, Giersdorf S, Hartmann O, Suttrop N, Welte T. Cardiovascular and inflammatory biomarkers to predict short- and long-term survival in community-acquired pneumonia: Results from the German Competence Network, CAPNETZ. *Am J Respir Crit Care Med.* Dec 1 2010;182(11):1426-1434.
- [5] Renaud B, Schuetz P, Claessens YE, Labarere J, Albrich W, Mueller B. Proadrenomedullin improves Risk of Early Admission to ICU score for predicting early severe community-acquired pneumonia. *Chest.* Dec 2012;142(6):1447-1454.
- [6] Gouya G, Sturm G, Lamina C, et al. The association of mid-regional pro-adrenomedullin and mid-regional pro-atrial natriuretic peptide with mortality in an incident dialysis cohort. *PLoS One.* 2011;6(3):e17803.
- [7] Litke A, Bossart R, Regez K, et al. The potential impact of biomarkerguided triage decisions for patients with urinary tract infections. *Infection.* Aug 2013;41(4):799-809.
- [8] Abraham E, Laterre PF, Garg R, et al. Drotrecogin alfa (activated) for adults with severe sepsis and a low risk of death. *N Engl J Med.* Sep 29 2005;353(13):1332-1341.
- [9] Christ-Crain M, Morgenthaler NG, Struck J, Harbarth S, Bergmann A, Muller B. Mid-regional pro-adrenomedullin as a prognostic marker in sepsis: an observational study. *Crit Care.* 2005;9(6):R816-824.
- [10] Guignant C, Voirin N, Venet F, et al. Assessment of pro-vasopressin and pro-adrenomedullin as predictors of 28-day mortality in septic shock patients. *Intensive Care Med.* Nov 2009;35(11):1859-1867.
- [11] Hagag AA, Elmahdy HS, Ezzat AA. Prognostic value of plasma proadrenomedullin and antithrombin levels in neonatal sepsis. *Indian Pediatr.* Jun 2011;48(6):471-473.
- [12] Klip IT, Voors AA, Anker SD, et al. Prognostic value of mid-regional pro-adrenomedullin in patients with heart failure after an acute myocardial infarction. *Heart.* Jun 2011;97(11):892-898.
- [13] Adlbrecht C, Hulsmann M, Strunk G, et al. Prognostic value of plasma midregional pro-adrenomedullin and C-terminal-pro-endothelin-1 in chronic heart failure outpatients. *Eur J Heart Fail.* Apr 2009;11(4): 361-366.
- [14] Nickel CH, Messmer AS, Geigy N, et al. Stress markers predict mortality in patients with nonspecific complaints presenting to the emergency department and may be a useful risk stratification tool to support disposition planning. *Acad Emerg Med.* Jul 2013;20(7): 670-679.

Technology

- [A] Mathis G, Lehn JM.: Trace – Another Story of Time. *Isotopics* 1995; Vol. 9.
- [B] Mathis G. , *Clin. Chem.* 1993; 39:1953-9.
- [C] B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD User Manual.
- [D] B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/KRYPTOR GOLD Manual.

■ Symbols

Symbols used in Instruction for Use and Product Labelling of B·R·A·H·M·S KRYPTOR/KRYPTOR compact/KRYPTOR compact PLUS/ KRYPTOR GOLD products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
Intended Use	Reference to the Intended use of the Medical Device	IVD	In Vitro Diagnostic Medical Device	LOT	Batch Code
CONT	Contents	CAL	Calibrator	CONTROL	Control
BUF	Buffer	SOLN 1	B·R·A·H·M·S KRYPTOR SOLUTION 1/B·R·A·H·M·S KRYPTOR compact SOLUTION 1/ B·R·A·H·M·S Solution 1 KRYPTOR GOLD	SOLN 2	B·R·A·H·M·S KRYPTOR SOLUTION 2/B·R·A·H·M·S KRYPTOR compact SOLUTION 2/ B·R·A·H·M·S Solution 2 KRYPTOR GOLD
SOLN 3	B·R·A·H·M·S KRYPTOR SOLUTION 3/B·R·A·H·M·S KRYPTOR compact SOLUTION 3/ B·R·A·H·M·S Solution 3 KRYPTOR GOLD	SOLN 4	B·R·A·H·M·S KRYPTOR SOLUTION 4/B·R·A·H·M·S KRYPTOR compact SOLUTION 4/ B·R·A·H·M·S Solution 4 KRYPTOR GOLD	CONT BAGS	Bags contained
BAGS	Bags	CONT PLATES	Plates contained	PLATES	Plates
CONT VIALS	Vials contained	VIALS	Vials	VIAL	Vial
H2O	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL	LYOPH	Lyophilized, Freeze Dried	RCNS	Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
®	Registered Trade Mark	TM	Trade Mark	REF	Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		GHS05, Corrosion
	GHS07, Exclamation Mark	TRACE	Trade Mark for TRACE™-technology	CE	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices
CE 0483	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg.Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste