## Magtrace® skirtas naudoti su Sentimag®

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# **Magtrace® for use with Sentimag®**

The information contained in this manual is subject to change without notice. Please visit [www.endomag.com](http://www.endomag.com) for the latest version.

## Intended Purpose

Magtrace® is a magnetic tracer intended to mark and locate lymph nodes in cancer patients prior to their surgical removal as part of Sentinel Lymph Node Biopsy (SLNB) procedure in adult populations. Magtrace® can be used as an SLNB marker wherever the physician deems SLNB appropriate. Within the adult population the use of Magtrace® is not restricted to the type of cancer, nor is it restricted by age or ethnicity within the adult population.

## Intended Patient Target Group/s

Magtrace® is intended for all cancer patients in the adult population for whom SLNB is deemed appropriate by the physician.

## Clinical Benefits

Magtrace® is a magnetic tracer. The magnetic signal from Magtrace® is induced by the Sentimag® device. Magtrace® does not rely for its effectiveness on radioisotopes. This means that all hospitals, many of which do not have nuclear medicine departments can use Magtrace® to make SLNB available to the patient.

## Indications for Use

Magtrace® is a magnetic tracer intended and calibrated for use with the Sentimag® device ONLY as part of an overall system to mark and locate lymph nodes in adult cancer patients prior to their surgical removal. The use of Magtrace® is not restricted to the type of cancer, provided the use of the Sentimag probe is viable for the particular cancer type.

According to its Intended Use, Magtrace® must be detected within 30 days after injection. In Europe, Magtrace® is classified as a class IIa medical device by Medical Device Directive 93/42/EEC.

## Description

Magtrace® is a blackish-brown sterile aqueous suspension of superparamagnetic carboxydextran-coated iron oxide particles in injectable water, supplied aseptically in single-use vials containing circa 2.2 ml. Each millilitre of Magtrace® contains circa 28 milligrams of iron.

## Dosage and Administration

The maximum dose is 2.0ml, with an equivalent iron content of 55 mg +/- 4 mg per dose. Inspect the seal of the vial before use to ensure it is unbroken. Do not use if the vial cap is broken, the vial is leaking, or if the expiry date has passed.

Draw Magtrace® via a sterile needle into an appropriate sterile hypodermic syringe and check the quantity.

For Breast:

**If using intra-operatively, administer 2.0ml Magtrace®** by subcutaneous injection into either subareolar or peritumoral interstitial tissue and follow with 5 minutes vigorous massage at the injection site. For subareolar injection, surgeons should wait at least **20 minutes** before attempting transcutaneous measurement of the axilla. Peritumoral injection may require a longer wait.

**If injecting pre-operatively (i.e. at least the day before surgery), administer 2.0ml Magtrace® or 1.0ml Magtrace®** into either subareolar or peritumoral interstitial tissue. There is no need to massage with pre-operative injection. The Magtrace® vial must only be used for a single administration to the patient. This single administration may be performed in multiple small injection volumes totaling 2.0ml Magtrace® or 1.0ml Magtrace® from a single volume drawn from the Magtrace® vial without the needle leaving the sterile field. The Magtrace® vial must be disposed of after drawing the required volume from the vial (refer to Disposal section below)

Proceed with incision only after obtaining a **clear transcutaneous** signal with the Sentimag®. Migration may be slower in older or larger patients.

For other indications: Administer Magtrace® by injection into peritumoral tissue and, where appropriate, follow with 5 minutes massage at the injection site. Proceed as above.

## Demographics of the Clinical studies

Magtrace® has been studied in populations with a range of ages (27-89 years) and BMI (17.1-51.7) reported. No limitations in age or BMI were reported.

### 1.1.1. tech. parametrai

## **Contraindications**

1. Magtrace® is contraindicated in any patient with hypersensitivity to iron oxide or dextran compounds.
2. Do not administer to any patient with iron overload disease.
3. Do not administer to any patient with a metal implant close to the expected sentinel lymph node location.

## **Warnings and Precautions**

1. Magtrace® is intended ONLY for use with the Sentimag® device, and is therefore subject to the Warnings and Precautions of the Sentimag® device, including the precaution regarding the use in patients with pacemakers.
2. For use ONLY by appropriately trained clinicians for detection of sentinel lymph nodes during surgical procedures.
3. Magtrace® can alter magnetic resonance imaging (MRI) studies of the injection and drainage site. Some amount of alteration may be long-term. Because of this, surgeons should consider whether Magtrace® is appropriate for a patient on a case by case basis.
4. Magtrace® is not intended for treatment of iron deficiency anaemia in patients or any other medicinal applications.
5. If inadvertently administered intravenously, anaphylactoid or cardiovascular reactions may occur.
6. Some transient or long-term brownish skin coloration may occur.
7. Do not reuse, sterility cannot be guaranteed if the rubber seal on the vial has already been punctured. Reuse of an already opened vial can lead to cross contamination and cross infection between patients.
8. There have been no studies of Magtrace® in pregnant women, nursing mothers or paediatric patients.

## **Adverse Reactions**

When similar material to that used in Magtrace® has been injected **intravenously**, the following undesirable effects have been reported: Common (<2%) – pain at the injection site, vasodilation, paresthesia. Uncommon (≥0.1% to <1%) – asthenia, back pain, injection site reactions, chest pain, nausea, vomiting, headache, taste perversion, pruritus, rash. Rare (≥0.01% to <0.1%) – hypersensitivity & anaphylaxis, hypertension, phlebitis, hyperesthesia, anxiety, dizziness, convulsion, parosmia, dyspnea, increased cough, rhinitis, eczema, urticaria. There have been a

small number of reports of inflammatory and hypersensitivity response with intradermal injection. There is no evidence of adverse reaction following interstitial injection. **Magtrace® is not intended for intravenous injection.**

## Overdose

Overdose is unlikely if used as specified with a single 2 ml volume of Magtrace® is administered.

## Non-Clinical Toxicology and Clinical Safety

Magtrace® has been reviewed and tested as specified in EN 10993-1:2009 based on the specified site of injection and duration. Detailed non-clinical toxicology of materials similar to that used in Magtrace® including human pharmacodynamics and pharmacokinetics studies shows that the material and its potential breakdown products are safe.

## Mechanism of Action

The lymphatic system comprises a network of vessels that drain tissue fluid (lymph) including the injected Magtrace® into lymph nodes and then the specialized organs involved in the immune system. The lymph nodes act as a filter, removing invading organisms or abnormal cells and other materials such as bacteria and cancer cells from the lymph fluid. Magtrace® particles are also trapped in the filter. They therefore become physically trapped in the node, allowing them to be used as a superparamagnetic marker which can be identified by the Sentimag® device.

The Magtrace® particles are trapped by the natural physical filtration action of the lymph node and are not part of a medicinal action.

## End user Profile and Training

Magtrace® is for use ONLY by appropriately trained clinicians for detection of sentinel lymph nodes during surgical procedures. Additionally, clinicians MUST be trained in the use of Magtrace® and Sentimag® system before use.

## Magnetic Fields

When the Magtrace® material is exposed to the excitation field of the Sentimag® the Magtrace® material responds with a temporarily induced magnetic field.

## Storage

Store between +2°C and +30°C and such that temperature variations will be minimised.

DO NOT FREEZE.

Do not use after expiry date specified on vial.

## Disposal

Magtrace® should not be disposed via public recycling.

- Incineration of contaminated parts is recommended.
- local legal requirements should be followed at all times.

## Patient Counselling

1. Question patient regarding any prior history of reactions to iron oxide or dextran products.
2. Inform patient to report any signs and symptoms of hypersensitivity that may develop during and/or following administration, such as rash, itching, dizziness, lightheadedness.
3. Advise patient that some long-term brownish skin coloration may occur.
4. Advise patient that Magtrace® may alter post-operative magnetic resonance imaging (MRI) scans. Some alteration may be long-term.

Any serious incident that has occurred in relation to the device should be reported to Endomagnetics Ltd using the following email address:

- [info@endomag.com](mailto:info@endomag.com)

## Symbols applicable to the Magtrace® medical device

Single Use

Expiry Date specified on vial

Lot or batch number specified on vial

Catalogue Number

Aseptically Filled


Read Instructions

Warnings and Cautions specified in instructions

Do not use if vial is open or damaged

Store between temperatures indicated

Manufacturer

CE mark for Medical Device as specified by the Medical Device Directive 93/42/EEC. 2797 signifies the approving Notified Body.

For Use by, or on the Order of, a Physician

MR Conditional

European Authorized Representative

Endomagnetics Ltd  
The Jeffreys Building  
Cowley Road, Cambridge  
CB4 0WS, United Kingdom

Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

Sentimag® is a registered European Union Trademark of Endomagnetics Ltd.  
Magtrace® is a registered Trademark of Endomagnetics Ltd in the United Kingdom.