
# Endomagnetics SentiMag sistema

## Aptarnavimo vadovas

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Reason For Return:

Please tick the appropriate option(s) below:

☐ I certify that I have decontaminated the product as per above.

Decontaminant Used:

☐ I certify that the product has not been exposed to any chemical or biological materials.

Title:

Name:

Signature:

Date:

Telephone:

Email:

**NOTE:** Please include a copy of this form with the product being returned.


Endomagnetics Ltd. incorporated in England and Wales – No: 06227698, VAT: 947 7709 68

Correspondence: Endomagnetics Ltd., 23 Cambridge Science Park, Milton Road, Cambridge, CB4 0EY, United Kingdom

Registered Office: Endomagnetics Ltd., Tower 42, 33rd Floor, 25 Old Broad Street, London, EC2N 1HQ, United Kingdom

43 psl

# Endomagnetics SentiMag System Service Manual

© March 2013. The information contained in this manual is the property of Endomagnetics, 325 Cambridge Science Park, Milton Road, Cambridge, CB4 0WG, UK, and may not be reproduced or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, nor stored in any retrieval system of any nature without the express written authority of Endomagnetics.

## **IMPORTANT**

IT IS ESSENTIAL THAT THE USER OF THIS MANUAL IS AWARE OF THE POTENTIAL HAZARDS ASSOCIATED WITH THE UNIT.

ALL SERVICE OPERATORS SHOULD BE FAMILIAR WITH THE SAFETY PRECAUTIONS AND WARNINGS GIVEN IN THE TRAINING COURSE PRIOR TO ATTEMPTING TO SERVICE THE UNIT.

IF THE UNIT IS SERVICED IN A MANNER NOT SPECIFIED BY THIS MANUAL, THE CORRECT FUNCTION MAY BE IMPAIRED.

### **INTRODUCTION**

This manual is intended for use by an authorised, competent service engineer or technician. It contains information on fault diagnosis, spares kits and the serving process for the Endomagnetics SentiMag System.

Removal and assembly of parts is supported by individual part lists and instructions.

Further technical support can be obtained directly from:

Integrated Technologies Ltd  
Viking House  
Ellingham Way  
Ashford  
Kent  
TN23 6NF  
UK

Tel: +44 (0)1233 638383

Fax: +44 (0)1233 639401

## **Table of Contents**

<b>1. SERVICE TOOLS REQUIRED</b>	<b>1</b>
<b>2. FAULT DIAGNOSIS OVERVIEW</b>	<b>2</b>
2.1 Cosmetic damage; casework, knobs/button, feet.	3
2.2 Probe related faults; not detected, abnormal readings, signal drift	4
2.3 Power related faults	6
2.4 Display related faults	7
2.5 Sound related faults	8
2.6 IV Pole Boss	9
2.7 Chassis	9
<b>3. SERVICING</b>	<b>10</b>
3.1 Servicing a unit with damaged volume and/or sense knobs	10
3.1.1 Spares kit required	10
3.1.2 Servicing process	10
3.2 Servicing a unit with damaged rear casework	11
3.2.1 Spares kit required	11
3.2.2 Servicing process	12
3.3 Servicing a unit with damaged feet	13
3.3.1 Spares kit required	13
3.3.2 Servicing process	13
3.4 Servicing a unit with damaged front casework	14
3.4.1 Spares kit required	14
3.4.2 Servicing process	15
3.5 Servicing a unit with damaged connector bracket, including mains inlet	15
3.5.1 Spares kit required	15
3.5.2 Servicing process	16
3.6 Servicing a unit with damaged IV pole boss	18
3.6.1 Spares kit required	18
3.6.2 Servicing process	19
3.7 Servicing a unit with damaged PSU PCB	20
3.7.1 Spares kit required	20
3.7.2 Servicing process	20
3.8 Servicing a unit with damaged Main Arm PCB	22

3.8.1 Spares kit required	22
3.8.2 Servicing process	23
<b>3.9 Servicing a unit with damaged Pressure Switch</b>	<b>24</b>
3.9.1 Spares kit required	24
3.9.2 Servicing process	25
<b>3.10 Servicing a unit with damaged Coil Interface PCB</b>	<b>26</b>
3.10.1 Spares kit required	26
3.10.2 Servicing process	27
<b>3.11 Servicing a unit with damaged speaker</b>	<b>28</b>
3.11.1 Spares kit required	28
3.11.2 Servicing process	29
<b>3.12 Servicing a unit with damaged LCD ARM PCB</b>	<b>30</b>
3.12.1 Spares kit required	30
3.12.2 Servicing process	31
<b>3.13 Servicing a unit with damaged LCD</b>	<b>33</b>
3.13.1 Spares kit required	33
3.13.2 Servicing process	33
<b>3.14 Servicing a unit with damaged Chassis</b>	<b>34</b>
3.14.1 Spares kit required	34
3.14.2 Servicing process	35
<b>4. DOCUMENT HISTORY</b>	<b>36</b>
<b>5. APPENDICES</b>	<b>37</b>
<b>5.1 Appendix 1 – Wiring schematic for the SentiMag System</b>	<b>37</b>
<b>5.2 Appendix 2 – Wiring diagrams</b>	<b>38</b>
5.2.1 775 Cable PSU PCB Mains power	38
5.2.2 775 Cable LCD PCB digital	38
5.2.3 775 Cable LCD PCB backlit	38
5.2.4 775 Cable LCD PCB audio	38
5.2.5 775 Cable Main PCB digital	39
5.2.6 775 Cable Main PCB analog	39
5.2.7 775 Cable Coil PCB driver	39
5.2.8 775 Cable LCD PCB download USB	40
5.2.9 775 Cable LCD PCB speaker	40
5.2.10 775 Cable LCD PCB comms RXD	40
5.2.11 775 Cable LCD PCB comms TXD	41
5.2.12 775 Cable LCD PCB footswitch	41
5.2.13 775 Cable Main PCB debug USB	41
5.2.14 775 Cable Coil PCB sense	42
5.2.15 775 Cable Coil PCB probe ok	42
5.2.16 775 Cable Coil PCB driver sig	42
5.2.17 775 Cable Coil PCB sense sig	42
<b>5.3 Appendix 3 – Decontamination</b>	<b>43</b>

## 1. Service Tools Required

For a complete service of the Endomagnetics SentiMag System the following tools are required:

For removing and fitting parts:

- Pozidrive screwdriver (Long shaft)
- Pozidrive screwdriver (Short shaft)
- Small flat head screwdriver
- M10 nut runner
- 19mm Spanner
- 16mm socket with extension, or nut runner
- Adjustable mole grips
- Anti-static wrist strap and bonding plug

For testing parts:

- DVM (Digital Volt Meter)
- Spare mains power lead
- Spare standard fuse

**Note:**

Units and probes which are used in the field may be biohazardous and must be correctly cleaned before handled or moved into a clean area. Before servicing a SentiMag System ensure that it is safe to work on and has been correctly decontaminated. See Appendix 3, section 5.3.


During factory manufacture, each SentiMag System is insulation and continuity tested to ensure that the unit is electrically insulated and correctly earthed. When serving a unit in the field it is not always appropriate to perform insulation and continuity testing. This service manual does not give instructions on how to carry out such testing but the service engineer/technician should ensure that the unit is correctly earthed and reassembled after servicing a unit.

An anti-static wrist band should be worn when handling PCBs and other sensitive electronic equipment (LCD & PSU).

## 2. Fault Diagnosis Overview

Top level overview of areas where faults may occur and servicing may be required:

Note: Pins in the probe connectors which plug into the unit may bend if a user has tried to plug the connector in at an angle. If pins are slightly bent out of position then the probe may not plug into a unit correctly or make good connection. A flat head screwdriver can be used with a small amount of force to straighten them. For all other probe related issues, contact the manufacturer as return to base for repair may be required.

**2.1 Cosmetic damage; casework, knobs/button, feet.**

The casework, knobs/button, and/or feet may be aesthetically damaged or damaged in a way which affects function:

**Front casework:**

- Scratched or broken moulding
- Scratched or damaged window
- Damaged or missing balance push button

Replace with front casework spares kit from ITL - Part No. S126197 (775 SPARE FRONT CASEWORK).

**Rear casework:**

- Scratched or broken moulding

Replace with rear casework spares kit from ITL - Part No. S126196 (775 SPARE REAR CASEWORK).

**Volume and sense knobs:**

- Damaged or missing knobs


Replace with volume and sense knobs spares kit from ITL - Part No. S126195 (775 SPARE VOL & SEN KNOBS).


**Feet:**

- Damaged or missing feet

Replace with feet spares kit from ITL - Part No. S126198 (775 SPARE FEET).


## 2.2 Probe related faults; not detected, abnormal readings, signal drift



## 2.3 Power related faults

## 2.4 Display related faults

## 2.5 Sound related faults

## 2.6 IV Pole Boss

The IV Pole Boss may be damaged in a way which prevents a unit being attached to an IV Pole:

- Cross threaded.

Replace with IV pole boss spares kit from ITL - Part No. S126200 (775 SPARE IV POLE BOSS).

## 2.7 Chassis

The chassis may be damaged in a way which prevents components in the unit being held together correctly. This may occur if a unit which has been dropped and as a result the chassis has twisted/buckled:

- Twisted or buckled chassis.
- Cross threaded holes caused by incorrect servicing.

Replace with chassis spares kit from ITL - Part No. S126206 (775 SPARE CHASSIS).

### 3. Servicing

#### 3.1 Servicing a unit with damaged volume and/or sense knobs

##### 3.1.1 Spares kit required

Order the volume and sense knobs spares kit from ITL - Part No. S126195 (775 SPARE VOL & SEN KNOBS)

The spares kit should consist of:

QTY	DESCRIPTION
1	KNOB 20MM (SENSITIVITY) (INCLUDES COVER, BRASS PART, WASHER AND NUT)
1	KNOB 40MM (VOLUME) (INCLUDES COVER, BRASS PART, WASHER AND NUT)

##### 3.1.2 Servicing process

1. Using a flat head screwdriver, remove the cover over the front of the knob(s) that need to be replaced.

2. Using a 10mm nut runner, undo the nut which holds the plastic knob in place.

3. Slide the knob(s) off.

4. If no further servicing of the unit is required, replace with the volume and sense knobs spares kit from ITL - Part No. S126195 (775 SPARE VOL & SEN KNOBS).

5. Ensure that the three positions on the sensitivity knob match the turning positions of the rotary switch. The grey "pointer" on the sensitivity knob must stop at positions '1', '2', and '3'.

## 3.2 Servicing a unit with damaged rear casework

### 3.2.1 Spares kit required

Order the rear casework spares kit from ITL - Part No. S126196 (775 SPARE REAR CASEWORK)

The spares kit should consist of:

QTY	DESCRIPTION
1	C_CASE_REAR_INJECTION_MOULDING
5	M3 SCREW PAN RECD SS
1	POREX MEMBRANE

### 3.2.2 Servicing process

1. Using a pozidrive screwdriver, unscrew and remove the 5 screws (M3 SCREW PAN RECD SS) which attach the rear casework to the front casework.

2. Remove the rear casework from the unit, by pulling it away from the rest of the unit and over the feet.

3. If no further servicing of the unit is required, replace with the rear casework spares kit from ITL - Part No. S126196 (775 SPARE REAR CASEWORK).

### 3.3 Servicing a unit with damaged feet

#### 3.3.1 Spares kit required

Order the feet spares kit from ITL - Part No. S126198 (775 SPARE FEET)


The spares kit should consist of:

QTY	DESCRIPTION
4	FOOT (20 DIA. X 15 MM)

#### 3.3.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.
2. Lay the unit on its side so that it is supported and the feet are accessible.

3. Using adjustable mole grips, grip the feet which need replacing and gently pull them away from the base plate.

4. Replace with the feet spares kit from ITL - Part No. S126198 (775 SPARE FEET).
5. If no further servicing of the unit is required, reattach the rear casework.

### 3.4 Servicing a unit with damaged front casework

#### 3.4.1 Spares kit required

Order the front casework spares kit from ITL - Part No. S126197 (775 SPARE FRONT CASEWORK)

The spares kit should consist of:

QTY	DESCRIPTION
1	C_CASE_FRONT_INJECTION_MOULDING
1	C_CASE_WINDOW_INJECTION_MOULDING
6	M3 X 6 SCREW PAN HD. REC'D
1	PUSH BUTTON TOP (RS 157-3650)
1	M3 X 20 SCREW (IN PUSH BUTTON)
1	SPINDLE BEARING (VOLUME BUTTON)
1	SPINDLE BEARING LARGE DIAMETER (SENSITIVITY BUTTON)

### 3.4.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.	
2. Remove the volume and sensitivity knobs, as shown in 3.1.2.	
3. Using a pozidrive screwdriver, unscrew and remove the 6 screws (M3 SCREW PAN HD. REC'D SS) which attach the front casework to the chassis.	
4. Remove the front casework from the unit by pulling it away from the rest of the unit.	
5. If no further servicing of the unit is required, replace with the front casework spares kit from ITL - Part No. S126197 (775 SPARE FRONT CASEWORK)	
6. Reattach the rear casework.	

## 3.5 Servicing a unit with damaged connector bracket, including mains inlet

### 3.5.1 Spares kit required

Order the connector bracket spares kit from ITL - Part No. S126194 (775 SPARE CONNECT BRACKET)

The spares kit should consist of:


QTY	DESCRIPTION
6	M3 X 6 SCREW PAN HD. REC'D
1	CONNECTOR BRACKET
1	POWER INPUT MODULE
1	UNIT LABEL
1	USB CONNECTOR COVER
2	SCREW, M3X6, SKT BTN, ST STL
1	USB CONNECTOR PCB
1	S123419 4.2.1 775 Cable PSU PCB Mains power
1	S123426 4.2.8 775 Cable LCD PCB download USB
1	S123431 4.2.13 775 Cable Main PCB debug USB
1	M4 SCREW PAN REC (EARTH)
2	WASHER SPROOF SS (EARTH)
1	FUSE

### 3.5.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.
2. Using a pozidrive screwdriver, unscrew and remove the 6 screws (M3 SCREW PAN HD. REC'D) which attach the connector bracket to the chassis back panel.


3. Disconnect the wiring connectors at PL1 and PL2 on the USB connector PCB.


4. Disconnect wiring connector at PL1 on the PSU PCB.

5. Using a small pozidrive screwdriver, unscrew and remove the earth screw (M4 SCREW PAN REC) which attaches the earth wire to the base plate.

6. Remove the connector bracket assembly from the unit.

<p>7. Assign the unit 'ref' and 'serial number' to the serial number label: place a suitable sized piece of paper with serial number printed on it under the window of the label. Ensure it is clearly visible.</p>	
<p>8. Stick the serial number label to the connector bracket in the correct place.</p>	
<p>9. If no further servicing of the unit is required, replace with the connector bracket spares kit from ITL - Part No. S126194 (775 SPARE CONNECT BRACKET).</p>	
<p>10. Ensure that the earth wire is correctly fixed to the base plate using the M4 SCREW PAN REC and two WASHER SPROOF SS. Ensure that the connector to PL1 on the PSU PCB is correctly connected.</p>	
<p>11. Reassemble the rest of the unit.</p>	

### 3.6 Servicing a unit with damaged IV pole boss

#### 3.6.1 Spares kit required

Order the IV pole boss spares kit from ITL - Part No. S126200 (775 SPARE IV POLE BOSS)

The spares kit should consist of:

QTY	DESCRIPTION
1	IV POLE BOSS
4	M3 X 10 CSK HA. RECD.SCREW SLDPRT

### 3.6.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.	
2. Remove the Connector bracket, as shown in 3.5.2.	
3. Lay the unit on its side so that it is supported and the M3 X 10 CSK HA. RECD.SCREW SLDPRT are accessible.	A photograph of the unit on its side. Red arrows point to the four screws that hold the IV pole boss to the base plate.
4. Using a pozidrive screwdriver, unscrew and remove the 4 screws which attaches the IV pole boss to the base plate.	
5. If no further servicing of the unit is required, replace with the IV pole boss spares kit from ITL - Part No. S126200 (775 SPARE IV POLE BOSS).	
6. Reassemble the rest of the unit.	

### 3.7 Servicing a unit with damaged PSU PCB

#### 3.7.1 Spares kit required

Order the PSU PCB spares kit from ITL - Part No. S126203 (775 SPARE PSU PCB)

The spares kit should consist of:


QTY	DESCRIPTION
1	TDK LAMDA PSU
1	PSU PCB
4	M3 X 10 HEX SPACER
8	M3 X 6 PAN HD. REC'D SCREW
1	S123420 775 Cable LCD PCB digital
1	S123421 775 Cable LCD PCB backlit
1	S123422 775 Cable LCD PCB audio
1	S123423 775 Cable Main PCB digital
1	S123424 775 Cable Main PCB analog
1	S123425 Cable Coil PCB driver

#### 3.7.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.	
2. Remove the Connector bracket, as shown in 3.5.2.	

3. Using a pozidrive screwdriver, unscrew and remove the 6 screws (M3 SCREW PAN HD. REC'D) which attach the back panel to the chassis. Remove the back panel.

4. Disconnect the wiring connectors at PL5, PL4, PL6, PL8, PL7 and PL2 on the PSU PCB.


5. Using a pozidrive screwdriver, unscrew and remove the 4 screws (M3 SCREW PAN HD. REC'D) which attach the spacers on the back of the PCB to the chassis.

6. Remove the PSU PCB from the unit.

7. If no further servicing of the unit is required, replace with the PSU PCB spares kit from ITL - Part No. S126203 (775 SPARE PSU PCB).

The wiring schematic and diagrams included in the Appendices, sections 5.1 and 5.2 respectively, may help to reconnect wires.

8. Reassemble the rest of the unit.

### 3.8 Servicing a unit with damaged Main Arm PCB

#### 3.8.1 Spares kit required

Order the Main Arm PCB spares kit from ITL - Part No. S126204 (775 SPARE MAIN ARM PCB)

The spares kit should consist of:

QTY	DESCRIPTION
1	MAIN ARM PCB
5	M3 X 10 HEX SPACER
10	M3 X 6 PAN HD. REC'D SCREW
1	S123423 775 Cable Main PCB digital
1	S123424 775 Cable Main PCB analog


1	S123428 775 Cable LCD PCB comms RXD
1	S123429 775 Cable LCD PCB comms TXD
1	S123431 775 Cable Main PCB debug USB
1	S123432 775 Cable Coil PCB sense
1	S123433 775 Cable Coil PCB probe ok
1	S123434 775 Cable Coil PCB driver sig
1	S123435 775 Cable Coil PCB sense sig

### 3.8.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.
2. Remove the Connector bracket, as shown in 3.5.2.
3. Remove the back panel, as shown in 3.7.2.
4. Disconnect the wiring connectors at PL3, PL6, PL5, PL2, PL7, PL9 and PL10 on the Main Arm PCB.
5. Using a pozidrive screwdriver, unscrew and remove the 5 screws (M3 SCREW PAN HD. REC'D) which attach the spacers on the back of the PCB to the chassis.

6. Remove the Main Arm PCB from the unit.

7. If no further servicing of the unit is required, replace with the Main Arm PCB spares kit from ITL - Part No. S126204 (775 SPARE MAIN ARM PCB).

The wiring schematic and diagrams included in the Appendices, sections 5.1 and 5.2 respectively, may help to reconnect wires.

8. Reassemble the rest of the unit.

### 3.9 Servicing a unit with damaged Pressure Switch

#### 3.9.1 Spares kit required

Order the pressure switch spares kit from ITL - Part No. S126201 (775 SPARE PRESSURE SWITCH)

The spares kit should consist of:

QTY	DESCRIPTION
1	PRESSURE SWITCH
1	S123430 775 CABLE LCD PCB FOOTSWITCH
2	SCREW NO6 X 1-2 ST PH REC ZCP
1	TYGON TUBING

1	COUPLING BODY PMCD (INCLUDES NUT)
1	SPACER A

### 3.9.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.
2. Remove the Connector bracket, as shown in 3.5.2.
3. Remove the back panel, as shown in 3.7.2.
4. Using a pozidrive screwdriver, unscrew and remove the 2 screws (SCREW NO6 X 1-2 ST PH REC ZCP) which attach the pressure switch to the chassis.

5. Disconnect the wiring connectors at PL2 on the LCD Arm PCB and the tygon tubing from the fitting on the back of the coupling body connector.

6. Using a 16mm socket with extension, or nut runner, unscrew the nut from the back the coupling body connector and remove it from the chassis.

7. If no further servicing of the unit is required, replace with the pressure switch spares kit from ITL - Part No. S126201 (775 SPARE PRESSURE SWITCH).

8. Reassemble the rest of the unit.

### 3.10 Servicing a unit with damaged Coil Interface PCB

#### 3.10.1 Spares kit required

Order the Coil Interface PCB spares kit from ITL - Part No. S126202 (775 SPARE COIL INT PCB)

The spares kit should consist of:

QTY	DESCRIPTION
1	COIL INTERFACE PCB
2	PA MODULE
2	M2 X 10 CSK HD. SLOTTED SCREW
2	M2 PLAIN WASHER
2	M2 FULL NUT
2	INSULATION SLEEVE
2	TO 220 MICA WASHER

1	PCB HEATSINK
1	HEATPAD
2	M3 X 6 PAN HD. REC'D SCREW
2	M3 FULL NUT
2	M3 NYLON SCREW
2	ODU SIZE 2 CONNECTOR (WITH FRONT FIXING)
1	S123425 775 Cable Coil PCB driver
1	S123432 775 Cable Coil PCB sense
1	S123433 775 Cable Coil PCB probe ok
1	S123434 775 Cable Coil PCB driver sig
1	S123435 775 Cable Coil PCB sense sig



### 3.10.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.
2. Remove the Connector bracket, as shown in 3.5.2.
3. Remove the back panel, as shown in 3.7.2.
4. Remove the IV Pole Boss, as shown in 3.6.2
5. Using a 19 mm spanner, unscrew and remove the 2 fixings on the front of the ODU size 2 connectors.

6. Disconnect the wiring connectors at PL3, PL4, PL5, PL6, and PL7, on the Coil Interface PCB.

<p>7. Using a small flathead screwdriver, unscrew and remove the 2 screws (M3 NYLON SCREW) which attach the PCB heatsink to the base plate.</p>	
<p>8. Remove the Coil Interface PCB from the unit. The heatpad may also need to be removed.</p>	
<p>9. If no further servicing of the unit is required, replace with the Coil Interface PCB spares kit from ITL - Part No. S126202 (775 SPARE COIL INT PCB).</p> <p>The wiring schematic and diagrams included in the Appendices, sections 5.1 and 5.2 respectively, may help to reconnect wires.</p>	
<p>10. Reassemble the rest of the unit.</p>	

### 3.11 Servicing a unit with damaged speaker

#### 3.11.1 Spares kit required

Order the speaker spares kit from ITL - Part No. S126199 (775 SPARE SPEAKER)


The spares kit should consist of:

QTY	DESCRIPTION
1	SPEAKER
1	S123427 775 CABLE LCD PCB SPEAKER
4	M3 X 12 PAN HEAD, RECESSED SCREW
4	M3 FULL NUT
6	M3 X 6 SCREW PAN HD. REC'D
1	SPEAKER PANEL


### 3.11.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.
2. Remove the Connector bracket, as shown in 3.5.2.
3. Remove the back panel, as shown in 3.7.2.
4. Disconnect the wiring connectors at PL14 on the LCD Arm PCB.


5. Using a pozidrive screwdriver, unscrew and remove the 6 screws (M3 SCREW PAN HD. REC'D) which attach the speaker panel to the chassis.

6. Remove the speaker assembly from the unit.

7. If no further servicing of the unit is required, replace with the speaker spares kit from ITL - Part No. S126199 (775 SPARE SPEAKER)

8. Reassemble the rest of the unit.

### 3.12 Servicing a unit with damaged LCD ARM PCB

#### 3.12.1 Spares kit required

Order the LCD ARM PCB spares kit from ITL - Part No. S126205 (775 SPARE LCD ARM PCB)

The spares kit should consist of:


QTY	DESCRIPTION
1	LCD ARM PCB
6	M3 X 6 SCREW PAN HD. REC'D
1	775 SWITCH ENDSTOP COLLECT
1	ROTARY SWITCH
1	ENCODER ROTARY SWITCH
2	M2.5 SCREW PAN RECD SS
1	S123426 775 Cable LCD PCB download USB
1	S123428 775 Cable LCD PCB comms RXD
1	S123429 775 Cable LCD PCB comms TXD
1	S123420 775 Cable LCD PCB digital
1	S123421 775 Cable LCD PCB backlit
1	S123422 775 Cable LCD PCB audio

### 3.12.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.	
2. Remove the Connector bracket, as shown in 3.5.2.	
3. Remove the back panel, as shown in 3.7.2.	
4. Remove the volume and sensitivity knobs, as shown in 3.1.2.	
5. Remove the front casework, as shown in 3.4.2.	

6. Using a pozidrive screwdriver, unscrew and remove the 2 screws (M2.5 SCREW PAN RECD SS) which attach the 775 end stop collect to the rotary switch. Remove the end stop.


7. Remove the PSU PCB, as shown in 3.7.2.


8. Remove the Main Arm PCB, as shown in 3.8.2.

9. Disconnect the wiring connectors at PL2, PL14, PL10 and SK4 on the LCD Arm PCB.

10. Using a pozidrive screwdriver, unscrew and remove the 6 screws (M3 SCREW PAN HD. REC'D) which attach the LCD Arm PCB to the chassis and display mounting plate.

11. Remove the LCD Arm PCB, and if no further servicing of the unit is required, replace with the LCD ARM PCB spares kit from ITL - Part No. S126205 (775 SPARE LCD ARM PCB).

12. Reassemble the rest of the unit.

### 3.13 Servicing a unit with damaged LCD

#### 3.13.1 Spares kit required

Order the LCD spares kit from ITL - Part No. S126207 (775 SPARE LCD)

The spares kit should consist of:

QTY	DESCRIPTION
1	LCD DISPLAY
1	DISPLAY MOUNTING PLATE
6	M3 X 6 SCREW PAN HD. REC'D
4	NO.4 X 6.4 PAN HD. RECD PLASTITE
4	M3 X 10 HEX SPACER
2	M3 X 18 HEX SPACER

#### 3.13.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.	
2. Remove the Connector bracket, as shown in 3.5.2.	
3. Remove the back panel, as shown in 3.7.2.	
4. Remove the volume and sensitivity knobs, as shown in 3.1.2.	
5. Remove the front casework, as shown in 3.4.2.	
6. Remove the PSU PCB, as shown in 3.7.2.	
7. Remove the Main Arm PCB, as shown in 3.8.2.	
8. Remove the LCD Arm PCB, as shown in 3.12.2.	

9. Using a pozidrive screwdriver, unscrew and remove the 4 screws (M3 SCREW PAN HD. REC'D) which attach the display mounting plate to the chassis.

10. Remove the LCD, and if no further servicing of the unit is required, replace with the LCD spares kit from ITL - Part No. S126207 (775 SPARE LCD)

11. Reassemble the rest of the unit.

### 3.14 Servicing a unit with damaged Chassis

#### 3.14.1 Spares kit required

Order the Chassis spares kit from ITL - Part No. S126206 (775 SPARE Chassis)

The spares kit should consist of:

QTY	DESCRIPTION
1	CHASSIS
1	BASE PLATE
10	M3 X 6 SCREW PAN HD. REC'D
1	EARTH LABEL
1	CHASSIS FRONT LABEL
2	M3 X 6 SCREW PAN HD. REC'D
2	M3 FULL NUT
1	BACK PANEL

### 3.14.2 Servicing process

1. Remove the rear casework, as shown in 3.2.2.	
2. Remove the Connector bracket, as shown in 3.5.2.	
3. Remove the back panel, as shown in 3.7.2.	
4. Remove the volume and sensitivity knobs, as shown in 3.1.2.	
5. Remove the front casework, as shown in 3.4.2.	
6. Remove the speaker, as shown in 3.11.2.	
7. Remove the IV Pole Boss, as shown in 3.6.2.	
8. Remove the pressure switch, as shown in 3.9.2.	
9. Remove the Coil Interface PCB, as shown in 3.10.2	
10. Remove the PSU PCB, as shown in 3.7.2.	
11. Remove the Main Arm PCB, as shown in 3.8.2.	
12. Remove the LCD Arm PCB, as shown in 3.12.2.	
13. Remove the LCD, as shown in 3.13.2	
14. Remove the feet, as shown in 3.3.2.	
15. Reassemble all removed parts onto LCD	

---


spares kit from ITL - Part No. S126206 (775 SPARE Chassis)	
16. Reassemble the rest of the unit.	

#### 4. Document History

Date	By	Issue	Reason
06/03/2013	Thomas Jull	1d1	Initial draft
07/03/2013	Thomas Jull	1d2	Internal review, minor grammar changes made
18/03/2013	Thomas Jull	1d3	To reflect feedback from Endomagnetics
25/03/2013	Thomas Jull	1	Initial issue

## 5. Appendices

### 5.1 Appendix 1 – Wiring schematic for the SentiMag System

## 5.2 Appendix 2 – Wiring diagrams

### 5.2.1 775 Cable PSU PCB Mains power

ITL Part No. S123419

FROM		USING		TO	
INLET,IEC,PNL MTG,SWITCHED,2A		-		PSU PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
Mains Inlet	L	LIVE	BRN	PL1	1
	N	NEUTRAL	BLU		5
	E	EARTH	GRN/YEL	CHASSIS	RING

### 5.2.2 775 Cable LCD PCB digital

ITL Part No. S123420

FROM		USING		TO	
PSU PCB ASSY		TWISTED PAIR [1,2]		LCD ARM PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL6	1	+3V3	RED	PL3	1
	2	0V	BLK		2
	3	N/C	N/C		3

### 5.2.3 775 Cable LCD PCB backlit

ITL Part No. S123421

FROM		USING		TO	
PSU PCB ASSY		TWISTED PAIR [1,2]		LCD ARM PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL4	1	N/C	N/C	PL12	1
	2	0V	BLK		2
	3	-5V	VIO		3

### 5.2.4 775 Cable LCD PCB audio

ITL Part No. S123422

FROM		USING		TO	
PSU PCB ASSY		TWISTED PAIR [1,2]		LCD ARM PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL5	1	+5V	YEL	PL13	1
	2	0V	BLK		2
	3	N/C	N/C		3

**5.2.5 775 Cable Main PCB digital**

ITL Part No. S123423

FROM		USING		TO	
PSU PCB ASSY		TWISTED [1,2,3]		MAIN PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL8	1	+3V3	RED	PL3	1
	2	0V	BLK		2
	3	+5V	YEL		3

**5.2.6 775 Cable Main PCB analog**

ITL Part No. S123424

FROM		USING		TO	
PSU PCB ASSY		TWISTED PAIRS [1,2] [3,4] [5,6]		MAIN PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL7	1	+3V3	RED	PL8	1
	2	0V	BLK		2
	3	+5V	YEL		3
	4	0V	BLK		4
	5	-5V	VIO		5
	6	0V	BLK		6

**5.2.7 775 Cable Coil PCB driver**

ITL Part No. S123425

FROM		USING		TO	
PSU PCB ASSY		TWISTED [1,2,3,4]		COIL INTERFACE PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL2	1	+5V	YEL	PL7	1
	2	0V	BLK		2
	3	-5V	BIO		3
	4	D-SHLD	BLK		4

**5.2.8 775 Cable LCD PCB download USB**

ITL Part No. S123426

FROM		USING		TO	
LCD ARM PCB ASSY		TWISTED PAIRS [1,4] [2,3]		775 USB CONNECTOR PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL1	1	+5V	YEL	PL1	1
	2	D-	BLU		2
	3	D+	WHT		3
	4	0V	BLK		4
	5	N/C	N/C		5

**5.2.9 775 Cable LCD PCB speaker**

ITL Part No. S123427

FROM		USING		TO	
LCD ARM PCB ASSY		TWISTED PAIR [1,2]		SPEAKER,50mm,8ohm	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL14	1	OP+	WHT	SOLDER	1
	2	OP-	BLU		2

**5.2.10 775 Cable LCD PCB comms RXD**

ITL Part No. S123428

FROM		USING		TO	
LCD ARM PCB ASSY		TWISTED PAIRS [1,2] [3,4] [5,6] [7,8]		MAIN PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL6	1	RX1-	BLU	PL6	3
	2	RX1+	WHT		4
	3	RX2-	BLU		1
	4	RX2+	WHT		2
	5	SP_RX1-	BLU		5
	6	SP_RX1+	WHT		6
	7	SP_RX2-	BLU		7
	8	SP_RX2+	WHT		8

**5.2.11 775 Cable LCD PCB comms TXD**

ITL Part No. S123429

FROM		USING		TO	
LCD ARM PCB ASSY		TWISTED PAIRS [1,2] [3,4] [5,6] [7,8]		MAIN PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL7	1	TX1-	BLU	PL5	3
	2	TX1+	WHT		4
	3	TX2-	BLU		1
	4	TX2+	WHT		2
	5	SP_TX1-	BLU		5
	6	SP_TX1+	WHT		6
	7	SP_TX2-	BLU		7
	8	SP_TX2+	WHT		8

**5.2.12 775 Cable LCD PCB footswitch**

ITL Part No. S123430

FROM		USING		TO	
LCD ARM PCB ASSY		TWISTED PAIR [1,2]		FOOTSWITCH ASSEMBLY PARTS LIST	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL2	1	N.O	WHT	HERGA -	4
	2	COM	BLK	SWITCH	3

**5.2.13 775 Cable Main PCB debug USB**

ITL Part No. S123431

FROM		USING		TO	
MAIN PCB ASSY		TWISTED PAIR [1,4] [2,3]		775 USB CONNECTOR PCB	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL1	1	+5V	YEL	PL2	1
	2	D-	BLU		2
	3	D+	WHT		3
	4	0V	BLK		4
	5	N/C	N/C		5

**5.2.14 775 Cable Coil PCB sense**

ITL Part No. S123432

FROM		USING		TO	
COIL INTERFACE PCB ASSY		TWISTED PAIR [1,2,3,4]		MAIN PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL4	1	+5V	YEL	PL7	1
	2	0V	BLK		2
	3	-5V	VIO		3
	4	S-SHLD	BLK		4

**5.2.15 775 Cable Coil PCB probe ok**

ITL Part No. S123433

FROM		USING		TO	
COIL INTERFACE PCB ASSY		TWISTED [1,2,3]		MAIN PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL5	1	/DRVEN	WHT	PL2	1
	2	OVS	BLK		2
	3	/PRBOK	WHT		3

**5.2.16 775 Cable Coil PCB driver sig**

ITL Part No. S123434

FROM		USING		TO	
MAIN PCB ASSY		CABLE,1PAIR,7/0.2,SCREEN,BLACK		COIL INTERFACE PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL9	1	DRV	WHT	PL3	1
	2	OVA	BLU		2
	3	SHLD	TIN	N/C	N/C

**5.2.17 775 Cable Coil PCB sense sig**

ITL Part No. S123435

FROM		USING		TO	
MAIN PCB ASSY		CABLE,1PAIR,7/0.2,SCREEN,BLACK		COIL INTERFACE PCB ASSY	
CONNECTOR	PIN NO.	SIG. NAME	WIRE COL.	CONNECTOR	PIN NO.
PL10	1	SIG. NAME	WHT	PL6	1
	2	GND	BLU		2
	3	SHLD	TIN	N/C	N/C

Wires which may have come out of connectors can be re-inserted into their plug pin holes and made secure by raising the barb on the back of the pin.

### 5.3 Appendix 3 – Decontamination

Probes and units used in the field should be correctly decontaminated before handled for servicing. All surfaces on the unit and probe should be cleaned with a decontaminate such as Clinell Alcoholic 2% Chlorhexidine disposable disinfection wipes.

A decontamination form should be filled out so that the service engineer undertaking the work knows that the unit and/or probe is/are safe to work on. The decontamination form may look like, and include, the same information shown in the below example:

ENDOMAGNETICS

**Decontamination Declaration**

Hospital  
or Clinic Name:

Address:

Product Code:

Serial Number:

Reason For Return:

Please tick the appropriate option(s) below:

☐ I certify that I have decontaminated the product as per above.

Decontaminant Used:

☐ I certify that the product has not been exposed to any chemical or biological materials.

Title:Name:

Signature:Date:

Telephone:Email:

**NOTE:** Please include a copy of this form with the product being returned.

Endomagnetics Ltd. incorporated in England and Wales – No: 06227698, VAT: 947 7709 68  
Correspondence: Endomagnetics Ltd., 23 Cambridge Science Park, Milton Road, Cambridge, CB4 0EY, United Kingdom  
Registered Office: Endomagnetics Ltd., Tower 42, 33rd Floor, 25 Old Broad Street, London, EC2N 1HQ, United Kingdom