


## Gliding nails; set no. 19.290.10



19.290.10  
Gliding nails

Item no.	Designation	Quantity	Figure
Implants - Nails			
<b>Gliding nail; titanium 4 poz. Titaninés elastinés vinyś</b>			
7.281.15	Gliding nail; Ø 1.5 mm x 450 mm	2	
7.281.20	Gliding nail; Ø 2.0 mm x 450 mm	2	
7.281.25	Gliding nail; Ø 2.5 mm x 450 mm	2	
7.281.30	Gliding nail; Ø 3.0 mm x 450 mm	2	
7.281.32	Gliding nail; Ø 3.2 mm x 450 mm	2	
7.281.40	Gliding nail; Ø 4.0 mm x 450 mm	2	
7.281.45	Gliding nail; Ø 4.5 mm x 450 mm	2	
7.281.50	Gliding nail; Ø 5.0 mm x 450 mm	2	
Instruments			
10.421.01	Extraction forceps for wires; flat	1	
10.509.00	Awl; Ø 8.0 mm	1	
6.903.02	Insertion instrument with quick-action chuck; cannulated, inner diameter 5.2 mm	1	
6.906.03	Power cutter for wires to maximum Ø 3.0 mm (steel) and to maximum Ø 4.0 mm (titanium)	1	
Container			
19.291.00	Perforated autoclave container with insert for instruments and gliding nails	1	




# OP-Set gliding nails, length 450 mm

4. pozicija,  
pagaminta  
iš titano

titanium Set-No. 19.290.10

Vinies proksimalinis galas be angų, distalinis  
galas plokščias ir lenktas




Item description	Item no.	quantity
<b>Implants - nails</b>		
Gliding nails; length 450 mm; titanium 4 pozicija		
diameter 1.5 mm	7.281.15	2
diameter 2.0 mm	7.281.20	2
diameter 2.5 mm	7.281.25	2
diameter 3.2 mm	7.281.32	2
diameter 3.5 mm	7.281.35	2
diameter 4.0 mm	7.281.40	2
diameter 4.5 mm	7.281.45	2
diameter 5.0 mm	7.281.50	2

## Depiction of the instruments in set



Item no. 6.903.02



Item no. 10.421.01



Item no. 10.509.00



Item no. 6.906.03

<b>Instruments</b>		
impactor with quick-action chuck; cannulated; inside diameter 5.2 mm	6.903.02	1
power cutter for wires to maximum Ø 3.0 mm (stainless steel) and Ø 4.0 mm (titanium)	6.906.03	1
extraction forceps for wires; flat	10.421.01	1
awl; Ø 8.0 mm	10.509.00	1

<b>Container</b>		
perforated autoclavable container with inset for instruments and gliding nails	19.291.00	1



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany  
Fon + 49 36738 498-0  
Fax + 49 36738 498-559  
E-Mail: export@koenigsee-implantate.de

03/09/2014/KaMo

For further information please visit:  
[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)

## Advantages of treatment with gliding nails

- The gliding nails enable rapid healing while maintaining the full scope of movement and function.
- Specially suitable for children because the application prevents growth impairments.
- Optimal stabilisation of the fracture via three-point support of the nail in the bone.
- Early functional post-operative measures are possible.
- Simple, minimally invasive implantation.
- Good cosmetic result.
- Intervention away from the fracture.

4 poz. Naudojimas: Ilgųjų kaulų diafizės, metafizės ir epimetafizės lūžių gydymui

### Indications

Long bone fractures for which three-point support is possible.

In children:

- Diaphyseal fractures of longer long bones
- Select metapyseal and epimetaphyseal fractures of long bones
- Complex fractures of the clavicle

In adults:

- Diaphyseal fractures of long bones of the upper extremities
- Clavicle shaft fractures

## Indications and contraindications

### Contraindications

- Disorders of vascularisation, which can lead to bone necrosis or infections of the surrounding soft tissues and an impaired healing process after the operation.
- Extensive soft tissue injury in the area for surgery.
- Inflammation in the area for surgery.
- Poor patient compliance.
- Allergies to the implant material.
- Implants made of implant steel or chrome-cobalt alloy or titanium already present in the vicinity of the area for treatment that are not compatible with the implant material to be inserted.
- The diameter of the medullary cavity is too low to permit insertion of a gliding nail.
- Unstable fractures for which no sufficient cortical support is possible.
- Joint fractures that should be treated open.
- The implantation is not recommended for overweight children with BMI >30.

# Table of contents

4 poz. Naudojimo instrukcijos:  
 proksimalinio žastikaulio,  
 žastikaulio diafizė, distalinio  
 žastikaulio, dilbio kaulų,  
 proksimalinio, distalinio  
 šlaunikaulio, šlaunikaulio  
 diafizės, blauzdikaulio kaulų,  
 raktikaulio sintezei


4	In detail
6	Introduction
9	Surgical instructions
9	General preliminary considerations
10	Proximal humerus and humerus shaft fractures
4	Distal humerus fractures
6	Forearm shaft fractures
8	Proximal femoral and femoral shaft fractures
20	Distal femoral fractures
22	Tibia and lower leg fractures
24	Clavicle fractures
26	Post-operative treatment
27	Implant removal
28	Ordering information
35	References and index of images





## Implantation of the nails


- Radiological check of the correct position of the humeral head relative to the shaft.
- The nail is finally anchored in the proximal epiphysis. The runners must be oriented divergently.



Positioned nails Fig. 8

- It must be ensured that the full range of motion of the humerus is restored.
- The nails are shortened with the power cutter (item no. 6.906.03 or 6.906.06 for the respective nail diameter) so that the ends lie deep within the subcutaneous tissue.

4 poz. Elastinės vinys, esant poreikiui, patrumpinamos, naudojant vielų kirpimo žnyples (priklausomai nuo vinies skersmens).



Implanted nails Fig. 9

## Conclusion and wound closure

- X-ray documentation
- Flushing of the surgical area
- Insertion of a Redon drain, if necessary
- Wound closure in layers
- Application of a sterile bandage

# PROXIMAL HUMERUS AND HUMERUS SHAFT FRACTURES

## Surgical instructions

### Preliminary considerations

- Gliding nails are used for proximal humerus fractures and humerus shaft fractures in both children and adults. The depicted operative technique is the same.
- Fractures on the proximal humerus or the humerus diaphysis are generally treated with double nail osteosynthesis.
- The nail insertion is retrograde anterolateral.

### Positioning

- The operation is carried out in the supine position.
- The fractured humerus is laid on an arm table; the shoulder and elbow must be freely accessible.

### Access and reduction




Fig. 3 Access to the distal humerus

- A closed reduction is performed first.
- Access takes place at the anterolateral (radial) aspect of the distal humerus, approximately 1 cm above the lateral epicondyle, proximally over a length of 3 to 4 cm.

**Note** Medial access should be avoided in order to avoid injury of the ulnar nerve.

# DISTAL HUMERUS FRACTURES

## Surgical instructions

### Preliminary considerations

- Gliding nails are used for distal humerus fractures in both children and adults. The depicted operative technique is the same.
- Fractures on the distal humerus or supracondylar humerus fractures are generally treated with a double nail osteosynthesis.
- The nail insertion is antegrade lateral.

### Positioning

- The operation is carried out in the supine position.
- The fractured humerus is laid on an arm table; the shoulder and elbow must be freely accessible.

### Access and reduction

- A closed reduction is performed first.
- The access takes place laterally at the height of the caudal end of the deltoid muscle distally over a length of 3 to 4 cm.

**Note** The incision may not be made too far distally to avoid injuring the radial nerve.



# FOREARM SHAFT FRACTURES

## Surgical instructions

### Preliminary considerations

- Gliding nails are used for forearm fractures in both children and adults. The depicted operative technique is the same.
- Forearm fractures are generally treated with one nail per bone.
- The nails are implanted in opposing directions (radius ascending, ulna descending) or both ascending.
- Isolated fractures of ulna or radius are treated with a single nail in the corresponding bone.

**Attention** The nail diameter to be selected is generally approximately 60% of the medullary cavity, measured at the narrowest point.

### Positioning

- The operation is carried out in the supine position.
- The fractured forearm is laid on an arm table.

### Access and reduction

- A closed reduction is performed first.
- Access at the distal radius:
  - Incision approximately 1 - 2 cm proximal of the radius epiphyseal plate.
  - At the height of the dorsal tubercle of radius, a 2 - 3 cm long transversal or longitudinal skin incision is made.
- Access at the proximal ulna:
  - At the height of the proximal lateral aspect of the olecranon, a 1.5 - 2 cm long transversal incision is made 3 cm distal of the apophysis.

**Note** Care must be taken to spare the radial nerve during access at the distal radius.



# PROXIMAL FEMORAL AND FEMORAL SHAFT FRACTURES

## Surgical instructions

### Preliminary considerations

- Use of the gliding nails for proximal femoral fractures or fractures of the femoral diaphysis in children generally consists of a double nail osteosynthesis.
- The insertion for both nails is retrograde.

### Positioning

- The operation is carried out in the supine position.
- Use of a fracture table, if necessary.

### Access and reduction

- A closed reduction is performed first.
- Access takes place symmetrically on both sides (medial and lateral) at the distal femur at the height of the proximal patella pole (above the epiphyseal plate) proximally over a length of 2 to 3 cm.

### Insertion of the nails

- Opening of the cortical bone according to the previously described procedure with the awl (item no. 10.509.00) with a perforation (2 - 3 cm proximal of the superior patella pole) at the cranial end of the incision.
- Prebend the nails according to the previously described procedure to achieve optimal three-point support with two C-shaped nails.
- The nails clamped into the insertion instrument (item no. 6.903.02) are inserted one after the other and driven proximally up to the fracture zone.
- The nails clamped into the insertion instrument (item no. 6.903.02) are inserted one after the other and driven proximally up to the fracture zone.

# DISTAL FEMORAL FRACTURES

## Surgical instructions

### Preliminary considerations

- Use of the gliding nails for distal femoral fractures in children generally consists of a double nail osteosynthesis.
- The insertion for both nails is antegrade.

### Positioning

- The operation is carried out in the supine position.
- Use of a fracture table, if necessary.

### Access and reduction

- A closed reduction is performed first.
- The access takes place laterally below the trochanter major over a distal length of 3 - 5 cm.

# TIBIA AND LOWER LEG FRACTURES

## Surgical instructions

### Preliminary considerations

- Use of the gliding nails for isolated tibia shaft or lower leg fractures in children generally consists of a double nail osteosynthesis.
- The insertion for both nails is antegrade.

### Positioning

- The operation is carried out in the supine position.
- Use of a fracture table, if necessary.
- The knee should be bent with a support.

### Access and reduction


- A closed reduction is performed first.
- The access takes place symmetrically on both sides (medial and lateral) at the proximal tibia at the height of the tuberosity distally over a length of 2 cm.

### Insertion of the nails

- Opening of the cortical bone according to the previously described procedure with the awl (item no. 10.509.00) with two symmetrical perforations at the caudal end of the incision.
- Prebend the nails according to the previously described procedure to achieve optimal three-point support with two C-shaped nails.



## E. Guided Growth System (Plus)



5. pozicija

**1**  
T80212 Guided Growth Plate System Plus  
TI eight-Plate L 12mm

**2**  
T80216 Guided Growth Plate System Plus  
TI eight-Plate L 16mm

**3**  
T80220 Guided Growth Plate System Plus  
TI eight-Plate L 20mm

**4**  
T80416 Guided Growth Plate System Plus  
TI quad-Plate L 16mm

**5**  
T80422 Guided Growth Plate System Plus  
TI quad-Plate L 22mm

Note: All the Plates are also available packaged sterile. They can be ordered using the above code numbers preceded by 99- (e.g. 99-T80212).

See Operative Technique EP-1701-OPT  
"Guided Growth System Plus"



Manufactured by:  
**ORTHOFIX Srl**  
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy  
Telephone +39 045 6719000, Fax +39 045 6719380



# 1. Screws

1



Guided Growth Plate System Plus TI Solid Screw Ø3.5

Length (mm)	Code
12	T80312
14	T80314
16	T80316

2



Guided Growth Plate System Plus TI Solid Screw Ø4.5

Length (mm)	Code
24	T80024
32	T80032
36	T80036

2



Guided Growth Plate System Plus TI Cannulated Screw Ø4.5

Length (mm)	Code
16	T80116
24	T80124
32	T80132
36	T80136

5 pozicija:  
4,5 mm  
kaniuluoti  
sraigai

1

Ø3.5 Solid Screw eight-Plate Plus Color Coded - Yellow

2

Ø4.5 Solid Screw eight-Plate Plus Color Coded - Blue

3

Ø4.5 Cannulated Screw eight-Plate Plus Color Coded - Green

Note: All the Screws are also available packaged sterile. They can be ordered using the above code numbers preceded by 99- (e.g. 99-T80312).



GUIDED GROWTH SYSTEM (PLUS)





## Features and benefits

- Dedicated for redirecting the growth of long bones in growing children
- Color coded plates and screws
- Sterile and non-sterile implants
- eight-Plate (for two screws placement) and quad-Plate (for four screws placement)
- Unique screw hole design for a screw angulation up to 60 degrees

5 pozicija. Aštuoniukės tipo plokštelė, skirta ilgųjų kaulų augimo linijos stabdymui

Aštuoniukės tipo plokštelė, fiksuojama dvejais sraigtais



## Benefits to Surgeons

- Plates sizes: 12, 16 and 20mm
- Cannulated and solid screw options
- Titanium alloy implants
- Optimized and lean instrumentation
- Central "dome" designed to aid application and removal across the growth plate

5 pozicija

Plokštelė pagaminta iš titano lydinio, Plokštelės centre - viena anga, skirta nukreipiklio fiksacijai



## Benefits to Patients

- Minimally invasive
- Early weight bearing as tolerated by the patient and under surgeon discretion
- Titanium alloy implants to avoid allergic reaction to nickel
- Low profile plates

