




# Broncho-Fiberscopes

| Broncho-Fiberscopes                                 |   | Set. no.   |   | Deflection up/down  | Direction of view | Angle of view | Working length |
|---|---|------------|---|---|-------------------|---------------|----------------|
| 6.4 mm x 54 cm<br>Broncho-Fiberscope<br>11004 BC1   |   | 11004 BCK1 |    |    | 0°                | 110°          | 54 cm          |
| ■ 5.6 mm x 54 cm<br>Broncho-Fiberscope<br>11009 BC1 | ■ | 11009 BCK1 |    |    | 0°                | 110°          | 54 cm          |
| ■ 5.2 mm x 54 cm<br>Broncho-Fiberscope<br>11001 BN1 | ■ | 11001 BNK1 |   |   | 0°                | 110°          | 54 cm          |
| 3.7 mm x 54 cm<br>Broncho-Fiberscope<br>11002 BD1   |   | 11002 BDK1 |  |  | 0°                | 90°           | 54 cm          |
| 2.8 mm x 54 cm<br>Broncho-Fiberscope<br>11003 BC3   |   | 11003 BCK3 |  |  | 0°                | 90°           | 54 cm          |

2.2 Lenkimo aukštyn/žemyn kampai: 180 /100

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2.5 Darbinis ilgis 540 mm

1. 2 Lenkimo aukštyn/žemyn kampai: 180 /100

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1.5 Darbinis ilgis 540 mm

Components/Spare Parts see chapter 13


| Length | Working channel inner diameter | Distal tip outer diameter | Accessories (included in delivery) |   |                               |                               |                               |                               |                               |
|--------|--------------------------------|---------------------------|------------------------------------|---|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|
|        |                                |                           |                                    |   | Broncho-Fiberscope 11004 BCK1 | Broncho-Fiberscope 11009 BCK1 | Broncho-Fiberscope 11001 BNK1 | Broncho-Fiberscope 11002 BDK1 | Broncho-Fiberscope 11003 BCK3 |
| 82 cm  | 2.8 mm                         | 6.4 mm                    | 6011590                            | Plug, for LUER-Lock connector             | ●                             | ●                             | ●                             | ●                             | ●                             |
| 82 cm  | 2.6 mm                         | 5.6 mm                    | 11025 E                            | Pressure Compensation Cap                 | ●                             | ●                             | ●                             | ●                             | ●                             |
| 82 cm  | 2.6 mm                         | 5.6 mm                    | 13242 XL                           | Leakage Tester                            | ●                             | ●                             | ●                             | ●                             | ●                             |
| 82 cm  | 2.6 mm                         | 5.6 mm                    | 11301 CD                           | Irrigation Adaptor                        | ●                             | ●                             | ●                             | ●                             | ●                             |
| 82 cm  | 2.3 mm                         | 5.2 mm                    | 11301 CE*                          | Suction Valve, unsterile, for single use  | ●                             | ●                             | ●                             | ●                             | ●                             |
| 82 cm  | 2.3 mm                         | 5.2 mm                    | 110950-50*                         | Cleaning Brush, unsterile, for single use | ●                             | ●                             | —                             | —                             | —                             |
| 82 cm  | 1.5 mm                         | 3.7 mm                    | 110941-50*                         | Cleaning Brush, unsterile, for single use | —                             | —                             | ●                             | —                             | —                             |
| 82 cm  | 1.5 mm                         | 3.7 mm                    | 110931-50*                         | Cleaning Brush, unsterile, for single use | —                             | —                             | —                             | ●                             | ●                             |
| 82 cm  | 1.2 mm                         | 2.8 mm                    | 110380-10*                         | Biopsy Forceps, sterile, for single use   | ●                             | ●                             | ●                             | —                             | —                             |
| 82 cm  | 1.2 mm                         | 2.8 mm                    | 600007                             | LUER-Lock Tube Connector                  | ●                             | ●                             | ●                             | ●                             | ●                             |

2.4 Kanalo diametras  
2,6 mm

2.3 Įvedamo vamzdelio  
diametras 5,6 mm

1.4 Kanalo diametras  
2,3 mm

1.3 Įvedamo vamzdelio  
diametras 5,2 mm

### Special Features:

- Available as battery and rechargeable versions
- LED service life of up to 40,000 hours
- Burning time > 120 min

1.7/2.7 Komplekte pateikiamas  
miniatiūrinis LED šviesos šaltinis,  
suderinamas su fibrobronchoskopu

11301 DE

11301 DG

11301 D4

**Battery Light Source LED for Endoscopes**, with fast screw thread, boost mode for temporary increase in brightness, burning time > 120 min, brightness > 110 lm / > 150 klx, weight approx. 150 g ready for use, **suitable for wipe disinfection**, for use with KARL STORZ endoscopes

11301 DE

**Battery Light Source LED for Endoscopes**, rechargeable, with click connection, boost mode for temporary increase in brightness, lithium-ion batteries, charging time 60 min, burning time at 100% brightness 40 min, weight approx. 150 g, **suitable for wipe disinfection**

11301 DF

**Same**, with coarse screw thread

11301 DG

**Charging Unit**, for two LED battery light sources, with fixed integrated power supply and adaptor for EU, UK, USA and Australia, power supply 110 – 240 VAC, 50/60 Hz, **suitable for wipe disinfection**, for use with Battery Light Source LED 11301 DE/DF

094129

**Battery Charger Li-Ion**, for charging Rechargeable Battery Box 094124 or Battery Light Source LED 11301 DE/DF, power supply 100 – 240 VAC, 50/60 Hz, for use with Mains Cord 094127

094127

**Mains Cord**, for Battery Charger 094129, length 150 cm

### Note:

Please order charging unit or battery charger with mains cord for Battery Light Sources LED 11301 DE/DF.

**Components/Spare Parts** see chapter 13

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ. This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of fiberoscopes manufactured by KARL STORZ. The fiberoscopes are sensitive precision instruments which require particular care during reprocessing. To prevent damage during reprocessing, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ fiberoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

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
This instruction manual describes three fiberolescope versions, distinguished from one another by the imaging system. If a section refers exclusively to one of these versions, it is thus appropriately labeled.

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Telescope with special filters for PDD/AF

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DCI®-Fiberskope (Direct-Coupled-Interface)

DCI® fiberoscopes (Direct Coupled Interface)

Fibroscopios DCI® (Direct-Coupled-Interface)

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




KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope for supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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- ① Focusing ring (if available)
- ② Suction valve insert (if available)  
Press = suction activated  
Release = suction deactivated
- ③ Valve nozzle (if available)
- ④ Up/down control
- ⑤ Valve for pressure compensation and leakage  
tester connection
- ⑥ Suction outlet
- ⑦ Light inlet piece
- ⑧ Instrument channel
- ⑨ Objective lens
- ⑩ Light outlet
- ⑪ Instrument in instrument channel

- ⑫ PDD filter changer:  
ALA for special tumor marker based on  
aminolevulinic acid  
W for white light  
AF autofluorescence mode

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
**2 Controls**


**Connectors, accessories and their uses**

|         |  |
|---------|--|
| 495F    | Light cable adaptor for Wolf   |
| 495G    | Light cable adaptor for KARL STORZ, Olympus, Winter & Ibe              |
| 11025E  | Pressure compensation cap (red)  |
| 11301CA | Lip valve, disposable version, pack of 20                              |
| 13242XA | Adaptor for leakage test in cleaning and disinfection machines (green) |
| 29100   | Plug for instrument channel, pack of 10, reusable                      |
| 11301CB | Suction valve insert, reusable version                                 |
| 5917800 | Seal ring (suction valve insert)                                       |
| 6885391 | Cap (suction valve insert)   |
| 5923110 | Spring (suction valve insert)  |
| 6927791 | Housing (suction valve insert)   |
| 8567290 | Plunger (suction valve insert)   |
| 11301CE | Suction valve insert, disposable version, pack of 20                   |
| 11301D1 | LED battery-powered light source with fine screw thread                |
| 11301D3 | LED battery-powered source with quick-release thread                   |
| 6927691 | Three-way double stopcock adaptor with LUER-Lock                       |
| 6011590 | LUER cap   |
| 600007  | Barbed tube, thin, with LUER-Lock                                      |

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
20260030 DCI® adaptor, for connecting endoscopes with DIN eyepiece to DCI® camera heads

20260031 DCI® eyepiece adaptor for connecting DCI® telescopes to camera heads with mounts for DIN eyepieces

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39107ALK Cleaning adaptor for use with small LUER stopcocks

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### 3 Intended Use

Fiberscopes are used in diagnostic and surgical endoscopy. They must only be used by appropriately trained and qualified technical staff.

**Indications:** Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in various medical fields.

**Contraindications:** No contraindications are known which directly relate to the product. The fibroscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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#### 3.1 Safety precautions when using the fiberscopes

Fiberscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

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#### **4 User qualification**

KARL STORZ fiberscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique.

The information given in these instructions only serves to instruct in the correct handling, cleaning, disinfection and sterilization of the fiberscopes. It is not intended as an introduction to the endoscopic technique.

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## 5 Warnings and cautions

Please read this manual and follow its instructions carefully. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the text should be carefully reviewed to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the product.

**NOTE:** A Note indicates special information about operating the product, or clarifies important information.

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**WARNING:** Before using the fiberscope, read the following safety instructions carefully, to avoid putting your patients, personnel or yourself at risk.

**WARNING:** Inspect the fiberscope prior to each surgical procedure to ensure that it is functioning correctly and that it has been correctly cleaned, disinfected and (where applicable) sterilized.

**WARNING:** This equipment may not be used during the discharge of a cardiac defibrillator. Remove all endoscopic equipment from patient contact and from the surgical site.

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**WARNING:** Verify the sound condition of the fibroscope as well as of accessories used in combination before every application. Damaged fibscopes or damaged accessories must not be used.

**WARNING:** Combinations of medical devices are only then assured to be safe if

- they are identified as such in the respective instruction manual or
- the intended purpose and the interface specifications of the devices used in combination permit this.

Pay careful attention to the instructions and interface specifications of medical products used in combination.

**WARNING:** Lengthy direct contact between the endoscope tip and surrounding tissue may cause tissue damage due to the potentially high temperatures on the tip of the endoscope and must thus be avoided.

**WARNING:** The fibroscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

**ATTENTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

**WARNING:** Before the fibroscope is removed from the patient, the locking mechanism 15 must be released and the distal tip of the fibroscope put in the neutral, non-angled position. Otherwise the instrument may be damaged or the patient injured.

**WARNING:** If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the fibroscope in the neutral, non-angled position and remove the fibroscope slowly and carefully from the patient.

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**CAUTION:** Do not steam sterilize (autoclave) the fiberscopes.

**CAUTION:** The maximum permissible temperature for cleaning, disinfection, sterilization, transport and storage is 65 °C.

**CAUTION:** Never immerse fiberscopes in normal saline solution, since even only brief contact can lead to pitting and corrosion.

**CAUTION:** Material damage: if the outer casing of the fiberscope is cleaned with alcohol, this will damage the material.

**CAUTION:** Fiberscopes must not be cleaned in an ultrasound bath under any circumstances.

**CAUTION:** Only use KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

**CAUTION:** Only use guide probes with flexible ends. Other types of guide probe may cause considerable damage to fiberscopes.

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**NOTE:** Any damage to the fiberscope resulting from incorrect operation is not covered by the guarantee.

**NOTE:** Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

**NOTE:** If the instrument is in the working channel, its additional rigidity may limit angling.

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**NOTE:** Federal law restricts this device to sale by or on the order of a physician (USA only).

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**6 Unpacking**

Check for missing items and evidence of shipping damage. File any complaints to KARL STORZ or the supplier immediately.

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**6.1 Basic equipment**

|  | Art. no.   |
|--|------------|
| 1 Fiberscope   |            |
| 1 Suction valve insert, disposable version, pack of 20 inserts   | 11301 CE   |
| 1 Leakage tester   | 13242 XL   |
| 1 Pressure compensation cap  | 11025 E    |
| 1 Plug for instrument channel, pack of 10  | 29100      |
| 2 Adaptor for cleaning the working channel in cleaning machines (green)  | 11301 CD   |
| 1 Cleaning brush for inner diameter 1.2–1.5 mm   | 11093x-50  |
| for inner diameter 2.0–2.3 mm  | 11094x-50  |
| for inner diameter 2.6–3.0 mm  | 11095x-50  |
| for 2.8 x 50 and for 2.8 x 70  | 27651 A*)  |
| for 3.7 x 65   | 27651 AL*) |
| for 5.0 x 65   | 27651 B    |
| for 5.6 x 65   | 27651 C    |
| for 6.4 x 65   | 27651 C    |
| 1 Biopsy forceps (only bronchofiberscopes)   | 11001 KL   |
| 1 Grasping forceps (only bronchofiberscopes)   | 11002 KS   |
| 1 Three-way double stopcock adaptor with LUER-Lock for O <sub>2</sub> connection (only intubation fiberscopes) | 6927691    |
| 1 LIPP tube holder (only intubation fiberscopes)   | 11301 CF   |
| 1 Mouth piece (only gastro-, intubation and bronchofiberscopes)  | 13272      |

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**6.1.1 Accessories**

- |  |                          |
|--|--------------------------|
| 1 Suction valve insert,<br>reusable version                                    | 11301 CB                 |
| 1 Battery light source   | 11301 D                  |
| 1 Adaptor for leakage test in<br>cleaning and disinfection<br>machines (green) | 13242 XA                 |
| 4 Cleaning adaptors for cleaning<br>machines                                   | 11301 GG2/<br>FF2/HH/KK2 |

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**7 Initial operation**

**WARNING:** The fiberscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

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**7.1 Fitting the suction valve**

Place the suction valve insert ② on the valve nozzle ③ and turn until the projections on the side of the insert lie over the nozzle recesses.

Push gently until the suction valve insert clicks into the valve nozzle.

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**7.2 Fitting the tube holder (only intubation fiberscopes)**

Push the tube holder onto the fiberscope from distal and secure it by turning the two parts in opposite directions.

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**7.3 Fitting the lip valve (only bronchofiberscopes)**

Place the disposable lip valve onto the instrument channel ⑧ and tighten gently by turning it in a clockwise direction (plastic part).

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## 8 Handling the fiberscopes

The optical section of KARL STORZ fiberscopes consists of a fiberoptic bundle in a flexible and watertight casing, a light guide, as well as objective and eyepiece lenses. Like all optical systems, these parts are very delicate. For this reason, careful handling of the fiberscopes is essential for a long service life and a consistent quality. Therefore be aware of the following points:

- Never kink or tightly coil the fiberscopes and light guides, so that no optical fibers break off. It is therefore best to suspend the fiberscopes for storage or keep in large containers.
- Do not jerk the up/down control ④ or move forcefully.
- Do not insert the instruments or cleaning brushes forcefully through the channels. This may cause damage to the casing, especially if the tip of the instrument is deflected.
- Carefully protect the objective lens ⑨ at the distal end from damage. Do not hit. Never remove soiling with sharp-edged objects.
- Never store or disinfect the fiberscope together with sharp-edged objects to prevent damage to the exterior.
- Never place heavy objects on top of the fiberscope.

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**NOTE:** Every fiberscope supplied by KARL STORZ includes a photo showing the condition of the fiber bundle upon delivery. A copy of this photo is kept on file in our quality control office. The wedge-shaped marking (at 12 o'clock) serves as a guide for the user.

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## 9 Using the fiberscopes

The following instructions are recommended for the proper operation of the fiberscope and accessories. This is not a reference for endoscopic techniques.

Consult the appropriate medical literature for information on endoscopic techniques.

You can find accessories for your fiberscope (catheters, flexible scissors and forceps, guide wires, stone baskets, lithotripter probes, laser fibers, electrodes etc.) in the appropriate KARL STORZ catalog.

**WARNING:** Always adjust the light source to the minimum illumination intensity necessary to achieve optimum visualization of the operating field. The higher the intensity setting of the light source, the greater the heat energy that will be generated at the tip of the fiberscope.

**WARNING:** Never lay the end of the fiber optic cable or the distal tip of the fiberscope on or under a sterile cover cloth. The light intensity may lead to burns on the patient and/or sterile cover cloth.

- When connecting and operating light sources and suction/irrigation equipment, be certain to read the appropriate instruction manuals for the equipment. The irrigation equipment recommended by KARL STORZ is set to a maximum pressure of 54 kpa (400 mmHg/0,54 bar/7.8 psi).
- Endoscopic examinations must be carried out under aseptic conditions.
- When introducing the fiberscope into the patient, the distal tip of the fiberscope must be in a straight, nondeflected position (control lock-in can be felt and heard). The fiberscope must be handled according to the recognized medical rules and procedures for endoscopic methods.

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- If no instrument is inserted, close the instrument channel proximally with a plug (Art. no. 29100).
- If an instrument is to be inserted, position the disposable lip valve (Art. no. 11301 CA) to seal off the instrument channel proximally.
- Instruments especially designed for this may be inserted through the instrument channel of the fiberscopes. Make certain the accessories have a suitable length and diameter. Unsuitable accessories may damage the fiberscopes. For questions on accessories, please contact KARL STORZ directly or an authorized dealer.

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**NOTE:** When the lip valve is first used, the instrument is subjected to increased resistance. This resistance will decrease the more often the lip valve is used.

- When inserting an instrument into the instrument channel, the fiberscope must be in a straight, nondeflected position.

**CAUTION:** When inserting instruments, increasing resistance is felt with increasing deflection. Forced introduction of instruments can lead to damage to the fiberscope or the instrument channel.

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**CAUTION:** If undue resistance is encountered when inserting an instrument, do not attempt to force the instrument. Withdraw the instrument and apply a water soluble lubricant to the tip of the instrument and attempt to insert again. Prior to removal of the instrument from the patient, the distal tip must be in a straight position.

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- When using forceps, scissors or other accessories that can be opened or closed, ensure that the instrument is completely closed prior to insertion into the instrument channel. Do not attempt to open the jaws until the distal tip is completely free of the channel. Opening the jaws inside the instrument channel can lead to damage to the fiberscope.

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**WARNING:** When using HF electrodes or laser probes, make sure that the active parts are clear of the instrument channel and in full view of the operating surgeon before power is applied. Activation within the instrument channel may damage the fiberscope and/or injure the patient.

**WARNING:** Make certain to read the instructions for the laser unit and/or HF generator.

**CAUTION:** If accessories with sharp tips, such as EHL probes or laser probes, are inserted into the instrument channel of a deflected fiberscope, this can damage the fiberscope.

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**CAUTION:** When using a stone basket, the stone must not be retracted through the instrument channel but instead the entire fiberscope must be removed. Once the stone is removed, the stone basket may then be taken out of the instrument channel.

**CAUTION:** The fiberscopes must be cleaned and then reprocessed immediately after the intervention in order to prevent incrustations. Dried-on incrustations may jeopardize successful decontamination and disinfection.

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**CAUTION:** Use only KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

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
**NOTE:** If the instrument is in the working channel, its additional rigidity may limit maneuverability.

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DCI®

## 10 Preparation of the fiberscope

Before using the fiberscope on a patient, exterior damage and operational problems must be eliminated.

Check the outer sheath of the fiberscope. Be especially aware of visible damage to the casing, particularly in the area of the fiberscope tip.

Look through the eyepiece and make certain that you have a clear image. The grid (Moiré effect) should be easily recognizable. If the image is not clear but cloudy, clean the eyepiece and lens with a cotton applicator soaked in alcohol, maintaining hygienic conditions at all times.

DCI® fiberscope image quality can be controlled on screen. If necessary, clean image waveguides, light cables and eyepieces as described above.

If it is still cloudy, it is possible that the objective lens is scratched or that there is moisture inside the telescope. The fiberscope must not be used in this case, but instead sent in for repairs.

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### 10.1 Connect these supply lines

- If necessary, remove the adaptor 495 G or 495 F from the light inlet piece and connect the fiberoptic light cable or battery-powered light source to the light inlet piece ⑦.

**NOTE:** Not necessary for battery-powered light sources with coarse threads (11301 D3).

**WARNING:** The instruction manuals of the supply units used must be observed.

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- suction to suction outlet ⑥
- or O<sub>2</sub> supply (under anesthesia) to instrument channel ⑧ via the appropriate adaptor

**CAUTION:** The instruction manuals of the supply units used must be observed.

Switch on the respective supply units.

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Insert the DCI® fiberscope into the DCI® camera head until it locks into place.

Insert the light cable (Art. no. 495 DV) into the DCI® camera head as far as it will go.

**WARNING:** The instruction manuals of the supply units used must be observed.

Connect the light cable and the video signal cable up to the appropriate supply units.

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
## 10.2 Controlling the distal tip

The up/down movement of the distal tip on the fiberscope can be adjusted at the control ④: move smoothly and slowly.

### Some fiberscopes have:

- contrapositive deflection direction (urology);
- a locking mechanism ⑮ for securing the fiberscope tip (gynecology).

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### Models with locking mechanism ⑮:

To secure the deflection position of the distal tip, turn the locking mechanism ⑮ slowly to distal until resistance can be felt.

**CAUTION:** Forcing the locking mechanism ⑮ further in the same direction may damage the locking mechanism ⑮.

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**! WARNING:** Before the fiberscope is removed from the patient, the locking mechanism 15 must be released and the distal tip of the fiberscope put in the neutral, non-deflected position. Otherwise the instrument may be damaged or the patient injured.

To release the locking mechanism 15 from the secure position, turn it slowly to proximal.

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### 10.3 Suction

Press the suction valve 2 to start suction (when the suction pump has been connected and switched on).

**! CAUTION:** There should be no instrument in the suction channel during suction. The instrument channel inlet should be closed so that no secondary air is drawn in.

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### 10.4 Insufflation

Under anesthesia, connect the O<sub>2</sub> hose to the instrument channel:

- using a barbed tube (Art. no. 600007) on the instrument channel or
- if you are using the three-way double stopcock adapter with LUER-Lock (Art. no. 6927691), attach the hose to the tube holder on the angled part of the adapter. This ensures that the straight-line passage of the adapter remains free for instrument insertion.

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### 10.5 Focusing

The image is focused with the focusing ring 1.

**i NOTE:** When connecting photo or video systems recommended by KARL STORZ, the fiberscope focusing ring must be at the zero position (mark).


DCI® fiberscopes are focused using the adjustment wheel of the DCI® camera (see page 28).

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### 10.6 Instrument channel

Flexible instruments can be inserted through the instrument channel ⑧ as required.

If the instrument channel is not being used: replace the sealing cap.

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### 10.7 PDD version

Fiberscopes intended for use in Photodynamic Diagnosis (PDD) or autofluorescence (AF) are illuminated with a special diagnosis light (D-Light C). Adjust the PDD filter changer ⑫ for special filters designed for this light and possibly for marker substances:

|     |   |
|-----|---|
| ALA | for a special tumor marker based on aminolevulinic acid |
| W   | for white light   |
| AF  | autofluorescence mode                                   |

**NOTE:** If the filter changer is set to 'ALA' or 'AF', a yellowish image is produced when viewing through the integrated special filter.

**NOTE:** No filter changer is available on PDD fiberscope art. no. 11272 CI. The filter (ALA) is permanently installed.

For optimal illumination, use the fluid light cable 495 FP (2.5 m) or 495 FO (1.8 m) between the light source and PDD fiberscope.

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## 11 DCI® fiberscope with camera

### 11.1 Assembly

Slide the DCI® into the snap-in DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

Connecting the special light cable 495 DC or 495 DV to the DCI® camera head:  
Insert the light cable plug into the DCI® camera head until it locks into place.  
While doing this, hold the light cable by the handle piece. Never pull on the cable!

**NOTE:** The entrance and exit faces of the light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Clean the surfaces with a cotton applicator soaked in cleaning solution, then disinfect it with a cotton applicator soaked in 70 % isopropyl alcohol and dry.

DCI® fiberscopes are focused using the adjustment wheel on the DCI® camera.  
The camera control unit is controlled using the two programmable control keys (silver, blue) on the DCI® camera head.

**NOTE:** Further controls are described in the instruction manual of the camera control unit used.

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### 11.2 Disassembly

Hold the DCI® camera and the DCI® fiberscope still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released.  
Pull the light cable out of the DCI® camera head.

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### 12 DIN eyepiece on DCI® fiberscope

The DCI® eyepiece adaptor 20260031 ensures the connection of DCI® telescopes to camera heads with DIN eyepiece mounts.

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### 12.1 Assembly

Slide the DCI® fiberscope into the 'snap-in' DCI® coupling mechanism of the DCI® eyepiece adaptor until the coupling snaps in automatically.  
Connect a light cable or a battery-powered light source.

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### 12.2 Disassembly

Remove the light cable or battery-powered light source.

Hold the DCI® eyepiece adaptor and the DCI® fiberscope. Press the release button on the DCI® eyepiece adaptor. The 'snap-in' DCI® coupling mechanism is released.

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### 13 DCI® camera on DIN eyepiece

#### 13.1 Assembly

Put the focusing ring (if available) of the fiberscope in the zero position (mark).

Slide the DCI® adaptor 20260030 into the 'snap-in' DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

Connect the standard light cable at the light inlet piece ⑦ of the fiberscope.

**NOTE:** The special light cable 495 DC or 495 DV cannot be used and may have to be removed from the DCI® camera.

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Opening the instrument coupling: turn the outer ring against the spring force, insert the fiberscope DIN eyepiece and let go of the outer ring. The instrument coupling closes automatically as a result of the spring force.

**CAUTION:** No further tightening of the outer ring to position the eyepiece is required. This may, among other things, cause damage to the instrument coupler!

#### 13.2 Disassembly

Hold the DCI® camera and the fiberscope still. Open the instrument coupler for the DCI® adaptor and remove the fiberscope.

Hold the DCI® camera and the DCI® adaptor still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released. Remove the DCI® adaptor.

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## 14 Reprocessing

### General warnings

-  **WARNING:** Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from [www.karlstorz.com](http://www.karlstorz.com).
-  **WARNING:** Risk of infection from non-sterile delivery. Reprocess the medical device before each use. Check the medical device for visible soiling before use.
-  **WARNING:** Before reusing the device, replace any defective parts or send the flexible endoscope for repair.
-  **ATTENTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.
-  **ATTENTION:** Never clean flexible endoscopes with metal brushes.
-  **ATTENTION:** Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning devices, for example, is possible.
-  **ATTENTION:** Flexible endoscopes must never be steam sterilized (autoclaved).
-  **ATTENTION:** To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

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- ① **NOTE:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.
- ① **NOTE:** Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.
- ① **NOTE:** National laws and regulations must be observed.
- ① **NOTE:** We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.
- ① **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

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#### 14.1 Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester (different adaptor sizes)

| Application site | Art. number for<br>brush/adaptor          |
|------------------|---|
| Surface          | 27652, dry, low-lint cloth                |
| Working channel  | Use a suitable brush, adaptor and syringe |
| Suction channel  | Use a suitable brush, adaptor and syringe |
| Vent port        | Use a suitable adaptor                    |

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#### 14.2 Precleaning

- Remove the entire instrument set from the instrument channel.
- Remove the suction valve set from the valve nozzle: press the projections on the side together and pull out the suction valve insert.
- Dispose of the removed single-use parts.

Dismantle other parts if necessary and reprocess them separately.

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use. Wipe-clean the exterior surfaces of the endoscope with a damp cloth.

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Rinse out the working channel in order to prevent drying.

**i** **NOTE:** To avoid damaging the sheath, do not use pushing and pulling movements.

**i** **NOTE:** Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

### 14.3 Leakage test

**i** **NOTE:** The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

#### 14.3.1 Phase 1: Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry.

1. Attach the adaptor of the leakage tester securely to the product's vent port and lock in place using the bayonet mount.
2. Using the leakage tester, pump air into the product until a pressure of 160 – 180 mmHg (blue area) is reached.
- ☛ The angled rubber at the distal end can easily expand while the air is pumped into the product.
3. Using the angling mechanism, move the distal end on the product's control body in all directions as far as possible while observing the pressure display of the leakage tester.
- ☛ The air pressure in the product can drop by up to 10 mmHg when moving the distal end.

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4. After stabilizing the air pressure in the product, wait at least 2 minutes and check on the leakage tester whether the air pressure in the product drops.
5. If the air pressure drops, the product is leaking and may no longer be used, immersed in liquid or reprocessed by machine. Complete the repairs form. Clean the product by wiping to remove organic soiling/residue and mark it as not disinfected. Contact KARL STORZ and send the product to KARL STORZ for repair. Line the carrying case for the product with the transport protective film 13990 SFN.
- ☛ If the air pressure does not drop, the product is leakproof. The leakage test: Phase 2 can be performed with the product.

**i** **NOTE:** To increase safety, the leakage test: Phase 1 can also be carried out directly after using the product.

#### 14.3.2 Phase 2: Performing the leakage test – in liquid

1. Lay the product in liquid and ensure that the leakage tester is not moist.
2. Fully cover the product with liquid and wet any lumens present without bubbles.
3. Using the leakage tester, pump air into the product until a pressure of 160 – 180 mmHg (blue area) is reached.
- ☛ The angled rubber at the distal end can easily expand while the air is pumped into the product.
4. Using the angling mechanism, move the distal end on the product's control body in all directions as far as possible while observing the pressure display of the leakage tester.
- ☛ The air pressure in the product can drop by up to 10 mmHg when moving the distal end.

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
5. After stabilizing the air pressure in the product, wait at least 3 minutes and check on the leakage tester whether the air pressure in the product drops or whether bubbles are continuously escaping along the product.
- ☛ If the air pressure drops or bubbles emerge from the product, the product is leaking and may no longer be used, immersed in liquid or reprocessed by machine. Complete the repairs form. Clean the product by wiping to remove organic soiling/residue and mark it as not disinfected. Contact KARL STORZ and send the product to KARL STORZ for repair. Line the carrying case for the product with the transport protective film 13990 SFN.
- ☛ If the air pressure does not drop and no bubbles emerge from the product, the product is leakproof. Remove the product from the liquid and then press the vent button on the pressure tester. Wait until the air pressure has dropped to 0 mmHg. The product can now be cleaned.

#### 14.4 Disassembly

Before cleaning and disinfection, remove the adaptor and other detachable parts from the light inlet piece, the instrument channel and the stopcock in their parts.

#### 14.5 Manual reprocessing

##### 14.5.1 Manual cleaning

1. If relevant, connect the adaptor.
2. Completely immerse the product in the cleaning solution.
3. If relevant, fill the lumen with cleaning solution using the syringe and ensure bubble-free wetting.
4. If relevant, remove the adaptor.
5. Observe the exposure time as specified by the chemical manufacturer during the cleaning process.
6. Brush the surfaces of the handle at least 3 times until no residue can be seen anymore.

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7. The mobile parts must be brushed.
8. If relevant, clean the valve seats by brushing them and turning them 360° 3 times until no residue can be seen anymore.
9. Clean the surface of the sheath with a low-lint cloth until no residue can be seen anymore.
10. Brush the lumen at least once as necessary until the brush head is visibly clean. Insert the brush into the proximal end of the channel until it protrudes from the distal end, then pull the brush out completely from the distal end. Clean the brush heads between each stroke of the brush.
11. If relevant, connect the adaptor.
12. If relevant, irrigate the lumen with cleaning solution 3 times.
13. After the exposure time, remove the product from the cleaning solution.
14. If relevant, irrigate the lumen with air.
15. If relevant, remove the adaptor.
16. Dry the surface with a soft, low-lint cloth.
17. If relevant, connect the adaptor.
18. Completely immerse the product in water of the best possible quality.
19. If relevant, fill the lumen with water using the syringe and ensure bubble-free wetting.
20. If relevant, remove the adaptor.
21. Remove the product from the water.
22. Dry the product.

**ATTENTION:** The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

#### 14.5.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. The working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

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
#### 14.5.3 Manual High-Level Disinfection (HLD)

The endoscope must be completely immersed in a high-level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

#### 14.5.4 Drying

Finally, the surfaces and channels/lumens are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

**ATTENTION:** The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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#### 14.6 Machine reprocessing

The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'.

##### 14.6.1 Manual preparation

1. If relevant, connect the adaptor.
2. Completely immerse the product in the cleaning solution.
3. If relevant, fill the lumen with cleaning solution using the syringe and ensure bubble-free wetting.
4. If relevant, remove the adaptor.
5. Observe the exposure time as specified by the chemical manufacturer during the cleaning process.
6. Brush the surfaces of the handle at least 3 times until no residue can be seen anymore.
7. The mobile parts must be brushed.
8. If relevant, clean the valve seats by brushing them and turning them 360° 3 times until no residue can be seen anymore.
9. Clean the surface of the sheath with a low-lint cloth until no residue can be seen anymore.
10. Brush the lumen at least once as necessary until the brush head is visibly clean. Insert the brush into the proximal end of the channel until it protrudes from the distal end, then pull the brush out completely from the distal end. Clean the brush heads between each stroke of the brush.
11. If relevant, connect the adaptor.
12. If relevant, irrigate the lumen with cleaning solution 3 times.
13. After the exposure time, remove the product from the cleaning solution.

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14. If relevant, irrigate the lumen with air.
15. If relevant, remove the adaptor.
16. Dry the surface with a soft, low-lint cloth.
17. If relevant, connect the adaptor.
18. Completely immerse the product in water of the best possible quality.
19. If relevant, fill the lumen with water using the syringe and ensure bubble-free wetting.
20. If relevant, remove the adaptor.
21. Remove the product from the water.

#### 14.6.2 Connecting

In order to ensure effective machine cleaning and disinfection, the endoscope's working and suction channel must be connected to the washer and disinfectant using the recommended adaptors.

**CAUTION:** The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

#### 14.6.3 Machine cleaning/ chemical-thermal disinfection

A maximum temperature of 60°C and a maximum pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

**NOTE:** If necessary, the instrument must be dried off afterwards by hand.

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#### 14.7 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

**i NOTE:** Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

**i NOTE:** Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connections) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

**! ATTENTION:** The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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**NOTE:** After treating with grease, the stopcocks must be sterilized separately using the fractionated steam sterilization process. Low-temperature sterilization is not possible on surfaces which have been lubricated and oiled.

- Afterwards, a functional check must be carried out.

**NOTE:** Check for a clear view and steerability.

- Clean the eyepiece and lens with a soft cloth, sponge or cotton-tipped applicator soaked in 70% alcohol.

**CAUTION:** Cleaning with the special cleaning paste should only be carried out if the image as viewed through the endoscope is visibly cloudy or if there are visible deposits on the optical surfaces, it should certainly not be used routinely during reprocessing.

Stubborn soiling and deposits can be removed with cleaning paste (Art. no. 27661). For this purpose, some cleaning paste is put on a clean, damp cleaning bud and distributed on the glass surfaces. Polishing ensures that even stubborn deposits are removed.

**NOTE:** The instrument must be reprocessed again following treatment with the cleaning paste!

Finally, dry the surfaces of the fibers and the optical end faces with a cotton-tipped applicator soaked in 70% isopropyl alcohol.

**NOTE:** The fiber area in the light inlet piece in particular ⑦ must be dried carefully with alcohol. Disinfectant and cleaning agent residue in the light inlet piece can become burnt on when the light cable is connected, thereby seriously impairing light transmission.

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#### 14.9 Sterilization

**WARNING:** Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

**WARNING:** Risk of infection: Sterilize the medical device in a disassembled state.

**ATTENTION:** Flexible endoscopes must never be steam sterilized (autoclaved).

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**NOTE:** All deviations from the recommended sterilization parameters must be validated by the user.

**NOTE:** When using reprocessing devices, please observe the device manufacturers' instructions.

**NOTE:** Sterilization is not possible on surfaces which have been lubricated and oiled.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

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**ATTENTION:** Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

**NOTE:** Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle with/without booster
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

#### 14.9.2 Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) sterilization – STERIS® V-PRO®

**NOTE:** The 'STERIS Low-Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® 1
- V-PRO® 1 Plus Lumen Cycle
- V-PRO® maX Flexible Cycle
- V-PRO® 60 Flexible Cycle

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#### 14.9.4 Chemical low-temperature sterilization with peracetic acid – STERIS System 1E®

Detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

**ATTENTION:** The pressure compensation cap must not be in place on the vent port.

- The glide surfaces of the stopcock must be lubricated locally with special grease.

#### 14.10 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

#### 14.11 Storage

**The storage method must be selected taking into account the respective applicable national requirements.**

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

**ATTENTION:** Do not expose flexible endoscopes to direct sunlight or excessive heat.

**ATTENTION:** The distal tip of the flexible endoscope must be in the neutral position during storage.

**ATTENTION:** Do not store the flexible endoscope in the case.

**ATTENTION:** The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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**14.12 Reprocessing of accessories**

Extensive accompanying documentation is available for the accessories.  
It is essential that you refer to these documents for further information on handling and reprocessing.

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## 15 Standard compliance

### 15.1 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

### 15.2 Responsibility

As the supplier of this instrument, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions:

assembly, system expansion, readjustment, modification, or repair have been performed by persons authorized by KARL STORZ and that the instrument is used in accordance with the instruction manual.

### 15.3 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see chapter 'Subsidiaries'), even during the warranty period.

Opening the device or performance of any repairs or modifications to the device by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair or modification performed during the warranty period shall void all warranty.

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#### 15.4 Directive compliance

This medical device bears the CE mark in accordance with the Medical Devices Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

The products listed here may not yet be available in all countries due to differences in approval requirements.

#### 16 Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new device and are only charged the repair costs for the defective device. For telescopes, you receive a warranty of 1 year, for instruments 1/2 year.

For fiberscopes and devices, individual repair is necessary. To bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In Germany, you can refer repairs direct to:

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Straße 34  
78532 Tuttlingen

In other countries, please contact your local KARL STORZ subsidiary or authorized dealer.

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### **16.1 Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among personnel). We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no warranty on the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga, Ontario  
L5N 3R3 Canada  
Phone: +1 905 816-4500, Fax: +1 905 816-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8525  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
1 South Los Carneros Road  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopy Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.  
Av. Ejercito Nacional No. 453 Piso 2,  
Colonia Granada, Alcaldia Miguel Hidalgo  
C.P. 11520 Ciudad de México  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Stamveien1  
1483 Hagan, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Sörstärågränd 14  
127 39 Skärholmen, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Taivaltie 5  
01610 Vantaa, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)96824775  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG  
Representative Office  
Zalgirio St. 94, LT9300  
Vilnius, Lithuania  
Phone: +370 5 272 0448  
Mobile: +370 685 67 000  
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B.V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.  
Phone: +31 (0)33 4545890  
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S.A.S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S.A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info-hrv@karlstorz.com

KARL STORZ Endoskopija d.o.o.  
Cesta v Gorice 34b  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gliwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited  
11<sup>th</sup> Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG  
Interchange 21 Tower, Level 32, Unit 3230  
399 Sukhumvit Road,  
North Klongtoey, Wattana,  
10110 Bangkok, Thailand  
Phone: +66 2 660 3669  
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG  
Resident Representative Office  
14<sup>th</sup> Floor, MPlaza Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
100022, Beijing, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch  
Room 701A Building 5 & Room 501 Building 7,  
No. 3000 Longdong Avenue, Pilot Free Trade Zone,  
201203, Shanghai, People's Republic of China  
Phone: +86 21 60339888, Fax: +86 21 60339808  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch  
Room 803-805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
610041, Chengdu, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
110014, Shenyang, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch  
Room 02B & 03 & 04A, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
510620, Guangzhou, People's Republic of China

Phone: +86 20 87321281, Fax: +86 20 87321286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG  
Representative Office Indonesia  
Sinarmas MSIG Tower Level 37  
Jl. Jend. Sudirman No. Kav. 21  
Jakarta Selatan  
DKI Jakarta 12920  
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.  
9F Hyowon-Building  
97, Jungdae-ro, Songpa-gu  
05719 Seoul, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: info-tw@karlstorz.com

KARL STORZ Endoscopy Philippines Inc.  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy New Zealand Ltd.  
31 Morningside Drive Mt Albert  
Auckland, 1025, New Zealand  
PO Box 56 511, Dominion Rd  
Auckland, 1446, New Zealand  
Phone: +64 9 846 6044  
Toll free: +64 508 84 84 84 (New Zealand only)  
E-Mail: sales-nz@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.  
68 Waterloo Road, Macquarie Park NSW 2113, Australia  
PO Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: karlstorz@karlstorz.com.au

www.karlstorz.com

## **KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 7461 708-0  
Telefax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)