# COULTER DxH Retic Pack

**REF** 628021

## TABLE OF CONTENTS

English . . . . .	2
Français (FR) . . . . .	4
Deutsche (DE) . . . . .	6
Italiano (IT) . . . . .	8
Español (ES) . . . . .	10
Português Portugal (PT-PT) . . . . .	12
Dansk (DK) . . . . .	14
Svenska (SE) . . . . .	16
Norsk (NO) . . . . .	18
Ελληνικά (GR) . . . . .	20
日本語 (JP) . . . . .	22
Lietuvių (LT) . . . . .	24
Magyar (HU) . . . . .	26
Polski (PL) . . . . .	28
Čeština (CZ) . . . . .	30
Slovák (SK) . . . . .	32
한국어 (KR) . . . . .	34
Türkçe (TR) . . . . .	36
Русский (RU) . . . . .	38
Hrvatski (HR) . . . . .	40
Български (BG) . . . . .	42
中文 (台灣) (ZH-TW) . . . . .	44
Română (RO) . . . . .	46
Srpski (Latinski) (SP) . . . . .	48
Українська (UA) . . . . .	50
Português Brasil (PT-BR) . . . . .	52
Nederlands (NL) . . . . .	54
Tiếng Việt (VN) . . . . .	56
REFERENCES . . . . .	58

# COULTER DxH Retic Pack

**REF** 628021

**For *In Vitro* Diagnostic Use**  
**Rx Only in the U.S.A.**

## INTENDED USE

For use on UniCel DxH Coulter Cellular Analysis Systems for clearing red cells and staining reticulocytes.

## INTENDED USER

This device is intended to be used by a laboratory professional.

## SUMMARY

Reticulocytes are immature erythrocytes (red blood cells) that retain a small amount of RNA as the nucleus is extruded during the maturation process. To detect the presence of this RNA, the red blood cells are stained while they are still viable in a process called supravital staining. Once a whole-blood sample is prepared, reticulocytes are differentiated from mature red blood cells and other cell populations using VCSn and reticulocyte counting technologies.

## PRINCIPLE

New Methylene Blue (Reagent A), a supravital dye, is incubated with a whole-blood sample. The dye precipitates the basophilic RNA network found in reticulocytes.<sup>1</sup> A clearing reagent (Reagent B) removes hemoglobin and unbound stain from the erythrocytes (red blood cells or RBCs) while keeping the cell and its membranes intact and without removing the precipitated dye - RNA complex. This process results in clear spherical mature RBCs and darkly stained reticulocytes that are ready for analysis. The total RBC population is analyzed on a hematology analyzer with reticulocyte counting technology using impedance, radio frequency, and VCSn technology then reticulocytes are differentiated from mature red blood cells and other cell populations.

## REACTIVE INGREDIENTS

The COULTER DxH Retic Pack contains the following:

### REAGENT A - Retic Stain

New Methylene Blue: 0.06% w/w

### REAGENT B - Retic Clear

Sulfuric Acid: 0.15% w/w

Potassium Thiocyanate: 0.3% w/w

## NON-REACTIVE INGREDIENTS

Stabilizers and Preservatives

## WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINERS.

GHS HAZARD CLASSIFICATION		
REAGENT B		H314: Causes severe skin burns and eye damage.
		P280: Wear protective gloves, protective clothing, and eye/face protection.
		P303+P361+P353: IF ON SKIN (or hair): Rinse skin with water.
		P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
		P310: Immediately call a POISON CENTER or doctor/physician.
		Sulfuric Acid 0.1-0.5%

## STORAGE, STABILITY, AND DISPOSAL

- Store COULTER DxH Retic Pack at 2-30°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

## MATERIALS PROVIDED

COULTER Retic Pack

## MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the UniCel DxH Coulter Cellular Analysis Systems IFU.
- System HELP or Instructions for Use.

## PREPARATION

Reagents in the COULTER DxH Retic Pack are ready to use. Replace the reagent kit as directed in your instrument product manuals and/or online help.

To avoid cross contamination when replacing the DxH Retic Pack, transfer pickup tubes one at a time beginning with Reagent A (Retic Stain) and then Reagent B (Retic Clear).

**IMPORTANT:** If product has been partially or completely frozen, allow product to warm to room temperature. Mix product by gentle inversion prior to placement on the instrument. Install and prime the reagent kit as directed in your instrument product manuals and/or online help. Verify background counts are acceptable before analyzing patient samples.

## PRODUCT AVAILABILITY

COULTER DxH Retic Pack

628021

(Contains 1 x 380 mL Reagent A and 1 x 1,900 mL Reagent B)

## TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

## ADDITIONAL INFORMATION

Glossary of Symbols is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs) by searching for part number C07163.

May be covered by one or more pat. - see [www.beckmancoulter.com/patents](http://www.beckmancoulter.com/patents)

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

Safety Data Sheet is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs)

## NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

## REVISION HISTORY

Revision AA, 04/2022

- Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes to product or instructions for use.

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Phiên bản AB, 03/2023
Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.