COULTER S-CAL Calibrator Kit

REF 628026

TABLE OF CONTENTS

English	2
Français (FR)	5
Deutsche (DE)	8
Italiano (IT)	11
Español (ES)	14
Português Portugal (PT-PT)	17
Dansk (DK)	20
Svenska (SE)	23
Norsk (NO)	26
Ελληνικά (GR)	29
日本語 (JP)	32
Lietuvių (LT)	35
Magyar (HU)	38
Polski (PL)	41
Čeština (CZ)	44
Slovák (SK)	47
한국어 (KR)	50
Türkçe (TR)	53
Русский (RU)	56
Hrvatski (HR)	59
中文 (台灣) (ZH-TW)	62
Română (RO)	65
Srpski (Latinski) (SP)	68
Українська (UA)	71
Português Brasil (PT-BR)	74
Nederlands (NL)	77
Tiếng Việt (VN)	80
REFERENCES	83

COULTER S-CAL Calibrator Kit

REF 628026

**For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.**

INTENDED USE

S-CAL Calibrator is intended for the determination of calibration factors for UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents. Refer to the Instructions for Use or System HELP.

SUMMARY AND PRINCIPLE

Coulter Cellular Analysis Systems, for which this kit was designed, require a calibrator to convert electronic measurements of each sample into accurate results expressed in clinical terms.

The parameters (WBC, RBC, HGB, MCV, PLT and MPV) are calibrated using S-CAL Calibrator, a stabilized human-blood preparation. S-CAL Calibrator is an acceptable alternative to whole-blood calibration.

The calibration procedure uses replicate measurements of S-CAL Calibrator. The average result for each parameter is divided into the S-CAL Calibrator Assigned Value to give a calibration factor.

HCT, MCH, MCHC, RDW, NRBC, Retic and DIFF parameters do not require calibration.

REAGENTS

S-CAL Calibrator consists of treated, stabilized, human erythrocytes and platelet sized components in an isotonic medium. Fixed erythrocytes are added to simulate leukocytes.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced materials: Human, Avian, Reptile and Ungulate.

Each human donor unit used in preparation of the material was tested by an FDA approved method for the presence of antibodies to Human Immunodeficiency Virus (HIV-1 and HIV-2) and Hepatitis C Virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative (were not repeatedly reactive.)

Because no test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV-1 and HIV-2), or other infectious agents are absent, this specimen/reagent should be handled at Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 1988.

Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (8/16/76).

STORAGE, STABILITY, AND DISPOSAL

The S-CAL Calibrator kit is shipped in a thermally insulated container designed to keep it cool. When stored at 2-8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. When the tubes are opened, the product is stable for 1 hour.

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (darkly-colored supernatant) is indicative of product deterioration. However, a slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

MATERIALS PROVIDED

- COULTER S-CAL Calibrator.
- Table of Expected Results.

MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient COULTER reagents to complete calibration.
- Instructions for Use or System HELP.
- 6C Cell Control Kit.

INSTRUCTIONS FOR USE


NOTE: Keep S-CAL Calibrator in the refrigerator until ready to begin Section II of INSTRUCTIONS FOR USE.

1. Before calibrating the instrument, refer to the appropriate Instructions for Use or System HELP, then go on to Section II.
2. DO NOT PROCEED IF INSTRUMENT PERFORMANCE DOES NOT MEET SPECIFICATIONS. CALL BECKMAN COULTER CUSTOMER SERVICE OR CONTACT YOUR LOCAL BECKMAN COULTER REPRESENTATIVE.
3. Calibration errors can be caused by:
 - Inadequate mixing
 - Excessive handling
 - Exceeding time limits
 - Failure to analyze a sufficient number of samples

IMPORTANT: Perform the following procedures exactly as described. Complete the entire procedure within 1 hour.

4. Remove S-CAL calibrator kit from the refrigerator and warm at AMBIENT TEMPERATURE for 15 minutes.
5. While the calibrator is equilibrating to ambient temperature, prepare the analyzer as follows:
 - A. Select Calibration Setup from CBC Calibration Screen in QA Menu.
 - B. Scan 2D barcode from the ASSAY SHEET to upload the assigned values.
6. After warming the S-CAL Calibrator, mix by hand as follows:

NOTE: Do **NOT** use a mechanical mixer.

- A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.
- 
- D. Repeat steps 6A through 6C.
 - E. Inspect the tube contents to determine if all cells have been uniformly distributed.
 - F. Repeat the mixing procedure if tube contents have not been uniformly distributed.
7. Analyze S-CAL Calibrator 10 times in Cassette Presentation Mode, following Instructions for Use or System HELP.

CHECK CALIBRATION CRITERIA

The DxH systems check the following criteria automatically and flag the calibration results with a yellow background if you need to calibrate a parameter or with a red background if there may be a problem.

Check calibration results. Calibration is required if the Difference is greater than DIFFERENCE LOWER LIMIT MAX and the Factor % Diff is greater than the FACTOR % DIFFERENCE LOWER LIMIT MAX.

PARAMETER	CALIBRATE IF DIFFERENCE LOWER LIMIT MAX	CALIBRATE IF FACTOR % DIFF LOWER LIMIT MAX
WBC	>0.1	>1.25%
RBC	>0.03	>0.7%
Hgb	>0.1	>0.78%
MCV	>1.0	>1.18%
PLT	>6.0	>2.70%
MPV	>0.5	>5.0%

VERIFY CALIBRATION

1. Analyze 6C Cell Controls according to package insert instructions.
2. Your results for the cell control levels should be within the expected range as stated on the TABLE OF EXPECTED RESULTS in the 6C Cell Control kit.
3. If your cell control results are not within the expected range, repeat the sample. If the results of the second sample are not within the expected range, follow the Troubleshooting Procedure on the 6C Cell Control package insert and call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ASSIGNED VALUES

The assigned values for WBC, RBC, HGB, MCV, PLT and MPV parameters are derived from sufficient samples to recover mean values with a high degree of statistical validity. Measurements are made on each of several COULTER instruments.

All instruments are traceable to whole-blood calibration and operated according to the appropriate Instructions for Use. Calibration stability is verified by comparison of instrument results and Reference Method results for fresh whole-blood specimens. System performance is monitored daily using 6C Cell Control at multiple levels.

Whole-blood calibration consists of K₃EDTA anticoagulated whole blood analyzed using the following reference methods:

WBC and RBC: A single-aperture impedance cell counter such as a COULTER Z Series cell counter and the manufacturer's recommended reagents. Macro dilutions are made using Class A glassware. Both WBC and RBC data are corrected for coincidence. The erythrocytes and leukocytes measurands (analytes) in this calibrator are traceable to Clin Lab Haemat, 1994, 16(2):131-138.

HGB: The hemoglobin measurand (analyte) in this calibrator is traceable to CLSI Standard H15-A3.¹ This hemiglobincyanide spectrophotometric procedure employs modified Drabkins (Ziljstra) Reagent and is referenced to NIST certified filters.

MCV: MCV is calculated: PCV/RBC x 10. Packed-cell volume (PCV) is measured by a hematocrit procedure according to CLSI Standard H7-A3.²

PLT: Phase contrast microscopy.

MPV: Latex particles (PN 7531820).

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the applicable Beckman Coulter reagent systems indicated in the calibrator Table of Expected Results (Assay Sheet).

Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

QUALITY CONTROL

After completing the calibration and verification procedures, the instrument is within the accuracy limit stated in your Product Manuals. Good laboratory practices recommend that each series of patient samples be preceded and followed by a quality-control check. Use COULTER cell controls to check the performance of the instrument

PRODUCT AVAILABILITY

COULTER S-CAL Calibrator Kit

628026 - 1 x 3.3 mL

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

Safety Data Sheet is available at beckmancoulter.com/techdocs

REVISION HISTORY

Revision AJ, 01/2018

Changes were made to:

- Added new Serbian language
- Added new Vietnamese language

Revision AK, 09/2018

Changes were made to:

- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

Revision AL, 05/2023

- Updated EC REP Address
- Updated copyright date

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REFERENCES

1. Clinical and Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood- Third Edition; Approved Standard. CLSI publication H15-A3. Wayne, PA, 2000.
2. Clinical and Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhematocrit Method- Third Edition; Approved Standard. CLSI publication H7-A3. Wayne, PA, 2000.

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.