# Shiley™




## Neonatal, Pediatric and Long Pediatric Tracheostomy Tube Cuffless

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de uso
ru	Инструкция по применению
zh	使用说明
pl	Instrukcja użycia
cs	Návod k použití
sk	Návod na použitie
sl	Navodila za uporabo
hu	Használati utasítás
el	Οδηγίες χρήσης
tr	Kullanma Talimatı
ar	تعليمات الاستخدام
bg	Инструкции за употреба
ro	Instrucțiuni de utilizare
ko	사용 설명서

**REF** NEF

**REF** PEF

**REF** PELF

REF			B = 0 $\theta = 120^\circ$ (105°–120°)							
(1) Part Number	(2) ID (mm)	(3) OD (mm)	(4) Length (mm)					(5) Cannula (mm) A + B + C	(6) Obturator	(7) Neck Strap
			A	C	D	E				
2.5NEF	2.5	4.2	8.8	19.2	16.5	16.9		28	X	X
3.0NEF	3.0	4.8	9.6	20.4	17.5	18.1		30	X	X
3.5NEF	3.5	5.4	10.3	21.7	18.5	19.2		32	X	X
4.0NEF	4.0	6.0	11.1	22.9	19.4	20.4		34	X	X
4.5NEF	4.5	6.7	11.9	24.1	20.4	21.6		36	X	X
2.5PEF	2.5	4.2	11.9	26.1	23.1	23.4		38	X	X
3.0PEF	3.0	4.8	12.4	26.6	23.4	23.9		39	X	X
3.5PEF	3.5	5.4	12.9	27.1	23.7	24.4		40	X	X
4.0PEF	4.0	6.0	13.4	27.6	24.0	25.0		41	X	X
4.5PEF	4.5	6.7	13.9	28.1	24.3	25.5		42	X	X
5.0PEF	5.0	7.3	14.8	29.2	25.2	26.7		44	X	X
5.5PEF	5.5	7.9	15.8	30.2	26.1	28.0		46	X	X
5.0PELF	5.0	7.3	20.8	29.2	27.5	32.3		50	X	X
5.5PELF	5.5	7.9	21.8	30.2	28.4	33.5		52	X	X
6.0PELF	6.0	8.5	22.8	31.2	29.3	34.8		54	X	X
6.5PELF	6.5	9.0	23.8	32.2	30.2	36.1		56	X	X

(1)	<p><b>en:</b> Part Number <b>fr:</b> Référence <b>de:</b> Artikelnummer <b>nl:</b> Onderdeelnummer <b>it:</b> Numero codice <b>es:</b> Número de pieza <b>sv:</b> Artikelnummer <b>da:</b> Delnummer <b>no:</b> Delnummer <b>fi:</b> Osanumero <b>pt:</b> Número da peça <b>ru:</b> Каталогный номер <b>zh:</b> 产品编号 <b>pl:</b> Numer katalogowy <b>cs:</b> Číslo dílu <b>sk:</b> Číslo dielu <b>sl:</b> Številka dela <b>hu:</b> Cikkszám <b>el:</b> Κωδικός είδους <b>tr:</b> Parça Numarası <b>ar:</b> رقم الجزء <b>bg:</b> Каталоген номер <b>ro:</b> Cod piesă <b>ko:</b> 품번</p>
(2)	<p><b>en:</b> ID – Inside Diameter (mm) <b>fr:</b> DI – Diamètre intérieur (mm) <b>de:</b> ID – Innendurchmesser (mm) <b>nl:</b> ID – binnendiameter (mm) <b>it:</b> ID – Diametro interno (mm) <b>es:</b> DI – Diámetro interno (mm) <b>sv:</b> ID – Innerdiameter (mm) <b>da:</b> ID – Indvendig diameter (mm) <b>no:</b> ID – Indre diameter (mm) <b>fi:</b> ID – sisäläpimitta (mm) <b>pt:</b> DI – Diâmetro interno (mm) <b>ru:</b> ID — внутренний диаметр (мм) <b>zh:</b> ID - 内径 ( mm ) <b>pl:</b> ID – średnica wewnętrzna (mm) <b>cs:</b> ID – vnitřní průměr (mm) <b>sk:</b> ID – vnútorný priemer (mm) <b>sl:</b> ID – notranji premer (mm) <b>hu:</b> ID – belső átmérő (mm) <b>el:</b> ID – Εσωτερική διάμετρος (mm) <b>tr:</b> I.Ç. – İç Çap (mm) <b>ar:</b> القطر الداخلي - ID <b>bg:</b> ID – Вътрешен диаметър (mm) <b>ro:</b> DI – Diametru interior (mm) <b>ko:</b> ID – 내경(mm)</p>
(3)	<p><b>en:</b> OD – Outside Diameter (mm) <b>fr:</b> DE – Diamètre extérieur (mm) <b>de:</b> AD – Außendurchmesser (mm) <b>nl:</b> OD – buitendiameter (mm) <b>it:</b> OD – Diametro esterno (mm) <b>es:</b> DE – Diámetro externo (mm) <b>sv:</b> YD - Ytterdiameter (mm) <b>da:</b> OD – Udvendig diameter (mm) <b>no:</b> YD – Ytre diameter (mm) <b>fi:</b> OD – ulkoläpimitta (mm) <b>pt:</b> DE – Diâmetro externo (mm) <b>ru:</b> OD — наружный диаметр (мм) <b>zh:</b> OD - 外径 ( mm ) <b>pl:</b> OD – średnica zewnętrzna (mm) <b>cs:</b> OD – vnější průměr (mm) <b>sk:</b> OD – vonkajší priemer (mm) <b>sl:</b> OD – zunanji premer (mm) <b>hu:</b> OD – külső átmérő (mm) <b>el:</b> OD – Εξωτερική διάμετρος (mm) <b>tr:</b> D.Ç. – Dış Çap (mm) <b>ar:</b> القطر الخارجي - OD <b>bg:</b> OD – Външен диаметър (mm) <b>ro:</b> DE – Diametru exterior (mm) <b>ko:</b> OD – 외경(mm)</p>

(4)	<p><b>en:</b> Length (mm) <b>fr:</b> Longueur (mm) <b>de:</b> Länge (mm) <b>nl:</b> Lengte (mm) <b>it:</b> Lunghezza (mm) <b>es:</b> Longitud (mm) <b>sv:</b> Längd (mm), <b>da:</b> Længde (mm) <b>no:</b> Lengde (mm) <b>fi:</b> Pituus (mm) <b>pt:</b> Comprimento (mm) <b>ru:</b> Длина (мм) <b>zh:</b> 长度 ( mm ) <b>pl:</b> Długość (mm) <b>cs:</b> Délka (mm) <b>sk:</b> Dĺžka (mm) <b>sl:</b> Dolžina (mm) <b>hu:</b> Hossz (mm) <b>el:</b> Μήκος (mm) <b>tr:</b> Uzunluk (mm) <b>ar:</b> الطول (ملم) <b>bg:</b> Дължина (mm) <b>ro:</b> Lungime (mm) <b>ko:</b> 길이(mm)</p>
(5)	<p><b>en:</b> Cannula (mm) <b>fr:</b> Canule (mm) <b>de:</b> Kanüle (mm) <b>nl:</b> Canule (mm) <b>it:</b> Cannula (mm) <b>es:</b> Cánula (mm) <b>sv:</b> Kanyl (mm), <b>da:</b> Kanyle (mm) <b>no:</b> Kanyle (mm) <b>fi:</b> Kanyyli (mm) <b>pt:</b> Cânula (mm) <b>ru:</b> Канюля (mm) <b>zh:</b> 套管 ( mm ) <b>pl:</b> Kaniula (mm) <b>cs:</b> Kanyla (mm) <b>sk:</b> Kanyla (mm) <b>sl:</b> Kanila (mm) <b>hu:</b> Kanül (mm) <b>el:</b> Σωληνίσκος (mm) <b>tr:</b> Kanül (mm) <b>ar:</b> الكانيولا (ملم) <b>bg:</b> Канюла (mm) <b>ro:</b> Canulă (mm) <b>ko:</b> 캐놀라(mm)</p>
(6)	<p><b>en:</b> Obturator <b>fr:</b> Obturateur <b>de:</b> Obturator <b>nl:</b> Obturator <b>it:</b> Otturatore <b>es:</b> Obturador <b>sv:</b> Obturator <b>da:</b> Obturator <b>no:</b> Obturator <b>fi:</b> Obturaattori <b>pt:</b> Obturador <b>ru:</b> Обтуратор <b>zh:</b> 闭孔器 <b>pl:</b> Obturator <b>cs:</b> Obturátor <b>sk:</b> Obturátor <b>sl:</b> Obturator <b>hu:</b> Obturátor <b>el:</b> Επιπωματικό <b>tr:</b> Obtüratör <b>ar:</b> المُغَسِّد <b>bg:</b> Обтуратор <b>ro:</b> Obturator <b>ko:</b> 닫개</p>
(7)	<p><b>en:</b> Neck Strap <b>fr:</b> Attache cervicale <b>de:</b> Halteband <b>nl:</b> Nekband <b>it:</b> Fettuccia di fissaggio <b>es:</b> Cinta para el cuello <b>sv:</b> Halskrage <b>da:</b> Halsrem <b>no:</b> Halsstropp <b>fi:</b> Kiinnityshihna <b>pt:</b> Fita para o pescoço <b>ru:</b> Шейная лента <b>zh:</b> 颈圈 <b>pl:</b> Pasek szyjny <b>cs:</b> Krční popruh <b>sk:</b> Krčný popruh <b>sl:</b> Trak za pritrnitev na vrat <b>hu:</b> Nyakpánt <b>el:</b> Ιμάντας λαιμού <b>tr:</b> Boyun Şeridi <b>ar:</b> شريط الرقبة <b>bg:</b> Лента за шия <b>ro:</b> Suport pentru gât <b>ko:</b> 목 스트랩</p>

<div data-bbox="298 1437 336 1507"><b>REF</b></div> <div data-bbox="303 1380 329 1429">NEF</div>	<p><b>en:</b> Neonatal Tracheostomy Tube Cuffless <b>fr:</b> Tube de trachéostomie néonatal sans ballonnet  <b>de:</b> Tracheostomiekannüle für Neugeborene, ohne Cuff <b>nl:</b> Neonatale tracheostomietube zonder cuff  <b>it:</b> Cannula tracheostomica neonatale non cuffiata <b>es:</b> Tubo de traqueotomía neonatal sin manguito  <b>sv:</b> Neonatal trakeostomitub utan kuff <b>da:</b> Neonatal trakeostomitube, uden manchete <b>no:</b> Neonatal trakeostomitube uten manskjett <b>fi:</b> Vastasyntyneen trakeostoomakanyyli, ilman mansettia <b>pt:</b> Tubo de traqueostomia neonatal sem manguito <b>ru:</b> Трахеостомическая трубка для новорожденных, без манжеты <b>zh:</b> 新生儿气管造口导管 (无袖带) <b>pl:</b> Rurka tracheostomijna noworodkowa, bez mankietu <b>cs:</b> Neonatální tracheostomická trubice, bez manžety <b>sk:</b> Novorodenecká tracheostomická kanyla, bez manžety <b>sl:</b> Neonatalna cevka za traheostomijo brez manšete <b>hu:</b> Újszülöttkekhez való használatra szolgáló tracheosztómias tubus mandzsetta nélkül <b>el:</b> Νεογνικός σωλήνας τραχειοστομίας χωρίς αερωθάλαμο <b>tr:</b> Neonatal Trakeostomi Tüpü, Kafsız <b>ar:</b> دون كفة <b>bg:</b> Трахеостомна тръба за новородени пациенти без маншет <b>ro:</b> Tub de traheostomie neonatal, neprevăzut cu manșon <b>ko:</b> 신생아용 기관절개 튜브, 커프 없음</p>
<div data-bbox="717 1437 756 1507"><b>REF</b></div> <div data-bbox="723 1380 749 1429">PEF</div>	<p><b>en:</b> Pediatric Tracheostomy Tube Cuffless <b>fr:</b> Tube de trachéostomie pédiatrique sans ballonnet  <b>de:</b> Tracheostomiekannüle für Kinder, ohne Cuff <b>nl:</b> Pediatrische tracheostomietube zonder cuff <b>it:</b> Cannula tracheostomica pediatrica non cuffiata <b>es:</b> Tubo de traqueotomía pediátrico sin manguito <b>sv:</b> Pediatrisk trakeostomitub utan kuff <b>da:</b> Pædiatrisk trakeostomitube, uden manchete <b>no:</b> Pediatrisk trakeostomitube uten manskjett <b>fi:</b> Lapsen trakeostoomakanyyli, ilman mansettia <b>pt:</b> Tubo de traqueostomia pediátrico sem manguito <b>ru:</b> Трахеостомическая трубка для детей, без манжеты <b>zh:</b> 儿童气管造口导管 (无袖带) <b>pl:</b> Rurka tracheostomijna pediatryczna, bez mankietu <b>cs:</b> Pediatrická tracheostomická trubice, bez manžety <b>sk:</b> Pediatrická tracheostomická kanyla, bez manžety <b>sl:</b> Pediatrična cevka za traheostomijo brez manšete <b>hu:</b> Gyermkekhez való használatra szolgáló tracheosztómias tubus mandzsetta nélkül <b>el:</b> Παιδιατρικός σωλήνας τραχειοστομίας, χωρίς αερωθάλαμο <b>tr:</b> Pediatrik Trakeostomi Tüpü, Kafsız <b>ar:</b> دون كفة <b>bg:</b> Трахеостомна тръба за педиатрични пациенти без маншет <b>ro:</b> Tub de traheostomie pediatric, neprevăzut cu manșon <b>ko:</b> 소아용 기관절개 튜브, 커프 없음</p>

<div data-bbox="502 1377 539 1513" data-label="Text"> <div data-bbox="502 1445 539 1513" data-label="Text">REF</div> <div data-bbox="502 1377 539 1442" data-label="Text">PELF</div> </div>	<div data-bbox="298 126 737 1351" data-label="Text"> <p><b>en:</b> Long Pediatric Tracheostomy Tube Cuffless <b>fr:</b> Tube de trachéostomie pédiatrique long sans ballonnet  <b>de:</b> Lange Tracheostomiekannüle für Kinder, ohne Cuff <b>nl:</b> Lange pediatrische tracheostomietube zonder cuff <b>it:</b> Cannula tracheostomica pediatrica lunga non cuffiata <b>es:</b> Tubo de traqueotomía pediátrico largo, sin manguito <b>sv:</b> Lång pediatriisk trakeostomitub utan kuff <b>da:</b> Lang pediatriisk trakeostomitube, uden manchet <b>no:</b> Lang pediatriisk trakeostomitube uten mansjett <b>fi:</b> Lapsen pitkä trakeostoomakanyyli, ilman mansettia <b>pt:</b> Tubo de traqueostomia pediátrico longo, sem manguito <b>ru:</b> Трахеостомиическая трубка для детей длинная, без манжеты <b>zh:</b> 儿童气管造口导管 (无袖带) <b>pl:</b> Rurka tracheostomijna pediatriczna długa, bez mankietu <b>cs:</b> Dlouhá pediatrická tracheostomická trubice, bez manžety <b>sk:</b> Dlhá pediatrická tracheostomická kanyla, bez manžety <b>sl:</b> Dolga pediatrična cevka za traheostomijo brez manšete <b>hu:</b> Gyermekekhez való használatra szolgáló, hosszú tracheosztómiás tubus mandzsetta nélkül <b>el:</b> Μακρὸς παιδιατρικὸς σωλήνας τραχειοστομίας, χωρίς αερωθάκη <b>tr:</b> Uzun Pediatrik Trakeostomi Tüpü, Kafsız <b>ar:</b> دون كفة ، طويل للأطفال الصغار <b>bg:</b> Дълга трахеостомна тръба за педиатрични пациенти без маншет <b>ro:</b> Tub de traheostomie pediatric lung, neprevăzut cu manșon <b>ko:</b> 긴 소아용 기관절개 튜브, 커프 없음</p> </div>
---	---

**en:** Identification of a substance that is not contained or present within the product or packaging  
**fr:** Identification d'une substance qui n'est ni contenue ni présente dans le produit ou dans l'emballage  
**de:** Identifikation einer Substanz, die nicht im Produkt oder in der Verpackung enthalten oder vorhanden ist  
**nl:** Aanduiding voor een stof die geen deel uitmaakt van en niet aanwezig is in het product en de verpakking  
**it:** Identificazione di una sostanza non contenuta né presente all'interno del prodotto o della confezione  
**es:** Identificación de una sustancia no contenida ni presente en el producto o en el envase  
**sv:** Identifiering av ett ämne som inte finns i produkten eller förpackningen  
**da:** Identifikation af et stof, der ikke er indeholdt eller til stede i produktet eller emballagen  
**no:** Identifikasjon av et stoff som ikke finnes i produktet eller pakningen  
**fi:** Tunnistaminen aineelle, joka ei sisälly tuotteen tai sen pakkaukseen  
**pt:** Identificação de uma substância não contida ou não presente no produto ou na embalagem  
**ru:** Обозначение вещества, не содержащегося в изделии или упаковке и не входящего в их состав  
**zh:** 表示产品或包装内不含或不存在的物质的标识  
**pl:** Identyfikacja substancji niezawartej i nieobecnej w produkcie i opakowaniu  
**cs:** Identifikace látky, která není obsažena nebo přítomna ve výrobku nebo obalu  
**sk:** Identifikácia látky, ktorá nie je obsiahnutá alebo prítomná vo výrobku alebo balení  
**sl:** Navedba snovi, ki jo izdelek ali embalaža ne vsebuje  
**hu:** Olyan anyag jelzése, amelyet sem a termék, sem a csomagolás nem tartalmaz  
**el:** Ταυτοποίηση ουσίας που δεν περιέχεται ή δεν υπάρχει στο προϊόν ή στη συσκευασία  
**tr:** Ürün veya ambalajda içerilmeyen veya bulunmayan madde işaretleri  
**ar:** التعرف على المادة غير المحتواة أو غير الموجودة داخل المنتج أو العبوة  
**bg:** Идентификация на вещество, което не се съдържа или не е налично в продукта или опаковката  
**ro:** Identificarea unei substanțe care nu este inclusă sau prezentă în produs sau ambalaj  
**ko:** 제품 또는 포장에 포함되어 있지 않거나 존재하지 않는 성분의 표시

Page Left Intentionally Blank

# Neonatal, Pediatric Tracheostomy Tubes

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## Description

The products have a radiopaque tube with a clear flange and standard 15mm connector for direct connection to standard ventilation and anesthesia equipment. Each tube is packaged with an obturator and cotton twill tie. The products are not made with natural rubber latex or DEHP.

## Indications For Use

This device is intended for use in providing tracheal access for airway management.

**Contraindications:** None known

## WARNINGS:

- **Carefully read these instructions prior to use.**
- **Insertion of a new tracheostomy tube or the replacement of an existing tracheostomy tube can be associated with patient discomfort, desaturation, or acute respiratory distress. During tube insertions and changes, patients should be in an appropriately monitored environment with trained care providers and the availability of appropriate airway devices.**
- **Sterile only if the protective package is not open, damaged or broken.**
- **Handle with due care to prevent user contamination prior to use.**
- **Avoid contact of tube with laser surgery beams or electrosurgery electrodes especially in the presence of oxygen enriched gas or nitrous oxide because this could result in combustion of the tube with harmful burns and emission of corrosive and toxic combustion products, including hydrochloric acid (HCl).**
- **Replacement tube and obturator should be kept at the patient's bedside.**
- **Do not use excessive force to connect accessories to the 15mm connector. This may result in difficult disconnection.**

- **Capping or plugging the tube is not recommended because there may not be sufficient airway for the patient to breathe.**
- **Excessive amounts of lubricant can dry on the inner surface of the tracheostomy tube resulting in either a plug or film that may partially or totally block the airway.**
- **Use expert clinical judgement in the selection of a tube for a patient with consideration to the inner diameter, outer diameter, and length:**
  - **An inner diameter too small increases resistance through the tube and airway clearance is more difficult.**
  - **An outer diameter too small increases the cuff pressure required to seal the trachea leading to tracheal damage or persistent air leak and hypoventilation.**
  - **An outer diameter too large will be difficult to pass through the stoma.**
  - **A tube too short may fall out or press against the posterior tracheal wall and obstruct the distal tip or cause tracheal stenosis.**
  - **A tube too long may cause excessive coughing or become obstructed by the anterior wall of the trachea.**
- **The tube and accessories are supplied sterile and cannot be adequately resterilized for safe reuse and are intended for single patient use. The tube and accessories are for single use and single patient use only.**
- **Attempts to resterilize the tube and accessories may result in product failure and increased risk to the patient.**

#### **CAUTIONS:**

- **Federal law (USA) restricts this device to sale by or on the order of a physician.**
- **The tube and obturator are for single patient use. Duration should not exceed twenty-nine (29) days. Covidien has not substantiated use of these devices beyond 29 days. Decisions about tracheostomy tube changes should be made by the responsible physician or designate using accepted medical techniques and judgement.**
- **Use only with equipment having standard 15mm connectors.**
- **Prior to tube insertion, check that the obturator can be inserted and removed without difficulty.**
- **The device is MRI safe.**

## Instructions

---

1. The physician should determine whether a Shiley™ neonatal or pediatric tube should be used.
2. Insert the obturator into the proximal end of the tracheostomy tube. When completely seated it will protrude slightly from the distal tip.
3. A thin film of water soluble lubricant should be applied to the tube and protruding portion of the obturator to facilitate ease of insertion. Insert the tube with the obturator fully seated in place.
4. After the tube is properly in place, immediately remove the obturator and confirm proper ventilation. Store replacement tube and obturator in an accessible location near the patient should it be needed in the event of dislodgement.
5. Attach the neck strap to the flange eyelet. Flex the patient's neck forward and tie the neck strap. When properly adjusted, one finger space between the neck strap and the patient's neck should exist.
6. Follow current medical accepted procedures regarding the suctioning to remove secretions.
7. Prior to removing the tube from the patient, hold on to the tracheostomy tube's connector, then gently twist and pull the attached accessories to disconnect them from the tracheostomy tube.
8. Remove the tube from the patient following currently accepted medical techniques.
9. Discard the tracheostomy tube and accessories in accordance with applicable national regulations.

## Adverse Events

Adverse events associated with the use of tracheostomy tubes during the tracheotomy procedure, during cannulation, or after decannulation in descending order of severity and frequency are: obstruction of the tube, accidental decannulation, granuloma, infection, post-decannulation tracheocutaneous fistula, paratracheal air leak, stenosis, bleeding, trauma to surrounding tissue, suprastomal collapse, hyperventilation, hypercarbia, hypoxia, respiratory/cardiac arrest, toxicity, tracheal fistulas, cough, dysphagia, aspiration, and foreign body aspiration.

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**Rx**  
ONLY

**DEHP**  
Not made  
with DEHP

Do not use if  
package is opened  
or damaged

**CE**  
0123

Part No. PT00103700 Rev A 2019-05

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG.

<sup>TM</sup>® brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

© 2017 Covidien.

 Covidien llc, 15 Hampshire Street, Mansfield, MA 02048 USA.

 Covidien Ireland Limited, IDA Business & Technology Park, Tullamore, Ireland.

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