

## Vista 120 S Patient Monitoring Solution

Dräger understands the growing need for a patient monitor with built-in connectivity that provides essential monitoring at a good value. The Vista 120 S supports adult, pediatric and neonatal patients and can be used on its own or with a Dräger therapy device as a fully integrated workstation.


**305 mm (12.1") TFT color touchscreen** 5.1.1  
High resolution display (800 x 600) is bright and easy to read, even from a distance

**Medibus/Medubus-X connectivity**  
Offers true integrated workstation functionality across care segments

**Configurable layout**  
Lets you see the information you want, the way you want to see it

**Enhanced trending**

- Stores up to 120 hours of trend data for all parameters in tabular and graphic formats
- Stores up to 1,200 NIBP measurements and 200 alarm events



**Core set of essential parameters** 5.1.2  
3/5 lead ECG, SpO<sub>2</sub>, non-invasive blood pressure, respiration and dual temperature

**Networking capabilities**  
Enables central monitoring

**Alarms**  
Alarm indicator and alarm pause/off

**Shortcut keys**  
Fast access to main functions

## Benefits

---

### Part of a complete department solution

The Vista 120 S supports adult, pediatric and neonatal patients in a variety of care environments – including Intensive Care Units, Operating Rooms, Emergency Departments and Neonatal Intensive Care Units. Use the Vista 120 S on its own or with a complementary Dräger device, such as a ventilator, or anesthesia machine for a fully integrated workstation.

---

### Essential monitoring capabilities, exceptional value

The Vista 120 S displays up to eleven waveforms in an easy-to-configure layout and offers a core set of essential parameters including 3/5 lead ECG, non-invasive blood pressure, respiration and dual temperature comes standard. Advanced parameters including three invasive blood pressures, mainstream etCO<sub>2</sub>, cardiac output and sidestream etCO<sub>2</sub> are also available.

Users are free to add external parameter modules including SCIO for all models and CO<sub>2</sub> on model C/C+ after initial device purchase.

---

### Supports workflow efficiently and cost-effectively

The Vista 120 S is easy to learn and easy to use. You can configure the display to see the information you want to see, the way you want to see it. Fast access keys and simplified menus put the data you need right at your fingertips. It's light, portable and ready to go with an integrated bed hook for easy patient transport.

---

### Built-in recorder

The Vista 120 S has an integrated recorder that prints out up to three channels of information – saving time by providing documentation when and where you need it.

---

### Clear view of patient data

The Vista 120 S has a 305 mm (12.1") TFT color touchscreen that is bright and clear. Displaying comprehensive patient data with ease.

---

### Dräger heritage of quality

At Dräger, every life is unique. Protecting, supporting and saving lives is the foundation of our company philosophy. Our goal is to provide product and solutions that support acute care, help improve patient outcomes, reduce costs and achieve greater overall patient satisfaction.

## Related Products

D-68804-2012



### Vista 120 Central Monitoring System

The easy-to-use Vista 120 Central Monitoring System (CMS) lets you centrally monitor the vital signs of up to 64 patients connected to Vista 120/Vista 120 S bedside monitors. This central surveillance streamlines workflow for clinicians, while significantly increasing patient safety.

D-6829-2014



### Vista 120

Hospitals around the world share a common challenge – to provide the best possible care in locations with growing populations, stricter financial regulations and caregivers that are increasingly overloaded. The Vista 120 was engineered to meet your clinical needs and stay within your budget, allowing you to deliver efficient and high-quality patient care.

D-11590-2019



### Vista 120 SC

Reduce clinicians' workload with an easy-to-use and intuitive user interface. The Vista 120 SC is designed for spot check and continuous vital signs monitoring to complete Dräger's hospital-wide solution offerings.

## Technical Data

### Classification

Protection class	Class I equipment and internal powered equipment
Degree of protection against electric shock	CF: ECG (RESP), TEMP, IBP, C.O. BF: SpO <sub>2</sub> , NIBP, CO <sub>2</sub> , AG, BIS
Defibrillation protection	Yes
Liquid ingress protection	IPX1
Disinfection/sterilization method	Refer to chapter "Care and Cleaning" for details.
Mode of operation	Continuous
Compliant with standards	IEC 60601-1: 2005+A1: 2012; IEC 60601-1-2: 2014; EN 60601-1: 2006+A1: 2013; EN 60601-2: 2015; IEC 60601-2-49: 2011

### Supported Parameters

<b>ECG</b> <b>5.2.1</b>	
Lead mode	3-lead wire: I, II, III 5-lead wire: I, II, III, aVR, aVL, aVF, V
Waveform	3-lead wire: 1-channel waveform 5-lead wire: 2-channel waveform, max. seven waveforms
Lead naming style	AHA, IEC
Display sensitivity	1.25 mm/mV (x0.125), 2.5 mm/mV (x0.25), 5 mm/mV (x0.5), 10 mm/mV (x1), 20 mm/mV (x2), 40 mm/mV (x4), AUTO gain
Sweep	6.25, 12.5, 25, 50 mm/s
Bandwidth (-3 dB)	Diagnosis: 0.05 to 150 Hz Monitor: 0.5 to 40 Hz Surgery: 1 to 20 Hz
CMRR (Common mode rejection ratio)	Diagnostic: > 95 dB Monitor: > 105 dB Surgery: > 105 dB
Notch	50 Hz/60 Hz (Notch filter can be selected manually)
Differential input impedance	> 5 MΩ
Input signal range	±10 mVPP
Electrode offset potential tolerance	±800 mV
Auxiliary current (Leads off detection)	Active electrode: < 100 nA Reference electrode: < 900 nA
Recovery time after defibrillation	< 5 s (measured without electrodes as IEC60601-2-27:2011, Sect. 201.8.5.5.1 requires.)
Leakage current of patient	< 10 μA
Scale signal	1 mV <sub>PP</sub> , accuracy is ±5
System noise	< 30 μV <sub>PP</sub>
ESU protection	Cut mode: 300 W Coagulation mode: 100 W Recovery time: ≤ 10 s
ESU noise suppression	Tested according to the test method in ANSI/AAMI EC13-2002: Sect. 5.2.9.14, it accords with the standard
Minimum input slew rate (lead II)	> 2.5 V/s
<b>Pace Pulse</b>	
Pulse indicator	Pulse is marked if the requirements of IEC 60601-2-27: 2011, Sect. 201.12.1.101.12 are met: Amplitude: ±2 mV to ±700 mV Width: 0.1 ms to 2.0 ms Ascending time: 10 μs to 100 μs

## Technical Data

### Pulse rejection

Pulse is rejected if the requirements of IEC 60601-2-27: 2011, Sect. 201.12.1.101.13 are met: Amplitude:  $\pm 2$  mV to  $\pm 700$  mV  
Width: 0.1 ms to 2.0 ms Ascending time: 10  $\mu$ s to 100  $\mu$ s

### Heart Rate 5.2.2

Range	ADU: 15 to 300 bpm PED/NEO: 15 to 350 bpm
Accuracy	$\pm 1\%$ or $\pm 1$ bpm, whichever is greater
Resolution	1 bpm
Sensitivity	$\geq 300 \mu V_{PP}$

### PVC

Range	ADU: 0 to 300 PVCs/min PED/NEO: 0 to 350 PVCs/min
Accuracy	1 PVCs/min or 2% of measurement, whichever is greater
Resolution	1 PVCs/min

### ST Value

Range	-2.0 to 2.0 mV
Accuracy	-0.8 mV to +0.8 mV: $\pm 0.02$ mV or 10%, whichever is greater
Resolution	0.01 mV

### HR Averaging Method

Method 1	Heart rate is computed by excluding the minimum and maximum values from the 12 most recent RR intervals and averaging the residual 10 RR intervals.
Method 2	If each of three consecutive RR intervals is greater than 1,200 ms, then the four most recent RR intervals are averaged to compute the HR.

### Range of Sinus and SV Rhythm

Tachycardia	Adult: RR interval for 5 consecutive QRS complex $\leq 0.5$ s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex $\leq 0.375$ s.
Normal	Adult: $0.5$ s < RR interval for 5 consecutive QRS complex < $1.5$ s. Pediatric/neonatal: $0.375$ s < RR interval for 5 consecutive QRS complex < $1$ s.
Bradycardia	Adult: RR interval for 5 consecutive QRS complex $\geq 1.5$ s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex $\geq 1$ s.

### Range of Ventricular Rhythm

Ventricular tachycardia	The interval of 5 consecutive ventricular complexes is less than 600 ms
Ventricular rhythm	The interval of 5 consecutive ventricular complexes ranges from 600 ms to 1,000 ms
Ventricular bradycardia	The interval of 5 consecutive ventricular complexes is higher than 1,000 ms

### Startup Time for Tachycardia

Ventricular tachycardia 1 mV 206 bpm	Gain 0.5: 10 s Gain 1.0: 10 s Gain 2.0: 10 s
Ventricular tachycardia 2 mV 195 bpm	Gain 0.5: 10 s Gain 1.0: 10 s

## Technical Data

Response time of heart rate meter to change in HR	Gain 2.0: 10 s HR range: 80 to 120 bpm Range: Within 11 s HR range: 80 to 40 bpm Range: Within 11 s
Tall T-wave rejection	Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.17 minimum recommended 1.2 mV T-wave amplitude
Accuracy of heart rate meter and response to irregular rhythm	Complied with IEC 60601-2-27: 2011, Sect. 201.7.9.2.9.101 b) 4). The HR value after 20 s: Ventricular bigeminy: 80 ±1 bpm Slow alternating ventricular bigeminy: 60 ±1 bpm Rapid alternating ventricular bigeminy: 120 ±1 bpm Bidirectional systoles: 91 ±1 bpm
<b>Respiration</b>	
Method	Impedance between RA-LL, RA-LA
Baseline impedance range	200 Ω to 2,500 Ω (with ECG cables of 1 KΩ resistance)
Measuring sensitivity	Within the baseline impedance range: 0.3 Ω
Waveform bandwidth	0.2 to 2.5 Hz (-3 dB)
RR measuring and alarm range:	Adult: 0 to 120 rpm Neo/Ped: 0 to 150 rpm
Resolution	1 rpm
Accuracy	Adult: 6 rpm to 120 rpm: ±2 rpm 0 rpm to 5 rpm: not specified Neo/Ped: 6 rpm to 150 rpm: ±2 rpm 0 rpm to 5 rpm: not specified
Gain selection	x0.25, x0.5, x1, x2, x3, x4, x5
<b>NIBP 5.2.3</b>	
Method	Oscillometric
Mode	Manual, auto, continuous
Measuring interval in auto mode	1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240, and 480 min
Continuous	5 min, interval is 5 s
Measuring type	Systolic pressure, diastolic pressure, mean pressure
Alarm type	SYS, DIA, MAP
<b>Measuring and Alarm Range</b>	
Adult mode	SYS: 40 to 270 mmHg DIA: 10 to 215 mmHg MAP: 20 to 235 mmHg
Pediatric mode	SYS: 40 to 230 mmHg DIA: 10 to 180 mmHg MAP: 20 to 195 mmHg
Neonatal mode	SYS: 40 to 135 mmHg DIA: 10 to 100 mmHg MAP: 20 to 110 mmHg
Cuff pressure measuring range	0 to 300 mmHg
Pressure resolution	1 mmHg
Maximum standard deviation	8 mmHg
<b>Maximum Measuring Period</b>	
Adult/pediatric	120 s
Neonate	90 s

## Technical Data

Typical measuring period	20 to 35 s (depend on HR/motion disturbance)
<b>Overpressure Protection</b>	
Adult	297 ±3 mmHg
Pediatric	245 ±3 mmHg
Neonatal	147 ±3 mmHg
<b>PR</b>	
Measuring range	40 to 240 bpm
Accuracy	±3 bpm or 3.5%, whichever is larger
<b>SpO<sub>2</sub> 5.2.4</b>	
Measuring range	0 to 100%
Resolution	1%
<b>Accuracy</b>	
Adult (including pediatric)	±2% (70 to 100% SpO <sub>2</sub> ) Undefined (0 to 69% SpO <sub>2</sub> )
Neonate	±3% (70 to 100% SpO <sub>2</sub> ) Undefined (0 to 69% SpO <sub>2</sub> )
<b>PI</b>	
Measuring range	0–10, invalid PI value is 0
Resolution	1
<b>Pulse Rate</b>	
Pulse rate measuring range	25 to 300 bpm
Alarm range	30 to 300 bpm
Accuracy	±2 bpm
<b>Nellcor Module</b>	
Measuring range	1% to 100%
Alarm range	20% to 100%
Resolution	1%
Data update period	1 s
Accuracy (70% to 100% SpO <sub>2</sub> ):	
DS-100A, OXI-A/N (adult)	±3 %
OXI-A/N (neonate)	±4%
D-YS (infant to adult)	±3%
D-YS (neonate)	±4%
D-YS with D-YSE ear clip	±3.5%
MAX-FAST	±2%
<b>Pulse Rate</b>	
Measuring range	20 to 300 bpm
Resolution	1 bpm
Accuracy	3 bpm (20 to 250 bpm)
Sensor wavelength	Approximately 660 and 900 nm
Emitted light energy	<15 mW
<b>NOTE:</b> Information about the wavelength range can be especially useful to clinicians (for instance, when photodynamic therapy is performed).	
<b>Temperature</b>	
Channels	2
Measuring and alarm range	0 to 50°C

## Technical Data

		(32 to 122°F)	
Sensor type		YSI 10 k	
Resolution		±0.1°C (0.1°F)	
Accuracy (without sensor)		±0.1°C	
Refresh time		Every 1 to 2 s	
IBP			
Channels		3	
Accuracy		±2% or ±1 mmHg, whichever is greater	
Resolution		1 mmHg	
Pressure Sensor			
Sensitivity		5 (µV/V/mmHg)	
Impedance range		300 Ω to 3,000 Ω	
Filter		DC~ 12.5 Hz; DC~ 40 Hz	
Zero		Range: ±200 mmHg	
Measuring and Alarm Range			
Art		0 to 300 mmHg	
PA		-6 to 120 mmHg	
CVP/RAP/LAP/ICP		-10 to 40 mmHg	
P1/P2		-50 to 300 mmHg	
CO <sub>2</sub>			
Complies with ISO 80601-2-55: 2011			
G2 Module			
Intended patient		Adult, pediatric, neonatal	
Measure parameters		EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR	
Unit		mmHg, %, kPa	
Measuring range		CO <sub>2</sub>	0 mmHg to 150 mmHg (0% to 20%)
		AwRR	2 rpm to 150 rpm
Resolution		EtCO <sub>2</sub>	1 mmHg
		FiCO <sub>2</sub>	1 mmHg
		AwRR	1 rpm
Accuracy			
EtCO <sub>2</sub>	±2 mmHg, 0 mmHg to 40 mmHg	Respiratory rate ≤ 60 rpm	Typical conditions: <ul style="list-style-type: none"><li>– Ambient temperature: (25±3)°C</li><li>– Barometric pressure: (760±10) mmHg</li><li>– Balance gas: N<sub>2</sub></li><li>– Sample gas flow rate: 100 ml/min</li></ul>
	±5% of reading, 41 mmHg to 70 mmHg		
	±8% of reading, 71 mmHg to 100 mmHg		
	±10% of reading, 101 mmHg to 150 mmHg		
AwRR	±12% or ±4 mmHg of reading, whichever is greater	Respiratory rate > 60 rpm	All conditions
	±1 rpm		
Drift of measure accuracy		Meets the requirements of the measure accuracy	
Sample gas flow rate		70 ml/min or 100 ml/min (default), accuracy: ±15 ml/min	
Warm-up time		Display reading within 20 s; reach to the designed accuracy within 2 minutes	



## Technical Data

Rise time	< 400 ms (water trap with 2 m gas sampling tube, sample gas flowrate: 100 ml/min)
Response time	< 4 s (water trap with 2 m gas sampling tube, sample gas flowrate: 100 ml/min)
Work mode	Standby, measure
O <sub>2</sub> compensation	Range: 0% to 100% Resolution: 1% Default: 16%
N <sub>2</sub> O compensation	Range: 0% to 100% Resolution: 1% Default: 0%
AG compensation	Range: 0% to 20% Resolution: 0.1% Default: 0%
Humidity compensation method	ATPD (default), BTPS
Barometric pressure compensation	Automatic (The change of barometric pressure will not add additional errors to the measurement values.)
Zero calibration	Support
Calibration	Support
Alarm	EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR
Apnea alarm delay	10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s, 60 s; default value is 20 s
Data sample rate	100 Hz
EtCO <sub>2</sub> change <sup>1</sup>	AwRR >80 rpm, EtCO <sub>2</sub> descending 8% AwRR >120 rpm, EtCO <sub>2</sub> descending 10%

**NOTE:** Use a test device equivalent to EN ISO 80601-2-55 fig 201.101 to measure at 1:2 I/E ratio. Respiration rate accuracy is determined by frequency of device, and ET READING change refers to the nominal value.

### Interfering Gas Effects

Gas	Gas Level (%)	Quantitative Effect/Comments
Nitrous oxide	60	The interfering gas will have no effect on the measurement value if compensation of O <sub>2</sub> , N <sub>2</sub> O, anesthetic agents has been correctly set.
Halothane	4	
Enflurane	5	
Isoflurane	5	
Sevoflurane	5	
Desflurane	15	

### Respironics Module

Applicable patient type	Adult, pediatric and neonatal patients
Technique	Infra-red absorption technique
Measure parameters	EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR
Unit	mmHg, %, Kpa

### Measuring Range

EtCO <sub>2</sub>	0 mmHg to 150 mmHg	
FiCO <sub>2</sub>	3 mmHg to 50 mmHg	
AwRR	0 rpm to 150 rpm (mainstream)	
	2 rpm to 150 rpm (sidestream)	
Resolution	EtCO <sub>2</sub>	1 mmHg
	FiCO <sub>2</sub>	1 mmHg
	AwRR	1 rpm
EtCO <sub>2</sub> accuracy	±2 mmHg, 0 mmHg to 40 mmHg	
	±5% of reading, 41 mmHg to 70 mmHg	

## Technical Data

	±8% of reading, 71 mmHg to 100 mmHg
	±10% of reading, 101 mmHg to 150 mmHg
	±12% of reading, RR is over 80 rpm (sidestream). There will be no degradation in performance due to respiration rate (mainstream)
AwRR accuracy	±1 rpm
Operation mode	Measure, standby
Sample gas flow rate (sidestream)	(50 ±10) ml/min

### O<sub>2</sub> Compensation

Range	0% to 100%
Resolution	1%
Default	16%
Barometric pressure compensation	User setup

### Anesthetic Gas Compensation

Range	0% to 20%
Resolution	0.1%
Default	0.0%
Balance gas compensation	Room air, N <sub>2</sub> O, helium

### Stability

Short-term drift	Drift over 4 hours < 0.8 mmHg
Long-term drift	120 hours
Zero calibration	Support
Alarm type	EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR
Apnea alarm delay	10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s; default value is 20 s
Data sample rate	100 Hz
CO <sub>2</sub> rise time/response time (mainstream)	Less than 60 ms
Sensor response time (sidestream)	< 3 seconds, including transport time and rise time

### Interfering Gas and Vapor Effects on EtCO<sub>2</sub> Measurement Values

Gas or Vapor	Gas Level (%)	Quantitative Effect/Comments
Nitrous oxide	60	Dry and saturated gas
Halothane	4	(0 ~ 40) mmHg: ± 1 mmHg additional error
Enflurane	5	(41 ~ 70) mmHg: ± 2.5% additional error
Isoflurane	5	(71 ~ 100) mmHg: ± 4% additional error
Sevoflurane	5	(101 ~ 150) mmHg: ± 5% additional error
Xenon	80	*Additional worst case error when compensation for PB, O <sub>2</sub> , N <sub>2</sub> O, anesthetic agents, or helium is correctly selected for the actual fractional gas constituents present.
Helium	50	Desflurane:
Desflurane	15	The presence of desflurane in the exhaled breath at concentrations greater than 5% will positively bias carbon dioxide values by up to an additional 3 mmHg at 38 mmHg.
		Xenon:
		The presence of xenon in the exhaled breath will negatively bias carbon dioxide values by up to an additional 5 mmHg at

## Technical Data

38 mmHg.

### Barometric Pressure on EtCO<sub>2</sub> Measurement Values

#### Quantitative Effect

Ambient barometric, operational

(0 ~ 40) mmHg:  $\pm 1$  mmHg additional error

(41 ~ 70) mmHg:  $\pm 2.5\%$  additional error

(71 ~ 100) mmHg:  $\pm 4\%$  additional error

(101 ~ 150) mmHg:  $\pm 5\%$  additional error

\*Additional worst case error when compensation for PB, O<sub>2</sub>, N<sub>2</sub>O, anesthetic agents, or helium is correctly selected for the actual fractional gas constituents present.

**NOTE:** Respiration rate accuracy was verified by using a solenoid test setup to deliver a square wave of known CO<sub>2</sub> concentration to the device. 5% and 10% CO<sub>2</sub> concentrations were used. Respiration rate was varied over the range of the device. Pass/fail criteria was comparison of the respiratory rate output from the sensor to the frequency of the square wave.

#### C.O.

Intended patient

Adult

Measurement method

Thermodilution technique

#### Measuring Range

C.O.

0.1 l/min ~ 20 l/min

TB

23°C ~ 43°C

TI

-1°C ~ 27°C

#### Resolution

C.O.

0.1 l/min

TB, TI

0.1°C

#### Accuracy

C.O.

$\pm 5\%$  or  $\pm 0.2$  l/min, whichever is greater

TB

$\pm 0.1^\circ\text{C}$  (without sensor)

TI

$\pm 0.1^\circ\text{C}$  (without sensor)

#### Trend Review

Short

1 hr, 1 s. resolution

Long

120 hrs, 1 min. resolution

Review

1,200 sets NIBP measurement data

**NOTE:** Regarding the AG specifications, refer to the supplement Scio Four modules.

#### Recorder

Record width

48 mm (1.9 inch)

Paper width

50 mm

Paper speed

12.5, 25, 50 mm/s

Trace

Up to 3 waveforms

Recording types

- Continuous real-time recording
- 8 seconds real-time recording
- Time recording
- Alarm recording
- Trend graph recording
- Trend table recording
- NIBP review recording
- Arrhythmia review recording
- Alarm review recording

## Technical Data

- C.O. measurement recording
- Frozen waveform recording
- Drug calculation titration recording
- Hemodynamic calculation result recording

### Display Specifications

Display screen	305 mm (12.1 inch) color touch screen TFT
Resolution	800 x 600
Maximum number of waveforms	11
Indicator LEDs	1 power, 2 alarm, 1 charge

### Physical Specification

Size (H x W x D)	266 x 344 x 145 mm
Weight	5 kg

### Electrical Specification

Power supply	100 V – 240 V~, 50 Hz/60 Hz <b>5.3.1</b>
Pmax	110 VA
FUSE	T 3.15 AH, 250 V

### Classification

Anti-electroshock type/protection class	Class I equipment and internal powered equipment
EMC type	Class A
Anti-electroshock degree	CF: ECG (RESP), TEMP, IBP, CO BF: SpO <sub>2</sub> , NIBP, CO <sub>2</sub> , AG
Liquid ingress protection	IPX1
Disinfection/sterilization method	Refer to Instructions for Use: Care and cleaning
Mode of operation	Continuous running equipment
Power supply	100 V to 240 V~, 50 Hz/60 Hz Pmax = 110 VA FUSE T 3.15 AH, 250 VP

### Lithium-ion Battery (optional) **5.3.2**

Quantity	1
Capacity	5,000 mAh
Battery life	>= 350 min (At 25±2°C, with (a) new fully charged battery/ batteries, continuous SpO <sub>2</sub> measurement and NIBP automatic measurement mode at interval of 15 minutes, Dräger ECG/TEMP module connected, recording at interval of 10 minutes, brightness set to "1")
Battery charge time	≤ 390 min, 100% charge ≤ 351 min, 90% charge (monitor is off)

### Environmental Requirements

The monitor may not meet the performance specifications given here if stored or used outside the specified temperature, humidity and altitude ranges.

#### Temperature Range

Operating	0 to 40°C (32 to 104°F)
Transport and storage	-20 to 55°C (-4 to 131°F)

#### Relative Humidity

Operating	15% RH ~ 95% RH (non-condensing)
-----------	----------------------------------

## Technical Data

Transport and storage	15% RH ~ 95% RH (non-condensing)
-----------------------	----------------------------------

### Atmospheric Pressure

Operating	86 kPa ~ 106 kPa
Transport and storage	70 kPa ~ 106 kPa

### Standards

IEC 60601-1: 2005+A1 :2012; IEC 60601-1-2: 2007; EN 60601-1: 2006+A1 :2013; EN 60601-1-2: 2007; IEC 60601-2-49: 2011

The Vista 120 S monitors comply with the Medical Device Directive (MDD) 93/42/EEC.

Vista 120 S	Model A MS32996	Model A+ MS32998	Model C MS32997	Model C+ MS32999
3/5 lead ECG	X	X	X	X
Proprietary SpO <sub>2</sub>	X		X	
Nellcor SpO <sub>2</sub>		X		X
NBP	X	X	X	X
Dual temps	X	X	X	X
3IBP			X	X
CO			X	X
etCO <sub>2</sub>			X	X
Built-in recorder		X	X	X
Gas bench	X	X	X	X
LAN	X	X	X	X
Wireless		X	X	X

Vista 120 S patient monitors are available in select markets only.

Please contact your local sales office for more information.

## Notes

Not all products, features, or services are for sale in all countries.  
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to [www.draeger.com/trademarks](http://www.draeger.com/trademarks) to find the current status.

### **CORPORATE HEADQUARTERS**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
[www.draeger.com](http://www.draeger.com)

### **Manufacturer:**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23542 Lübeck, Germany

### **REGION EUROPE**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### **REGION MIDDLE EAST, AFRICA**

Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

### **REGION ASIA PACIFIC**

Draeger Singapore Pte. Ltd.  
61 Science Park Road  
The Galen #04-01  
Singapore 117525  
Tel: +65 6872 9288  
Fax: +65 6259 0398  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

### **REGION CENTRAL AND SOUTH AMERICA**

Dräger Indústria e Comércio Ltda.  
Al. Pucurui - 51 - Tamboré  
06460-100 - Barueri - São Paulo  
Tel. +55 (11) 4689-4900  
[relacionamento@draeger.com](mailto:relacionamento@draeger.com)  
Draeger  
Panamá Comercial

Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

