# **AQT90 FLEX**

**TnT CAL Cartridge**

**RADIOMETER** 

# Content

en	TnT CAL Cartridge .....	4
bg	Калибрираща касета за TnT .....	7
cs	Kazeta TnT CAL .....	10
da	TnT KAL-Kassette .....	13
de	TnT KAL-Kassette .....	16
el	Κασέτα βαθμονόμησης TnT .....	19
es	Cartucho CAL de TnT .....	22
et	TnT CAL-kassett .....	25
fi	TnT CAL -kasetti .....	28
fr	Cartouche CAL TnT .....	31
hr	TnT kalibracijska kaseta .....	34
hu	TnT CAL kazetta .....	37
it	Cartuccia CAL TnT .....	40
ja	TnT CALカートリッジ .....	43
lt	TnT CAL kasetė .....	46
lv	TnT CAL kārtidžs .....	49
nl	TnT CAL-cassette .....	52
no	TnT KAL-kassett .....	55
pl	Kaseta TnT CAL .....	58
pt	Cartucho CAL TnT .....	61
ro	Cartușul CAL pentru TnT .....	64
ru	CAL-картридж для TnT .....	67
sk	Kalibračná kazeta CAL pre TnT .....	70
sv	CAL-kassett för TnT .....	73
tr	TnT Kalibrasyon kartuşu .....	76
zh	TnT 定标测试卡 .....	79

en

### **TnT CAL Cartridge**

Proprietary name and common name: Cardiac troponin T (TnT)

Catalog no. 944-268

### **Intended use**

For *in vitro* diagnostic use.

The TnT **CAL** Cartridges are *in vitro* diagnostic reagents intended for the calibration adjustment of the TnT Test on the AQT90 FLEX analyzer by establishing points of reference to estimate TnT values.

**NOTICE:** The cartridges are ready for use.

### **Summary**

Calibration adjustment must be performed every time the Test Cartridge lot changes. The assay has been factory-calibrated using an 9-point calibration and the factory-defined calibration data is provided in the form of a barcode. This data is used as a lot-specific reference curve, which is adjusted to the correct instrument-specific signal level by using a CAL cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and instrument-specific calibration curve for the assay. When a CAL cartridge with a new lot number is loaded, the analyzer requests the factory-defined calibration data so that it can perform the calibration adjustment. If the user does not load a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after performing a calibration adjustment with a CAL Cartridge with the same lot number.

**NOTICE:** Liquid quality controls must be run after every calibration adjustment.

See the *AQT90 FLEX operator's manual* for details.

### **Contents:**

- One TnT CAL Cartridge for performing one calibration adjustment
- One "TnT Factory-defined Calibration Data" sheet with barcode

### **Components:**

The CAL Cartridge contains eight analyte-specific background cups and eight cups with added antigen.

### **Storage:**

Store unopened at 2-8 °C

**Shelf life:** If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the Expiration date shown on the pouch label.

**On-board stability:** See label on the CAL Cartridge pouch.

**Chemical composition:**

Reactive ingredients:

- Mouse monoclonal TnT capturing and tracing antibodies: 850 ng/cup
- Recombinant TnT (only in the eight CAL cups with added antigen): approx. 0.12 ng/cup

Other ingredients:

- Bovine serum albumin
- Bovine  $\gamma$ -globulin
- Mouse IgG (blocker substances for heterophilic interferences)

Main buffer component: TRIS (Tris Hydroxymethylaminomethane)

Toxic ingredient: Sodium azide < 0.1 %

**Data in the cartridge barcodes:**

• Checksum	• Cartridge ID
• Barcode type	• Expiration date
• Parameter code	• Barcode version
• Lot no.	• On-board stability

**NOTICE:** The barcode is lot-specific.

**CAUTION – Sensitive to humidity**

Do not use any cartridge that has not been kept according to storage specifications or that has been damaged or taken out of its protective pouch.

**Procedures**

**NOTICE:** To keep the analyzer ready for measurement at all times, the power of the analyzer must always be on.

**Calibration adjustment procedure**

**NOTICE:** Calibration adjustment must be performed every time a Test Cartridge with a new lot number is used.

1. Find the CAL Cartridge in the CAL Cartridge box with the new lot number.
2. Load it into the cartridge inlet with the barcode facing downwards.  
**NOTICE:** The analyzer does not accept cartridges after the expiration date printed on its pouch label.

**WARNING – Risk of incorrect results**

Do not use a cartridge that has already been used on another analyzer, otherwise you risk obtaining incorrect results.

3. Scan the barcode from the "TnT Factory-defined Calibration Data" sheet provided in the CAL Cartridge box.
4. Start the calibration adjustment procedure.  
Results are marked on the screen as successful or failed.

**Patents**

Manufactured under license from Roche Diagnostics. US Patent No. US6,376,206.

In addition, one or more of the following patents and patent applications may apply:

- US Patent No.: US6,429,026
- European Patent No.: EP892927
- Japanese Patent No.: JP3537829
- German Patent No.: DE69729186
- US Patent Application No.: US2008-0226508
- European Patent Application No.: EP2130050
- International Patent Application No.: WO2008/113352

**Explanation of graphical symbols**

See the *Explanation of graphical symbols* table at the end of this insert.

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Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQURE, PICO, CLINITUBES and QUALICHECK are trademarks of Radiometer Medical ApS.
All rights reserved. 995-640. 201112F.