EasyQC CE

MEDICA

Blood Gas and Electrolyte Quality Control, Level 1

REF 6303

Intended Use

To monitor the performance of the EasyStat and EasyBloodGas analyzers. Refer to Analyze QC in the Operator's Manuals for detailed instructions. For in vitro diagnostic use.

Active Ingredients

EasyQC controls contain buffered bicarbonate and electrolyte solution equilibrated with precise mixtures of O₂, CO₂, and N₂ gases. This control contains no preservatives or human-based materials.

Product Description

EasyQC is available for monitoring analyzer performance at different points in the clinical range. It is packaged in sealed glass ampules containing 1.7 mL of solution. Ampules are packaged thirty (30) per box.

Storage and Stability

Store control materials at 4–25°C. Changes in values may occur if stored at temperatures greater than 30°C for prolonged periods. DO NOT FREEZE. Discard any control material which exceeds the expiration date on the product label.

Recommended Use

Run quality control material once per shift before analyzing patient specimens. Run the controls at any time to verify analyzer performance. Use control material immediately after opening, following steps 1–5 in the instructions section of this sheet.

Expected Values:

Assigned values for each level of a specified lot of control material are based on the results of multiple tests run at Medica's laboratories. Approximately 95% of all EasyStat and EasyBloodGas analyzers, operating according to specifications, are expected to fall within the ranges listed on this sheet.

Instructions:

- 1 Controls must be brought to room temperature (23°C +/- 1) prior to use. If the control material has not been stored at room temperature, remove the ampules from the box and allow the material to equilibrate to room temperature for four hours. Refer to limitations 3 and 4.
- 2 Before use, hold the ampules at the tip and bottom (to minimize increasing the temperature) and shake 15–20 times (10 seconds) to mix the gas and liquid phase. Tap the ampule to restore liquid to the bottom.
- 3 Carefully snap off the top of the ampule at the score. Avoid injury by protecting fingers with tissue or gloves.
- 4 Immediately introduce control material to the analyzer using direct aspiration. Run control material directly from the ampule within one minute after opening.
- 5 Repeat steps 1–4 for remaining control levels.

Limitations:

- 1 This control is for use as an aid in the evaluation of performance of the EasyStat and EasyBloodGas analyzers, and should not be considered as a substitute for other aspects of total quality control, such as calibration, maintenance and record keeping.
- 2 This control does not contain red blood cells and, therefore, may not detect certain analyzer malfunctions which would affect the testing of whole blood.
- 3 The gas/liquid equilibrium in each ampule is dependent upon temperature. A controlled temperature of 23°C will produce the most accurate measurements of pH, PCO₂ and PO₂. (NOTE: PO₂ will vary inversely by about 1% per degree C that the temperature of the ampule varies from 23°C.)
- 4 In establishing the values shown in Expected Values, samples are analyzed at approximately 760 mmHg. The PO₂ value will decrease by approximately 1.0% per 1000 feet (305 meters) above sea level, and PCO₂ values decrease by approximately 0.5% under the same conditions.

EasyQC CE

MEDICA

Blood Gas and Electrolyte Quality Control, Level 2

REF 6304

Intended Use

To monitor the performance of the EasyStat and EasyBloodGas analyzers. Refer to Analyze QC in the Operator's Manuals for detailed instructions. For in vitro diagnostic use.

Active Ingredients

EasyQC controls contain buffered bicarbonate and electrolyte solution equilibrated with precise mixtures of O₂, CO₂, and N₂ gases. This control contains no preservatives or human-based materials.

Product Description

EasyQC is available for monitoring analyzer performance at different points in the clinical range. It is packaged in sealed glass ampules containing 1.7 mL of solution. Ampules are packaged thirty (30) per box.

Storage and Stability

Store control materials at 4–25°C. Changes in values may occur if stored at temperatures greater than 30°C for prolonged periods. DO NOT FREEZE. Discard any control material which exceeds the expiration date on the product label.

Recommended Use

Run quality control material once per shift before analyzing patient specimens. Run the controls at any time to verify analyzer performance. Use control material immediately after opening, following steps 1–5 in the instructions section of this sheet.

Expected Values

Assigned values for each level of a specified lot of control material are based on the results of multiple tests run at Medica's laboratories. Approximately 95% of all EasyStat and EasyBloodGas analyzers, operating according to specifications, are expected to fall within the ranges listed on this sheet.

Instructions

- 1 Controls must be brought to room temperature (23°C +/- 1) prior to use. If the control material has not been stored at room temperature, remove the ampules from the box and allow the material to equilibrate to room temperature for four hours. Refer to limitations 3 and 4.
- 2 Before use, hold the ampules at the tip and bottom (to minimize increasing the temperature) and shake 15–20 times (10 seconds) to mix the gas and liquid phase. Tap the ampule to restore liquid to the bottom.
- 3 Carefully snap off the top of the ampule at the score. Avoid injury by protecting fingers with tissue or gloves.
- 4 Immediately introduce control material to the analyzer using direct aspiration. Run control material directly from the ampule within one minute after opening.
- 5 Repeat steps 1–4 for remaining control levels.

Limitations

- 1 This control is for use as an aid in the evaluation of performance of the EasyStat and EasyBloodGas analyzers, and should not be considered as a substitute for other aspects of total quality control, such as calibration, maintenance and record keeping.
- 2 This control does not contain red blood cells and, therefore, may not detect certain analyzer malfunctions which would affect the testing of whole blood.
- 3 The gas/liquid equilibrium in each ampule is dependent upon temperature. A controlled temperature of 23°C will produce the most accurate measurements of pH, PCO₂ and PO₂. (NOTE: PO₂ will vary inversely by about 1% per degree C that the temperature of the ampule varies from 23°C.)
- 4 In establishing the values shown in Expected Values, samples are analyzed at approximately 760 mmHg. The PO₂ value will decrease by approximately 1.0% per 1000 feet (305 meters) above sea level, and PCO₂ values decrease by approximately 0.5% under the same conditions.

EasyQC CE

MEDICA

Blood Gas and Electrolyte Quality Control, Level 3

REF 6305

Intended Use

To monitor the performance of the EasyStat and EasyBloodGas analyzers. Refer to Analyze QC in the Operator's Manuals for detailed instructions. For in vitro diagnostic use.

Active Ingredients

EasyQC controls contain buffered bicarbonate and electrolyte solution equilibrated with precise mixtures of O₂, CO₂, and N₂ gases. This control contains no preservatives or human-based materials.

Product Description

EasyQC is available for monitoring analyzer performance at different points in the clinical range. It is packaged in sealed glass ampules containing 1.7 mL of solution. Ampules are packaged thirty (30) per box.

Storage and Stability

Store control materials at 4–25°C. Changes in values may occur if stored at temperatures greater than 30°C for prolonged periods. DO NOT FREEZE. Discard any control material which exceeds the expiration date on the product label.

Recommended Use

Run quality control material once per shift before analyzing patient specimens. Run the controls at any time to verify analyzer performance. Use control material immediately after opening, following steps 1–5 in the instructions section of this sheet.

Expected Values:

Assigned values for each level of a specified lot of control material are based on the results of multiple tests run at Medica's laboratories. Approximately 95% of all EasyStat and EasyBloodGas analyzers, operating according to specifications, are expected to fall within the ranges listed on this sheet.

Instructions:

- 1 Controls must be brought to room temperature (23°C +/- 1) prior to use. If the control material has not been stored at room temperature, remove the ampules from the box and allow the material to equilibrate to room temperature for four hours. Refer to limitations 3 and 4.
- 2 Before use, hold the ampules at the tip and bottom (to minimize increasing the temperature) and shake 15–20 times (10 seconds) to mix the gas and liquid phase. Tap the ampule to restore liquid to the bottom.
- 3 Carefully snap off the top of the ampule at the score. Avoid injury by protecting fingers with tissue or gloves.
- 4 Immediately introduce control material to the analyzer using direct aspiration. Run control material directly from the ampule within one minute after opening.
- 5 Repeat steps 1–4 for remaining control levels.

Limitations:

- 1 This control is for use as an aid in the evaluation of performance of the EasyStat and EasyBloodGas analyzers, and should not be considered as a substitute for other aspects of total quality control, such as calibration, maintenance and record keeping.
- 2 This control does not contain red blood cells and, therefore, may not detect certain analyzer malfunctions which would affect the testing of whole blood.
- 3 The gas/liquid equilibrium in each ampule is dependent upon temperature. A controlled temperature of 23°C will produce the most accurate measurements of pH, PCO₂ and PO₂. (NOTE: PO₂ will vary inversely by about 1% per degree C that the temperature of the ampule varies from 23°C.)
- 4 In establishing the values shown in Expected Values, samples are analyzed at approximately 760 mmHg. The PO₂ value will decrease by approximately 1.0% per 1000 feet (305 meters) above sea level, and PCO₂ values decrease by approximately 0.5% under the same conditions.

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EasyStat Material Safety Data Sheets


Attention: Safety Officer / MSDS Requester

Following is a list of chemicals/solutions required to operate and maintain an EasyStat analyzer. This list identifies which solutions require and which do not require a Material safety Data Sheet (MSDS).

According to OSHA regulations, substances containing more than 1% of a hazardous material or more than 0.1% of an extremely hazardous material, require the availability of an MSDS.

Products in this list which do not require an MSDS, should be handled with care and avoid contact with skin. If contact occurs, wash exposed area with plenty of water. As a reference, refer to the labels for contents of active ingredients.

<u>Product Description</u>	<u>MSDS</u>	<u>Substance, (content)</u>
Reagent Module	No	Non-hazardous
Daily Cleaning Solution Kit	Yes	Pepsin, (powder)
Quality Control Kits (all levels)	No	Saline solution & buffer
Red Test Dye	No	Non-hazardous

4410 SOUTH 102ND STREET

OMAHA, NEBRASKA 68127

PHONE: 402-339-2494

FAX: 402-339-0801

EMAIL: ali@americanlaboratories.com

NOTE: SEE LAST PAGE FOR CONDITIONS UNDER WHICH THE INFORMATION AND RECOMMENDATIONS OF THIS DATA IS FURNISHED.

MATERIAL SAFETY DATA

DATE: December 1, 1996

AMERICAN LABORATORIES, INC.
4410 SOUTH 102nd STREET
OMAHA, NE 68127
(402-339-2494)

NA- Not Applicable

Page 1 of 4

SECTION I - IDENTIFICATION

Trade Name: K-ZYME AP LIQUID

Name and/or Family or Description: AQUEOUS ENZYME MIXTURE Formula: NA
(No IUC name for the mixture is known to exist)

Chemical Name: Proteolytic Enzyme DOT HAZARD CLASS: NA

CAS No. : See information listed below in Section II DOT LABEL(S): NA

DOT PROPER SHIPPING NAME: Chemicals Not Otherwise Indexed (NOI), Non Hazardous

SECTION II - (Hazardous & other)

INGREDIENT	CAS-NO.	PERCENT	OCCUPATIONAL EXPOSURE LIMITS TLV	VAPOR PRESSURE
------------	---------	---------	----------------------------------------	-------------------

BACTERIAL SUBTILISIN PROTEASE	9014-01-1	NA	$0.10 \times 10^{-6} \text{ AU/m}^3$	As water
-------------------------------------	-----------	----	--------------------------------------	----------

AGGIH = TLV (as pure crystalline enzyme) 0.06 ug/cu meter or as above

These compounds are not considered carcinogens by the National Toxicology Program, the International Agency for research on Cancer, or the Occupational Safety and Health Administration nor are they known to be mutagenic.

Identity	OTHER INGREDIENTS	CAS NO.	WT%
Water		(7732-18-5	@
Propylene Glycol		(57-55-6)	@

@ - among 5 most predominant compounds NA- not applicable

SECTION III - PHYSICAL DATA

STATE: Liquid APPEARANCE: Brown to amber ODOR: characteristic
BOILING PT: 100-105C MELTING PT: NA EVAPORATION RATE: as water
STORAGE CONDITIONS: Store tightly closed in cool, dry place
SOLUBILITY IN WATER: NA (Aqueous Solution) 100% Miscible

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

FLASH POINT * IGNITION TEMPERATURE *

FLAMMABLE LIMITS: NA(Aqueous solution) * Information not available

No special requirements, however, self-contained breathing apparatus and protective clothing should be worn in fighting fires involving chemicals.

Extinguishing Media: No special requirements.

SECTION V - HEALTH HAZARD DATA

EFFECTS OF OVEREXPOSURE: Possible irritation or allergic response in affected individuals to skin contact, eye contact, or inhalation of dust from product.

CHRONIC EFFECTS OF OVEREXPOSURE: May cause allergic type sensitivity in susceptible individuals resulting in delayed hyperallergenic response.

MEDICAL CONDITIONS PRONE TO AGGRAVATION BY EXPOSURE: Individuals with known allergies, particularly to enzymes, as well as those individuals with chronic respiratory disease or particular susceptibility thereto and those with asthmatic conditions should avoid contact with dusty conditions.

PRIMARY ROUTES OF ENTRY: Ingestion and Inhalation, no appreciable entry known through intact skin although irritation possible, see above.

EMERGENCY AND FIRST AID PROCEDURES

Rinse mouth/throat with water if ingested. Wash affected area of skin with soap and remove contaminated clothing. In case of inhalation, remove from area of exposure; if allergic symptoms occur in 0-48 hours, seek medical attention. If eye contact occurs rinse thoroughly with water for several minutes, seek medical attention if irritation occurs. Supplementary information available upon request concerning concentrated enzyme exposure.

WARNING STATEMENTS

CAUTION! MAY CAUSE SKIN OR EYE IRRITATION. STORE IN COOL DRY PLACE. AVOID CONTACT WITH SKIN AND CLOTHING AND INHALATION OF AIRBORNE PARTICLES. WASH HANDS, CLOTHING, AND EQUIPMENT AFTER USE.

SECTION VI - REACTIVITY DATA

MATERIALS TO AVOID: None known

HAZARDOUS DECOMPOSITION PRODUCTS: None Known

SECTION VII - SPILL, LEAK, AND DISPOSAL INFORMATION

WASTE DISPOSAL - Dispose of this material and any adsorber material (see below) in closed containers in a sanitary chemical landfill which complies with all local, state, and federal regulations.

SPILL OR LEAKAGE PROCEDURES - Dike and adsorb spill with inert materials - small quantities can be flushed thoroughly with water to appropriate chemical waste sewer (avoid high pressure washing).

SECTION VIII - SAFE HANDLING AND USE INFORMATION

Respiratory Protection: Wear suitable NIOSH-approved dust masks or respiratory protector where local ventilation exhaust is not available or situations where inhalation of mist is likely to result.

Ventilation: Local ventilation required to prevent ambient air entrainment of fugitive aerosol material.

OTHER PROTECTIVE EQUIPMENT: Eye wash facility, whatever equipment sufficient to prevent skin contact and splash protection.

Storage and Handling: Maintain good housekeeping, avoid creating aerosols. Keep out of direct sunlight and in cool dry place as deterioration of product and loss of activity may occur under above conditions.

Eye Protection: Liquid tight chemical goggles

Skin Protection: Long sleeved shirt, trousers, safety shoes, and rubber/polyethylene gloves.

Hygienic Practices: Wash hands after handling, as well as any other affected skin areas. Avoid contact with food or food preparatory surfaces. IF this occurs wash the area thoroughly with suitable detergent and water. Remove and wash contaminated clothing.

SECTION IX - SPECIAL PRECAUTIONS

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORING - Maintain good housekeeping, avoid creating aerosols. For maintenance of product quality, store in dry, cool (under 110 deg. F.) warehouse conditions.

KEEP OUT OF REACH OF CHILDREN

ALTHOUGH THE INFORMATION AND RECOMMENDATIONS SET FORTH IN THIS SHEET ARE BELIEVED TO BE CORRECT AS OF THE DATE HEREOF, AMERICAN LABORATORIES, INC. MAKES NO REPRESENTATION AS TO THE COMPLETENESS OR ACCURACY OF SUCH INFORMATION AND RECOMMENDATIONS FOR SPECIFIC APPLICATIONS. AMERICAN LABORATORIES, INC. SHALL IN NO EVENT BE RESPONSIBLE FOR ANY DAMAGES OF WHATSOEVER NATURE DIRECTLY OR INDIRECTLY RESULTING FROM THE PUBLICATION OR USE OF OR SOLE RELIANCE UPON SUCH INFORMATION AND RECOMMENDATIONS. YOU ARE ENCOURAGED TO ADVISE ANYONE WORKING WITH OR EXPOSED TO SUCH PRODUCTS OF THE INFORMATION CONTAINED HEREIN.

NO WARRANTY, EITHER EXPRESS OR IMPLIED, OF MERCHANTABILITY OR FITNESS OR OF ANY OTHER NATURE WITH RESPECT TO THE PRODUCT OR TO THE INFORMATION AND RECOMMENDATIONS HEREIN IS MADE HEREUNDER.