ASO (StrepA) test kit

for quantitative in vitro determination of Anti Streptolysin O in Finger blood, EDTA venous blood or serum on smart or CUBE laboratory photometer.

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Kit size

16 tests/kit
2 x 1ml decision range

Test kit preparation: Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The ASO test is used to determine post streptococcal complications including rheumatic fever and glomerulonephritis as well as recent streptococcal infection.

The presence and level of ASO antibodies in human blood/serum directly reflects the extent and degree of infection. Elevated levels of ASO may also be present in other conditions including tonsillitis, acute rheumatoid arthritis, scarlet fever, and various other streptococcal infections.

Method / Measurement Range

Latex Enhanced Immunturbidimetric assay:
100-2000 IU/ml

Principle

Latex particles coated with Streptolysin O react with ASO in the patient sample resulting in an increase in turbidity, the absorbance is measured at 700nm and is directly proportional to the concentration of ASO.

Test Kit

ERS Cuvette pre-filled with buffer reagent
ERS Cap pre-filled with latex reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

Warnings and Precautions

The reagents contain sodium azide (0.95g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Freshly drawn whole blood from finger tip is preferred. Alternatively EDTA venous blood or serum can be used.

Reference Range

≤100 IU/ml for preschool children
≤ 250 IU/ml for school age children
≤ 200 IU/ml for adults

It is recommended that each laboratory should establish its own expected range.

Quality Control

For internal quality control the ASO control kit should be used.

Order number: ST2500

Precision

following ...

References

1. Klein G.C., Baker C.N., Jones W.L., Upper Limits of normal Antistreptolysin-O and Antideoxyribonuclease B titres. Applied Microbiology. 1971; 21: 999-1001
2. Thomas L. ed., Clinical Laboratory Diagnosis. Streptococcus pyogenes Infection. In: Thomas L.ed., Clinical Laboratory Diagnosis. Use and assessment of Clinical Laboratory results. 1st Edition (1998). TH-Books, Frankfurt/Main Germany.

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ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (> 20 °C)!
Use the integrated capillary ONLY for aspirating whole blood from finger tip! Add EDTA venous blood or serum ONLY with pipette!

Correct use of the integrated capillary

- A CORRECT: Capillary is completely filled with blood
- B WRONG: Capillary is partially filled with blood
- C WRONG: Too much blood around the capillary and/or ERS cap

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

2. Select whole blood or serum


Please consult the laboratory photometer user manual for more information.

- 2.1 Sampletype menu: Select the sample type
- 2.2 Hematocrit correction (if desired)
- 2.3 Instrument ready for start of analysis

3. Sample preparation with whole blood from finger tip


- 3.1a Aspirate 5µl whole blood with integrated capillary from finger tip
- 3.2a Apply ERS cap firmly onto ERS cuvette
- 3.3a Shake ERS cartridge 3 - 4x gently until the finger blood is fully dispensed out of capillary into ERS liquid

or with EDTA venous blood or serum


- 3.1b Aspirate 5µl EDTA venous blood or serum with pipette
- 3.2b Dispense 5µl EDTA venous blood or serum INTO LIQUID in ERS cuvette
- 3.3b Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the  start button on the smart laboratory photometer or by closing the door of the CUBE laboratory photometer.

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
Performing an ASO (anti-Streptolysin-O) test

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



CUBE smart

1. 1.1 → 1.2 → 1.3 → 1.4

1.1: Inserting the cartridge into the device.

1.2: The device is shown with the cartridge inserted.

1.3: The device is shown with the cartridge inserted.

1.4: The device is shown with the cartridge inserted.

CUBE smart

2. 2.1 → 2.2 → 2.3

2.1: The device screen shows patient information: Name: 29, ID: 29, Sex: Male, Sample type: Blood, Hematocrit: 40. There are checkmarks for 'Serum' and 'Blood'.

2.2: The device screen shows a keypad with numbers 1-9, 0, and a 'Done' button.

2.3: The device screen shows the instruction 'Insert cartridge, close door' and a red icon of the device.

CUBE smart

3. 3.1a → 3.2a → 3.3a → 3.4

3.1a: The device is shown with the cartridge inserted.

3.2a: The device is shown with the cartridge inserted.

3.3a: The device is shown with the cartridge inserted.

3.4: The device is shown with the cartridge inserted.