# CRP test kit

## with integrated capillary

for quantitative in vitro determination of c-reactive protein in fingerblood, venous blood or serum on the smart or CUBE laboratory photometer.

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### English

**Order Information**  
Order number: ST0102  
Order number: ST1000

**Indication**  
CRP test kit with integrated capillary  
CRP control kit

**Kit size**  
32 tests/pack  
2 x 1 ml (low high range)

**Test kit preparation:** Allow single test min. 10 minutes to warm up to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

### Summary

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. It increases already after 6 hours in acute inflammatory processes associated with bacterial infections, post operative conditions or tissue damage. Many studies have demonstrated that in apparently healthy subjects there is a direct correlation between CRP concentrations and the risk for developing coronary heart disease (CHD).

### Method / Measurement Range

Immunturbidimetric Assay:

Serum: 0,5 - 120mg/l at 546nm  
1,0 - 120mg/l at 700nm  
Whole blood: 2,0 - 240mg/l at 546 & 700nm

### Principle

Kinetic determination of the concentration of CRP by photometric measurement at 546nm or 700nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

### Test Kit

ERS Cuvette pre-filled with buffer reagent  
ERS Cap pre-filled with latex reagent with rabbit anti human CRP antibodies

### Storage Instructions and Stability of the Test Kit

The reagents are stable until the indicated expiry date, if stored at 2-8°C.

### Warnings and Precautions

The reagents contain sodium azide (0.95g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

### Waste Management

Please refer to local legal requirements.

### Specimen Sample Material

Freshly drawn fingertip blood is preferred.  
Alternatively EDTA venous blood or serum can be used.

### Reference Range

Values <10 mg/l are expected as normal. Any condition that results in sudden or severe inflammation may increase CRP level.  
It is recommended that each laboratory should establish its own expected range.

### Quality Control

For internal quality control the CRP control kit should be used.  
Order number: ST1000

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## English

### ATTENTION!

Use the integrated capillary ONLY for aspirating whole blood from finger tip! Add EDTA venous blood or serum ONLY with pipette! Do NOT wipe the capillary respectively pipette tip after aspirating the sample! Allow single test min. 10 minutes to warm up to room temperature before use!

### Correct use of the integrated capillary

- A CORRECT: Capillary is completely filled with blood
- B WRONG: Capillary is partially filled with blood
- C WRONG: Too much blood around the capillary and/or ERS cap

### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

### 2. Select whole blood or serum

- Please consult the laboratory photometer user manual for more information.
- 2.1 Sampletype menu: Select the sample type
- 2.2 Hematocrit correction (if desired)
- 2.3 Instrument ready for start of analysis

### 3. Sample preparation

#### with whole blood from finger tip

- 3.1a Aspirate 5µl whole blood with integrated capillary from finger tip
- 3.2a Apply ERS cap firmly onto ERS cuvette
- 3.3a Shake ERS cartridge 3 - 4x gently until the finger blood is fully dispensed out of capillary into ERS liquid

#### or with EDTA venous blood or serum

- 3.1b Aspirate 5µl EDTA venous blood or serum with pipette
- 3.2b Dispense 5µl EDTA venous blood or serum INTO LIQUID in ERS cuvette
- 3.3b Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

Technical details subject to change without notice.

Doc ID: 2adba42-4198-4309-91b0-efba5d30ed1a ; latest revision DP 5.0 ; date of print: 2013-12-02

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