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EN Instructions for use

Patient tubing for ulrich medical CT/MRI injectors

About this document

- These instructions for use are part of the product and they describe its safe use as intended.
- Read the instructions for use before using the product.
- Keep the instructions for use with the product so that they are accessible at all times.

Highlighting and symbols

Meaning

- ⚠ WARNING** Hazards to persons. Non-compliance can cause death or serious injuries to persons and damage to the product.
- Handling instructions.

Symbols on the product and packaging

Icon	Description	Icon	Description
	Caution		Sterilized using ethylene oxide
	Manufacturer's CE mark		Do not use if package is damaged
	Product number		Do not use tubing if the protective caps are not properly attached
	Lot code		Pyrogen-free
	Use by date		Keep dry
	Manufacturer		Temperature limit
	Consult the instructions for use		Relative humidity limit
	Do not re-use		Protect from sunlight

Intended use

The patient tubing for ulrich medical CT/MRI injectors is used for venous-side administration of contrast agents and physiological saline solutions (NAC) during a CT/MRI examination. The patient tubing creates the connection from the ulrich medical pump tubing (device side) to the patient.

Product description

The patient tubing consists of tubing with two protective caps and has a male Luer lock connector with integrated check valve on the patient side. The pump tubing is connected to the female Luer lock connector. A second check valve is integrated in the tubing. The check valves are self-closing. They only open when a medium flows in the defined direction. If back pressure is exerted or if there is no media flow, the valves remain closed.

- | | |
|---|---------------------------------------|
| A | Device side with pump tubing |
| B | Patient side |
| 1 | Female Luer lock |
| 2 | Protective caps for tubing ends |
| 3 | Male Luer lock with first check valve |
| 4 | Second check valve |
| 5 | Media flow direction |

Delivery contents and scope of application

Art. no.	Designation of the patient tubing
XD 2030	Patient hose 1,20 m
XD 2035	Patient hose 1,50 m
XD 2040	Patient hose 2,50 m
XD 2045	Patient hose 3,20 m

Sterili, be lateko, su
dvim apsauginiais
vožtuvais, 250 cm.
ilgio.

Product for single use

- The product is EO sterilized and delivered in sterile packaging.
- Do not reuse product.
- Do not resterilize product.
- Do not use product from open or damaged sterile packaging.
- Discontinue using product if the use-by date has expired.
- Do not use product that was not stored under the recommended storage conditions.

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- Before using the product, check its operability and proper condition.
- To avoid damage due to improper use, use product according to these instructions for use and the instructions for use of the injector.
- Product and accessories should only be used by persons who have the required knowledge, training or experience.

Risk of injury!

The use of tubing combinations or systems that have not been approved by ulrich medical jeopardizes the safety of the patient and/or the user as well as the proper functioning of the injector.

- Only use original products approved by ulrich medical.

Risk of injury! Air injection, overdose and flow reduction due to pressure regulation if tubing extensions are used.

Tubing extensions jeopardize the safety of the patient and risk contaminating the tubing system.

- Do not use tubing extensions (e.g. "Heidelberger extension").

Risk of injury! Air injection and spraying of contrast agent. Risk of damage to property.

If screwed in too tightly, the thread of the Luer lock connector may break and lead to open cracks.

- Close the Luer lock connector properly.

Risk of injury due to infections and contaminations!

Improper handling jeopardizes the safety of the patient.

- Handle all sterile components carefully in order to ensure sterility.

→ Do not use product if the protective caps are missing or are not attached securely.

→ Replace patient tubing for each new patient.

→ If new patient tubing is not connected immediately after disconnecting the patient tubing from the pump tubing: Place a sterile protective cap on the pump tubing.

- Dispose of patient tubing upon disconnecting it from the patient in order to prevent contamination of the pump tubing.

Use

Preparation

- Remove the product from the packaging.
- Check product before use for damaged parts.
- Do not use product with damaged parts.

Connect product to pump tubing

- Connect product only as directed in the instructions for use of the injector.

Performance testing

- Conduct performance testing only as directed in the instructions for use of the injector.

Operation

→ Operate product only as directed in the instructions for use of the injector.

Validated processing method

Single-use product

WARNUNG! Patient-/user infection due to re-use and/or processing.

If re-used and/or processed, there is a risk of contamination and material fatigue. Soiling, contamination, impaired function may lead to injury, illness or death.
→ Do not process product.

Technical service

ulrich GmbH & Co. KG
Injector Systems Customer Service
Buchbrunnenweg 12
89081 Ulm / Germany
Telephone: +49 (0)731 9654-0
Fax: +49 (0)731 9654-2805
Email: injectors@ulrichmedical.com

Technical Data

Parameter	Value
Pressure tightness	20 bar max.
Material	PVC (DEHP-free, latex free)
Max. permissible duty	12 hours for one patient (single-use item)
Length	Filling volume
1,2 m	6,0 ml
1,5 m	7,5 ml
2,5 m	12,5 ml
3,2 m	16,0 ml

Storage


- Store product in sterile packaging in a dry, dark and temperature-controlled room.
- Protect product in sterile packaging from environmental influences such as soil, dust, moisture, heat and sunlight.

Ambient conditions


	Use	Storage
Temperature	+5°C ... +40°C	+15°C ... +25°C
Relative humidity	10% ... 90%, non-condensing	30% ... 60%, non-condensing

Disposal

- Dispose of product after use according to local regulations.

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Technical specifications – tubing system for 3T tennessee™, mississippi™, missouri™, ohio M™ and ohio tandem™

Consumption	Length(s)	Specials
Pump tubing	120cm	2 check valves Proven CA compatibility Non-pyrogenic Latex-free DEHP-free
Patient tubing	150cm	3 accesses for reservoir containers, each with air filter Particle filter Integrated pressure control unit Proven resistance to pressure Proven CA compatibility Non-pyrogenic Latex-free DEHP-free
1 tubing per patient (max. 24 h)	250cm	
1 tubing per working day	320cm	

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Ulrich GmbH & Co. KG