Nutritub® ENFit®

GB	Gastral/Intestinal polyurethane tube, Instruction for use	10
DE	Gastrische/intestinale Sonde aus Polyurethan, Gebrauchsanweisung	14
BG	Гастрална/интестинална полиуретанова сонда, Инструкции за употреба	18
CN	聚氨酯胃/肠管, 使用说明	23
CZ	Žaludeční/střevní polyuretanová sonda, Návod k použití	26
EE	Poliüretaanist mao-/soolesond, Kasutusjuhised	30
ES	Sonda gástrica/intestinal de poliuretano, Instrucciones de Uso	34
FI	Polyuretaanista valmistettu maha- ja suoliletku, Käyttöohjeet	38
FR	Sonde gastrique/intestinale en polyuréthane, Directives pour l'utilisation	42
HR	Gastrički/intestinalni poliuretanski tubus, Upute za uporabu	47
HU	Poliuretán gyomor-/ bélsonda, Használati útmutató	51
ID	Slang Gastrik/Intestinal berbahan poliuretan, Petunjuk Penggunaan	55
IT	Sonda gastrica/intestinale in poliuretano, Indicazioni d'uso	60
KR	위장/장내 폴리우레탄 튜브, 사용 지침	64
LV	Gastrointestinālā poliuretāna zonde	68
NL	Maag-darmsonde van polyurethaan, Gebruiksaanwijzing	72
PL	Zgłębnik poliuretanowy do żywienia dożołądkowego/dojelitowego, Instrukcje użycia	77
PT	Sonda gástrica/intestinal de poliuretano, Instruções de uso	81
RO	Tub gastric/intestinal din poliuretan, Instrucțiuni de utilizare	85
RU	Желудочный/кишечный полиуретановый зонд, Инструкция по применению	90
SE	Gastriskt/intestinalt polyuretanrör, Bruksanvisning	94
SI	Gastrična/intestinalna sonda iz poliuretana, Navodila za uporabo	98
SK	Gastričná/intestinálna polyuretánová sonda, Návod na používanie	102
UA	Шлунково-кишковий поліуретановий зонд	106

B | BRAUN

grunewald 200990

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











L

GB General Information

Carefully read all the instructions before using Nutritub® ENFit®

Transnasal insertion. EtO Sterile.

Intended use:

Polyurethane feeding tube for nasogastral / nasointestinal enteral nutrition ONLY.

Transnasal/Oral insertion.

The user must be a trained medical professional according to local practices.

Indications:

- Enteral nutrition.
- Enteral medication.

Features:

Tubes are made of biocompatible polyurethane and PVC DEHP-free with radiopaque stripe and equipped with ENFit® Male connector according to ISO 80369-3.

A transparent polyurethane tube CH 5 (50 cm), CH 6 (60 cm), CH 8 (100 cm or 125 cm), CH 10 (110 cm) CH 12 (110 cm or 125 cm), CH 14 (110 cm), CH 16 (110 cm) or CH 18 (110 cm) with:

- Radio-opaque stripe
- Closed round shaped tip or rigid polyurethane white end cap
- Stainless steel stylet with ENFit®-lock lift (stylet is used for all Intestinal Nutritub ENFit® and in CH 10, CH 12, CH 14, CH 16, CH 18 Gastral).
- Length markers at 2 cm intervals for CH 5 and CH 6, and at 5 cm intervals for all the others CH.
- ENFit® capsule
- ENFit® connectors

Only for Intestinal Nutritub ENFit® CH 8 and CH 12 the tube used is 125 cm long.

Gastral Nutritub ENFit® without white PU end cap

have a series of holes on both side. In particular, in CH 5, CH 6 and CH 8 4 holes are present, while for CH 14, CH 16 and CH 18 3 holes can be found. Intestinal Nutritub ENFit® have a ballast containing tungsten attached to the end cap.

Nutritub ENFit® are equipped with ENFit® connectors and they must be connected ONLY to giving sets provided with ENFit® connectors.

NOTE: if the stylet is present, its connector is hollow and allows the insufflation of air if this is required by medical center procedures for the determination of the enteral tube position after the complete insertion. The stylet **SHALL** be removed **ONLY** as the position of the tube has been confirmed.

Contraindications:

- Abnormalities of the oesophageal and intestinal tract (e.g. stenosis, perforation, neoplasm etc.).
- Uncooperative patients, these patients have a greater risk of aspiration.
- Coagulopathy / thrombocytopenia.
- Known sensitivities or allergies to the materials and substances used in the product.
- Anatomical anomalies or serious alterations in the nasal, tracheal, oesophageal, pharyngeal or gastric regions.
- Threatening haemorrhage of oesophageal varices.
- Pre-existing contraindications to enteral nutrition.
- Long-term enteral nutrition (percutaneous tube techniques are preferable for this indication).

Side effects:

- Inadvertent malposition.
- Pulmonary aspiration or aspiration pneumonia.
- Inadvertent tube removal.
- Tube obstruction (clogging).
- Diarrhoea.

- Knotting.
- Ulceration.
- Breakdown of the skin surrounding the feeding tube.
- Ab Ingestis Pneumonia caused by vomiting, gastric reflux or malpositioning.
- Nostril bedsores.
- Clogged or dislodged feeding tube.

Warnings:

- If a stylet was present during the insertion, it shall be removed after the positioning, since it can obstruct the flow of nutrients and it can create interferences with NMR diagnostic techniques.
- In case of presence of stylet, never re-insert it when the tube is in patient. Inserting the stylet into an indwelling tube can result in perforation of the tube, leading to injury of the surrounding organs and haemorrhage.
- Before inserting the feeding tube, carefully check that the tip is smooth and with no burrs in order to avoid patient irritation/damage. Discard defective tubes.
- Start therapy after confirmation of the correct positioning of the feeding tube is recommended according to national/hospital guidelines.
- If resistance is noted during tube insertion, stop procedure, notify physician (if not present during the insertion), retract the tube and re-attempt insertion once patient is comfortable. Cases have been reported in which inadvertent placement of a feeding tube in the lung has resulted in pneumothorax in neurologically compromised or intubated patients.
- This is a disposable device intended for single use. Do not reuse. Re-use of single use devices creates a potential risk to patient or user. It may lead to contamination and/ or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

- Not for administration and admixture of liquids and drugs which mean to be infused via IV-access (e.g. IV solutions/drugs).
- Do not use Nutritub® after the expiry date printed on the packaging. The expiry date only applies if the product is correctly stored.
- Confirmation of the correct positioning of the feeding tube with X-ray technique is recommended.
- There is a risk that the tube may slip into the trachea. Insertion of the tube and the eventual stylet into the bronchial tree, can result in damage to the lung with severe pneumonia.
- Vigorous syringe force should not be used to irrigate, administer liquids or unblock the tube.

Precautions:

- Feeding tubes should be irrigated regularly in accordance with accepted medical practice with an ENFit® syringe to keep their patency:
 - a) Before and after each feeding as well as after any interruption of feeding;
 - b) Before and after each administration of medication (Important: medication should be administered only in liquid form to rule out the danger of tube occlusion);
 - c) Once a day when the tube is not in use during continuous feeding.
- Do not use if packaging or product is damaged.
- Place and replace feeding tubes according to medical standards techniques described in medical text books and according to national and/or local hospital guide lines.
- Remove the tube after maximum four weeks.
- Apply hygienic requirements according to accepted medical practice.
- Use only non-carbonated bottle water at room temperature to irrigate the tube and clean the connector.
- Do not attempt to unclog a blocked tube with force; in this case, replace the tube.

- Perform periodical checks of position.
- Pay attention when handling used devices to avoid patient/operator/environment contamination.
- It is under the physician's responsibility to determine the correct feeding tube size depending on the patient's characteristics.
- Furthermore, the tube may be displaced during surgical procedures, patient nursing and vomiting; it is necessary to assure the correct positioning before starting enteral nutrition.
- The doctor in charge is responsible for choosing the suitable drug to be administered to the patient via feeding tube basing on: tube characteristics, medical and analytical findings, general health condition of the patient, individual risk-benefit ratio.

Storage conditions:

- Keep dry.
- Keep away from sunlight.
- Store between 0°C and 30°C.

Instructions for Use

Read all instructions before using the enteral feeding tube.

Please refer to the pictograms displayed at the beginning of this IFU.

The operative actions described in the pictures are explained below:

- Place the patient in a beach chair position, with a degree between 30° and 45° (fig A), leaning slightly backwards. Fully explain the procedure and the associated side effects to the patient (if collaborative).
- Inside the package, pre-lubricated stylet is already inserted in the feeding tube. Make sure that it is firmly secured inside the tube. If, for any reason, it has been removed, insert it back

inside the tube before positioning (fig B) and verify its patency. Stretch the tube to remove any memory retained from the packaging.

- Identify the insertion depth (fig C): calculate the distance from the auricular lobe to the nose and from this to the xiphoid process (NEX measurement). NEX measurement is used as an estimation of the insertion depth, more precise methods to confirm the correct positioning shall be used when the tube is inserted. For intestinal tubes add 30 – 40 cm. For gastric tubes add 5 cm. Take note of the tube length. Before proceeding with the insertion, use a water soluble anaesthetizing gel/spray and cover the tube tip with water soluble lubricant that, in case of erroneous tube positioning in the airways, is easily absorbed without the risk of obstruction. Identify nostril with best ventilation. Ask the patient, if collaborative, which is the nostril with the best ventilation.
- Insert the tube in the nostril with best ventilation following hospital procedures. Observe nasal hygiene, verifying that there are no ulcers and rashes. Ask the patient to keep the head slightly bent backward to align the pharynx and the larynx, simplifying the insertion (fig D).
- If the patient is collaborative, ask him to swallow in order to occlude the tracheal hollow and prevent the hazardous tube insertion into the trachea and simultaneously drink small amounts of water (fig E); visually check the presence of the tube in the back of the throat. Direct tube posteriorly, aiming tip parallel to nasal septum and superior surface of hard palate. Advance tube to nasopharynx, allowing tip to seek its own passage.
- In case of a gastric tube, maintaining the position described in figure D, advance tube

through oesophagus into stomach with gentle motion, and proceed with the insertion of the tube until the chosen depth is met (fig F).

- G. In case of an intestinal tube, lay down the patient on the right side till the position required is reached (fig G).
- H. Check and verify the position of the tube, without removing the stylet (if present) (fig H). The position of the tube shall be checked by: X-ray technique, suction of gastric content, auscultation during air insufflation. Note: the syringe can be directly connected to the stylet connector.
- I. Position range of the tube in case of gastric insertion (fig I).
- J. Position range of the tube in case of intestinal insertion (fig J).
- K. If the stylet was present in the version, unscrew and remove it, maintaining the tube in place (fig K).
- L. Flush the tube lumen using 20 ml of water or as stated in the local medical practice (fig L).
- M. Fix the tube in place near the nostril and on the cheek, using appropriate devices and according to accepted medical practice (fig M). Begin feeding according to physician's order and usual institution protocol. At the end of the administration, disconnect the administration set and flush/irrigate the feeding tubes. Clean the ENFit® connector and close the feeding tube with the dedicated capsule.

Feeding tube removal:

- Before feeding tube removal it is suggested to introduce 30 mL air volume in order to remove eventual gastric secretions that could be irritating in contact with different tissue than the gastric ones.
- Remove the tube with gentle motion.
- During feeding tube removal suggest the patient (if collaborative) to take a deep breath and keep the breath in order to promote the glottis closure and reduce the risk of gastric material suction in the airways
- After use, this product might be a potential hazard. Handle and dispose it in accordance with accepted medical practice and according to local and federal laws and regulations.

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B.Braun Avitum Italy S.p.A.,