## Produkto informacija

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## Product fact sheet

# Sentimag® for lymph node and lesion localisation

Version 3.0

## Product name

Sentimag®

## Manufacturer information

Endomagnetics Ltd.

[www.endomag.com](http://www.endomag.com)

## Summary

The Sentimag® instrument is used during cancer treatment to localise with Magtrace® marked lymph nodes and with Magseed® marked soft tissue for their surgical removal. The instrument uses magnetic susceptometry and generates a magnetic field that temporarily magnetises the iron oxide in the Magtrace tracer and the Magseed® marker. With the probe the magnetic signatures generated by Magtrace® and Magseed® are detected. The device can be used pre- and post-incision and is simple to use after a short time of familiarisation.

The Sentimag® instrument comprises a portable mains-powered base unit, a detachable hand-held probe that is connected to the base unit with a flexible cable, an optional probe holder and a detachable air-operated footswitch that is connected to the base unit with a flexible hose.

The sensing of Magtrace® or Magseed® is indicated by a change in frequency of an audio signal from the base unit. A visible numerical representation of the detected signal level is simultaneously displayed on the base unit's liquid crystal display.

## Sentimag® product benefits for SLNB

- The Sentimag® - Magtrace® system provides the best standard of cancer staging, regardless of the hospital setting
- Magtrace® is injected when it is convenient - from 20 min up to 30 days prior to surgery
- Eliminate issues with radioactive materials; reach equivalent clinical outcomes [1-3]
- Highly accurate, proximity-based detection for intuitive node localisation
- Gamma system users will have no problems adapting to the Sentimag® procedure as the probe handling is remarkably similar

[1] Alvarado et al. (2019): Ann Surg Oncol. 26(11):3510–6.

[2] Karakatsanis et al. (2016): Breast Cancer Res Treat. 157(2):281–94.

[3] Teshome et al. (2016): Ann Surg Oncol. 23(5):1508–14.

## Sentimag® product benefits for soft tissue localisation

- Highly accurate, distance-based detection for intuitive magnetic soft tissue localisation incl. lesions and lymph nodes
- 360-degree, depth- and side-sensing enable distance calculation which makes the surgery particularly accurate
- The world's only system for magnetic lesion and lymph node localisation

## Article numbers

Article no.	EAN / GTIN	Item	Q
SE000203	Bill of Material (BOM)	Sentimag System SMS02 BOM consists out of SMB02, SMP02, SMA, SMC02, SMH02	1
SMB02	15060391210008	Sentimag Base Unit	1
SMP02	15060391210015	Sentimag Probe Assembly	1
SMA	15060391210022	Sentimag Footswitch Assembly	1
SMC02	15060391210039	Sentimag Probe Test Phantom	1
SMH02	15060391210012	Sentimag Probe Holder	1

Article no.	EAN / GTIN	Item	Q
SPC02UK	25060391210043	Sentimag Power Cord - UK	1
SPC02EU	45060391210047	Sentimag Power Cord - EU	1
SPC02CH	45060391210054	Sentimag Power Cord - CH	1
SPC02DK	25060391210050	Sentimag Power Cord - DK	1
SPC02IL	15060391210220	Sentimag Power Cord - IL	1
SPC02SA	15060391210237	Sentimag Power Cord - SA	1

## Features

Feature	Description
<b>Regulatory info</b>	<ul style="list-style-type: none"> <li>• CE-marked as medical device class IIa</li> <li>• FDA cleared</li> </ul>
<b>Working principle</b>	<ul style="list-style-type: none"> <li>• The Sentimag® is a magnetic susceptometer. It delivers small amplitude, time-varying magnetic field via a hand-held probe, and electronically detects the presence of any magnetic materials in the vicinity of the probe head.</li> <li>• The clinically introduced tracer Magtrace® and the marker Magseed® have a high magnetic susceptibility and are suited for the detection with the Sentimag®.</li> </ul>
<b>Instrument dimensions</b>	<ul style="list-style-type: none"> <li>• (W x H x D):</li> <li>• 240 mm x 370 mm x 210 mm</li> </ul>
<b>Instrument weight</b>	<ul style="list-style-type: none"> <li>• 3.8 kg (4.5 kg in case)</li> </ul>
<b>Probe weight</b>	<ul style="list-style-type: none"> <li>• 1.1 kg</li> </ul>
<b>Probe diameter</b>	<ul style="list-style-type: none"> <li>• At the tip: 18.5 mm</li> <li>• At thickest point: 20.5 mm</li> </ul>
<b>Probe head length</b>	<ul style="list-style-type: none"> <li>• 235 mm</li> </ul>
<b>Cable lengths</b>	<ul style="list-style-type: none"> <li>• Power cord: 1.8 m</li> <li>• Probe cable: 2.75 m</li> <li>• Footswitch cable: 3.06 m</li> </ul>

Feature	Description
<b>Shown on the LCD-display</b>	<ul style="list-style-type: none"> <li>• Numerical identification of the signal strength (ferrous signals: yellow; non-ferrous signals: red)</li> <li>• Probe generation 1 or 2</li> <li>• Volume adjustment</li> <li>• Sensitivity setting adjustment</li> <li>• During balancing a scale symbol</li> <li>• Error Codes</li> </ul>
<b>Default sensitivity setting</b>	<ul style="list-style-type: none"> <li>• Magtrace®: setting 1</li> <li>• Magseed®: setting 2</li> </ul>
<b>Cleaning and disinfection</b>	<ul style="list-style-type: none"> <li>• The Sentimag® probe (including cable) must be cleaned both immediately before and immediately after use.</li> <li>• Only use legally marketed wipes and solutions intended for medical devices and indicated as suitable for plastics.</li> <li>• Please contact your local representative in case of any questions.</li> </ul>
<b>Operating temperature range</b>	<ul style="list-style-type: none"> <li>• Base unit: 18 °C to 32 °C</li> <li>• Probe: 18 °C to 39 °C</li> </ul>
<b>Storage temperature range</b>	<ul style="list-style-type: none"> <li>• 0 °C to 40 °C</li> <li>• Prevent extreme temperatures</li> </ul>
<b>Transportation temperature range</b>	<ul style="list-style-type: none"> <li>• -10 °C to 50 °C</li> <li>• Prevent extreme temperatures</li> </ul>
<b>Atmospheric pressure range</b>	<ul style="list-style-type: none"> <li>• Operating, storage and transportation:</li> <li>• 80 kPa to 105 kPa</li> </ul>
<b>Use of device</b>	<ul style="list-style-type: none"> <li>• Less than 2000 metres</li> </ul>
<b>Power lead connector</b>	<ul style="list-style-type: none"> <li>• IEC 60320 C13</li> </ul>
<b>Power supply voltage</b>	<ul style="list-style-type: none"> <li>• 110 V to 230 V AC</li> <li>• 5 Hz to 60 Hz nominal</li> </ul>
<b>Fuse type and rating</b>	<ul style="list-style-type: none"> <li>• T1AH 250 V Ø 5x20 mm</li> </ul>
<b>Protection against electrical shock</b>	<ul style="list-style-type: none"> <li>• Class I Protectively Earthed</li> </ul>
<b>Applied part isolation</b>	<ul style="list-style-type: none"> <li>• Type B</li> </ul>
<b>Ingress protection rating</b>	<ul style="list-style-type: none"> <li>• IPX0 (not protected)</li> <li>• Parts rated for 5 years lifetime</li> </ul>

Feature	Description
<b>Sentimag® heat emission</b>	<ul style="list-style-type: none"> <li>• The Sentimag® system has fully passed IEC 60601-1 3rd Edition.</li> <li>• Under normal use, the hottest Sentimag® part is the power supply unit, with a temperature of 59.4 °C.</li> <li>• PCB boards reach a temperature of 34.8 – 52.5 °C.</li> <li>• The outer handle achieves a temperature of 29.3 °C.</li> <li>• The probe tip achieves 25.6 °C and the probe handle 25.2 °C.</li> <li>• These values are all significantly below the maximum allowed.</li> <li>• There has been no study on what the overall radiated heat emission is from a system.</li> </ul>

<b>Accuracy</b>	<ul style="list-style-type: none"> <li>• The device provides a qualitative measurement of the presence of magnetic material in the vicinity of the probe.</li> </ul>
-----------------	--

<b>Performance</b>	<ul style="list-style-type: none"> <li>• In the longitudinal direction the signal decreases approx. exponentially with distance.</li> <li>• In the transverse direction the signal decreases in a Gaussian-like manner with the distance from the tip of the probe.</li> </ul>
--------------------	--

<b>Magnetic field type</b>	<ul style="list-style-type: none"> <li>• The magnetic field is generated by passing a 10 kHz sinusoidally varying alternating current through a wire coil mounted in the probe.</li> </ul>
----------------------------	--

<b>Magnetic field intensity</b>	<ul style="list-style-type: none"> <li>• The maximum field intensity at any point on the probe casing is no greater than 251 µT.</li> </ul>
---------------------------------	---

<b>Magnetic field distribution</b>	<ul style="list-style-type: none"> <li>• The magnetic field distribution around the probe head is linearly proportional to the spatial variation in signal.</li> </ul>
------------------------------------	--

<b>Electromagnetic emissions</b>	<ul style="list-style-type: none"> <li>• The Sentimag® uses RF energy only for its internal function.</li> <li>• Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</li> <li>• There is no known risk of reciprocal interference posed by Sentimag® with any other equipment during specific investigations or treatments.</li> <li>• The Sentimag® is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</li> </ul>
----------------------------------	--

Feature	Description
<b>Electromagnetic Immunity</b>	<ul style="list-style-type: none"> <li>• Floor should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.</li> <li>• Mains power quality should be that of a typical commercial and/or hospital environment.</li> </ul>

<b>Disposal</b>	<ul style="list-style-type: none"> <li>• The Sentimag® instrument should be disposed of in accordance with the European Union WEEE Directive 2002/96/EC, on Waste Electrical and Electronic Equipment.</li> </ul>
-----------------	---