# Endo-SPONGE®

<b>(DE)</b>	<b>Gebrauchsanweisung</b>	<b>– 3</b>	<b>(CZ)</b>	<b>Návod k použití</b>	<b>– 42</b>
<b>(GB)</b>	<b>Instructions for use</b>	<b>– 6</b>	<b>(SK)</b>	<b>Návod na použitie</b>	<b>– 45</b>
<b>(FR)</b>	<b>Mode d'emploi</b>	<b>– 9</b>	<b>(SI)</b>	<b>Navodila za uporabo</b>	<b>– 48</b>
<b>(ES)</b>	<b>Instrucciones de empleo</b>	<b>– 12</b>	<b>(PL)</b>	<b>Instrukcja użycia</b>	<b>– 51</b>
<b>(IT)</b>	<b>Istruzioni per l'uso</b>	<b>– 15</b>	<b>(LV)</b>	<b>Lietošanas instrukcija</b>	<b>– 54</b>
<b>(PT)</b>	<b>Instruções de utilização</b>	<b>– 18</b>	<b>(LT)</b>	<b>Naudojimo nurodymai</b>	<b>– 57</b>
<b>(SE)</b>	<b>Bruksanvisning</b>	<b>– 21</b>	<b>(EE)</b>	<b>Kasutusjuhend</b>	<b>– 60</b>
<b>(NL)</b>	<b>Gebruiksaanwijzing</b>	<b>– 24</b>	<b>(TR)</b>	<b>Kullanma talimatları</b>	<b>– 63</b>
<b>(DK)</b>	<b>Brugsanvisning</b>	<b>– 27</b>	<b>(RO)</b>	<b>Instrucţiuni de folosire</b>	<b>– 66</b>
<b>(NO)</b>	<b>Bruksanvisning</b>	<b>– 30</b>	<b>(BG)</b>	<b>Инструкции за употреба</b>	<b>– 69</b>
<b>(FI)</b>	<b>Käyttöohjeet</b>	<b>– 33</b>	<b>(SG)</b>	<b>Instructions for use</b>	<b>– 72</b>
<b>(GR)</b>	<b>Οδηγίες χρήσης</b>	<b>– 36</b>	<b>(USA)</b>	<b>Instructions for use</b>	<b>– 75</b>
<b>(HU)</b>	<b>Használati utasítás</b>	<b>– 39</b>			

**CE** 0123

01-09/21 506999

**B | BRAUN**

SHARING EXPERTISE

B. Braun Surgical, S.A.  
Carretera de Terrassa, 121  
08191 Rubí, Spain

Phone +34 93 586 62 00  
Fax +34 93 588 10 96

[www.bbraun.com](http://www.bbraun.com)

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# Endo-SPONGE®

### Produktbeschreibung

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#### a) Insertion of Endo-SPONGE®

1. Prepare the patient according to the usual standards for an endoscopy and explain the treatment method to him or her.
2. Before use, the leak cavity should be sterilised using appropriate methods, such as irrigating with Ringer's solution.
3. Measure the leak cavity (length and diameter) visually using the endoscope in order to determine the length and diameter of the sponge. The size of the sponge should be chosen so that the sponge lies only within the leak cavity, but does not touch intact mucosa. After measurement, the endoscope is removed (Fig. 1).
4. Size the sponge according to the size of the leak cavity (several sponges may be inserted, depending on the size):
  - It is advisable that the sponge fills the leak cavity without touching the intact mucosa.
  - When cutting the sponge, round off any edges and points.
  - After cutting, remove all sponge residues and particles from the sponge surface.
  - Shorten the sponge length by up to 1 cm by cutting off the sponge at the end.
  - When shortening the sponge by more than 1 cm, shorten the drainage tube as well so that the sponge protrudes at least 3 mm from the end of the tube; if necessary, suture the edges of the sponge together over the end of the drainage tube.
5. Select the correct overtube size for the endoscope and sponge size. Check that the endoscope can be pushed through the overtube. Only sponges that have been re-sized (in length and width) can be used with the small overtube. Use sterile hydrogel based on glycerol (not supplied in the set), if necessary to insert the endoscope into the overtube.
6. When a rigid rectoscope or proctoscope is used, it will serve as the overtube as well. In this case, the length of the pusher must be properly adapted to the length of the rectoscope or proctoscope used; the pusher must be no more than 1 cm longer than the rectoscope or proctoscope).
7. First, push the non-tapered end of the overtube over the endoscope until the overtube tip (tapered end) is about 3 mm distance of the end of the endoscope.
8. Insert the endoscope and overtube (or rectoscope or proctoscope) carefully into the anal canal using suitable lubricant, and bring it to the level of the leak cavity (Fig. 2).
9. Push the endoscope with the overtube (or rectoscope or proctoscope) forward to the end of the leak cavity.
10. While viewing, push the overtube forward over the stationary endoscope until the overtube tip (tapered end) reaches the end of the leak cavity.
11. Fix the overtube in place and remove the endoscope.
12. Wet the overtube lumen (or rectoscope or proctoscope) and the re-sized sponge sufficiently with sterile hydrogel based on glycerol (not supplied in the set). Apply the sponge without delay; otherwise the sterile hydrogel based on glycerol may dry out (no more lubrication).
13. Compress the sponge and insert it fully into the overtube (or rectoscope or proctoscope). Continue to hold the overtube (or rectoscope or proctoscope) firmly in place to be sure to avoid further penetration.
14. Guide the pusher via the drainage tube and, using the pusher, push the sponge forward to the end of the overtube (or rectoscope or proctoscope). The tip of the sponge is now positioned at the deepest point in the leak cavity (Fig. 3). Meanwhile, continue to hold the overtube (or rectoscope or proctoscope) firmly in place.
15. Now keeping the pusher fixed in place, withdraw the overtube (or

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- A dwell time of 48 hours is recommended; more than 72 hours must be ruled out due to the risk of granulation tissue overgrowing into the sponge, as a result the sponge could get broken during removal, leaving part of it into the application area and embedded in granulation tissue. In case this happened, an endoscopic loop must snare or forceps could be optionally used to detach the sponge from the surrounding tissue for its removal.

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Do not reuse

Use-by date Year+Month+Date

Sterile unless package is opened or damaged. Method of sterilization: Ethylene Oxide

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# Endo-SPONGE®

### Description of the device

Endo-SPONGE, an endoluminal vacuum therapy (EVT) device, is a minimally invasive method for the treatment of anastomotic or Hartmann's stump leakages in the lower pelvic area.

The Endo-SPONGE system consists of a drainage tube with an attached open-pore polyurethane sponge, a connection system (Y-piece with an extensor tube), an application system (overtube, pusher) and a rinsing set (syringe, tip). The drainage consists of a Redon-Drainage (12 CH) with the open-pore sponge attached on the end. The area of the drainage tube in the sponge has perforated holes on its side. The size of the sponge on the drainage tube can be cut to size for the respective application.

The connection system consists of a Y-piece with a double-drain connector on one end where can be connected up to two redon- drainage systems and a hole on the other end for the junction with the extensor tube. The Y-piece-extensor tube ensemble permits the connection with the vacuum source (Redyrob® Trans Plus bottle) by means of luer lock connector.

The Endo-SPONGE application system consists of one coaxially arranged tube, the insertion tube made of silicone (overtube) and the inner more rigid tube with a handle (pusher). The lumen of the overtube is slightly larger than the outside diameter of the applied endoscope and is used as a guide for the insertion of the sponge system. The pusher is used to push forward and position the sponge.

### Contents

1. Endo-SPONGE®, open-pore PUR sponge (ø 3.3 x 7.5 cm) with Redon drain CH12, med. PVC, 40 cm long
2. Pusher, ABS + PVC, CH 30, 30 cm long
3. Overtubes in 2 sizes, depending on device and sponge size
  - Silicon tube, each 29 cm long
  - Tapered rounded tip
  - Size 1: inner diameter 13 mm, outer diameter 17 mm
  - Size 2: inner diameter 15 mm, outer diameter 19 mm
4. Irrigation set comprised of 20 ml syringe + cap + slide clamp
5. Y connecting tube with Luer lock fitting for Redyrob® Trans Plus bottle (purchased separately – Art. No. 5526604)

### NOT INCLUDED

1. Redyrob® Trans Plus – controllable wound drainage system, Art. No. 5526604
2. Sterile hydrogel based on glycerol
3. Standard or therapeutic flexible endoscope
4. Additional endoscopic accessories to be used at the physician discretion

### Materials used

Acrylonitrile Butadiene Styrene (ABS), Polyvinyl chloride (PVC), Polyurethane (PUR), Polyethyleneterephthalate (PET), Hydrogel, Silicone, Polyethylene (PE), Polypropylene (PP), Isoprene rubber (IR).

### Indications for use

- Treatment of anastomotic or Hartmann's stump leakages following colorectal surgery in the lower pelvic area (extraperitoneal position) by means of negative pressure.

The leak must have created a drainable cavity with or without a local infection.

### Mode of action

Endo-SPONGE consists of an open-pored sponge connected to a drainage tube. After endoscopic insertion of the sponge into the leakage cavity the drainage tube is routed out through the anus and connected to a vacuum system. The application of the vacuum leads to a continuous drainage of the fluid and the sponge in the cavity promotes the cleaning of the surface. To achieve an effective treatment the size of the sponge is cut to fit the cavity. Depending on the size of the leakage cavity it can be necessary to place more than one sponge into the cavity. The sponge system is changed every 48-72 hours. To change the sponge the vacuum is disconnected. Removal of the sponge must be done with prior irrigation with 0.9 % saline solution to remove the granulated tissue from the surface of the sponge. The sponge is removed through the anus and size of the new sponge is adapted to the size of the leakage cavity. Endo-SPONGE treatment is stopped as far as the cavity reaches a size of 2 x 1 cm, because no further reduction of the sponge size is technically possible.

### Contraindications

- Malignant tumor wound
- Necrotic tissue/gangrene
- Untreated osteomyelitis
- Sponge position directly adjacent to vessels, urinary bladder or small intestinal loops
- Non-drainable septic focus
- Clotting disorders
- Treatment with a therapeutic dose of anticoagulant or platelets aggregation inhibitor drugs
- Generalised peritonitis or sepsis
- Patients with known sensitivities or allergies to its components (refer to **Materials used**, for details about components name).

### Materials used, for details about components name).

- There is no clinical evidence from the use of Endo-SPONGE in paediatric population. As a precautionary measure, Endo-SPONGE is not recommended in paediatric population.
- There is no clinical evidence from the use of Endo-SPONGE in pregnant and breastfeeding women. As a precautionary measure, Endo-SPONGE is not indicated for use during pregnancy or breastfeeding.

### Mode of application

The following steps are required for the application of Endo-SPONGE®.


Fig. 1

Situs


Fig. 2

Insertion of endoscope and overtube.


Fig. 3

Correct positioning, PUR sponge is released in the wound cavity.


Fig. 4

Correct position of the PUR sponge in the wound cavity after removal of the pusher and overtube. Suction is applied.

#### a) Insertion of Endo-SPONGE®

1. Prepare the patient according to the usual standards for an endoscopy and explain the treatment method to him or her.
2. Before use, the leak cavity should be sterilised using appropriate methods, such as irrigating with Ringer's solution.
3. Measure the leak cavity (length and diameter) visually using the endoscope in order to determine the length and diameter of the sponge. The size of the sponge should be chosen so that the sponge lies only within the leak cavity, but does not touch intact mucosa. After measurement, the endoscope is removed (Fig. 1).
4. Size the sponge according to the size of the leak cavity (several sponges may be inserted, depending on the size):
  - It is advisable that the sponge fills the leak cavity without touching the intact mucosa.
  - When cutting the sponge, round off any edges and points.
  - After cutting, remove all sponge residues and particles from the sponge surface.
  - Shorten the sponge length by up to 1 cm by cutting off the sponge at the end.
  - When shortening the sponge by more than 1 cm, shorten the drainage tube as well so that the sponge protrudes at least 3 mm from the end of the tube; if necessary, suture the edges of the sponge together over the end of the drainage tube.
5. Select the correct overtube size for the endoscope and sponge size. Check that the endoscope can be pushed through the overtube. Only sponges that have been re-sized (in length and width) can be used with the small overtube. Use sterile hydrogel based on glycerol (not supplied in the set), if necessary to insert the endoscope into the overtube.
6. When a rigid rectoscope or proctoscope is used, it will serve as the overtube as well. In this case, the length of the pusher must be properly adapted to the length of the rectoscope or proctoscope unless the pusher must be no more than 1 cm longer than the rectoscope or proctoscope).
7. First, push the non-tapered end of the overtube over the endoscope until the overtube tip (tapered end) is about 3 mm distance of the end of the endoscope.
8. Insert the endoscope and overtube (or rectoscope or proctoscope) carefully into the anal canal using suitable lubricant, and bring it to the level of the leak cavity (Fig. 2).
9. Push the endoscope with the overtube (or rectoscope or proctoscope) forward to the end of the leak cavity.
10. While viewing, push the overtube forward over the stationary endoscope until the overtube tip (tapered end) reaches the end of the leak cavity.
11. Fix the overtube in place and remove the endoscope.
12. Wet the overtube lumen (or rectoscope or proctoscope) and the re-sized sponge sufficiently with sterile hydrogel based on glycerol (not supplied in the set). Apply the sponge without delay; otherwise the sterile hydrogel based on glycerol may dry out (no more lubrication).
13. Compress the sponge and insert it fully into the overtube (or rectoscope or proctoscope). Continue to hold the overtube (or rectoscope or proctoscope) firmly in place to be sure to avoid further penetration.
14. Guide the pusher via the drainage tube and, using the pusher, push the sponge forward to the end of the overtube (or rectoscope or proctoscope). The tip of the sponge is now positioned at the deepest point in the leak cavity (Fig. 3). Meanwhile, continue to hold the overtube (or rectoscope or proctoscope) firmly in place.
15. Now keeping the pusher fixed in place, withdraw the overtube (or

- rectoscope or proctoscope), but no further than the funnel of the pusher (the sponge will now expand in the leak cavity).
16. Remove the pusher and overtube (or rectoscope).
17. Check the position of the sponge using the endoscope. (Fig. 4)
18. Connect the enclosed Y connecting tube with Luer lock fitting and Redyrob® Trans Plus bottle.
19. Set the control knob on the bottle to suction power 1 (see Redyrob® Trans Plus instruction leaflet). Do not use conventional high vacuum or medium vacuum wound drainage bottles, as they produce too much suction. Observe the flow of secretion and check the quantity. Abduction of secretion normally begins immediately. If it does not, the connection between the Endo-SPONGE® and the Redyrob® Trans Plus bottle must be checked, as well as the setting of the control knob.
20. Additional sealing is not necessary, since the patient's sphincter apparatus will create a tight closure.
21. Additional fixing of the drainage, which is sticking out per anum, is likewise unnecessary, since the sponge is kept in the leak cavity by the application of suction.
22. If needed, several Endo-SPONGE® can be placed consecutively in a leak cavity.
23. Maximum use of Endo-SPONGE®: 72 hours. By that point, it must be replaced, since the growth of granulation tissue into the sponge pores can make it difficult to remove the sponge.
24. It is advisable to inspect the system regularly.

#### b) Removal of Endo-SPONGE®

25. Disconnect Endo-SPONGE® from Redyrob® Trans Plus.
26. Irrigate the leak cavity with about 20 ml of Ringer's solution or 0.9% saline solution using the enclosed irrigation set via the Redon drain in place and close the remaining Redon drainages with the slide clamp.
27. Removal of Endo-SPONGE®:
  - Pull off the drainage tube carefully, steadily and with increasing force until the sponge detaches from the leak cavity (loss of resistance).
  - After the sponge has detached, continue to pull with a small amount of force in order to remove the tube per anum.
28. Check the integrity of the sponge.
29. Using the endoscope, inspect the treated leak cavity for sponge residues, and document the outcome of treatment. If it is necessary to insert another Endo-SPONGE®, follow the steps described in point 3.
30. Remove sponge residue using suitable grasping endoscopy.
31. For the effective closure of the leak cavity, it is advisable to reduce the diameter of the sponge, while initially keeping the same length, in order to close the remaining canal from back to front in the subsequent treatment cycles.
32. We recommend using Endo-SPONGE® until the length of the leak cavity is less than 2 cm and the diameter is less than 1 cm.
33. After the end of treatment (no later than one week after), the persistent result of treatment must be documented.

#### Warnings

- The sponge could suffer damages when reshaping and/or removing, generating residual sponge particles.
- The residual sponge particles could cause fistula formation, foreign body reactions. Possible need of surgical removal.
- Endo-SPONGE® must not be used in body openings other than indicated.
- Due to the underlying disease most patients have a localized infection

which can lead to sepsis (i.e. peritonitis, necrosis...).

### Precautions

- This treatment should only be performed by experienced doctors with practice both in the interventional treatment of the lower gastrointestinal tract using flexible endoscopy and in negative pressure wound therapy in general.

- Particles and sponge residues are produced when the sponge is cut to size. The sponge must be cut (e.g. with scissors or a scalpel) at an appropriate distance from the patient and in a suitable environment in which particles are permitted.

- After cutting the sponge, all residues and particles must be removed from the surface by tapping the sponge, collected and disposed of in the usual manner.

- When reshaping the sponge, ensure that there are no sharp edges or points and round it off, otherwise these can easily break off during the removal of the sponge. After cutting no cut in the sponge or the drainage tube should exist.

- When shortening the sponge length also shorten the drainage tube. The sponge must protrude at least 3 mm from the end of the drainage tube.

- The dwell time depends on the local clinical situation. If applicable, it is recommended to adjust the time of sponge change depending on the amount of debris and growth of granulation tissue.

- A dwell time of 48 hours is recommended; more than 72 hours must be ruled out due to the risk of granulation tissue overgrowing into the sponge, as a result the sponge could get broken during removal, leaving part of it into the application area and embedded in granulation tissue. In case this happened, an endoscopic snare or forceps could be optionally used to detach the sponge from the surrounding tissue for its removal.

- The flow of secretion must be monitored and the quantity checked. The drainage of secretion usually begins immediately. If this is not the case, the connection between Endo-SPONGE and the Redyrob® Trans Plus bottle must be checked, as well as the setting of the control knob.

- All items are for single use.

- Do not use any items unless the packaging is intact.

- Be careful not to damage to the sheathing of endoscope

- Damage to cables of endoscope can happen in case of excess bending of flexible distal end within the overtube

### Side effects

- Erosion of structures adjacent to sponge (vessels, urinary bladder, small intestine, colon, etc.)

- Injury to intestinal wall and perforation has been reported in some cases.

- Bleeding, that depending on patient condition could lead to severe bleeding.

- Sponge dislocation

- Post-interventional stricture / stenosis

- Pain

### Sterilization

The Endo-SPONGE kit is sterilized by ethylene oxide gas. Opened or damaged packages must not be used.

### Storage

Endo-SPONGE® should be stored at room temperature. Do not expose to extremes of temperatures for long periods.

Use the properly stored product by the date stated on the packaging.

### Symbols used on labelling

Do not reuse!

Use until Year + Month

Sterile unless package is opened and damaged. Method of sterilization: Ethylene oxide

See I nstructions for use!

Batch number

Cat. No.

CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

Date of information: 06/2021

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost."

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
Title: FOLLETO ENDO-SPONGE V.10 Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)  
Title: SUT-PR-ES04-Técnico Sistemas Impresión  
Date: Friday, 10 September 2021, 08:41 W. Europe Daylight Time  
Meaning: Document signed as Author

---

---

UserName: García, Guillermo (garcgues)  
Title: SU-PR-ES04-Technical support manager  
Date: Friday, 10 September 2021, 09:18 W. Europe Daylight Time  
Meaning: Precheck of Document

---

---

UserName: Piera, Maria Eva (piermaes)  
Title: Senior Regulatory Affairs  
Date: Monday, 13 September 2021, 08:08 W. Europe Daylight Time  
Meaning: Approve Document

---

---