

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<b>AFI/GPeakSysSL_Endo(A2C) Alias: GPeakSysSL Endo(A2C)</b>	<b>(GEU-106-0001, 99GEMS, “Global Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0129, 99GEMS, “AFI on Endocardium”) (111031, DCM, “Image View”) = (G-A19B, SRT, “Apical two chamber”)</b>
<b>AFI/GPeakSysSL_Endo(A4C) Alias: GPeakSysSL Endo(A4C)</b>	<b>(GEU-106-0001, 99GEMS, “Global Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0129, 99GEMS, “AFI on Endocardium”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)</b>
<b>AFI/GPeakSysSL_Endo(APLAX) Alias: GPeakSysSL Endo(APLAX)</b>	<b>(GEU-106-0001, 99GEMS, “Global Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0129, 99GEMS, “AFI on Endocardium”) (111031, DCM, “Image View”) = (G-0395, SRT, “Apical long axis”)</b>
<b>AFI/GPeakSysSL_Endo(Avg) Alias: GPeakSysSL Endo(Avg)</b>	<b>(GEU-106-0001, 99GEMS, “Global Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0129, 99GEMS, “AFI on Endocardium”)</b>
<b>AFI/PSD Alias: PSD</b>	<b>(GEU-106-0131, 99GEMS, “Peak Strain Dispersion”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”)</b>
<b>AFI/PSD_ASE18 Alias: PSD ASE18</b>	<b>(GEU-106-0131, 99GEMS, “Peak Strain Dispersion”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”)</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		= (F-32020, SRT, “Systole”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0128, 99GEMS, “AFI with 18 segments following 2015 ASE recommendations”)
<b>AFI/PSD_ENDO Alias: PSD ENDO</b>	<b>(GEU-106-0131, 99GEMS, “Peak Strain Dispersion”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0129, 99GEMS, “AFI on endocardium”)</b>
<b>AFI/BA PeakSysSL Alias: BA PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (T-32619, SRT, “left ventricle basal anterior segment”)</b>
<b>AFI/BAS PeakSysSL Alias: BAS PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R-10075, SRT, “left</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>ventricle basal anteroseptal segment”)</b>
<b>AFI/BS PeakSysSL Alias: BS PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R- 10076, SRT, “left ventricle basal inferoseptal segment”)</b>
<b>AFI/BI PeakSysSL Alias: BI PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (T- 32615, SRT, “left ventricle basal inferior segment”)</b>
<b>AFI/BP PeakSysSL Alias: BP PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R- 10079, SRT, “left ventricle basal inferolateral segment”)</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<b>AFI/BL PeakSysSL Alias: BL PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R- 1007A, SRT, “left ventricle basal anterolateral segment”)</b>
<b>AFI/MA PeakSysSL Alias: MA PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (T- 32617, SRT, “left ventricle mid anterior segment”)</b>
<b>AFI/MAS PeakSysSL Alias: MAS PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R- 10077, SRT, “left ventricle mid anteroseptal segment”)</b>
<b>AFI/MS PeakSysSL Alias: MS PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0018, 99GEMS, “AFI”)</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R-10078, SRT, “left ventricle mid inferoseptal segment”)
<b>AFI/MI PeakSysSL Alias: MI PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (T-32616, SRT, “left ventricle mid inferior segment”)
<b>AFI/MP PeakSysSL Alias: MP PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (R-1007B, SRT, “left ventricle mid inferolateral segment”)
<b>AFI/ML PeakSysSL Alias: ML PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		(G-C0E3, SRT, “Finding Site”) = (R-1007C, SRT, “left ventricle mid anterolateral segment”)
<b>AFI/AA PeakSysSL Alias: AA PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (T-32613, SRT, “left ventricle apical anterior segment”)
<b>AFI/AAS PeakSysSL Alias: AAS PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (GEU-106-0026, 99GEMS, “left ventricle apical anteroseptal segment”)
<b>AFI/AS PeakSysSL Alias: AS PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, “Peak Longitudinal Strain”)</b>	(G-C036, SRT, “Measurement Method”) = (GEU-106-0018, 99GEMS, “AFI”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (G-C0E3, SRT, “Finding Site”) = (T-32614, SRT, “left

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		ventricle apical septal segment")
<b>AFI/AI PeakSysSL Alias: AI PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C036, SRT, "Measurement Method") = (GEU-106- 0018, 99GEMS, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") (G-C0E3, SRT, "Finding Site") = (T- 32618, SRT, "left ventricle apical inferior segment")</b>
<b>AFI/AP PeakSysSL Alias: AP PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C036, SRT, "Measurement Method") = (GEU-106- 0018, 99GEMS, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") (G-C0E3, SRT, "Finding Site") = (GEU-106-0025, 99GEMS, "left ventricle apical posterior segment")</b>
<b>AFI/AL PeakSysSL Alias: AL PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C036, SRT, "Measurement Method") = (GEU-106- 0018, 99GEMS, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") (G-C0E3, SRT, "Finding Site") = (T- 3261C, SRT, "left ventricle apical lateral segment")</b>

**DIRECTION DOC2652554 REV 3**

<b>4DAutoLAQ/Vmin</b> Alias: 4DAutoLAQ Vmin	(GEU-106-0123, 99GEMS, “Left Atrium minimal volume by 4D LA quantification tool”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/Vmax</b> Alias: 4DAutoLAQ Vmax	(G-0383, SRT, “Left Atrium Systolic Volume”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/VpreA</b> Alias: 4DAutoLAQ VpreA	(GEU-106-0124, 99GEMS, “Left Atrium volume at preA time by 4D LA quantification tool”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/GLS_R</b> Alias: 4DAutoLAQ GLS_R	(GEU-106-0115, 99GEMS, “Left Atrium global longitudinal strain in reservoir phase”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/GLS_CD</b> Alias: 4DAutoLAQ GLS_CD	(GEU-106-0116, 99GEMS, “Left Atrium global longitudinal strain in conduit phase”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/GLS_CT</b> Alias: 4DAutoLAQ GLS_CT	(GEU-106-0117, 99GEMS, “Left Atrium global longitudinal strain in contractile phase”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/GCS_R</b> Alias: 4DAutoLAQ GCS_R	(GEU-106-0118, 99GEMS, “Left Atrium global circumferential strain in reservoir phase”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)



**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<b>4DAutoLAQ/GCS_CD</b> Alias: 4DAutoLAQ GCS_CD	(GEU-106-0119, 99GEMS, “Left Atrium global circumferential strain in conduit phase”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/GCS_CT</b> Alias: 4DAutoLAQ GCS_CT	(GEU-106-0120, 99GEMS, “Left Atrium global circumferential strain in contractile phase”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/EV</b> Alias: 4DAutoLAQ EV	(GEU-106-0121, 99GEMS, “Left Atrium Emptying Volume by 4D LA quantification tool”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>4DAutoLAQ/EF</b> Alias: 4DAutoLAQ EF	(GEU-106-0122, 99GEMS, “Left Atrium Emptying Fraction by 4D LA quantification tool”)	(G-C036, SRT, “Measurement Method”) = (GEU-106- 0114, 99GEMS, “4D auto LA quantification tool”)
<b>CardiacWork/CW(GWI)</b> Alias: CW(GWI)	(GEU-106-0062, 99GEMS, “Global Work Index”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (G-0387, SRT, “Mitral Valve Closure to Opening Time”) (G-C036, SRT, “Measurement Method”) = (GEU-106- 0061, 99GEMS, “MyocardialWork Quantification Tool”)
<b>CardiacWork/CW(GCW)</b> Alias: CW(GCW)	(GEU-106-0063, 99GEMS, “Global Constructive Work”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (G-0387, SRT, “Mitral Valve Closure to Opening Time”) (G-C036, SRT, “Measurement Method”) = (GEU-106-

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>0061, 99GEMS, “MyocardialWork Quantification Tool”)</b>
<b>CardiacWork/CW(GWW) Alias: CW(GWW)</b>	<b>(GEU-106-0064, 99GEMS, “Global Wasted Work”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”) = (G-0387, SRT, “Mitral Valve Closure to Opening Time”) (G-C036, SRT, “Measurement Method”) = (GEU-106- 0061, 99GEMS, “MyocardialWork Quantification Tool”)</b>
<b>CardiacWork/CW(GWE) Alias: CW(GWE)</b>	<b>(GEU-106-0065, 99GEMS, “Global Work Efficiency”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”) = (G-0387, SRT, “Mitral Valve Closure to Opening Time”) (G-C036, SRT, “Measurement Method”) = (GEU-106- 0061, 99GEMS, “MyocardialWork Quantification Tool”)</b>
<b>CardiacWork/CW(BPS) Alias: CW(BPS)</b>	<b>(F-008EC, SRT, “Systolic Blood Pressure”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0061, 99GEMS, “MyocardialWork Quantification Tool”)</b>
<b>CardiacWork/CW(BPD) Alias: CW(BPD)</b>	<b>(F-008ED, SRT, “Diastolic Blood Pressure”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106- 0061, 99GEMS, “MyocardialWork Quantification Tool”)</b>
<b>PV A Dur-MV A Dur Alias: PV A- PV D</b>	<b>(GEU-106-0070, 99GEMS, “Difference between the Pulmonary Vein A- wave duration and the</b>	<b>(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)</b>

**DIRECTION DOC2652554 REV 3**

	mitral inflow A-wave duration")	
<b>PV A Dur/MV VTI Alias: PV A Dur/MV VTI</b>	(GEU-106-0071, 99GEMS, "Ratio of Pulmonary Vein A-wave duration by Mitral Valve VTI")	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)
<b>PV A Dur/MV A Dur Alias: PV A Dur/MV A Dur</b>	(GEU-106-0072, 99GEMS, "Ratio of Pulmonary Vein A-wave duration by Mitral inflow A-wave duration")	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)
<b>MV Eprime/Aprime Septal Ratio/Calc Alias : E'/A' Septal</b>	(G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave")	(G-C0E3, SRT, "Finding Site") = (SRT, G-0391, "Medial Mitral Annulus")
<b>MV Eprime/Aprime Lateral Ratio/Calc Alias : E'/A' Lateral</b>	(G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave")	(G-C0E3, SRT, "Finding Site") = (SRT, G-0392, "Lateral Mitral Annulus")
<b>MV Aprime Septal Velocity Alias : A' Septal</b>	(G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole")	(G-C0E3, SRT, "Finding Site") = (SRT, G-0391, "Medial Mitral Annulus")
<b>MV Aprime Lateral Velocity Alias : A' Lateral</b>	(G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole")	(G-C0E3, SRT, "Finding Site") = (SRT, G-0392, "Lateral Mitral Annulus")
<b>MV Sprime Septal Velocity Alias : S' Septal</b>	(G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity")	(G-C0E3, SRT, "Finding Site") = (SRT, G-0391, "Medial Mitral Annulus")
<b>MV Sprime Lateral Velocity Alias : S' Lateral</b>	(G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity")	(G-C0E3, SRT, "Finding Site") = (SRT, G-0392, "Lateral Mitral Annulus")

**Section Right Ventricle**

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<b>GEU Parameter ID (and corresponding alias)</b>	<b>Base Measurement Concept Name</b>	<b>Concept or Acquisition Context Modifier</b>
<b>TomTec/RVFunction/EDV</b>  <b>Alias: RVEDV(TomTec)</b>	<b>(8822-3, LN, “Right Ventricular ED Volume”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106-0022, 99GEMS, “4D Right Ventricle Volume”)</b>
<b>TomTec/RVFunction/ESV</b>  <b>Alias: RVESV(TomTec)</b>	<b>(8824-5, LN, “Right Ventricular ES Volume”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106-0022, 99GEMS, “4D Right Ventricle Volume”)</b>
<b>TomTec/RVFunction/SV</b>  <b>Alias: RVSV(TomTec)</b>	<b>(F-32120, SRT, “Stroke Volume”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106-0022, 99GEMS, “4D Right Ventricle Volume”)</b>
<b>TomTec/RVFunction/EF</b>  <b>Alias: RVEF(TomTec)</b>	<b>(10231-9, LN, “RV Ejection Fraction”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU-106-0022, 99GEMS, “4D Right Ventricle Volume”)</b>
<b>MM/RVOT</b>  <b>Alias: RVOT</b>	<b>(G-038F, SRT, “Cardiovascular Orifice Diameter”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”) (G-0373, SRT, “Image Mode”) = (G-0394, SRT, “M mode”)</b>
<b>2D/RVEF(A-L A4C)</b> <b>Alias: RVEF A-L A4C</b>	<b>(F-02268, SRT, “Right Ventricular Ejection Fraction”)</b>	<b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		(G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area-Length Single Plane”)
<b>2D/RVEF(MOD A4C)</b> <b>Alias: RVEF MOD A4C</b>	(F-02268, SRT, “Right Ventricular Ejection Fraction”)	(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125208, DCM, “Method of Disks, Single Plane”)
<b>2D/RVOT Diam</b> <b>Alias: RVOT Diam</b>	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
<b>2D/RVAWd</b> <b>Alias: RVAWd</b>	(18153-7, LN, “Right Ventricular Anterior Wall Diastolic Thickness”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
<b>2D/RVAWs</b> <b>Alias: RVAWs</b>	(18157-8, LN, “Right Ventricular Anterior Wall	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT,

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

	<b>Systolic Thickness”)</b>	<b>“Doppler Color Flow”) depending on scan mode</b>
<b>2D/RVIDd</b> <b>Alias: RVIDd</b>	<b>(20304-2, LN, “Right Ventricular Internal Diastolic Dimension”)</b>	
<b>2D/RVIDs</b> <b>Alias: RVIDs</b>	<b>(20305-9, LN, “Right Ventricular Internal Systolic Dimension”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode</b>
<b>2D/RVD Major</b> <b>Alias: RV length</b>	<b>(G-A193, SRT, “Major Axis”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode</b>
<b>2D/RVD Minor</b> <b>Alias: RV mid diameter</b>	<b>(G-A194, SRT, “Minor Axis”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode</b>
<b>2D/RVD Minor Base</b> <b>Alias: RV base diameter</b>	<b>(80080-5, LN, “Right ventricular base Minor Axis at end diastole, 4-chamber view”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32011, SRT, “End Diastole”) (111031, DCM, “Image View”) =</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		(G-A19C, SRT, “Apical four chamber”)
<b>2D/RVA Diastole</b> <b>Alias: RVA Diastole</b>	(GEU-106-0054, 99GEMS, “Right Ventricular Diastolic Area”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32011, SRT, “End Diastole”)
<b>2D/RVA Systole</b> <b>Alias: RVA Systole</b>	(GEU-106-0055, 99GEMS, “Right Ventricular Systolic Area”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109070, DCM, “End Systole”)
<b>2D/RVOT Area</b> <b>Alias: RVOT Area</b>	(G-038E, SRT, “Cardiovascular Orifice Area”)	(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
<b>2D/RV FAC</b> <b>Alias: RV FAC</b>	(78175-7, LN, “Right ventricular Fraction area change by US 2D”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”)
<b>4DAutoRVQ/EDV</b> <b>Alias: RV EDV</b>	(8822-3, LN, “Right Ventricular ED Volume”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”)

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		(G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109022, DCM, “End Diastole”)
<b>4DAutoRVQ/ESV</b> Alias: RV ESV	(8824-5, LN, “Right Ventricular ES Volume”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
<b>4DAutoRVQ/EF</b> Alias: RV EF	(F-02268, SRT, “Right Ventricular Ejection Fraction”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”)
<b>4DAutoRVQ/SV</b> Alias: RV SV	(F-02268, SRT, “Right Ventricular Stroke Volume”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”)



**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”)
4DAutoRVQ/Dd_base Alias: RV Dd base	(GEU-106-0049, 99GEMS, “RV basal minor axis at end diastole on A4C from 4D image”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109022, DCM, “End Diastole”)
4DAutoRVQ/Dd_mid Alias: RV Dd mid	(GEU-106-0050, 99GEMS, “RV mid-cavity minor axis at end diastole on A4C from 4D image”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (R-4089A, SRT, “Cardiac Cycle

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		Point”) = (109022, DCM, “End Diastole”)
4DAutoRVQ/Ld Alias: RV Ld	(GEU-106-0051, 99GEMS, “RV major axis at end diastole on A4C from 4D image”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109022, DCM, “End Diastole”)
4DAutoRVQ/FAC Alias: RV FAC	(GEU-106-0053, 99GEMS, “RV FAC on user-selected from 4D image”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)
4DAutoRVQ/EDV_index Alias: RV EDV index	(8822-3, LN, “Right Ventricular ED Volume”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”)

DIRECTION DOC2652554 REV 3

		(G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (R-4089A, SRT, “Cardiac Cycle Point”) = (109022, DCM, “End Diastole”) (121425, DCM, “Index”) = (8277-6, LN, “Body Surface Area”)
4DAutoRVQ/ESV_index Alias: RV ESV index	(8824-5, LN, “Right Ventricular ES Volume”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV quantification tool”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”) (121425, DCM, “Index”) = (8277-6, LN, “Body Surface Area”)
4DAutoRVQ/SV_index Alias: RV SV index	(F-02268, SRT, “Right Ventricular Stroke Volume”)	(G-0373, SRT, “Image Mode”) = (125231, DCM, “3D mode”) (G-C036, SRT, “Measurement Method”) = (GEU-106-0048, 99GEMS, “4D auto RV

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		quantification tool") (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area")
<b>RVLd(A4C)</b>  <b>Alias: RVLd(A4C)</b>	<b>(18078-6, LN, "Right Ventricular Major Axis Diastolic Dimension")</b>	<b>(111031, DCM, "Image View") = (G-A19C, SRT, "Apical four chamber")</b>
<b>RVAd(A4C)</b>  <b>Alias: RVAd(A4C)</b>	<b>(G-A166, SRT, "Area")</b>	<b>(111031, DCM, "Image View") = (G-A19C, SRT, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")</b>
<b>RVEDV(A-L A4C)</b>  <b>Alias: RVEDV A-L A4C</b>	<b>(8822-3, LN, "Right Ventricular ED Volume")</b>	<b>(111031, DCM, "Image View") = (G-A19C, SRT, "Apical four chamber") (G-C036, SRT, "Measurement Method") = (125205, DCM, "Area-Length Single Plane")</b>
<b>RVEDV(MOD A4C)</b>  <b>Alias: RVEDV(MOD A4C)</b>	<b>(8822-3, LN, "Right Ventricular ED Volume")</b>	<b>(111031, DCM, "Image View") = (G-A19C, SRT, "Apical four chamber") (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method of Disks, Single Plane")</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<p><b>RVLs(A4C)</b>  <b>Alias: RVLs(A4C)</b></p>	<p><b>(18079-4, LN, “Right Ventricular Major Axis Systolic Dimension”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”)</b></p>
<p><b>RVAs(A4C)</b>  <b>Alias: RVAs(A4C)</b></p>	<p><b>(G-A166, SRT, “Area”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)</b></p>
<p><b>RVESV(A-L A4C)</b>  <b>Alias: RVESV A-L A4C</b></p>	<p><b>(8824-5, LN, “Right Ventricular ES Volume”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area-Length Single Plane”)</b></p>
<p><b>RVESV(MOD A4C)</b>  <b>Alias: RVESV(MOD A4C)</b></p>	<p><b>(8824-5, LN, “Right Ventricular ES Volume”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125208, DCM, “Method of Disks, Single Plane”)</b></p>
<p><b>MM/RVIDd</b>  <b>Alias: RVIDd</b></p>	<p><b>(20304-2, LN, “Right Ventricular Internal Diastolic Dimension”)</b></p>	<p><b>(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)</b></p>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<b>MM/RVIDs</b> <b>Alias: RVIDs</b>	(20305-9, LN, “Right Ventricular Internal Systolic Dimension”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/RVAWd</b> <b>Alias: RVAWd</b>	(18153-7, LN, “Right Ventricular Anterior Wall Diastolic Thickness”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/RVAWs</b> <b>Alias: RVAWs</b>	(18157-8, LN, “Right Ventricular Anterior Wall Systolic Thickness”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/RVPEP</b> <b>Alias: RVPEP</b>	(20301-8, LN, “Right Ventricle Pre Ejection Period”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/RVET</b> <b>Alias: RVET</b>	(20222-6, LN, “Ejection Time”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/RVPEP/ET Ratio</b> <b>Alias: RVPEP ET Ratio</b>	(59088-5, LN, “Right Ventricular Pre- ejection time/Ejection Time”)	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>Est RVSP</b> <b>Alias: RVSP</b>	(G-0380, SRT, “Right Ventricular Peak Systolic Pressure”)	
<b>RVOT Vmax</b> <b>Alias: RVOT Vmax</b>	(11726-7, LN, “Peak Velocity”)	(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<b>RVOT Vmax P</b> <b>Alias: RVOT Vmax</b>	<b>(11726-7, LN, “Peak Velocity”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT maxPG</b> <b>Alias: RVOT maxPG</b>	<b>(20247-3, LN, “Peak Gradient”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT Vmean</b> <b>Alias: RVOT Vmean</b>	<b>(20352-1, LN, “Mean Velocity”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT meanPG</b> <b>Alias: RVOT meanPG</b>	<b>(20256-4, LN, “Mean Gradient”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT VTI</b> <b>Alias: RVOT VTI</b>	<b>(20354-7, LN, “Velocity Time Integral”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT Env.Ti</b> <b>Alias: RVOT Env.Ti</b>	<b>(GEU-106-0085, 99GEMS, “Time duration of the VTI trace on RVOT”)</b>	<b>(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)</b>
<b>RVOT HR</b> <b>Alias: HR</b>	<b>(8867-4, LN, “Heart rate”)</b>	
<b>RVOT SV</b> <b>Alias: RVOT SV</b>	<b>(F-32120, SRT, “Stroke Volume”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT SI</b> <b>Alias: RVOT SI</b>	<b>(F-00078, SRT, “Stroke Index”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>“Right Ventricle Outflow Tract”)</b>
<b>RVOT CO</b> <b>Alias: RVOT CO</b>	<b>(F-32100, SRT, “Cardiac Output”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVOT CI</b> <b>Alias: RVOT CI</b>	<b>(F-32110, SRT, “Cardiac Index”)</b>	<b>(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVET</b> <b>Alias: RVET</b>	<b>(79929-6, LN, “Right Ventricular Ejection Time by US doppler”)</b>	<b>(G-0373, SRT, “Image Mode”) = (R-409E4, SRT, “Doppler Pulsed”) or (G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”) (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVPEP</b> <b>Alias: RVPEP</b>	<b>(20301-8, LN, “Right Ventricle Pre Ejection Period”)</b>	<b>(G-0373, SRT, “Image Mode”) = (R-409E4, SRT, “Doppler Pulsed”) or (G-0373, SRT, “Image Mode”) = (R-409E3, SRT, “Doppler Continuous Wave”) (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)</b>
<b>RVPEP/ET Ratio</b> <b>Alias: RVPEP/ET Ratio</b>	<b>(59088-5, LN, “Ventricular Pre Ejection</b>	<b>(G-0373, SRT, “Image Mode”) =</b>



**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

	Time/Ejection time by US")	(R-409E4, SRT, "Doppler Pulsed") (G-C0E3, SRT, "Finding Site") = (T-32550, SRT, "Right Ventricle Outflow Tract")
<b>RIMP</b>  Alias: RIMP	(G-0381, SRT, "Right Ventricular Index of Myocardial Performance")	(G-C0E3, SRT, "Finding Site") = (T-32500, SRT, "Right Ventricle") (G-0373, SRT, "Image Mode") = (R-409E4, SRT, "Doppler Pulsed")
<b>S' (Doppler)</b>  Alias: RV S'	(59133-9, LN, "Peak Tissue Velocity")	(R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole")
<b>AFIRV/TAPSE</b> Alias : TAPSE	(GEU-106-0030, 99GEMS, "Tricuspid Annular Plane Systolic Excursion (TAPSE)")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
<b>AFIRV/GPeakSysSL(A4C_RV)</b> Alias : GPeakSysSL(A4C_RV)	(GEU-106-0001, 99GEMS, "Global Peak Longitudinal Strain")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
<b>AFIRV/GPeakSysSL(A4C_RVFW)</b> Alias : GPeakSysSL(A4C_RVFW)	(GEU-106-0153, 99GEMS, "Longitudinal Strain on the	(G-C036, SRT, "Measurement Method") = (GEU,

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

	<b>Freewall of the Right Ventricle"))</b>	<b>GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber"))</b>
<b>AFIRV/Basal FW PeakSysSL Alias : Basal FW PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106- 0147, "Basal segment of the Freewall of the Right Ventricle") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")) (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F- 32020, "Systole")</b>
<b>AFIRV/Mid FW PeakSysSL Alias : Mid FW PeakSysSL</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106- 0148, "Mid segment of the Freewall of the Right Ventricle") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F- 32020, "Systole")
AFIRV/Apical FW PeakSysSL Alias : Apical FW PeakSysSL	(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")	(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106- 0149, "Apical segment of the Freewall of the Right Ventricle") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F- 32020, "Systole")
AFIRV/Apical IVS PeakSysSL Alias : Apical IVS PeakSysSL	(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")	(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106- 0152, "Apical segment of the Intraventricular Septum") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT,

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		"Cardiac Cycle Point") = (SRT, F-32020, "Systole")
AFIRV/Mid IVS PeakSysSL Alias : Mid IVS PeakSysSL	(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")	(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0151, "Mid segment of the Intraventricular Septum") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F-32020, "Systole")
AFIRV/Basal IVS PeakSysSL Alias : Basal IVS PeakSysSL	(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")	(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0150, "Basal segment of the Intraventricular Septum") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>Point") = (SRT, F-32020, "Systole")</b>
<b>AFIRV/GPeakSysSL_Endo(A4C_RV) Alias : GPeakSysSL Endo(A4C_RV)</b>	<b>(GEU-106-0001, 99GEMS, "Global Peak Longitudinal Strain")</b>	<b>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")</b>
<b>AFIRV/GPeakSysSL_Endo(A4C_RVFW) Alias : GPeakSysSL Endo(A4C_RVFW)</b>	<b>(GEU-106-0153, 99GEMS, "Longitudinal Strain on the Freewall of the Right Ventricle")</b>	<b>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")</b>
<b>AFIRV/Basal FW PeakSysSL_Endo Alias : Basal FW PeakSysSL Endo</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0147, "Basal segment of the Freewall of the Right Ventricle") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>Point") = (SRT, F-32020, "Systole")</b>
<b>AFIRV/Mid FW PeakSysSL_Endo Alias : Mid FW PeakSysSL Endo</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0148, "Mid segment of the Freewall of the Right Ventricle") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F-32020, "Systole")</b>
<b>AFIRV/Apical FW PeakSysSL_Endo Alias : Apical FW PeakSysSL Endo</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0149, "Apical segment of the Freewall of the Right Ventricle") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>Point") = (SRT, F-32020, "Systole")</b>
<b>AFIRV/Apical IVS PeakSysSL_Endo Alias : Apical IVS PeakSysSL Endo</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0152, "Apical segment of the Intraventricular Septum") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F-32020, "Systole")</b>
<b>AFIRV/Mid IVS PeakSysSL_Endo Alias : Mid IVS PeakSysSL Endo</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0151, "Mid segment of the Intraventricular Septum") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>Point") = (SRT, F-32020, "Systole")</b>
<b>AFIRV/Basal IVS PeakSysSL_Endo</b> <b>Alias : Basal IVS PeakSysSL Endo</b>	<b>(GEU-106-0002, 99GEMS, "Peak Longitudinal Strain")</b>	<b>(G-C0E3, SRT, "Finding Site") = (GEU, GEU-106-0150, "Basal segment of the Intraventricular Septum") (G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0129, "AFI on endocardium") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (SRT, F-32020, "Systole")</b>
<b>2D/RV FAC General</b> <b>Alias : RV FAC</b>	<b>(79936-1, LN, "Right ventricular Fractional Area Change")</b>	
<b>2D/RV/LV Basal Ratio</b> <b>Alias: RV/LV Basal Ratio</b>	<b>(GEU-106-0180, 99Gems, "RV/LV basal diameter ratio")</b>	<b>(G-0373, SRT, "Image Mode")= (G-03A2, SRT, "2D mode") (R-4089A, SRT, "Cardiac Cycle Point") = (F-32011, SRT, "End Diastole") (G-C0E3, SRT, "Finding Site") = (GEU-106-0182, 99Gems, "LV base level") (111031, DCM,</b>



DIRECTION DOC2652554 REV 3

		"Image View") = (G-A19C, SRT, "Apical four chamber")
--	--	---------------------------------------------------------------

Section Left Atrium

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
<b>2D/Ao/LA</b>  <b>Alias: Ao/LA</b>	(17985-3, LN, "Left Atrium to Aortic Root Ratio")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R- 409E2, SRT, "Doppler Color Flow") depending on scan mode
<b>MM/Ao/LA</b>  <b>Alias: Ao/LA</b>	(17985-3, LN, "Left Atrium to Aortic Root Ratio")	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/LAAo/LA/Ao</b>  <b>Alias: LA/Ao</b>	(17985-3, LN, "Left Atrium to Aortic Root Ratio")	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>MM/LAAo/Ao/LA</b>  <b>Alias: Ao/LA</b>	(17985-3, LN, "Left Atrium to Aortic Root Ratio")	(G-0373, SRT, «Image Mode») = (G-0394, SRT, «M mode»)
<b>2D/LA Major</b>  <b>Alias: LA Major</b>	(G-A193, SRT, "Major Axis")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R- 409E2, SRT, "Doppler Color Flow")

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		depending on scan mode
<p><b>2D/LA Minor</b> <b>Alias: LA Minor</b></p>	<p><b>(G-A194, SRT, “Minor Axis”)</b></p>	<p><b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode</b></p>
<p><b>2D/LA</b> <b>Alias: LA Diam</b></p>	<p><b>(M-02550, SRT, “Diameter”)</b></p>	<p><b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode</b></p>
<p><b>2D/LA D1</b> <b>Alias: LA Diam 1</b></p>	<p><b>(GEU-106-0092, 99GEMS, “Left atrium diameter 1 in LA Ellipsoid Volume”)</b></p>	
<p><b>2D/LA D2</b> <b>Alias: LA Diam 2</b></p>	<p><b>(GEU-106-0093, 99GEMS, “Left atrium diameter 2 in LA Ellipsoid Volume”)</b></p>	
<p><b>2D/LA D3</b> <b>Alias: LA Diam 3</b></p>	<p><b>(GEU-106-0094, 99GEMS, “Left atrium diameter 3 in LA Ellipsoid Volume”)</b></p>	
<p><b>2D/LA SupInf D</b> <b>Alias: LA SupInf D</b></p>	<p><b>(GEU-106-0095, 99GEMS, “Left Atrium</b></p>	

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

	<b>Superior-Inferior Diameter”)</b>	
<b>2D/LA MedLat D Alias: LA MedLat D</b>	<b>(GEU-106-0096, 99GEMS, “Left Atrium Medial-Lateral Diameter”)</b>	
<b>2D/LA AntPost D Alias: LA AntPost D</b>	<b>(29469-4, LN, “Left Atrium Antero-posterior Systolic Dimension”)</b>	
<b>2D/LAA Systole Alias: LAA Systole</b>	<b>(17977-0, LN, “Left Atrium Systolic Area”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)</b>
<b>2D/LAA Diastole Alias: LAA Diastole</b>	<b>(GEU-106-0058, 99GEMS, “Left Atrium Area at diastole”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32010, SRT, “Diastole”)</b>
<b>2D/LA/Ao Alias: LA/Ao</b>	<b>(17985-3, LN, “Left Atrium to Aortic Root Ratio”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		<b>Flow”) depending on scan mode</b>
<b>2D/LA Area</b>  <b>Alias: LA Area</b>	<b>(G-A166, SRT, “Area”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode</b>
<b>2D/LAvol/Ravol</b> <b>Alias: Lavol Ravol Ratio</b>	<b>(59131-3, LN, “Left Atrium volume/Right Atrium volume”)</b>	
<b>2D/LA Volume</b> <b>Alias: LA Volume</b>	<b>(GEU-106- 0097, 99GEMS, “Left Atrium volume by Ellipsoid Method”)</b>	
<b>2D/LA Ellipsoid Volume</b> <b>Alias: LA Ellipsoid Volume</b>	<b>(GEU-106- 0097, 99GEMS, “Left Atrium volume by Ellipsoid Method”)</b>	
<b>LALd(A4C)</b>  <b>Alias: LALd A4C</b>	<b>(29467-8, LN, “Left Atrium Superior- Inferior Dimension, 4- chamber view”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”) = (F- 32010, SRT, “Diastole”)</b>
<b>LAAd(A4C)</b>  <b>Alias: LAAd A4C</b>	<b>(17977-0, LN, “Left Atrium Area A4C view”)</b>	<b>(R-4089A, SRT, “Cardiac Cycle Point”) = (F- 32010, SRT, “Diastole”)</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<p><b>LAEDV(A-L A4C)</b>  <b>Alias: LAEDV A-L A4C</b></p>	<p><b>(122407, DCM, “Left Atrial End Diastolic Volume”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area-Length Single Plane”)</b></p>
<p><b>LAESV(A-L A4C)</b>  <b>Alias: LAESV A-L A4C</b></p>	<p><b>(G-0383, SRT, “Left Atrium Systolic Volume”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area-Length Single Plane”)</b></p>
<p><b>LAEDV(MOD A4C)</b>  <b>Alias: LAEDV MOD A4C</b></p>	<p><b>(122407, DCM, “Left Atrial End Diastolic Volume”)</b></p>	<p><b>(111031, DCM, “Image View”) = (G-A19C, SRT, “Apical four chamber”) (G-C036, SRT, “Measurement Method”) = (125208, DCM, “Method of Disks, Single Plane”)</b></p>
<p><b>LALs(A4C)</b>  <b>Alias: LALs A4C</b></p>	<p><b>(29467-8, LN, “Left Atrium Superior- Inferior Dimension, 4- chamber view”)</b></p>	<p><b>(R-4089A, SRT, “Cardiac Cycle Point”) = (F- 32020, SRT, “Systole”)</b></p>
<p><b>LAAs(A4C)</b>  <b>Alias: LAAs A4C</b></p>	<p><b>(17977-0, LN, “Left Atrium</b></p>	<p><b>(R-4089A, SRT, “Cardiac Cycle Point”) = (F-</b></p>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

	<b>Area A4C view"))</b>	<b>32020, SRT, "Systole")</b>
<b>LAESV(MOD A4C)</b>  <b>Alias: LAESV MOD A4C</b>	<b>(G-0383, SRT, "Left Atrium Systolic Volume")</b>	<b>(111031, DCM, "Image View") = (G-A19C, SRT, "Apical four chamber") (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method of Disks, Single Plane")</b>
<b>LALd(A2C)</b> <b>Alias: LALd A2C</b>	<b>(GEU-106- 0102, 99GEMS, "Left Atrium Length on Apical two chamber view")</b>	<b>(R-4089A, SRT, "Cardiac Cycle Point") = (F- 32011, SRT, "End Diastole") (111031, DCM, "Image View") = (G-A19B, SRT, "Apical two chamber")</b>
<b>LALs(A2C)</b> <b>Alias: LALs A2C</b>	<b>(GEU-106- 0102, 99GEMS, "Left Atrium Length on Apical two chamber view")</b>	<b>(R-4089A, SRT, "Cardiac Cycle Point") = (109070, DCM, "End Systole") (111031, DCM, "Image View") = (G-A19B, SRT, "Apical two chamber")</b>
<b>LAAd(A2C)</b> <b>Alias: LAAd A2C</b>	<b>(GEU-106- 0103, 99GEMS, "Left Atrium Area on Apical two chamger view")</b>	<b>(R-4089A, SRT, "Cardiac Cycle Point") = (F- 32011, SRT, "End Diastole") (111031, DCM, "Image View") = (G-A19B,</b>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		SRT, “Apical two chamber”)
<p><b>LAA(A2C)</b> Alias: LAA A2C</p>	<p>(GEU-106-0103, 99GEMS, “Left Atrium Area on Apical two chamber view”)</p>	<p>(R-4089A, SRT, “Cardiac Cycle Point”) = (109070, DCM, “End Systole”) (111031, DCM, “Image View”) = (G-A19B, SRT, “Apical two chamber”)</p>
<p><b>LAEDV(A-L A2C)</b> Alias: LAEDV A-L A2C</p>	<p>(122407, DCM, “Left Atrial End Diastolic Volume”)</p>	<p>(111031, DCM, “Image View”) = (G-A19B, SRT, “Apical two chamber”) (G-C036, SRT, “Measurement Method”) = (125205, DCM, “Area-Length Single Plane”)</p>
<p><b>LAEDV(MOD A2C)</b> Alias: LAEDV MOD A2C</p>	<p>(122407, DCM, “Left Atrial End Diastolic Volume”)</p>	<p>(111031, DCM, “Image View”) = (G-A19B, SRT, “Apical two chamber”) (G-C036, SRT, “Measurement Method”) = (125208, DCM, “Method of Disks, Single Plane”)</p>
<p><b>LAESV(A-L A2C)</b> Alias: LAESV A-L A2C</p>	<p>(G-0383, SRT, “Left Atrium Systolic Volume”)</p>	<p>(111031, DCM, “Image View”) = (G-A19B, SRT, “Apical two chamber”) (G-C036, SRT, “Measurement Method”) =</p>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		(125205, DCM, “Area-Length Single Plane”)
<p><b>LAESV(MOD A2C)</b></p> <p><b>Alias: LAESV MOD A2C</b></p>	(G-0383, SRT, “Left Atrium Systolic Volume”)	<p>(111031, DCM, “Image View”)</p> <p>= (G-A19B, SRT, “Apical two chamber”)</p> <p>(G-C036, SRT, “Measurement Method”) =</p> <p>(125208, DCM, “Method of Disks, Single Plane”)</p>
<p><b>MM/LA/Ao</b></p> <p><b>Alias: LA/Ao</b></p>	(17985-3, LN, “Left Atrium to Aortic Root Ratio”)	<p>(G-0373, SRT, «Image Mode»)</p> <p>= (G-0394, SRT, «M mode»)</p>
<p><b>MM/LA</b></p> <p><b>Alias: LA Diam</b></p>	(29469-4, LN, “Left Atrium Antero-posterior Systolic Dimension”)	<p>(G-0373, SRT, «Image Mode»)</p> <p>= (G-0394, SRT, «M mode»)</p>
<p><b>2D/LAEDV(A-L)</b></p> <p><b>Alias: LAEDV(A-L)</b></p>	(122407, DCM, “Left Atrial End Diastolic Volume”)	<p>(G-0373, SRT, “Image Mode”)</p> <p>= (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”)</p> <p>depending on scan mode</p> <p>(G-C036, SRT, “Measurement Method”) =</p> <p>(125204, DCM, “Area-Length Biplane”)</p>



**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<p><b>2D/LAEDVI(A-L)</b></p> <p><b>Alias: LAEDV Index (A-L)</b></p>	<p><b>(GEU-106-0027, 99GEMS, “Left Atrial End Diastolic Volume Index”)</b></p>	<p><b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125204, DCM, “Area-Length Biplane”)</b></p>
<p><b>2D/LAESV(A-L)</b></p> <p><b>Alias: LAESV(A-L)</b></p>	<p><b>(G-0383, SRT, “Left Atrium Systolic Volume”)</b></p>	<p><b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125204, DCM, “Area-Length Biplane”)</b></p>
<p><b>LAEDV(MOD BP)</b></p> <p><b>Alias: LAEDV(MOD BP)</b></p>	<p><b>(122407, DCM, “Left Atrial End Diastolic Volume”)</b></p>	<p><b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement</b></p>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		Method”) = (125207, DCM, “Method of Disks, Biplane”)
<b>LAESV(MOD BP)</b>  <b>Alias: LAESV(MOD BP)</b>	<b>(G-0383, SRT, “Left Atrium Systolic Volume”)</b>	<b>(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R- 409E2, SRT, “Doppler Color Flow”) depending on scan mode (G-C036, SRT, “Measurement Method”) = (125207, DCM, “Method of Disks, Biplane”)</b>
<b>SD/Laappendix Vmax</b> <b>Alias: Laappendix Vmax</b>	<b>(29486-8, LN, “Left Atrial Appendage Peak Velocity”)</b>	
<b>AFILA/2DLA_ReservoirStrain_R_Wave(A4C)</b> <b>Alias : 2DLA_ReservoirStrain_R_Wave(A4C)</b>	<b>(GEU-106- 0161, 99GEMS, “Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave”)</b>	<b>(G-C036, SRT, “Measurement Method”) = (GEU, GEU- 106-0018, “AFI”) (111031, DCM, “Image View”) = (SRT, G- A19C, “Apical four chamber”) (R-4089A, SRT, “Cardiac Cycle Point”) = (GEU, GEU-106-0158, “Reservoir phase of the left atrium”)</b>

DIRECTION DOC2652554 REV 3

<p>AFILA/2DLA_ConduitStrain_R_Wave(A4C) Alias : 2DLA_ConduitStrain_R_Wave(A4C)</p>	<p>(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")</p>
<p>AFILA/2DLA_ContractileStrain_R_Wave(A4C) Alias : 2DLA_ContractileStrain_R_Wave(A4C)</p>	<p>(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")</p>
<p>AFILA/2DLA_ReservoirStrain_P_Wave(A4C) Alias : 2DLA_ReservoirStrain_P_Wave(A4C)</p>	<p>(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View")</p>

DIRECTION DOC2652554 REV 3

		= (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_P_Wave(A4C) Alias : 2DLA_ConduitStrain_P_Wave(A4C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_P_Wave(A4C) Alias : 2DLA_ContractileStrain_P_Wave(A4C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		phase of the left atrium")
AFILA/2DLA_ReservoirStrain_R_Wave(A2C) Alias : 2DLA_ReservoirStrain_R_Wave(A2C)	(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_R_Wave(A2C) Alias : 2DLA_ConduitStrain_R_Wave(A2C)	(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_R_Wave(A2C) Alias : 2DLA_ContractileStrain_R_Wave(A2C)	(GEU-106-0161, 99GEMS, "Longitudinal Strain of the Left Atrium	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI")

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

	with 0 strain at ECG R-wave"))	(111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium"))
AFILA/2DLA_ReservoirStrain_P_Wave(A2C) Alias : 2DLA_ReservoirStrain_P_Wave(A2C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave"))	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium"))
AFILA/2DLA_ConduitStrain_P_Wave(A2C) Alias : 2DLA_ConduitStrain_P_Wave(A2C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave"))	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU,

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_P_Wave(A2C) Alias : 2DLA_ContractileStrain_P_Wave(A2C)	(GEU-106-0162, 99GEMS, "Longitudinal Strain of the Left Atrium with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")
AFILA/2DLA_ReservoirStrain_R_Wave(BiP) Alias : 2DLA_ReservoirStrain_R_Wave(BiP)	(GEU-106-0163, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_R_Wave(BiP) Alias : 2DLA_ConduitStrain_R_Wave(BiP)	(GEU-106-0163, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU,

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		GEU-106-0159, "Conduit phase of the left atrium")
AFILA/2DLA_ContractileStrain_R_Wave(BiP) Alias : 2DLA_ContractileStrain_R_Wave(BiP)	(GEU-106-0163, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG R-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")
AFILA/2DLA_ReservoirStrain_P_Wave(BiP) Alias : 2DLA_ReservoirStrain_P_Wave(BiP)	(GEU-106-0164, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0158, "Reservoir phase of the left atrium")
AFILA/2DLA_ConduitStrain_P_Wave(BiP) Alias : 2DLA_ConduitStrain_P_Wave(BiP)	(GEU-106-0164, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG P-wave")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0159, "Conduit phase of the left atrium")



**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

<p>AFILA/2DLA_ContractileStrain_P_Wave(BiP) Alias : 2DLA_ContractileStrain_P_Wave(BiP)</p>	<p>(GEU-106-0164, 99GEMS, "Longitudinal Strain of the Left Atrium from biplane measurements with 0 strain at ECG P-wave")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (R-4089A, SRT, "Cardiac Cycle Point") = (GEU, GEU-106-0160, "Contractile phase of the left atrium")</p>
<p>AFILA/2DLA_EF(A2C) Alias : 2DLA_EF(A2C)</p>	<p>(GEU-106-0165, 99GEMS, "Left Atrium Emptying Fraction by speckle tracking")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")</p>
<p>AFILA/2DLA_EV(A2C) Alias : 2DLA_EV(A2C)</p>	<p>(GEU-106-0166, 99GEMS, "Left Atrium Emptying Volume by speckle tracking")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")</p>
<p>AFILA/2DLA_Vmax(A2C) Alias : 2DLA_Vmax(A2C)</p>	<p>(GEU-106-0167, 99GEMS, "Left Atrium maximal volume by speckle tracking")</p>	<p>(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-</p>

**GE HEALTHCARE  
STATEMENT**

**VIVID AND ECHOPAC V206  
CONFORMANCE**

**DIRECTION DOC2652554 REV 3**

		A19B, "Apical two chamber")
AFILA/2DLA_Vmin(A2C) Alias : 2DLA_Vmin(A2C)	(GEU-106-0168, 99GEMS, "Left Atrium minimal volume by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")
AFILA/2DLA_VpreA(A2C) Alias : 2DLA_VpreA(A2C)	(GEU-106-0169, 99GEMS, "Left Atrium volume at preA time by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19B, "Apical two chamber")
AFILA/2DLA_EF(A4C) Alias : 2DLA_EF(A4C)	(GEU-106-0165, 99GEMS, "Left Atrium Emptying Fraction by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View") = (SRT, G-A19C, "Apical four chamber")
AFILA/2DLA_EV(A4C) Alias : 2DLA_EV(A4C)	(GEU-106-0166, 99GEMS, "Left Atrium Emptying Volume by speckle tracking")	(G-C036, SRT, "Measurement Method") = (GEU, GEU-106-0018, "AFI") (111031, DCM, "Image View")