

DIRECTION DOC2652554 REV 3

2D/IVC Alias: IVC	(18006-7, LN, “Inferior Vena Cava Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/SVC Diam Ins Alias: SVC Diam Ins	(18007-5, LN, “Superior Vena Cava Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-40899, SRT, “Respiratory Cycle Point”) = (F-20010, SRT, “During Inspiration”)
2D/SVC DIAM Exp Alias: AVC Diam Exp	(18007-5, LN, “Superior Vena Cava Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (R-40899, SRT, “Respiratory Cycle Point”) = (F-20020, SRT, “During Expiration”)
2D/IVC Collapse Index Alias: IVC Collapse Index	(18050-5, LN, “Inferior Vena Cava % Collapse”)	
2D/SVC Collapse Index Alias: SVC Collapse Index	(GEU-106-0133, 99GEMS, “Superior Vena Cava % Collapse”)	

Section Cardiac Shunt Study

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
Qp/Qs Alias: Qp/Qs	(29462-9, LN, “Pulmonary-to-Systemic Shunt Flow Ratio”)	
Systemic VTI Alias: Systemic VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (F-32330, SRT, “Left to right cardiovascular shunt”)
Pulmonic VTI Alias: Pulmonic VTI	(20354-7, LN, “Velocity Time Integral”)	(G-C048, SRT, “Direction of Flow”) = (F-32340, SRT, “Right to left cardiovascular shunt”)

Section Congenital Anomaly of Cardiovascular System

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/ASD Diam Alias: ASD Diam	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C0E3, SRT, “Finding Site”) = (D4-31220, SRT, “Atrial Septal Defect”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/VSD Diam Alias: VSD Diam	(G-038F, SRT, “Cardiovascular Orifice Diameter”)	(G-C0E3, SRT, “Finding Site”) = (D4-31150, SRT, “Ventricular Septal Defect”) (G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) or (R-409E2, SRT, “Doppler Color Flow”) depending on scan mode
2D/Pre Ductal Alias: Pre Ductal	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (G-C0E3, SRT, “Finding Site”) = (T-42340, SRT, “Preductal region of aortic arch”)
2D/Post Ductal Alias: Post Ductal	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (G-C0E3, SRT, “Finding Site”) = (T-42350, SRT, “Postductal region of aortic arch”)
2D/Systemic Diam Alias: Systemic Diam	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

2D/Pulmonic Diam Alias: Pulmonic Diam	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”) (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
VSD Vmax Alias: VSD Vmax	(11726-7, LN, “Peak Velocity”)	(G-C0E3, SRT, “Finding Site”) = (D4-31150, SRT, “Ventricular Septal Defect”)
VSD maxPG Alias: VSD maxPG	(20247-3, LN, “Peak Gradient”)	(G-C0E3, SRT, “Finding Site”) = (D4-31150, SRT, “Ventricular Septal Defect”)
ASD Vmax Alias: ASD Vmax	(11726-7, LN, “Peak Velocity”)	(G-C0E3, SRT, “Finding Site”) = (D4-31220, SRT, “Atrial Septal Defect”)
ASD maxPG Alias: ASD maxPG	(20247-3, LN, “Peak Gradient”)	(G-C0E3, SRT, “Finding Site”) = (D4-31220, SRT, “Atrial Septal Defect”)
Systemic HR Alias: Systemic HR	(8867-4, LN, “Heart rate”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic Vmax Alias: Systemic Vmax	(11726-7, LN, “Peak Velocity on systemic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic Vmean Alias: Systemic Vmean	(20352-1, LN, “Mean Velocity on systemic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic maxPG Alias: Systemic maxPG	(20247-3, LN, “MPeak Gradient on systemic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		(G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic meanPG Alias: Systemic meanPG	(20256-4, LN, “Mean Gradient on systemic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic SV Alias: Systemic SV	(F-32120, SRT, “Stroke Volume on systemic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic CO Alias: Systemic CO	(F-32100, SRT, “Cardiac Output on systemic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Systemic Env.Ti Alias: Systemic Env.Ti	(GEU-106-0090, 99GEMS, “Time duration of the VTI trace on Systemic side flow”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32550, SRT, “Right Ventricle Outflow Tract”)
Pulmonic HR Alias: Pulmonic HR	(8867-4, LN, “Heart rate”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic Vmax Alias: Pulmonic Vmax	(11726-7, LN, “Peak Velocity on pulmonic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed»)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		(G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic Vmean Alias: Pulmonic Vmean	(20352-1, LN, “Mean Velocity on pulmonic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic maxPG Alias: Pulmonic maxPG	(20247-3, LN, “Peak Gradient on pulmonic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic meanPG Alias: Pulmonic meanPG	(20256-4, LN, “Mean Gradient on pulmonic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic SV Alias: Pulmonic SV	(F-32120, SRT, “Stroke Volume on pulmonic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic CO Alias: Pulmonic CO	(F32100, SRT, “Cardiac Output on pulmonic side”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Pulmonic Env.Ti Alias: Pulmonic Env.Ti	(GEU-106-0091, 99GEMS, “Time duration of the VTI trace on Pulmonic side flow”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (G-C0E3, SRT, “Finding Site”) = (T-32650, SRT, “Left Ventricle Outflow Tract”)
Coarctation/Post-Ductal	(17995-2, LN, “Thoracic Aorta Coarctation Systolic	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

Alias: Coarctation Post Ductal	Peak Instantaneous Gradient")	
Coarctation/Post-Ductal PG Alias: Coarctation Post Ductal PG	(20256-4, LN, "Mean Gradient")	(G-C0E3, SRT, "Finding Site") = (D4-32030, SRT, "Thoracic Aortic Coarctation")
Coarctation/Pre-Ductal Alias: Pre-Ductal	(GEU-106-0107, 99GEMS, "Peak Velocity in the Pre-Ductal area of the Aortic Coarctation")	(G-C0E3, SRT, "Finding Site") = (D4-32030, SRT, "Thoracic Aortic Coarctation")
Coarctation/Pre-Ductal PG Alias: Pre-Ductal PG	(GEU-106-0108, 99GEMS, "Maximum Pressure Gradient in the Pre-Ductal area of the Aortic Coarctation")	(G-C0E3, SRT, "Finding Site") = (D4-32030, SRT, "Thoracic Aortic Coarctation")
ASD Vmean Alias: ASD Vmean	(20352-1, LN, "Mean Velocity")	(G-C0E3, SRT, "Finding Site") = (D4-31220, SRT, "Atrial Septal Defect")
ASD meanPG Alias: ASD meanPG	(20256-4, LN, "Mean Gradient")	(G-C0E3, SRT, "Finding Site") = (D4-31220, SRT, "Atrial Septal Defect")
ASD Env. Ti Alias: ASD Env. Ti	(GEU-106-0132, 99GEMS, "Time duration of the VTI trace")	(G-C0E3, SRT, "Finding Site") = (D4-31220, SRT, "Atrial Septal Defect")
ASD VTI Alias: ASD VTI	(20354-7, LN, "Velocity Time Integral")	(G-C0E3, SRT, "Finding Site") = (D4-31220, SRT, "Atrial Septal Defect")
VSD Vmean Alias: VSD Vmean	(20352-1, LN, "Mean Velocity")	(G-C0E3, SRT, "Finding Site") = (D4-31150, SRT, "Ventricular Septal Defect")
VSD meanPG Alias: VSD meanPG	(20256-4, LN, "Mean Gradient")	(G-C0E3, SRT, "Finding Site") = (D4-31150, SRT, "Ventricular Septal Defect")
VSD Env. Ti Alias: VSD Env. Ti	(GEU-106-0132, 99GEMS, "Time duration of the VTI trace")	(G-C0E3, SRT, "Finding Site") = (D4-31150, SRT, "Ventricular Septal Defect")
VSD VTI Alias: VSD VTI	(20354-7, LN, "Velocity Time Integral")	(G-C0E3, SRT, "Finding Site") = (D4-31150, SRT, "Ventricular Septal Defect")

Section Pericardial cavity

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/Pes Alias: Pes	(121206, DCM, "Distance")	(G-C0E3, SRT, "Finding Site") = (D3-90008, SRT, "Pericardial effusion") (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole")
2D/Ped Alias: Ped	(121206, DCM, "Distance")	(G-C0E3, SRT, "Finding Site") = (D3-90008, SRT, "Pericardial effusion") (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") or (R-409E2, SRT, "Doppler Color Flow") depending on scan mode (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
MM/Ped Alias: Ped	(121206, DCM, "Distance")	(G-C0E3, SRT, "Finding Site") = (D3-90008, SRT, "Pericardial effusion") (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")

Section Aortic Sinotubular Junction

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/Ao st junct Alias: Ao st junct	(M-02550, SRT, "Diameter")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode")
2D/Ao st junct/Ao Alias: Ao st junct/Ao	(59116-4, LN, "Aortic sinotubular junction diameter/Aortic root diameter by US")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode")

DIRECTION DOC2652554 REV 3

2D/Ao Diam Stub Alias: Ao Diam Stub	(GEU-106-0068, 99GEMS, “Aortic diameter at sinotubular transition”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”)
--	--	---

Section Sinus Valsalva

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/Ao Diam Svals Alias: Ao Diam Svals	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”)
2D/SinusesOfValsalva Alias: Sinuses of Val	(M-02550, SRT, “Diameter”)	

Section Patent Ductus Arteriosus

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/PDA Diam Alias: PDA Diam	(M-02550, SRT, “Diameter”)	(G-0373, SRT, “Image Mode”) = (G-03A2, SRT, “2D mode”)
PDA Systolic Alias: PDA Systolic	(11726-7, LN, “Peak Velocity”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
PDA Systolic PG Alias: PDA Systolic PG	(20247-3, LN, “Peak Gradient”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
PDA Systolic Vmean Alias: PDA Systolic Vmean	(20352-1, LN, “Mean Velocity”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

PDA Systolic meanPG Alias: PDA Systolic meanPG	(20256-4, LN, “Mean Gradient”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
PDA Systolic Env. Ti Alias: PDA Systolic Env. Ti	(GEU-106-0132, 99GEMS, “Time duration of the VTI trace”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
PDA Systolic VTI Alias: PDA Systolic VTI	(20354-7, LN, “Velocity Time Integral”)	(R-4089A, SRT, “Cardiac Cycle Point”) = (F-32020, SRT, “Systole”)
PDA Diastolic Alias: PDA Diastolic	(11726-7, LN, “Peak Velocity”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
Alias: PDA Diastolic PG	(20247-3, LN, “Peak Gradient”)	(G-0373, SRT, «Image Mode») = (R-409E4, SRT, «Doppler Pulsed») (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
PDA Diastolic Vmean Alias: PDA Diastolic Vmean	(20352-1, LN, “Mean Velocity”)	(R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
PDA Diastolic meanPG Alias: PDA Diastolic meanPG	(20256-4, LN, “Mean Gradient”)	(R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
PDA Diastolic Env. Ti Alias: PDA Diastolic Env. Ti	(GEU-106-0132, 99GEMS, “Time duration of the VTI trace”)	(R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")
PDA Diastolic VTI Alias: PDA Diastolic VTI	(20354-7, LN, “Velocity Time Integral”)	(R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole")

Section Patent Foramen Ovale

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
---	--------------------------------------	--

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

2D/PEs Alias: PEs	(M-02550, SRT, "Diameter")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode")
--	---------------------------------------	---

Section Coronary Artery

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/LCA Alias: LCA	(M-02550, SRT, "Diameter")	(G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") (G-C0E3, SRT, "Finding Site") = (T-43107, SRT, "Left Main Coronary Artery")
2D/RCA Alias: RCA	(M-02550, SRT, "Diameter")	(G-C0E3, SRT, "Finding Site") = (T-43203, SRT, "Right Coronary Artery")
2D/LCX Alias: LCX	(M-02550, SRT, "Diameter")	(G-C0E3, SRT, "Finding Site") = (T-43120, SRT, "Circumflex Coronary Artery")

Section Mitral Valve (prosthetics)

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
MP VTI Alias: MP VTI	(20354-7, LN, "Velocity Time Integral")	(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")
MP Area Alias: MP Area	(G-038E, SRT, "Cardiovascular Orifice Area")	(125215, DCM, "Continuity Equation by Velocity Time Integral")

Section Aortic Valve (prosthetics)

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
---	--	--

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

AP VTI Alias: AP VTI	(20354-7, LN, "Velocity Time Integral")	(G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow")
AP Area Alias: AP Area	(G-038E, SRT, "Cardiovascular Orifice Area")	(125215, DCM, "Continuity Equation by Velocity Time Integral")

Section Aortic Arch

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/ProxAoArch Alias : ProxAoArch	(18011-7, LN, "Aortic Arch Diameter")	(G-C0E3, SRT, "Finding Site") = (SRT, G-A118, "Proximal")

Section LCA Descending Branch

GEU Parameter ID (and corresponding alias)	Base Measurement Concept Name	Concept or Acquisition Context Modifier
2D/LAD Alias: LAD	(M-02550, SRT, "Diameter")	

DIRECTION DOC2652554 REV 3

15. VASCULAR ULTRASOUND PROCEDURE REPORT

This section describes the contents of the Vascular Ultrasound Procedure Report (TID 5100) SR.

Note: If “Use older SR version” is enabled (see 2.6, 3.6 and 7) the corresponding section present in the DICOM Conformance Statement of the selected version should be used.

15.1 USAGE AND EXTENSION OF TID 5100 VASCULAR ULTRASOUND REPORT

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	EV (125100, DCM, “Vascular Ultrasound Procedure Report”)	1	M		
	>	HAS OBS CONTEXT	INCLUDE	DTID (1001) Observation Context	1	M		
	>	CONTAINS	INCLUDE	DTID (5101) Vascular Patient Characteristics	1	U		
	>	CONTAINS	INCLUDE	DTID (5102) Vascular Procedure Summary Section	1	U		
	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1-n	U		

15.2 TID 5101 VASCULAR PATIENT CHARACTERISTICS

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	EV (121118, DCM, “Patient Characteristics”)	1	M		
	>	CONTAINS	NUM	EV (121033, DCM, “Subject Age”)	1	U		Units = DCID (7456) Units of Measure for Age
	>	CONTAINS	CODE	EV (121032, DCM, “Subject Sex”)	1	U		DCID (7455) Sex
	>	CONTAINS	NUM	EV (8867-4, LN, “Heart Rate”)	1	U		
	>	CONTAINS	NUM	EV (F-008EC, SRT, “Systolic Blood Pressure”)	1	U		
	>	CONTAINS	NUM	EV (F-008ED, SRT, “Diastolic Blood Pressure”)	1	U		

15.3 TID 5102 VASCULAR PROCEDURE SUMMARY SECTION

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (121111, DCM, “Summary”)	1	M		

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	>	CONTAINS	TEXT	EV (121106, DCM, “Comment”)	1	M		
--	---	----------	------	--------------------------------	---	---	--	--

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

15.4 TID 5103 VASCULAR ULTRASOUND SECTION (EXTENDED)

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
			CONTAINER	DT (121070, DCM, "Findings")	1	M		
	>	HAS CONCEPT MOD	CODE	EV (G-C0E3, SRT, "Finding Site")	1	M		See 15.6. GEU Applications and Extensions - \$SectionScope
	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, "Laterality")	1	U		See 15.9 GE Ultrasound Sidedness and Vessel Location
	>	HAS CONCEPT MOD	CODE	EV (G-0373, SRT, "Image Mode")	1	M		See 15.8 GE Ultrasound modes.
	>	CONTAINS	INCLUDE	DTID (5104) Vascular Measurement Group	1-n	M		See 15.5 TID 5104 Vascular Ultrasound Measurement Group
	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Measurement = \$AnatomyRatio

* This template is extended with the Image Mode row.

15.5 TID 5104 VASCULAR ULTRASOUND MEASUREMENT GROUP

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$Anatomy GEU Parameters	1	M		See 15.6 GEU Applications and Extensions – Anatomy GEU Parameter
	>	HAS CONCEPT MOD	CODE	EV (G-A1F8, SRT, "Topographical Modifier")	1	U		See 15.9 GE Ultrasound Sidedness and Vessel Location
	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Measurement = See 15.10 Vascular Base Measurement \$Derivation = DCID (3626) Measurement Type

GE HEALTHCARE
STATEMENT

DIRECTION DOC2652554 REV 3

15.6 GEU APPLICATIONS AND EXTENSIONS

Section Scope	Section Laterality	Anatomy	Anatomy Ratio	GEU Parameters Base Measurement Concept Name	
DT (121070, DCM, “Findings”)	EV (G-C171, SRT, “Laterality”)				
(T-40501, SRT, “Blood Vessel of Head”)	(G-A101, SRT, “Left”) for Left, (G-A100, SRT, “Right”) for Right. or (G-A103, SRT, “Unilateral”)	DCID 12105 Intracranial Cerebral Vessels or DCID 12106 Intracranial Cerebral Vessels (Unilateral)		Anatomy GEU parameter	Code and Description
				ICA	(T-45300, SRT, “Internal Carotid Artery”)
				MCA	G (T-45600, SRT, “Middle Cerebral Artery”)
				ACA	(T-45540, SRT, “Anterior Cerebral Artery”)
				PCA	(T-45900, SRT, “Posterior Cerebral Artery”)
				PCoMA	(T-45320, SRT, “Posterior Communicating Artery”)
				ACoMA	(T-45530, SRT, “Anterior Communicating Artery”)
				VERT	(T-45700, SRT, “Vertebral Artery”)
				BA	(T-45800, SRT, “Basilar Artery”)
				TABLE 15.6.1 TCD Study Folder Code Maps	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

(T-45005, SRT, “Artery of neck)	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right.	DCID 12104 Extracranial Arteries	DCID 12123 Carotid Ratios	Anatomy GEU parameter	Code and Description
				VERT	(T-45700, SRT, “Vertebral Artery”)
				CCA	(T-45100, SRT, “Common Carotid Artery”)
				ICA, ICA1	(T-45300, SRT, “Internal Carotid Artery”)
				Innominate	(T-46010, SRT, ‘Innominate Artery’)
				BULB	(T-45170, SRT, “Carotid Bulb”)
				ECA	(T-45200, SRT, “External Carotid Artery”)
				SUBC	(T-46100, SRT, “Subclavian Artery”)
				BIF	(SRT, T-45160, “Carotid Bifurcation”)
				Stent	(A-25500, SRT, 'Stent')
				Pre-Stent	(GEU-1004-71, 99GEMS, 'Pre-Stent')
				Post-Stent	(GEU-1004-72, 99GEMS, 'Post-Stent')
				TABLE 15.6.2 Carotid Study Folder Code Maps	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

(T-47040, SRT, “Artery of Lower Extremity”)	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right. or (G-A103, SRT, “Unilateral”)	DCID 12109 Lower Extremity Arteries or DCID 12112 Abdominal Arteries (unilateral)		Anatomy GEU parameter	Code and Description
				ComIliac	(T-46710, SRT, “Common Iliac Artery”)
				ExtIliac (EIA)	(T-46910, SRT, “External Iliac Artery”)
				ComFemoral (CFA)	(T-47400, SRT, “Common Femoral Artery”)
				SupFemoral (SFA)	(T-47403, SRT, “Superficial Femoral Artery”)
				Popliteal (Pop A)	(T-47500, SRT, “Popliteal Artery”)
				AntTibial (ATA)	(T-47700, SRT, “Anterior Tibial Artery”)
				PostTibial (PTA)	T-47600, SRT, “Posterior Tibial Artery”)
				Peroneal (Peron A)	(T-47630, SRT, “Peroneal Artery”)
				DorsPedis (DPA)	(T-47741, SRT, “Dorsalis Pedis Artery”)
				DeepFemoral (DFA)	(T-47440, SRT, “Profunda Femoris Artery”)
				Profunda (Pro)	(T-47440, SRT, “Profunda Femoris Artery”)
				Aorta	(T-4200, SRT, ”Aorta”)
				Stent	(A-25500, SRT, 'Stent')
				Pre-Stent	(GEU-1004-71, 99GEMS, 'Pre-Stent')
				Post-Stent	(GEU-1004-72, 99GEMS, 'Post-Stent')
TABLE 15.6.3 LEA Study Folder Code Maps					

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

(T-49403, SRT, “Vein of Lower Extremity”)	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right. Or (G-A103, SRT, “Unilateral”)	DCID 12110 Lower Extremity of Veins or DCID 12114 Abdominal Veins (unilateral)		Anatomy GEU parameter	Code and Description
				Popliteal	(T-49640, SRT, “Popliteal Vein”)
				LSaphenous	(T-49550, SRT, “Lesser Saphenous Vein”)
				AntTibial	(T-49630, SRT, “Anterior Tibial Vein”)
				PostTibial	(T-49620, SRT, “Posterior Tibial Vein”)
				Peroneal	(T-49650, SRT, “Peroneal Vein”)
				Profunda	(T-49660, SRT, “Profunda Femoris Vein”)
				ExtIliac	(T-48930, SRT, “External Iliac Vein”)
				ComFemoral	(G-035B, SRT, “Common Femoral Vein”)
				ComIliac	(T-48920, SRT, “Common Iliac Vein”)
				Great saphenous	(T-49530, SRT, “Great Saphenous Vein”)
				Femoral	(G-035B, SRT, “Femoral Vein”)
				IVC	(T-48710, SRT, “Inferior Vena Cava”)
				DeepFemoral	(T-49660, SRT, “Profunda Femoris Vein”)
				Profunda	(T-49660, SRT, “Profunda Femoris Vein”)
				SaphFemJunc	(T-D930A, SRT, ‘Saphenofemoral Junction’)
				GreatSaphCalf	(R-1025A, SRT, ‘Great Saphenous Vein of Calf’)
				GreatSaphAccess	(GEU-1004-73, 99GEMS, ‘Great Saphenous Vein of Accessory’)
				Perforator	(GEU-1005-6, 99GEMS, ‘User Vessel Anatomy’)
				SaphPopJunc	(T-4941A, SRT, ‘Saphenopopliteal junction’)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

(T-49403, SRT, “Vein of Lower Extremity”) Continued	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right. Or (G-A103, SRT, “Unilateral”)	DCID 12110 Lower Extremity of Veins or DCID 12114 Abdominal Veins (unilateral)		Anatomy GEU parameter	Code and Description
				FemoralPopJunc	(GEU-1005-6, 99GEMS, ‘User Vessel Anatomy’)
				PopTibialJunc	(GEU-1005-6, 99GEMS, ‘User Vessel Anatomy’)
				VaricoseVein	(GEU-1005-6, 99GEMS, ‘User Vessel Anatomy’)
				AntAccessSaphV	(GEU-1005-6, 99GEMS, ‘User Vessel Anatomy’)
				PostAccessSaphV	(GEU-1005-6, 99GEMS, ‘User Vessel Anatomy’)
				GreatSaphCalf	(R-10259, SRT, 'Great Saphenous Vein of Thigh')
				Pseudo	(M-32390, SRT, 'Pseudo Aneurysm')
				TABLE 15.6.4 LEV Study Folder Code Maps	
(T-47020, SRT, “Artery of Upper Extremity”)	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right.	DCID (12107) Upper Extremity Arteries		Anatomy GEU parameter	Code and Description
				SUBC	(T-46100, SRT, “Subclavian artery”)
				Axill	(T-47100, SRT, “Axillary artery”)
				BrachialA	(T-47160, SRT, “Brachial artery”)
				RadialA	(T-47300, SRT, “Radial artery”)
				UlnarA	(T-47200, SRT, “Ulnar artery”)
				Palmar	(T-47340, SRT, “Deep Palmar Arch of Radial Artery”)
				Innominate	(T-46010, SRT, “Innominate Artery”)
				Pseudo	(M-32390, SRT, 'Pseudo Aneurysm')
				Stent	(A-25500, SRT, 'Stent')
				Pre-Stent	(GEU-1004-71, 99GEMS, 'Pre-Stent')
				Post-Stent	(GEU-1004-72, 99GEMS, 'Post-Stent')
				TABLE 15.6.5 UEA Study Folder Code Maps	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

(T-49103, SRT, “Vein of Upper Extremity”)	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right.	DCID 12108 Upper Extremity Veins		<table><tr><th>Anatomy GEU parameter</th><th>Code and Description</th></tr><tr><td>JugularV</td><td>(T-48170, SRT, “Internal Jugular vein”)</td></tr><tr><td>InnoV</td><td>(T-48620, SRT, “Innominate vein”)</td></tr><tr><td>SUBCV</td><td>(T-48330, SRT, “Subclavian vein”)</td></tr><tr><td>AxillV</td><td>(T-49110, SRT, “Axillary vein”)</td></tr><tr><td>CephV</td><td>(T-49240, SRT, “Cephalic vein”)</td></tr><tr><td>BasilV</td><td>(T-48052, SRT, “Basilic vein”)</td></tr><tr><td>BracV</td><td>(T-49350, SRT, “Brachial vein”)</td></tr><tr><td>McubV</td><td>(T-49250, SRT, “Median Cubital vein”)</td></tr><tr><td>RadialV</td><td>(T-49340, SRT, “Radial vein”)</td></tr><tr><td>UlnarV</td><td>(T-49330, SRT, “Ulnar vein”)</td></tr><tr><td>Pseudo</td><td>(M-32390, SRT, 'Pseudo Aneurysm')</td></tr><tr><td>AVF</td><td>(M-39390, SRT, 'AV Fistula')</td></tr><tr><td>Axill</td><td>(T-49110, SRT, “Axillary vein”)</td></tr></table>	Anatomy GEU parameter	Code and Description	JugularV	(T-48170, SRT, “Internal Jugular vein”)	InnoV	(T-48620, SRT, “Innominate vein”)	SUBCV	(T-48330, SRT, “Subclavian vein”)	AxillV	(T-49110, SRT, “Axillary vein”)	CephV	(T-49240, SRT, “Cephalic vein”)	BasilV	(T-48052, SRT, “Basilic vein”)	BracV	(T-49350, SRT, “Brachial vein”)	McubV	(T-49250, SRT, “Median Cubital vein”)	RadialV	(T-49340, SRT, “Radial vein”)	UlnarV	(T-49330, SRT, “Ulnar vein”)	Pseudo	(M-32390, SRT, 'Pseudo Aneurysm')	AVF	(M-39390, SRT, 'AV Fistula')	Axill	(T-49110, SRT, “Axillary vein”)
Anatomy GEU parameter	Code and Description																															
JugularV	(T-48170, SRT, “Internal Jugular vein”)																															
InnoV	(T-48620, SRT, “Innominate vein”)																															
SUBCV	(T-48330, SRT, “Subclavian vein”)																															
AxillV	(T-49110, SRT, “Axillary vein”)																															
CephV	(T-49240, SRT, “Cephalic vein”)																															
BasilV	(T-48052, SRT, “Basilic vein”)																															
BracV	(T-49350, SRT, “Brachial vein”)																															
McubV	(T-49250, SRT, “Median Cubital vein”)																															
RadialV	(T-49340, SRT, “Radial vein”)																															
UlnarV	(T-49330, SRT, “Ulnar vein”)																															
Pseudo	(M-32390, SRT, 'Pseudo Aneurysm')																															
AVF	(M-39390, SRT, 'AV Fistula')																															
Axill	(T-49110, SRT, “Axillary vein”)																															

TABLE 15.6.6 UEV Study Folder Code Maps

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

(T-71019, SRT, “Vascular Structure of Kidney”)	(G-A101, SRT, “Left”) for Left, or (G-A100, SRT, “Right”) for Right.	DCID 12115 Renal Vessels Or SRT T-71000 Kidney Or SRT T74000 Bladder	DCID 12124 Renal Ratios	<table><tr><th>Anatomy GEU parameter</th><th>Code and Description</th></tr><tr><td>MRenalA</td><td>(T-46600, SRT, “Renal Artery”)</td></tr><tr><td>RenalV</td><td>(T-48740, SRT, “Renal Vein”)</td></tr><tr><td>SegmentalA</td><td>(T-46659, SRT, “Segmental Artery”)</td></tr><tr><td>InterlobarA</td><td>(T-4667D, SRT, “Interlobar Artery of Kidney”)</td></tr><tr><td>ArcurateA</td><td>(T-4668A, SRT, “Arcuate Artery of the Kidney”)</td></tr><tr><td>Aorta</td><td>(T-4200, SRT, ”Aorta”)</td></tr><tr><td>RenalVolume</td><td>(G-D705, SRT, “Renal Volume”)</td></tr><tr><td>NewBladderDistL</td><td>(G-A22A, SRT, “Bladder Length”)</td></tr><tr><td>NewBladderDistH</td><td>(121207, DCM, “Bladder Height”)</td></tr><tr><td>NewBladderDistW</td><td>(G-A220, SRT, “Bladder Width”)</td></tr><tr><td>NewBladderVolume</td><td>(G-D705, SRT, “Bladder Volume”)</td></tr><tr><td>NewPostBladderDistL</td><td>(GEU-1004-32, GEU, “Post Void Bladder Length”)</td></tr><tr><td>NewPostBladderDistH</td><td>(GEU-1004-33, GEU, “Post Void Bladder Height”)</td></tr><tr><td>NewPostBladderDistW</td><td>(GEU-1004-34, GEU, “Post Void Bladder Width”)</td></tr><tr><td>NewPostBladderVolume</td><td>(GEU-1004-35, GEU, “Post Void Bladder Volume”)</td></tr><tr><td>CorticalThickness</td><td>(GEU-1007-13, GEU, “CorticalThickness”)</td></tr></table>	Anatomy GEU parameter	Code and Description	MRenalA	(T-46600, SRT, “Renal Artery”)	RenalV	(T-48740, SRT, “Renal Vein”)	SegmentalA	(T-46659, SRT, “Segmental Artery”)	InterlobarA	(T-4667D, SRT, “Interlobar Artery of Kidney”)	ArcurateA	(T-4668A, SRT, “Arcuate Artery of the Kidney”)	Aorta	(T-4200, SRT, ”Aorta”)	RenalVolume	(G-D705, SRT, “Renal Volume”)	NewBladderDistL	(G-A22A, SRT, “Bladder Length”)	NewBladderDistH	(121207, DCM, “Bladder Height”)	NewBladderDistW	(G-A220, SRT, “Bladder Width”)	NewBladderVolume	(G-D705, SRT, “Bladder Volume”)	NewPostBladderDistL	(GEU-1004-32, GEU, “Post Void Bladder Length”)	NewPostBladderDistH	(GEU-1004-33, GEU, “Post Void Bladder Height”)	NewPostBladderDistW	(GEU-1004-34, GEU, “Post Void Bladder Width”)	NewPostBladderVolume	(GEU-1004-35, GEU, “Post Void Bladder Volume”)	CorticalThickness	(GEU-1007-13, GEU, “CorticalThickness”)
Anatomy GEU parameter	Code and Description																																					
MRenalA	(T-46600, SRT, “Renal Artery”)																																					
RenalV	(T-48740, SRT, “Renal Vein”)																																					
SegmentalA	(T-46659, SRT, “Segmental Artery”)																																					
InterlobarA	(T-4667D, SRT, “Interlobar Artery of Kidney”)																																					
ArcurateA	(T-4668A, SRT, “Arcuate Artery of the Kidney”)																																					
Aorta	(T-4200, SRT, ”Aorta”)																																					
RenalVolume	(G-D705, SRT, “Renal Volume”)																																					
NewBladderDistL	(G-A22A, SRT, “Bladder Length”)																																					
NewBladderDistH	(121207, DCM, “Bladder Height”)																																					
NewBladderDistW	(G-A220, SRT, “Bladder Width”)																																					
NewBladderVolume	(G-D705, SRT, “Bladder Volume”)																																					
NewPostBladderDistL	(GEU-1004-32, GEU, “Post Void Bladder Length”)																																					
NewPostBladderDistH	(GEU-1004-33, GEU, “Post Void Bladder Height”)																																					
NewPostBladderDistW	(GEU-1004-34, GEU, “Post Void Bladder Width”)																																					
NewPostBladderVolume	(GEU-1004-35, GEU, “Post Void Bladder Volume”)																																					
CorticalThickness	(GEU-1007-13, GEU, “CorticalThickness”)																																					

TABLE 15.6.7. Renal Study Folder Code Maps

**VIVID AND ECHOPAC V206
CONFORMANCE**

**GE HEALTHCARE
STATEMENT**

DIRECTION DOC2652554 REV 3

(T-46002, SRT, "Artery of Abdomen")	(G-A101, SRT, "Left") for Left, (G-A100, SRT, "Right") for Right or (G-A103, SRT, "Unilateral") Or (G-C0E3, SRT, "Finding")	DCID 12111 or 12112 Abdominal Arteries (lateral or unilateral). DCID 12113 or 12114 Abdominal Veins (lateral or unilateral) or DCID 12115 Renal Vessels Or SRT T-C3000 Spleen Or SRT T-71000 Kidney		Anatomy GEU parameter	Code and Description
				Aorta	(T-42000, SRT, "Aorta")
				Celiac	(T-46400, SRT, "Celiac Axis")
				CHA	(T-46421, SRT, "Common Hepatic Artery")
				Splenic A	(T-46460, SRT, "Splenic Artery")
				SMA	(T-46510, SRT, "Superior Mesenteric Artery")
				IMA	(T-46520, SRT, "Inferior Mesenteric Artery")
				MRenalA	(T-46600, SRT, "Renal Artery")
				RenalV	(T-48740, SRT, "Renal Vein")
				SegmentalA	(T-46659, SRT, "Segmental Artery")
				InterlobarA	(T-4667D, SRT, "Interlobar Artery of Kidney")
				ArcuateA	(T-4668A, SRT, "Arcuate Artery of the Kidney")
				CIA	(T-46710, SRT, "Common Iliac Artery")
				PrHepatic	(T-46422, SRT, "Proper Hepatic Artery")
				GDA	(T-46440, SRT, Gastroduodenal Artery)
				IVC	(T-48710, SRT, "Inferior Vena Cava")
				Splenic V	(T-48890, SRT, Splenic Vein")
				Hepatic V	(T-48720, SRT, Hepatic Vein")
				MHV	(T-48726, SRT, Middle Hepatic Vein")
				MPV	(GEU-1004-65, 99GEMS, "Main Branch of Portal Vein")
				Portal V	(T-48810, SRT, "Portal Vein")
				SMV	(T-48840, SRT, "Superior Mesenteric Vein")
				TIPS	(G-036C, SRT, "Transjugular Intrahepatic Portosystemic Shunt")
				CIV	(T-48920, SRT, "Common Iliac Vein")

TABLE 15.6.8 Abdomen Study Folder Code Maps

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

[illegible]

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

15.7 TID 300 MEASUREMENT

	NL	Relation with Parent	Value Type	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			NUM	\$Measurement	1	M		Units = \$Units
3	>	HAS CONCEPT MOD	CODE	EV(GEU-1005-5, 99GEMS, Measurement Label)	1	UC	Only for AP or Trans	- AP (122675, Anterior-Posterior, DCM,) - Trans (G-A117. Transverse, SRT)
4	>	HAS CONCEPT MOD	CODE	EV(G-A1F8, SRT, "Topographical modifier")	1	U		See 15.9 GE Ultrasound Sidedness and Vessel Location
5	>	HAS CONCEPT MOD	CODE	EV(121401, DCM, "Derivation")	1	U		See 15.11 Derivation and Selection
6	>	HAS PROPERTIES	CODE	EV(121404, DCM, "Selection Status")	1	U		See 15.11 Derivation and Selection

15.8 GE ULTRASOUND MODES

GE Ultrasound Modes	Code Value
2D	(G-03A2, SRT, "2D mode")
CF	(R-409E2, SRT, "Doppler Color Flow")
PW	(R-409E4, SRT, "Doppler Pulsed")
MM	(G-0394, SRT, "M mode")
CW	(R-409E3, SRT, "Doppler Continuous Wave")

15.9 GE ULTRASOUND SIDEDNESS AND VESSEL LOCATION

Side	Code Value
Rt	(G-A100, SRT, "Right")
Lt	(G-A101, SRT, "Left")

Vessel Location	Code Value
Prox	(G-A118, SRT, "Proximal")
Mid	(G-A188, SRT, "Mid-longitudinal")
Dist	(G-A119, SRT, "Distal")

Note (*) when there is no Sidedness or Locations, the SR nodes are not populated.

15.10 SR MAPPING TABLE FOR VASCULAR BASE MEASUREMENT CONCEPT

GEU Measurement Parameter	Standard Measurement Concept Name
---------------------------	-----------------------------------

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

PS	(11726-7, LN, "Peak Systolic Velocity")
ED	(11653-3, LN, "End Diastolic Velocity")
MD	(11665-7, LN, "Minimum Diastolic Velocity")
Tamax	(11692-1, LN, "Time averaged peak velocity")
PI	(12008-9, LN, "Pulsatility Index")
RI	(12023-8, LN, "Resistivity Index")
PV	(11726-7, LN, Peak Velocity)
SD Ratio	(12144-2, LN, "Systolic to Diastolic Velocity Ratio")
DS Ratio	(122218, DCM, Diastolic/Systolic velocity ratio)
Accel	(20167-3, LN, "Acceleration Index")
AT	(20168-1, LN, "Acceleration Time")
TAMEAN	(20352-1, LN, "Time averaged mean velocity")
VOLFLOW	(33878-0, LN, "Volume flow")
ICACCA Ratio (PS)	(33868-1, LN, "ICA/CCA velocity ratio")
HR (Heart Rate)	(8867-4, LN, Heart Rate)
AC	(GEU-1004-9, 99GEMS, "Angular Correction")
RAR	(33869-9, LN, "Renal Artery/Aorta velocity ratio")
PS/Hz	(GEU-1004-13, 99GEMS, 'Peak Systolic Frequency')
ED/Hz	(GEU-1004-14, 99GEMS, 'End Diastolic Frequency')
MD/Hz	(GEU-1004-16, 99GEMS, 'Minimum Diastolic Frequency')
PV/Hz	(GEU-1004-15, 99GEMS, 'Peak Velocity Frequency')
RefluxTime	(GEU-1004-22, 99GEMS, 'Reflux Time')
GreatSaphAccess	(GEU-1004-73, 99GEMS, 'Great Saphenous Vein of Accessory')

TABLE 15.10.1 PWD-MODE MEASUREMENTS

GEU Measurement Parameter	Standard Measurement Concept Name
DiamStenD1/D2	(G-0364 , SRT, " Vessel Lumen Diameter")
AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Sectional Area")
StenosisD	(R-101BB, SRT, " Lumen Diameter Stenosis")
StenosisA	(R-101BA, SRT, "Lumen Area Stenosis")
IMT Ant Avg	(GEU-1005-20, 99GEMS, "IMT Anterior Average")
IMT Ant Max	(GEU-1005-21, 99GEMS, "IMT Anterior Max")
IMT Ant Min	(GEU-1005-22, 99GEMS, "IMT Anterior Min")
IMT Ant SD	(GEU-1005-23, 99GEMS, "IMT Anterior SD")
IMT Ant nMeas	(GEU-1005-24, 99GEMS, "IMT Anterior nMeas")
IMT Post Avg	(GEU-1005-26, 99GEMS, "IMT Posterior Average")
IMT Post Max	(GEU-1005-27, 99GEMS, "IMT Posterior Max")
IMT Post Min	(GEU-1005-28, 99GEMS, "IMT Posterior Min")
IMT Post SD	(GEU-1005-29, 99GEMS, "IMT Posterior SD")
IMT Post nMeas	(GEU-1005-30, 99GEMS, "IMT Posterior nMeas")
Aorta AP	(Diameter, SRT, M-02550), with Measurement Label ("Anterior-Posterior", DCM 122675)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

Aorta Trans	(Diameter, SRT, M-02550), with Measurement Label (“Transverse”, SRT, G-A117)
ILA AP	(Diameter, SRT, M-02550) with Measurement Label (“Anterior-Posterior”, DCM, 122675)
ILA Trans	(Diameter, SRT, M-02440), with Measurement Label (“Transverse”, SRT, G-A117)
DiamRatioD1, DiamRatioD2	(121206, DCM, “Distance”)
ABDiamRatio	(GEU-1004-55, 99GEMS, “A/B Diameter Ratio”)
AreaRatioA1, AreaRatioA2	(121056, DCM, “Area Outline”)
ABAreaRatio	(GEU-1004-66, 99GEMS, “A/B Area Ratio”)
ComIliac AP	(M-02550, SRT, “Common Iliac Artery Diameter”) With Measurement Label (122675, DCM, “Anterior-Posterior”) Inside (T-46710, SRT, “Common Iliac Artery”)
ComIliacTrans	(M-02550, SRT, “Common Iliac Artery Diameter”) With Measurement Label (G-A117, SRT, “Transverse”) Inside (T-46710, SRT, “Common Iliac Artery”)
ComIliac DiamStenD1/D2	(G-0364, SRT, “Vessel Lumen Diameter”) Inside (T-46710, SRT, “Common Iliac Artery”)
ComIliac StenosisD	(R-101BB, SRT, “Lumen Diameter Stenosis”) Inside (T-46710, SRT, “Common Iliac Artery”)
ComIliac AreaStenA1/A2	(G-0366, SRT, “Vessel Lumen Cross-Sectional Area”) Inside (T-46710, SRT, “Common Iliac Artery”)
ComIliac StenosisA	(R-101BA, SRT, “Lumen Area Stenosis”) Inside (T-46710, SRT, “Common Iliac Artery”)
ExtIliac AP	(M-02550, SRT, “External Iliac Artery Diameter”) with Measurement Label (122675, DCM, “Anterior-Posterior”) Inside (T-46910, SRT, “External Iliac Artery”)
ExtIliac Trans	(M-02550, SRT, “External Iliac Artery Diameter”) with Measurement Label (G-A117, SRT, “Transverse”) Inside (T-46910, SRT, “External Iliac Artery”)
ExtIliac DiamStenD1/D2	(G-0364, SRT, “Vessel Lumen Diameter”) Inside (T-46910, SRT, “External Iliac Artery”)
ExtIliac StenosisD	(R-101BB, SRT, “Lumen Diameter Stenosis”) Inside (T-46910, SRT, “External Iliac Artery”)
ExtIliac AreaStenA1/A2	(G-0366, SRT, “Vessel Lumen Cross-Sectional Area”) Inside (T-46910, SRT, “External Iliac Artery”)
ExtIliac StenosisA	(R-101BA, SRT, “Lumen Area Stenosis”) Inside (T-46910, SRT, “External Iliac Artery”)
ComFemoral AP	(M-02550, SRT, “Common Femoral Artery Diameter”) with Measurement Label (122675, SRT, “Anterior-Posterior”) Inside (T-47400, SRT, “Common Femoral Artery”)
ComFemoral Trans	(M-02550, SRT, “Common Femoral Artery Diameter”) with Measurement Label (G-A117, SRT, “Transverse”) Inside (T-47400, SRT, “Common Femoral Artery”)
ComFemoral Diam StenD1/D2	(G-0364, SRT, “Vessel Lumen Diameter”) Inside (T-47400, SRT, “Common Femoral Artery”)
ComFemoral StenosisD	(R-101BB, SRT, “Lumen Diameter Stenosis”) Inside (T-47400, SRT, “Common Femoral Artery”)
ComFemoral AreaStenA1/A2	(G-0366, SRT, “Vessel Lumen Cross-Sectional Area”) Inside (T-47400, SRT, “Common Femoral Artery”)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

ComFemoral StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") Inside (T-47400, SRT, "Common Femoral Artery")
AP	(M-02550, SRT, 'Diameter') with Measurement Labels (122675, DCM, 'Anterior-Posterior')
2D AP	(M-02550, SRT, "Aorta") with Measurement Label (122675, SRT, "Anterior-Posterior") inside (T-42000, SRT, Aorta)
Trans	(M-02550, SRT, 'Diameter') with Measurement Labels (G-A117, SRT, 'Transverse')
2D Trans	(M-02550, SRT, "Aorta") with Measurement Label (G-A117, SRT, "Transverse") inside (T-42000, SRT, Aorta)
VFDiam	(GEU-1004-49, 99GEMS, 'Volume Flow Diameter')

TABLE 15.10.2 VASCULAR B-MODE MEASUREMENTS

GEU Measurement Parameter	Standard Measurement Concept Name
LiverDistL	(G-A22A, SRT, "Liver Length") Inside (T-62002, SRT, Liver) container
LiverDistH	(121207, DCM, "Liver Height") Inside (T-62002, SRT, Liver) container
LiverDistW	(G-A220, SRT, "Liver Width") Inside (T-62002, SRT, Liver) container
LiverVolume	(G-D705, SRT, "Liver Volume") Inside (T-62002, SRT, Liver) container
MassDistL	(G-A22A, SRT, "Mass Length") Inside (M-03000, SRT, Mass) container
MassDistH	(121207, DCM, "Mass Height") Inside (M-03000, SRT, Mass) container
MassDistW	(G-A220, SRT, "Mass Width") Inside (M-03000, SRT, Mass) container
MassVolume	(G-D705, SRT, "Mass Volume") Inside (M-03000, SRT, Mass) container
CystDistL	(G-A22A, SRT, "Cyst Length") Inside (M-3340A, SRT, Cyst) container
CystDistH	(121207, DCM, "Cyst Height") Inside (M-3340A, SRT, Cyst) container
CystDiswW	(G-A220, SRT, "Cyst Width") Inside (M-3340A, SRT, Cyst) container
CystDistVolume	(G-A220, SRT, "Cyst Volume") Inside (M-3340A, SRT, Cyst) container
PancHead	(GEU-1004-59, GEU, "Pancreas Head Diameter") Inside (T-D4034, SRT, Pancreas) container
PancDuct	(GEU-1004-61, GEU, "Pancreas Duct Diameter") Inside (T-D4034, SRT, Pancreas) container
PancBody	(GEU-1004-60, GEU, "Pancreas Body Diameter") Inside (T-D4034, SRT, Pancreas) container
CBDporta	(GEU-1004-67, GEU, "Common bile duct porta") Inside (T-60610, SRT, Bile Duct) container
CBDpanc	(GEU-1004-68, GEU, "Common bile duct pancreas")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	Inside (T-60610, SRT, Bile Duct) container
GBW	(GEU-1004-38, GEU, "Thickness of Gall Bladder Wall") Inside (T-63000, SRT, Gall Bladder) container
GBL	(G-A22A, SRT, "Gall Bladder Length") Inside (T-63000, SRT, Gall Bladder) container
AoProxAP	(M-02550, SRT, "Aorta Diameter") with (122675, DCM, "Anterior-Posterior") inside (T-42000, SRT, Aorta) with (Proximal, SRT,G-A118) as Topographical Modifier
AoProxTrans	(M-02550, SRT, "Aorta Diameter") with (G-A117, SRT, "Transverse") inside (T-42000, SRT, Aorta) with (Proximal, SRT,G-A118) as Topographical Modifier
AoMidAP	(M-02550, SRT, "Aorta Diameter") with (122675, DCM, "Anterior-Posterior") inside (T-42000, SRT, Aorta) with (Mid-longitudinal, SRT,G-A188) as Topographical Modifier
AoMidTrans	(M-02550, SRT, "Aorta Diameter") with (122675, DCM, "Anterior-Posterior") inside (T-42000, SRT, Aorta) with (Mid-longitudinal, SRT,G-A188) as Topographical Modifier
AoDistAP	(M-02550, SRT, "Aorta Diameter") with (122675, DCM, "Anterior-Posterior") inside (T-42000, SRT, Aorta) with (Distal, SRT,G-A119) as Topographical Modifier
AoDistTrans	(M-02550, SRT, "Aorta Diameter") with (G-A117, SRT, "Transverse") inside (T-42000, SRT, Aorta) with (Distal, SRT,G-A119) as Topographical Modifier
IliacAP	(M-02550, SRT, "Iliac artery Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (T-46710, SRT, Common Iliac Artery)
IliacTrans	(M-02550, SRT, "Iliac artery Diameter") with Measurement Label (G-A117, SRT, "Transverse") inside (T-46710, SRT, Common Iliac Artery)
BladderDistL	(G-A22A, SRT, "Bladder Length") Inside (T-74000, SRT, Bladder) container
BladderDistH	(121207, DCM, "Bladder Height") Inside (T-74000, SRT, Bladder) container
BladderDistW	(G-A220, SRT, "Bladder Width") Inside (T-74000, SRT, Bladder) container
BladderVolume	(G-D705, SRT, "Bladder Volume") Inside (T-74000, SRT, Bladder) container
CeliacA AP	(M-02550, SRT, "Celiac Axis Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (T-46400, SRT, Celiac Axis)
CeliacA Trans	(M-02550, SRT, "Celiac Axis Diameter") with Measurement Label (G-A117, SRT, "Transverse") inside (T-46400, SRT, Celiac Axis)
CeliacA DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (T-46400, SRT, Celiac Axis)
CeliacA Stenosis D	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (T-46400, SRT, Celiac Axis)
CeliacA AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (T-46400, SRT, Celiac Axis)
CeliacA StenosisA	(R-101BA, SRT, "Lumen Area Stenosis")

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

	inside (T-46400, SRT, Celiac Axis)
CeliacA SMA AP	(R-101BA, SRT, "Lumen Area Stenosis") inside (T-46400, SRT, Celiac Axis)
SMA AP	(M-02550, SRT, "Superior Mesenteric Artery Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (T-46510, SRT, Superior Mesenteric Artery)
SMA Trans	(M-02550, SRT, "Superior Mesenteric Artery Diameter") with Measurement Label (G-A117, SRT, "Transverse") inside (T-46510, SRT, Superior Mesenteric Artery)
SMA DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (T-46510, SRT, Superior Mesenteric Artery)
SMA StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (T-46510, SRT, Superior Mesenteric Artery)
SMA AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (T-46510, SRT, Superior Mesenteric Artery)
SMA StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (T-46510, SRT, Superior Mesenteric Artery)
IMA AP	(M-02550, SRT, "Inferior Mesenteric Artery Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (T-46520, SRT, Inferior Mesenteric Artery)
IMA Trans	(M-02550, SRT, "Inferior Mesenteric Artery Diameter") with Measurement Label (G-A117, SRT, "Transverse") inside (T-46520, SRT, Inferior Mesenteric Artery)
IMA DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (T-46520, SRT, Inferior Mesenteric Artery)
IMA StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (T-46520, SRT, Inferior Mesenteric Artery)
IMA AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (T-46520, SRT, Inferior Mesenteric Artery)
IMA StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (T-46520, SRT, Inferior Mesenteric Artery)
Stent AP	(M-02550, SRT, "Vessel Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (A-25500, SRT, Stent)
Stent Trans	(M-02550, SRT, "Vessel Diameter") with Measurement Label (G-A117, SRT, "Transverse") inside (A-25500, SRT, Stent)
Stent DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (A-25500, SRT, Stent)
Stent StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (A-25500, SRT, Stent)
Stent AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (A-25500, SRT, Stent)
Stent StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (A-25500, SRT, Stent)
PreStent AP	(M-02550, SRT, "Vessel Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (GEU-1004-71, GEU, Pre-Stent)
PreStent Trans	(M-02550, SRT, "Stent") with Measurement Label (G-A117, SRT, "Transverse") inside (GEU-1004-71, GEU, Pre-Stent)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

PreStent DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (GEU-1004-71, GEU, Pre-Stent)
PreStent StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (GEU-1004-71, GEU, Pre-Stent)
PreStent AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (GEU-1004-71, GEU, Pre-Stent)
PreStent StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (GEU-1004-71, GEU, Pre-Stent)
PostStent AP	(M-02550, SRT, "Vessel Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (GEU-1004-72, GEU, Post-Stent)
PostStent Trans	(M-02550, SRT, "Stent") with Measurement Label (G-A117, SRT, "Transverse") inside (GEU-1004-72, GEU, Post-Stent)
PostStent DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (GEU-1004-72, GEU, Post-Stent)
PostStent StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (GEU-1004-72, GEU, Post-Stent)
PostStent AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (GEU-1004-72, GEU, Post-Stent)
PostStent StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (GEU-1004-72, GEU, Post-Stent)
AortaAP	(M-02550, SRT, "Aorta Diameter") with Measurement Label (122675, DCM, "Anterior-Posterior") inside (T-42000, SRT, Aorta)
AortaTrans	(M-02550, SRT, "Aorta Diameter") with Measurement Label (G-A117, SRT, "Transverse") inside (T-42000, SRT, Aorta)
Aorta Sagittal	(M-02550, SRT, "Aorta Diameter") with Measurement Label (G-A117, SRT, "Sagittal") inside (T-42000, SRT, Aorta)
IIA DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (T-46740, SRT, Internal Iliac Artery)
IIA StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (T-46740, SRT, Internal Iliac Artery)
IIA AreaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (T-46740, SRT, Internal Iliac Artery)
IIA StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (T-46740, SRT, Internal Iliac Artery)
SpleenDistL	(G-A22A, SRT, Spleen Length) Inside (T-C3000, SRT, Spleen)
SpleenDistH	(121207, DCM, Spleen Height) Inside (T-C3000, SRT, Spleen)
SpleenDistW	(G-A220, SRT, Spleen Width) Inside (T-C3000, SRT, Spleen)
SpleenVolume	(G-D705, SRT, Spleen Volume) Inside (T-C3000, SRT, Spleen)
Aorta DiamStenD1/D2	(G-0364, SRT, "Vessel Lumen Diameter") inside (T-42000, SRT, Aorta)
Aorta StenosisD	(R-101BB, SRT, "Lumen Diameter Stenosis") inside (T-42000, SRT, Aorta)

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

AortaStenA1/A2	(G-0366, SRT, "Vessel Lumen Cross-Section Area") inside (T-42000, SRT, Aorta)
Aorta StenosisA	(R-101BA, SRT, "Lumen Area Stenosis") inside (T-42000, SRT, Aorta)

TABLE 15.10.3 ABDOMEN B-MODE MEASUREMENTS

15.11 DERIVATION AND SELECTION

GEU Name	Derivation	Selection
Av	Mean	User chosen value
Mx	Maximum	User chosen value
Mn	Minimum	User chosen value
Lt	Most recent value chosen	User chosen value
* (decided by another parameter)	Best Value	User chosen value

DIRECTION DOC2652554 REV 3


16. DICOM STRUCTURE REPORTS – USER DEFINED OBJECTS

The Vivid Ultrasound Scanner and EchoPAC provide the operator the ability to define measurements and export them in a manner compliant with established DICOM templates. This functionality is defined explicitly within the user manual for the Vivid scanner and EchoPAC application. Please consult the User Manual for more details and how to use this functionality.

16.1 GE DEFAULT DICOM EXPORT FORMAT

The GE Vivid Ultrasound Scanner and EchoPAC application allow the user to define measurements without specifying codes or meanings for the measurements. This default value is used when creating a user defined measurement within the Measurement menu. The measurements defined in this manner use a common structure.

Example of a user defined measurement format:



GE HEALTHCARE
STATEMENT

DIRECTION DOC2652554 REV 3

16.2 TEMPLATE EXTENSIONS

16.2.1 TID5100: Vascular Ultrasound Report Extensions

	NL	Rel with Parent	VT	Concept Name	V M	Reg Type	Condition	Value Set Constraint
...
1 2	>	CONTAINS	INCLUDES	TID(9900) User-defined concepts	1	U		

16.2.2 TID5200: Adult Echo Template Report Extensions

	NL	Rel with Parent	VT	Concept Name	V M	Reg Type	Condition	Value Set Constraint
...
3 1	>	CONTAINS	INCLUDES	TID(9900) User-defined concepts	1	U		

16.2.3 TID5220: Pediatric Template Report Extensions

	NL	Rel with Parent	VT	Concept Name	V M	Reg Type	Condition	Value Set Constraint
...
3 1	>	CONTAINS	INCLUDES	TID(9900) User-defined concepts	1	U		

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

16.2.4 TID9900: User-defined concepts

	NL	Rel with Parent	VT	Concept Name	V M	Reg Type	Condition	Value Set Constraint
1			CONTAINER	DT(T9900-01, 99GEMS, “User-defined concepts”)	1	M		
2	>	CONTAINS	INCLUDE	TID(9901) User-defined concept	1-n	MC	One of row 2 and 3 must be present	\$Type = DT (T9900-02, 99GEMS, “Measurement”)
3	>	CONTAINS	INCLUDE	TID(9901) User-defined concept	1-n	MC	One of row 2 and 3 must be present	\$Type = DT (T9900-03, 99GEMS, “Calculation”)

16.2.5 TID9901: User-defined concept

	NL	Rel with Parent	VT	Concept Name	V M	Reg Type	Condition	Value Set Constraint
1			CONTAINER	\$Type	1	M		
2	>	HAS CONCEPT MOD	TEXT	DT(T9900-04,99GEMS, “Label”)	1	M		
3	>	CONTAINS	INCLUDE	TID (300)	1	1-n	IFF \$Type = “Measurement”	\$Measurement = DCID (99008) Results \$Derivation – (3627) Measurement Type
4	>	CONTAINS	INCLUDE	TID (300)	1	1	IFF \$Type = “Calculation”	\$Measurement = DT(T9900-05, 99GEMS, “Value”)
5	>	INFERRED FROM	TEXT	DCID(228)	1	U		

DIRECTION DOC2652554 REV 3

16.3 USER DEFINED DICOM MAPPINGS

The Vivid Ultrasound Scanner and EchoPAC application do allow for more specific DICOM mappings for user defined measurements.

Within the DICOM Mapping configuration Interface, the user will be presented with the Category and Parameter to map. The desired measurement must be chosen from these two fields to initiate the mapping process.

When defining the measurement in the Mappings Interface, the user is required to specify the following DICOM fields:

- Finding Site
- Coding Scheme Value
- Coding Scheme Designator
- Coding Scheme Meaning

There are other DICOM attributes available within the Mapping Interface which are optional to define. If defined, they may help describe the measurement more accurately.

These include:

- View
- Phase
- Method
- Target
- Direction
- Respiratory Cycle Point

The values defined for the measurement will be presented to the user within a DICOM Encoding dialogue box within the DICOM Mapping configuration page. This visual representation shows the DICOM fields mapped when exporting the user defined object within the DICOM Structured Report.

17. SECURITY

17.1 INTRODUCTION

The security section describes security features implemented by this product. It includes description of non-DICOM network protocols, information to configure firewalls and application whitelists, list of supported DICOM security profiles as well as Web Security features. Additionally, secured media storage, VPN, etc. are also specified in this security section.

17.2 EXTERNAL NETWORK REQUIREMENTS

TABLE 17.2.1
EXTERNAL NETWORK REQUIREMENTS

Profile	Actor	Transaction	Protocol Used	RFCs	Security support	Reference
Basic Network Address Management	DHCP Client	Find and Use DHCP Server	DHCP	RFC2131 RFC2132 RFC2563		C.1.2
		Maintain Lease	DHCP	RFC2131 RFC2132		C.1.2
	DNS Client	Resolve Hostname	DNS	RFC1035 RFC2181		C.1.2
Application Configuration Management	LDAP Client	Find LDAP Server	LDAP	RFC2181 RFC2219 RFC2782		C.1.3
		Query LDAP Server	LDAP	RFC2251		C.1.3

17.3 TCP PORT CONFIGURATION

Ports and firewall configuration needed for interconnections supported by the Vivid scanner are described in Vivid™ Ultrasound Systems and EchoPAC™ Privacy and Security Manual - GD092163-1EN.

17.4 DICOM® Security PROFILE AVAILABILITY

17.4.1 Secure Use and User Identity Profiles

TABLE 17.4.1
SECURE USE AND USER IDENTITY PROFILES

Profile	Creator/Sender	Consumer/Receiver	Reference
---------	----------------	-------------------	-----------

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Online Electronic Storage Secure Use	N	N	C.2.1
Audit Trail Message Format	N	N	C.2.2
Audit Trail Message Transmission Profile - SYSLOG-TLS	N	N	C.2.3
Audit Trail Message Transmission Profile - SYSLOG-UDP	N	N	C.2.4
Basic User Identity Association	N	N	8.5
User Identity Plus Passcode Association	N	N	8.5
Kerberos Identity Negotiation Association	N	N	8.5
Generic SAML Assertion Identity Negotiation Association	N	N	8.5

17.4.2 Secure Transport Connection Profiles

TABLE 17.4.2
SECURE TRANSPORT CONNECTION PROFILES

Profile	Creator/Sender	Consumer/Receiver	Reference
BCP195 TLS Secure Transport Connection	Y	N	C.2.5
Non-Downgrading BCP195 TLS Secure Transport Connection	N	N	C.2.5
CRYPTREC TLS	N	N	C.2.5

17.4.3 Media Storage Security Profiles

Not Applicable.

17.4.4 Digital Signature Profiles

Not Applicable.

17.4.5 Additional DICOM® Security Profiles supported

Not Applicable.

17.5 USER IDENTITY NEGOTIATION SUPPORTED

Not Applicable.

17.6 WEB SERVICES SECURITY FEATURES

Not Applicable.

17.7 ADDITIONAL SECURITY FEATURES

Not Applicable.

**GE HEALTHCARE
STATEMENT**

DIRECTION DOC2652554 REV 3

17.7.1 Media storage security

Not Applicable.

17.7.2 Network security

Not Applicable.

17.7.3 Other security features

The Vivid system has multiple layers of security controls, including the following:

- user authentication with configurable user and password policies
- application whitelisting
- OS hardening
- embedded firewall
- encryption of patient data stored on the system.

See Vivid™ Ultrasound Systems and EchoPAC™ Privacy and Security Manual - GD092163-1EN for more details.

APPENDICES

A. Information Object Definitions (IODs)

This section provides the detailed content of the IODs natively created by the Vivid scanner, EchoPAC Software Only and EchoPAC Plug-in application, e.g., images created by an acquisition modality or evidence documents created on a review workstation application.

Throughout the tables listed in Annex A the following codes are used for the Source and Presence columns.

In the Source Column, the following values are supported:

FIXED: the value is pre-defined and cannot be modified.

GENERATED: the value is generated by the system.

CONFIGURATION: the value is copied from system configuration.

MWL: the value is copied from modality worklist.

USER: the value is entered by the user.

SCANNED: the value is read from a barcode scanner or similar device.

EMPTY: the attribute is sent without value.

SRC_INSTANCE: the value is copied from previously created instances.

The Presence columns reflect the usage of the module, functional group macro, attributes or value in the Vivid scanner and EchoPAC Software Only application implementation and is not necessarily the same as defined in the DICOM standard. For the Presence column the following values are supported:

ALWAYS: the module, functional group macro, attributes or value is always present

CONDITIONAL: the presence of the module, functional group macro, attributes or value is dependent on a condition. The condition must be listed in the Conditions column

EMPTY: The attribute is present but without a value (zero length)

A.1 Information Shared across multiple IODs

A.1.1 Shared Modules

All IODs generated by the system use the following common modules or a subset of them, as defined in the IOD specific subsections below.

TABLE A.1.1
SHARED MODULES

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Patient							
Referenced Patient Sequence	(0008,1120)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there. Not used in SR	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

						Documents	
>Referenced SOP Class UID	(0008,1150)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	Not used in SR Documents
>Referenced SOP Instance UID	(0008,1155)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	Not used in SR Documents
Patient's Name	(0010,0010)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Patient ID	(0010,0020)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Issuer of Patient ID	(0010,0021)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	
Issuer of Patient ID Qualifiers Sequence	(0010,0024)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Universal Entity ID	(0040,0032)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Universal Entity ID Type	(0040,0033)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Identifier Type Code	(0040,0035)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
Patient's Birth Date	(0010,0030)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Patient's Birth Time	(0010,0032)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface.	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

						Taken from worklist if it is there.	
Patient's Sex	(0010,0040)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Other Patient IDs	(0010,1000)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Other Patient IDs Sequence	(0010,1002)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
>Patient ID	(0010,0020)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
>Issuer of Patient ID	(0010,0021)	MWL	CONDITIONAL	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
>Type of Patient ID	(0010,0022)	FIXED	ALWAYS	"TEXT"			
>Issuer of Patient ID Qualifiers Sequence	(0010,0024)	MWL	CONDITIONAL	ALWAYS		Taken from worklist if it is there.	
>>Universal Entity ID	(0040,0032)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>>Universal Entity ID Type	(0040,0033)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

>>Identifier Type Code	(0040,0035)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
Ethnic Group	(0010,2160)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	
Patient Comments	(0010,4000)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	
Patient Identity Removed	(0012,0062)	GENERATED	CONDITIONAL	CONDITIONAL		Present and set to "YES" if patient is anonymized.	
De-identification Method	(0012,0063)	GENERATED	CONDITIONAL	CONDITIONAL		Present and set to "GEVU anonymization" if patient is anonymized.	
General Study							
Study Date	(0008,0020)	GENERATED	ALWAYS	ALWAYS			Is set to examination date
Study Time	(0008,0030)	GENERATED	ALWAYS	ALWAYS			Is set to examination time
Accession Number	(0008,0050)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Referring Physician's Name	(0008,0090)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Study Description	(0008,1030)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

						(from Requested Procedure Description).	
Physician(s) of Record	(0008,1048)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there (from Names of Intended Recipients of Result)	Not used in SR Documents
Referenced Study Sequence	(0008,1110)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there. Not used in SR Documents .	
>Referenced SOP Class UID	(0008,1150)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	
>Referenced SOP Instance UID	(0008,1155)	MWL	ALWAYS	CONDITIONAL		Taken from worklist if it is there.	
Study Instance UID	(0020,000D)	GENERATED/MWL	ALWAYS	CONDITIONAL		Uniquely generated by the equipment. Taken from worklist if it is there.	
Study ID	(0020,0010)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there (from Requested Procedure Id)	
Patient Study							
Patient's Size	(0010,1020)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Patient's Weight	(0010,1030)	USER/MWL	ALWAYS	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there.	
Additional Patient History	(0010,21B0)	USER/MWL	CONDITIONAL	CONDITIONAL		May be entered from User Interface (in Referral reason). Taken from worklist if it is there.	
Admission ID	(0038,0010)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
General Series							
Series Date	(0008,0021)	GENERATED	ALWAYS	ALWAYS			Is set to Series date
Series Time	(0008,0031)	GENERATED	ALWAYS	ALWAYS			Is set to Series time
Modality	(0008,0060)	GENERATED	ALWAYS	ALWAYS	"US" or "SR"		
Series Description	(0008,103E)	USER	ALWAYS	CONDITIONAL		May be entered from User Interface (in Diagnosis).	
Performing Physician's Name	(0008,1050)	USER/MWL	CONDITIONAL	CONDITIONAL		May be entered from User Interface. Taken from worklist if it is there (from Scheduled Performing Physician's Name)	
Operators' Name	(0008,1070)	USER/CONFIGURATION	ALWAYS	CONDITIONAL		May be entered from User Interface. Default is login id.	
Referenced Performed	(0008,1111)	SRC_INSTANCE	CONDITIONAL	ALWAYS		Used if Modality	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Procedure Step Sequence						Performed Procedure Step is enabled.	
>Referenced SOP Class UID	(0008,1150)	SRC_INSTANCE	ALWAYS	ALWAYS		Used if Modality Performed Procedure Step is enabled.	
>Referenced SOP Instance UID	(0008,1155)	SRC_INSTANCE	ALWAYS	ALWAYS		Used if Modality Performed Procedure Step is enabled.	
Protocol Name	(0018,1030)	CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if the image is acquired in a stress protocol.	
Series Instance UID	(0020,000E)	GENERATED	ALWAYS	ALWAYS			Uniquely generated by the equipment
Series Number	(0020,0011)	GENERATED	ALWAYS	ALWAYS			Internal number which is incremented for each new series within a study.
Request Attributes Sequence	(0040,0275)		ALWAYS	ALWAYS			
>Accession Number	(0008,0050)	USER/MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Referenced Study Sequence	(0008,1110)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Study Instance UID	(0020,000D)	GENERATED/MWL	ALWAYS	ALWAYS		Taken from worklist if it is there.	
>Requested Procedure Description	(0032,1060)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Requested Procedure Code Sequence	(0032,1064)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Scheduled Procedure Step Description	(0040,0007)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

>Scheduled Protocol Code Sequence	(0040,0008)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Scheduled Procedure Step ID	(0040,0009)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
>Requested Procedure ID	(0040,1001)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	
Performed Procedure Step Start Date	(0040,0244)	GENERATED	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	
Performed Procedure Step Start Time	(0040,0245)	GENERATED	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	
Performed Procedure Step ID	(0040,0253)	GENERATED	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	
Performed Procedure Step Description	(0040,0254)	GENERATED	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	
Performed Protocol Code Sequence	(0040,0260)	MWL/CONFIGURATION	CONDITIONAL	CONDITIONAL		Taken from worklist or selected protocol.	
General Equipment							
Manufacturer	(0008,0070)	FIXED	ALWAYS	ALWAYS	See Table A.1.3		
Institution Name	(0008,0080)	CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if the configured value is not empty.	
Station Name	(0008,1010)	CONFIGURATION	ALWAYS	ALWAYS			
Institutional Department Name	(0008,1040)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		The default value is the configured Department name.	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

						May be overridden by the value entered in the Patient info screen.	
Manufacturer's Model Name	(0008,1090)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Device Serial Number	(0018,1000)	GENERATED	ALWAYS	ALWAYS		The value is the serial number of the scanner. Not used in EchoPAC	
Software Versions	(0018,1020)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
General Image							
Content Date	(0008,0023)	GENERATED	ALWAYS	ALWAYS			Set from Image date.
Content Time	(0008,0033)	GENERATED	ALWAYS	ALWAYS			Set from Image time.
Image Type	(0008,0008)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Acquisition DateTime	(0008,002A)	GENERATED	ALWAYS	ALWAYS			
Referenced Image Sequence	(0008,1140)	GENERATED	CONDITIONAL	CONDITIONAL		Written if the content of the object has been created based on other object(s)	
>Referenced SOP Class UID	(0008,1150)	SRC_INSTANCE	ALWAYS	ALWAYS			
>Referenced SOP Instance UID	(0008,1155)	SRC_INSTANCE	ALWAYS	ALWAYS			
Derivation Description	(0008,2111)	GENERATED	CONDITIONAL	CONDITIONAL		May be written and contain additional derivation information if the Image Type is DERIVED.	

VIVID AND ECHOPAC V206 CONFORMANCE

GE HEALTHCARE STATEMENT

DIRECTION DOC2652554 REV 3

Source Image Sequence	(0008,2112)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Used if the image is derived from another image	
>Referenced SOP Class UID	(0008,1150)	SRC_INSTANCE	ALWAYS	ALWAYS			May be used for source images.
>Referenced SOP Instance UID	(0008,1155)	SRC_INSTANCE	ALWAYS	ALWAYS			May be used for source images.
Instance Number	(0020,0013)	GENERATED	ALWAYS	ALWAYS			Internal value which is incremented for each image within a series.
Patient Orientation	(0020,0020)	FIXED	ALWAYS	ALWAYS	EMPTY		
Image Comments	(0020,4000)	GENERATED	CONDITIONAL	CONDITIONAL		Sent if the image is acquired in a stress protocol.	
Burned In Annotation	(0028,0301)	FIXED	CONDITIONAL	ALWAYS	YES	Only sent if the image contains burned in annotations	
Lossy Image Compression	(0028,2110)	GENERATED	CONDITIONAL	ALWAYS	01	If the image is lossy compressed	
					00	Otherwise	
Lossy Image Compression Ratio	(0028,2112)	GENERATED	CONDITIONAL	ALWAYS		Only sent if the image is lossy compressed	
Image Pixel							
Samples per Pixel	(0028,0002)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Photometric Interpretation	(0028,0004)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Rows	(0028,0010)	GENERATED	ALWAYS	ALWAYS			Value depends on scanning mode and configuration setup.

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Columns	(0028,0011)	GENERATED	ALWAYS	ALWAYS			Value depends on scanning mode and configuration setup.
Bits Allocated	(0028,0100)	FIXED	ALWAYS	ALWAYS	8		
Bits Stored	(0028,0101)	FIXED	ALWAYS	ALWAYS	8		
High Bit	(0028,0102)	FIXED	ALWAYS	ALWAYS	7		
Pixel Representation	(0028,0103)	FIXED	ALWAYS	ALWAYS	0		0 – unsigned integer
Planar Configuration	(0028,0006)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Red Palette Color Lookup Table Descriptor	(0028,1101)	GENERATED	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Green Palette Color Lookup Table Descriptor	(0028,1102)	GENERATED	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Blue Palette Color Lookup Table Descriptor	(0028,1103)	GENERATED	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Red Palette Color Lookup Table Data	(0028,1201)	GENERATED	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Green Palette Color Lookup Table Data	(0028,1202)	GENERATED	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Blue Palette Color Lookup Table Data	(0028,1203)	GENERATED	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Pixel Data	(7FE0,0010)	GENERATED	ALWAYS	ALWAYS			Pixel Data of image.
Palette Color Lookup Table							
Red Palette Color Lookup Table Descriptor	(0028,1101)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Green Palette Color Lookup Table Descriptor	(0028,1102)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Blue Palette Color Lookup Table Descriptor	(0028,1103)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Red Palette Color Lookup Table Data	(0028,1201)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Green Palette Color Lookup Table Data	(0028,1202)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
Blue Palette Color Lookup Table Data	(0028,1203)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading Palette images.	
US Region Calibration							
Sequence of Ultrasound Regions	(0018,6011)	GENERATED	ALWAYS	ALWAYS			
>Region Spatial Format	(0018,6012)	GENERATED	ALWAYS	ALWAYS			
>Region Data Type	(0018,6014)	GENERATED	ALWAYS	ALWAYS			
>Region Flags	(0018,6016)	GENERATED	ALWAYS	ALWAYS	Bit 0: 0 (Opaque) Bit 1: 0 (Not Protected because there may be other regions within the image) Bit 2: 0 (Velocity)		
>Region Location Min X0	(0018,6018)	GENERATED	ALWAYS	ALWAYS			Varies with scanning mode.
>Region Location Min Y0	(0018,601A)	GENERATED	ALWAYS	ALWAYS			Varies with scanning mode.
>Region Location Max X1	(0018,601C)	GENERATED	ALWAYS	ALWAYS			Value is image width-1.

VIVID AND ECHOPAC V206 CONFORMANCE

GE HEALTHCARE STATEMENT

DIRECTION DOC2652554 REV 3

>Region Location Max Y1	(0018,601E)	GENERATED	ALWAYS	ALWAYS			Value is image height-1.
>Reference Pixel X0	(0018,6020)	GENERATED	CONDITIONAL	ALWAYS		Varies with scanning mode.	
>Reference Pixel Y0	(0018,6022)	GENERATED	CONDITIONAL	ALWAYS		Varies with scanning mode.	
>Physical Units X Direction	(0018,6024)	GENERATED	ALWAYS	ALWAYS			Values supported: 3 (cm) 4 (seconds)
>Physical Units Y Direction	(0018,6026)	GENERATED	ALWAYS	ALWAYS			Values supported: 3 (cm) 4 (seconds) 7 (cm/sec)
>Reference Pixel Physical Value X	(0018,6028)	GENERATED	CONDITIONAL	ALWAYS		Varies with scanning mode.	
>Reference Pixel Physical Value Y	(0018,602A)	GENERATED	CONDITIONAL	ALWAYS		Varies with scanning mode.	
>Physical Delta X	(0018,602C)	GENERATED	ALWAYS	ALWAYS			Varies with scanning mode.
>Physical Delta Y	(0018,602E)	GENERATED	ALWAYS	ALWAYS			Varies with scanning mode.
>Transducer Frequency	(0018,6030)	GENERATED	CONDITIONAL	ALWAYS		Sent if relevant for the image	
>Pulse Repetition Frequency	(0018,6032)	GENERATED	CONDITIONAL	ALWAYS		Sent if relevant for the image	
>Pixel Component Organization	(0018,6044)	GENERATED	CONDITIONAL	ALWAYS		Pixel component calibration data does not exist for any region.	
US Image							
Image Type	(0008,0008)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Stage Name	(0008,2120)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in	Name of stage of stress test.

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

						a stress test.	
Stage Number	(0008,2122)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in a stress test.	The stress test stage number when acquiring the image; starting at value 1.
Number of Stages	(0008,2124)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in a stress test.	
View Name	(0008,2127)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in a stress test.	The stress test view name when acquiring the image. The name is defined in the User Interface.
View Number	(0008,2128)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in a stress test.	The stress test view number when acquiring the image; starting at value 1.
Number of Event Timers	(0008,2129)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading images.	
Number of Views in Stage	(0008,212A)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in a stress test.	
Event Elapsed Time(s)	(0008,2130)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading images.	
Event Timer Name(s)	(0008,2132)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only used when reading images.	
Heart Rate	(0018,1088)	GENERATED	ALWAYS	ALWAYS			Set to the detected heart rate.
R Wave Time Vector	(0018,6060)	GENERATED	CONDITIONAL	CONDITIONAL		Used if relevant for the image	May be filled in with timing information.
Samples per Pixel	(0028,0002)	GENERATED	ALWAYS	ALWAYS	See		

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

					Table A.1.3		
Photometric Interpretation	(0028,0004)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Planar Configuration	(0028,0006)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
Frame Increment Pointer	(0028,0009)	GENERATED	CONDITIONAL	CONDITIONAL			
Ultrasound Color Data Present	(0028,0014)	GENERATED	ALWAYS	ALWAYS			
Bits Allocated	(0028,0100)	FIXED	ALWAYS	ALWAYS	8		
Bits Stored	(0028,0101)	FIXED	ALWAYS	ALWAYS	8		
High Bit	(0028,0102)	FIXED	ALWAYS	ALWAYS	7		
Pixel Representation	(0028,0103)	FIXED	ALWAYS	ALWAYS	0		
Lossy Image Compression	(0028,2110)	GENERATED	ALWAYS	CONDITIONAL		Set to 1 if the image is compressed using JPEG Baseline compression. Otherwise, set to 0.	
Stage Code Sequence	(0040,000A)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if image is acquired in a stress test.	Coded stage name of stress test. The name is defined in the User Interface.
View Code Sequence	(0054,0220)	USER/CONFIGURATION	CONDITIONAL	CONDITIONAL		Sent if the image is acquired in a stress test.	The coded view name of the stress test. The name is defined in the User Interface.
VOI LUT							
Window Center	(0028,1050)	FIXED	CONDITIONAL	CONDITIONAL	Value set to 127 if Photometric Interpretation has value MONOCHROME2	Only used for images created from CT data.	
Window Width	(0028,1051)	FIXED	CONDITIONAL	CONDITIONAL	Value set to 256 if Photometric Interpretation has value MONOCHROME2	Only used for images created from CT data.	
SOP Common							

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

SOP Class UID	(0008,0016)	GENERATED	ALWAYS	ALWAYS	See Table A.1.3		
SOP Instance UID	(0008,0018)	GENERATED	ALWAYS	ALWAYS			Uniquely generated by the equipment.

A.1.2 Common Functional Group Macros

N/A

A.1.3 Shared Private Modules

Table A.1.2 lists private attributes that are used in multiple IODs generated by the system. For documentation convenience and readability, they are organized in modules, although the concept of modules does not exist in the standard for private attributes.

TABLE A.1.2
SHARED PRIVATE MODULES

Attribute Name	Tag	VR	VM	Con- tains PHI	Presence (Attribute)	Presence (Value)	Value	Condition s	Description
GEMS_Ultrasound_MovieGroup_001									
Private Creator	7FE1,00xx	LO	1	No	CONDITIO NAL	ALWAYS	GEMS_Ultr asound_M ovieGroup _001	Used in DICOM US modality objects if Allow Raw Data is enabled	If so configured, the product will send ultrasound raw data information in private data elements designated by this Private Creator element. All private tags starting with 7FE1,xx will belong to the GEMS_Ultrasoun d_MovieGroup_ 001.
GEMS_Ultrasound_ImageGroup_001									
Private Creator	6003,00xx	LO	1	Yes	CONDITIO NAL	ALWAYS	GEMS_Ultr asound_Im ageGroup_ 001	Used in DICOM US modality objects	The product will include ultrasound preview image in private data elements designated by the Private Creator element.

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Attribute Name	Tag	VR	VM	Contains PHI	Presence (Attribute)	Presence (Value)	Value	Conditions	Description
									All private tags starting with 6003,xx will belong to the GEMS_Ultrasound_ImageGroup_001.
GEMS_Ultrasound_ExamGroup_001									
Private Creator	6005,00xx	LO	1	Yes	CONDITIONAL	ALWAYS	GEMS_Ultrasound_ExamGroup_001	Used in DICOM SR modality objects with Allow Private Data enabled	The product will send exam information in private data elements designated by the Private Creator element All private tags starting with 6005,00xx will belong to the GEMS_Ultrasound_ExamGroup_001.

A.1.4 Shared Values and Code Sets

The following Shared Values and Code Sets are used in multiple IODs generated by the system.

TABLE A.1.3
SHARED VALUES AND CODE SETS

Attribute Name	Tag	Value/Code	Condition	Comments
Image Type	(0008,0008)	The first two values contain "ORIGINAL\PRIMARY" or "DERIVED\PRIMARY". Value 3 is always empty. Value 4 is a description of the mode. Values beyond this may be used for private data.		
Manufacturer	(0008,0070)	GE Vingmed Ultrasound		
SOP Class UID	(0008,0016)	1.2.840.10008.5.1.4.1.1.3.1	Ultrasound Multiframe Image Storage	
		1.2.840.10008.5.1.4.1.1.3	Ultrasound Multiframe Image Storage (retired)	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

		1.2.840.10008.5.1.4.1.1.6.1	Ultrasound Image Storage	
		1.2.840.10008.5.1.4.1.1.6	Ultrasound Image Storage (retired)	
		1.2.840.10008.5.1.4.1.1.7	Secondary Capture Image Storage	
		1.2.840.10008.5.1.4.1.1.88.33	Comprehensive SR	
Manufacturer's Model Name	(0008,1090)	Vivid E80	Based on the product, the value is set to any of these	
		Vivid E90		
		Vivid E95		
		Vivid iq		
		Vivid S60		
		Vivid S70		
		Vivid T8		
		Vivid T9		
		EchoPAC		
		EchoPAC PC Integrated		
Software Versions	(0018,1020)	<Model Name>:206.<x>	The <Model Name> has the same value as element (0008,1090), whereas <x> denotes the current minor/patch/build version of the SW.	
Samples per Pixel	(0028,0002)	1	If the Photometric Interpretation element value is any of: <ul style="list-style-type: none"> • MONOCHROME 2 • PALETTE COLOR (only for read) 	
		3	If the Photometric Interpretation element value is any of: <ul style="list-style-type: none"> • RGB • YBR_FULL • YBR_FULL_422 	
Photometric Interpretation	(0028,0004)	MONOCHROME2	Grayscale Images	
		RGB	When compression is any of: <ul style="list-style-type: none"> • None • JPEG Lossless 	
		YBR_FULL	When compression is: <ul style="list-style-type: none"> • RLE 	
		YBR_FULL_422	When compression is:	

DIRECTION DOC2652554 REV 3

Planar Configuration	(0028,0006)		<ul style="list-style-type: none"> JPEG Lossy 	
		PALETTE COLOR	Only supported for read	
		0 (color by pixel)	When the Photometric Interpretation element value is any of: <ul style="list-style-type: none"> RGB YBR_FULL_422 	
		1 (color by plane)	When the Photometric Interpretation element value is: <ul style="list-style-type: none"> YBR_FULL 	

A.2 Ultrasound Image IOD

The following table defines the structure of Ultrasound Image IOD.

TABLE A.2.1
ULTRASOUND IMAGE IOD MODULES

Module Name	Presence (Module)	Condition	Reference
Patient	ALWAYS		Table A.1.1
General Study	ALWAYS		Table A.1.1
Patient Study	ALWAYS		Table A.1.1
General Series	ALWAYS		Table A.1.1
General Equipment	ALWAYS		Table A.1.1
General Image	ALWAYS		Table A.1.1
Image Pixel	ALWAYS		Table A.1.1
Contrast/Bolus	CONDITIONAL	Only used when filled in by the user	Table A.2.1
Palette Color Lookup Table	CONDITIONAL	Only used when reading Palette images.	Table A.1.1
US Region Calibration	CONDITIONAL	Only used when the characteristics of the image (depth, scale, etc.) are constant throughout the whole recording interval.	Table A.1.1
US Image	ALWAYS		Table A.1.1
VOI LUT	CONDITIONAL	Only used for images created from CT data.	Table A.1.1
SOP Common	ALWAYS		Table A.1.1

A.2.1 Ultrasound Image IOD Specific Modules

See Table A.1.1

Shared modules”. None of the modules are specific to the Ultrasound Image IOD.

DIRECTION DOC2652554 REV 3

A.2.2 Ultrasound Image IOD Specific Functional Group Macros

N/A

A.2.3 Ultrasound Image IOD Specific Private Modules

N/A

A.2.4 Ultrasound Image IOD Specific Values and Code Sets

N/A

A.3 *Ultrasound Multi-frame Image IOD*

The following table defines the structure of Ultrasound Multi-frame Image IOD.

TABLE A.3.1
ULTRASOUND MULTI-FRAME IMAGE IOD

Module Name	Presence (Module)	Condition	Reference
Patient	ALWAYS		Table A.1.1
General Study	ALWAYS		Table A.1.1
Patient Study	ALWAYS		Table A.1.1
General Series	ALWAYS		Table A.1.1
General Equipment	ALWAYS		Table A.1.1
General Image	ALWAYS		Table A.1.1
Image Pixel	ALWAYS		Table A.1.1
Contrast/Bolus	CONDITIONAL	Only used when filled in by the user	Table A.1.1
Cine	ALWAYS		Table A.3.2
Multi-frame	ALWAYS		Table A.3.2
Palette Color Lookup Table	CONDITIONAL	Only used when reading Palette images	Table A.1.1
US Region Calibration	CONDITIONAL	Only used when the characteristics of the image (depth, scale, etc.) are constant throughout the whole recording interval.	Table A.1.1
US Image	ALWAYS		Table A.1.1
VOI LUT	CONDITIONAL	Only used for images created from CT data.	Table A.1.1
SOP Common	ALWAYS		Table A.1.1

A.3.1 Ultrasound Multi-frame Image IOD Specific Modules

TABLE A.3.2
ULTRASOUND MULTI-FRAME IMAGE IOD SPECIFIC MODULES

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Cine							
Start Trim	(0008,2142)	GENERATED	ALWAYS	ALWAYS			
Stop Trim	(0008,2143)	GENERATED	ALWAYS	ALWAYS			
Recommended Display Frame Rate	(0008,2144)	GENERATED	ALWAYS	ALWAYS			
Cine Rate	(0018,0040)	GENERATED	ALWAYS	ALWAYS			
Effective Duration	(0018,0072)	GENERATED	ALWAYS	ALWAYS			
Frame Time	(0018,1063)	GENERATED	ALWAYS	ALWAYS			Is set to the interframe time
Frame Time Vector	(0018,1065)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Only reading is supported.	The average frame time is calculated and written to the Frame Time attribute.
Frame Delay	(0018,1066)	GENERATED	ALWAYS	ALWAYS			
Actual Frame Duration	(0018,1242)	GENERATED	ALWAYS	ALWAYS			
Preferred Playback Sequencing	(0018,1244)	FIXED	ALWAYS	ALWAYS	0 (looping)		
Multiframe							
Number of Frames	(0028,0008)	GENERATED	ALWAYS	ALWAYS			Is set to the number of frames in image
Frame Increment Pointer	(0028,0009)	GENERATED	ALWAYS	ALWAYS			Is set to Frame Time (0018,1063) or Frame Time Vector (0018,1065)

A.3.2 Ultrasound Multi-frame Image IOD Specific Functional Group Macros

N/A

DIRECTION DOC2652554 REV 3

A.3.3 Ultrasound Multi-frame Image IOD Specific Private Modules

N/A

A.3.4 Ultrasound Multi-frame Image IOD Specific Values and Code Sets

N/A

A.4 Secondary Capture Image IOD

The following table defines the structure of Secondary Capture Image IOD.

TABLE A.4.1
SECONDARY CAPTURE IMAGE IOD

Module Name	Presence (Module)	Condition	Reference
Patient	ALWAYS		Table A.1.1
General Study	ALWAYS		Table A.1.1
Patient Study	ALWAYS		Table A.1.1
General Series	ALWAYS		Table A.1.1
General Equipment	ALWAYS		Table A.1.1
SC Equipment	ALWAYS		Table A.4.2
General Image	ALWAYS		Table A.1.1
Image Pixel	ALWAYS		Table A.1.1
SC Image	ALWAYS		Table A.4.2
VOI LUT	ALWAYS		Table A.1.1
SOP Common	ALWAYS		Table A.1.1

A.4.1 Secondary Capture Image IOD Specific Modules

TABLE A.4.2
SECONDARY CAPTURE IMAGE IOD SPECIFIC MODULES

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
SC Equipment							
Modality	(0008,0060)	FIXED	ALWAYS	ALWAYS	"US"		
Conversion Type	(0008,0064)	FIXED	ALWAYS	ALWAYS	"WSD"		
Secondary Capture Device ID	(0018,1010)	FIXED	ALWAYS	ALWAYS	Same as the value of (0008,1090) "Manufacturer's Model Name"		

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Secondary Capture Device Manufacturer	(0018,1016)	FIXED	ALWAYS	ALWAYS	Same as the value of (0008,0070) "Manufacturer"		
Secondary Capture Device Manufacturer's Model Name	(0018,1018)	FIXED	ALWAYS	ALWAYS	Same as the value of (0008,1090) "Manufacturer's Model Name"		
Secondary Capture Device Software Versions	(0018,1019)	FIXED	ALWAYS	ALWAYS	Same as the value of (0018,1020) "Software Versions"		
SC Image							
Date of Secondary Capture	(0018,1012)	GENERATED	ALWAYS	ALWAYS			
Time of Secondary Capture	(0018,1014)	GENERATED	ALWAYS	ALWAYS			

A.4.2 Secondary Capture Image IOD Specific Functional Group Macros

N/A

A.4.3 Secondary Capture Image IOD Specific Private Modules

N/A

A.4.4 Secondary Capture Image IOD Specific Values and Code Sets

N/A

A.5 Basic Directory IOD

The following table defines the structure of the Basic Directory IOD:

TABLE A.5.1
BASIC DIRECTORY IOD MODULES

Module Name	Presence (Module)	Condition	Reference
File Set Identification	ALWAYS		Table A.5.1
Directory Information	ALWAYS		Table A.5.1

A.5.1 Basic Directory IOD Specific Modules

TABLE A.5.2
BASIC DIRECTORY IOD SPECIFIC MODULES

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
File Set Identification Module							
File-set ID	(0004,1130)	FIXED	ALWAYS	ALWAYS	Empty		
Directory Information Module							
Offset of the First Directory Record of the Root Directory Entity	(0004,1200)	GENERATED	ALWAYS	ALWAYS			
Offset of the Last Directory Record of the Root Directory Entity	(0004,1202)	GENERATED	ALWAYS	ALWAYS			
File-set Consistency Flag	(0004,1212)	GENERATED	ALWAYS	ALWAYS	0	If FSC/FSU. Ignored when reading (if FSR)	0 = no known inconsistencies
Directory Record Sequence	(0004,1220)	GENERATED	ALWAYS	ALWAYS			Is created by FSC
>Offset of the Next Directory Record	(0004,1400)	GENERATED	ALWAYS	ALWAYS			
>Record In-use Flag	(0004,1410)	GENERATED	ALWAYS	ALWAYS	FFFFH (65535)	If FSC/FSU	The whole record is skipped during reading (if FSR) if the value is 0.
>Offset of Referenced Lower-Level Directory Entity	(0004,1420)	GENERATED	ALWAYS	ALWAYS			
>Directory Record Type	(0004,1430)	SRC_INSTANCE, GENERATED	ALWAYS	ALWAYS	PATIENT, STUDY, SERIES, IMAGE, SR DOCUMENT		These are the supported values when writing

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

>Referenced File ID	(0004,1500)	GENERATED	CONDITIONAL	CONDITIONAL	File path consisting of 5 elements: 1.GEMS_IMG or GEMS_REP 2.Month of exam 3.Day of exam 4.Patient initials and time of exam 5.Time stamp	Set only if Record Type is IMAGE or SR DOCUMENT	GEMS_IMG" Is used if IMAGE. GEMS_REP is used if SR DOCUMENT.
>Referenced SOP Class UID in File	(0004,1510)	GENERATED	CONDITIONAL	ALWAYS		Set only if Record Type is IMAGE or SR DOCUMENT	
>Referenced SOP Instance UID in File	(0004,1511)	GENERATED	CONDITIONAL	ALWAYS		Set only if Record Type is IMAGE or SR DOCUMENT	
>Referenced Transfer Syntax UID in File	(0004,1512)	GENERATED	CONDITIONAL	ALWAYS		Set only if Record Type is IMAGE or SR DOCUMENT	
>Record Selection Keys							See Table A.5.3 – Table A.5.7

A.5.2 Basic Directory IOD Specific Functional Group Macros

The records in the Directory Records Sequence are filled with additional attributes depending on the Record Type, as described in the below tables.

The following comment is referenced in the below tables:

- 1) Is filled in by FSC and FSU as in chapters 4, 5, 6 and 7.

TABLE A.5.3
ATTRIBUTES SPECIFIC TO RECORD TYPE PATIENT

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Specific Character Set	(0008,0005)	GENERATED	CONDITIONAL	CONDITIONAL			See 1)
Patient's Name	(0010,0010)	USER/MWL	ALWAYS	CONDITIONAL			
Patient ID	(0010,0020)	USER/MWL	ALWAYS	CONDITIONAL			See 1). If empty, a Patient Id is created by the equipment.
Patient's Birth Date	(0010,0030)	USER/MWL	ALWAYS	CONDITIONAL			See 1)

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Patient's Sex	(0010,0040)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Referenced Patient Sequence	(0008,1120)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Taken from worklist if it is there.	See 1)
>Referenced SOP Class UID	(0008,1150)	MWL	ALWAYS	CONDITIONAL			See 1)
>Referenced SOP Instance UID	(0008,1155)	MWL	ALWAYS	CONDITIONAL			See 1)
Patient's Birth Time	(0010,0032)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Other Patient Ids	(0010,1000)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Other Patient Names	(0010,1001)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Taken from source image if there	See 1)
Ethnic Group	(0010,2160)	MWL	ALWAYS	CONDITIONAL			See 1)
Patient Comments	(0010,4000)	MWL	ALWAYS	CONDITIONAL			See 1)

TABLE A.5.4
ATTRIBUTES SPECIFIC TO RECORD TYPE STUDY

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Specific Character Set	(0008,0005)	GENERATED	CONDITIONAL	CONDITIONAL			See 1)
Study Date	(0008,0020)	GENERATED	ALWAYS	ALWAYS			See 1). If empty, a Study Date is created by the equipment.
Study Time	(0008,0030)	GENERATED	ALWAYS	ALWAYS			See 1). If empty, a Study Time is created by the equipment.
Accession Number	(0008,0050)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Referring Physician's Name	(0008,0090)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Study Description	(0008,1030)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Physician(s) of Record	(0008,1048)	MWL	CONDITIONAL	CONDITIONAL			See 1)
Name of Physician(s) Reading Study	(0008,1060)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Taken from source image if there	See 1)
Referenced Study Sequence	(0008,1110)	FIXED	ALWAYS	ALWAYS	Empty		See 1)

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

>Referenced SOP Class UID	(0008,1150)	MWL	ALWAYS	CONDITIONAL			See 1)
>Referenced SOP Instance UID	(0008,1155)	MWL	ALWAYS	CONDITIONAL			See 1)
Admitting Diagnoses Description	(0008,1080)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Taken from source image if there	See 1)
Patient's Age	(0010,1010)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Taken from source image if there	See 1)
Patient's Size	(0010,1020)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Patient's Weight	(0010,1030)	USER/MWL	ALWAYS	CONDITIONAL			See 1)
Occupation	(0010,2180)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Taken from source image if there	See 1)
Additional Patient's History	(0010,2180)	USER/MWL	CONDITIONAL	CONDITIONAL			See 1)
Study Instance UID	(0020,000D)	GENERATED/MWL	ALWAYS	CONDITIONAL			See 1)
Study ID	(0020,0010)	USER/MWL	ALWAYS	CONDITIONAL			See 1). If empty, a Study Id is created by the equipment.

TABLE A.5.5
ATTRIBUTES SPECIFIC TO RECORD TYPE SERIES

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Specific Character Set	(0008,0005)	GENERATED	CONDITIONAL	CONDITIONAL			Is filled in by FSC or FSU as contained in the image or SR document message, if one of the tags contains extended characters
Series Date	(0008,0021)	GENERATED	ALWAYS	ALWAYS			As in 1) if the instance is IMAGE.
Series Time	(0008,0031)	GENERATED	ALWAYS	ALWAYS			As in 1) if the instance is IMAGE.
Modality	(0008,0060)	GENERATED	ALWAYS	ALWAYS	"US" or "SR"		See 1)
Manufacturer	(0008,0070)	FIXED	ALWAYS	ALWAYS	See Table A.1.3		See 1)

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Institution Name	(0008,0080)	CONFIGURATIO N	CONDITIONAL	CONDITIONAL			See 1)
Station Name	(0008,1010)	CONFIGURATIO N	ALWAYS	ALWAYS			See 1)
Series Description	(0008,103E)	USER	ALWAYS	CONDITIONAL			As in 1) if the instance is IMAGE.
Institutional Department Name	(0008,1040)	USER/ CONFIGURATIO N	CONDITIONAL	CONDITIONAL			See 1)
Performing Physicians' Name	(0008,1050)	USER/MWL	CONDITIONAL	CONDITIONAL			As in 1) if the instance is IMAGE.
Operator's Name	(0008,1070)	USER/ CONFIGURATIO N	ALWAYS	CONDITIONAL			As in 1) if the instance is IMAGE.
Manufacturer's Model Name	(0008,1090)	GENERATED	ALWAYS	ALWAYS			See 1)
Referenced Performed Procedure Step Sequence	(0008,1111)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	See 1)
>Referenced SOP Class UID	(0008,1150)	SRC_INSTANCE	ALWAYS	ALWAYS			See 1)
>Referenced SOP Instance UID	(0008,1155)	SRC_INSTANCE	ALWAYS	ALWAYS			See 1)
Software Versions	(0018,1020)	GENERATED	ALWAYS	ALWAYS			See 1)
Protocol Name	(0018,1030)	CONFIGURATIO N	CONDITIONAL	CONDITIONAL			As in 1) if the instance is IMAGE.
Series Instance UID	(0020,000E)	GENERATED	ALWAYS	ALWAYS			See 1)
Series Number	(0020,0011)	GENERATED	ALWAYS	ALWAYS			See 1). If empty, a Series Number is created by the equipment.
Performed Procedure Step Start Date	(0040,0244)	GENERATED	CONDITIONAL	CONDITIONAL			As in 1) if the instance is IMAGE.
Performed Procedure Step Start Time	(0040,0245)	GENERATED	CONDITIONAL	CONDITIONAL			As in 1) if the instance is IMAGE.
Performed Procedure Step ID	(0040,0253)	SRC_INSTANCE	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	As in 1) if the instance is IMAGE.

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Performed Procedure Step Description	(0040,0254)	GENERATED	CONDITIONAL	CONDITIONAL		Used if Modality Performed Procedure Step is enabled.	As in 1) if the instance is IMAGE.
Performed Protocol Code Sequence	(0040,0260)	MWL/ CONFIGURATIO N	CONDITIONAL	CONDITIONAL		Taken from worklist or selected protocol.	
Request Attributes Sequence	(0040,0275)		ALWAYS	CONDITIONAL		The content depends on the available values	As in 1) if the instance is IMAGE.
>Scheduled Procedure Step Description	(0040,0007)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist	As in 1) if the instance is IMAGE.
>Scheduled Protocol Code Sequence	(0040,0008)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist	As in 1) if the instance is IMAGE.
>Scheduled Procedure Step ID	(0040,0009)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist	As in 1) if the instance is IMAGE.
>Requested Procedure ID	(0040,1001)	MWL	CONDITIONAL	CONDITIONAL		Taken from worklist	As in 1) if the instance is IMAGE.

TABLE A.5.6
ATTRIBUTES SPECIFIC TO RECORD TYPE IMAGE

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Specific Character Set	(0008,0005)	GENERATED	CONDITIONAL	CONDITIONAL			See 1)
Image Type	(0008,0008)	GENERATED	ALWAYS	ALWAYS			See 1)
Content Date	(0008,0023)	GENERATED	ALWAYS	ALWAYS			See 1)
Content Time	(0008,0033)	GENERATED	ALWAYS	ALWAYS			See 1)
Contrast/Bolus Agent	(0018,0010)	USER	CONDITIONAL	CONDITIONAL			See 1)
Instance Number	(0020,0013)	GENERATED	ALWAYS	ALWAYS			See 1). If empty, an Instance Number is created by the equipment.
Photometric Interpretation	(0028,0004)	GENERATED	ALWAYS	ALWAYS			See 1)
Number Of Frames	(0028,0008)	GENERATED	ALWAYS	ALWAYS			See 1)
Rows	(0028,0010)	GENERATED	ALWAYS	ALWAYS			See 1)

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Columns	(0028,0011)	GENERATED	ALWAYS	ALWAYS			See 1)
Lossy Image Compression	(0028,2110)	GENERATED	CONDITIONAL	ALWAYS			See 1)
Lossy Image Compression Ratio	(0028,2112)	GENERATED	CONDITIONAL	ALWAYS			See 1)

TABLE A.5.7
ATTRIBUTES SPECIFIC TO RECORD TYPE SR DOCUMENT

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
Specific Character Set	(0008,0005)	GENERATED	CONDITIONAL	CONDITIONAL			See 1)
Content Date	(0008,0023)	GENERATED	ALWAYS	ALWAYS			See 1)
Content Time	(0008,0033)	GENERATED	ALWAYS	ALWAYS			See 1)
Instance Number	(0020,0013)	GENERATED	ALWAYS	ALWAYS			See 1). If empty, an Instance Number is created by the equipment.
Concept Name Code Sequence	(0040,A043)	GENERATED	ALWAYS	ALWAYS			See 1)
Completion Flag	(0040,A491)	FIXED	ALWAYS	ALWAYS	PARTIAL		See 1)
Verification Flag	(0040,A493)	FIXED	ALWAYS	ALWAYS	UNVERIFIED		See 1)

A.5.3 Basic Directory IOD Specific Private Modules

N/A

A.5.4 Basic Directory IOD Specific Values and Code Sets

N/A

A.6 Comprehensive SR IOD

The following table defines the structure of the Comprehensive Structured Report IOD:

TABLE A.6.1
COMPREHENSIVE SR IOD MODULES

Module Name	Presence (Module)	Condition	Reference
Patient	ALWAYS		Table A.1.1
General Study	ALWAYS		Table A.1.1
Patient Study	ALWAYS		Table A.1.1
SR Document Series	ALWAYS		Table A.6.2
General Equipment	ALWAYS		Table A.1.1
SR Document General	ALWAYS		Table A.6.2

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

SR Document Content	ALWAYS		Table A.6.2
SOP Common	ALWAYS		Table A.1.1

A.6.1 Comprehensive SR IOD Specific Modules

TABLE A.6.2
COMPREHENSIVE SR SPECIFIC MODULES

Attribute Name	Tag	Source	Presence (Attribute)	Presence (Value)	Value	Conditions	Comment
SR Document Series							
Modality	(0008,0060)	FIXED	ALWAYS	ALWAYS	SR		
Referenced Performed Procedure Step Sequence	(0008,1111)	GENERATED	ALWAYS	ALWAYS	Performed Procedure Step SOP Instance.	If MPPS enabled.	
					Empty.	If MPPS disabled.	
>Referenced SOP Class UID	(0008,1150)	GENERATED	CONDITIONAL	ALWAYS		Used if Modality Performed Procedure Step is enabled.	
>Referenced SOP Instance UID	(0008,1155)	GENERATED	CONDITIONAL	ALWAYS		Used if Modality Performed Procedure Step is enabled.	
Series Instance UID	(0020,000E)	GENERATED	ALWAYS	ALWAYS			Uniquely generated by the equipment
Series Number	(0020,0011)	GENERATED	ALWAYS	ALWAYS			Internal number which is incremented for each new series within a study.
SR Document General							
Content Date	(0008,0023)	GENERATED	ALWAYS	ALWAYS			Date of creation
Content Time	(0008,0033)	GENERATED	ALWAYS	ALWAYS			Time of creation
Instance Number	(0020,0013)	GENERATED	ALWAYS	ALWAYS			Internal number which is incremented for each new SR document within a series.
Referenced Request Sequence	(0040,A370)	GENERATED	ALWAYS	ALWAYS			

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

>Study Instance UID	(0020,000D)	GENERATED	ALWAYS	ALWAYS			Taken from Study Instance UID in General Study Module
>Referenced Study Sequence	(0008,1110)	FIXED	ALWAYS	ALWAYS	Empty		
>Accession Number	(0008,0050)	GENERATED	ALWAYS	ALWAYS			Taken from Accession Number in General Study Module
>Placer Order Number/Imaging Service Request	(0040,2016)	FIXED	ALWAYS	ALWAYS	Empty		
>Filler Order Number/Imaging Service Request	(0040,2017)	FIXED	ALWAYS	ALWAYS	Empty		
>Requested Procedure ID	(0040,1001)	GENERATED	ALWAYS	ALWAYS	Taken from worklist	If worklist is the source	
					Empty	If worklist is not the source	
>Requested Procedure Description	(0032,1060)				Taken from worklist	If worklist is the source	
					Empty	If worklist is not the source	
>Requested Procedure Code Sequence	(0032,1064)	FIXED	ALWAYS	ALWAYS	Empty		
Performed Procedure Code Sequence	(0040,A372)	FIXED	ALWAYS	ALWAYS	Empty		
Completion Flag	(0040,A491)	FIXED	ALWAYS	ALWAYS	PARTIAL		
Verification Flag	(0040,A493)	FIXED	ALWAYS	ALWAYS	UNVERIFIED		
SR Document Content							
Value Type	(0040,A040)	FIXED	ALWAYS	ALWAYS	CONTAINER		
Concept Name Code Sequence	(0040,A043)	GENERATED	ALWAYS	ALWAYS			
>Code Value	(0008,0100)	GENERATED	ALWAYS	ALWAYS	125200	If TID 5200	
					125195	If TID 5220	
					125100	If TID 5100	
>Coding Scheme Designator	(0008,0102)	GENERATED	ALWAYS	ALWAYS	DCM		
>Code Meaning	(0008,0104)	GENERATED	ALWAYS	ALWAYS	Adult Echocardiograph	If TID 5200	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

					y Procedure Report		
					Pediatric Cardiac Ultrasound Report	If TID 5220	
					Vascular Ultrasound Procedure Report	If TID 5100	
Continuity of Content	(0040,A050)	FIXED	ALWAYS	ALWAYS	SEPARATE		
Content Template Sequence	(0040,A504)	GENERATED	ALWAYS	ALWAYS			
>Mapping Resource	(0008,0105)	FIXED	ALWAYS	ALWAYS	DCMR		
>Template Identifier	(0040,DB00)	GENERATED	ALWAYS	ALWAYS	5200 or 5100 or 5220		
Content Sequence	(0040,A730)	GENERATED	ALWAYS	ALWAYS			See template "Adult Echocardiography Procedure Report" or "Pediatric Cardiac Ultrasound Report" or "Vascular Ultrasound Procedure Report".
> Relationship Type	(0040,A010)	GENERATED	ALWAYS	ALWAYS			See template "Adult Echocardiography Procedure Report" or "Pediatric Cardiac Ultrasound Report" or "Vascular Ultrasound Procedure Report".
> SR Document Content Module		GENERATED	ALWAYS	ALWAYS			See template "Adult Echocardiography Procedure Report" or "Pediatric Cardiac Ultrasound Report" or "Vascular Ultrasound Procedure Report".

GE HEALTHCARE
STATEMENT

DIRECTION DOC2652554 REV 3

**A.6.2 Comprehensive SR IOD Specific Functional Group
Macros**

N/A

A.6.3 Comprehensive SR IOD Specific Private Modules

N/A

A.6.4 Comprehensive SR IOD Specific Values and Code Sets

N/A

B. Structured Report Content Encoding

B.1 Vascular Ultrasound Procedure Report (TID 5100)

See chapter 15 “Vascular Ultrasound Procedure Report”.

B.2 Adult Echocardiography Procedure Report (TID 5200)

See chapters 14 “Adult and Pediatric Echocardiography Procedure Reports”.

B.3 Pediatric Cardiac Ultrasound Reports (TID 5220)

See chapters 14 “Adult and Pediatric Echocardiography Procedure Reports”.

B.4 User Defined Template Structure

See chapter 16 “DICOM Structure Reports – User Defined objects”.

C. Security Details

C.1 External Network Requirements details

C.1.1 Basic Time Synchronization

Not Applicable.

C.1.2 Basic Network Address Management

Standard functionality supported. Some of the configurable parameters are listed in the table below.

Parameters
DHCP Enabled/disabled
IP-address
Subnet mask
Default gateway
Obtain DNS address automatically
Preferred and alternative DNS server IP-address
List of DNS servers
NetBIOS settings
Speed and Duplex

C.1.3 Application Configuration Management

Actor	LDAP Security Pattern	Supported	Comments
LDAP Client	Anonymous	Y	
	Authenticated	Y	
	TLS - Anonymous	Y	
	TLS – Authenticated	Y	

C.1.4 DNS Service Discovery

Standard functionality supported.

C.2 DICOM® Security Profile Details

C.2.1 Online Electronic Storage Secure Use

Not Applicable.

DIRECTION DOC2652554 REV 3

C.2.2 Audit Trail Messages

Not Applicable.

C.2.3 Audit Trail Message Transmission Profile – Syslog – TLS

Not Applicable.

C.2.4 Audit Trail Message Transmission Profile - Syslog-UDP

Not Applicable.

C.2.5 Secure Transport Connection Details

The scanner supports import of public certificates and private certificate key through configuration user interface. The product uses these certificate/key for authentication and encryption when it acts as a TLS server and when it acts as a TLS client.

SECURE TRANSPORT CONNECTION PROFILES AND CIPHER SUITES

Profile	Cipher Suite	Default Preference Order (from 1=preferred to n=less preferred)
BCP195 TLS Secure Transport Connection	TLS_RSA_WITH_AES_128_CBC_SHA	1
	TLS_RSA_WITH_3DES_EDE_CBC_SHA	2
Non-Downgrading BCP195 TLS Secure Transport Connection	NA	
CRYPTREC TLS	NA	

The following table describes the Secure Transport Connection configuration parameters supported by this product:

Local Secure Transport Connection configuration parameters			
Parameter	Configurable	Default Value	Comment
Port	See Section Error! Reference source not found. Error! Reference source not found.		
Remote Secure Transport Connection configuration parameters			
Parameter	Configurable	Default Value	Comment
Port	See Section Error! Reference source not found. Error! Reference source not found.		

DIRECTION DOC2652554 REV 3

C.2.6 Attribute Confidentiality Details

Not Applicable.

C.2.7 Digital Signature details

Not Applicable.

C.2.8 Additional DICOM® Security Profiles supported

Not Applicable.

D. Mapping of Attributes

D.1 Modality Worklist, Instances and MPPS messages

The table below describes the mapping of attributes between Modality Worklist, Instances and MPPS messages.

In the Scenarios column the following values are supported:

SCHEDULED: the image acquisition was scheduled at the RIS and procedure details have been communicated in the MWL query)

UNSCHEDULED: the image acquisition was performed without Modality Worklist information

APPEND: instances acquired are added to an existing study after the initial procedure was finalized

GROUP: multiple requested procedures are grouped into one study.]

In the Value Source columns, the following values are supported:

GENERATED: the value is generated by the system.

SRC_INSTANCE: the value is copied from previously created instances.

MWL: the value is copied from modality worklist.

USER: the value is entered by the user.

SCANNED: the value is read from a barcode scanner or similar device.

EMPTY: the attribute is sent without value.

The Destination column uses either ROOT, if the attribute is added to the root of the instance or list the attribute Tag of the Sequence, the attribute will be added to. The comment column can be used to provide additional information regarding the values added to the IOD.

TABLE D.1.1
ATTRIBUTE MAPPING

Attribute Name	Tag	Scenario	Image		MPPS		Comments
			Value Source	Destination Tag	Value Source	Destination Tag	
Study Instance UID	(0020,000D)	SCHEDULED	MWL	ROOT	SRC_INSTANCE	(0040,0270)	
		UNSCHEDULED	GENERATED	ROOT	EMPTY	(0040,0270)	
		APPEND	SRC_INSTANCE	ROOT	SRC_INSTANCE	(0040,0270)	
Accession Number	(0008,0050)	SCHEDULED	MWL	ROOT	SRC_INSTANCE	(0040,0270)	
		UNSCHEDULED	EMPTY	ROOT	EMPTY	(0040,0270)	
		APPEND	SRC_INSTANCE	ROOT	SRC_INSTANCE	(0040,0270)	
Modality	(0008,0060)	SCHEDULED	MWL	ROOT	SRC_INSTANCE	ROOT	US or SR, from Scheduled Procedure Step Sequence (0040,0100), MWL from (0008,0060)

GE HEALTHCARE STATEMENT

VIVID AND ECHOPAC V206 CONFORMANCE

DIRECTION DOC2652554 REV 3

Attribute Name	Tag	Scenario	Image		MPPS		Comments
			Value Source	Destination Tag	Value Source	Destination Tag	
		UNSCHEDULED	GENERATED	ROOT	SRC_INSTANCE	ROOT	US or SR, from Scheduled Procedure Step Sequence (0040,0100), MWL from (0008,0060)
		APPEND	SRC_INSTANCE	ROOT	SRC_INSTANCE	ROOT	US or SR, from Scheduled Procedure Step Sequence (0040,0100), MWL from (0008,0060)
Referenced Study Sequence	(0008,1110)	SCHEDULED	MWL	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INSTANCE	(0040,0270)	
		UNSCHEDULED	N/A	N/A	EMPTY	(0040,0270)	
		APPEND	SRC_INSTANCE	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INSTANCE	(0040,0270)	
Referenced SOP Class UID	(0008,1150)	SCHEDULED	MWL	ROOT	SRC_INSTANCE	(0008,1120)	
		UNSCHEDULED	EMPTY	ROOT	SRC_INSTANCE	(0008,1120)	
		APPEND	SRC_INSTANCE	ROOT	SRC_INSTANCE	(0008,1120)	
Referenced SOP Instance UID	(0008,1155)	SCHEDULED	MWL	ROOT	SRC_INSTANCE	(0008,1120)	
		UNSCHEDULED	EMPTY	ROOT	SRC_INSTANCE	(0008,1120)	
		APPEND	SRC_INSTANCE	ROOT	SRC_INSTANCE	(0008,1120)	
Study ID	(0020,0010)	SCHEDULED	GENERATED	ROOT	SRC_INSTANCE	ROOT	Copied from Requested Procedure ID (0040,1001)
		UNSCHEDULED	USER	ROOT	SRC_INSTANCE	ROOT	Copied from Requested Procedure ID (0040,1001)
		APPEND	GENERATED	ROOT	SRC_INSTANCE	(0040,0270)	Copied from Requested Procedure ID (0040,1001)
Requested Procedure Description	(0032,1060)	SCHEDULED	MWL	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INSTANCE	(0040,0270)	To Performed Procedure Step Description (0040,0254)
		UNSCHEDULED	N/A	N/A	EMPTY	(0040,0270)	To Performed Procedure Step Description (0040,0254)
		APPEND	SRC_INSTANCE	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INSTANCE	(0040,0270)	To Performed Procedure Step Description (0040,0254)
Requested Procedure Code Sequence	(0032,1064)	SCHEDULED	MWL	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INSTANCE	(0008,1032)	To Procedure Code Sequence
		UNSCHEDULED	N/A	N/A	EMPTY	(0008,1032)	To Procedure Code Sequence
		APPEND	SRC_INSTANCE	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INSTANCE	(0008,1032)	To Procedure Code Sequence

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

Attribute Name	Tag	Scenario	Image		MPPS		Comments
			Value Source	Destination Tag	Value Source	Destination Tag	
Requested Procedure ID	(0040,1001)	SCHEDULED	MWL	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INST ANCE	(0040,0270)	^(a) for use in Image IODs ^(b) for use in Evidence Documents
		UNSCHEDULED	N/A	N/A	EMPTY	(0040,0270)	
		APPEND	SRC_INST ANCE	(0040,0275) ^(a) (0040,A370) ^(b)	SRC_INST ANCE	(0040,0270)	^(a) for use in Image IODs ^(b) for use in Evidence Documents
Scheduled Performing Physician's Name	(0040,0006)	SCHEDULED	MWL	ROOT	SRC_INST ANCE	ROOT	Maps to the attribute "Performing Physician's Name"
		UNSCHEDULED	EMPTY	ROOT	SRC_INST ANCE	ROOT	Maps to the attribute "Performing Physician's Name"
		APPEND	SRC_INST ANCE	ROOT	SRC_INST ANCE	ROOT	Maps to the attribute "Performing Physician's Name"
Scheduled Procedure Step Description	(0040,0007)	SCHEDULED	MWL	(0040,0275) ^(a) N/A ^(b)	SRC_INST ANCE	(0040,0270)	
		UNSCHEDULED	N/A	N/A	EMPTY	(0040,0270)	
		APPEND	SRC_INST ANCE	(0040,0275) ^(a) N/A ^(b)	SRC_INST ANCE	(0040,0270)	
Scheduled Protocol Code Sequence	(0040,0008)	SCHEDULED	MWL	(0040,0275) ^(a) N/A ^(b)	SRC_INST ANCE	(0040,0270)	
		UNSCHEDULED	N/A	N/A	EMPTY	(0040,0270)	
		APPEND	SRC_INST ANCE	(0040,0275) ^(a) N/A ^(b)	SRC_INST ANCE	(0040,0270)	
Scheduled Procedure Step ID	(0040,0009)	SCHEDULED	MWL	(0040,0275) ^(a) N/A ^(b)	SRC_INST ANCE	(0040,0270)	
		UNSCHEDULED	N/A	N/A	EMPTY	(0040,0270)	
		APPEND	SRC_INST ANCE	(0040,0275) ^(a) N/A ^(b)	SRC_INST ANCE	(0040,0270)	
Patient Name	(0010,0010)	SCHEDULED	MWL	ROOT	SRC_INST ANCE	ROOT	
		UNSCHEDULED	USER	ROOT	SRC_INST ANCE	ROOT	
		APPEND	SRC_INST ANCE	ROOT	SRC_INST ANCE	ROOT	
Patient ID	(0010,0020)	SCHEDULED	MWL	ROOT	SRC_INST ANCE	ROOT	
		UNSCHEDULED	USER	ROOT	SRC_INST ANCE	ROOT	

**GE HEALTHCARE
STATEMENT**

**VIVID AND ECHOPAC V206
CONFORMANCE**

DIRECTION DOC2652554 REV 3

Attribute Name	Tag	Scenario	Image		MPPS		Comments
			Value Source	Destination Tag	Value Source	Destination Tag	
		APPEND	SRC_INST ANCE	ROOT	SRC_INST ANCE	ROOT	
Patient Birth Date	(0010,0030)	SCHEDULED	MWL	ROOT	SRC_INST ANCE	ROOT	
		UNSCHEDULED	USER	ROOT	SRC_INST ANCE	ROOT	
		APPEND	SRC_INST ANCE	ROOT	SRC_INST ANCE	ROOT	
Patient Sex	(0010,0040)	SCHEDULED	MWL	ROOT	SRC_INST ANCE	ROOT	
		UNSCHEDULED	USER	ROOT	SRC_INST ANCE	ROOT	
		APPEND	SRC_INST ANCE	ROOT	SRC_INST ANCE	ROOT	
Admission ID	(0038,0010)	SCHEDULED	MWL	ROOT	SRC_INST ANCE	ROOT	
		UNSCHEDULED	EMPTY	ROOT	SRC_INST ANCE	ROOT	
		APPEND	SRC_INST ANCE	ROOT	SRC_INST ANCE	ROOT	

^(a) for use in Image IODs

^(b) for use in Evidence Documents

E. Index of Tables

Table 0.1 Network Services.....	i
Table 0.2 Media Services.....	i
Table 2.7.1 Supported Specific Character Sets.....	23
Table 3.3.1 Exam Save – Supported SOP Classes and Transfer syntaxes	27
Table 3.3.2 DICOM Browse – SUPPORTED SOP CLASSES AND TRANSFER SYNTAXES.....	28
Table 3.3.3 Exam Read – Supported SOP Classes and Transfer syntaxes	29
Table 4.3.1 Mapping of DICOM entities to Vivid v206 entities	32
Table 5.3.1 Mapping of DICOM Entities to Vivid v206 Entities.....	37
Table 6.3.1 Mapping of DICOM Entities to Vivid v206 Entities.....	40
Table 7.3.1 Mapping of DICOM Entities to Vivid v206 Entities.....	43
Table 7.5.1 SR Root Templates	44
Table 8.3.1 Mapping of DICOM Entities to Vivid v206 Entities.....	48
Table 9.3.1 Mapping of DICOM Entities to Vivid v206 Entities.....	53
Table 9.4.1 Supported C-FIND Query Parameters for Modality Worklist-SCU	54
Table 11.2.1 STORAGE COMMITMENT MODULE	63
Table 11.2.2 STORAGE COMMITMENT REQUEST - ACTION INFORMATION... ..	64
Table 11.2.3 STORAGE COMMITMENT RESULT - EVENT INFORMATION.....	65
Table 12.2.1 Basic Grayscale Print Management Meta SOP Class.....	67
Table 12.2.2 Basic Color Print Management Meta sop class	67
Table 12.3.1 DIMSE Service Group.....	68
Table 12.3.2 DIMSE Service Group.....	69
Table 12.3.3 DIMSE Service Group.....	69
Table 12.3.4 DIMSE Service Group.....	70
Table 12.3.5 DIMSE Service Group.....	70
Table 12.4.1 Film Session IOD Modules.....	71
Table 12.4.2 Basic Film Box IOD Modules	71
Table 12.4.3 Basic Image Box IOD Modules.....	71
Table 12.4.4 Printer IOD Modules	72
Table 12.5.1 SOP Common Module Attributes.....	72
Table 12.5.2 Basic Film Session Presentation Module Attributes.....	73
Table 12.5.3 Basic Film Session Relationship Module Attributes	73
Table 12.5.4 Basic Film Box Presentation Module Attributes	74
Table 12.5.5 Basic Film Box Relationship Module Attributes.....	75
Table 12.5.6 Image Box Pixel Presentation Module Attributes	76
Table 12.5.7 Printer Module Attributes	77
Table 13.2.1 Mapping of DICOM Entities to Vivid v206 Entities.....	79
Table 13.3.1 Supported STUDY level C-FIND Matching Keys for Study Root Q/R Model -SCU	80
Table 17.2.1 External Network Requirements.....	325
Table 17.4.1 Secure Use and User Identity Profiles.....	325
Table 17.4.2 Secure Transport Connection Profiles	326

**GE HEALTHCARE
STATEMENT**

DIRECTION DOC2652554 REV 3

Table A.1.1 Shared modules	329
Table A.1.2 Shared private modules	344
Table A.1.3 Shared values and code sets	345
Table A.2.1 Ultrasound Image IOD Modules	347
Table A.3.1 Ultrasound Multi-frame Image IOD	348
Table A.3.2 Ultrasound Multi-frame Image IOD Specific Modules	348
Table A.4.1 Secondary Capture Image IOD	349
Table A.4.2 Secondary Capture Image IOD Specific Modules	350
Table A.5.1 Basic Directory IOD Modules	351
Table A.5.2 Basic Directory IOD Specific Modules	351
Table A.5.3 Attributes specific to record type PATIENT	353
Table A.5.4 Attributes specific to record type STUDY	353
Table A.5.5 Attributes specific to record type SERIES	355
Table A.5.6 Attributes specific to record type IMAGE	357
Table A.5.7 Attributes specific to record type SR DOCUMENT	357
Table A.6.1 Comprehensive SR IOD modules	358
Table A.6.2 Comprehensive SR Specific Modules	358
Table D.1.1 Attribute mapping	366