
To

HYDROPHOBIC

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# enVista® TORIC

ONE-PIECE HYDROPHOBIC ACRYLIC  
TORIC IOL

Ref **MX60TPxxxx**

## MATERIAL

Glistening-Free Hydrophobic Acrylic

4 % water content

UV-blocker

Refractive index: 1.54

## DESIGN

One-Piece, Aberration-Free Aspheric Optic

Step-vaulted haptics; Modified C-loop haptics

360° posterior square edge

Fenestrated haptics

Optic diameter: 6.0 mm

Overall diameter: 12.5 mm

## DIOPTER RANGE

+6.0 D to +30.0 D in 0.5 D increments

Cylinder powers-IOL plane: +1.25 D / +2.00 D /

+2.75 D / +3.50 D / +4.25 D / +5.00 D / +5.75 D

Cylinder powers-corneal plane: +0.90 D / +1.40 D /

+1.93 D / +2.45 D / +2.98 D / +3.50 D / +4.03 D

## INJECTOR

Reusable BLIS-R

with single-use cartridge BLIS-X0 from 0.0 D to +34.0 D (10/box)  
and **BLIS-X1** from +10.0 D to +34.0 D (10/box)

Recommended incision size: 2.2 mm WAT

INJ100 (10/box)

Recommended incision size: 2.2 mm WAT

## CONSTANTS\*

Immersion A-scan  
and IOL Master

A Constant SRK/T: 119.1

ACD: 5.61

Surgeon Factor: 1.46

Haigis Constant:  $a_0$ : 1.923 /  $a_1$ : 0.4 /  $a_2$ : 0.1

Applanation  
A-Scan

A Constant: 118.7

ACD: 5.37

Surgeon Factor: 1.62

\* Constants are estimates only.  
It is recommended that each surgeon  
develops his / her own values.  
Latest update: July 2015

TORIC

HYDROPHOBIC

MINI-INCISION


## enVista® TORIC INTRAOCULAR LENS

Lock in superior rotational stability<sup>1</sup>

Unique haptics are designed to secure  
a predictable astigmatism correction

- ✓ Glistening-Free Hydrophobic Acrylic
- ✓ The ideal combination of stable performance and predictability
- ✓ Aberration-Free Aspheric
- ✓ Fenestrated, step-vaulted haptics with 56° Contact angle and square posterior edge optic are designed to optimize 360° Capsular contact<sup>2</sup>
- ✓ 360° posterior square edge with haptic-optic junction designed to minimise PCO
- ✓ Polished for a smooth optic surface


Unique fenestrated, step-vaulted haptics with 56° contact angle are designed to maximise stability

- 91% of patients had  $\leq 5^\circ$  rotation from day of surgery to 6 months<sup>1</sup>
- 3° absolute mean rotation at 6 months<sup>1</sup>
- 0.28 mm mean decentration<sup>1</sup>

1. Packer M and al. Safety and effectiveness of a glistening-free single-piece hydrophobic acrylic intraocular lens (enVista). Clinical Ophthalmology 2013;7:1905-1912

2. Nishi O, Nishi K, Osakabe Y. Effect of intraocular lenses on preventing posterior capsule opacification: design versus material. J Cataract Refract Surg. 2004;30(10):2170-2176

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- Hoffer KJ. The Hoffer Q formula: a comparison of theoretic and regression formulas. Journal of Cataract and Refractive Surgery, Vol. 19, pp. 700-712, 1993; ERRATA, Vol. 20, pp. 677, 1994.

- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY, Ruiz RS. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery, Vol. 14, pp. 14-24, 1988.
- Norby NES. Unfortunate Discrepancies, Letter to the Editor and Reply by Holladay JT. Journal of Cataract and Refractive Surgery, Vol. 24, pp. 433-434, 1998.
- Olsen T, Olesen H, Thim K, and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery, Vol. 18, pp. 280-285, 1992.
- Retzlaff JA, Sanders DR, Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery, Vol. 16, pp. 333-340, 1990; ERRATA, Vol. 16, pp. 528, 1990.
- Haigis W: The Haigis Formula. In: Intraocular lens power calculations. H. John Shammas (eds), Slack Incorporated, Thorofare, NJ, USA, pp. 39-57, 2004.

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Bausch & Lomb Incorporated  
Rochester, NY 14609 USA

Bausch & Lomb Incorporated  
106 London Road, Kingston-upon-Thames, KT2 6TN, UK

Gamybos vieta: Bausch & Lomb Incorporated, 21 Park Place Blvd. North Clearwater, FL 33759 JAV

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

The enVista™ Toric intraocular lens (IOL) is a single-piece ultra-violet absorbing posterior chamber intraocular lens designed to replace the natural crystalline lens in adult patients in whom the cataractous lens has been removed. The enVista Toric IOL has an aspheric optic and is designed to be free of spherical aberration. Clinical studies have not been conducted with the enVista Toric IOL to assess the effect of the aspheric surface on spherical aberration, visual acuity, or contrast sensitivity.

PHYSICAL CHARACTERISTICS OF ENVISTA™ TORIC MODEL MX60T

Lens / Haptic Material	Hydrophobic acrylic (hydroxyethyl methacrylate (HEMA)-polyethylene glycol phenyl ether acrylate (poly(EG)PEA)-styrene copolymer, crosslinked with ethylene glycol dimethacrylate)
Material Characteristics	Index Of Refraction: 1.54 @ 35°C; Specific Gravity: 1.19 g/ml
Optic Type / Powers	Aspheric / +6.0 to +30.0 Diopters in 0.5 Diopter increments (SE – Spherical Equivalent)
IOL Cylinder Powers	1.25, 2.00, 2.75, 3.50, 4.25, 5.00, and 5.75 Diopters (CYL)
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 12.5 mm; Haptic Angle: 0°
Spectral Transmittance	<b>Ultraviolet:</b> 10% transmittance at 363 nanometers for +20.0 diopter IOL. See figure 1 with chart's X value = Wavelength (nm) and Y value = % Transmittance; chart compares the transmittance curve of an enVista Toric MX60T Lens to a 53 Year Old Human Lens. <b>A</b> = +20 Diopter enVista Toric MX60T Lens and <b>B</b> = 53 Year Old Human Lens. <b>NOTE:</b> Light transmittance values for an IOL material may vary slightly depending on the method of measurement. <b>Reference:</b> 53 year old human lens data from Boettner, E.A. and Welter, J. R., "Transmission of the Ocular Media," Investigative Ophthalmology, 1:776-783, 1962

**INDICATIONS**  
Indicated for primary implantation for the visual correction of aphakia and pre-existing corneal astigmatism in adult patients in whom the cataractous lens has been removed. The lens is intended for placement in the capsular bag.

**WARNINGS**  
Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

1. Recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the intraocular lens may affect the ability to observe, diagnose, or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
4. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage to the endothelium during implantation.
6. Suspected microbial infection.
7. Children under the age of 2 years are not suitable candidates for intraocular lenses.
8. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.
9. Rotation of enVista Toric away from the intended axis can reduce their astigmatic correction. Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, lens positioning should occur prior to capsule fibrosis and lens encapsulation.

**PRECAUTIONS**

1. Do not attempt to resterilize the lens as this can produce undesirable side effects.
2. Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage (such as the loss of saline storage solution, or the presence of salt crystallization).
3. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
4. Do not store the lens at a temperature greater than 43°C (110°F). DO NOT FREEZE. Do not autoclave the intraocular lens.
5. Do not reuse the lens. It is intended for permanent implantation. If explanted, sterility and proper function cannot be assured.
6. The safety and effectiveness of the enVista Toric IOL have not been substantiated in patients with preexisting ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting a lens in a patient with one or more of these conditions. Physicians considering lens implantation in such patients should explore the use of alternative methods of aphakic correction and consider lens implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.

- Before Surgery**
- Retinal conditions or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this lens.
  - Amblyopia
  - Clinically severe corneal dystrophy (e.g., Fuchs’)
  - Rubella, congenital, traumatic or complicated cataracts
  - Extremely shallow anterior chamber, not due to swollen cataract
  - Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye (e.g. iritis or uveitis).
  - Aniridia
  - Iris neovascularization
  - Glaucoma (uncontrolled or controlled with medication)
  - Microphthalmos or macrophthalmos
  - Optic nerve atrophy
  - Previous corneal transplant
  - Pre-existing ocular conditions which may negatively impact stability of the implant.
- During Surgery**
- Mechanical or surgical manipulation required to enlarge the pupil
  - Vitreous loss (significant)
  - Anterior chamber bleeding (significant)
  - Uncontrollable positive intraocular pressure
  - Complications in which the IOL stability could be compromised
7. Patients with preoperative problems such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, and chronic drug misuse may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from lens implantation when such conditions exist.
  8. A high level of surgical skill is required for intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.
  9. As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block, cyclitic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: lens repositioning, lens replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.
  10. Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens, to minimize the possibility of lens rotation causing misalignment of enVista Toric from the intended axis placement.

**MODE OF ACTION**  
enVista Toric IOLs are intended to be folded and placed in the posterior chamber of the eye after the cataractous crystalline lens has been removed. These IOLs have biconvex aspheric toric optics with cylinder axis marks to denote the flat meridian, which should be aligned with the postoperative steep corneal meridian to correct the astigmatism. The axis marks (lines) are on the posterior surface of the enVista Toric lens, denoting the flat meridian of the IOL. Optic.

- CALCULATION OF LENS POWER**  
The recommended A-constant listed on the lens carton is intended for use with axial length measurements obtained by optical biometry. Use of axial length measurements by other techniques (e.g. Applanation A-scan) will normally require a different lens constant. This number is a guideline only and is based on an evaluation of clinical data obtained using the IOL Master. The physician should determine preoperatively the power of the lens to be implanted. Lens power calculation methods are described in the following references:
- Hoffer K.J. The Hoffer Q formula: a comparison of theoretic and regression formulas, Journal of Cataract and Refractive Surgery Vol. 19, pp. 700-712, 1993; ERRATA, Vol. 20, pp. 677, 1994.
  - Holladay J.T, Musgrove KH, Prager T.C, Lewis J.W, Chandler T.V, Ruiz R.S. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery, Vol. 14, pp. 17-24, 1988.
  - Norrby NES. Unfortunate Discrepancies, Letter to the Editor and Reply by Holladay J.T. Journal of Cataract and Refractive Surgery, Vol. 24, pp. 433-434, 1998.
  - Olsen T, Olesen H, Thim K, and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery, Vol. 18, pp. 280-285, 1992.
  - Rietveld J.A, Sanders DR, Koffa M.C. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery, Vol. 16, pp. 333-340, 1990; ERRATA, Vol. 16, pp. 528, 1990.
  - Haigis W. The Haigis Formula. In: Intraocular lens power calculations. H. John Shammas (eds), Slack Incorporated, Thorofore, NJ, USA, pp. 39-57, 2004.

- DIRECTIONS FOR USE**
1. Prior to implanting, examine the lens package for type, power, and proper configuration.
  2. Open the peel pouch and remove the vial in a sterile environment.
  3. Remove the lid from the vial.
  4. With a pair of smooth forceps, remove the lens from the vial by gently grasping the lens haptic.
  5. Rinse the entire lens with sterile balanced salt solution or sterile normal saline.
  6. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.
  7. The lens may be soaked in sterile balanced salt solution until ready for implantation.
  8. Armix, Armix Plus, or OcuCoat viscoelastic should be used for lubrication of the delivery system when inserting the lens.
  9. Bausch + Lomb recommends using a Bausch + Lomb approved delivery system.
  10. There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.

**OVERVIEW OF CLINICAL STUDIES**  
Clinical studies have been conducted on the enVista single-piece IOL (model MX60) and the parent exact X-60 three-piece IOL (model X-60). The results of clinical investigations provide reasonable assurance that the Model MX60 IOL is safe and effective for the visual correction of aphakia following cataract extraction. For the enVista MX60 clinical study, all subjects in the safety analysis set were evaluated for IOL glisterings at Form 3 and Form 4 visits. IOL glisterings were evaluated via retroillumination slit lamp examination utilizing a photographic grading scale provided in the protocol. The grading scale consisted of (in order of severity), "none, grade 0 (trace), grade 1, 2, 3, or 4." No glisterings of any grade were reported for any subject at any visit in the clinical study. Also, rotational stability between Form 3 (30 to 60 days post-op) and Form 4 (120 to 180 days post-op) was evaluated on 113 eyes. 100% percent of eyes reported a rotation of ≤ 5°, which exceeds the American National Standards Institute (ANSI) standard of 90 percent of subjects with ≤ 5° of rotation at 2 consecutive visits at least 3 months apart.

**HOW SUPPLIED**  
The lens is individually packaged in a sterile vial (containing a 0.9% saline solution), within a peel pouch, and should only be opened under sterile conditions. A patient card and self-adhesive labels are supplied to provide traceability of the lens. The package is sterilized by gamma irradiation.

**EXPIRATION DATE**  
Sterility is guaranteed unless the pouch is damaged or opened. The expiration date on the lens package is the sterility expiration date. This lens should not be implanted after the indicated sterility expiration date.

**ADVERSE EVENT REPORTING**  
Adverse events and/or potentially sight threatening complications that may be regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be reported within five (5) days to Bausch + Lomb Incorporated. This information is being requested from all surgeons in order to document potential long-term effects of intraocular lens implantation.

Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with intraocular lenses. These problems may be related to a specific lot of lenses or may be indicative of long-term effects associated with these lenses or with IOLs in general. If you wish to report a problem, please call Bausch + Lomb at 1-800-338-2020.

**PATIENT REGISTRATION INSTRUCTIONS AND REPORTING REGISTRATION**  
Each patient who receives an enVista IOL must be registered with Bausch + Lomb at the time of lens implantation. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to Bausch + Lomb. Patient registration is essential and will assist Bausch + Lomb in responding to adverse reaction reports and/or potentially sight-threatening complications. An implant identification card is supplied in the lens package and must be given to the patient.

**RETURNED GOODS POLICY**  
All lenses being returned must be accompanied by a returned goods authorization number issued by Bausch + Lomb Customer Service. Call 1-800-338-2020 for return authorization and full policy information.

**WARRANTY**  
Bausch & Lomb Incorporated warrants that the intraocular lens, when delivered, will conform to all applicable laws and the manufacturer's then current version of the published specifications for such intraocular lens in all material respects and will be free from defects in material and workmanship.  
BAUSCH & LOMB INCORPORATED DISCLAIMS ALL OTHER WARRANTIES, EXPRESS, IMPLIED OR BY OPERATION OF LAW, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. BAUSCH & LOMB INCORPORATED SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, INDIRECT OR EXEMPLARY DAMAGES OF ANY KIND, DIRECTLY OR INDIRECTLY ARISING FROM THE PURCHASE OR USE OF THIS PRODUCT EVEN IF BAUSCH & LOMB INCORPORATED HAD BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, DAMAGE OR EXPENSE.

Symbol	Description	Symbol	Description
PC	Posterior Chamber	Ø <sub>b</sub>	Body Diameter (Optic Diameter)
UV	Ultraviolet	Ø <sub>l</sub>	Overall Diameter (Overall Length)
D	Diopter	SN	Serial Number
	Member Green Dot Scheme		Caution: Federal (US) law restricts this device to sale by or on the order of a physician

Bausch & Lomb Incorporated  
Rochester, NY 14609 USA  
 Bausch & Lomb Incorporated  
106 London Road, Kingston-upon-Thames, KT2 6TN, UK  
Manufacturing Site: Bausch & Lomb Incorporated, 21 Park Place Blvd. North, Clearwater, FL 33759 USA  
U.S. Patents: 6281319 and 6635731. "i"/"m" are trademarks of Bausch & Lomb Incorporated or its affiliates.  
All other brand/product names are trademarks of their respective owners. © Bausch & Lomb Incorporated.