# COULTER DxH Diluent

**REF** 628017

## TABLE OF CONTENTS

English . . . . .	2
Français (FR) . . . . .	4
Deutsch (DE) . . . . .	6
Italiano (IT) . . . . .	8
Español (ES) . . . . .	10
Português Portugal (PT-PT) . . . . .	12
Dansk (DA) . . . . .	14
Svenska (SV) . . . . .	16
Norsk (NO) . . . . .	18
Ελληνικά (EL) . . . . .	20
日本語 (JP) . . . . .	22
Lietuviškai (LT) . . . . .	24
Magyar (HU) . . . . .	26
Polski (PL) . . . . .	28
Čeština (CZ) . . . . .	30
Slovenčina (SK) . . . . .	32
한국어 (KO) . . . . .	34
Türkçe (TR) . . . . .	36
Русский (RU) . . . . .	38
Hrvatski (HR) . . . . .	40
Български (BG) . . . . .	42
中文 (ZH-TW) . . . . .	44
Română (RO) . . . . .	46
Србија (SR) . . . . .	48
Українська (UK) . . . . .	50
Português Brasil (PT-BR) . . . . .	52
Nederlands (NL) . . . . .	54
Tiếng Việt (VI) . . . . .	56

# COULTER DxH Diluent

**REF** 628017

**For *In Vitro* Diagnostic Use**  
**Rx Only in the U.S.A.**

## INTENDED USE

For use as an isotonic buffered diluent in conjunction with a cyanide-free lytic agent for counting and sizing blood cells on all UniCel DxH Coulter Cellular Analysis Systems. Use only with COULTER DxH Cell Lyse. Refer to your instrument product manuals and/or online help, as applicable. For use as a rinsing agent on the UniCel DxH Slidemaker Stainers.

## INTENDED USER

This device is intended to be used by a laboratory professional.

## SUMMARY

Blood cell analysis comprises diluting a whole-blood sample with a solution that functions as a diluent. The diluent provides the ability to analyze portions of the diluted blood sample for different blood cell types and the measurement of hemoglobin. COULTER DxH Diluent also acts as a rinsing agent on all UniCel DxH Coulter Cellular Analysis Systems.

## PRINCIPLE

The chemical composition of COULTER DxH Diluent provides an isotonic solution with buffering capacity for cell component stabilization and measurement of hemoglobin. Biocides are included for product preservation against microbial growth. Alkaline salts provide a suitable osmolality so as not to adversely affect cell volume (iso-osmotic).

When used on the UniCel DxH Coulter Cellular Analysis Systems in conjunction with COULTER DxH Cell Lyse, COULTER DxH Diff Pack, and COULTER DxH Retic Pack, the diluent is used for erythrocyte, hemoglobin, platelet, NRBC and reticulocyte measurements, as well as, for the measurement and differentiation of leukocytes into five subpopulations – neutrophils, lymphocytes, monocytes, eosinophils, and basophils.

## REACTIVE INGREDIENTS

Imidazole ..... 0.1 - <0.3% w/w

## NON-REACTIVE INGREDIENTS

Buffer, Stabilizers, Preservatives

## WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINER.

GHS HAZARD CLASSIFICATION		
COULTER DxH Diluent	<b>DANGER</b> 	EUH208: May produce an allergic reaction.
		H316: Causes mild skin irritation.
		H360: May damage fertility or the unborn child.
		P201: Obtain special instructions before use.
		P280: Wear protective gloves, protective clothing and eye/face protection.
		P308+P313: IF exposed or concerned: Get medical advice/attention.
		Imidazole 0.1 - <0.3% w/w
		reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) < 0.0015% w/w

## STORAGE, STABILITY, AND DISPOSAL

- Store Coulter DxH Diluent at 2-40°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

## MATERIALS PROVIDED

- COULTER DxH Diluent

## MATERIALS REQUIRED BUT NOT PROVIDED

- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the UniCel DxH Coulter Cellular Analysis Systems IFU.
- System HELP or Instructions for Use.

## PREPARATION

Coulter DxH Diluent is ready to use. Replace the reagent as directed in your instrument product manuals and/or online help.

**IMPORTANT:** If product has been partially or completely frozen, allow product to warm to room temperature. Invert the container 16 times to ensure complete mixing prior to placement on the instrument. Install and prime the diluent as directed in your instrument product manuals and/or online help. Verify background counts are acceptable before analyzing patient samples.

**PRODUCT AVAILABILITY**

COULTER DxH Diluent - 1 x 10 L

**REF** 628017

**TRADEMARKS**

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

**ADDITIONAL INFORMATION**

Glossary of Symbols is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs) by searching for part number C07163.

May be covered by one or more pat. - see [www.beckmancoulter.com/patents](http://www.beckmancoulter.com/patents)

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

**SDS** Safety Data Sheet is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs)

**NOTICE TO USER**

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

**REVISION HISTORY**

Revision AA, 03/2022

- Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes to product or instructions for use.

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Beckman Coulter Ireland Inc.,
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.