COULTER DxH Diff Pack

REF 628020

TABLE OF CONTENTS

English	2
Français (France) (fr-fr)	4
Deutsche (de-de)	6
Italiano (it-it)	8
Español (es-es)	10
Português (Portugal) (pt-pt)	12
Dansk (da-dk)	14
Svenska (sv-se)	16
Norsk (nb-no)	18
Suomi (fi-fi)	20
Ελληνικά (el-gr)	22
日本語 (ja-jp)	24
Lietuvių (lt-lt)	26
Magyar (hu-hu)	28
Polski (pl-pl)	30
Čeština (cs-cz)	32
Slovák (sk-sk)	34
한국어 (ko-kr)	36
Türkçe (tr-tr)	38
Русский (ru-ru)	40
Eestlane (et-ee)	43
Hrvatski (hr-hr)	45
Български (bg-bg)	47
中文 (台灣) (zh-tw)	49
Română (ro-ro)	51
Slovenščina (sl-si)	53
Srpski (sr-sp)	55
Latviešu (lv-lv)	57
Українська (uk-ua)	59
Português (Brasil) (pt-br)	61
Nederlands (nl-nl)	63
Tiếng Việt (vi-vn)	65
Қазақша (kk-kz)	67

COULTER DxH Diff Pack

REF 628020

**For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.**

INTENDED USE

For use on UniCel DxH Coulter Cellular Analysis Systems to lyse erythrocytes and preserve leukocytes. Refer to your instrument product manuals and/or online help, as applicable.

INTENDED USER

This device is intended to be used by a laboratory professional.

SUMMARY

An erythrolytic reagent and a leukocyte preservative produce a sample for performing a five-part WBC differential analysis using VCSn technology.

PRINCIPLE

WBC differential analysis begins by diluting a whole-blood sample with a lytic solution of formic acid that rapidly lyses the erythrocytes and reduces the cellular debris to an insignificant level without altering the leukocytes. The leukocyte preservative, consisting of an aqueous solution containing soluble salts, is added to stop and/or neutralize the lytic activity and restore the ionic balance to the blood, thus preserving the leukocytes in their near-native state for differentiation into five subpopulations (neutrophils, lymphocytes, monocytes, eosinophils, and basophils) using impedance, radio frequency, and VCSn technology.

REACTIVE INGREDIENTS

Coulter DxH Diff Pack contains the following:

COULTER Erythrolyse II Reagent

Wetting Agent 0.3-1.5g/L

Formic Acid 1.2mL/L

COULTER StabiLyse Reagent

Sodium Carbonate 6.0g/L

NON-REACTIVE INGREDIENTS

Stabilizers, Preservatives

WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINERS.

GHS HAZARD CLASSIFICATION		
Erythrolyse II Reagent	WARNING	EUH 208: May produce an allergic reaction.
		H316: Causes mild skin irritation.
		P332+P313: If skin irritation occurs: Get medical advice/attention.
		Formic Acid < 1% w/w
		reaction mass of: 5-chloro-2-methyl-4-isothiazolin -3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) < 0.0015% w/w
StabilYse Reagent		H402: Harmful to aquatic life.
		P273: Avoid release to the environment.
		Pyridine-2-thiol 1-oxide, Sodium Salt < 0.1% w/w

STORAGE, STABILITY, AND DISPOSAL

- Store COULTER DxH Diff Pack at 2-25°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.
- Closed container shelf life per stability study: 9 months

MATERIALS PROVIDED

- COULTER Diff Pack

MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the UniCel DxH Coulter Cellular Analysis Systems IFU.
- System HELP or Instructions for Use.

PREPARATION

Reagents in the COULTER DxH Diff Pack are ready to use. Replace the reagent kit as directed in your instrument product manuals and/or online help.

To avoid cross contamination when replacing the DxH Diff Pack, transfer Erythrolyse II and StabiLyse pickup tubes one at a time.

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Mix product by gentle inversion prior to placement on the instrument. Install and prime the reagent kit as directed in your instrument product manuals and/or online help. Verify background counts are acceptable before analyzing patient samples.

PRODUCT AVAILABILITY

COULTER DxH Diff Pack

628020

(Contains 1 x 1,900 mL Erythrolyse II and 1 x 850 mL StabiLyse reagents)

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs by searching for part number C07163.

May be covered by one or more pat. - see www.beckmancoulter.com/patents

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

Safety Data Sheet is available at beckmancoulter.com/techdocs

NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

REVISION HISTORY

Revision AA, 01/2022

- Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes to product or instructions for use.

Revision AB, 08/2023

- Updated GHS Hazard Classification table
- Updated STORAGE, STABILITY AND DISPOSAL section
- Added PRODUCT DESCRIPTION to Russian Translation
- Inserted Legal Manufacturer Phone Number
- Updated Legal Manufacturer Address Brazilian Translation

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.