COULTER Retic-X Cell Control

REF 628028

TABLE OF CONTENTS

English	2
Français (FR)	5
Deutsche (DE)	8
Italiano (IT)	11
Español (ES)	14
Português Portugal (PT-PT)	17
Dansk (DK)	20
Svenska (SE)	23
Norsk (NO)	26
Ελληνικά (GR)	29
日本語 (JP)	32
Lietuvių (LT)	35
Magyar (HU)	38
Polski (PL)	41
Čeština (CZ)	44
Slovák (SK)	47
한국어 (KR)	50
Türkçe (TR)	53
Русский (RU)	56
Hrvatski (HR)	59
Български (BG)	62
中文 (台灣) (ZH-TW)	65
Română (RO)	68
Srpski (Latinski) (SP)	71
Українська (UA)	74
Português Brasil (PT-BR)	77
Nederlands (NL)	80
Tiếng Việt (VN)	83

COULTER Retic-X Cell Control

REF 628028

For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.

INTENDED USE

Retic-X Cell Control is a hematology quality control material used to monitor the performance of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents.

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor instrument performance. This product can also be used to establish your own laboratory mean.

SUMMARY AND PRINCIPLE

Retic-X Cell Control is a reference product prepared from stabilized red blood cells. Retic-X Cell Control confirms and monitors instrument accuracy and precision performance by providing measurements for reticulocyte identification using VCSn technology.

REAGENTS

Retic-X Cell Control consists of treated, stabilized human erythrocytes in an isotonic medium. Retic-X Cell Control also contains a stabilized, reticulocyte-like component.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced materials: Human and Avian.

Each human donor unit used in preparation of this material was tested by an FDA approved method for the presence of the antibodies to Human Immunodeficiency Virus (HIV-1 and HIV-2) and Hepatitis C Virus (HCV) as well as for hepatitis B virus surface antigen and found to be negative (were not repeatedly reactive).

Because no test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV-1 and HIV-2), or other infectious agents are absent, this specimen/reagent should be handled at Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 1988.

Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (8/16/76).

GHS HAZARD CLASSIFICATION	
COULTER Retic-X Cell Control	May produce an allergic reaction.
	Neomycin Sulfate <0.1%
	Penicillin-G, sodium salt <0.1%

STORAGE, STABILITY, AND DISPOSAL

Retic-X Cell Control is shipped in a thermally insulated container designed to keep it cool. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. For opened vial stability, refer to TABLE OF EXPECTED RESULTS.

Storage of the control product in a cap down (inverted) position might require additional mixing for complete resuspension of cellular components. Over mixing causes hemolysis.

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

System performance should be considered acceptable when 95% of recovered values fall within the Expected Range. Slight hemolysis occurs over time through routine use and does not interfere with control performance. Failure to recover values within the Expected Range might indicate product deterioration in the absence of known instrument problems or gross hemolysis. A pink to light red color of the supernatant is normal and should not be considered deterioration of the product.

PERFORMANCE LIMITS

Refer to your Instructions for Use or System HELP for performance specifications for your instrument.

MATERIALS PROVIDED

- Retic-X Cell Control. (Refer to the Product Availability section of this insert for quantities and fill volumes).
- Table of Expected Results containing Assigned Values and Expected Ranges.


MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient COULTER reagents to complete analysis.
- Instructions for Use or System HELP.

INSTRUCTIONS FOR USE

1. Remove the COULTER Retic-X Cell Control tubes from the refrigerator and warm at AMBIENT TEMPERATURE for 15 minutes.
2. Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or Instructions for Use when a new lot of controls is received.
3. After warming, mix BY HAND as follows:
 - A. Roll tube slowly between the palms of the hands eight times in an upright position.

- B. Invert the tube and slowly roll it between the palms eight times.
- C. Gently invert the tube eight times.

4. Inspect the tube contents for complete resuspension. If necessary, repeat steps 3A and 3C until tube contents are uniformly suspended. **DO NOT OVERMIX. DO NOT** use any type of mechanical mixer.
5. Prepare and analyze COULTER Retic-X Cell Control in the same manner as a patient sample. Refer to the Instructions for Use or System HELP for sample preparation and instrument operation instructions.
6. Return the control tubes to the refrigerator **within 30 minutes**.
7. Compare instrument values to those given in the TABLE OF EXPECTED RESULTS.

ASSIGNED VALUES AND EXPECTED RESULTS

Refer to TABLE OF EXPECTED RESULTS.

ASSIGNED VALUES are determined on validated systems using specific COULTER reagents. The values are confirmed by multiple analysis of the control.

EXPECTED RANGES include variation between lots and between instruments, and represent 95% confidence limits for well maintained instrument systems.

Whether you use the Beckman Coulter Assigned Value or your own laboratory mean, the instrument is considered well maintained and operating correctly if:

- $\geq 95\%$ of the recovered values fall within the EXPECTED RANGE of the Beckman Coulter Assigned Value.
- Recovered values do not trend (five or more consecutive values drift in a single direction) OUTSIDE the EXPECTED RANGE.
- No more than three consecutive values fall outside the EXPECTED RANGE.

Before your current cell control lot(s) expire, perform the following on your new lot(s):

- Confirm that recovered values are within the TABLE OF EXPECTED RESULTS.
- OR
- Establish your own laboratory mean.

QUALITY CONTROL

This product is designed to monitor the reticulocyte capabilities of the operating system. Assuring the quality of the complete system, however, requires monitoring all reagents, controls and system indicators. Refer to the inserts for COULTER 6C Cell Control, COULTER LATRON CP-X Control, and Instructions for Use or System HELP. Inability to obtain any assigned value and/or electronic measurement could indicate poor system performance.

TROUBLESHOOTING PROCEDURE

The Instructions for Use or System HELP contain troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (US and Canada) or contact your local Beckman Coulter Representative.

To provide faster handling of your inquiry, please have the following information available when you call:

- Instrument model and serial number.
- Lot numbers, expiration dates and supporting data for all reagents, for control(s) in question and for other cell control lots currently in use.
- Previous cell control lot numbers and the data for these previous lots.
- Data history for COULTER LATRON CP-X Control.
- Data from your last IQAP Survey Report.

PRODUCT AVAILABILITY

COULTER Retic-X Cell Control:

628028

4 x 3.5 mL Level 1

4 x 3.5 mL Level 2

4 x 3.5 mL Level 3

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (US or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

Safety Data Sheet is available at beckmancoulter.com/techdocs

REVISION HISTORY

Revision AJ, 01/2018

- Added new Serbian language
- Added new Vietnamese language

Revision AK, 09/2018

- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

Revision AL, 05/2023

- Updated EC REP Address
 - Updated copyright date
-

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.