STORZ

KARL STORZ—ENDOSKOPE

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Flexibles Intubations-Videoendoskop – FIVE Serie 11302 BDx,
11301 BNx, 11303 BNx, 11304 BCx

INSTRUCTION MANUAL

Flexible intubation videoscope – FIVE Series 11302 BDx,
11301 BNx, 11303 BNx, 11304 BCx

MANUAL DE INSTRUCCIONES

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1 Important information

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper use, cleaning, disinfection and, if necessary, sterilization of the intubation video-scope. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

WARNING:

Risk of infection:

KARL STORZ instruments are not sterile when delivered, and must be cleaned, disinfected and, if necessary, sterilized prior to initial use and each subsequent reuse.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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Controls, displays, connectors and their uses

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3 Controls, displays, connectors and their uses

- ① Connecting cable
- ② Suction valve
- ③ Connection for leakage tester or pressure compensation cap (pressure compensation valve)
- ④ Connection for suction valve/irrigation adaptor 11301 CB1/11301 CE1/11301 CD1/091011-01 (only for endoscopes with working channel)
- ⑤ Working channel with lip valve 29100 (only for endoscopes with working channel)
- ⑥ Tube holder
- ⑦ Sheath
- ⑧ Distal tip
- ⑨ Deflection lever
- ⑩ Video recording Start/Stop
- ⑪ White balance
- ⑫ Single image capture
- ⑬ Monitor 8403 ZX
- ⑭ Second connection for videoscope
- ⑮ Primary connection for videoscope
- ⑯ Mains connection
- ⑰ Reset button
- ⑱ Holder for SD memory card
- ⑲ Status bar
- ⑳ Pilot lamp
- ㉑ On/Off switch button
- ㉒ Charging pilot lamp (orange)
- ㉓ Softkeys
- ㉔ Function symbols
- ㉕ HDMI port
- ㉖ USB port
- ㉗ Leakage tester 13242 XL
- ㉘ Pressure compensation cap 11025 E

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5 Intended use

5.1 Intended use Intubation endoscopes, flexible, CMOS with working channel

Flexible intubation endoscopes with working channel (CMOS) are designed for both endotracheal and nasal intubation. Flexible intubation endoscopes with working channel (CMOS) are designed for endoscopic visualization and image transmission as well as for creating an access to introduce additional instruments for diagnosis and treatment of the upper and lower airways. Flexible intubation endoscopes with working channel (CMOS) are intended for temporary use in interventions performed through natural orifices.

NOTE: The flexible intubation videoscope 11301 BNX can also be used with the C-MAC monitor 8402 ZX. All other endoscopes in the FIVE series are only compatible with the C-MAC monitor 8403 ZX.

Flexible intubation videoscopes are applied parts, which must only be used with KARL STORZ supply units which are intended for this purpose.

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5.2 Indications

The use of intubation endoscopes is indicated if, in the opinion of the attending physician, a fiberoptic endotracheal intubation or inspection of the oropharyngeal, supra- or subglottic region is indicated or PCT is performed.

This includes, amongst other things:

- Endotracheal intubation with a difficult respiratory tract
- Endoscopic inspection and monitoring of the mouth, nose and throat region.
- Monitoring of percutaneous tracheostomy
- Anesthesiological bronchoscopy

5.3 Contraindications

The use of the intubation endoscope is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Intubation endoscopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

No other contraindications are known.

5.4 Patient population

Men and women as well as children and adolescents. Use is not limited to a certain patient profile (sex, age, weight, etc.).

NOTE: The attending physician decides which patients to use the product on.

5.5 User qualification

KARL STORZ flexible intubation videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible laryngoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible intubation videoSCOPE. It is not intended as an introduction to the techniques of flexible laryngoscopy.

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6 Safety instructions

Please read this manual carefully and follow the instructions exactly. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this instrument. To further emphasize the words **Warning**, **Caution** and **Note**, they are preceded by a pictogram.

⚠️ WARNING:

‘Warning’ indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

⚠️ ATTENTION:

‘Attention’ indicates that particular service procedures or safety precautions must be followed to avoid any damage to the product.

NOTE:

‘Notes’ indicate special information about operating the instrument, or clarify important information.

6.1 Safety precautions when using flexible intubation videoscopes

Flexible intubation videoscopes must only be used according to the medical rules and procedures of endoscopy recognized for this method.

NOTE: The flexible intubation videoscopes do not offer any electric insulation against high frequency (HF) currents. The use of HF instrument sets can result in image faults and damage to the endoscope.

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⚠ ATTENTION:

Damage to product: Whenever the flexible intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.

- Damaged flexible intubation videoscopes or damaged accessories must not be used.

⚠ WARNING:

Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this.

The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.

⚠ WARNING:

Risk of infection:

- The flexible intubation videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.

⚠ WARNING:

The pressure compensation cap ② must be removed from the pressure compensation valve ③ during:

- endoscopic interventions
- immersion in liquid, e.g., disinfectant solution and for machine reprocessing

⚠ WARNING:

Risk from electromagnetic incompatibility: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the stated EMC instructions during installation and operation.

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NOTE: Damage to the flexible intubation videoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

⚠️ WARNING:

Always unplug the flexible intubation videoscope before carrying out any maintenance and cleaning procedures on it.

⚠️ WARNING:

Risk of injury: The flexible intubation videoscope must not be used during the discharge of a cardiac defibrillator.

- Remove the flexible intubation videoscope from the operation area before discharge takes place.

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⚠️ WARNING:

Risk of infection: In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices.

- To avoid contact and aerogen infections, medical devices must be decontaminated on-site and wrapped in protective film.

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⚠ WARNING:

Risk of infection:

- The flexible intubation videoscope is not sterile when delivered, and must therefore be properly cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

7.2 Basic equipment

The 11302 BDXX, 11301 BNXX, 11303 BNXX and 11304 BCXX sets consist of:

- 1x Flexible intubation videoscope
 - 11302 BDXX 4.0 x 65 cm or
 - 11301 BNX 5.5 x 65 cm or
 - 11303 BNX 5.5 x 65 cm or
 - 11304 BCX 6.5 x 65 cm
- 1x Case
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025E
- 1x Suction valve insert for single use, pack of 20, 11301 CE1
- 1x Cap for LUER lock irrigation port, black, pack of 10, autoclavable, 29100
- 2x Adaptors for cleaning the working channel in cleaning machines (black), 11301 CD1
- 50x Cleaning brush
- 1x Tube holder for video endoscopes 11301 CFX
- 1x Bronchoscope insertion tube, size 4, 10309
- 1x Bronchoscope insertion tube, size 2, 10310
- 1x Protective film 13990 SFN
- 1x Protective cap 8403 YZ

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7.4 External system components

Further information on external system components can be found in the instruction manual for the C-MAC monitor 8403 ZX (96076008D) and C-HUB II (96206529D).

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7.5 General information on use

The flexible intubation videoscope must be used according to the medical rules and procedures recognized for laryngoscopy.

NOTE: Any damage to the flexible intubation videoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

⚠ WARNING:

Risk of infection: The flexible intubation videoscope is not sterile when delivered, and must therefore be properly reprocessed prior to initial use and each subsequent reuse.

⚠ ATTENTION:

Damage to product:

- In order to avoid damage to the exterior, do not store or disinfect the flexible intubation videoscope together with sharp objects.

The following instructions outline the correct operation of the flexible intubation videoscope and the corresponding accessories. This is not an introduction to the techniques of laryngoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

The flexible intubation videoscope is used for endotracheal intubation and inspection of the oropharynx.

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7.6 Assembly before use

7.6.1 Fitting the sealing cap

Place the sealing cap (Art. no. 29100) on the proximal end of the working channel.

7.6.2 Fitting the suction valve

NOTE: The suction valve only needs to be fitted if the flexible videoendoscope has a working channel.

1. Place the suction valve (Art. no. 11301 CE1) onto the valve connecting piece in such a way that the tube connection is positioned above the function buttons and the unit's side nubs engage with the recesses on the connecting piece.

2. Using gentle pressure, rotate the suction valve in a clockwise or counter-clockwise direction until it noticeably engages.

For additional information on handling the suction valve, see the instruction manual for the suction valve 97000282.

7.6.3 Fitting the tube holder

Push the tube holder ⑥ 11301 CFX from distal onto the cone of the flexible intubation videoendoscope until the tube holder is securely in place.

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7.7 Checking before use

NOTE: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.

⚠ ATTENTION:

Damage to product: Damaged flexible intubation videoscopes or damaged accessories must not be used.

- Whenever the flexible intubation video scope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.

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7.7.1 Checking the suction

Suction can be triggered by pressing the suction valve ② (when the suction pump is connected and switched on).

NOTE: There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed (Art. no. 29100) to prevent additional air from being drawn in.

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7.7.2 Checking the insufflation

Connect the O₂ tube to the instrument channel:

- Use either the olive connector (Art. no. 600007) on the instrument channel or,
- if you are using the three-way double stopcock with LUER lock (Art. no. 6927691), connect the tube to the tube holder on the angled part of the adaptor. This keeps the straight opening of the adaptor free for instrumentation.

⚠ WARNING:

Observe the instruction manuals of the supply units used.

- Switch on the respective supply units.

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
⚠ ATTENTION:

Damage to product:

- The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

⚠ ATTENTION:

Damage to product:

- Do not squeeze the distal tip of the flexible intubation videoendoscope too tightly.

Every time the flexible intubation videoendoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Intubation videoendoscopes with sharp edges must not be used.

They may cause serious injury to the patient or user.

Also check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for loose connections.

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
7.8 Connecting and switching on

The flexible intubation videoscopes (FIVE) can be operated with the following monitors:

- C-MAC monitor 8403 ZX (11301 BNX also with C-MAC monitor 8402 ZX)
- also with additional monitors via the camera control unit C-HUB II 20 2903 01

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7.8.1 Connection to the C-MAC monitor

The flexible intubation videoscopes 11302 BDx, 11303 BNX and 11304 BCX can be used together with the C-MAC monitor 8403 ZX. The instruction manual for the monitor must be observed (see manual 96076008D).

NOTE: Older monitor models (8402 ZX/8401 ZX) are not compatible.

The flexible intubation videolescope 11301 BNX can be connected to both the C-MAC monitor 8402 ZX and the C-MAC monitor 8403 ZX.

⚠ WARNING:

Risk of injury: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.

⚠ WARNING:

Risk of injury: Optical radiation poses a risk of injury to eyes.

- Never look directly into the light output at the distal tip.

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Plug the connecting cable ① of the flexible intubation videoendoscope into one of the connection sockets ⑭ or ⑮ on the monitor.

The connecting cable ① can be connected or disconnected during operation. If necessary, connect the cord of the power supply up to the connection socket ⑯ on the C-MAC monitor and plug it into the power supply. The power cord can be fixed to the serial port using the union nut.

The battery is charged automatically as soon as the C-MAC monitor is connected to the power supply. This is indicated by the lighting up of the orange LED ⑳.

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7.8.2 Switching the flexible intubation videoscope On/Off

Switching On/Off: Pressing the on/off button ② on the C-MAC monitor switches on the flexible intubation videoscope with the integrated LED light source. Pressing the button a second time switches the unit off.

⚠ WARNING:

Danger of glare:

- Do not look directly into the light beam.

⚠ WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible intubation videoscope.

- Switch off the flexible intubation videoscope before putting it down, and do not place it in the vicinity of the patient (as per IEC 60601-1-1).

NOTE: The C-MAC monitor can be switched on or off by pressing the On/Off button ②. The flexible intubation video endoscope itself does not have an On/Off button.

Pressing and holding the On/Off button switches off the monitor together with the connected endoscope. The green indicator light ② goes out. If the On/Off button is actuated only briefly during operation, the monitor will switch to standby mode. This is indicated by the flashing indicator light on the monitor.

NOTE: When the monitor is switched off, the integrated light source in the flexible intubation videoscope also switches off.

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7.8.3 Functions of the C-MAC monitor 8403 ZX

A more detailed description of the operation and functions of the monitor 8403 ZX can be found in the corresponding instruction manual supplied with every monitor (see manual 96076008D).

If you lose the instruction manual, you can also download it from our homepage at www.karlstorz.com.

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7.8.4 Connection to the C-HUB II

The C-HUB II (Art. no. 20290320) is a plug & play device for processing and transmitting image and control signals to external devices in medical video endoscopy. The flexible intubation videoscope can be used in combination with the C-HUB II and an external monitor.

To do so, insert the connector of the flexible intubation videoscope into the serial port of the C-HUB II.

A more detailed description of the operation and functions of the C-HUB II can be found in the corresponding manual (manual 96206529D).

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NOTE: The flexible intubation videoscopes 11302 BDx, 11303 BNX and 11304 BCX can only be connected to the C-HUB II 20290301. Older models are not compatible. The flexible videoendoscope 11301BNX can be used with both the C-HUB II (Art. no. 20290101) and the C-HUB II (Art. no. 20290320).

NOTE: Observe the instruction manual for the C-HUB II (Art. no. 96206529D).

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1. Connect the power supply unit to the C-HUB II.


2. Connect the video connecting cable from the flexible intubation videoendoscope to the video input on the front of the C-HUB II.

3. Connect the monitor to the S-Video or HDMI output on the rear of the C-HUB II.

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Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

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
8 Handling

8.1 Holding the device correctly

Take the handle of the flexible intubation video-scope in your right or left hand and wrap your fingers around it to secure, as shown in the Figure. The thumb is located on the deflection lever ①. With this kind of control grip, the flexible intubation video-scope can be operated with one hand. (Fig. 1 Intubation video-scope with working channel; Fig. 2 Intubation video-scope without working channel).

8.2 Handling the sheath

⚠ ATTENTION:

Damage to product:

- The flexible intubation video-scope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the flexible intubation video-scope.

⚠ ATTENTION:

Damage to product:

- The sheath of the flexible intubation video-scope must never be clamped or coiled up too tightly. It is best to store the flexible intubation video-scope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406 AS).

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8.3 Handling the deflection

mechanism

The deflection lever ⑨ can be used to angle the distal tip of the flexible intubation videoendoscope. To do so, the deflection lever needs to be moved upwards (tip angled downwards) or downwards (tip angled upwards) using the thumb. When the distal tip is moved back into the neutral position, there will be a click. This ensures that the user is aware that the distal tip is back in the neutral position.

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

⚠️ WARNING:

Risk of injury: When inserting the flexible intubation videoendoscope into the patient, the distal tip of the flexible intubation videoendoscope must be in the straight, non-deflected position.

⚠️ WARNING:

Risk of injury: The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device!

- Actuate the deflection lever gently and carefully as far as the limit stop.
- Move the deflection lever smoothly and slowly!

⚠️ WARNING:

Risk of injury:

- Before the endotracheal tube is pushed via the sheath of the flexible intubation videoendoscope, the distal tip must be in the neutral position as the endoscope could otherwise be damaged.

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8.4 Handling the working channel

Do not force instruments or cleaning brushes through the channels. This can lead to damage to the working channel, especially if the tip of the sheath is deflected.

NOTE: Check instruments, cleaning brushes and other accessories for compatibility before inserting these in the working channel.

⚠ ATTENTION:

Damage to product: Unsuitable or damaged instruments may damage the working channel or cause the endoscope to malfunction.

- Check instruments, cleaning brushes and other accessories for sharp edges and corners before inserting them in the working channel.

When connecting and operating suction/irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kpa (400 mmHg/0.54 bar/7.8 psi).

- If no instrument is inserted, seal the proximal end of the working channel with a plug (Art. no. 29100).
- If an instrument is to be inserted, attach the lip valve (Art. no. 11301 CA) for single use to seal the proximal end of the instrument channel.
- The working channel of the flexible intubation videoscope means that specially developed instruments can be introduced. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage fiberscopes and flexible intubation videoscopes. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

NOTE: When using the lip valve, the instrument is subjected to increased resistance.

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- The flexible intubation videoscope must be in a straight, undeflected position when inserting accessories into the working channel.

⚠ ATTENTION:

Damage to product: The resistance that can be felt when inserting the instruments increases as the angle of deflection increases.

- *If the instruments are inserted with force, this can damage the flexible intubation videoscope or the instrument channel.*

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the flexible endoscope.

⚠ ATTENTION:

Damage to product: Accessories from other manufacturers may damage the flexible endoscope.

- *Only use accessories (guide probes, etc.) from KARL STORZ.*

8.5 Handling the suction valve

When the suction valve is assembled, a suction pump must be connected to the suction valve via a tube. Suction is activated upon pressing the valve button. If the suction valve is released, the suction is interrupted.

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8.6 Function buttons of the endoscope

The endoscope is equipped with three buttons. These are used to activate the functions displayed as symbols on the endoscope. The following functions are available:

	Video recording ⑩: Starts or stops the recording of a video stream on the SD memory card inserted in the receptacle ⑬ of the C-MAC monitor. The video recording symbol appears in the top left-hand corner of the live image during video recording.
	White balance ⑪: Enables the user to carry out a manual white balance. To do this, point the tip of the videoscope at a white surface and press the button for 1 second. NOTE: With the flexible intubation videoscopes, the performance of the white balance procedure is visualized by a progress bar in the menu bar. The white balance can also be performed via the C-MAC monitor menu.
	Single image capture ⑫: Saves the live image on the screen to the SD memory card of the C-MAC monitor.

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9 Reprocessing

NOTE: Remove the connector from the C-MAC monitor or the C-HUB II. Remove the tube holder, suction valves, lip valves and other accessories from the endoscope, and reprocess them separately.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

NOTE: Empty fields are not part of the validated process.

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⚠️ WARNING:

Risk of infection:

- Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

⚠️ WARNING:

Risk of infection from non-sterile delivery.

- Reprocess the device before each use. Check the medical device for visible soiling before use.

⚠️ WARNING:

Damage to product:

- Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

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⚠️ ATTENTION:

Damage to product:

- Never clean flexible endoscopes with metal brushes.

⚠️ ATTENTION:

Damage to product: Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

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Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Clean surface	27652
Working channel	27651 H or 110950-01 (11304 BCX), 110941-01 (11303 BNX), 110931-01 (11302 BDX), 110941-01 (11301 BNX)
Suction channel	11301 CD1, 27651 H or 110950-01 (11304 BCX), 110941-01 (11303 BNX), 110931-01 (11302 BDX), 110941-01 (11301 BNX)
Pressure compensation valve	13242 XA

9.2 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe-clean the exterior surfaces of the endoscope with a damp cloth.

Rinse out the working channel in order to prevent drying.

NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

9.3 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

9.3.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

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NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

9.3.2 Performing the leakage test – in liquid

NOTE: Damage to product: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

NOTE: Any lumina must be covered with no air bubbles by filling in a targeted manner.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid. Cleaning can then begin.

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9.4.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. The working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time (as specified by the manufacturer), the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

9.4.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time (as specified by the manufacturer), the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Di-aldehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

9.4.4 Drying

Finally, the surfaces and channels/lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

⚠ ATTENTION:

Damage to product:

- The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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9.5 Machine cleaning and disinfection

The following method for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

9.5.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channel is brushed until the brush head is visibly clean. Never brush back and forth within the channel.

Finally, it must be rinsed with cold water to ensure neutralization.

9.5.2 Connecting

In order to ensure effective machine cleaning and disinfection, the instrument must be connected to the suction and working channel of the endoscope in the washer and disinfectant.

⚠ ATTENTION:

Damage to product:

- During cleaning and disinfection the pressure compensation cap must not be in place on the pressure compensation valve.

9.5.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded.

The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

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9.6 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connection) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

⚠ ATTENTION:

Damage to product:

- The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.
- Afterwards, a functional check must be carried out.

NOTE: Check for a clear view and steerability.

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⚠️ WARNING:

Risk of infection:

- Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

⚠️ WARNING:

Risk of infection:

- Medical devices must be sterilized in dismantled condition.

⚠️ ATTENTION:

Damage to product:

- Flexible endoscopes must never be steam sterilized (autoclaved).

⚠️ ATTENTION:

Damage to product:

- The pressure compensation cap must be in place on the pressure compensation valve for sterilization and transport to ensure pressure compensation.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

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9.8.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

⚠ ATTENTION:

Damage to product: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: For detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the ‘STERRAD® Sterility Guide’ to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle (without channel)
- STERRAD® 100S Long Cycle with Boosters (with channel)
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

9.8.2 Hydrogen peroxide (H₂O₂) sterilization – STERIS® V-PRO™

NOTE: For detailed information, please consult the user handbook of the respective device.

NOTE: The ‘STERIS Low Temperature Reprocessing Device Compatibility Matrix’ can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® max Flexible Cycle
- V-PRO® 60 Flexible Cycle
- V-PRO® 1 Plus Lumen Cycle
- V-PRO® 1

9.8.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 60 minutes as per ANSI/AAMI/ISO Norm 10993-7.

EO sterilization is validated with the following parameters:

Pre-conditioning:

- 55°C
- 70% relative humidity
- 30 minutes

Cycle:

- 55°C
- 70% relative humidity
- 725 mg/L conc.
- 120 minutes
- 100% EO

Aeration:

- 12 hours
- 55 +/- 4 °C

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

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9.9 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

9.10 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

⚠ ATTENTION:

Damage to product:

- Do not expose flexible endoscopes to direct sunlight or excessive heat.

⚠ ATTENTION:

Damage to product:

- The distal tip of the flexible endoscope must be in the neutral position during storage.

⚠ ATTENTION:

Damage to product:

- Do not store the flexible endoscope in the case.

⚠ ATTENTION:

Damage to product:

- The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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⚠️ WARNING:

Risk of infection:

- Improperly reprocessed medical devices may pose a risk of infection for patients, users and third parties and may cause malfunction of the medical device. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

⚠️ WARNING:

Risk of infection:

- These medical devices are not supplied in sterile condition. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

⚠️ WARNING:

Risk of infection:

- When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

⚠️ ATTENTION:

Damage to product:

- When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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9.11.4 Manual precleaning

Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

Brushing the lumina

The working channels, lumina and hollow spaces must be precleaned under cold running water with the aid of a suitable brush.

9.11.5 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national regulations and the chemical manufacturer's instructions must be taken into account when using this method.

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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9.11.6 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, if necessary, the instrument must be connected up to the washer and disinfectant.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Disassembled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

9.11.8 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

9.11.9 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

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The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization in the fractionated prevacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C – 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

Or:

- 134 °C – 137 °C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

⚠ WARNING:

Risk of infection:

- Medical devices must be sterilized in disassembled condition.

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The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 50, 100S, 200 without booster
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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10.2 Flexible intubation videoscope 11301 BNX/11303 BNX

Direction of view:	0°
Angle of view:	100°
Deflection:	up 140°, down 140°
Working length:	65 cm
Outer diameter distal tip:	5.5 mm
Working channel:	2.1 mm
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight	400 g
Drip water protection:	IPX8
Classification:	Ila

Applied part: Provides Type BF protection when used in conjunction with a KARL STORZ camera control unit with a BF symbol on the connection port.

Operating/storage conditions

Operating conditions:	15°C ... +40°C 15–80% rel. humidity
Storage conditions:	-0°C ... +55°C 15–90% rel. humidity
Atmospheric pressure:	700–1080 hPa

ATTENTION:

The specifications in relation to atmospheric pressure can only be guaranteed if the pressure compensation cap is fitted to the pressure compensation valve during transportation.

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10.3 Flexible intubation videoscope 11302 BDX

Direction of view:	0°
Angle of view:	100°
Deflection:	up 140°, down 140°
Working length:	65 cm
Outer diameter distal tip:	4.0 mm
Working channel:	1.5 mm
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight	385 g
Drip water protection:	IPX8
Classification:	Ila

Applied part: Provides Type BF protection when used in conjunction with a KARL STORZ camera control unit with a BF symbol on the connection port.

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Operating/storage conditions

Operating conditions:	15°C ... +40°C
Storage conditions:	15–80% rel. humidity
	-0°C ... +65°C
	15–90% rel. humidity

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11 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

12 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

13 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

14 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see Chapter 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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15 Directive compliance

In accordance with the Medical Devices Directive
(MDD):

11302 BDX, 11301 BNX, 11303 BNX and
11304 BCX are Class II a medical devices.

This medical device bears a CE mark in accor-
dance with the Medical Device Directive (MDD)
93/42/EEC. A code number after the CE mark
indicates the responsible notified body.

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16 Electromagnetic Compatibility

The applied part described in this instruction
manual has been tested as a system with
the following basic devices. For relevant
electromagnetic compatibility (EMC) information,
refer to 'Electromagnetic Compatibility (EMC)
information' in the instruction manuals for these
basic devices:

Basic device	Instruction manual
C-MAC monitor 8403 ZX	96076008 D
C-HUB II 20 2903 20	96206529 D

The EMC warning statements, precautions,
notes and Emission/immunity limits specified
in the instruction manuals for the basic devices
also apply to the applied part described in this
instruction manual.

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17 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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18 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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
11025 E

20 Appendix Accessories

29100 Sealing cap,
black sealing cap, for LUER lock irrigation port,
black, pack of 10 autoclavable

11301 CD1 Irrigation adaptor,
for machine cleaning, reusable, for use with
flexible intubation videoscopes 11301 BNX,
11302 BDX, 11303 BNX and 11304 BCX

11301CE1 Suction valve,
for single use, pack of 20, for use with flexible
intubation videoscopes 11301 BNX, 11302 BDX,
11303 BNX and 11304 BCX

091011-20 Suction Valve,
sterile, for single use, package of 20, for use
with flexible Intubation Video Endoscopes
11301 BNX, 11303 BNX, 11302 BDX, 11304 BCX
and flexible Video Cystoscopes 11272 VH/VHU

10309/10310 Bronchoscope insertion tube,
size 4, with integrated mouthpiece, for single
use, sterile, insertion length 85 mm, made from
EVA, pack of 10, for use with flexible intubation
endoscopes

11301 CFX Tube holder,
for use with flexible intubation videoscopes
11302 BDX, 11301 BNX and 11303 BNX

11025 E Pressure compensation cap,
for venting during gas and plasma sterilization

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110284-10 Biopsy Forceps,
with Oval Jaws, coated, working length 120 cm,
sterile, for single use, package of 10, for use
with flexible endoscopes with working channels
as of diameter 2.0 mm

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110345-10* Grasping Forceps Type Pelican,
with large opening angle, diameter 2.4 mm,
working length 170 cm, sterile, for single use,
package of 10, for use with flexible endoscopes
with working channels as of diameter 2.8 mm

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110352-10* Grasping Forceps Type Alligator,
with ribbing, diameter 1.8 mm, working length
120 cm, sterile, for single use, package of 10,
for use with flexible endoscopes with working
channels as of diameter 2.0 mm

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110354-10* Grasping Forceps Type Alligator,
with ribbing, diameter 2.4 mm, working length
160 cm, sterile, for single use, package of 10,
for use with flexible endoscopes with working
channels as of diameter 2.8 mm

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KARL STORZ Endoscopy Canada Ltd.
17171 Millcreek Drive, Mississauga, Ontario
L5M 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksaa.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopy México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcadia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº 413, 2º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopy Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Slamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storstjärnsgård 14
12739 Skärholmén, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Tavallie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representative Office
Vilnius, Lithuania
Zalgirio St. 94, LT19300
Mobile: +370 685 07 048
Phone: +370 5 272 0400
E-Mail: info-lt@karlstorz.com
KARL STORZ Endoskopi Danmark AS
Skovlyoftvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TG, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com
KARL STORZ Endoscopie Nederland B.V.
Displayweg 2
3821 BT Amerfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com
KARL STORZ Endoscopie France S.A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraber Hauptstr. 148/7/61
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com
KARL STORZ Endoscopy Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopy Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-it@karlstorz.com
KARL STORZ Croatia d.o.o.
Capraska 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info-hrv@karlstorz.com

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com
KARL STORZ Endoskóp Magyarország Kft.
Toberek utca 2. sz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopy Romania srl
Str. Prof. Dr. Anton Colanin, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com
KARL STORZ Endoskope Greece M.E.P.E.*
Patriarhou Grigorios E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
E-Mail: info-gr@karlstorz.com

*Repair & Service Subsidiary
KARL STORZ Industriest**
Gedik is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy
000 KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
TOV LLC KARL STORZ Ukraine
Avenue Geroyiv Stalingrada Str. 2D, office 717
Kyiv, 04210/Ukraine
Phone: +38 095 000-895-0, +38-097-000-895-0,
+38 073 000-895-0
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG
Representation Office
Sabit Orudschow 1184, apt. 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: info-az@karlstorz.com
KARL STORZ ENDOSKOPE –
East Mediterranean and Gulf (Offshore) S.A.L.
Spark Tower 1st floor
Charles Helou St., Horch Tabet – Sin El Fil
Beirut, Lebanon
Phone: +961 1 501105, Fax: +961 1 501950
E-Mail: info@karlstorz-emg.com
KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TÜO KARL STORZ Endoskopi Kasachstan
Saryarka, b. BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415892
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in
KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 32, Unit 3230
399 Sukhumvit Road,
North Kongroey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 660 3669
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPIaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com
KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339868, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China

Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415892
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonweath Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ Endoscopy Philippines Inc.
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy New Zealand Ltd.
31 Morningside Drive Mt Albert
Auckland, 1025, New Zealand
PO Box 56 511, Dominion Rd
Auckland, 1446, New Zealand
Phone: +64 9 846 6044
Toll free: +64 508 84 84 84 (New Zealand only)
E-Mail: sales-nz@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113, Australia
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen
Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com