DIAPATH

Constant Innovation in Anatomic Pathology

INSTRUCTION FOR USE

010229

EN

IT

FR

ES

DE

May-Grunwald Giemsa acc. to Romanowsky for sections

Kit for special stains

Manufacturer

Diapath S.p.A.

Via Savoldini, 71 - 24057 Martinengo (BG) - Italy

Tel. 0363 986411

Fax 0363 948000

www.diapath.com

info@diapath.com

EMDN code: W01030708

UDI code: 08033895000932

Code	Test	Reagent	Reagent code	Format
010229	100	May Grunwald acc. to Pappenheim	C050AA	1x100 ml
		Giemsa acc. to Pappenheim	C045AB	1x30 ml
		Acetic water	G066AA	1x30 ml

Intended use

Single-use, non-sterile in vitro diagnostic device, for clinical use, indicated for typing cells of the haemolymphopoietic system, for the identification of parasites and of the endoplasmic reticulum.

The device allows expert personnel to make a diagnosis on the biological sample sectioned for research and diagnostic purposes.

The device has been designed for professional use only, in medical and veterinary laboratories and medical facilities, in manual mode, by direct contact or by immersion, or in combination with other devices.

For further information, warnings and precautions for use, please read the safety data sheet on the website.

Description of the reagents and any limitations of use

Device	Stability of closed product	Stability after opening	Ready for use	Storage temperature
010229	2 years	2 years	YES	15°C – 30°C
May Grunwald acc. to Pappenheim	2 years	2 years	YES	15°C – 30°C
Giemsa acc. to Pappenheim	2 years	2 years	YES	15°C – 30°C
Acetic water	2 years	2 years	YES	15°C – 30°C

Information relating to interfering substances that may affect the performance of the device

Not applicable

List of materials provided and specific materials required but not provided

The device is packaged, ready to use, in individual units to ensure product integrity and cleanliness.

The end result is only guaranteed by combining the above-mentioned components.

In terms of staining, the result obtained from the combination of the different components is intended for inspection under a microscope and interpretation of the result by experienced and qualified personnel.

The materials required for the analysis, but not provided, are:

- Distilled water
- Running water
- Reagents for tissue hydration and dehydration (xylene and substitutes, graded alcohols)
- Mounting medium
- Glass coverslip or substitutes (at user's discretion)
- Positive control slides

Additional equipment, if provided for in the protocol of use:

- Humidity chamber
- Oven

It is advisable to use Diapath's intended reagents. Any other reagent can be used, subject to the laboratory's approval.

Information on the combined use of the device

The kit is for manual use.

The device is compatible with the other devices with which it is used in combination and mentioned above.

Functional characteristics

Type of preparation	Section in paraffin - Cytological
Recommended fixative	Formalin
Control tissue	Bone marrow biopsies
Storage temperature	+15°/+30°C
Indicative time	1 hour and 5 minutes
Critical phase	Dilution of reagents
Additional equipment	Oven, humidity chamber

Test principle

The May Grunwald solution is composed of methylene blue eosinate and stains the nuclei blue and the basophilic cytoplasm red-pink.

The Giemsa acc. to Pappenheim solution is a complex mix, consisting of methylene blue salts and azure II, that increase the intensity of the nuclear staining and the ability to selectively highlight the cell elements.

The intensity of the stain may vary both depending on differentiation times and on the temperature used in the staining method.

Instructions for use

It is good practice to use a humidity chamber to avoid excessive section drying.

Diluted working solutions cannot be stored for a prolonged time and must be prepared before use. Formalin-fixed tissues are more basophilic and may require longer immersion times in the dyes or higher dye concentrations in the working solutions.

For paraffin-embedded histological samples start at step 1; For spray-fixed cytological samples, heat-fixed or air-fixed cytological samples or fixed by immersion in 95% alcohol start at step 2.

Staining protocol:

- 1. De-wax and immerse the sections in distilled water*
- 2. Dilute **May Grunwald acc. to Pappenheim** with distilled water in the ratio 1:4
- 3. Cover the sections with the prepared solution (step 2) for 20 minutes at 37°C
- 4. Wash with distilled water
- 5. Dilute 20 drops of **Giemsa acc. to Pappenheim** with 10 ml of distilled water
- 6. Cover the sections with the prepared solution (step 5) for 40 minutes at 37°C
- 7. Wash with distilled water
- 8. Cover the sections with **Acetic Water** for 30 seconds
- 9. Drain the slide and dry it with filter paper
- 10. Dehydrate with 1:1 solution of **absolute ethyl alcohol** and **acetone**
- 11. Clear and mount**

*Before proceeding with the staining protocol, it is necessary to deparaffinize and hydrate the sections . The preliminary hydration of the section must be carried out using specific reagents to remove the paraffin layer from the tissue and prepare the tissue for subsequent staining. Example of protocol:		**On completion of staining, the sections must be dehydrated, cleared, and mounted . Example of protocol:	
Oven at 60° for 20 min		70% alcohol: 1 min	
Xylene: 3 min		95% alcohol: 1 min	
Xylene: 3 min		99% alcohol: 2 min	
Xylene: 3 min		99% alcohol: 2 min	
99% alcohol: 3 min		Xylene: 3 min	
99% alcohol: 3 min		Xylene: 3 min	
95% alcohol: 3 min		Xylene: 3 min	
		Mounting medium compatible with histological slides	

It is good practice to follow all the manufacturer’s instructions in order to obtain the desired result, as an incorrect procedure may hinder the quality of the staining.


All operational procedures, prior to staining, i.e. fixation, dehydration, clearing, paraffin impregnation and finally cutting, are supervised by the operator and must be validated according to normal routine operations within the reference laboratory.

Please refer to the Bibliography section

It is neither a self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel, since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified and trained personnel.

Interpretation of the results

Basophilic cytoplasm:	Variants of Blue
Acidophilic cytoplasm:	Pink-Red
Nuclei:	Blue
Parasites:	Blue

Quality control

It is necessary to include a positive control slide with each staining run in order to verify the validity of the test performed. It is important that the positive control is treated similarly to the test samples, i.e. tested in accordance with the operational techniques of the test samples.

Analytical performance characteristics

Not applicable

Conditions for collection, handling and preparation of the sample

The preparation of the sample follows the procedures set forth in an anatomical pathology laboratory.
Do not handle until all safety precautions have been read and understood.

Storage

Store the device at a temperature between +15°C and +30°C.

Do not expose the kit to environmental conditions other than those prescribed, as failure to comply with storage conditions could lead to a change in the chemical and physical characteristics of the device, resulting in underperforming staining.

If properly stored and in unopened packaging, the product is stable as stated in the section *Description of the reagents and any limitations of use*. Do not use after the expiry date.

If the product is not stored as instructed, its yield may vary and it will have to be tested by the operator.

Stability when handling the device

The use of the device is under the absolute control of the skilled and qualified operator.

The environmental conditions in the laboratory must be controlled:

- Temperature: +15°C / +30°C
- Relative humidity ≤80%

Disposal instructions / Safe final disposal of the device

For specific instructions on how to safely dispose of the device after use, please refer to the local regulations of the country in which the device is used.

It is good practice to dispose of the kit as special waste in dedicated containers and in accordance with the internal disposal procedures of the use environment.

For information on the content of the reagents, please refer to the safety data sheet. If the product has expired and/or is unused, it must be disposed of in accordance with local waste regulations, based on the hazard classification on the label, after checking for any contamination. In some cases, analytical evaluation may be required so as to correctly classify the waste and determine the hazard characteristics to be attributed to it.

Warnings, precautions and restrictions on the use of the device

The device must only be used for the intended purpose.

For proper handling, it is recommended to:

- Wear personal protective equipment, such as a lab coat, gloves and goggles. Protect your eyes and face
- Handle the device under a fume hood
- In order to obtain the expected result, it is advisable to follow the instructions given in the protocol of use and the validated components within the kit
- Maintain a clean work environment and follow the cleaning protocols in the laboratory of use

Device malfunction or degradation

The device retains its properties throughout its estimated useful life if the storage requirements are guaranteed and complied with.

The device does not generate risks from any debris due to wear and tear, degradation products and processing residues.

In case of malfunction, isolate the device, replace it with a new one and dispose of the non-functioning kit. Should this occur, please contact the manufacturer by sending an e-mail to claims@diapath.com

Error	Possible cause	Possible corrective action
Ineffective staining and artefact generation causing a delay or missed diagnosis	Fraudulent use beyond the expiry date	Including appropriate warnings and markings on the label
	Inadequate storage of the device	Appropriate symbols and adequate device storage warnings on the label

Warnings from the user, with the need to report any serious accident that occurred in relation to the device

In case of a serious incident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Key of the symbols used in the labelling and instructions for use

Symbol	Symbol name	Explanatory notes
	Manufacturer	Indicates the manufacturer of the device
	Expiration date	Indicates the date after which the device must not be used
	Lot number	Indicates the manufacturer's batch necessary for batch identification
	Catalogue code	Indicates the code assigned by the manufacturer allowing its identification
	Distributor	Indicates the company distributing the medical device
	Not sterile	Indicates a device that has not undergone a sterilization process

	Keep away from light	Indicates a device that needs protection from light sources.
	Temperature limit	Indicates the temperature limits to which the device can be safely exposed
	Not reusable	Indicates a device that is to be used once, for single use only
	Refer to instructions	Indicates the need for the user to consult the instructions
	Attention	Indicates the need for the user to be cautious when operating the device or that the situation requires greater awareness to avoid undesirable consequences
	In vitro diagnostic medical device	Indicates a device intended to be used as an in vitro diagnostic medical device
	CE Marking	The CE mark certifies the conformity of the product to the applicable directives of the European Union
	Content sufficient for <n> tests	Indicates the number of tests that can be performed with the device
 XX.XX.XXXX	Manufacturing date	Indicates the date of manufacture of the medical device
	Unique identifier of the device	Indicates a vector containing the fixed and variable information related to the UDI

Identification of the substance/mixture and identification of the hazards

Substance or mixture	Hazard pictograms	Hazard statements	Cautionary advice
MAY GRUNWALD ACC. TO PAPPENHEIM		H225 Highly flammable liquid and vapour H301+H311+H331 Toxic if swallowed, in contact with skin or if inhaled. H370 Causes damage to organs.	P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection/... P301+P310: IF SWALLOWED: immediately call a POISON CENTRE or doctor/physician/... P308+P311 IF exposed or concerned: call a POISON CENTRE or doctor/physician/... P370+P378 In case of fire: use ... to extinguish.

			<p>P403+P233 Keep the container well closed in a well-ventilated place.</p> <p>P403+P235 Store in a well-ventilated place. Keep cool.</p>
GIEMSA ACC. TO PAPPENHEIM		<p>H225 Highly flammable liquid and vapour.</p> <p>H301+H311+H331 Toxic if swallowed, in contact with skin or if inhaled.</p> <p>H370 Causes damage to organs.</p>	<p>P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources.</p> <p>No smoking.</p> <p>P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection/...</p> <p>P301+P310: IF SWALLOWED: immediately call a POISON CENTRE or doctor/physician/...</p> <p>P308+P311 IF exposed or concerned: call a POISON CENTRE or doctor/physician/...</p> <p>P370+P378 In case of fire: use ... to extinguish.</p> <p>P403+P233 Keep the container well closed in a well-ventilated place.</p> <p>P403+P235 Store in a well-ventilated place. Keep cool.</p>
ACETIC WATER		<p>The product is not considered hazardous according to EC Regulation 1272/2008 (CLP).</p>	

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“Laboratory methods in histotechnology. Stain methods of the Armed Forces Institute of Pathology”. Prophet, Mills, Arrinton, Sobin. Editor McGRAW Hill. Washington D.C. 1968

“Theory and practice of Histological Techniques”. Edited by Bancroft JD and Gamble M. Churchill Livingstone. New York. 2002

Revision table

Rev.	Date	Purpose of the revision
1st issue	27/05/2024	Compliance with the new Regulation 746/17

(*) If the request is due to a change, please provide a detailed description of the reason for revising the document and the changes that have been made since the last revision; if available, please enclose recalls, non-conformities, NC CA PA statistics, etc.

You can find this document in the reserved area of the www.diapath.com website. If you want to download the document and keep up to date on the latest available version, you will need to log in after registering.

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