DIAPATH

Constant Innovation in Anatomic Pathology

INSTRUCTION FOR USE

C045X

EN

IT

FR

ES

DE

Giemsa acc. to Pappenheim

Nuclear stain

Manufacturer

Diapath S.p.A.

Via Savoldini, 71 - 24057 Martinengo (BG) - Italy

Tel. 0363 986411

Fax 0363 948000

www.diapath.com

info@diapath.com

EMDN code: W01030708

Code	Format	UDI code
C0450	100 ml	08033895003476
C0452	500 ml	08033895003971
C0453	1 l	08033895003988
C0454	2 l	08033895003506
C0455	2.5 l	08033895004008
C0456	5 l	08033895003520

Intended use

Disposable, non-sterile, in vitro diagnostic device for clinical use to identify parasites (e.g., *Helicobacter Pylori*) and mast cells in paraffin-embedded histological sections.

When used in association with the May Grünwald acc. to Pappenheim staining solution it can highlight the various component of blood on blood smears.

In terms of staining, the result obtained from the combination of the different components can be intended for inspection under a microscope and interpretation of the result by experienced and qualified personnel.

The device is intended for manual use or in combination with an automatic stainer, to be used on preparations previously prepared according to the methods set forth by the standard sample processing technique used by the laboratory.

Description of the reagents and any limitations of use

Glycerol	CAS No. 56-81-5	EC No. 200-289-5	
Methyl alcohol	CAS No. 67-56-1	EC No. 200-659-6	
Azure B	CAS No. 531-53-3	EC No. 208-510-7	C.I. 52010
Eosin Y	CAS No. 17372-87-1	EC No. 241-409-6	C.I. 45380
Methylene blue	CAS No. 61-73-4	EC No. 200-515-2	C.I. 42780

Device	Stability of closed product	Stability after opening	Ready for use	Storage temperature
C0450	2 years	2 years	YES	15°C – 30°C
C0452	2 years	2 years	YES	15°C – 30°C
C0453	2 years	2 years	YES	15°C – 30°C
C0454	2 years	2 years	YES	15°C – 30°C
C0455	2 years	2 years	YES	15°C – 30°C
C0456	2 years	2 years	YES	15°C – 30°C

Diapath S.p.A.

Via Pietro Savoldini, 71 - 24057 Martinengo (BG) Italy

T +39 0363 986411 | F +39 0363 948000

info@diapath.com | diapath.com

Rev.00 27/05/2024

Information relating to interfering substances that may affect the performance of the device

Not applicable

List of materials provided and specific materials required but not provided

The materials required for the analysis, but not provided, are:

- Distilled water
- Running water
- Reagents for tissue dehydration (xylene and substitutes, graded alcohols)
- May Grünwald acc. to Pappenheim solution
- Phosphate buffer pH 6.8
- Mounting medium
- Glass coverslip or substitutes (at user's discretion)
- Positive control slides

Additional equipment, if provided for in the protocol of use:

- Automatic colouring machine

Information on the combined use of the device

The device is intended for use within a staining protocol, in combination with specific reagents, previously listed in the section on materials that are necessary but not supplied.

Test Principle

The Giemsa stain is a complex mix composed of methylene blue and azure that can bind to the nuclei and micro-organisms, turning them dark blue when viewed through the microscope. This solution increases the intensity of the nuclear staining and can selectively highlight the cell elements.

It can be used in combination with the May Grünwald acc. to Pappenheim solution, which consists of methylene blue eosinate and which stains the nuclei blue and the basophilic cytoplasm pink-red, with the granules taking on a stain ranging from red to purple, depending on their nature, that can be acidophilic or basophilic.

Instructions for use

It is good practice to use a humidity chamber to avoid excessive section drying.

Diluted working solutions must be prepared before use.

For paraffin-embedded histological specimens start at step 1; for fixed cytological specimens, start at step 2.

May Grünwald - Giemsa staining protocol on blood smears

1. Air-dry the slides
2. Dilute **Phosphate Buffer pH 6.8** with distilled water in a ratio of 1:10.
Example: to prepare 500 ml phosphate buffer solution pH 6.8, dilute 50 ml Phosphate Buffer pH 6.8 in 450 ml distilled water.
3. Dilute **Giemsa acc. to Pappenheim** with Phosphate Buffer solution pH 6.8 in a 1:10 ratio.
Example: to prepare 500 ml Giemsa working solution, dilute 50 ml Giemsa acc. to Pappenheim in 450 ml phosphate buffer solution pH 6.8 (filter if necessary).
4. Immerse the slides in **May Grünwald acc. to Pappenheim** for 5 minutes
5. Wash with running water for 1 minute
6. Immerse the slides in **Giemsa working solution**¹ for 20-30 minutes
7. Wash with running water for 1 minute
8. Air-dry the slides
9. Clear and mount**

Notes: ¹The Giemsa working solution must be freshly prepared before use and disposed of after each working day.

Example of Giemsa staining protocol to search for *Helicobacter Pylori* on histological sections (Lennert staining method)

- | | |
|--|------------|
| 1. De-wax and hydrate the sections with water* | |
| 2. Giemsa solution (undiluted) | 30 minutes |
| 3. Distilled water | 15 seconds |
| 4. Acidified Water 0.5% | 10 seconds |
| 5. Distilled water | 15 seconds |
| 6. 95% alcohol | 30 seconds |
| 7. 99% alcohol | 1 minute |
| 8. 99% alcohol | 1 minute |
| 9. Xylene | 3 minutes |
| 10. Xylene | 3 minutes |
| 11. Xylene | 3 minutes |
| 12. Mount the slide ** | |

*Before proceeding with the staining protocol, it is necessary to deparaffinize and hydrate the sections . The preliminary hydration of the section must be carried out using specific reagents to remove the paraffin layer from the tissue and prepare the tissue for subsequent staining. Example of protocol:	**On completion of staining, the sections must be dehydrated, cleared, and mounted . Example of protocol:
Oven at 60° for 20 min	
Xylene: 3 min	70% alcohol: 1 min
Xylene: 3 min	95% alcohol: 1 min
Xylene: 3 min	99% alcohol: 2 min
99% alcohol: 3 min	99% alcohol: 2 min
99% alcohol: 3 min	Xylene: 3 min
95% alcohol: 3 min	Xylene: 3 min
70% alcohol: 3 min	Xylene: 3 min
Running water: 1 min	Mounting medium compatible with histological slides
Immersion in distilled water	

It is good practice to follow all the manufacturer’s instructions in order to obtain the desired result, as an incorrect procedure may hinder the quality of the staining.

All operational procedures, prior to staining, i.e. fixation, dehydration, clearing, paraffin impregnation and finally cutting, are supervised by the operator and must be validated according to normal routine operations within the reference laboratory.

Used in combination with automatic instruments, it requires validation of the staining protocol.

It is neither self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel, since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified and trained personnel.

Interpretation of the results

on blood smears:

Basophilic cytoplasm:	Variants of Blue
Acidophilic cytoplasm:	Pink-Red
Nuclei:	Purple

on histological sections:

Nuclei:	Blue
Background:	Azure - Pink
H. Pylori:	Deep blue
Mast cells:	Reddish-Purple

Quality control

It is necessary to include a positive control slide with each staining run in order to verify the validity of the test performed. It is important that the positive control is treated similarly to the test samples, i.e. tested in accordance with the operational techniques of the test samples.

Analytical performance characteristics

Not applicable

Conditions for collection, handling and preparation of the sample

The preparation of the cytological and histological sample follows the activities expected in the normal workflow of a pathology laboratory.

Do not handle until all safety precautions have been read and understood.

Storage

Store the device at a temperature between +15°C and +30°C.

Do not expose the reagent to environmental conditions other than those prescribed, as failure to comply with storage conditions would lead to a change in the chemical and physical characteristics of the device, resulting in underperforming staining.

If properly stored and in unopened packaging, the product is stable as stated in the section *Description of the reagents and any limitations of use*. Do not use after the expiry date.

If the product is not stored as instructed, its yield may vary and it will have to be tested by the operator.

Stability when handling the device

The use of the device is under the absolute control of the skilled and qualified operator.

The environmental conditions in the laboratory must be controlled:

- Temperature: +15°C /+30°C
- Relative humidity ≤80%

Disposal instructions / Safe final disposal of the device

For specific instructions on how to safely dispose of the device after use, please refer to the local regulations of the country in which the device is used.

It is good practice to dispose of the device as special waste in dedicated containers and in accordance with internal disposal procedures.

For information on the content of the reagents, please refer to the safety data sheet. If the product has expired and/or is unused, it must be disposed of in accordance with local waste regulations, based on the hazard classification on the label, after checking for any contamination. In some cases, analytical evaluation may be required so as to correctly classify the waste and determine the hazard characteristics to be attributed to it

Warnings, precautions and restrictions on the use of the device

The device is to be used only for its intended use.

For proper handling, it is recommended to:

- Wear personal protective equipment, such as a lab coat, gloves and goggles. Protect your eyes and face
- Work with the device under a fume hood
- Avoid inhaling the fumes of vapours while handling the device
- In order to obtain the expected result, it is advisable to follow the indications given in the protocol of use
- Keep the work area clean in order to avoid compromising the diagnostic results

Device malfunction or degradation

The device retains its properties throughout its estimated useful life if the storage requirements are guaranteed and complied with.

The device does not generate risks from any debris due to wear and tear, degradation products and processing residues.

In case of malfunction, isolate the device, replace it with a new one and dispose of the non-functioning one. Should this occur, please contact the manufacturer by sending an e-mail to claims@diapath.com

Error	Possible cause	Possible corrective action
Loss of staining specificity	Failure to observe the instructions on correct preservation of the device after opening or unsuitable storage conditions	Consult the storage instructions and temperatures beforehand for suitable maintenance of the device, also provided on the label
Unsuitable handling during use	Non-performing device	Staining protocol provided on the IFU given to personnel

Warnings from the user, with the need to report any serious accident that occurred in relation to the device

In case of a serious incident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Key of the symbols used in the labelling and instructions for use

Symbol	Symbol name	Explanatory notes
	Manufacturer	Indicates the manufacturer of the device
	Expiration date	Indicates the date after which the device must not be used
	Lot number	Indicates the manufacturer's batch necessary for batch identification
	Catalogue code	Indicates the code assigned by the manufacturer allowing its identification
	Distributor	Indicates the company distributing the medical device
	Not sterile	Indicates a device that has not undergone a sterilization process

	Keep away from light	Indicates a device that needs protection from light sources.
	Temperature limit	Indicates the temperature limits to which the device can be safely exposed
	Not reusable	Indicates a device that is to be used once, for single use only
	Refer to instructions	Indicates the need for the user to consult the instructions
	Attention	Indicates the need for the user to be cautious when operating the device or that the situation requires greater awareness to avoid undesirable consequences
	In vitro diagnostic medical device	Indicates a device intended to be used as an in vitro diagnostic medical device
	CE Marking	The CE mark certifies the conformity of the product to the applicable directives of the European Union
	Content sufficient for <n> tests	Indicates the number of tests that can be performed with the device
 XX.XX.XXXX	Manufacturing date	Indicates the date of manufacture of the medical device
	Unique identifier of the device	Indicates a vector containing the fixed and variable information related to the UDI

Identification of the substance/mixture and identification of the hazards

Substance or mixture	Hazard pictograms	Hazard statements	Cautionary advice
GIEMSA ACC. TO PAPPENHEIM	  	<p>H225 Highly flammable liquid and vapour.</p> <p>H301+H311+H331 Toxic if swallowed, in contact with skin or if inhaled.</p> <p>H370 Causes damage to organs.</p>	<p>P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.</p> <p>P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection/...</p> <p>P301+P310 IF SWALLOWED: immediately call a POISON CENTRE or doctor/physician/...</p> <p>P308+P311 If exposed or concerned: call a POISON CENTRE or doctor/physician/...</p>

			P370+P378 In case of fire: use ... to extinguish. P403+P233 Keep the container well closed in a well-ventilated place. P403+P235 Store in a well ventilated place. Keep cool
--	--	--	--

Bibliography

- “Laboratory methods in histotechnology. Stain methods of the Armed Forces Institute of Pathology”. Prophet, Mills, Arrinton, Sobin. Editor McGRAW Hill. Washington D.C. 1968
- “Tecnica di anatomia patologica”. M.Melis, F. Carpino, U. Di Tondo. Edi Ermes. 1989
- “Theory and practice of Histological Techniques”. Edited by Bancroft JD and Gamble M. Churchill Livingstone. New York. 2002

Revision table

Rev.	Date	Purpose of the revision
1st issue	27/05/2024	Compliance with the new Regulation 746/17

(*) If the request is due to a change, please provide a detailed description of the reason for revising the document and the changes that have been made since the last revision; if available, please enclose recalls, non-conformities, NC CA PA statistics, etc.

You can find this document in the reserved area of the www.diapath.com website. If you want to download the document and keep up to date on the latest available version, you will need to log in after registering.

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