DIAPATH

Constant Innovation in Anatomic Pathology

INSTRUCTION FOR USE

T0212

EN

IT

FR

ES

DE

Phosphate buffer pH 6.8

Buffers for stain dilution

EN

Manufacturer

Diapath S.p.A.

Via Savoldini, 71 - 24057 Martinengo (BG) - Italy

Tel. +39 0363 986411

Fax +39 0363 948000

www.diapath.com

info@diapath.com

EMDN code: W01030799

Code	Format	UDI code
T0212	500 ml	08033895012119

Intended use

Single-use, non-sterile, in vitro diagnostic device, for clinical use, used for the dilution of Giemsa stain according to Pappenheim. Device involved in the treatment process for sample preparation and intended for professional use only, in medical and veterinary laboratories and medical facilities, in manual mode or in combination with automatic stainers.

For further information, warnings and precautions for use, please read the safety data sheet and consult the website.

Description of the reagents and any limitations of use

Sodium dibasic phosphate	CAS No. 7558-79-4	EC No. 231-448-7
Sodium phosphate monobasic	CAS No. 10049-21-5	EC No. 231-449-2
Distilled water	CAS No. 7732-18-5	EC No. 231-791-2

Device	Stability of closed product	Stability after opening	Ready for use	Storage temperature
T0212	2 years	2 years	YES	15°C – 30°C

Information relating to interfering substances that may affect the performance of the device

The device may be contaminated by reagents found in the work environment.

List of materials provided and specific materials required but not provided

The materials required for the analysis, but not provided, are:

- Giemsa stain acc. to Pappenheim, used in combination with May Grunwald acc. to Pappenheim for May-Grunwald Giemsa staining.
- Reagents for histological and cytological staining such as alcohols, xylenes or substitutes, distilled water
- Mounting medium
- Coverglass or substitutes
- Pipettes

Instrumentation required:

- Automatic, semi-automatic and manual staining trays

Information on the combined use of the device

The device has been designed and manufactured to be compatible with tissues and with other devices it is used with (reagents, automatic stainer, slides), without hindering performance.

Test principle

Phosphate buffer solution pH 6.8 is used as a diluent for Giemsa to maintain a stable pH level when a small amount of acid or alkali is added to the solution.

Instructions for use

The device is concentrated, therefore it must be used after dilution.

1. Dilute **Phosphate Buffer pH 6.8** with **distilled water** in a ratio of **1:10**.
Example: to prepare 500 ml phosphate buffer solution pH 6.8, dilute 50 ml Phosphate Buffer pH 6.8 in 450 ml distilled water.
2. Dilute **Giemsa acc. to Pappenheim** with **Phosphate Buffer** solution **pH 6.8** in a **1:10** ratio.
Example: to prepare 500 ml Giemsa working solution, dilute 50 ml Giemsa acc. to Pappenheim in 450 ml phosphate buffer solution pH 6.8 (filter if necessary).

Recommended staining protocol May-Grunwald Giemsa for air-dried cytological smears

1. Air-dry the slides
2. Immerse the slides in **May Grunwald acc. to Pappenheim** for 5 minutes
3. Wash with running water for 1 minute
4. Immerse the slides in **Giemsa working solution**¹ (prepared according to steps 1,2) for 20-30 minutes
5. Wash with running water for 1 minute
6. Air-dry the slides
7. Clear and mount

Notes: ¹ The Giemsa working solution must be freshly prepared before use and disposed of after each working day.

It is neither a self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel, since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified and trained personnel.

Interpretation of the results

Nuclei:	Purple
Basophilic cytoplasm:	Different shades of Blue
Acidophilic cytoplasm:	Light Red - Pink

Analytical performance characteristics

Not applicable

Conditions for sampling, handling and preparation of the sample

The device is designed for processing cytological samples, previously collected and handled by experienced and qualified personnel.

Storage

Store the device at a temperature between +15°C and +30°C.

Keep the container well closed after each use, away from sources of light and heat.

The product is stable when properly stored in unopened packaging. Do not use after the expiry date.

If the product is not stored as instructed, its yield may vary and it will have to be tested by the operator.

Stability when handling the device

If handled incorrectly, the device may come into superficial contact with the operator.

It is necessary to wear the PPE required in the workplace.

Do not use the device after the expiry date.

Disposal instructions / Safe final disposal of the device

Specific instructions on how to safely dispose of the device are to be found in the local regulations of the country where the device is used.

If the product has expired and/or is unused, it must be disposed of in accordance with local waste regulations, based on the hazard classification indicated on the label and after checking for possible contamination. In some cases, analytical evaluation may be required so as to correctly classify the waste and determine the hazard characteristics attributed to it.

Warnings, precautions and restrictions on the use of the device

- Device used exclusively for the purpose specified in *Intended Use*
- Use the device only after wearing gloves and other PPE
- Comply with all safety requirements set out for the work environment where the device is used
- Maintain a clean work environment

Device malfunction or degradation

The device retains its properties throughout its estimated useful life if the storage requirements are guaranteed and complied with.

The device does not generate risks from any debris due to wear and tear, degradation products and processing residues.

In case of malfunction, isolate the device, replace it with a new one and dispose of the non-functioning one. Should this occur, please contact the manufacturer by sending an e-mail to claims@diapath.com

Error	Possible cause	Possible corrective action
Failure to maintain a clean work environment	Contamination of the device and noticeable traces of dirt particles	It is recommended to maintain a clean work environment so as to prevent possible contamination of the sample
	Artefact generation	
Delayed diagnosis or non-diagnosis	Device stored under unsuitable temperature conditions	Symbols and appropriate warnings on the device label
Deliberate use past the expiry date or intended use restrictions	Non-performing device	Symbols and appropriate warnings on the device label

Warnings from the user, with the need to report any serious accident that occurred in relation to the device

In case of a serious incident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Key of the symbols used in the labelling and instructions for use

Symbol	Symbol name	Explanatory notes
	Manufacturer	Indicates the manufacturer of the device
	Expiration date	Indicates the date after which the device must not be used
	Lot number	Indicates the manufacturer's batch necessary for batch identification
	Catalogue code	Indicates the code assigned by the manufacturer allowing its identification
	Distributor	Indicates the company distributing the medical device
	Not sterile	Indicates a device that has not undergone a sterilization process
	Keep away from light	Indicates a device that needs protection from light sources.
	Temperature limit	Indicates the temperature limits to which the device can be safely exposed
	Not reusable	Indicates a device that is to be used once, for single use only
	Refer to instructions	Indicates the need for the user to consult the instructions
	Attention	Indicates the need for the user to be cautious when operating the device or that the situation requires greater awareness to avoid undesirable consequences
	In vitro diagnostic medical device	Indicates a device intended to be used as an in vitro diagnostic medical device
	CE Marking	The CE mark certifies the conformity of the product to the applicable directives of the European Union
	Content sufficient for <n> tests	Indicates the number of tests that can be performed with the device
	Manufacturing date	Indicates the date of manufacture of the medical device

	Country of origin	Identification of the country where the device was manufactured
	Unique identifier of the device	Indicates a vector containing the fixed and variable information related to the UDI

Identification of the substance/ mixture and identification of the hazards

The product is not considered hazardous according to EC Regulation 1272/2008 (CLP).

Bibliography

- Theory and Practice of Histological Techniques, John d. Bancroft, Marilyn Gamble, 2008, Churchill Livingstone Elsevier, Sixth Edition.
- Routine Cytological Staining Techniques: Theoretical Background and Practice, Mathilde E. Boon, Johanna S. Drijver, 1986, Elsevier Science Publishing Company
- Laboratory methods in histotechnology. Stain methods of the Armed Forces Institute of Pathology". Prophet, Mills, Arrinton, Sobin. Editor McGRAW Hill. Washington D.C. 1968

Revision table

Rev.	Date	Purpose of the revision
1st issue	20/07/2023	First issue

(*) If the request is due to a change, please provide a detailed description of the reason for revising the document and the changes that have been made since the last revision; if available, please enclose recalls, non-conformities, NC AC AP statistics, etc.

You can find this document in the reserved area of the www.diapath.com website. If you want to download the document and keep up to date on the latest available version, you will need to log in after registering.

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Bibliografia

- Theory and Practice of Histological Techniques, John d. Bancroft, Marilyn Gamble, 2008, Churchill Livingstone Elsevier, Sixth Edition.
- Routine Cytological Staining Techniques: Theoretical Background and Practice, Mathilde E. Boon, Johanna S. Drijver, 1986, Elsevier Science Publishing Company
- Laboratory methods in histotechnology. Stain methods of the Arm Forces Institute of Pathology". Prophet, Mills, Arrinton, Sobin. Editor McGRAW Hil. Washington D.C. 1968

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