CALCIDOR

10. P.D.

CALCIDOR

CALCIDOR is a radiopaque, self-curing, hard setting calcium hydroxide material recommended for pulp capping and for cavity lining. CALCIDOR forms a protective base under permanent dental filling materials, including amalgam and composites. CALCIDOR will not inhibit the polymerisation of acrylic or composite restorations. CALCIDOR is resistant to phosphoric acid used in acid etching technique. CALCIDOR stimulates the formation of secondary dentin.

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
CALCIDOR

Kalcio hidroksidas

Radiopaque Calcium-Hydroxide composition

Instruction for use for the dentist

Product description:

Calcidor is a rigid, self-setting Calciumhydroxide used for pulp capping and as a protective cavity liner under dental filling materials, cements and composites. Calcidor will not inhibit the polymerization of acrylic or composite restorations. Calcidor serves as a protective barrier between dentin/pulp and filling materials. Calcidor also isolates and seals damaged dentinal tubules. Calcidor is resistant to phosphoric acid used in etching technique. Calcidor stimulates the formation of secondary dentin. Calcidor has very good flow characteristics, is easy to mix and to apply. Calcidor has high compressive strength and low solubility and is radiopaque.

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© 2010 DoriDent Dr. Hirschberg GmbH