DiaSorin Italia S.p.A.UK Branch,
Central Rd, Dartford DA1 5LR, UK
www.diasorin.com
Tel. +44 1322 317949

en

REF 7G79-09 / 11

GE41/42

Revised August, 2023

Murex HIV Ag/Ab Combination

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Enzyme immunoassay for improved detection of seroconversion to human immunodeficiency virus types 1 (HIV-1, HIV-1 group O) and detection of anti-HIV-2 antibodies.

The assay is intended to screen individual human donors for the presence of HIV p24 antigen and antibodies to HIV-1, including group O, and HIV-2 or as an aid to the diagnosis of HIV infection.

Customer Service

For additional product information, please contact your local customer service organization.

These instructions for use must be read carefully prior to use. The instructions for use must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions for use.

IVD

Key to symbols used			
	List Number		<i>In Vitro</i> Diagnostic Medical Device
	Lot Number		Store at 2-8°C
	Expiration Date		CAUTION: Consult accompanying documents
	Manufacturer		Consult instructions for use
			Keep away from sunlight

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

INTENDED PURPOSE

The Murex HIV Ag/Ab Combination is an enzyme immunoassay for the qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 (HIV-1, HIV-1 group O) and type 2 (HIV-2) in human serum or plasma samples. The assay is intended to screen individual human donors for the presence of HIV p24 antigen and antibodies to HIV-1, including group O, and HIV-2. It is also intended as an aid to the diagnosis of HIV infection in individuals with or without the symptoms of HIV. This assay is intended for *in-vitro* diagnostic use, on either a semi-automated or automated platform.

SUMMARY AND EXPLANATION OF THE TEST

HIV stands for human immunodeficiency virus and it can lead to acquired immunodeficiency syndrome or AIDS if not treated. HIV is a retrovirus, with two subtypes: HIV-1 and HIV-2¹. The HIV-1 subtype is the most common and responsible for AIDS throughout most of the world, while HIV-2 is found primarily in Western Africa.

HIV has a worldwide and pandemic distribution² and, since its identification, 39 million people have died from HIV infection and more than 35 million people are currently living with HIV infection³. HIV transmission occurs through sexual contact between HIV-infected individuals, exposure to contaminated blood or blood products, and prenatal infection of a fetus or perinatal infection of a newborn from an infected mother⁴.

HIV has a general tropism for the immune system, specifically the CD4 cells (T cells), leading opportunistic infections or cancers development during the last stage of HIV infection. The disease is characterized by three different stages: acute HIV infection, within the first 2-4 weeks, clinical latency, which can last decades with low viral reproduction and the AIDS phase^{5, 6}.

During the course of infection, several immunological and viral markers appear and can be monitored and used to identify HIV infection: HIV RNA, HIV p24 antigen (coded by the gag gene), and antibodies to HIV antigens. The time interval before HIV-specific antibodies appear is known as the serological "window period". The detection of specific antibody to HIV defines the end of the window period and labels the individual as "seropositive"⁷.

The introduction of fourth generation serological assays, able to detect both HIV p24 antigen and HIV-1/2 antibodies, allows the earlier identification of infected individuals and can shorten the diagnostic window period⁸.

Transfusion-transmitted infection risk continue to be a big challenge and the screening of blood and blood products for the presence of HIV plays an essential role in the prevention of HIV infections. The implementation of nucleic acid testing (NAT) together with the existing antigen/antibody-based assays for donors screening has further reduced the residual risk of recipient's infection⁹⁻¹¹.

PRINCIPLE OF THE PROCEDURE

Murex HIV Ag/Ab Combination is based on microwells coated with synthetic peptide representing immunodominant regions of HIV-1 (O) and HIV-2, recombinant protein derived from the envelope regions of HIV-1 and HIV-2 and HIV pol protein, together with monoclonal antibodies raised against p24 of HIV-1. The Conjugate is a mixture of the same antigen epitopes, and different monoclonal antibodies, also raised against p24, all labelled with horseradish peroxidase.


Test specimens and control sera are incubated in the wells and reactive HIV core and/or antibodies to HIV in the sample or control sera bind to the antibodies and/or antigens on the microwell; sample and any excess antibodies are then washed away. In a subsequent step, Conjugate is added which in turn binds to any reactive HIV core and/or specific antibody already bound to the reagents on the well. Samples not containing reactive core antigen or specific antibody will not cause the Conjugate to bind to the well.

Unbound Conjugate is washed away and a solution containing 3, 3', 5, 5'-tetramethylbenzidine (TMB) and hydrogen peroxide is added to the wells. Wells with bound Conjugate develop a blue green colour which is converted to an orange colour which may be read at 450 nm after the reaction has been stopped with sulphuric acid.

REAGENTS

DESCRIPTION, PREPARATION FOR USE AND RECOMMENDED STORAGE CONDITIONS

See also **Warnings and Precautions**.

All components must be stored at 2 to 8°C, unless otherwise stated, under which condition they will retain activity until the expiry date of the kit.

COATED	WELLS
--------	-------

1. Coated Wells

One plate (7G79-09) or five plates (7G79-11) of 96 microwells coated with HIV antigens (0.02 to 0.25 µg/mL) and monoclonal antibodies (0.6 µg/mL).

Allow the wells to reach room temperature (18 to 30°C) before removal from the bag.

Place unused wells in the sealable storage bag provided and return to 2 to 8°C.

SAMPLE	DIL
--------	-----

2. Sample Diluent

One bottle containing 8 mL (7G79-09) or 18 mL (7G79-11) of a green/brown buffer solution, bovine and murine protein, detergent and saponin. Contains 0.05% ProClin[®] 300 preservative.

CONJUGATE

3. Conjugate

One bottle (7G79-09) or three bottles (7G79-11) containing 1.1 mL of HIV antigens (0.2 to 1.8 µg/mL) and monoclonal antibodies (0.6 to 2.5 µg/mL) conjugated to horseradish peroxidase and freeze dried. When reconstituted each bottle is sufficient for up to two plates.

CONJUGATE	DIL
-----------	-----

4. Conjugate Diluent

One bottle (7G79-09) or three bottles (7G79-11) containing 22 mL of a yellow solution consisting of buffer, bovine protein, saponin and detergent, sufficient to reconstitute one bottle of Conjugate. Contains 0.1% ProClin[®] 300 preservative.

Reconstitution of Conjugate

Tap the bottle of Conjugate gently on the bench to remove any material adhering to the rubber stopper. Pour the whole contents of the bottle of conjugate diluent into the bottle of conjugate, recap the latter and mix by gentle inversion. Allow to rehydrate for at least 30 minutes with occasional swirling. The reconstituted conjugate will be red in colour. Reconstituted conjugates may be returned to and pooled in the plastic conjugate diluent bottles if required.

After reconstitution the Conjugate may be stored at 2 to 8°C for up to four weeks.

5. Anti-HIV-1 Positive Control

One bottle containing 1.7 mL of inactivated human serum in a buffer containing bovine protein. Contains 0.05% Bronidox[®] preservative.

6. Anti-HIV-2 Positive Control

One bottle containing 1.7 mL of inactivated human serum in a buffer containing bovine protein. Contains 0.05% Bronidox[®] preservative.

7. HIV-1 p24 Positive Control

One bottle containing 1.7 mL of p24 (recombinant antigen) in a buffer containing bovine protein. Contains 0.05% Bronidox[®] preservative.

8. Negative Control

Two bottles containing 2.5 mL of normal human serum diluted in a bovine protein buffer. Contains 0.05% Bronidox[®] preservative.

9. Substrate Diluent

One bottle containing 35 mL of a colourless solution of tri-sodium citrate and hydrogen peroxide.

10. Substrate Concentrate

One bottle containing 35 mL of 3, 3', 5, 5'-tetramethylbenzidine (TMB) and stabilisers in an orange solution.

Substrate Solution

To prepare the Substrate Solution add a volume of colourless Substrate Diluent to an equal volume of orange Substrate Concentrate in either a clean glass or plastic vessel. **It is important that this order of addition is followed and that any pipettes and glassware used to prepare Substrate Solution are clean.** Alternatively, the Substrate Solution may be made by pouring the entire contents of the bottle of Substrate Diluent into the bottle of Substrate Concentrate. One bottle of Substrate Solution provides sufficient reagent for at least five plates - see **Table 1**.


Table 1
Volume of Substrate Concentrate and Substrate Diluent Required

Number of Wells	Number of Plates
8 16 24 32 40 48 56 64 72 80 96	1 2 3 4
Substrate Concentrate (mL)	
0.5 1.0 2.0 2.5 2.5 3.0 3.5 4.0 4.5 4.5 6.0	6 12 18 22
Substrate Diluent (mL)	
0.5 1.0 2.0 2.5 2.5 3.0 3.5 4.0 4.5 4.5 6.0	6 12 18 22

Additional reagent may be required for use with automated systems. Keep away from all natural and artificial light. The Substrate Solution should be pale yellow; if it is green before being used it should be discarded and fresh Substrate Solution prepared.

The prepared Substrate Solution from this kit may be used interchangeably with that from all other Murex kits which use orange coloured Substrate Concentrate. Ensure that the Substrate Solution is prepared from the Substrate Diluent and Substrate Concentrate provided together.

The prepared Substrate Solution is stable refrigerated (2 to 8°C) or at 15 to 25°C for up to two days but it must be discarded if crystals have formed.

11. Wash Fluid

One (7G79-09) or two (7G79-11) bottles containing 125 mL of 20 times working strength Glycine/ Borate Wash Fluid. Contains 0.2% Bronidox[®] preservative.

Add one volume of Wash Fluid Concentrate to 19 volumes of distilled or deionised water to give the required volume or dilute the entire contents of one bottle of Wash Fluid to a final volume of 2500 mL. Crystals may be observed in the Wash Fluid Concentrate but these crystals will dissolve when the Wash Fluid is diluted to working strength. When diluted the Wash Fluid contains 0.01% Bronidox[®] preservative.

The Wash Fluid from this kit may be used interchangeably with the Glycine/Borate Wash Fluid from any other Murex kit.

Store the working strength Wash Fluid at 18 to 30°C in a closed vessel under which conditions it will retain activity for one month.

NOTE: The Wash Fluid may develop a yellow colour on storage. This will have no effect on the performance of the assay providing the Wash Fluid is fully aspirated from the wells.

NOTE: Although the Substrate Solution and Wash Fluid are interchangeable, they must not be used beyond the expiry date printed on the component labels.

WARNINGS AND PRECAUTIONS

The reagents are for *in vitro* diagnostic use only.

For professional use only.

Please refer to the manufacturer's safety data sheet and the product labelling for information on potentially hazardous components.

Low levels of fibrin precipitate may be observed in the Kit Controls and product performance is not affected by this. This is a product of certain serum batches used to manufacture the controls.

HEALTH AND SAFETY INFORMATION

CAUTION: This kit contains components of human origin.

The human sera used for manufacture have been screened and found reactive or non-reactive for analytes as shown in **Table 2** below.

Table 2

Component	Reactive for	Non-reactive for
Negative Control	N/A	HBsAg, antibodies to HCV, HIV-1 and HIV-2
Positive Control 1	antibodies to HIV-1	HBsAg
Positive Control 2	antibodies to HIV-2	HBsAg

Additionally human sera used for positive controls are also tested for antibodies to HCV and may be reactive. All reactive serum used has been inactivated prior to use in reagent preparation. However, all material of human origin should be considered as potentially infectious and it is recommended that this kit and test specimens be handled using established good laboratory practice.

Pursuant to EC Regulation 1272/2008 (CLP) hazardous reagents are classified and labeled as follows:

Reagents:	CONJUGATE DIL	SAMPLE DIL	CONJUGATE*
Classification:	Skin sens. 1 H317		
Signal Word:	Warning		
Symbols / Pictograms:			
Hazard Statements:	H317 May cause an allergic skin reaction.		
Precautionary Statements:	P280 Wear protective gloves/protective clothing/eye protection/face protection. P363 Wash contaminated clothing before reuse. P333+P313 If skin irritation or rash occurs: Get medical advice / attention.		
Contains:	Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6](3:1).		
* The reconstituted Conjugate contains 0.1% ProClin® 300 which is classified hazardous per EC Regulation 1272/2008.			

Reagents:	SUBSTRATE	CONC
Classification:	Eye Irrit. 2 H319	
Signal Word:	Warning	
Symbols / Pictograms:		
Hazard Statements:	H319 Causes serious eye irritation	
Precautionary Statements:	P264 Wash hands thoroughly after handling. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.	

Pursuant to EC Regulation 1272/2008 (CLP), WASH FLUID is labeled as EUH210, safety data sheets available on request.

For additional information see Safety Data Sheets available on www.diasorin.com

- Potentially contaminated materials should be disposed of safely according to local requirements.
- Spillage of potentially infectious material should be removed immediately with absorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued¹². Sodium hypochlorite should not be used on acid containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- Neutralised acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.

- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- If any of the reagents come into contact with the skin or eyes wash the area extensively with water.
- Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If either come into contact with the skin or eyes, wash thoroughly with water.

ANALYTICAL PRECAUTIONS

- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the **Test Procedure** or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
- Allow all reagents and samples to come to 18 to 30°C before use. Immediately after use return reagents to the recommended storage temperature.
- Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionised water.
- Avoid the use of self-defrosting freezers for the storage of reagents and samples.
- Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
- Do not allow wells to become dry during the assay procedure.
- Do not cross-contaminate reagents. Dedicate a pipette for use with the Substrate Solution of Murex assays. A pipette should also be dedicated for use with the Conjugate.
- The Sample Diluent in this assay has the potential to cause false positive results in anti hepatitis B surface antigen (anti-HBs) assays if reagent cross contamination occurs.
If running Murex HIV Ag/Ab Combination in conjunction with an anti-HBs assay on a fixed tip instrument ensure that the possibility of cross contamination is excluded during the validation process.
- Do not touch or splash the rim of the well with Conjugate. Do not blow out from micropipettes; reverse pipetting is recommended whenever possible.
- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate.
- Do not contaminate microwells with dust from disposable gloves.
- When using fully automated processors
 - It is not necessary to use plate lids and tap dry the wells.
 - Do not allow system fluids to contaminate samples or reagents.
 - The possibility of cross contamination between assays needs to be excluded when validating assays on fully automated processors.
- Ensure the assay is run within the temperature limits defined in the assay protocol.
- Do not use CO₂ incubators.
- Do not store the Stop Solution in a shallow dish or return it to a stock bottle after use.
- The possibility of cross contamination between assays needs to be excluded when validating assay protocols on instrumentation.

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SPECIMEN COLLECTION, TRANSPORT AND STORAGE

SPECIMEN COLLECTION

Serum, EDTA plasma or citrate plasma samples may be used. Ensure that the serum samples are fully clotted. Remove any visible particulate matter from the sample by centrifugation. If samples are prepared using liquid anti-coagulants e.g. citrate plasma, the dilution effect should be considered.

SPECIMEN TRANSPORT AND STORAGE

Store samples at 2 to 8°C. Samples not required for assay within 72 hours should be removed from the clot or cell pellet and stored frozen (-15°C or colder). Avoid multiple freeze-thaw cycles. Samples may be stored frozen for up to 3 months. After thawing ensure samples are thoroughly mixed before testing.

PROCEDURE

MATERIALS REQUIRED BUT NOT PROVIDED

- Stop Solution (0.5M to 2M Sulphuric Acid)**, e.g. add between 3.0 mL (for 0.5M) and 11 mL (for 2.0M) of analytical grade concentrated sulphuric acid (18M) to about 80 mL of distilled or deionised water and then make up to 100 mL with more water. Alternatively, the following reagent can be used: 1N Sulphuric Acid (Code N0164 - 15 vial pack and N0165 - 1 vial pack).
- Freshly distilled or high quality deionised water** is required for dilution of Wash Fluid, for preparation of the Stop Solution and for use in conjunction with automated washers.
- Micropipettes and Multichannel micropipettes** of appropriate volume.
- Incubator** capable of maintaining the temperature limits defined in the assay protocol.
- Moulded Heating Block** (Code 5F09-02). For use in laboratory incubators. The moulded heating block should ideally be kept in the incubator used. If this is not possible it must be placed in the incubator at least four hours before beginning the assay.
- Instrumentation**
 - Automated microplate stripwasher,
 - Microplate reader, or
 - Fully automated microplate processor.

All instruments must be alidated before use.
Please contact your representative for details of recommended systems, software protocols for instrumentation and validation procedures.
- Disposable Reagent Troughs**. (Code 5F24-01).
- Sodium hypochlorite** for decontamination. (Refer to **Health and Safety Information**).
- Sodium hydroxide solution** (0.1M). (For Instrument decontamination).

TEST PROCEDURE

Please read **Analytical Precautions** carefully before performing the test.
Addition of the various components of the assay to the wells may be confirmed visually by examining the plate for the following colours:

Sample Diluent is green/brown in colour. On addition of Sample or Control the colour will change to blue/green. The colour change will vary from sample to sample but some change should always be visible. The addition of sample or control may be confirmed using a microplate reader at 570 nm or 620 nm with a reference of 690 nm.

Reconstituted Conjugate is red in colour. The addition of Conjugate may be confirmed using a microplate reader at 490 nm with a reference of 690 nm.

Substrate Solution is initially pale yellow with any reactive wells becoming blue green. On addition of Stop Solution the blue green colour of the reactives will change to orange, whilst the negatives will change to pink. The addition of Substrate Solutions may be confirmed using a microplate reader at 450 nm (no reference).

SEMI AUTOMATED PROCESSING

Step 1	Reconstitute and mix the Conjugate , prepare the Substrate Solution and Wash Fluid .	
Step 2	Use only the number of wells required for the test. Avoid touching the tops or bottoms of the wells.	
Step 3	Add 25 µL of Sample Diluent to each well.	25 µL
Step 4	Add 100 µL of Samples or 100 µL Controls to the wells. For each plate use the first column of wells for the assay Controls. Add the Controls to the designated wells after dispensing the samples. Pipette 100 µL of the Negative Control into each of three wells A1 to C1 and 100 µL of the p24, anti-HIV-1 and HIV-2 Positive Controls into wells D1 to F1 respectively. Use of a white background will aid visualization of sample addition.	100 µL
Step 5	Cover the wells with the lid and incubate for 60 minutes at 37°C ±1°C	60 mins
Step 6	At the end of the incubation time wash the plate as described under Wash Procedures .	
Step 7	Immediately after washing the plate, add 100 µL of Conjugate to each well	100 µL
Step 8	Cover the wells with the lid and incubate for 30 minutes at 37°C ±1°C.	30 mins
Step 9	At the end of the incubation time wash the plate as described under Wash Procedures .	
Step 10	Immediately after washing the plate, add 100 µL of Substrate Solution to each well.	100 µL
Step 11	Cover the wells with a lid and incubate for 30 minutes at 37°C ±1°C. Keep away from direct sunlight. A blue green colour should develop in wells containing reactive samples.	30 mins
Step 12	Add 50 µL of Stop Solution (0.5M to 2M sulphuric acid) to each well.	50 µL
Step 13	Within 15 minutes read the absorbance at 450 nm using 620 nm to 690 nm as the reference wavelength if available. Blank the instrument on air (no plate in the carriage).	A₄₅₀

WASH PROCEDURES

Protocols for recommended washers and procedures for verifying washers and analysers can be obtained from your representative. The following protocol is recommended:

a. Protocol for automated stripwasher

Perform 5 wash cycles using working strength Wash Fluid. Ensure, where possible, that:

- i) Flow-through washing with a volume of 500 µL/well is used with instrumentation supplied by DiaSorin. When using other instrumentation for which this is not possible, ensure that the well is completely filled.
- ii) The dispense height is set to completely fill the well, with a slight positive meniscus, without causing an overflow.
- iii) The time taken to complete one aspirate/wash/soak cycle is approximately 30 seconds.
- iv) Ensure that no liquid is left in the well (by use of a double aspirate step in the final cycle where possible).
- v) After washing is completed, invert the plate and tap out any residual Wash Fluid onto absorbant paper.

NOTE: Do not allow the wells to become dry during the assay procedure.

Washers must be rinsed with distilled or deionised water at the end of the test to avoid blockage and corrosion.

FULLY AUTOMATED PROCESSORS

Contact your representative for details of currently available validated protocols. For instrumentation without established validated protocols, the following guidelines are recommended:

1. Do not programme times shorter than specified in the procedure.
2. For each incubation at 37°C, programmed times may be increased by up to 5 minutes.
3. Ensure all **Analytical Precautions** are followed. Protocols written following these guidelines must be fully validated prior to use according to local procedures

RESULTS

CALCULATION OF RESULTS

Each plate must be considered separately when calculating and interpreting results of the assay.

Approved software may be used for calculation and interpretation of results.

Negative Control

Calculate the mean absorbance of the Negative Controls.

Example:

Well 1	=	0.084	
Well 2	=	0.086	
Well 3	=	0.070	
Total	=	0.240	
Mean Negative Control	=	0.240/3	= 0.080

If one of the Negative Control Wells has an absorbance more than 0.15 O.D. above the mean of all three, discard that value and calculate the new Negative Control mean from two remaining replicates.

Cut-off value

Calculate the Cut-off value by adding 0.150 to the mean of the Negative Control replicates (see above).

Mean Negative Control	=	0.080	
Cut-Off value	=	0.080+0.150	=0.230

QUALITY CONTROL

Results of an assay are valid if the following criteria for the controls are met:

Negative Control

The mean absorbance is less than 0.15.

Positive Controls

The absorbance of each of the Positive Controls is more than 0.8 above the mean absorbance of the Negative Control.

Assays which do not meet these criteria should be repeated.

In the unlikely event of the results repeatedly failing to meet either the Quality Control criteria or the expected performance of the test, please contact your representative.

INTERPRETATION OF RESULTS

Non-reactive Results

Samples giving an absorbance less than the Cut-off value are considered negative in the assay.

Reactive Results

Samples giving an absorbance equal to or greater than the Cut-off value are considered initially reactive in the assay (see **Limitations of the Procedure**).

Unless local procedures state otherwise, such samples **must** be retested in duplicate using the original source. Samples that are reactive in at least one of the duplicate retests are considered repeatedly reactive in Murex HIV Ag/Ab Combination and are presumed to contain reactive HIV core antigen and/or antibodies to HIV-1 or HIV-2. Such samples **must** be further investigated and the results of this assay considered with any other clinical and/or assay information. Samples that are non-reactive in both wells on retest are considered non-reactive for HIV core antigen and HIV antibodies.

No sample addition

Absorbance values significantly higher than the Negative Control may be obtained in wells where the sample has been omitted but all the reagents have been added.

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance of Murex HIV Ag/Ab Combination has been determined by testing samples from random blood donors, patients with AIDS diagnosed according to CDC criteria, patients with AIDS Related Complex (ARC), other patients with known antibody to HIV-1 (including group O), patients with confirmed HIV-2 infection and patients at risk of HIV infection or in other clinical categories. In addition, its performance on commercially available seroconversion panels has been evaluated.

Diagnostic Sensitivity

A total of 497 specimens from patients with confirmed HIV-1 infection were tested and found to be reactive with Murex HIV Ag/Ab Combination. The specimens were taken from patients at various stages of HIV infection and included 24 specimens from patients with HIV-1 subtype O infection and a further 139 specimens from patients infected with HIV-1 subtypes other than subtype B.

In addition a total of 100 specimens from patients with confirmed HIV-2 infection were also tested with Murex HIV Ag/Ab Combination and found to be reactive.

The diagnostic sensitivity of Murex HIV Ag/Ab Combination on this population of specimens is therefore estimated to be 100% (597/597) with a lower 95% confidence limit of 99.38% (593/597) by the binomial distribution.

A total of 32 commercial HIV-1 seroconversion panels were tested with Murex HIV Ag/Ab Combination. Using the presence of both core (p24) and an envelope (gp120/160) band on Western blot as the reference criteria, Murex HIV Ag/Ab Combination detected antibody to HIV earlier or in the same sample as Western blot in all of the panels.

A total of 68 cell culture supernatants including HIV-1 and HIV-2 and 80 HIV-1 antigen positive specimens were tested with the Murex HIV Ag/Ab Combination and correctly identified.

Analytical sensitivity

Sensitivity on AFSSAPS HIV Ag Standard

Sensitivity of Murex HIV Ag/Ab Combination on the AFSSAPS HIV Ag standard was determined at three testing centres.

Table 3
AFSSAPS HIV Ag standard

Centre	Sensitivity HIV Ag pg/mL
1	31
2	28
3	25
Mean	28

The data shown in Table 3 was obtained during this testing but may not be exactly reproducible on other testing occasions.

Sensitivity with the First International Reference Reagent HIV-1 p24 Antigen NIBSC code 90/636

The sensitivity of the First international reference reagent HIV-1 p24 Antigen NIBSC code 90/636 with Murex HIV Ag/Ab Combination was determined to be 0.977 IU/mL.

Diagnostic Specificity

The Murex HIV Ag/Ab Combination assay demonstrated a specificity of ≥99.5% in a study where specimens from a European blood donor population were tested. A total of 9,290 routine donor plasma specimens were screened with Murex HIV Ag/Ab Combination at three European blood transfusion centres. The results are summarised in Table 3. In the study, 99.77% (9269/9290) of specimens were non-reactive and 0.23% (21/9290) were repeatedly reactive. One of the repeatedly reactive specimens was weakly positive with the Murex HIV Antigen mAb (8E77). None of the remaining 20 specimens were confirmed as positive for the presence of HIV-1 antigen or antibody to HIV-1 or HIV-2.

The specificity of Murex HIV Ag/Ab Combination on presumed negative European blood donors is estimated to be 99.78% (9269/9289) with 95% confidence limits of 99.67% (9258/9289) to 99.87% (9277/9289) by the binomial distribution.*

A total of 272 specimens from patients with conditions unrelated to HIV infection were also tested with Murex HIV Ag/Ab Combination. These included specimens from pregnant women and patients suffering with autoimmune disease and other acute viral infections. A total of five specimens were reactive with Murex HIV Ag/Ab, four were reactive with two other commercially available screening assays. In Western blot studies four produced indeterminate results and one was negative.

In addition, 30 lipaemic, icteric and haemolysed specimens were also tested and found to be non-reactive. 10 HAMA positive specimens were tested and found to be non-reactive.

The overall diagnostic specificity of Murex HIV Ag/Ab Combination on confirmed negative specimens during this performance evaluation is estimated to be 99.78% (9561/9582) with 95% confidence limits of 99.67% (9550/9582) to 99.86% (9568/9582) by the binomial distribution.*

*Representative performance data are shown. Results obtained at individual laboratories and with different populations may vary.

Table 4
Reactivity of Murex HIV Ag/Ab Combination with presumed negative specimens from routine European blood donors

Centre	Number of presumed negative specimens tested	Number of repeatedly reactive specimens
A	3095	6 ^a (0.19%)
B	2803	9 (0.32%)
C	3392	6 (0.18%)
TOTAL	9290	21 (0.23%)

^a Includes one specimen which was weakly positive in Murex HIV Antigen mAb (8E77).

Analytical Specificity

Controlled studies of the potentially interfering substances showed no interference at the concentrations for each substance summarised in Table 5 for the Murex HIV Ag/Ab Combination.

Table 5
Performance with potentially interfering substances

Substance (Endogenous)	Concentration tested
Total protein – low	No Interference at 50 g/L
Triglycerides	No Interference at 3000 mg/dL
Vitamin H (Biotin)	No Interference at 3500 ng/mL
Unconjugated bilirubin	No Interference at 40 mg/dL
Conjugated bilirubin	No Interference at 40 mg/dL
Cholesterol	No Interference at 400 mg/dL

Substance (Exogenous)	Result
Adefovir dipivoxil	No Interference at 10 mg/L
Entecavir	No Interference at 0.5 mg/L
Lamivudine	No Interference at 300 mg/L
Tenofovir	No Interference at 0.0978 mg/dL
Acetaminophen	No Interference at 15.6 mg/dL
Emtricitabine	No Interference at 0.9 mg/dL
Interferon alpha 2a	No Interference at 6000 IU/mL
Interferon alpha 2b	No Interference at 6000 IU/mL
Interferon alpha 1b (IFN-alpha-1)	No Interference at 6000 IU/mL
Ibuprofen	No Interference at 21.9 mg/dL
Telbivudine	No Interference at 600 mg/L
Vitamin C (L-Ascorbic acid)	No Interference at 300 mg/L
Sofosbuvir	No Interference at 0.190 mg/dL

Assay Reproducibility

The reproducibility of Murex HIV Ag/Ab Combination was assessed by testing two of the assay controls and four quality assurance panel members as ten replicates on four separate occasions. The results from the testing are summarised in Table 6.

Table 6
Murex HIV Ag/Ab Combination - Assay Reproducibility

Specimen	Number of Assays	Number of Replicates	Mean Absorbance/Cut-off ratio	Intra-assay %CV	Inter-assay %CV
Negative Control	4	10	0.266	8.7	11.3
HIV-1 Positive Control	4	10	8.287	4.3	4.7
QA01	4	10	3.672	4.6	7.3
QA02	4	10	4.696	5.6	12.9
QA03	4	10	3.006	3.9	4.2
QA04	4	10	1.663	6.8	9.2

LIMITATIONS OF THE PROCEDURE

1. The **Test Procedure** and **Interpretation of Results** must be followed.
2. This test has only been evaluated for use with individual (unpooled) serum, EDTA plasma or citrate plasma samples.
3. A negative result with an antigen/antibody detection test does not preclude the possibility of infection with HIV.
4. A positive result with Murex HIV Ag/Ab Combination should be confirmed by at least one other test.
5. Non-repeatable reactive results may be obtained with any EIA procedure.
The most common sources of error are:
 - a) Imprecise delivery of Sample, Conjugate or Substrate into the wells.
 - b) Contamination of Substrate with Conjugate.
 - c) Contamination with conjugates from other assays.
 - d) Blocked or partially blocked washer probes.
 - e) Insufficient aspiration leaving a small volume of Wash Fluid in the wells.
 - f) Failure to ensure that the bottom surface of the wells is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before a plate is read.
 - g) Failure to read at the correct wavelength (450 nm) or use of an incorrect reference wavelength (not 620 nm to 690 nm).
6. The use of highly haemolysed samples, incompletely clotted sera, plasma samples containing fibrin or samples with microbial contamination may give rise to erroneous results.
7. Samples which contain higher than normal physiological concentrations of total protein may give rise to erroneous results.
8. This test has not been evaluated for use with samples from cadavers.

Summary of safety and performance is available on EUDAMED.

For EU only: Please be aware that any serious incident that has occurred in relation to this IVD medical device shall be reported to DiaSorin Italia S.p.A. UK Branch and the Competent Authority of the EU Member State in which the user and/or the patient is established.

BIBLIOGRAPHY

1. **Sharp PM** et al. Origins of HIV and the AIDS pandemic. Cold Spring Harb Perspect Med. 2011 Sep;1(1).
2. **Becerra** et al. Recent Insights into the HIV/AIDS Pandemic. Microb Cell. 2016 Sep 05;3(9):451-475.
3. **Poorolajal J.** et al. Survival rate of AIDS disease and mortality in HIV-infected patients: a meta-analysis. Public Health. 2016 Oct;139:3-12.
4. **Shaw G.M.** and Hunter E. HIV Transmission. Cold Spring Harb Perspect Med 2012;2.
5. **Parekh BS.** et al. Diagnosis of Human Immunodeficiency Virus Infection. Clin Microbiol Rev. 2018 Nov 28;32(1).
6. **Pandey D.** et al. HIV Infection: A Review of Their Inhibitors Progression. Biomed Pharmacol J 2017;10(2).
7. **Buttò S.** et al. Laboratory diagnostics for HIV infection. Ann Ist Super Sanita. 2010;46(1):24-33.
8. **Bangalee A.** et al. Evaluation of serological assays for the diagnosis of HIV infection in adults. S Afr Fam Pract (2004). 2021 Oct 25;63(1):e1-e5.
9. **Dwyre DM.** et al. Hepatitis B, hepatitis C and HIV transfusion-transmitted infections in the 21st century. Vox Sang. 2011 Jan;100 (1):92-8.
10. **Salles NA.** Detection of HIV-1 infections in blood donors during the pre-seroconversion window period in São Paulo, Brazil. Rev Soc Bras Med Trop. 2019 May 16;52.
11. **WHO.** (2017) Blood safety and availability. Fact sheet.
12. **Centres for Disease Control.** (1985). Recommendations for preventing transmission of infection with human T-lymphotropic virus type III/lymphadenopathy-associated virus in the workplace. *MMWR*, **34**, No. 45, 681.

Bronidox[®] and ProClin[®] are not trade marks of DiaSorin.

DiaSorin Italia S.p.A. UK Branch
Central Road,
Dartford DA1 5LR
UK

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Enzyme immunoassay for the detection of antibodies to human immunodeficiency virus types 1 (HIV-1, HIV-1 group O) and 2 (HIV-2) in human serum or plasma

The assay is intended to screen individual human donors for the presence of antibodies to HIV-1, including group O, and HIV-2 or as an aid to the diagnosis of HIV infection.

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IVD

Key to symbols used			
	List Number		<i>In Vitro</i> Diagnostic Medical Device
	Lot Number		Store at 2-8°C
	Expiration Date		CAUTION: Consult accompanying documents
	Manufacturer		Consult instructions for use

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INTENDED USE

Enzyme immunoassay for the detection of antibodies to human immunodeficiency virus types 1 (HIV-1, HIV-1 group O) and 2 (HIV-2) in human serum or plasma.

The assay is intended to screen individual human donors for the presence of antibodies to HIV-1, including group O, and HIV-2 or as an aid to the diagnosis of HIV infection.

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Two types of human immunodeficiency virus, HIV-1 and HIV-2, have been described and implicated as causative of the Acquired Immunodeficiency Syndrome (AIDS). Both are retroviruses which are transmitted by exposure to certain infected body fluids, primarily blood and genital secretions, and by transplacental passage. Infection by HIV-1 has been reported worldwide; HIV-2 infection has been reported as occurring mainly in West Africa and some European countries¹. The two types of virus show substantial antigenic cross reactivity in their core proteins, but the envelope glycoproteins are less cross reactive. It is necessary for screening purposes to use epitopes from the envelope proteins of both viruses in addition to the major cross reacting core proteins to ensure detection of antibodies against both types of virus at all stages following infection². Variants of HIV-1, classified together as group O, have been identified in samples from Cameroon and Europe^{3,4}. Group O is highly divergent from the originally known subtypes of HIV-1 (together classified as group M). Specific epitopes from the envelope region of this virus can be used to detect antibody to group O in infected individuals; reliance on cross reactions to the known subtypes of HIV is not satisfactory⁵. The earliest specific antibody response following infection by HIV may be of immunoglobulin M (IgM) followed by a response in immunoglobulin G (IgG)⁶. Maximum sensitivity for detection of anti-HIV seroconversion is achieved by assays which respond to both IgM and IgG.

Murex HIV-1.2.O is expected to detect IgG and IgM and additionally IgA to the envelope glycoproteins and the cross-reacting core proteins of HIV-1 and HIV-2. Consequently potentially infectious samples of serum, EDTA plasma or citrate plasma can be identified.

PRINCIPLE OF THE PROCEDURE

Murex HIV-1.2.O is based on microwells coated with a synthetic peptide representing an immunodominant region of HIV-1 (O), recombinant protein derived from the envelope proteins of HIV-1 and HIV-2 and an HIV core protein. The Conjugate is a mixture of the same epitopes all labelled with horseradish peroxidase.

Test specimens and control sera are incubated in the wells and antibodies to HIV in the sample or control sera bind to the antigens on the microwell; sample and any excess antibodies are then washed away. In a subsequent step, Conjugate is added which in turn binds to any specific antibody already bound to the antigen on the well. Samples not containing specific antibody will not cause the Conjugate to bind to the well. Unbound Conjugate is washed away and a solution containing 3,3',5,5'-tetramethylbenzidine (TMB) and hydrogen peroxide is added to the wells. Wells with bound Conjugate develop a purple colour which is converted to an orange colour when the reaction is stopped with sulphuric acid. After incubation the enzymic reactions are terminated with sulphuric acid and the colour is read spectrophotometrically at 450 nm. The amount of Conjugate, and hence colour, in the wells is directly related to the concentration of antibody to HIV in the sample.

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COATED WELLS

1. Coated Wells

One plate (9E25-01) or five plates (9E25-02) of 96 microwells coated with HIV antigens.

If less than the whole plate is being used allow the wells to reach room temperature (18 to 30°C) before removal from the bag. Place unused wells in the sealable storage bag provided and return to 2 to 8°C.

SAMPLE DIL

2. Sample Diluent

One bottle containing 36 ml of a green/brown buffer solution and detergents. Contains 0.05% ProClin® preservative.

CONJUGATE

3. Conjugate

One bottle (9E25-01) or two bottles (9E25-02) containing a red freeze dried pellet of HIV antigen Conjugated to HRP. When reconstituted each bottle is sufficient for up to three plates.

CONJUGATE DIL

4. Conjugate Diluent

One bottle (9E25-01) or two bottles (9E25-02) each containing 18 ml of a yellow solution consisting of buffer, bovine protein and detergent. Contains 0.1% ProClin® preservative.

Reconstitution of Conjugate

Tap the bottle of Conjugate gently on the bench to remove any material adhering to the rubber stopper. Pour the whole contents of a bottle of conjugate diluent into a bottle of conjugate, recap the latter and mix by gentle inversion. Allow to rehydrate for at least 30 minutes with occasional swirling. The reconstituted conjugate will be red in colour.

After reconstitution the Conjugate may be stored at 2 to 8°C for up to eight weeks or at -15°C or colder for up to four weeks with up to four freeze/thaw cycles.

CONTROL 1 + !

5. Anti-HIV-1 Positive Control Serum

One bottle containing 1.7 ml of inactivated human serum in a buffer containing protein. Contains 0.05% Bronidox® preservative.

CONTROL 2 + !

6. Anti-HIV-2 Positive Control Serum

One bottle containing 1.7 ml of inactivated human serum in a buffer containing protein. Contains 0.05% Bronidox® preservative.

CONTROL - !

7. Negative Control

One bottle containing 2.5 ml of normal human serum diluted in a bovine protein buffer. Contains 0.05% Bronidox® preservative.

SUBSTRATE DIL

8. Substrate Diluent

One bottle containing 35 ml of a colourless solution of tri-sodium citrate and hydrogen peroxide.

SUBSTRATE	CONC
-----------	------

9. Substrate Concentrate

One bottle containing 35 ml of 3,3',5,5'-tetramethylbenzidine (TMB) and stabilisers in a pink solution.

Substrate Solution

To prepare the Substrate Solution add a volume of colourless Substrate Diluent to an equal volume of pink Substrate Concentrate in either a clean glass or plastic vessel. **It is important that this order of addition is followed and that any pipettes and glassware used to prepare Substrate Solution are clean.**

Alternatively, the Substrate Solution may be made by pouring the entire contents of the bottle of Substrate Diluent into the bottle of Substrate Concentrate.

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Additional reagent may be required for use with automated systems. Keep away from sunlight. The Substrate Solution should be pink; if it is purple before being used it should be discarded and fresh Substrate Solution prepared.

The prepared Substrate Solution from this kit may be used interchangeably with that from all other Murex kits which use pink coloured Substrate Concentrate. Ensure that the Substrate Solution is prepared from the Substrate Diluent and Substrate Concentrate provided together.

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10. Wash Fluid

One or two bottles containing 125 ml of 20 times working strength Glycine/Borate Wash Fluid. Contains 0.2% Bronidox® preservative.

Add one volume of Wash Fluid Concentrate to 19 volumes of distilled or deionised water to give the required volume or dilute the entire contents of one bottle of Wash Fluid to a final volume of 2500 ml. Crystals may be observed in the Wash Fluid Concentrate but these crystals will dissolve when the Wash Fluid is diluted to working strength. When diluted the Wash Fluid contains 0.01% Bronidox® preservative.

The Wash Fluid from this kit may be used interchangeably with Glycine/Borate Wash Fluid from any other Murex kit.

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The human sera used for the manufacture have been screened and found reactive or non-reactive for analytes as shown in **Table 2** below.

Table 2

Component	Reactive for	Non-reactive for
Negative Control	N/A	HBsAg, Antibodies to HIV (types 1 and 2), HCV and HTLV
Positive Control 1	Antibodies to HIV-1	HBsAg
Positive Control 2	Antibodies to HIV-2	HBsAg

Additionally human sera used for positive controls are also tested for antibodies to HCV and may be reactive.

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Reagents:	CONJUGATE DIL	SAMPLE DIL	CONJUGATE *
Classification:	Skin sens. 1 H317		
Signal Word:	Warning		
Symbols / Pictograms:			
Hazard Statements:	H317 May cause an allergic skin reaction		
Precautionary Statements:	P280 Wear protective gloves/protective clothing/eye protection/face protection. P363 Wash contaminated clothing before reuse. P333+P313 If skin irritation or rash occurs: Get medical advice / attention		
Contains:	Reaction mass of: 5-chloro-2-methyl-4-isothiazol-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1).		
* The reconstituted Conjugate contains 0.1% ProClin® 300 which is classified hazardous per EC Regulation 1272/2008			

Reagents:	SUBSTRATE CONC
Classification:	Eye Irrit. 2 H319
Signal Word:	Warning
Symbols / Pictograms:	
Hazard Statements:	H319 Causes serious eye irritation
Precautionary Statements:	P264 Wash hands thoroughly after handling P280 Wear protective gloves/protective clothing/eye protection/face protection P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

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- Potentially contaminated materials should be disposed of safely according to local requirements.
- Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued⁷. Sodium hypochlorite should not be used on acid containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- Neutralised acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.
- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- The following reagents contain low concentrations of harmful or irritant substances:
 - The Conjugate Diluent and Sample Diluent contain ProClin® 300 which can be absorbed through the skin and is a sensitising agent.
- Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If either come into contact with the skin or eyes, wash thoroughly with water.
- If any of the reagents come into contact with the skin or eyes wash the area extensively with water.

ANALYTICAL PRECAUTIONS

- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the **Test Procedure** or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
- Allow all reagents and samples to come to 18 to 30°C before use. Immediately after use return reagents to the recommended storage temperature.
- Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionised water.
- Avoid the use of self-defrosting freezers for the storage of reagents and samples.
- Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
- Do not allow wells to become dry during the assay procedure.
- Do not cross-contaminate reagents. Dedicate a pipette for use with the Substrate Solution of Murex assays. A pipette should also be dedicated for use with the Conjugate.
- Do not touch or splash the rim of the well with Conjugate. Do not blow out from micropipettes; reverse pipetting is recommended whenever possible.

- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate.
- Do not contaminate microwells with dust from disposable gloves.
- When using fully automated processors
 - It is not necessary to use plate lids and to tap dry the wells.
 - Do not allow system fluids from fully automated microplate processors to contaminate the samples or reagents.
 - The possibility of cross contamination between assays needs to be excluded when validating assays on fully automated processors.
- Ensure the assay is run within the temperature limits defined in the assay protocol.
- Do not use CO₂ Incubators.
- Do not store the Stop Solution in a shallow dish or return it to a stock bottle after use.
- The possibility of cross contamination between assays needs to be excluded when validating assay protocols on instrumentation.

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SPECIMEN TRANSPORT AND STORAGE

Store the samples at 2 to 8°C. Samples not required for assay within 7 days should be removed from the clot or cell pellet and stored frozen (-15°C or colder). Avoid multiple freeze-thaw cycles. After thawing ensure samples are thoroughly mixed before testing.

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- Stop Solution (0.5M to 2M Sulphuric Acid)**. e.g add between 3 ml (for 0.5M) and 11 ml (for 2.0M) of analytical grade concentrated sulphuric acid (18.0M) to about 80 ml of distilled or deionised water and then make up to 100 ml with more water. Alternatively, the following reagent can be used: 1N Sulphuric Acid (Code N0164 - 15 vial pack and N0165 - 1 vial pack).
- Freshly distilled or high quality deionized water** is required for dilution of Wash Fluid, for preparation of the Stop Solution and for use in conjunction with automated washers.
- Micropipettes and Multichannel micropipettes** of appropriate volume.
- Incubator** capable of maintaining the temperature limits defined in the assay protocol.
- Moulded Heating Block** (Code 5F09-02). For use in laboratory incubators. The moulded heating block should ideally be kept in the incubator used. If this is not possible it must be placed in the incubator at least four hours before beginning the assay.
- Instrumentation**
 - Automated microplate strip washer.
 - Microplate reader.
or
 - Fully automated microplate processor.
All instruments must be validated before use. Please contact your representative for details of recommended systems, software protocols for instrumentation and validation procedures.
- Disposable Reagent Troughs**. (Code 5F24-01).
- Sodium hypochlorite** for decontamination (Refer to **Health and Safety Information**).
- Sodium hydroxide solution** (0.1M) (for instrument decontamination)

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Addition of the various components of the assay to the wells may be confirmed visually by examining the plate for the following colours.

Sample Diluent is green/brown in colour. On addition of Sample or Control the colour will change to blue/green. The colour change will vary from sample to sample but some change should always be visible.

Reconstituted Conjugate is red in colour.

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The addition of sample or reagent can be confirmed using a microplate reader as follows: Sample Diluent plus Sample read at 570 or 620 nm with a reference at 690 nm, Conjugate at 490 nm with a reference at 690 nm, Substrate Solutions at 490 nm (no reference).

Wells containing Sample Diluent may be left for up to 30 minutes at 18 to 30°C prior to addition of Sample, and for up to 15 minutes after the addition of Sample and Controls before starting **Step 5**

SEMI AUTOMATED PROCESSING

Step 1	Reconstitute and mix the Conjugate , prepare the Substrate Solution and Wash Fluid .	
Step 2	Use only the number of wells required for the test. Avoid touching the tops or bottoms of the wells.	
Step 3	Add 50 µl of Sample Diluent to each well.	50 µl
Step 4	Add 50 µl of Samples or 50 µl Controls to the wells. For each plate use the first column of wells for the assay Controls. Add the Controls to the designated wells after dispensing the samples. Pipette 50 µl of the Negative Control into each of three wells A1 to C1 and 50 µl of the anti-HIV-1 and HIV-2 Positive Controls into wells D1 and E1 respectively. Use of a white background will aid visualisation of sample addition.	50 µl
Step 5	Cover the wells with the lid and incubate for 30 minutes at 37°C ±1°C.	30 mins
Step 6	At the end of the incubation time wash the plate as described under Wash Procedures .	
Step 7	Immediately after washing the plate, add 50 µl of Conjugate to each well.	50 µl
Step 8	Cover the wells with the lid and incubate for 30 minutes at 37°C ±1°C.	30 mins
Step 9	At the end of the incubation time wash the plate as described under Wash Procedures .	
Step 10	Immediately after washing the plate, add 100 µl of Substrate Solution to each well.	100 µl
Step 11	Cover the wells with a lid and incubate for 30 minutes at 37°C ±1°C. Keep away from direct sunlight. A purple colour should develop in wells containing reactive samples.	30 mins
Step 12	Add 50 µl of Stop Solution (0.5M to 2M sulphuric acid) to each well.	50 µl
Step 13	Within 15 minutes read the absorbance at 450 nm using 620 nm to 690 nm as the reference wavelength if available. Blank the instrument on air (no plate in the carriage).	A ₄₅₀

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a) Protocol for automated microplate stripwasher.

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- (i) Flow-through washing with a fill volume of 500 µl/well is used with instrumentation supplied by DiaSorin. When using other instrumentation for which this is not possible, ensure that the well is completely filled.

- (ii) The dispense height is set to completely fill the well with a slight positive meniscus, without causing an overflow.
- (iii) The time taken to complete one aspirate/wash/soak cycle is approximately 30 seconds.
- (iv) Ensure that no liquid is left in the well (by use of a double aspirate step in the final cycle where possible).
- (v) After washing is completed, invert the plate and tap out any residual Wash Fluid onto absorbent paper.

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Contact your representative for details of currently available validated protocols. For instrumentation without established validated protocols, the following guidelines are recommended.

1. Do not programme times shorter than specified in the procedure.
2. Incubation times between 30 and 35 minutes (32.5 ± 2.5 minutes) may be programmed.
3. Wells containing Sample Diluent may be left for up to 30 minutes at 18 to 30°C prior to addition of Sample, and for up to 15 minutes after the addition of Sample and Controls before starting **Step 5**.
4. Ensure all **Analytical Precautions** are followed.
Protocols written following these guidelines must be fully validated prior to use according to local procedures.

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Example:

Well 1	=	0.084	, Well 2 =	0.086	, Well 3 =	0.070
Total	=	0.240				
Mean Negative Control	=	0.240/3				
	=	0.080				

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Cut-off value

Calculate the Cut-off value by adding 0.2 to the mean of the Negative Control replicates (see above).

Mean Negative Control	=	0.080				
Cut-Off Value	=	0.080	+	0.200	=	0.280

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Negative Control

The mean absorbance must be less than 0.3.

Positive Control

The absorbance of each of the Positive Controls should be more than 0.8 above the mean absorbance of the Negative Control. Assays which do not meet these criteria should be repeated.

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Such samples should be retested in duplicate using the original source. Samples that are reactive in at least one of the duplicate retests are considered repeatedly reactive in Murex HIV-1.2.O and are presumed to contain antibodies to HIV-1 or HIV-2. Such samples should be further investigated and the presence of antibodies against HIV confirmed by other tests. Samples that are non-reactive i.e. with an absorbance less than that of the Cut-off value, should be considered non-reactive for HIV antibodies.

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The performance of Murex HIV-1.2.O has been determined by testing samples from routine blood donors, patients with AIDS diagnosed according to CDC criteria, patients with AIDS Related Complex (ARC), other patients with known antibody to HIV-1 (including group O), patients with confirmed HIV-2 infection and patients at risk of HIV infection or in other clinical categories. In addition, its performance on commercially available seroconversion panels has been evaluated.

1. Donor Samples

The Murex HIV-1.2.O assay demonstrated a specificity of $\geq 99.5\%$ in a study where a total of 8,225 routine donor samples were screened with Murex HIV-1.2.O at four blood transfusion centres. In the study, 99.91% (8218/8225) of samples were non-reactive and 0.09% (7/8225) were repeatedly reactive. None of the repeatedly reactive samples were confirmed as positive for the presence of antibody to HIV-1 or HIV-2.

The specificity of Murex HIV-1.2.O on presumed negative samples from European donors is estimated to be 99.91% (8218/8225) with a 95% confidence interval of 99.82% (8211/8225) to 99.97% (8223/8225) by the binomial distribution.*

2. Clinical Samples

A total of 1,895 specimens were tested with Murex HIV-1.2.O. The samples were collected at a number of European retrovirology laboratories and included:

- 1,455 samples where antibody was demonstrated by alternative format enzyme immunoassays and/or Western blot. The samples comprised 1,254 samples from patients at various stages of HIV-1 infection (including 58 samples believed to be from patients with group O infection) and 201 samples from patients infected with HIV-2.
- 440 samples from patients with other diseases unrelated to HIV.

In this limited study, Murex HIV-1.2.O detected antibody to HIV in all 1,455 samples in group a) and none of the 440 samples in group b). The relative sensitivity of Murex HIV-1.2.O on these populations of samples was 100% (1455/1455) with a 95% confidence interval of 99.75% (1451/1455) to 100% by the binomial distribution.

A total of 18 HIV-1 seroconversion panels were tested with Murex HIV-1.2.O. Using the presence of both core (p24) and an envelope (gp120/160) band on Western blot as the reference criteria, Murex HIV-1.2.O detected antibody to HIV earlier or in the same sample as Western blot in all of the panels.

3. Assay Reproducibility

The intra-lot and inter-lot reproducibility of Murex HIV-1.2.O has been evaluated, using the assay controls and quality control panel members, at three laboratories. The results of the study are summarised in Table 3 and Table 4.

*Representative performance data are shown: results obtained at individual laboratories and with different populations may vary.

Table 3

Murex HIV-1.2.O - Intra-lot Assay Reproducibility (Mean Absorbance and Coefficient of Variation)

Centre	Kit Negative Control	Kit HIV-1 Positive Control	Kit HIV-2 Positive Control	QC1362	QC1366
Site 1	0.157 (CV 17.5%)	1.465 (CV 17.3%)	1.399 (CV 17.1%)	0.987 (CV 19.9%)	1.300 (CV 23.4%)
Site 2	0.115 (CV 6.7%)	1.372 (CV 9.7%)	1.344 (CV 13.4%)	0.934 (CV 15.3%)	1.240 (CV 14.7%)
Site 3	0.125 (CV 6.8%)	1.203 (CV 5.2%)	1.069 (CV 6.3%)	0.902 (CV 6.3%)	0.971 (CV 7.0%)

Table 4
Murex HIV-1.2.O - Inter-lot Assay Reproducibility

	QC1362	QC1363	QC1364	QC1376
Mean	0.79	0.55	0.73	0.22
CV	10.1%	8.8%	11.7%	5.0%

LIMITATIONS OF PROCEDURE

- The **Test Procedure** and **Interpretation of Results** must be followed.
- This test has only been evaluated for use with individual (unpooled) serum, EDTA plasma or citrate plasma samples. Murex HIV-1.2.O has not been evaluated for any other purpose.
- A negative result with an antibody detection test does not preclude the possibility of infection.
- Non-repeatable reactive results may be obtained with any EIA procedure.
- The most common sources of error are:
 - Imprecise delivery of Sample, Conjugate or Substrate into the wells.
 - Contamination of Substrate with Conjugate.
 - Contamination with conjugates from other assays.
 - Blocked or partially blocked washer probes.
 - Insufficient aspiration leaving a small volume of Wash Fluid in the wells.
 - Failure to ensure that the bottom surface of the wells is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before a plate is read.
 - Failure to read at the correct wavelength or use of an incorrect reference wavelength.
- A specimen collected in citrate may give a raised absorbance.
- The use of highly haemolysed samples, incompletely clotted sera, plasma samples containing fibrin or samples with microbial contamination may give rise to erroneous results.
- This test has not been evaluated for use with samples from cadavers.

BIBLIOGRAPHY

- Clavel, F. (1987). HIV-2: the West African Aids virus. *AIDS*, **1**, 135.
- Denis, F., Leonard, G., et al., (1988). Comparison of 10 enzyme immunoassays for the detection of antibody to human immunodeficiency virus type 2 in West African sera. *J. Clin. Microbiol.* **26**, 1000.
- Vanden Haesevelde, M., Decourt, J., et al., (1994). Genomic cloning and complete sequence analysis of a highly divergent African human immunodeficiency virus isolate. *J. Virol.* **68**, 1586.
- Gürtler, L.G., Hauser, P.H., et al., (1994). A new subtype of human immunodeficiency virus type 1 (MVP-5180) from Cameroon. *J. Virol.* **68**, 1581.
- Loussert-Ajaka, I., Ly, T.D., et al., (1994). HIV-1/HIV-2 seronegativity in HIV subtype O infected patients. *Lancet*, **343**, 1393.
- Gains, H., von Sydons, H., et al., (1988). Detection of immunoglobulin M antibody in primary human immunodeficiency virus infection. *AIDS*, **2**, 11.
- Centres for Disease Control. (1985). Recommendations for preventing transmission of infection with human T-lymphotropic virus type III/lymphadenopathy-associated virus in the workplace. *MMWR*, **34**, No. 45, 681.

Bronidox® and ProClin® are not trade marks of DiaSorin.

DiaSorin S.p.A. UK Branch
Central Road,
Dartford DA1 5LR
UK

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Revised July, 2022

Murex anti-HCV (version 4.0)

An enzyme immunoassay for the detection of antibodies to hepatitis C virus (HCV) in human serum or plasma

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INTENDED USE

Murex anti-HCV (version 4.0) is an enzyme immunoassay for the detection of antibodies to hepatitis C virus (HCV) in human serum or plasma. It is intended to provide improved early detection of HCV infection in blood screening and clinical settings.

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Hepatitis C (HCV) is now recognised as the primary cause of transfusion associated non A non B hepatitis¹. HCV is a single stranded positive-sense RNA virus sharing similarity with flaviviruses and pestiviruses^{2,3} and is global in distribution. Although the acute presentation of HCV infection is generally mild, often clinically asymptomatic, with only 10 to 25% of patients developing jaundice, greater than 50% of infected individuals go on to develop chronic hepatitis with serious and possibly life threatening sequelae such as cirrhosis and hepatocellular carcinoma^{4,5}. Estimates of HCV prevalence in blood donors world-wide range from 0.5 to 8%⁶, although with improved sensitivity and specificity of the diagnostic tests and donor selection, declining incidence has been reported. Data suggests that the incidence is from 0.1 to 1.5% in Europe and 0.6% in the USA⁷.

Diagnosis of HCV is dependent on the direct detection of viral RNA by PCR or by detection of anti-HCV antibodies. Recombinant DNA techniques have been used to develop structural and non-structural proteins derived from HCV RNA with utility for testing of antibodies to HCV. Anti-HCV assays have evolved from first generation products incorporating NS4 proteins only through to third generation assays incorporating core (structural), NS3 protease/helicase (non-structural), NS4 (non-structural) and NS5 replicase (non-structural) proteins. Studies report that the third generation assays demonstrate significant improvements in sensitivity, particularly with regard to increased reactivity with the NS3 antigen and earlier detection of seroconversion⁸.

Murex anti-HCV (version 4.0) utilises antigens from the core, NS3, NS4 and NS5 regions of the virus. Antigens have been carefully developed and selected to provide a sensitive and specific diagnostic test.

PRINCIPLE OF THE PROCEDURE

In Murex anti-HCV (version 4.0) diluted specimen is incubated in microwells coated with highly purified antigens which contain sequences from the core, NS3, NS4 and NS5 regions of HCV. During the course of the first incubation any anti-HCV antibodies in the specimen will bind to the immobilised antigens. Following washing to remove unbound material, the captured anti-HCV antibodies are incubated with peroxidase conjugated monoclonal anti-human IgG. During the course of the second incubation the conjugate will bind to antibody immobilised in the first step. After removal of excess conjugate, bound enzyme is detected by the addition of a solution containing 3,3',5,5'-tetramethylbenzidine (TMB) and hydrogen peroxide. A purple colour will develop in the wells which contained anti-HCV reactive specimens.

The enzyme reaction is terminated with sulphuric acid to give an orange colour which is read photometrically. The amount of Conjugate bound, and hence colour, in the wells, is directly related to the concentration of antibody in the specimen.

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COATED WELLS

1. Coated Wells

One plate (7F5156) or five plates (7F5157) each made up of 96 microwells coated with purified HCV antigens.

Allow the wells to reach 18 to 30°C before removing from the bag. If less than the whole plate is being used, reseal unused strips in the original foil storage bag or in the plastic snap seal bag provided along with the desiccant sachet and return to 2 to 8°C for up to 6 months after the first opening.

SAMPLE DIL

2. Sample Diluent

One bottle of 20 ml (7F5156), or one bottle of 100 ml (7F5157) of buffer containing proteins of bovine origin. Contains 0.05% Bronidox[®] and 0.1% Sodium azide as preservatives.

CONTROL - ⚠

3. Negative Control

One bottle containing 0.8 ml of normal human serum diluted in a buffered solution containing protein of bovine origin. Contains 0.05% Bronidox[®] preservative.

CONTROL + ⚠

4. Anti-HCV Positive Control

One bottle containing 0.6 ml of inactivated human serum consisting of antibodies to HCV diluted in a buffered solution containing protein of bovine origin. Contains 0.05% Bronidox[®] preservative.

CONJUGATE DIL

5. Conjugate Diluent

One bottle (7F5156) or three bottles (7F5157) containing 20 ml of buffered solution consisting of inorganic salts and bovine protein with 0.05% Bronidox[®] preservative.

CONJUGATE

6. Conjugate

One bottle (7F5156) or three bottles (7F5157) each containing enough freeze dried horseradish peroxidase-labelled mouse monoclonal antibody to human IgG in a bovine protein base to perform 192 tests.

Reconstitute at least 15 minutes prior to use to ensure complete dissolution. Bring a bottle of Conjugate Diluent to room temperature. Tap the bottle of Conjugate gently on the bench to remove any material adhering to the rubber stopper. Carefully remove the stopper and pour the Conjugate Diluent into the bottle. Recap and allow to stand with occasional swirling and inversion.

After reconstitution, store at 2 to 8°C for up to seven days or frozen (-15°C or colder) in aliquots for up to six months. The reconstituted Conjugate can be freeze/thawed up to four times. Keep away from sunlight.

SUBSTRATE	DIL
-----------	-----

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Alternatively, the Substrate Solution may be made by pouring the entire contents of the bottle of Substrate Diluent into the bottle of Substrate Concentrate. One bottle of Substrate Solution provides sufficient reagent for at least five plates – see Table 1.

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9. Wash Fluid

One (7F5156) or two bottles (7F5157) containing 125 ml of 20 times working strength Glycine/Borate Wash Fluid. Contains 0.2% Bronidox® preservative.

Add one volume of Wash Fluid concentrate to 19 volumes of distilled or deionised water to give the required volume or dilute the entire contents of one bottle of Wash Fluid to a final volume of 2500 ml. Crystals may be observed in the Wash Fluid Concentrate but these crystals will dissolve when the Wash Fluid is diluted to working strength. When diluted the Wash Fluid contains 0.01% Bronidox® preservative.

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The human sera used for manufacture have been screened and found reactive or non-reactive for analytes as shown in Table 2 below.

Table 2

Component	Reactive for	Non-reactive for
Negative Control	N/A	antibodies to HCV, HIV (types 1 and 2), and HBsAg
Positive Control	antibodies to HCV	HBsAg and antibodies to HIV (types 1 and 2)

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Pursuant to EC Regulation 1272/2008 (CLP) hazardous reagents are classified and labelled as follow:

REAGENTS:	WASH FLUID
SIGNAL WORD:	EUH210 Safety Data Sheet available on request.
SYMBOLS / PICTOGRAMS:	---
HAZARD STATEMENTS:	---
PRECAUTIONARY STATEMENTS:	---
CONTAINS:	---

REAGENTS:	SUBSTRATE CONC
SIGNAL WORD:	Warning
SYMBOLS / PICTOGRAMS:	
HAZARD STATEMENTS:	H319 Causes serious eye irritation.
PRECAUTIONARY STATEMENTS:	---
CONTAINS:	P264 Wash hands thoroughly after handling. P280 Wear protective gloves / protective clothing / eye protection / face protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

REAGENTS:	SAMPLE DIL
SIGNAL WORD:	---
SYMBOLS / PICTOGRAMS:	---
HAZARD STATEMENTS:	EUH032 Contact with acids liberates very toxic gas.
PRECAUTIONARY STATEMENTS:	---
CONTAINS:	---

REAGENTS:	CONTROL + 
SIGNAL WORD:	EUH210 safety data sheet available on request.
SYMBOLS / PICTOGRAMS:	---
HAZARD STATEMENTS:	---
PRECAUTIONARY	---
STATEMENTS:	---
CONTAINS:	---

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- Potentially contaminated materials should be disposed of safely according to local requirement.
- Spillage of potentially infectious material should be removed immediately with absorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued⁹. Sodium hypochlorite should not be used on acid containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- Neutralised acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.
- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- The following reagents contain low concentrations of harmful substances:
 - The Sample Diluent contains saponin and detergents.
- Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If either come into contact with the skin or eyes, wash thoroughly with water.
- If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

ANALYTICAL PRECAUTIONS

- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the **Test Procedure** or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
- Allow all reagents and specimens to come to 18 to 30°C before use. Immediately after use return all reagents to the recommended storage temperature.
- Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionised water.
- Avoid the use of self-defrosting freezers for the storage of reagents and specimens.
- Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
- Do not allow wells to become dry during the assay procedure.
- Do not cross-contaminate reagents. Dedicate a pipette for use with the Substrate Solution of Murex assays. A pipette should also be dedicated for use with the Conjugate.
- Do not touch or splash the rim of the well with Conjugate. Do not blow out from micropipettes; reverse pipetting is recommended wherever possible.
- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate.
- Do not contaminate microwells with dust from disposable gloves.
- When using fully automated microplate processors:
 - It is not necessary to use plate lids and to tap dry the wells.
 - Do not allow system fluids from fully automated microplate processors to contaminate the specimens or reagents.
 - The possibility of cross contamination between assays needs to be excluded when validating assays on fully automated processors.

- Cross contamination between assays needs to be considered when validating assays protocols on automated microplate processors.
- Ensure the assay is run within the temperature limits defined in the assay protocol.
- Do not use CO₂ Incubators.
- Do not store the Stop Solution in a shallow dish nor return it to a stock bottle after use.
- It is important that specimens and controls are thoroughly mixed with the Sample Diluent. Failure to do this may cause erroneous results. Adequate mixing can be achieved by:
 - Manual mixing by pipetting up and down at least four times when adding specimens or controls.
 - Placing the plate on a microplate shaker at a speed of 800 rpm for 30 seconds.

Users of automatic sample dispensers may find it convenient to add 90 µl of Sample Diluent to the wells, followed by 20 µl of Specimen then the remaining 90 µl of Sample Diluent. This procedure will ensure adequate mixing.

- The possibility of cross contamination between assays needs to be excluded when validating assay protocols on instrumentation.

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SPECIMEN COLLECTION

Serum, EDTA plasma, citrate plasma or heparin plasma specimens may be used. Blood collected by venepuncture should be allowed to clot naturally. Ensure that the serum specimens are fully clotted. Remove any visible particulate matter from the specimen by centrifugation.

If specimens are prepared using liquid anti-coagulant e.g. citrate plasma, the dilution effect should be considered.

SPECIMEN TRANSPORT AND STORAGE

Store specimens at 2 to 8°C. Specimens not required for assay within seven days should be removed from the clot or cell pellet and stored frozen (-15°C or colder). Avoid multiple freeze-thaw cycles. After thawing, ensure specimens are thoroughly mixed before testing.

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- Stop Solution (0.5M to 2M Sulphuric Acid)** e.g. add between 3 ml (for 0.5M) and to 11 ml (for 2.0M) of analytical grade concentrated sulphuric acid (18.0 M) to about 80 ml of distilled or deionised water and then make up to 100 ml with more water. Alternatively, the following reagent can be used: 1N Sulphuric Acid Code N0165 (1 vial pack) or code N0164 (15 vial pack).
- Freshly distilled or high quality deionised water** is required for dilution of Wash Fluid, for preparation of the Stop Solution and for use in conjunction with automated washers.
- Micropipettes and Multichannel micropipettes** of appropriate volume.
- Incubator** capable of maintaining the temperature limits defined in the assay protocol.
- Moulded Heating Block** (code 5F09-02). For use in laboratory incubators. The moulded heating block should ideally be kept in the incubator used. If this is not possible it must be placed in the incubator at least four hours before beginning the assay.
- Instrumentation**
 - Automated microplate strip washer.
 - Microplate reader.

or

 - Fully automated microplate processor.

All instruments must be validated before use.

Please contact your representative for details of recommended systems, software protocols for instrumentation and validation procedures.
- Disposable Reagent Troughs.** (Code No. 5F24-01).
- Sodium hypochlorite** for decontamination. (**Refer to Health and Safety Information**).

TEST PROCEDURE

Please read '**Analytical Precautions**' carefully before performing the test. Addition of the various components of the assay to the wells may be confirmed visually by examining the plate for the following colours.

Sample Diluent is green/brown in colour. On addition of the Specimen or Control the colour will change to blue/green. The colour change may vary from specimen to specimen but some change should always be visible.

Conjugate is brown in colour.

Substrate Solution is initially pink with any positive wells becoming purple. On addition of **Stop Solution** the purple colour of the positives will change to orange, whilst the negatives remain pink.

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The addition of specimen or reagent can be confirmed using a microplate reader as follows: Sample Diluent plus Specimen read at 620 nm or 570 nm with a reference at 690 nm, Conjugate at 410 nm with a reference at 690 nm, Substrate Solution at 490 nm (no reference).

SEMI AUTOMATED PROCESSING

Step 1	Reconstitute the Conjugate with Conjugate Diluent and prepare the Substrate Solution and dilute the Wash Fluid .	
Step 2	Use only the number of strips required for the test.	
Step 3	Add 180 µl of Sample Diluent to each well.	180 µl
Step 4	Add 20 µl of Specimens or Controls to the wells. If using two strips or less pipette one Negative Control to well A1 and one Positive Control to well B1. <u>If using more than two strips it is recommended that two wells of Negative Control are used for additional security.</u> Add the Controls to the designated wells on each plate after dispensing the specimens. The use of a white background will aid visualisation of specimen addition. Specimens and Controls must be thoroughly mixed with the Sample Diluent.	20 µl
Step 5	Cover the wells with the lid and incubate for 1 hour at 37°C ± 1°C.	1 hour
Step 6	At the end of the incubation period wash the plate as described under Wash Procedures .	
Step 7	Immediately after washing the plate, add 100 µl of Conjugate to each well.	100 µl
Step 8	Cover the wells with the lid and incubate for 30 minutes at 37°C ± 1°C.	30 mins
Step 9	At the end of the incubation period wash the plate as described under Wash Procedures .	
Step 10	Immediately after washing the plate, add 100 µl of Substrate Solution to each well.	100 µl
Step 11	Cover the wells with a lid and incubate for exactly 30 minutes at 37°C ± 1°C while colour develops. Keep away from direct sunlight. A purple colour should develop in wells containing reactive specimens.	30 mins
Step 12	Add 50 µl of Stop Solution (0.5M to 2M sulphuric acid) to each well.	50 µl
Step 13	Within 15 minutes read the absorbance at 450 nm using 620 nm to 690 nm as the reference wavelength if available. Blank the instrument on air (no plate in the carriage).	450 nm

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a) **Protocol for automated microplate stripwasher**

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- Flow-through washing with a fill volume of 500 µl/well is used with instrumentation supplied by DiaSorin. When using other instrumentation, for which this is not possible, ensure that the well is completely filled.
- The dispense height is set to completely fill the well with a slight positive meniscus without causing an overflow.
- The time taken to complete one aspirate/wash/soak cycle is approximately 30 seconds.
- Ensure that no liquid is left in the well (by use of a double aspirate step in the final cycle where possible).
- After washing is completed, invert the plate and tap out any residual Wash Fluid onto absorbent paper.

NOTE: Do not allow the wells to become dry during the assay procedure. Washers must be rinsed with distilled water at the end of the test to avoid blockage and corrosion.

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- For the first incubation, incubation times between 60 and 70 minutes (or 65 ± 5 minutes) may be programmed.

- For the 30 minute incubations, incubation times between 30 and 35 minutes (or 32.5 ± 2.5 minutes) may be programmed.
- Ensure all '**Analytical Precautions**' are followed. Protocols written following these guidelines must be fully validated prior to use according to local procedure.

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Discard any Negative Control value which is >0.25.

Cut-off Value

Calculate the Cut-off value by adding 0.6 either to the Negative Control or to the mean of the Negative Control replicates (see above).

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The mean absorbance is less than 0.25.

Positive Control

The absorbance is more than 0.8 above the mean absorbance of the Negative Control.

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INTERPRETATION OF THE RESULTS

Non-Reactive Results

Specimens giving an absorbance less than the Cut-off value are considered negative in the assay.

Reactive Results

Specimens giving an absorbance equal to or greater than the Cut-off value are considered initially reactive in the assay. Such specimens must be repeated in duplicate using the original source. Specimens that are reactive in at least one of the repeat tests are considered repeatedly reactive and are presumed to contain antibody to HCV. Unless local procedures state otherwise these specimens must be further investigated.

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The performance of Murex anti-HCV (version 4.0) has been determined by testing specimens from random blood donors, patients with known antibody to HCV, patients with diseases related to HCV and patients with diseases unrelated to HCV.

In addition, its performance on commercially available seroconversion panels has been evaluated.

1. **Donor Specimens**

A total of 8835 routine donor specimens were screened with Murex anti-HCV (version 4.0) in Europe and Australia.

In the study, 99.82% (8819/8835) of specimens were non-reactive, 0.18% (16/8835) were initially reactive and 0.12% (11/8835) were repeatedly reactive. None of the repeatedly reactive specimens have been confirmed as positive for the presence of antibody to HCV.

The specificity of Murex anti-HCV (version 4.0) on this population of presumed negative specimens is estimated to be 99.88% (8824/8835) with a 95% confidence interval of 99.77% (8815/8835) to 99.94% (8830/8835) by the binomial distribution.

2. Clinical Specimens

A total of 69 specimens with antibody to HCV, confirmed with an alternative immunoassay, Ortho RIBA3 and/or Western blot, were reactive with Murex anti-HCV (version 4.0).

A total of 27 commercially available HCV seroconversion panels were also tested with Murex anti-HCV (version 4.0). Comparison with an alternative commercially available immunoassay for the detection of antibody to HCV showed that Murex anti-HCV (version 4.0) detected antibody two bleeds earlier in four panels, one bleed earlier in three panels, one bleed later in three panels and at the same bleed in 17 panels.

In addition, 873 potentially cross-reactive specimens from patients with conditions unrelated to HCV infection, including other acute viral infections, antenatal, lipaemic, icteric and haemolysed specimens, were tested with Murex anti-HCV (version 4.0). A total of 869 of these specimens were non-reactive with Murex anti-HCV (version 4.0), the remaining four specimens included two which gave indeterminate results with Ortho RIBA3.

3. Assay Reproducibility

Five replicates of each of four specimens were tested on ten separate occasions with two separate lots to ascertain the reproducibility of Murex anti-HCV (version 4.0). The results of the study are summarised in Table 3 and Table 4

Table 3
Murex anti-HCV (version 4.0) - Assay Reproducibility (Lot 1)

Specimen	Number of Assays	Mean Absorbance/ Cut-off Value	Intra-assay %CV	Inter-assay %CV
1	10	4.03	2.9	3.9
2	10	3.9	4.6	5.2
3	10	2.33	3.1	6.9
4	10	0.14	3.9	10.2

Table 4
Murex anti-HCV (version 4.0) - Assay Reproducibility (Lot 2)

Specimen	Number of assays	Mean Absorbance/ Cut-off value	Intra-assay % CV	Inter-assay % CV
1	10	3.22	6.5	9.1
2	10	3.14	7.2	9.2
3	10	1.69	5.2	9.5
4	10	0.13	4.0	8.4

LIMITATIONS OF THE PROCEDURE

- The **Test Procedure** and **Interpretation of Results** must be followed.
 - This test has only been evaluated for use with individual (unpooled) serum, EDTA plasma, heparin plasma or citrate plasma specimens.
 - A non-reactive specimen result with an antibody detection test does not preclude the possibility of infection.
 - Non-repeatable reactive results may be obtained with any EIA procedure.
 - The most common sources of error are:
 - Imprecise delivery of Specimen, Conjugate or Substrate into the wells.
 - Contamination of Substrate with Conjugate.
 - Contamination with conjugates from other assays.
 - Blocked or partially blocked washer probes.
 - Insufficient aspiration leaving a small volume of Wash Fluid in the wells.
 - Failure to ensure that the bottom surface of the wells is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before a plate is read.
 - Failure to read at the correct wavelength or use of an incorrect reference wavelength.
 - The use of highly haemolysed specimens, incompletely clotted sera, plasma specimens containing fibrin or specimens with microbial contamination may give rise to erroneous results.
 - This test has not been evaluated for use with specimens from cadavers.
- See also '**Analytical Precautions**'.

BIBLIOGRAPHY

- Choo, Q.L., Kuo, G., et al., (1989).** Isolation of a cDNA clone derived from a blood-borne non-A non-B viral hepatitis genome. *Science* 244, 359.
- Miller, R.H. and Purcell, R.H. (1990).** Hepatitis C virus shares amino acid sequence similarity with pestiviruses and flaviviruses as well as two plant virus supergroups. *Proc Natl Acad Sci* 87, 2057.
- Weiner, A.J., Brauer, M.J. et al., (1991).** Variable and hypervariable domains are found in the regions of HCV corresponding to the flavivirus envelope and NS1 proteins and the pestivirus envelope glycoproteins. *Virology* 180, 842.
- Jove, J., Sanchez-Tapias, J.M. et al., (1990).** Post-transfusional vs. sporadic non-A, non-B chronic hepatitis. A clinicopathological and evolutive study. *Liver*, 8, 42.
- Hopf, U., Möller, B. et al., (1990).** Long term follow up of post transfusion and sporadic chronic hepatitis non-A, non-B and frequency of circulating antibodies to hepatitis C virus (HCV). *Hepatology*, 10, 69.
- Conroy-Cantilena, C. (1997).** Hepatitis C virus diagnostics: technology, clinical applications and impacts. *Trends Biotech* 15, 71.
- Botte, C., Janot, C. (1996).** Epidemiology of HCV infection in the general population and in blood transfusion. *Nephrol Dial Transplant* 11 (Suppl 4), 19.
- Courouce, A.M. (1998).** Development of screening and confirmation tests for antibodies to Hepatitis C virus. Reesink, H.W. (ed) Hepatitis C virus – *Curr Stud Hematol Blood Transfus.* 62, 64, Basel, Karger.
- Centres for Disease Control. (1985).** Recommendations for preventing transmission of infection with human T-lymphotropic virus type III / lymphadenopathy-associated virus in the workplace. *MMWR*, 34, No. 45, 681.

Bronidox® is not a trade mark of DiaSorin

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LITERATŪRA

1. Choo, Q.L., Kuo, G., et al., (1989). Isolation of a cDNA clone derived from a blood-borne non-A non-B viral hepatitis genome. *Science* 244, 359.
2. Miller, R.H. and Purcell, R.H. (1990). Hepatitis C virus shares amino acid sequence similarity with pestiviruses and flaviviruses as well as two plant virus supergroups. *Proc Natl Acad Sci* 87, 2057.
3. Weiner, A.J., Brauer, M.J. et al., (1991). Variable and hypervariable domains are found in the regions of HCV corresponding to the flavivirus envelope and NS1 proteins and the pestivirus envelope glycoproteins. *Virology* 180, 842.
4. Jove, J., Sanchez-Tapias, J.M. et al., (1990). Post-transfusional vs. sporadic non-A, non-B chronic hepatitis. A clinicopathological and evolutive study. *Liver*, 8, 42.
5. Hopf, U., Möller, B. et al., (1990). Long term follow up of post transfusion and sporadic chronic hepatitis non-A, non-B and frequency of circulating antibodies to hepatitis C virus (HCV). *Hepatology*, 10, 69.
6. Conroy-Cantilena, C. (1997). Hepatitis C virus diagnostics: technology, clinical applications and impacts. *Trends Biotech* 15, 71.
7. Botte, C. Janot, C. (1996). Epidemiology of HCV infection in the general population and in blood transfusion. *Nephrol Dial Transplant* 11 (Suppl 4), 19.
8. Courouce, A.M. (1998). Development of screening and confirmation tests for antibodies to Hepatitis C virus. Reesink, H.W. (ed) Hepatitis C virus – *Curr Stud Hematol Blood Transfus.* 62, 64, Basel, Karger.
9. Centres for Disease Control. (1985). Recommendations for preventing transmission of infection with human T-lymphotropic virus type III / lymphadenopathy-associated virus in the workplace. *MMWR*, 34, No. 45, 681.

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 DiaSorin Italia S.p.A. UK Branch
Central Road,
Dartford, DA1 5LR
UK

 DiaSorin Italia S.p.A.
Via Crescentino snc
13040 Saluggia (VC) Italy

D13DS47LT
2022 m. liepos mėn.

Revised July, 2022

Murex anti-HBc (total)

Enzyme immunoassay for the detection of antibodies to hepatitis B core antigen (anti-HBc) in human serum or plasma

The assay is intended to screen individual human donors for the presence of anti-HBc or as an aid to the diagnosis of HBV infection.

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INTENDED USE

Murex anti-HBc (total) is intended for the detection of antibodies to hepatitis B core antigen (anti-HBc) in human serum or plasma.

The assay is intended to screen individual human donors for the presence of anti-HBc or as an aid to the diagnosis of HBV infection.

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Murex anti-HBc (total) is a competitive enzyme immunoassay for the detection of antibodies against hepatitis B core antigen (anti-HBc). Anti-HBc titres increase rapidly after exposure to HBV¹ and can persist for many years at lower titres after clearance of hepatitis B surface antigen (HBsAg) and resolution of infection². The presence of anti-HBc is therefore indicative both of past exposure to HBV and infection in the acute/chronic period. During active infection both immunoglobulin M (IgM) and immunoglobulin G (IgG) anti-HBc are usually present and these may be the only serological markers of HBV infection during the "window period" when HBsAg has cleared but before antibodies to HBsAg are detectable³. Murex anti-HBc (total) detects IgG and IgM antibodies to hepatitis B core antigen.


PRINCIPLE OF THE PROCEDURE

Murex anti-HBc (total) is based on microwells coated with recombinant hepatitis B core antigen (HBcAg). Samples and controls are incubated in the wells and any anti-HBc present in the sample or control binds to the HBcAg. Any excess antibody is then removed by washing. Conjugate (monoclonal anti-HBc conjugated to horseradish peroxidase) is then added to the wells. A second incubation is carried out during which the conjugate binds to any HBcAg on the well surface not blocked by anti-HBc from the test sample. After washing to remove any unbound conjugate, a solution containing 3, 3', 5, 5'-tetramethylbenzidine (TMB) and hydrogen peroxide is added to the wells. Wells which do not contain anti-HBc and therefore bind conjugate, will develop a blue/green colour which is converted to orange when the enzyme reaction is stopped with sulphuric acid. The intensity of colour can be determined spectrophotometrically. The intensity of colour is greatest in the absence of anti-HBc and falls with increasing concentrations of anti-HBc in the sample.

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COATED WELLS 1. Coated Wells

One plate (8G21-01) or 5 plates (8G21-02) of 96 wells coated with recombinant HBcAg.

Allow the wells to reach 18 to 30°C before removing from the bag. place any unused wells in the sealable storage bag provided and return to 2 to 8°C.

SAMPLE DIL 2. Sample Diluent

Five bottles (7 mL per bottle) containing a green/yellow buffered solution.

The Sample Diluent is susceptible to oxidation so the bottle must be kept closed with a bung and a cap whenever possible.

After first opening, the Sample Diluent is stable at 2 to 8°C for up to one week. The bottle must not be opened for use on more than three occasions.

CONJUGATE 3. Conjugate

One bottle (8G21-01) or five bottles (8G21-02) (7 mL) containing monoclonal anti-HBc conjugated to horseradish peroxidase in a red buffered solution. Contains 0.05% ProClin[®] 300 preservative.

CONTROL + ⚠ 4. Positive Control

One bottle containing 2 mL of anti-HBc diluted in a buffered diluent containing a blue dye. Contains 0.05% Bronidox[®] preservative.

CONTROL - ⚠ 5. Negative Control

One bottle (2 mL) of human serum containing a green dye. Contains 0.05% Bronidox[®] preservative.

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Substrate Solution

To prepare the Substrate Solution add a volume of colourless Substrate Diluent to an equal volume of orange Substrate Concentrate in either a clean glass or plastic vessel. **It is important that this order of addition is followed and that any pipettes and glassware used to prepare Substrate Solution are clean.** Alternatively, the Substrate Solution may be made by pouring the entire contents of the bottle of Substrate Diluent into the bottle of Substrate Concentrate. On addition of the Diluent, the Concentrate will change colour from orange to yellow. One bottle of Substrate Solution provides sufficient reagent for at least five plates - see **Table 1**.

Tabela 1
Volume of Substrate Concentrate and Substrate Diluent Required

Number of Wells	Number of Plates
8 16 24 32 40 48 56 64 72 80 88	1 2 3 4
Substrate Concentrate (mL)	
1.0 1.5 2.0 2.5 2.5 3.0 3.5 4.0 4.5 4.5 5.0	6 12 18 22
Substrate Diluent (mL)	
1.0 1.5 2.0 2.5 2.5 3.0 3.5 4.0 4.5 4.5 5.0	6 12 18 22

Additional reagent may be required for use with automated systems. Keep away from all natural and artificial light. The Substrate Solution should be yellow; if it is green before being used, it should be discarded and fresh Substrate Solution prepared.

The prepared Substrate Solution from this kit may be used interchangeably with that from all other Murex kits which use orange coloured Substrate Concentrate. Ensure that the Substrate Solution is prepared from Substrate Diluent and Substrate Concentrate provided together.

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8. Wash Fluid

One bottle (8G21-01) or two bottles (8G21-02) containing 125 mL of 20 times working strength Glycine/Borate Wash fluid. Contains 0.2% Bronidox® preservative.

Add one volume of Wash fluid concentrate to 19 volumes of distilled or deionized water to give the required volume or dilute the entire contents of one bottle of Wash fluid to a final volume of 2500 mL. Crystals may be observed in the Wash fluid concentrate but these crystals will dissolve when the Wash fluid is diluted to working strength. When diluted, the Wash fluid contains 0.01% Bronidox® preservative.

The Wash fluid from this kit may be used interchangeably with the Glycine/Borate Wash fluid from any other Murex kit.

Store the working strength Wash fluid at 18 to 30°C in a closed vessel under which conditions it will retain activity for one month.

NOTE: The Wash fluid may develop a yellow colour on storage. This will have no effect on the performance of the assay providing the Wash fluid is fully aspirated from the wells.

NOTE: Although the Substrate Solution and Wash fluid are interchangeable, they must not be used beyond the expiry date printed on the component labels.

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HEALTH AND SAFETY INFORMATION

CAUTION: This kit contains components of human origin. The human sera used for manufacture have been screened and found reactive or non-reactive for analytes as shown in **Table 2** below.

Table 2

Component	Reactive for	Non-reactive for
Negative Control	N/A	HBsAg, and antibodies to HIV (types 1 and 2), and HCV
Positive Control	HBsAg	Antibodies to HIV (types 1 and 2), and HCV

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Reagents:	CONJUGATE
Classification:	Skin sens. 1 H317
Signal Word:	Warning
Symbols / Pictograms:	
Hazard Statements:	H317 May cause an allergic skin reaction.
Precautionary Statements:	P280 Wear protective gloves/protective clothing/eye protection/face protection. P363 Wash contaminated clothing before reuse. P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
Contains:	Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1).
Reagents:	SAMPLE DIL
Classification:	Eye Irrit. 2 H319
Signal Word:	Warning
Symbols / Pictograms:	
Hazard Statements:	H319 Causes serious eye irritation; EUH031 Contact with acids liberates toxic gas.
Precautionary Statements:	P280 Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337+P313 If eye irritation persists: Get medical advice/attention.
Reagents:	SUBSTRATE CONC
Classification:	Eye Irrit. 2 H319
Signal Word:	Warning
Symbols / Pictograms:	
Hazard Statements:	H319 Causes serious eye irritation
Precautionary Statements:	P264 Wash hands thoroughly after handling. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

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1. Potentially contaminated materials should be disposed of safely according to local requirement.
2. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

3. Neutralised acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
5. The following reagent contains low concentrations of harmful or irritant substances.
 - a) The Sample Diluent contains Saponin.
6. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If either come into contact with the skin or eyes, wash thoroughly with water.
7. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

ANALYTICAL PRECAUTIONS

1. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
2. Do not modify the **Test procedure** or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
3. Allow all reagents and samples to come to 18 to 30°C before use. Immediately after use return all reagents to the recommended storage temperature.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionised water.
5. Avoid the use of self-defrosting freezers for the storage of reagents and samples.
6. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
7. Do not allow wells to become dry during the assay procedure.
8. Do not cross-contaminate reagents. Dedicate a pipette for use with the Substrate Solution of Murex assays. A pipette should also be dedicated for use with the Conjugate.
9. Do not touch or splash the rim of the well with Conjugate. Do not blow out from micropipettes; reverse pipetting is recommended wherever possible.
10. Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate.
11. Do not contaminate microwells with dust from disposable gloves.
12. When using fully automated microplate processors
 - i) It is not necessary to use plate lids and to tap dry the wells.
 - ii) Do not allow system fluids from fully automated microplate processors to contaminate the samples or reagents.
 - iii) The possibility of cross contamination between assays needs to be excluded when validating assays on fully automated processors.
13. Ensure the assay is run within the temperature limits defined in the assay protocol.
14. Do not use CO₂ incubators.
15. Do not store the Stop Solution in a shallow dish or return it to a stock bottle after use.
16. The possibility of cross contamination between assays needs to be excluded when validating assay protocols on instrumentation.

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SPECIMEN TRANSPORT AND STORAGE

Store samples at 2 to 8°C. Samples not required for assay within 7 days should be removed from the clot or cell pellet and stored frozen (-15°C or colder). Avoid multiple freeze-thaw cycles. After thawing, ensure samples are thoroughly mixed before testing.

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1. **Stop Solution (0.5M to 2M Sulphuric Acid)**. Add 3 to 11 mL of analytical grade concentrated sulphuric acid (18.0M) to about 80 mL of distilled or deionized water and then make up to 100 mL with more water. Alternatively, the following reagent can be used: 1N Sulphuric Acid (Code N0164 - 15 vial pack and N0165 - 1 vial pac).
2. **Freshly distilled or high quality deionised water** is required for dilution of Wash fluid, for preparation of the Stop Solution and for use in conjunction with automated washers.
3. **Micropipettes and Multichannel micropipettes** of appropriate volume.
4. **Incubator** capable of maintaining the temperature limits defined in the assay protocol.
5. **Moulded Heating Block** (Code 5f09-02) for use in laboratory incubators. The moulded heating block should ideally be kept in the incubator used. If this is not possible it must be placed in the incubator at least four hours before beginning the assay.
6. **Instrumentation**
 - a) Automated microplate stripwasher,
 - b) Microplate reader, or
 - c) Fully automated microplate processor.
 All instruments must be validated before use. Please contact your representative for details of recommended systems, software protocols for instrumentation and validation procedures.
7. **Disposable Reagent Troughs**. (Code 5F24-01).
8. **Sodium hypochlorite** for decontamination. (Refer to **Health and Safety Information**).
9. **Sodium hydroxide solution** (0.1M). (For instrument decontamination).

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Addition of the various components of the assay to the wells may be confirmed visually by examining the plate for the following colours:

Sample diluent is green/yellow in colour. On addition of Sample or Control the colour will change to blue/green. The colour change will vary from sample to sample but some change should always be visible.

Conjugate is red in colour.

Substrate Solution is initially pale yellow with any non-reactive wells becoming blue/green. On addition of Stop Solution the blue/green colour of the non-reactives will change to orange/yellow, whilst the reactives will change to a pink colour.

The addition of sample or reagent can be confirmed using a microplate reader as follows.

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Sample Diluent plus Sample read at 570 or 620 nm with a reference at 690 nm.
 Conjugate at 490 nm with a reference at 690 nm.
 Substrate Solution at 450 nm (no reference).

SEMI AUTOMATED PROCESSING

Step 1	Use only the number of Wells required for the test.	
Step 2	Prepare the Wash Fluid .	
Step 3	Add 50 µL of Sample Diluent to each well.	50 µL
Step 4	Add 50 µL of Sample or Control to the wells. For each run add Negative Control to Wells A1 and B1 and the Positive Control to Wells C1 and D1.	50 µL
Step 5	Add the Controls to the designated Wells after dispensing the samples. Cover the wells with the lid and incubate for 30 minutes at 37°C ±1°C.	30 mins
Step 6	At the end of the incubation period wash the plate as described under Wash Procedures . After washing is completed, invert the plate and tap out any residual Wash Fluid onto absorbant paper.	
Step 7	Add 50 µL of Conjugate to each Well.	50 µL
Step 8	Cover the wells with the lid and incubate for 30 minutes at 37°C ±1°C.	30 mins
Step 9	Prepare the Substrate Solution .	
Step 10	At the end of the incubation period wash the plate as described under Wash Procedures . After washing is completed invert the plate and tap out any residual Wash Fluid onto absorbent paper.	
Step 11	Immediately after washing the plate, add 100 µL of Substrate Solution to each well.	100 µL
Step 12	Cover the wells with a lid and incubate for 30 minutes at 37°C ±1°C. Keep away from direct sunlight.	30 mins
Step 13	Add 50 µL of Stop Solution .	50 µL
Step 14	Within 15 minutes read the absorbance at 450 nm using 620 nm to 690 nm as the reference wavelength if available. Blank the instrument on air (no plate in the carriage).	450 mm

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- i) Flow-through washing with a fill volume of 500 µL/well is used with instrumentation supplied by DiaSorin. When using other instrumentation for which this is not possible, ensure that the well is completely filled.
- ii) The dispense height is set to completely fill the well with a slight positive meniscus, without causing an overflow.
- iii) The time taken to complete one aspirate/wash/soak cycle is approximately 30 seconds.
- iv) Ensure that no liquid is left in the well (by use of a double aspirate step in the final cycle where possible).

- v) After washing is completed, invert the plate and tap out any residual Wash fluid onto absorbent paper.
- b) Protocol for manual washer
 - i) Aspirate the first row of wells.
 - ii) Completely fill this row with working strength Wash fluid.
 - iii) Repeat this procedure for each row of wells in turn.
 - iv) Ensure that each row of wells is left to soak for 30 seconds.
 - v) Repeat (i) to (iv) a further 4 times.
 - vi) Aspirate the contents of the wells. It is recommended that the wells are inverted and tapped dry on paper towel or tissue after the last wash.

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1. Do not programme times shorter than specified in the procedure.
2. For each incubation at 37°C, programmed times may be increased by up to 5 minutes.
3. Wells containing either Sample Diluent or Sample Diluent and sample/ control may be left at 18 to 30°C (room temperature) for up to 2 hours before starting **Step 4 and 5 respectively**.
4. Ensure all "**Analytical precautions**" are followed. Protocols written following these guidelines must be fully validated prior to use according to local procedures.

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Positive Control

Calculate the mean absorbance of the Positive Controls.

Cut-off value

Calculate the Cut-off value by adding the mean absorbance of the positive Control to the mean absorbance of the Negative Control and divide this value by two.

Example:

Negative Control Absorbance	
Well 1	= 1.652
Well 2	= 1.586
Mean Negative Control	= (1.652 + 1.586)/2 = 1.619
Positive Control absorbance	
Well 3	= 0.040
Well 4	= 0.042
Mean Positive Control	= (0.040 + 0.042)/2 = 0.041
Cut-Off value	= (1.619 + 0.041)/2 = 0.830

QUALITY CONTROL

Results of an assay are valid if the following criteria for the Controls are met:

Negative Control minus Positive Control

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INTERPRETATION OF RESULTS

Non-reactive Results

Samples giving an absorbance greater than the Cut-off value are considered non-reactive in the assay.

Reactive Results

Samples giving an absorbance equal to or less than the Cut-off value are considered reactive in Murex anti-HBc (total). Unless local procedures state otherwise, such samples must be retested in duplicate using the original source. Samples that are reactive in at least one of the duplicate retests are considered repeatedly reactive in Murex anti-HBc (total) and are presumed to contain antibodies to HBc. Such samples must be further investigated and the results from this assay considered with any other clinical and/or assay information. Samples which are non-reactive in both wells on retest must be considered non-reactive in this assay.

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The performance of Murex anti-HBc (total) has been determined by testing samples from random blood donors, patients with known antibody to hepatitis B core antigen and patients in other clinical categories.

Diagnostic Sensitivity

A total of 447 specimens from patients known to contain antibody to hepatitis B core antigen were tested and found to be reactive with Murex anti-HBc (total). The specimens were taken from patients at various stages of HBV infection and included 139 specimens with IgM antibody to hepatitis B core antigen.

The diagnostic sensitivity of Murex anti-HBc (total) on this population of specimens is therefore estimated to be 100% (447/447) with a lower 95% confidence limit of 99.18% by the binomial distribution.

Diagnostic Specificity

A total of 360 specimens from patients with conditions unrelated to HBV infection and confirmed as negative for antibody to hepatitis B core antigen were also tested with Murex anti-HBc (total). These included haemolysed specimens, specimens from pregnant women, patients suffering with autoimmune disease and patients with other acute viral infections.

In addition, a panel of lipaemic and icteric specimens were tested and found to be non-reactive.

The diagnostic specificity of Murex anti-HBc (total) on this population of clinical specimens is estimated to be 100% (360/360) with a lower 95% confidence limit of 98.98%.

Screening Specificity

A total of 5344 routine donor specimens from two European and two South American blood transfusion centres were screened with Murex anti-HBc (total). The results are summarised in Table 3. In the study, 99.14% (5298/5344) of specimens were non-reactive and 0.86%

(46/5344) were repeatedly reactive. Thirty two of the repeatedly reactive specimens were confirmed positive by testing with at least two other commercial assays for the detection of antibody to hepatitis B core antigen. None of the remaining fourteen specimens were confirmed as positive for the presence of antibody to hepatitis B core antigen.

The screening specificity of Murex anti-HBc (total) on this population of negative routine donor specimens is estimated to be 99.74% (5298/5312) with 95% confidence limits of 99.86% to 99.56% by the binomial distribution.

Tabela 3
Reactivity of Murex anti-HBc (total) with presumed negative specimens from routine blood donors

Centre	Number of routine donor specimens tested	Number of repeatedly reactive specimens	Number of confirmed positive specimens	Number of false reactive specimens
A	1975	9	2	7 (0.35%)
B	2081	10	7	1 (0.05%)
C	608	18	16	2 (0.33%)
D	680	11	7	4 (0.59%)
TOTAL	5344	46	32	14 (0.26%)

Assay Reproducibility

Five quality assurance panel members were tested as ten replicates on four separate occasions to assess the reproducibility of Murex anti-HBc (total). The results from the testing are summarised in Table 4.

Tabela 4
Murex anti-HBc (total) - Assay Reproducibility

Specimen	Number of assays	Number of replicates	Cut-off value/ Mean Absorbance	Intra-assay %CV	Inter-assay %CV
QA01	4	10	0.782	6.1	6.7
QA02	4	10	1.263	6.7	8.9
QA03	4	10	2.672	7.4	8.4
QA04	4	10	2.380	7.3	8.7
QA05	4	10	0.982	5.4	6.3

LIMITATIONS OF THE PROCEDURE

1. The **Test Procedure** and **Interpretation of Results** must be followed.
2. This test has only been evaluated for use with individual (unpooled) serum, EDTA plasma or citrate plasma samples.
3. A negative result with an antibody detection test does not preclude the possibility of infection.
4. Non-repeatable reactive results may be obtained with any EIA procedure.
5. The most common sources of error are:
 - a) Imprecise delivery of Sample, Conjugate or Substrate into the wells.
 - b) Contamination of Substrate with Conjugate.
 - c) Contamination with conjugates from other assays.
 - d) Blocked or partially blocked washer probes.
 - e) Insufficient aspiration leaving a small volume of Wash fluid in the wells.
 - f) Failure to ensure that the bottom surface of the wells is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before a plate is read.
 - g) Failure to read at the correct wavelength or use of an incorrect reference wavelength.

6. The use of highly haemolysed samples, incompletely clotted sera, plasma samples containing fibrin or samples with microbial contamination may give rise to erroneous results.
7. This test has not been evaluated for use with samples from cadavers.

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 DiaSorin Italia S.p.A. UK Branch
Central Road,
Dartford DA1 5LR
UK

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 **DiaSorin Italia S.p.A. UK Branch**
Central Road,
Dartford DA1 5LR
UK

0123

EC REP **DiaSorin Italia S.p.A.**
Via Crescentino snc
13040 Saluggia (VC) Italy

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DiaSorin Italia S.p.A.UK Branch,
Central Rd, Dartford DA1 5LR, UK
www.diasorin.com
Tel. +44 1322 317949

en

REF 9F80-01 / 05

GE34/36

Revised February, 2024

Murex HBsAg Version 3

Enzyme immunoassay for the detection of hepatitis B surface antigen in human serum or plasma.

The assay is intended to screen individual human donors for the presence of hepatitis B surface antigen or as an aid to the diagnosis of HBV infection.

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This instructions for use must be read carefully prior to use. The instructions for use must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions for use.

IVD

Key to symbols used	
 List Number	 <i>In Vitro</i> Diagnostic Medical Device
 Lot Number	 Store at 2-8°C
 Expiration Date	 CAUTION: Consult accompanying documents
 Manufacturer	 Consult instructions for use
	 Keep away from sunlight

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INTENDED PURPOSE

Murex HBsAg Version 3 is an enzyme immunoassay for the qualitative detection of hepatitis B surface antigen (HBsAg) in human serum or plasma samples. The assay is intended to screen individual human donors for the presence of HBsAg. It is also intended as an aid to the diagnosis of hepatitis B virus (HBV) infection in individuals with or without the symptoms of hepatitis. This assay is intended for *in vitro* diagnostic use, on either a semi-automated or automated platform.

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The HBV is a member of the *Hepadnaviridae*, a family of small, enveloped and primarily hepatotropic DNA viruses. Its small genome (3.2 kb) contains four open reading frames encoding 7 proteins: a secreted dimeric protein (HBV e antigen; HBeAg), a viral capsid protein (HBV core antigen; HBCAg), a reverse transcriptase (HBV Pol/RT), a transcriptional regulator (HBV x antigen, HBx) and, finally, three surface envelope glycoproteins (PreS1, PreS2 and hepatitis B surface antigen) originating from an overlapping open reading frame and representing large, medium and small surface envelope glycoproteins, respectively.¹

HBV replicates and assembles exclusively in hepatocytes with virion release through the cellular secretory pathway. Upon uptake into hepatocytes, its nucleocapsid is transported to the nucleus to release the relaxed-circular DNA (rcDNA) genome. The rcDNA is then converted into covalently closed circular DNA (cccDNA) which serves as a template for viral transcription and production of new rcDNA in a transcription/reverse transcription cycle. Hepatitis B surface (HBs) envelope proteins combine with mature nucleocapsids in the endoplasmic reticulum followed by secretion of the assembled virion into circulation. The production of these HBs proteins greatly exceeds the demand for virion formation. As a result, excessive HBs proteins bind each other and interact with lipids originating from the virion surface to form non-infectious filamentous and spherical subviral particles. These subviral particles are secreted into the systemic circulation and exceed the HBV virion concentrations by a factor of 10^2 to 10^5 . They are antigenic and, hence, are called hepatitis B surface antigen (HBsAg) particles. The HBs protein has the determinant antigenic component important for activating host immunity. HBsAg is found on both these particles and the circulating virions.²

The serum concentration of HBsAg depends on the amount of hepatitis B cccDNA genomes in hepatocytes, the transcriptional activity of these genomes and the host immune response against the virus.³ There are at least nine genotypes of HBV (A to I) on the criterion of at least 8% difference in genome sequences. Genotypes C, F and certain subtypes of genotype A are associated with higher rates of hepatocellular carcinoma. The efficacy of antiviral therapy and HBV vaccination is not affected by the genotype.⁵

There is universal consensus on the use of HBsAg as a marker in the initial assessment of hepatitis B, either for screening asymptomatic high-risk patients or diagnosing patients with symptoms of hepatitis B. HBsAg is the first serological marker to appear after an active infection⁷⁻¹⁰, even if HBV DNA is still undetectable. HBsAg represents the serological landmark of HBV infection and the following organizations recommend HBsAg for diagnosis or screening for HBV infection: the World Health Organization (WHO)⁵, European Association for the Study of the Liver (EASL)¹, Korean Association for the Study of the Liver (KASL)⁴, Asian-Pacific Association for the Study of the Liver (APASL)⁶.


PRINCIPLE OF THE PROCEDURE

In Murex HBsAg Version 3, the sample is pre-incubated in microwells coated with a mixture of mouse monoclonals specific for different epitopes on the 'a' determinant of HBsAg. Affinity purified goat antibody to HBsAg conjugated to horseradish peroxidase is then added to the sample in the well. During the two incubation steps any HBsAg present in the sample is bound to the well in an antibody-antigen-antibody-enzyme complex. In the absence of HBsAg no conjugate will be bound. After washing to remove sample and unbound Conjugate, a solution containing 3,3', 5,5'-tetramethylbenzidine (TMB) and hydrogen peroxide is added to the wells. Wells which contain HBsAg and hence bound Conjugate will develop a purple colour which is converted to orange when the enzyme reaction is terminated with sulphuric acid.

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1. Coated Wells

One plate (9F80-01) or five plates (9F80-05) of 96 wells coated with mouse monoclonal antibody to HBsAg (2.2 µg/mL).

Allow the wells to reach room temperature (18 to 30°C) before removal from the bag. place unused wells in the sealable storage bag provided and return to 2 to 8°C.

SAMPLE	DIL
--------	-----

2. Sample diluent

One bottle containing 16 mL of green/brown buffer containing detergents and proteins of goat and bovine origin. Mix by inversion before use. Contains 0.05% ProClin® 300 preservative.

CONTROL	-	⚠
---------	---	---

3. Negative Control

One bottle containing 2.5 mL of normal human serum. The serum is diluted in a buffer containing protein of bovine origin. Contains 0.05% Bronidox® preservative.

CONTROL	+	⚠
---------	---	---

4. Positive Control

One bottle containing 2 mL of inactivated human serum. The serum is diluted in a buffer containing protein of bovine origin. Contains 0.05% Bronidox® preservative.

CONJUGATE

5. Conjugate

One bottle containing 6 mL (9F80-01) or two bottles each containing 16 mL (9F80-05) of horseradish-peroxidase labelled goat antibody to HBsAg (0.01 mg/mL) in a red buffer containing proteins of bovine and goat origin. Mix by inversion before use. Contains 0.05% ProClin® 300 preservative.

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SUBSTRATE CONC

7. Substrate Concentrate

One bottle containing 35 mL of 3,3',5,5'- tetramethylbenzidine (TMB) and stabilisers in a pink solution.

Substrate Solution

To prepare the Substrate Solution add a volume of colourless Substrate Diluent to an equal volume of pink Substrate Concentrate in either a clean glass or plastic vessel. **It is important that this order of addition is followed and that any pipettes and glassware used to prepare Substrate Solution are clean.** Alternatively, the Substrate Solution may be made by pouring the entire contents of the bottle of Substrate Diluent into the bottle of Substrate Concentrate. One bottle of Substrate Solution provides sufficient reagent for at least five plates - see **Table 1**:

WASH FLUID

Additional reagent may be required for use with automated systems. Keep away from all natural and artificial light. The Substrate Solution should be pink; if it is purple before being used, it should be discarded and fresh Substrate Solution prepared.

The prepared Substrate Solution from this kit may be used interchangeably with that from all other Murex kits which use pink coloured Substrate Concentrate. Ensure that the Substrate Solution is prepared from Substrate Diluent and Substrate Concentrate provided together.

The prepared Substrate Solution is stable refrigerated (2 to 8°C) or at 15 to 25°C for up to two days but must be discarded if crystals have formed.

8. Wash Fluid

One bottle containing 125 mL of 20 times working strength glycine/Borate Wash Fluid. Contains 0.2% Bronidox® preservative.

Add one volume of Wash Fluid Concentrate to 19 volumes of distilled or deionised water to give the required volume or dilute the entire contents of one bottle of Wash Fluid to a final volume of 2500 mL. Crystals may be observed in the Wash Fluid Concentrate but these crystals will dissolve when the Wash Fluid is diluted to working strength. When diluted, the Wash Fluid contains 0.01% Bronidox® preservative.

The Wash Fluid from this kit may be used interchangeably with the glycine/Borate Wash Fluid from any other Murex kit. Store the working strength Wash Fluid at 18 to 30°C in a closed vessel under which conditions it will retain activity for one month.

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**Table 1
Volume of Substrate Concentrate and
Substrate diluent Required**

Number of Wells	No of plates
8 16 24 32 40 48 56 64 72 80 96	1 2 3 4
Substrate Concentrate (mL)	
0.5 1.0 2.0 2.5 2.5 3.0 3.5 4.0 4.5 4.5 6.0	6 12 18 22
Substrate Diluent (mL)	
0.5 1.0 2.0 2.5 2.5 3.0 3.5 4.0 4.5 4.5 6.0	6 12 18 22

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Table 2

Component	Reactive for	Non-reactive for
Negative Control	N/A	HBsAg, and antibodies to HIV-1 and 2, HCV and HTLV I + II
Positive Control	HBsAg	Antibodies to HIV-1 and HIV-2, and HCV

All reactive serum used has been inactivated prior to use in reagent preparation. However, all material of human origin should be considered as potentially infectious and it is recommended that this kit and test specimens be handled using established good laboratory practice.

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Reagents:	<input type="text" value="SAMPLE"/> <input type="text" value="DIL"/> · <input type="text" value="CONJUGATE"/>
Classification:	Skin sens. 1 H317
Signal Word:	Warning
Symbols / Pictograms:	
Hazard Statements:	H317 May cause an allergic skin reaction.
Precautionary Statements:	P280 Wear protective gloves/protective clothing/eye protection/face protection. P363 Wash contaminated clothing before reuse. P333+P313 If skin irritation or rash occurs: get medical advice/attention.
Contains:	reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1).
Reagents:	<input type="text" value="SUBSTRATE"/> <input type="text" value="CONC"/>
Classification:	Eye Irrit. 2 H319
Signal Word:	Warning
Symbols / Pictograms:	
Hazard Statements:	H319 Causes serious eye irritation
Precautionary Statements:	P264 Wash hands thoroughly after handling P280 Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338 IF IN EYES: rinse cautiously with water for several minutes. remove contact lenses, if present and easy to do. Continue rinsing.

Pursuant to EC regulation 1272/2008 (CLP), is labeled as EUH210, safety data sheets available on request.

For additional information see Safety Data Sheets available on www.diasorin.com.

- Potentially contaminated materials should be disposed of safely according to local requirement.
- Spillage of potentially infectious material should be removed immediately with absorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued¹⁰. Sodium hypochlorite should not be used on acid containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- Neutralised acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.
- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- The following reagents contain low concentrations of harmful substances.
 - The Conjugate and Sample Diluent contain detergents.
- Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If either come into contact with the skin or eyes, wash thoroughly with water.
- If any of the reagents come into contact with the skin or eyes wash the area extensively with water.

ANALYTICAL PRECAUTIONS

- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the Test procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
- Allow all reagents and samples to come to 18 to 30°C before use. Immediately after use return all reagents to the recommended storage temperature.
- Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionised water.
- Avoid the use of self-defrosting freezers for the storage of reagents and samples.
- Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
- Do not allow wells to become dry during the assay procedure.
- Do not cross-contaminate reagents. Dedicate a pipette for use with the Substrate Solution of Murex assays. A pipette should also be dedicated for use with the Conjugate.
- Do not touch or splash the rim of the well with Conjugate. Do not blow out from micropipettes; reverse pipetting is recommended wherever possible.
- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate.
- Do not contaminate microwells with the dust from disposable gloves.

12. When using fully automated microplate processors:
 - i) It is not necessary to use plate lids and to tap dry the wells.
 - ii) Do not allow system fluids from fully automated microplate processors to contaminate samples or reagents.
 - iii) The possibility of cross contamination between assays needs to be excluded when validating assays on fully automated processors.
13. Ensure the assay is run within the temperature limits defined in the assay protocol.
14. Do not use CO₂ incubators.
15. Do not store the Stop Solution in a shallow dish or return it to a stock bottle after use.
16. The possibility of cross contamination between assays needs to be excluded when validating assay protocols on instrumentation.

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SPECIMEN COLLECTION

Serum, EDTA plasma or citrate plasma samples may be used. Blood collected by venepuncture should be allowed to clot naturally. Ensure that the serum samples are fully clotted. Remove any visible particulate matter from the sample by centrifugation. If samples are prepared using liquid anti-coagulants e.g. citrate plasma, the dilution effect should be considered.

SPECIMEN TRANSPORT AND STORAGE

Store samples at 2 to 8°C. Samples not required for assay within 72 hours should be removed from the clot or cell pellet and stored frozen (-15°C or colder). Samples can be freeze-thawed up to four times. Samples may be stored frozen for up to 12 months.¹¹ After thawing, ensure samples are thoroughly mixed before testing.

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1. **Stop Solution (0.5M to 2M Sulphuric Acid)**. E.g. add between 3 mL (for 0.5M) and 11 mL (for 2.0M) of analytical grade concentrated sulphuric acid (18.0M) to about 80 mL of distilled or deionised water and then make up to 100 mL with more water. Alternatively, the following reagent can be used; 1N Sulphuric Acid (Code N0164, 15 vial pack and N0165, 1 vial pack).
2. **Freshly distilled or high quality deionised water** is required for dilution of Wash Fluid, for preparation of the Stop Solution and for use in conjunction with automated washers.
3. **Micropipettes and Multichannel micropipettes** of appropriate volume.
4. **Incubator** capable of maintaining the temperature limits defined in the assay protocol.
5. **Moulded Heating Block** (Code 5F09-02). For use in laboratory incubators. The moulded heating block should ideally be kept in the incubator used. If this is not possible it must be placed in the incubator at least four hours before beginning the assay.
6. **Instrumentation**
 - a) Automated microplate strip washer.
 - b) Microplate reader.

or

 - c) Fully automated microplate processor.

All instruments must be validated before use. Please contact your representative for details of recommended systems, software protocols for instrumentation and validation procedures.
7. **Disposable Reagent Troughs** (Code 5F24-01).
8. **Sodium hypochlorite** for decontamination (refer to **Health and Safety Information**).

9. **Sodium hydroxide solution (0.1M)** (for instrument decontamination).

TEST PROCEDURE

Please read **Analytical Precautions** carefully before performing the test.

Addition of the various components of the assay to the wells may be confirmed visually by examining the plate for the following colours: **Sample diluent** is green/brown in colour. On addition of the Sample or Control the colour will change to blue/green. The colour change will vary from sample to sample but some change should always be visible.

Conjugate is red in colour.

Substrate Solution is initially pink with any reactive wells becoming purple. On addition of Stop Solution the purple colour of the reactives will change to orange, whilst the negatives remain pink.

The addition of sample or reagent can be confirmed using a microplate reader as follows: Sample Diluent plus Sample read at 570 or 620 nm with a reference at 690 nm, Conjugate at 490 nm with a reference at 690 nm, Substrate Solution at 490 nm (no reference).

SEMI AUTOMATED PROCESSING

Step 1	Prepare Substrate Solution and Wash Fluid .	
Step 2	Use only the number of wells required for the test.	
Step 3	Add 25 µl of Sample diluent to each well.	25 µL
Step 4	Add 75 µl of Samples or Controls to the wells. To each plate add 75 µl of the Negative Control to wells A1 and B1 and 75 µl of positive Control into well C1. Add the Controls to the designated wells after dispensing the samples.	75 µL
Step 5	Cover the plate with a lid and incubate for 60 minutes at 37°C ±1°C.	60 min
Step 6	Add 50 µl of Conjugate to each well.	50 µL
Step 7	Shake the plate using a plate shaker for 10 seconds or manually agitate by gently tapping the sides for 10 seconds.	10 sec
Step 8	Cover the plate with the lid and incubate for 30 minutes at 37°C ±1°C.	30 min
Step 9	At the end of the incubation time wash the plate 5 times as described under Wash procedures . After washing is completed invert the plate and tap out any residual Wash Fluid onto absorbent paper.	5 washes
Step 10	Immediately after washing the plate, add 100 µl of Substrate Solution to each well.	100 µL
Step 11	Cover the plate with a lid and incubate for 30 minutes at 37°C ±1°C while colour develops. A purple colour should develop in wells containing reactive samples.	30 min
Step 12	Add 50 µl of Stop Solution to each well.	50 µL
Step 13	Within 15 minutes read the absorbance of each well at 450 nm using 620 nm to 690 nm as the reference wavelength if available. Blank the instrument on air (no plate in the carriage).	A_{450/Ref}

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- a. **Protocol for automated microplate stripwasher**
Perform 5 wash cycles using working strength Wash Fluid. Ensure, where possible, that:
- Flow-through washing with a fill volume of 500 µl/well is used with instrumentation supplied by DiaSorin. When using other instrumentation for which this is not possible, ensure that the well is completely filled.
 - The dispense height is set to completely fill the well with a slight positive meniscus, without causing an overflow.
 - The time taken to complete one aspirate/wash/soak cycle is approximately 30 seconds.
 - Ensure that no liquid is left in the well (by use of a double aspirate step in the final cycle where possible).
 - After washing is completed, invert the plate and tap out any residual Wash Fluid onto absorbent paper.

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Washers must be rinsed with distilled water at the end of the test to avoid blockage and corrosion.

FULLY AUTOMATED MICROPLATE PROCESSORS

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- Do not programme times shorter than specified in the procedure.
- For each incubation at 37°C, programmed times may be increased by up to 20%.
- Wells containing Sample Diluent may be left for up to 60 minutes at 18 to 30°C prior to addition of sample and controls and for up to 60 minutes after the addition of sample and controls before starting **Step 5**.
- Ensure all **Analytical Precautions** are followed. Protocols written following these guidelines must be fully validated prior to use according to local procedures.

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Negative Control

Calculate the mean absorbance of the replicates of the Negative Control. If one of the Negative Control wells has an absorbance more than 0.03 above the other discard the higher value.

Cut-off Value

Calculate the *Cut-off Value* by adding 0.05 to the mean of the Negative Control replicates.

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Negative Control

The mean $A_{450/ref}$ of the Negative Control is less than 0.15 or the mean A_{450} of the Negative Control is less than 0.2.

Positive Control

The $A_{450/ref}$ or A_{450} of the Positive Control is more than 0.8 above the mean $A_{450/ref}$ or A_{450} of the Negative Control.

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INTERPRETATION OF RESULTS

Non Reactive Result

Samples giving an absorbance less than the Cut-off Value are considered non-reactive in Murex HBsAg Version 3.

Reactive Result

Samples giving an absorbance equal to or greater than the Cut-off Value are considered initially reactive in the assay (see **Limitations of the Procedure**). Such samples should be retested in duplicate using the original sample source. Samples that are reactive in at least one of the re-tests are presumed to contain HBsAg and should be confirmed by using the Murex HBsAg Confirmatory Version 3 kit (2G27-01) and tests for other HBV markers. Samples that are non-reactive in both wells on retest should be considered non-reactive.

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The performance of Murex HBsAg Version 3 has been determined by testing samples from random blood donors, patients with acute and chronic hepatitis B infection, patients with mutant forms of hepatitis B infection and patients with diseases unrelated to hepatitis B infection.

In addition, its performance against the French A.F.S.S.A.P.S. panels and other commercially available seroconversion samples has been evaluated.

1. Diagnostic specificity

The Murex HBsAg assay demonstrated a specificity of $\geq 99.5\%$ in a study where a total of 12330 routine donor samples were screened with Murex HBsAg Version 3. In the study, 0.18% (22/12330) of samples were initially reactive and 0.03% (4/12330) were repeatedly reactive. None of the repeatedly reactive samples with Murex HBsAg version 3 and the alternative assays were confirmed as positive for the presence of hepatitis B surface antigen.

The specificity of Murex HBsAg Version 3 on presumed negative samples from donors is estimated to be 99.97% (12326/12330) with a 95% confidence interval of 99.92% (12320/12330) to 99.99% (12329/12330) by the binomial distribution.

2. Analytical specificity

Cross-reactivity

In addition, 798 potentially cross-reactive samples from patients with conditions unrelated to hepatitis B infection, including other acute viral infections, antenatal and autoimmune samples, were tested with Murex HBsAg Version 3. A total of 796 of these samples were non-reactive with Murex HBsAg Version 3. Of the two repeatedly reactive samples, one was false reactive and showed no other hepatitis markers, the remaining sample was anti-HBc positive.

Table 3
Performance of Murex HBsAg Version 3 with potentially cross-reacting samples

Clinical category	Total
Acute viral infections (hospitalised)	233 (1)
Autoimmune disease (hospitalised)	108
Clinically normal	27
Antenatal	417 (1)
Rheumatoid factor positive	10
Unknown	3
Total	798 (2)

(#) Indicates cross reactive samples

Interference

Controlled studies with potentially interfering substances showed no interference at the concentrations for each substance summarised in Table 4 for the Murex HBsAg version 3.

Table 4
Performance of Murex HBsAg Version 3 with potentially interfering samples

Substance	Result
Indirect bilirubin	No interference at 0.2mg/mL
Triglycerides	No interference at 30mg/mL
Heparin	No interference at 40 USP/mL
Cholesterol	No interference at 5mg/mL
Total protein	No interference between 40-120 mg/mL
Red blood cells	No interference
Interferon alpha 1	No interference at 6000 IU/mL
Interferon alpha 2a	No interference at 6000 IU/mL
Interferon alpha 2b	No interference at 6000 IU/mL
Acetaminophen	No interference at 207 mg/L
Ibuprofen	No interference at 512 mg/L
Adefovir dipivoxil	No interference at 10 mg/L
Telbivudine	No interference at 600 mg/L
Lamivudine	No interference at 1.05 mg/dL
Entecavir	No interference at 0.5 mg/L

NB*100 mg/dL=1 mg/mL

3. Diagnostic sensitivity

Samples from patients at various stages of hepatitis B infection and patients with conditions unrelated to hepatitis B were tested in three regional virus reference laboratories and at DiaSorin.

A total of 584 samples from patients suffering from acute and chronic hepatitis B infection were tested with Murex HBsAg version 3. All 584 samples were confirmed with an alternative immunoassay for HBsAg and found to be positive in both assays.

A panel of ten (10) recombinant mutants were tested with the Murex HBsAg version 3 and three comparator CE-marked assays to determine correct antigenic detection of the HBsAg structure. The mutants contained important epitope clusters within amino acids 100-160, that includes the "a determinant region" (amino acids 124-147), the most important target for serological diagnosis. Eight of ten recombinant mutants were recognised with Murex HBsAg, comparable with, or better than, the competitor assays.

Table 5
Performance Murex HBsAg Version 3 with recombinant mutants

Sample ID	Mutation	Comparator HBsAg assay ¹²			Murex, data from Poster 1 ¹²
		1	2	3	
Mutant-01	T123N	+	+	+	+
Mutant-02	T123N-T124S	-	-	-	-
Mutant-03	P142L-F/Y143H-D144E-G145R	+	-	-	-
Mutant-04	I110R-SS117I-G119R-T123N	+	+	-	+
Mutant-05	I22+DT	+	+	-	+
Mutant-06	I22+DT-G145R	-	+	-	+
Mutant-07	G145R	+	+	+	+
Mutant-08	D114A	+	+	+	+
Mutant-09	P142L-G145R	+	+	+	+
Mutant-10	P142S-G145R	+	+	+	+

Note: '+' and '-' indicate reactive and non-reactive respectively

A further five samples from patients infected with mutant forms of hepatitis B infection, confirmed by DNA sequencing, were also tested with Murex HBsAg version 3 and were all detected successfully.

Additionally, a total of 33 commercially available HBV seroconversion panels were tested with Murex HBsAg Version 3. Comparison with four other commercially available microplate based immunoassays for the detection of hepatitis B surface antigen showed that Murex HBsAg Version 3 detected HBsAg comparably to the other commercial available immunoassays.

4. Analytical sensitivity

Murex HBsAg Version 3 was tested with the WHO Third International Standard for HBsAg, (HBV genotype B4, subtypes ayw1/adw2), NIBSC code 12/226, and it was determined that it was detectable to a limit of 0.023IU/mL.

5. Assay Repeatability and Reproducibility

Ten replicates of each of five samples were tested on ten separate test occasions with two separate batches to ascertain the reproducibility of Murex HBsAg Version 3. The results of the study are summarised in Table 6 and Table 7.

Table 6
Murex HBsAg Version 3 - Assay Repeatability and Reproducibility (Kit 1)

Specimen	Number of assay	Mean S/CO	Intra assay %CV	Inter assay %CV
1	10	9.26	7.6	10.0
2	10	3.56	8.2	10.8
3	10	12.28	10.0	13.1
4	10	1.88	9.3	12.4
5	10	0.69	11.0	14.2

Table 7
Murex HBsAg Version 3 - Assay Repeatability and Reproducibility (Kit 2)

Specimen	Number of assay	Mean S/CO	Intra assay %CV	Inter assay %CV
1	10	10.46	6.8	8.1
2	10	3.89	7.2	9.2
3	10	12.22	7.7	9.5
4	10	1.80	9.3	13.1
5	10	0.58	7.4	7.9

Representative performance data are shown: results obtained at individual laboratories and with different populations may vary.

LIMITATIONS OF THE PROCEDURE

1. The **Test Procedure** and **Interpretation of Results** must be followed.
2. This test has only been evaluated for use with individual (unpooled) serum, EDTA plasma, or citrate plasma samples.
3. A negative result with an antigen detection test does not preclude the possibility of infection.
4. Non-repeatable reactive results may be obtained with any EIA procedure.
5. The most common sources of error are:
 - a) Imprecise delivery of Sample, Conjugate or Substrate into the wells.
 - b) Contamination of Substrate with Conjugate.
 - c) Contamination with conjugates from other assays
 - d) Blocked or partially blocked washer probes.
 - e) Insufficient aspiration leaving a small volume of Wash Fluid in the wells.
 - f) Failure to ensure that the bottom surface of the wells is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before a plate is read.
 - g) Failure to read at the correct wavelength or use of an incorrect reference wavelength.
6. The use of highly haemolysed samples, **samples containing human anti-mouse antibody**, incompletely clotted sera, plasma samples containing fibrin or samples with microbial contamination may give rise to erroneous results.
7. This test has not been evaluated for use with samples from cadavers.

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BIBLIOGRAPHY

1. EASL. Clinical Practice Guidelines on the management of hepatitis B virus infection. s.l.: **European Association for the Study of the Liver**. Electronic address, e.e.e.& European Association for the Study of the Liver, L. EASL, 2017. J Hepatol 67, 370-398. (2017).
2. Advances in the molecular diagnosis of hepatitis B infection: providing insight into the next generation of disease. **Bayliss, J. Nguyen, T., Lesmana, C.R., Bowden, S. & Revill, P.** 113-121, 2013, Vol. Semin Liver Dis. 33.
3. Management of chronic hepatitis B: Canadian Association for the Study of the Liver consensus guidelines. **Coffin, C.S., Fung, S.K. & Ma, M.M.** 917-938, s.l.: Can J. Gastroenterol, 2012, Vol. 26.
4. KASL, Clinical practice guidelines: management of chronic hepatitis B. (**KASL**), **Korea Association of the Study of the Liver**. 18-75, s.l.: Clin Mol Hepatol, 2016, Vol. 22.
5. Guidelines for the prevention, care and treatment of persons with chronic hepatitis B infection. **Organization, World Health**. s.l.: WHO, 2015.
6. Asian-Pacific clinical practice guidelines on the management of hepatitis B: a 2015 update. **Sarin, S.K. et al.** 1-98, s.l.: Hepatol. Int., 2016, Vol. 10.
7. SASLT practice guidelines for the management of hepatitis B virus. **Abaalkhail, F. et al.** 5-25, s.l.: Saudi J Gastroenterol, 2014, Vol. 20.
8. Updated CDC recommendations for the management of hepatitis B virus-infected health-care providers and students. (CDC), **Centers for Disease Control and Prevention**. 1-12, s.l.: MMWR Recomm Rep, 2012, Vol. 61.
9. South African guideline for the management of chronic hepatitis B. **Spearman, C.W. et al.** 337-349, s.l.: S Afr Med J, 2013, Vol. 103.
10. **Centres for disease Control**. (1985). recommendations for preventing transmission of infection with human T-lymphotropic virus type III / lymphadenopathy-associated virus in the workplace. MMWR, 34, No. 45, 681.
11. **Fung, J. et al.** (2011). Stability of hepatitis B surface antigen over time: Implications for studies using stored sera. Journal of Medical Virology, 83(11), 1900–1904.
12. Sarasso C. et al. A Novel Panel of Recombinant Hepatitis B Surface Antigens Mutants. Source: https://www.diasorin.com/sites/default/files/rd/publications/attachments/poster_lsn_xl_hbsag_q_escv_2010_sarasso.pdf.

Bronidox® and ProClin® are not trademarks of DiaSorin.

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LITERATŪRA

- EASL. Clinical Practice Guidelines on the management of hepatitis B virus infection. s.l.: **European Association for the Study of the Liver**. Electronic address, e.e.e.& European Association for the Study of the Liver, L. EASL, 2017. J Hepatol 67, 370-398. (2017).
- Advances in the molecular diagnosis of hepatitis B infection: providing insight into the next generation of disease. **Bayliss, J. Nguyen, T., Lesmana, C.R., Bowden, S. & Revill, P.** 113-121, 2013, Vol. Semin Liver Dis. 33.
- Management of chronic hepatitis B: Canadian Association for the Study of the Liver consensus guidelines. **Coffin, C.S., Fung, S.K. & Ma, M.M.** 917-938, s.l.: Can J. Gastroenterol, 2012, Vol. 26.
- KASL, Clinical practice guidelines: management of chronic hepatitis B. (**KASL**), **Korea Association of the Study of the Liver**. 18-75, s.l.: Clin Mol Hepatol, 2016, Vol. 22.
- Guidelines for the prevention, care and treatment of persons with chronic hepatitis B infection. **Organization, World Health**. s.l.: WHO, 2015.
- Asian-Pacific clinical practice guidelines on the management of hepatitis B: a 2015 update. **Sarin, S.K. et al.** 1-98, s.l.: Hepatol. Int., 2016, Vol. 10.
- SASLT practice guidelines for the management of hepatitis B virus. **Abalkhail, F. et al.** 5-25, s.l.: Saudi J Gastroenterol, 2014, Vol. 20.
- Updated CDC recommendations for the management of hepatitis B virus-infected health-care providers and students. (CDC), **Centers for Disease Control and Prevention**. 1-12, s.l.: MMWR Recomm Rep, 2012, Vol. 61.
- South African guideline for the management of chronic hepatitis B. **Spearman, C.W. et al.** 337-349, s.l.: S Afr Med J, 2013, Vol. 103.
- Centres for disease Control**. (1985). recommendations for preventing transmission of infection with human T-lymphotropic virus type III / lymphadenopathy-associated virus in the workplace. MMWR, 34, No. 45, 681.
- Fung, J. et al.** (2011). Stability of hepatitis B surface antigen over time: Implications for studies using stored sera. Journal of Medical Virology, 83(11), 1900–1904.
- Sarasso C. et al. A Novel Panel of Recombinant Hepatitis B Surface Antigens Mutants. Source: https://www.diasorin.com/sites/default/files/rd/publications/attachments/poster_Isn_xl_hbsag_q_escv_2_010_sarasso.pdf.

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0318

HBc IgM

**“Capture” Enzyme ImmunoAssay (ELISA)
for the quantitative/qualitative
determination of IgM class antibody to
Hepatitis B Virus core Antigen
in human plasma and sera**

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HBc IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgM class antibodies to Hepatitis B Virus core Antigen in human plasma and sera with the "capture" system.

The kit is intended for the classification of the viral agent and for the follow-up of chronic patients under therapy. For "in vitro" diagnostic use only.

B. INTRODUCTION

Hepatitis B core Antigen (or HBcAg) is the major component of the core particles of Hepatitis B virus (or HBV).

Particles have a size of 27nm and contain a circular double-stranded DNA molecule, a specific DNA-polymerase and HBcAg. HBcAg is composed of a single polypeptide of about 17 kD that is released upon disaggregation of the core particles ; the antigen contains at least one immunological determinant.

Upon primary infection, anti HBcAg IgM antibodies are one of the first markers of HBV hepatitis appearing in the serum of the patient, together or slightly later than HBsAg, the viral surface antigen.

Anti HBcAg IgM titers, very high during the acute phase, decrease along the illness, as IgG antibodies appear, down to undetectable levels in convalescent patients.

In chronic hepatitis, however, spikes of anti HBcAg IgM synthesis are present, confirming reactivation of HBV in hepatocytes and giving origin to permanent IgM low titers.

The determination of anti HBcAg IgM antibodies has become very important for the fast classification of the virus, of the phase of the illness and for the monitoring of patients under treatment with interferon.

C. PRINCIPLE OF THE TEST

The assay is based on the principle of "IgM capture" where IgM class antibodies in the sample are first captured by the solid phase coated with anti hIgM antibody.

After washing out all the other components of the sample and in particular IgG antibodies, the specific IgM captured on the solid phase are detected by the addition of a purified preparation of recombinant HBcAg, labelled with a monoclonal antibody conjugated with peroxidase (HRP).

After incubation, microwells are washed to remove unbound conjugate and then the chromogen/substrate is added.

In the presence of peroxidase the colourless substrate is hydrolysed to a coloured end-product, whose optical density may be detected and is proportional to the amount of IgM antibodies to HBcAg present in the sample.

D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

1. Microplate: MICROPLATE

8x12 microwell strips coated with purified anti human IgM specific mouse monoclonal antibody, post-coated with bovine serum proteins and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Curve: CAL N° ...

6x2.0 ml/vial. Ready to use and color coded standard curve calibrated on the HBcIgM reference preparation supplied by Paul Ehrlich Institute (HBc-Referenzserum-IgM 84), ranging: CAL1 = 0 U/ml // CAL2 = 5 U/ml // CAL3 = 10 U/ml // CAL4 = 20 U/ml // CAL 5 = 50 U/ml // CAL 6 = 100 U/ml.

It contains chemical inactivated HBcIgM positive human plasma, 100 mM Tris buffer pH 7.4+/-0.1, 0.5% Tween 20, 0.09% sodium azide and 0.045% ProClin 300 as preservatives. The Calibration Curve is coded with blue alimentary dye.

Important Note: Even if plasma has been chemically inactivated, handle this component as potentially infectious.

3. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle. 20x concentrated solution.

Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

4. Enzyme Conjugate (Immunocomplex) : CONJ

1x16.0 ml/vial. Ready-to-use solution. Contains an immunocomplex formed by a specific mouse monoclonal antibody, labelled with HRP, and a purified recombinant HBcAg. The reagent is dissolved into a buffer solution 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The component is red colour coded.

5. Specimen Diluent : DILSPE

2x60.0 ml/vial. Buffered solution for the dilution of samples; it contains 100 mM Tris buffer pH 7.4+/-0.1, 0.5% Tween 20, 2% Casein, 0.045% ProClin 300 and 0.09% sodium azide as preservatives. The component is blue color coded.

6. Control Serum : CONTROL ...ml

1 vial. Lyophilized. Contains fetal bovine serum, human HBcIgM positive human plasma calibrated at 20 ± 10% PEI U/ml. 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Important Notes

1. The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label .

2. Important Note: Even if plasma has been chemically inactivated, handle this component as potentially infectious.

7. Chromogen/Substrate : SUBS TMB

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂.

Note: To be stored protected from light as sensitive to strong illumination.

8. Sulphuric Acid: H₂SO₄ 0.3 M

1x15ml/vial. Contains 0.3 M H₂SO₄ solution.

Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

9. Plate sealing foils: n° 2

10. Package insert: n° 1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (150ul, 100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
14. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND RECOMMANDATIONS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been

observed in the preparation of the sample with citrate, EDTA and heparin.

2. Avoid any addition of preservatives; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results.
4. Haemolysed and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing. In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2°-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Calibration Curve:

Ready to use. Mix well on vortex before use.

Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bidistilled water up to 1200ml and mixed gently end-over-end before use.

During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: Once diluted, the wash solution is stable for 1 week at +2..8° C.

Enzyme conjugate:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

Specimen Diluent

Ready to use. Mix on vortex before use.

Control Serum

Dissolve the content of the vial with EIA grade water as reported in the label. Mix well on vortex before use. The dissolved control serum is ready to use.

Note: The control after dissolution is not stable. Store frozen in aliquots at -20°C.

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces.
If this component has to be transferred use only plastic, and if possible, sterile disposable container

Sulphuric Acid:

Ready to use. Mix well on vortex before use.
Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.
H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.
P332 + P313 – If skin irritation occurs: Get medical advice/attention.
P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313 – If eye irritation persists: Get medical advice/attention.
P362 + P363 – Take off contaminated clothing and wash it before reuse.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.
An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Blanking is carried out on the well identified in the section "Assay Procedure". The optical

system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate (TMB+H₂O₂) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Control Serum as described above and gently mix.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.

In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

Two procedures can be carried out with the device according to the request of the clinician.

M.1 Quantitative analysis

1. Place the required number of strips in the plastic holder and carefully identify the wells for standards and samples.

- Dilute samples **1:101** dispensing 1 ml Sample Diluent into a disposable tube and then 10 µl sample; mix on vortex before use. Do not dilute the Calibrators and the dissolved Control Serum as they are ready-to-use.
- Leave the A1+B1 wells empty for blanking purposes.
- Pipette 100 µl of the Calibrators in duplicate, 100 µl dissolved Control Serum in duplicate followed by 100 µl of diluted samples. The Control Serum is used to verify that the whole analytical system works as expected. Check that Calibrators, Control Serum and samples have been correctly added.
- Incubate the microplate **for 60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section I.3)
- In all the wells except A1+B1, pipette 100 µl Enzyme Conjugate. Incubate the microplate **for 60 min at +37°C**.

Important note: Be careful not to touch the inner surface of the well with the pipette tip and not to immerse the top of it into samples or controls. Contamination might occur.

- When the second incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1+B1 included.

Important note: Do not expose to strong direct light. as a high background might be generated.

- Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with positive samples, the control serum and the positive calibrators, as well, will turn from clear to blue.
- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9 to block the enzymatic reaction.. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

M.2 Qualitative analysis

- Place the required number of strips in the plastic holder and carefully identify the wells for standards and samples.
- Dilute samples **1:101** dispensing 1 ml Sample Diluent into a disposable tube and then 10 µl sample; mix on vortex before use. Do not dilute the Calibrators as they are ready-to-use.
- Leave the A1 well empty for blanking purposes.
- Pipette 100 µl Calibrator 0 U/ml in duplicate, 100 µl Calibrator 10 U/ml in duplicate and 100 µl Calibrator 100 U/ml in single. Then dispense 100 µl diluted samples in proper sample wells. Check that Calibrators and samples have been correctly added.
- Incubate the microplate **for 60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section I.3)
- In all the wells except A1, pipette 100 µl Enzyme Conjugate. Incubate the microplate **for 60 min at +37°C**.

Important note: Be careful not to touch the inner surface of the well with the pipette tip and not to immerse the top of it into samples or controls. Contamination might occur.

- When the second incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1 included.

Important note: Do not expose to strong direct light. as a high background might be generated.

- Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with positive samples, the control serum and the positive calibrators, as well, will turn from clear to blue.
- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9 to block the enzymatic reaction. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

Important notes:

- Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
- Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.
- The Control Serum (CS) does not affect the cut-off calculation and therefore the test results calculation. The Control Serum may be used only when a laboratory internal quality control is required by the management

N. ASSAY SCHEME

The assay protocol can be summarized in the table below:

Calibrators & diluted samples & dissolved Control Serum	100 ul
1st incubation	60 min
Temperature	+37°C
Washing steps	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing steps	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Chromogen/Substrate	100ul
3rd incubation	20 min
Temperature	room
Sulphuric Acid	100 ul
Reading OD	450nm /620-630nm

An example of dispensation scheme in quantitative assays is reported below:

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Legenda: BLK = Blank // CAL = Calibrators
CS = Control Serum // S = Sample

An example of dispensation scheme in qualitative assays is reported below:

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S 3	S 11										
B	CAL1	S 4	S 12										
C	CAL1	S 5	S 13										
D	CAL3	S 6	S 14										
E	CAL3	S 7	S 15										
F	CAL6	S 8	S 16										
G	S 1	S 9	S 17										
H	S 2	S 10	S 18										

Legenda: BLK = Blank // CAL = Calibrators// S = Sample

	with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator 5 U/ml < CAL 0 + 5SD or < CAL 0 + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Calibrator 10 U/ml < CAL 0 + 0.200	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Calibrator 100 U/ml < 1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the calibrator; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Control Serum Different from expected value	First verify that: 1. the procedure has been correctly performed; 2. no mistake has occurred during its distribution (ex.: dispensation of a wrong sample); 3. the washing procedure and the washer settings are correct; 4. no external contamination of the standard has occurred. 5. the Control Serum has been dissolved with the right volume reported on the label. If a mistake has been pointed out, the assay has to be repeated after eliminating the reason of this error. If no mistake has been found, proceed as follows: a) a value up to +/-20% is obtained: the overall Precision of the laboratory might not enable the test to match the expected value +/-10%. Report the problem to the Supervisor for acceptance or refusal of this result. b) a value higher than +/-20% is obtained: in this case the test is invalid and the DiaPro's customer service has to be called.

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as qualified.

Control that the following data are matched:

Parameter	Requirements
Blank well	< 0.100 OD450nm
Calibrator 0 PEI U/ml coefficient of variation	< 0.150 OD450nm after blanking < 30%
Calibrator 0 PEI U/ml	OD450nm > OD450nm Cal 0 U/ml + 5SD and anyway > OD450nm Cal 0 U/ml + 0.100
Calibrator 10 PEI U/ml	OD450nm > OD450nm Cal 0 U/ml + 0.200
Calibrator 100 PEI U/ml	> 1.000 OD450nm
Control Serum	OD450nm = OD450nm of the Calibrator 20 U/ml ± 10%

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Calibrator 0 U/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive calibrators instead of Cal 0); 4. that no contamination of the Cal 0, or of the wells where this was dispensed, has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated

If any of the above problems have occurred, report the problem to the supervisor for further actions.

Important note:

The analysis must be done proceeding as the reading step described in the section M, point 12.


P. RESULTS

P.1 Quantitative method

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm/620-630nm (4-parameters interpolation is suggested).

Then on the calibration curve calculate the concentration of anti HBc IgM antibody in samples.

An example of Calibration curve is reported below.

Important Note: Do not use this example to make real calculations on samples.

P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm/620-630nm values for the Calibrators 0 and 10 U/ml and then check that the assay is valid.

Example of calculation (data obtained proceeding as the the reading step described in the section M, point 12).

The following data must not be used instead of real figures obtained by the user.

Calibrator 0 U/ml: 0.020 – 0.024 OD450nm
 Mean Value: 0.022 OD450nm
 Lower than 0.150 – Accepted
 Calibrator 10 U/ml: 0.350 – 0.330 OD450nm
 Mean Value: 0.340 OD450nm
 Higher than Cal 0 + 0.200 – Accepted
 Calibrator 100 U/ml: 2.845 OD450nm
 Higher than 1.000 – Accepted

Q. INTERPRETATION OF RESULTS

Q.1 Qualitative results

For qualitative interpretations, the medical literature generally considers positive samples showing a concentration of HBc IgM ≥ 10 PEI U/ml.

Test results are therefore interpreted as a ratio of the sample OD450nm and the OD450nm/620-630nm of the Cal 10 PEI U/ml (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

Q.2 Quantitative results

The calibration curve is used to determine the concentration of IgM antibodies to HBcAg in samples.

Samples with a concentration lower than 5 PEI U/ml are considered negative for HBcIgM.

Samples with a concentration between 5 and 10 PEI U/ml are considered in a gray-zone.

In the follow up of chronic hepatitis, however, values higher of 5 PEI U/ml may be considered positive for HBcIgM, when in presence of other clinical signs.

Samples with a concentration higher than 10 PEI U/ml are considered positive for HBcIgM.

Important general notes:

- When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to produce the calibration curve, calculate sample concentration and generate the correct interpretation of results.
- Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
- A positive result is indicative of HBV infection and therefore the patient should be treated accordingly.
- When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
- Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

1. Limit of detection

The limit of detection of the assay has been calculated by means of :

- the HBcIgM reference preparation supplied by Paul Erlich Institute, Germany (HBc-Referenzserum-IgM 84), on which the Standard Curve has been calibrated.
- Accurun 113 (cat. N° A113-5001) supplied by Boston Biomedica Inc., USA

Results of Quality Control for three lots are given in the following tables:

BCM.CE	Lot #	0103	Lot #	0103/2	Lot #	0303
PEI U/ml	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
100	2.752	8.9	2.883	9.7	2.911	9.1
50	1.917	6.2	1.972	6.7	2.053	6.4
20	0.980	3.2	0.914	3.1	1.095	3.4
10	0.544	1.8	0.513	1.7	0.592	1.8
5	0.310	1.0	0.296	1.0	0.321	1.0
2.5	0.155	0.5	0.149	0.5	0.161	0.5
1.25	0.084	0.3	0.084	0.3	0.093	0.3
negative	0.040		0.035		0.044	

BBI Accurun # 113 lot # 48-9999-0621

BCM.CE	Lot #	0103	Lot #	0103/2	Lot #	0303
BBI 113	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
1 x	3.336	10.8	3.195	10.4	3.269	10.3
2 x	2.472	8.0	2.385	7.8	2.385	7.5
4 x	1.467	4.7	1.413	4.6	1.429	4.5
8 x	0.865	2.8	0.807	2.6	0.856	2.7
16 x	0.430	1.4	0.427	1.4	0.410	1.3
32 x	0.234	0.8	0.234	0.8	0.248	0.8
64 x	0.129	0.4	0.133	0.4	0.122	0.4
128 x	0.086	0.3	0.082	0.3	0.089	0.3
negative	0.040		0.040		0.052	

Moreover the BBI's panel # PHE 102 was also examined in three lots of product; data are reported below with reference to a European kit (BBI's results).

BBI – Panel code PHE 102

	Lot # 0103	Lot # 0103/2	Lot # 0303	Sorin EIA
Member	S/Co	S/Co	S/Co	S/Co
01	6.7	6.3	6.5	2.0
02	11.3	10.0	10.7	6.1
03	9.5	7.2	8.4	3.0
04	5.8	3.4	4.1	2.1
05	11.3	11.4	11.2	3.1
06	12.1	11.6	11.8	4.1
07	0.1	0.1	0.1	0.2
08	9.2	8.5	8.8	2.3
09	12.2	11.7	11.9	4.2
10	11.7	10.2	10.8	2.8
11	5.9	5.8	5.8	2.1
12	12.7	11.4	11.7	5.2
13	11.6	11.0	11.3	3.6
14	7.0	6.3	6.6	2.3
15	12.4	11.5	11.8	4.5

2. Diagnostic Sensitivity:

It is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

The diagnostic sensitivity has been tested internally and externally in a qualified Clinical Laboratory on panels of samples classified positive by a US FDA approved kit.

Positive samples were collected from different patients and from different HBV pathologies (acute and chronic hepatitis).

An overall value > 98% has been found in the study conducted on a total number of more than 200 samples.

A Seroconversion panel produced by BBI, USA, code # PHM 935A, has also been studied; results are reported below with reference to two commercial kits (BBI's results).

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3. Diagnostic Specificity:

It is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

The diagnostic specificity has been determined internally and externally in a qualified Clinical Laboratory on panels of negative samples from normal individuals and blood donors, classified negative with a US FDA approved kit.

A total number of more than 400 negative specimens were tested. A diagnostic specificity > 98% has been found.

Moreover, the diagnostic specificity was assessed by testing more than 50 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolyzed, lipemic, etc.).

No interference was observed in the study.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been

used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

Frozen specimens have also been tested to check whether this interferes with the performance of the test. No interference was observed on clean and particle free samples.

4. Precision:

It has been calculated on three samples examined in 16 replicate in three different runs, carried out on three different lots. The values found were as follows:

BCM.CE: lot # 0103

Cal 0 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.055	0.053	0.051	0.053
Std.Deviation	0.005	0.006	0.005	0.006
CV %	9.9	12.3	10.7	10.9

Cal 5 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.324	0.308	0.321	0.318
Std.Deviation	0.022	0.018	0.024	0.021
CV %	6.8	5.7	7.5	6.7

Cal 50 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.109	2.048	2.052	2.070
Std.Deviation	0.101	0.088	0.136	0.109
CV %	4.8	4.3	6.7	5.2

BCM.CE: lot # 0103/2

Cal 0 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.057	0.053	0.054	0.055
Std.Deviation	0.005	0.005	0.004	0.004
CV %	8.3	9.0	7.3	8.2

Cal 5 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.332	0.331	0.322	0.328
Std.Deviation	0.017	0.018	0.016	0.017
CV %	5.0	5.5	4.9	5.1

Cal 50 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.311	2.208	2.212	2.244
Std.Deviation	0.110	0.090	0.095	0.098
CV %	4.7	4.1	4.3	4.4

BCM.CE: lot # 0303

Cal 0 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.043	0.042	0.040	0.042
Std.Deviation	0.004	0.005	0.004	0.004
CV %	10.3	11.1	10.9	10.8

Cal 5 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.320	0.326	0.314	0.320
Std.Deviation	0.023	0.024	0.026	0.024
CV %	7.1	7.4	8.2	7.6

Cal 50 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.150	2.163	2.092	2.135
Std.Deviation	0.057	0.067	0.076	0.067
CV %	2.6	3.1	3.6	3.1

Important note:

The performance data have been obtained proceeding as the reading step described in the section M, point 12.

S. LIMITATIONS

Frozen samples containing fibrin particles or aggregates may generate false positive results.

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.

This test is suitable only for testing single samples and not pooled ones.

Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.

REFERENCES

1. Engvall E. and Perlmann P.. J. Immunochemistry, 8, 871-874, 1971
2. Engvall E. and Perlmann P.. J.Immunol. 109, 129-135, 1971
3. Remington J.S. and Klein J.O.. In "Infectious diseases of the fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
4. Volk W.A.. In "Essential of Medical Microbiology". 2nd ed. pp 729, G.B.Lippincott Company, Philadelphia, New York, S.Josè, Toronto.
5. Snyderman D.R., Bryan J.A. and Dixon R.E.. Ann.Int.Med., 83, pp 838, 1975.
6. Barker L.F., Gerety R.J., Lorenz D.E.. Viral Hepatitis. 581-587, 1978.
7. Cossart Y.. Brit.Med.Bull.. 28, pp 156, 1972
8. Lander J.J., Alter H. and Purcell R.. J.Immunol.. 106, pp 1066, 1971
9. Mushawar I.K., Dienstag J.L., Polesky H.F. et al.. Ann.J.Clin.Pathol.. 76, pp 773, 1981.
10. Grebenchtchikov N. et al.. J.Immunol. Methods, 15(2) :219-231, 2002
11. Schrijver RS and Kramps JA, Rev.Sci.Tech. 17(2):550-561, 1998

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

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#204.

PERFORMANCE EVALUATION KIT: PIPETTOR - WASHER

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30 plates **REF** 89897

KIT FOR THE EVALUATION OF PIPETTING SYSTEMS AND
MICROPLATE WASHERS USING A COLORIMETRIC METHOD

883529 - 2014/05

BIO-RAD

TABLE OF CONTENTS

1. INTENDED USE	3
2. PRINCIPLE OF THE TEST	3
3. KIT COMPONENTS	4
4. PRECAUTIONS.....	4
5. HEALTH AND SAFETY INSTRUCTIONS	4
6. MATERIALS REQUIRED BUT NOT PROVIDED	5
7. PREPARATION OF REAGENTS	5
8. STORAGE CONDITIONS - SHELF LIFE.....	5
9. ASSAY PROCEDURE	5
10. RESULTS CALCULATION AND INTERPRETATION	5

1. INTENDED USE

The Pipettor-Washer Performance Evaluation kit is a colorimetric assay kit for the evaluation of precision and accuracy, of the volume dispensed by pipettor and microplate washer, and of the residual volume after aspiration by microplate washers.

To ensure conformity to original manufacturer's conditions and specifications, and to guarantee reliable operation, instruments used for *in-vitro* diagnostics testing must be evaluated at installation, and at regular intervals. Both accuracy and precision must be tested. Any variation from the expected value should lead to proper corrective actions until acceptable values are obtained.

The test passed/failed criteria, and the corrective action guide, both depend upon the equipment being tested, and are detailed within each equipment verification procedure.

2. PRINCIPLE OF THE TEST

The Pipettor - Washer Performance Evaluation kit consists of dye solutions specifically adapted to measure various liquid volume dispenses (10-300 μl) or residual volumes after aspiration (0-15 μl). The kit contains 3 empty bottles to be filled with clean water plus 16 bottles of dye (for 10 plate kit) or 21 bottles of dye (for 30 plate kit) with different concentrations of dye solution. The dye is chosen to conform to the Beer-Lambert law to a high degree over a wide dynamic range, and to be totally safe for the user and environment.

To test dispensing performance, the equipment dispenses the appropriate dye volume in a microplate. After the microplate wells have been thoroughly mixed, a microplate reader measures the absorbance value, and calculates the amount dispensed. To test aspirating performance, the equipment aspirates the content of a microplate previously filled with dye. Water is then added to the residual dye left in the wells. After mixing, the microplate reader measures the absorbance value of the resulting solution, and calculates the residual volume with reference to a linear curve which is derived from wells manually pipetted by a qualified technician. The statistics of the run can then be calculated.

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3. KIT COMPONENTS

All reagents are exclusively for *in vitro* diagnostic use.

Identification on label		Description	Quantity	
			89894	89897
	Microplate	Microplate Empty, Greiner 8 x 12 microplate	10 plates	30 plates
	Clean Water	Empty 15-ml Vial Empty vial. (To be filled with clean water)	3 vials (Empty)	3 vials (Empty)
C0	Color Solution C0	Color solution C0 Dye solution in Tris NaCl buffer containing Tween® 20 Preservative: ProClin 300 ($\leq 0.1\%$)	4 vials (14 ml)	9 vials (14 ml)
C1	Color Solution C1	Color solution C1 Dye solution in water Preservative: ProClin 300 ($\leq 0.1\%$)	3 vials (2 ml)	3 vials (2 ml)
C2	Color Solution C2	Color solution C2 Dye solution in water Preservative: ProClin 300 ($\leq 0.1\%$)	3 vials (2 ml)	3 vials (2 ml)
C3	Color Solution C3	Color solution C3 Dye solution in water Preservative: ProClin 300 ($\leq 0.1\%$)	3 vials (10 ml)	3 vials (10 ml)
C4	Color Solution C4	Color solution C4 Dye solution in water Preservative: ProClin 300 ($\leq 0.1\%$)	3 vials (10 ml)	3 vials (10 ml)

4. PRECAUTIONS

The reliability of the results depends on correct implementation of the following Good Laboratory Practices:

- Do not use a kit beyond its expiration date.
- Do not interchange reagents between kit lots.
- Approximately 30 minutes before beginning the assay, remove the kit from the refrigerator (2°C-8°C) and allow the kit components to reach room temperature (18°C-30°C). Mix the reagents thoroughly by gently inverting the container several times before use. Once the testing is completed, the remaining assay materials should be immediately refrigerated (2°C-8°C).
- Use glassware thoroughly washed and rinsed with deionized water or, preferably, disposable material.

5. HEALTH AND SAFETY INSTRUCTIONS

- All the reagents included in the kit are intended for "*in vitro* use".
- Wear personal protective equipment while handling all reagents and while operating the washer and reader.
- Do not pipette by mouth.
- This product does not contain any biohazardous materials. However, it is used to check the correct operation of equipment which is regularly used with biohazardous material, and it is normal to assume that it comes in contact with contaminated parts (tubing, work surfaces, wastes tanks, manifold, absorbance reader, etc.) during the course of its use. As a consequence, waste material should be considered biohazardous when disposing or treating.

- Clean up all spills immediately and thoroughly. Rinse with water or a mild detergent. Disinfect the area in case of spills involving biohazardous materials. Dispose of all contaminated materials appropriately.
- Dispose of all waste in accordance with applicable national and/or local regulations.
- Chemical reagents should be handled in accordance with Good Laboratory Practices.

WARNING: Components C0, C1, C2, C3 and C4 contain $\leq 0.1\%$ ProClin 300

H317: May cause an allergic skin reaction.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352: IF ON SKIN: Wash with plenty of soap and water.

P333+P313: If skin irritation or rash occurs: Get medical advice/attention.

P501: Dispose of contents/container in accordance with local/regional/national/international regulations.

- For hazard and precaution recommendations related to some chemical components in this test kit, please refer to the pictogram(s) mentioned on the labels and the information supplied at the end of instructions for use. The Safety Data Sheet is available on www.bio-rad.com.

6. MATERIALS REQUIRED BUT NOT PROVIDED

- Clean water
- Calibrated single and multichannel, electronic or mechanical pipettors (10 μ l to 500 μ l)
- Absorbent paper.
- Calibrated microplate reader : dual-wavelength photometer (492/620 nm).

7. PREPARATION OF THE REAGENTS

Note: Before use, allow reagents to reach room temperature (18°C-30°C).

1) Ready to use reagents

Reagent 1 (C0): Color Solution C0

Reagent 2 (C1): Color Solution C1

Reagent 3 (C2): Color Solution C2

Reagent 4 (C3): Color Solution C3

Reagent 5 (C4): Color Solution C4

2) Reagents needing preparation

Water vials: Fill the required number of empty vials with 14 ml of clean water.

8. STORAGE CONDITIONS - SHELF LIFE

The kit should be stored in a refrigerator (at 2°C to 8°C). When kept at storage temperature, each reagent contained in the Pipettor – Washer PE Kit can be used until the expiration date on the package.

Once the testing is completed, remaining assay materials should be immediately refrigerated (at 2°C to 8°C).

9. ASSAY PROCEDURE

Strictly follow the manufacturer's procedure.

The assay procedure depends on the instrument you intend to test. Refer to each instrument's performance evaluation procedure.

10. RESULTS CALCULATION AND INTERPRETATION

Result calculations and interpretations depend on the tested instrument and on the specific verification procedure used. Refer to each instrument's performance evaluation procedure.

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Highlights of Immunomat™

- **SERION ELISA Analyzer** for medium to high throughput applications
- Validated for **antibody detection** in serum, plasma or cerebrospinal fluid (CSF), for **avidity determination** and **antigen detection** with **SERION ELISA immunoassays**
- **Short loading times** by **barcode identification** of samples, reagents and microtiter plates
- Processing of **up to 16 different assays per plate**
- **Reload function** for patient samples, reagents and microtiter plates
- **2D hand barcode scanner** for **parameters of quality control certificates**
- **Multishot dispenser function** and **memory function for tip racks**
- **Optimize function** for scheduling for an **efficient workload** of the instrument
- **Clot detection** and **bubble kill function** guarantee **walk-away functionality**
- **Level sensors for fluid containers** with **automated warning** in case of **lack of reagents or filled waste container**
- **Easy access to buffers and reagents**
- Fast and **quantitative evaluation** of SERION ELISA immunoassays
- **Parameter- or patient-orientated result reports**
- **Listing of test reagents**
- **Import/export function for reagents**
- Integrated **reports for quality controls** of standards and controls
- **Bi-directional connection to laboratory software systems** via ASTM interface

Order Information of Immunomat™ and Lab Ware

Immunomat™

Order Nr.: VT 020

SERION Clean, cleaning solution, 500 ml, 5 x conc.

Order Nr.: VT 125

Deep well microtiter plates, 50 pieces

Order Nr.: VT 124

Pipetting tips, 300 µl, 18 x 960 tips

Order Nr.: VT 111

Pipetting tips, 1100 µl, 10 x 960 tips

Order Nr.: VT 112

Plastic bottles, white, with cap, 35 ml, 100 pieces

Order Nr.: VT 113

Plastic bottles, yellow, with cap, 35 ml, 100 pieces

Order Nr.: VT 114

Plastic bottles, white, with cap, 50 ml, 100 pieces

Order Nr.: VT 115

Glas bottles, with cap, 2.5 ml, 323 pieces

Order Nr.: VT 116

Glas bottles, with cap, 2.5 ml, 30 pieces

Order Nr.: VT 116-30

Waste bags, 10 pieces

Order Nr.: VT 051