# superDimension™

## Aspirating Needle

<b>en</b>	Instructions for Use
<b>fr</b>	Mode d'emploi
<b>de</b>	Gebrauchsanleitung
<b>es</b>	Instrucciones de uso
<b>it</b>	Istruzioni per l'uso
<b>sv</b>	Bruksanvisning
<b>cs</b>	Návod k použití
<b>da</b>	Brugsanvisning
<b>no</b>	Bruksanvisning
<b>tr</b>	Kullanma Talimatı

1.12 pirkimo dalis

**REF** SDAN1000

**REF** SDAN2000

1.11 pirkimo dalis

**REF** AKI00101-01

**THE BACK OF COVER IS BLANK - DO NOT PRINT ON THIS PAGE**

## Symbols on Product and Packaging

Sterile using ethylene oxide

Single use

1.11 ir 1.12 pirkimo  
dalys

Federal law (USA) restricts this device to sale by or on the order of a physician

Do not resterilize

Do not use if package is opened or damaged

Consult instructions for use

Keep dry

Authorized representative in the European Community

Catalog number

Manufacturer

Use by

Batch code

Date of manufacture

Package quantity

Pouch quantity

CE mark. Designates that the labeled product is authorized for sale in European countries

Not made with natural rubber latex

en

## superDimension™

**REF** SDAN1000 **Pre-marked Aspiring Needle**  
21 G (0.813 mm)

**REF** SDAN2000 **Pre-marked Aspiring Needle**  
19 G (1.067 mm)

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1.11 pirkimo dalis.  
adatos dydis - 21 G;

### Indications for Use

To be utilized during flexible bronchoscopy by a trained physician to obtain a cell aspirate in the bronchial tree.

### Contraindications

Include but are not limited to:

- Patient whose general health and/or respiratory function are compromised to the point that the patient would not tolerate flexible bronchoscopy
- Uncooperative patient

### Description

The superDimension aspiring needle features a positive stop retraction/extension design. This design provides protection during endoscope passage while allowing accurate and quick release when the needle is at the targeted site. This device has been sterilized using ethylene oxide.

### Warning

Potential adverse effects include but are not limited to: fever, minor bleeding, pneumothorax, infection, bacteremia, and pneumomediastinum.

**Warning**

This device has been validated and warranted for single-use only. Any re-use, reprocessing, or re-sterilization could result in device failure causing patient or user injury, illness, or death. Re-use, reprocessing, or re-sterilization could cause contamination of the device, which could result in patient infection.

A thorough understanding of the technical principles, clinical applications, and risks associated with pulmonary endoscopy, and with the procedure for which this device is intended, is necessary.

**Caution**

Use of this device is restricted to devices or equipment with a biopsy channel minimum inside diameter of 2.08 mm (0.082 inch).

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Excessive pressure or force applied may cause damage to either the device or the biopsy channel.


**Instructions for Use**

1. Examine the package carefully before opening to confirm its integrity and that the use by date has not passed.
2. Remove the superDimension aspirating needle from the package. Remove the small tip protector from the cannula. Remove the protective tape from the proximal end of the device.
3. Inspect for any functional abnormality. If any irregularities are found, call Covidien for a return authorization number.
4. Ensure that the needle extends and retracts smoothly by sliding the proximal end gently. Note the position stops in the extended and retracted positions.

**Important**

No twisting or turning of the cannula hub handle is necessary when using this device.

- With the needle in the retracted position, securely attach a luer-lock syringe to the luer-lock fitting of the cannula hub handle. The cannula may be primed with sterile saline at this time. Insert the distal tip of the device, with the needle in the retracted position, into the superDimension extended working channel. Advance the device using short, 2 cm strokes until reaching the appropriate mark on the catheter body (see illustration). Note: Use standard endoscopic technique when using the unmarked aspirating needle.

- To aspirate a cell specimen, with the needle extended into the specimen site, quickly and fully retract the syringe plunger. To retain the specimen, maintain consistent negative pressure.
- Once the specimen has been obtained, retract the needle into the cannula. Remove the device from the biopsy channel with the needle remaining in the retracted position.
- To retrieve the specimen, release the syringe plunger to cease negative pressure. Prepare the specimen using the preferred method.
- Dispose of the superDimension aspirating needle in accordance with hospital or clinic protocol.

### Further Information

If further information is required, please contact Covidien.

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DLB00586-01 Rev. B  
PLB00460-01 Rev. B

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, <sup>TM</sup> brands are trademarks of their respective owner.

© 2014 Covidien.

Manufactured for:  
Covidien Inc,  
15 Hampshire Street  
Mansfield, MA 02048 USA.

Covidien Ireland Limited  
IDA Business and Technology Park,  
Tullamore.

Made in USA. Printed in USA.

[www.covidien.com](http://www.covidien.com)  
[www.superdimension.com](http://www.superdimension.com)

2014-12