COVIDIEN™

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REF SDMDK1000

superDimension™

Marker Delivery Kit

en Instructions for Use

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sv Bruksanvisning

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Sterile using
ethylene oxide

Does not contain
natural rubber
latex components

Single use

Federal law (USA)
restricts this
device to sale by
or on the order of
a physician

Do not resterilize

Do not use if
package is opened
or damaged

Consult
instructions for use

Warning or caution

Keep dry

CE mark.
Designates that
the product is
authorized for
sale in European
countries

Authorized
representative
in the European
Community

Catalog number

Manufacturer

Use-by

Batch code

Date of
manufacture

Package quantity

Pouch quantity

Length of catheter
or guidewire

French size of
catheter

Compatible
guidewire (inches)

Diameter of
guidewire

superDimension™ Marker Delivery Kit

Indications for Use

Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

Contraindications

Do not use with fiducial markers that exceed 0.9 mm (0.036 inch) outer diameter.

Description

The superDimension marker delivery kit contains (1) catheter and (1) guidewire. Both components are supplied sterile, non-pyrogenic, and are intended for single use.

Warning

This device has been validated and warranted for single use only. Any reuse, reprocessing, or resterilization could result in device failure causing patient or user injury, illness, or death. Reuse, reprocessing, or resterilization could cause contamination of the device, which could result in patient infection.

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Inspect Package Contents

Examine the package carefully before opening to confirm its integrity and that the use by date has not passed. The device is provided in a sterile package. Do not use the marker delivery kit if the package is damaged, opened, or if the use by date has passed. The marker delivery kit is sterilized by ethylene oxide. Do not resterilize.

Instructions for Use


1. **Unpack.** Remove the catheter and guidewire from the package.
2. **Insert the marker delivery catheter into the extended working channel.** After successful navigation to the desired site in the lungs, place the catheter into the proximal hub of the extended working channel and advance catheter using short, 2 cm strokes until reaching the appropriate mark on the catheter body (see illustration).
3. **Confirm the location.** Use fluoroscopy to confirm the location of catheter distal end.
4. **Load the marker.** Load a single compatible fiducial marker into the proximal hub of the marker delivery catheter.
5. **Place the marker into tissue.** Place the floppy end of the guidewire into the proximal hub of the catheter and advance the guidewire using short, 2 cm strokes until the guidewire reaches

the distal tip. Slowly advance the guidewire to push the fiducial marker into tissue.

6. **Remove the guidewire from the catheter.**
7. **Place additional markers.** If additional markers are desired, rotate the catheter to facilitate a different placement orientation and repeat steps 3 through 6.
8. **Remove the catheter from the extended working channel.**
9. **Disposal.** When the procedure is complete, dispose of the catheter and guidewire in accordance with hospital or clinic protocol.

Further Information

For information on the superDimension navigation system, refer to the *superDimension Navigation System User Manual*. For more information on the ILLUMISITE platform, refer to the *ILLUMISITE Platform User Manual*. If further information is required, please contact Medtronic.

* Second set of reference marks for EDGE tools are compatible for use with EDGE and ILLUMISITE tools.

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IFU SDI-0564-01 Rev. A 2020-05
P/N 490085003/001

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, TM* brands are trademarks of their respective owner.
©2020 Covidien.

May be covered by U.S. Patents:
[covidien.com/patents](https://www.covidien.com/patents)

Manufactured for:
Covidien LLC,
15 Hampshire Street
Mansfield, MA 02048 USA.

1-888-586-4767
[covidien.com](https://www.covidien.com)
[superdimension.com](https://www.superdimension.com)
[illumisite.com](https://www.illumisite.com)

Product of USA (from US and foreign materials). Printed in USA.

Covidien Ireland Limited
IDA Business and Technology Park,
Tullamore, Ireland.