# SuperLock™

Fiducial Marker

Nitinol Coil

4 mm x 7 mm

1.7 pirkimo dalis.  
Nitolono ritė

en	Instructions for Use
fr	Mode d'emploi
pt	Instruções de uso
de	Gebrauchsanleitung
es	Instrucciones de uso
it	Istruzioni per l'uso

nl	Gebruiksaanwijzing
sv	Bruksanvisning
ru	Инструкция по применению
zh	使用说明
cs	Návod k použití
da	Brugsanvisning

ko	사용 방법
no	Bruksanvisning
pl	Instrukcja stosowania
tr	Kullanma Talimatı
zh-t	使用說明

THE BACK OF COVER IS BLANK • DO NOT PRINT ON THIS PAGE

## Symbols on Product and Packaging

Sterilized using ethylene oxide

Not made with natural rubber latex

Do not re-use

Do not re-sterilize

Do not use if package is opened or damaged

Consult instructions for use

MR conditional

Keep dry

1.7 pirkimo dalis.  
Vienkartinis, sterilus

Authorized representative in the European Community

Catalog number

Manufacturer

Use-by

Batch code

Date of manufacture

Package quantity

Pouch quantity

CE mark. Designates that the product is authorized for sale in European countries

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### Indications for Use

The SuperLock fiducial marker is intended to be used to radiographically mark soft lung tissue for future surgical or therapeutic purposes.

### Contraindications

- No more than 10 devices should be used per patient
- Devices should not be used in an MRI machine with magnetic coils that exceed 3T

### Description

The SuperLock fiducial marker is constructed with a gold tube mechanically fastened to a nitinol wire that has been formed into a coil shape. This marker is designed to be delivered through a catheter with a minimum inner diameter of 0.038 in. (0.97 mm) using a 0.038 in. (0.97 mm) guidewire. The SuperLock fiducial marker comes preloaded in a delivery cartridge for transfer into the catheter.

### Warning

This device has been validated and warranted for single-use only. Any re-use, reprocessing, or re-sterilization could result in device failure causing patient or user injury, illness, or death. Re-use, reprocessing, or re-sterilization could cause contamination of the device, which could result in patient infection.

Care should be taken while advancing the guidewire near the pleura of the lungs to avoid a pneumothorax.

### Caution

To avoid risk of marker pre-deployment, do not advance the SuperLock fiducial marker before the delivery cartridge is mated with marker delivery catheter.

Do not use with fiducial marker delivery tools with an inside diameter of less than 0.038 in. (0.97 mm).


Do not use if package is opened or damaged, as product sterility may be compromised.

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
1.7 markeris įvedamas  
endoskopiškai


Fiducial marker as preloaded in delivery cartridge

Fiducial marker as deployed

Outside diameter  
4.1 mm

## MR Conditional

MR conditional

- Maximum spatial gradient magnetic field of 720-Gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning

## MRI-Related Heating

In non-clinical testing, the SuperLock fiducial marker produced a temperature rise of equal to or less than 1.5 degrees C at a maximum whole body absorption rate (SAR) of 2.7 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3-Tesla/128-MHz, Excite, HDx™\*, Software 14.XM5, General Electric Healthcare System MR scanner.

## Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the SuperLock fiducial marker. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this device. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5 mm relative to the size and shape of the SuperLock fiducial marker.

## Implant Procedure

Use this implant procedure with the superDimension™ marker delivery kit (model SDMDK1000 510(k) K110093, available from Covidien).

1. Inspect for any damage. If any irregularities are noted, contact Covidien for a return authorization number.
2. Verify that the marker delivery catheter is at the desired location.
3. Open the pouch and remove the delivery cartridge with the preloaded fiducial marker.
4. Place the distal end of the delivery cartridge into the proximal end of the delivery catheter hub and rotate the delivery cartridge until it is fully seated on the luer feature of the catheter hub.

### Notice

To avoid fiducial marker pre-deployment, ensure that the delivery cartridge is fully seated in the delivery catheter hub.

5. Ensuring that the delivery cartridge is fully seated in the delivery catheter hub, push the floppy end of the guidewire through the delivery cartridge to advance the fiducial marker into the catheter.

6. Advance the guidewire using short, 2 cm strokes until the fiducial marker reaches the distal tip. Slowly advance the guidewire until the fiducial marker is deployed at the desired location.
7. Remove the guidewire and delivery cartridge from the catheter.
8. Repeat steps 2 through 7 to place additional fiducial markers.

**Disposal:** When the procedure is complete, dispose of the delivery cartridge in accordance with hospital or clinic protocol.

**Storage:** The SuperLock fiducial marker is supplied sterile and non-pyrogenic in an unopened inner package. Store in a cool, dry place.

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COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, <sup>TM</sup>\* brands are trademarks of their respective owner.

© 2014, 2018 Covidien.

May be covered by U.S. Patents:  
[www.covidien.com/patents](http://www.covidien.com/patents)

Manufacturer:  
Covidien llc,  
15 Hampshire Street  
Mansfield, MA 02048 USA.  
1-888-586-4767  
[covidien.com](http://covidien.com)  
[superDimension.com](http://superDimension.com)

Made in USA. Printed in USA.  
2018-01

Covidien Ireland Limited  
IDA Business and Technology Park,  
Tullamore.

IFU No. SDI-0119-01 Rev. A