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BIBLIOGRAFIA

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9^a edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- CLSI H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture: Approved Standard, 6^a edition, 2007.
- CLSI H18-A4: Procedures for the handling and processing of blood specimens: Approved Guideline, 4^a edition, 2010.
- Technical Manual, 18^a edition, American Association of Blood Banks, Bethesda, Maryland, 2014.
- Judd WJ, Johnson ST, Story JR. Judd's Methods in immunohematology, 3^a edition, AABB Press Bethesda, 516, 2008.

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
Figure 1. Pattern of reaction grades

Stability of the results

After centrifuging the cards it is recommended that the results be read immediately. Do not leave processed cards in a horizontal position. If necessary, a delayed reading can be performed up to 24 hours after processing the cards if they are kept in an upright position, refrigerated (2 - 8 °C) and sealed with a laboratory covering film to avoid evaporation of the supernatant.

Interpretation of the results

ABO system. The expected reaction with microtubes A and B and its interpretation are shown in the following table (+ = positive, and 0 = negative).

Microtute A	Microtute B	Microtute CtL	Interpretation
0	0	0	O group
+	0	0	A group
0	+	0	B group
+	+	0	AB group

D antigen. A positive result in the corresponding microtube indicates the presence of antigen D (Rh system).

Notes:

- The acronym "CtL" means Control.
- The CtL microtube should be negative. If it is positive, due to the formation of rouleaux, to strong cold autoagglutinins or other causes, invalidate the test. Repeat the determination after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the CtL microtube of the repeat test is negative, the results of the test can be interpreted; if it is positive, invalidate the test.
- The ABO group obtained should be verified comparing the result obtained with a previously typed sample. If a previously typed sample is not available, it is recommended to determine the reverse group. Forward (cell grouping) discrepancies should be investigated before releasing the result.
- To verify a negative status or to ensure detection of weak and partial D, other reagents and techniques (e.g. Indirect Antiglobulin testing) which may detect different weak and partial D variants should be used.
- Precaution should be taken at interpretation of Double Population events. Not all mixed-cells situations are detected. Additional information on patient history and additional testing will be necessary for resolution. Transfused patients or those subjected to bone marrow transplant may present images of Double Population. Double Population is also observed in some ABO subgroups (A₁), in Tn cryptantigens, in blood chimerism in fraternal twins, and in the very rare case of mosaicism arising from dispemias¹.
- The observation of complete or partial hemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of collection and/or handling of the sample.

LIMITATIONS OF THE PROCEDURE

- Grossly hemolyzed, cloudy or contaminated samples or samples with presence of a clot, may cause false positive or false negative results.
- Aged or hemolyzed specimens may cause weaker reactions compared to those obtained with a fresh sample.
- Samples with high-potency antibodies may coat the red blood cells completely, causing spontaneous agglutination¹.
- Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the nonspecific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test¹.
- Red blood cells from individuals with A or B variants may present a weak expression of the antigens, and may not be detected.
- Antigen expression may be weakened in red blood cells of people with leukaemia or other malignant diseases¹.
- A very weak expression or variants of the D antigens may not be detected.
- The anti-A reagent contained in this card could react with Tn cryptantigens.
- On occasions, unagglutinated red blood cells may be retained somewhere in the gel column with the appearance of very minute red dot or flecks. However, this nonspecific retention should not interfere with the interpretation of the result.

IVD	In vitro diagnostic medical device
LOT	Batch code
	Use by
	Temperature limitation
	Consult instructions for use
REF	Catalogue number
	Cards
	Manufacturer
	This way up
	Fragile, handle with care
	Keep dry

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ABO/Rh system

The diagnostic sensitivity and specificity of the antibodies present in the DG Gel Confirm card, for the determination of the antigens of the ABO and Rh systems, have been studied in a representative number of positive and negative samples.

Antibody	No. of Samples	Sensitivity ^{a)}	Specificity ^{b)}
Anti-A	3041	100%	100%
Anti-B	3046	100%	100%
Anti-D ¹⁾	3057 ¹⁾	99,50%	100%

- (a) Sensitivity: (number of true positive results / (number of true positive results + number of false negative results)) x 100
(b) Specificity: (number of true negative results / (number of true negative results + number of false positive results)) x 100
(c) These samples include 52 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D¹⁾ reagent is 73,10%.

Precaution

The precision of the reagents present in the DG Gel Confirm card has been determined, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

BIBLIOGRAPHY

- Mollison P, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9^a edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- CLSI H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture: Approved Standard, 6^a edition, 2007.
- CLSI H18-A4: Procedures for the handling and processing of blood specimens: Approved Guideline, 4^a edition, 2010.
- Technical Manual, 18^a edition, American Association of Blood Banks, Bethesda, Maryland, 2014.
- Judd WJ, Johnson ST, Story JR. Judd's Methods in immunohematology, 3^a edition, AABB Press Bethesda, 516, 2008.

PRESENTATION

210339 DG Gel Confirm 50 Cards Profile: 2x(A/B/D¹⁾/CtL

Manufactured by:

Diagnostic Grifols, S.A.

Passeig Fluvial 24, 08150 Parets del Vallès (Barcelona), Spain

Date of last version: October 2016

This document is available in several languages. The translations have been made from the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

SYMBOLS KEY

One or more of these symbols may have been used in the labeling/packaging of this product.

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. CLSI H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
4. CLSI H18-A4: Procedures for the handling and processing of blood specimens; Approved Guideline, 4th edition, 2014.
5. Technical Manual, 18th edition, American Association of Blood Banks, Bethesda, Maryland, 2014.
5. External evaluation Diagnostic Grifols, S.A. (REGD-0002998)M Noviembre 2007.
6. Judd WJ, Johnson ST, Storry JR. Judd's Methods in immunohematology, 3rd edition, AABB Press Bethesda, 516, 2008.

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Read carefully before using the product. For use in “in vitro” diagnosis only

INTENDED USE

Determination of the antigens of the Rh and Kell system, in gel technique.

INTRODUCTION

The clinical importance of the Rh blood group stems from the fact that the D antigen is highly immunogenic. The D (RH1), C (RH2), E (RH3) , c (RH4) and e (RH5) antigens are the most important of the group of antigens that comprise this system¹.

The K (KEL1) antigen has a great clinical importance, since it is present in cases of haemolytic disease¹. The determination of the Rh phenotypes and K antigen is defined by the presence or absence of the antigens D (RH1), C (RH2), E (RH3), c (RH4), e (RH5) and C^w (RH8) and K (KEL1) in the red blood cells.

The anti-D^{W+}, anti-C, anti-E, anti-c, anti-e, anti-C^w and anti-K reagents are used to perform the Rh blood group and K antigen typing.

PRINCIPLE OF THE TEST

The principle of the test is based on the gel technique described by Y. Lapiere² for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

The DG Gel card is a plastic support consisting of 8 microtubes. Each microtube consists of a column and a dispensing/incubation chamber.

Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextrans are mixed with a reagent that contains specific antibodies or a buffer.

The microtubes that contain the specific antibodies into the gel solution act as a reaction medium and the red blood cells agglutinate in contact with the antibodies.

The microtubes without antibodies are used for controls.

During centrifugation and depending on their size, the red blood cell agglutinations are trapped in the surface or along the gel column. The non-agglutinated red blood cells descend to the bottom of the microtube.

REAGENTS

Each microtube of the DG Gel Rh Pheno+Kell card contains polymerised dextrans in buffered medium with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

– microtube D^{W+}: monoclonal anti-D (mixture of IgM antibodies of human origin, clones RUM-1 and ESD-1M).

This anti-D monoclonal reagent detects weak D and partial variants of the D antigen, including the D^W variant.

– microtube C: monoclonal anti-C (IgM antibodies of human origin, clone MS-24).

– microtube E: monoclonal anti-E (IgM antibodies of human origin, clone DEM-1).

– microtube c: monoclonal anti-c (IgM antibodies of human origin, clone H-48).

– microtube e: monoclonal anti-e (mixture of IgM antibodies of human origin, clones MS-21, MS-63 and MS-16).

– microtube C^w: monoclonal anti-C^w (IgM antibodies of human origin, clone MS-110).

– microtube Kell: anti-K monoclonal (IgM antibodies of human origin, clone MS-56).

– microtube Ctl.: buffered solution without antibodies (microtube control).

Reagent ready to use. Use the microtubes immediately once the seal has been opened.

STABILITY

DG Gel Rh Pheno+Kell is stable up to the expiry date stated on the label, with the seal intact and stored at 2–25 °C in the position indicated on the outer packaging. Do not freeze.

MATERIAL REQUIRED BUT NOT PROVIDED

– 10 µl, 50 µl and 1 ml automatic pipettes.

– Disposable pipette tips.

– Glass tubes.

– DG Gel Sol.

– Centrifuge for DG Gel cards.

SAMPLES

Blood samples recently collected with the routinely used anticoagulants (EDTA, citrate and heparin).

Do not use haemolysed, cloudy or contaminated samples or with clot presence.

The procedure of extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives^{3,4} and following the instructions of the manufacturer of the material used for collecting the sample.

• Determination of the antigens of the Rh/Kell system: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2–8 °C can be used up to 48 hours after their extraction.

Red blood cells from bags collected in CPD, CPDA or SAG–Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2–8 °C.

If red blood cells from the bag segment are used, it is recommended to wash them with physiological saline solution before preparing the suspension. Do not use if clots or haemolysis are observed.

TEST METHOD

DG Gel Rh Pheno+Kell can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Allow the samples and reagents to reach room temperature (18–25 °C).

Inspect the state of the cards before use (see section LIMITATIONS "In relation to the product").

Identify the cards and samples to be used.

Each card is individually identified by a bar code. Unless a bar code reader is used, identify them manually.

If manual method is used: carefully peel off the metal film that covers the microtubes to prevent cross-contaminations among them and carefully dispense the red blood cell suspension, avoiding the pipette tip to come into contact with the wall or the content of the microtubes.

Manual method:

1. Determination of the antigens of the Rh/Kell system (microtubes D^{W+}/C/E/c/e/C^w/Kell/Ctl.):

– Prepare a 5% red blood cell suspension in DG Gel Sol (50 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Ensure the resuspension of the red blood cells before use.

– Add 10 µl of 5% red blood cell suspension into each of the microtubes indicated.

2. Centrifuge in centrifuge for DG Gel cards.

3. Read the results.

RESULTS

Result reading:

Negative:	–	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
	+/–	Scarce small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
Positive:	2+	Small or medium-sized agglutinations throughout the column.
	3+	Upper band of medium-sized agglutinations in the upper half of the column.
	4+	Band of agglutinated red blood cells in the upper part of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column).

Stability of the results: it is recommended an immediate reading of the results after centrifuging the cards. Do not leave processed cards in horizontal position.

If necessary, a delayed reading can be made up until 24 hours after processing the cards if kept in vertical position, refrigerated (2–8 °C) and sealed with parafilm or similar material to avoid evaporation of the supernatant.

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+ = positive
0 = negative

Rh/Kell system (C, E, c, e, C^w and K antigens): A positive result in microtubes C, E, c, e, C^w and K indicates the presence of C, E, c, e, C^w or K antigens of the Rh/Kell system.

Notes:

1. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.

2. The microtube Ctl. must be negative. If it is positive, invalidate the test.

Repeat the test after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the microtube Ctl. of the repetition is negative, the results of the test can be interpreted; if it is positive, invalidate the test.

3. It is recommended to check/investigate any discrepant results obtained.

4. Rh system (D antigen): in the event of +/- to 3+ reactions, presence of a weak antigen must be investigated.

5. Rh system (C, E, c, e and C^w antigens): in the event of positive reactions below 3+, presence of a weak antigen must be investigated.

6. Kell System: positive reactions below 2+, may indicate presence of a weak antigen.

7. Negative reactions with the anti-D^{W+} reagent must be verified using other reagents and techniques which may detect different variants of the D antigen.

8. The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.

QUALITY CONTROL

1. It is recommendable to include known positive and negative controls in each series of tests.

2. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity:

The diagnostic sensitivity and diagnostic specificity of the antibodies present in the DG Gel Rh Pheno+Kell card to determine the antigens of the Rh and Kell system have been studied^{5,6} in a representative number of positive and negative samples.

Antibody	No. Of samples	Sensitivity ^(a)	Specificity ^(b)
Anti-D ^{VI+}	3057 ^{(c),5}	99.6%	100%
Anti-C	1044 ⁵	100%	100%
Anti-E	1048 ⁵	100%	100%
Anti-c	1048 ⁵	100%	100%
Anti-e	1047 ⁵	100%	100%
Anti-C ^w	1049 ⁵	100%	100%
Anti-K	317 ⁶	100%	100%

- (a) Sensitivity: (number of true positive results/(number of true positive results + number of false negative results)) x 100
(b) Specificity: (number of true negative results/(number of true negative results + number of false positive results)) x 100
(c) These samples include 52 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D^{VI+} reagent is 78.8%.

Precision:

The precision of the reagents present in the DG Gel Rh Pheno+Kell card has been determined^{7,8}, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

LIMITATIONS

In relation to the sample:

Do not use haemolysed, cloudy or contaminated blood samples or with clot presence.

Sample: red blood cells

- 1.Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test⁹.
- 2.Transfused patients or those subjected to bone marrow transplant may present images of double population⁹.
- 3.Patients with high-potency autoantibodies may coat the red blood cells completely, causing spontaneous agglutination⁹.

In relation to the product:

Before use, inspect the microtubes of the cards:

- 1.Do not use the card if microbiological contamination or change in microtube colour is observed.
- 2.If due to improper transport or storage, dispersed drops at the top of the microtubes, are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
- 3.If trapped bubbles in the gel, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.

WARNINGS AND PRECAUTIONS

The use of "in vitro" diagnostic reagents for professional use requires taking into account the following indications:

- The reagents of the DG Gel Rh Pheno+Kell card of human monoclonal origin are manufactured using a non-reactive material for the HB_s antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin will not transmit hepatitis or AIDS. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

- 1.Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- 2.Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- 3.NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
- 4.NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
- 5.External evaluation Diagnostic Grifols, S.A. (REGD-0002998), November 2007.
- 6.External evaluation Diagnostic Grifols, S.A. (REGD-0012604), February 2013.
- 7.Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
- 8.Internal study Diagnostic Grifols, S.A.(DG-RI013/07-3), May 2007.
9. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.

PRESENTATION

210350 DG Gel Rh Pheno+Kell 50 cards Profile: D^{VI+}/C/E/c/e/C^w/Kell/Ctl.

Revision date: February 2013

This document is available in several languages. The translations have been made from the master document in Spanish. In the event of doubts or discrepancies, the wording in the master document in Spanish shall take precedence.

	"In vitro" diagnostic medical device
	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
	Catalogue number
	Cards

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture;

Approved Standard, 6th edition, 2007.

4. NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.

5. External evaluation Diagnostic Grifols, S.A. (REGD-0002998), November 2007.

6. External evaluation Diagnostic Grifols, S.A. (REGD-0012604), February 2013.

7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.

8. Internal study Diagnostic Grifols, S.A.(DG-RI013/07-3), May 2007.



9. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.

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INTENDED USE

DG Gel Sol is a reagent for preparing red blood cell suspensions used in gel techniques (DG Gel).

PRINCIPLE

A reduction in the ionic strength of the medium of the red blood cell suspension enhances significantly the association of the antibody with the antigen. DG Gel Sol is a low ionic strength solution that facilitates the attachment of the antibody to the red blood cells by lowering the density of the cation cloud around them, thus improving sensitization.

COMPOSITION

Each DG Gel Sol kit contains:

DG Gel Sol:

- 2 vials of 100 ml, liquid. Each vial contains low ionic strength buffered solution, being its main components glycine 1.37% and glucose 0.85%.

Preparation of the reagents

DG Gel Sol:

Ready to use:

Stability once opened: at 2-8 °C, it has the same stability as when unopened, provided that it is handled and stored correctly, avoiding contamination of the contents and prolonged exposure to temperatures that are not suitable for storage.

Do not freeze.

STABILITY

All unopened reagents are stable up to the expiry date stated on the label if stored at 2-8 °C.

METHOD

DG Gel Sol can be used manually or with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Manual method (DG Gel):

Follow the instructions for use of the DG Gel card to be used, particularly the points that indicate the use of DG Gel Sol.

LIMITATIONS

Limitations in relation to the reagent:

Before use, inspect the reagent.

1. Do not use if the container is deteriorated and there is a loss of content.
2. Do not use if microbial or fungal contamination is observed.

WARNINGS AND PRECAUTIONS

The use of “in vitro” diagnosis reagents for professional use requires taking into account the following indications:

- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positives or negatives.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Follow the recommended procedures strictly, periodically checking the instruments used.
- Do not use the product beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have doubts or need further information on the use of this product, consult the authorised distributor in your country.

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IVD	In vitro diagnostic medical device
LOT	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
REF	Catalogue number

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Essential II Control

English 

Quality control for use with DG Gel technique, MDmulticard or conventional tube technique

SUMMARY AND PRINCIPLE

Blood Bank Quality Assurance includes a check of all reagents routinely employed in immunohematologic procedure on the day of use to ensure that the reagents are functioning properly and to verify techniques.

Essential II Control is intended to allow regular control of materials, work procedures and instrument procedures for (i) the determination of the ABO, Rhesus and K antigens, (ii) the determination of the appropriate ABO blood group antibodies and (iii) the detection of irregular antibodies in both antiglobulin and enzyme technique. The product is suitable for use with DG Gel technique, MDmulticard and conventional tube technique.

REAGENT

Warning: Do not use for grouping donors or recipients. Use only in the laboratory for quality assurance of blood bank reagents.

Essential II Control suspensions have been prepared from red blood cells collected from blood donors. Each tube may be prepared from a single donor or a pool of donors. The suspensions are composed of human erythrocytes (haematocrit $15 \pm 2\%$) in a buffered isotonic nutrient medium with preservatives (0.03% (w/v) neomycin and 0.05% (w/v) chloramphenicol) and added antibodies as appropriate which contain $< 0.01\%$ (w/v) sodium azide as a preservative. Frozen/thawed red blood cells may have been used in this product.

Each Essential II Control kit contains 2x2 tubes of 6 ml:

Tube	ABO	Rhesus	Other Typing	Antibodies
QCS1	A	ccddee	K pos	Anti-B, Anti-D (~0.05 IU/ml)
QCS2	B	CcD.Ee	K neg, Fy ^a neg	Anti-A, Anti-Fy ^a

Antibody	Clones used	Origin
Anti-A monoclonal	Birma-1	Murine
Anti-B monoclonal	LB2	Murine
Anti-D polyclonal	-	Human
Anti-Fy ^a monoclonal	P3TIM	Murine/human

Direct antiglobulin tests are negative on all red blood cell suspensions of Essential II Control.

For professional use only. Ready-to-use.

Caution: All blood products should be treated as potentially infectious. All human source material used in the manufacture of this product were tested and found negative for HBsAg and antibodies to HCV, HIV-1/HIV-2. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents, appropriate care should be taken in the use of this product. The absence of murine viruses has not been established.

Once used, the product must be disposed in special containers for biological waste.

Warning: This product contains sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. If discarded into sink, flush with a large volume of water to prevent azide build-up.

STABILITY

Do not expose to extreme temperatures. Store at 2-8°C when not in use. Turbidity may indicate reagent deterioration or microbial contamination; discard if present. Do not use beyond expiration date. If stored appropriately at 2-8°C the product is stable after the first opening until the indicated expiration date.

It is normal to observe a slight change of colour in the supernatant.

PROCEDURE

Reagent provided

Essential II Control, 2x2x6 ml, cat.no. 213287

Materials required but not provided

Consider the products and equipment related to the appropriate technique intended to be used.

The reagents must be brought to room temperature (18-25°C) prior testing.

PROCEDURE OUTLINE

Process the Essential II Control in the same way as patients' samples.

Manual method, DG Gel technique

Use the suspensions of Essential II Control as control samples with DG Gel reagents of Diagnostic Grifols S.A. and/or 0.8% Reagent Red Blood Cells of Medion Grifols Diagnostics AG or Diagnostic Grifols S.A.. Follow the instructions for use of the DG Gel card and 0.8% Reagent Red Blood Cells employed.

Manual method, MDmulticard

Use the suspensions of Essential II Control as control samples with MDmulticard reagents of Medion Grifols Diagnostics AG. Follow the instructions for use of the MDmulticard reagent employed.

Manual method, conventional tube technique

Use the suspensions of Essential II Control as control samples with tube technique reagents of Medion Grifols Diagnostics AG. Follow the instructions for use of the tube technique reagent and the Reagent Red Blood Cells employed.

Automated method

The suspensions of Essential II Control can be used in Diagnostic Grifols S.A. instruments for DG Gel. Follow the user's manual of the automatic instrument employed. When working with automated systems, it is recommended to perform a quality control after decontamination processes. Essential II Control is valid for this purpose.

QUALITY CONTROL

To ensure proper centrifugation, each individual centrifuge should be calibrated for the specific test procedure being performed.

RESULTS

Interpretation

If the test results differ from the awaited results listed in section "REAGENT", check immediately the working procedure, test method, automatic instrument/equipment and other related materials used. Repeat the test after having applied the appropriate corrections.

LIMITATIONS OF PROCEDURE

- Improper techniques may invalidate the results obtained with Essential II Control.
- False positive or false negative results can occur due to contamination of test material, improper reaction temperature and/or time, improper centrifugation, improper storage conditions of material and omission of test reagents.
- Do not use if evident changes of colour and/or obvious haemolysis is observed.
- Cells may lose antigen strength during the shelf life of the product giving weaker positive reactions than fresh cells.
- Some interaction between the Rh negative (D) red blood cells and the Anti-D reagent (clones: P3x290, P3x35, P3x61 and P3x21223B10) may be observed as a result of reagents' inherent properties.
- Take under consideration all limitation statements in the instructions for use of the used DG Gel cards, the 0.8% Reagent Red Blood Cells, the conventional tube reagents, the MDmulticard reagents and the automatic instruments.

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance of the product was validated according to the limited field of application (refer to intended use, section “SUMMARY AND PRINCIPLE”). The requirements of the Common Technical Specifications (CTS) 2009/886/EC are fulfilled as far as applicable.

According to the repeatability and reproducibility testing, the total number of tests (n) is displayed below and the sensitivity as well as the specificity was calculated for each technique for which an agreement of 100% for sensitivity and specificity was obtained.

Technique	Antigen Parameters: Control for Antibody Reagents							
	n determination							
	A	B	D	C	E	c	e	K
Gel	160	160	240	160	160	120	120	80
MDmulticard	80	80	160	80	80	80	80	80
Tube	40	40	40	40	40	40	40	40

Technique	Antibody Parameters: Control for Reverse Grouping Cells			
	n determination			
	Anti-A	Anti-B	Anti-AB	No Isoagglutinin
Gel	80	80	-	-
Tube	80	80	-	-

Technique	Antibody Parameters: Control for Antibody Screening Cells and Identification Cells	
	n determination	
	Anti-D	Anti-Fy ^a
Gel	120	120
Tube	60	60

Sensitivity: probability to obtain a positive result in presence of a positive sample.

Specificity: probability to obtain a negative result in presence of a negative sample.

WARRANTY

This product is warranted to perform as described in its labelling and in the product literature and Medion Grifols Diagnostics AG disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Medion Grifols Diagnostics AG be liable for any consequential damages arising out of the aforesaid express warranty.

IVD

CE

0123

Used symbols

IVD

"In Vitro" Diagnostic Medical Device

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Tel. 8 5 279 0080

Serigrup Diana A/B Serigrup Diana 4 Serigrup Diana A₂

Reagent Red Blood Cells for reverse grouping Instructions for use. For *in vitro* diagnostic use.

INTENDED USE

Serigrup Diana A/B, Serigrup Diana 4 and Serigrup Diana A₂ Reagent Red Blood Cells are used for the determination of the reverse ABO group in gel technique.

SUMMARY AND EXPLANATION

The ABO system is defined by the presence or absence of the A and/or B antigens in the red blood cells and by the presence of antibodies in the serum, corresponding to the antigen or antigens missing in the red blood cells.

The determination of the serum or reverse ABO group, is based on the presence or absence of anti-A and/or anti-B isagglutinins in the patient's serum, using red blood cells from a known blood group: A₁, A₂, B and O. Anti-A and anti-B isagglutinins are natural antibodies produced by individuals lacking the respective A and B antigens.

Some individuals with subgroups A₁ or A₂ may present anti-A antibodies, which interfere with the procedure for determining the reverse ABO group using A₂ Reagent Red Blood Cells. The Serigrup Diana A₂ reagent is used to investigate such discrepancies.

PRINCIPLE OF THE TEST

Anti-A and anti-B from human serum bind to red blood cells possessing the corresponding antigens, resulting in direct agglutination. Serigrup Diana Reagent Red Blood Cells are utilized in the gel technique to detect antibodies to human blood groups A and B.

REAGENTS

Observable indications

Inspect the conditions of the Reagent Red Blood Cells before use (see Warnings and Precautions). Do not use the product, if you detect clausiness and/or large-scale hemolysis, alterations or changes in color, microbiological contamination or other artifacts. Do not use the product, if the package, vials or caps are damaged or if contents have leaked.

Reagent Red Blood Cells

All Serigrup Diana Reagent Red Blood Cells consist of a 0.8% suspension ready-to-use of human red blood cells of A₁, B, A₂ and O blood group produced by a single donor for each vial. The cells come in a buffered solution with preservatives (neomycin 0.010% (w/v) and chloramphenicol 0.017% (w/v)). Each glass vial of Serigrup Diana contains 10 mL of human red blood cells. Vials come with a blue cap for Serigrup Diana A₂, a black cap for Serigrup Diana A₁, a yellow cap for Serigrup Diana B and a white cap for Serigrup Diana O.

Material required but not provided

Please refer to the Instructions for Use of Grifols gel cards.

STORAGE AND STABILITY

- Do not use beyond the expiration date.
- Store upright (as indicated by the two arrows on the outer packaging) at 2 - 8 °C.
- Do not freeze.

Notes:

- If you identify incorrect temperature conditions during storage or shipment do not use the product.
- To avoid contamination, close the vials with their caps when they are not being used. Ensure that the caps of the Reagent Red Blood Cell vials have not been swapped.
- Once a vial has been used, keep it at the indicated storage temperature.

WARNINGS AND PRECAUTIONS

- The results do not constitute a clinical diagnosis by themselves. Evaluate the results together with the patient's clinical information and other data.
- The product should only be used by qualified personnel.
- Any modifications of the test procedures described in these Instructions for Use require validation by the user.
- The use of volumes other than those indicated in the method, may modify the reaction and lead to incorrect test results, i.e. false-positive or false-negative results.
- Use disposable plastic pipette tips in order to avoid cross-contamination with samples and other reagents.
- Do not mix residues from different vials.
- Serigrup Diana Reagent Red Blood Cells are of human origin and are manufactured using materials that have been tested and found non-reactive for the HBs antigen, and for anti-HIV and anti-HCV antibodies. However, there is no known procedure to ensure that products of human origin will not transmit infectious agents. Human blood products and samples should be handled as if they were potentially capable of transmitting infectious diseases.

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LITERATŪRA

1. CLSI H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
2. CLSI H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
3. Technical Manual, 18th edition, American Association of Blood Banks, Bethesda, Maryland, 2014.

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/CE žyma/

Concentrated wash solution for Grifols analyzers. For laboratory use.

APPLICATION

DG Fluid B is used for the washing of the Grifols analyzers' fluidic system. The solution must be diluted prior to use.

REAGENT

Twelve 125 mL bottles of a surfactant concentrated solution with color added.

Preservative: 0.1% sodium azide in final diluted state.

Warning: The concentrated solution contains 1.6% sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azide. On disposal, flush with a large volume of water to prevent azide buildup.

STABILITY

The concentrated solution is stable until the expiration date indicated on the label when stored or used between 2°C and 25°C. In the final dilution state it is stable 4 weeks at room temperature (18°C to 25°C). However, the solution may suffer a loss of color during the product life without affecting its functionality.

PROCEDURE

Dilute the contents of one bottle of DG Fluid B to 2 litres with purified water for laboratory use or equivalent.

Diluted DG Fluid B should be poured into the appropriate container provided with the Grifols analyzers.

WARNINGS AND PRECAUTIONS

- 1. H302: Harmful if swallowed.
- 2. H312: Harmful in contact with skin.
- 3. P270: Do not eat, drink or smoke when using this product.
- 4. P264: Wash hands thoroughly after handling.
- 5. P280: Wear protective gloves/protective clothing/eye protection/face protection.
- 6. P301+P330+P312: IF SWALLOWED. Rinse mouth and immediately call a POISON CENTER or doctor / physician if you feel unwell.
- 7. P302+P352+P312: IF ON SKIN: Wash with plenty of soap and water and immediately call a POISON CENTER or doctor / physician if you feel unwell.
- 8. P363: Wash contaminated clothing before reuse.
- 9. P501: Dispose of contents / container according to local, state and national regulations.
- 10. Do not use beyond expiration date.

PRESENTATION

213678 DG Fluid B 12 x 125 mL

Revision date: September 2014

This document is available in several languages. The translations have been made from the original document in English. In the event of doubts or discrepancies, the wording in the original document in English shall prevail.

SYMBOLS KEY:

One or more of these symbols may have been used in the labeling/package of this product.

	Warning
	<i>In vitro</i> diagnostic medical device
	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
	Catalogue number

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APPLICATION

DG Fluid A is used for the washing of the Grifols analyzers' fluidic system. The solution must be diluted prior to use.

REAGENT

Twelve 125 mL bottles of a saline-based concentrated solution with color added.

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PROCEDURE

Dilute the contents of one bottle of DG Fluid A to 2 litres with purified water for laboratory use or equivalent.

Diluted DG Fluid A should be poured into the appropriate container provided with the Grifols analyzers.

WARNINGS AND PRECAUTIONS

1. H302:Harmful if swallowed.
2. H312: Harmful in contact with skin.
3. P270: Do not eat, drink or smoke when using this product.
4. P264: Wash hands thoroughly after handling.
5. P280: Wear protective gloves/protective clothing/eye protection/face protection.
6. P301+P330+P312: IF SWALLOWED. Rinse mouth and immediately call a POISON CENTER or doctor / physician if you feel unwell.
7. P302+P352+P312: IF ON SKIN: Wash with plenty of soap and water and immediately call a POISON CENTER or doctor / physician if you feel unwell.
8. P363: Wash contaminated clothing before reuse.
9. P501: Dispose of contents / container according to local, state and national regulations.
10. Do not use beyond expiration date.

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Diagnostic Grifols, S.A. Passeig Fluvial 24, 08150 Parets del Vallès (Barcelona), Ispanija