Neuro Plating System

P-NEU-IFU-N01G(Rev.12.1)2023.08

CE1639

ENGLISH

1. Indications for Use

Neuro Plating system is intended for use in selective trauma of the cranial skeleton, cranial surgery and reconstructive procedures.

2. Patient population

- Age: Adults and adolescents (age 12 and higher)
 - * People with other conditions depend on the surgeon's decision
- Gender: Not relevant
- Anatomy, physiology, possibly other aspects: Selective trauma of the cranial skeleton, cranial surgery and reconstructive procedure
- Possible other aspects: Pregnant woman is according to the surgeon's decision

3. Contraindication

- Patient who has active infection or potential infection
- Patient who has bad bone quality or has disease which may cause bone weakness.
- Patient who has problem with an immune system.
- Patient who has an allergy against titanium and titanium alloy (if an allergy is suspected, an allergy test is recommended.)
- Patient who has a blood poisoning.
- Patients who can't take care of them selves after the surgery (mentally or physically)

4. Side Effects

- Heavy physical exercise or pressure and incomplete treatment can damage an implant.
- Implant loosening or moving
- An allergy against titanium or titanium alloy
- Pain caused by implant
- Neural damage caused by surgical surgery
- Bone necrosis or bone absorption or incomplete treatment
- Primary and secondary infections.

5. Warning or Caution

- Do not use the products other than its purpose mentioned on IFU.
- Do not use it for those patients who suffer from a low healing ability or long healing period. Ex: Too old, mental illness, alcohol addicted patients.
- If an allergy against titanium or titanium alloy is suspected, an appropriate test should be done before using the products on the patient.
- Always follow appropriate precautions.
- Do not, under any circumstances, combine implant components from different manufacturers.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- The implant components were tested and approved in combination with Osteonic components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Use the most suitable size product for a patient.
- Its function can't be carried out properly if it is used in a wrong or improper purpose.
- Warn patients that side effect might be occur if the patient doesn't follow surgeon's treatment instruction after surgery.
- There are possibilities of implant loosening or implant moving caused by improper transplant. Therefore, informing patient about possibility of implant moving or loosening and possible side effect should be done.
- Excessive and repetitive bending may cause fracture.
- Dispose all risky products or wrongly handled products.
- Do not reuse surgically contaminated implants.
- Intraoperatively contaminated implants must be disposed.
- All the instruments used together with the implants should be reused after reprocessing.

6. MRI Safety Information

Non-clinical testing has demonstrated the Neuro Plating System is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial magnetic field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1 W/kg and head SAR of 1 W/kg for landmarks above the shoulder
- Normal Operating Mode for gradient output.

Under the scan conditions defined above, the Neuro Plating System is expected to produce a maximum temperature rise of less than 5.00°C after 7 minutes of continuous scanning. A cool down period of 5 minutes is needed after each 7 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 2 mm from the Neuro Plating System when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Neuro Plating System Instruments are MR Unsafe. Do not bring the Neuro Plating System Instruments into the MR environment.

7. Direction for use

1) Pre-operative preparations

- (1) Check if there is any damage of the product or its package.
- (2) Before surgery, select the screws and plates after checking patient's bone condition and any problem expected after surgery.
- (3) The operating surgeons should be properly informed about the instruments and their usage, clinical indications and contraindications.
- (4) Check if there is any biochemical and biological factors which take a bad effect on surgery.
- (5) Carefully read IFU before using the products.
- (6) Surgeons/physician should clearly inform the patient about risks related to surgery.
- (7) Since the product is supplied in a non-sterile condition, it must be sterilized before use.

2) Implantation of plates & screws

(1) Plate & Screw Selection

Select appropriate screws and plates according to the treatment area and the patient's bony anatomy.

(2) Plate handling

- Be careful not to deform the screw hole when bending the plate.
- Use plate bender for bending of plates.
- Use plate with different hole sizes for fixation of StarPore only (larger hole: StarPore, smaller hole: bone).
- Use forceps to place a plate at the bone fixation site.

(3) Screw handling

- To pick up screw from magazine, use manual or Power hand driver in combination with according shaft.

- Insert the screw at right angles to the bone surface plane and the implant.
- For use of Self-Tapping screws, drill a pilot hole for easy screw insertion if necessary.
- Use 1,95mm screw for fixation of StarPore bone flap only.

(4) Drilling pilot holes.

- Making a pilot hole before screw insertion can reduce bone loss and the risk of screw fracture.
- Use a minimum RPM to prevent bone necrosis during Pilot Drill.

(5) Plate & screw fixation.

- Insert screws by using screw driver and shaft.
- Shaft can be removed from Screw Driver.
- Use a shaft that corresponds to the size of the cross groove on the screw head.
- In case of high density bone, pilot drill needs to be done before screw insertion.

(6) Screw Removal

- Fix the screwdriver to the screw head firmly, then slowly rotate the screw opposite to the direction of screw insertion. Then, remove the screw completely from the bone.

8. Material

Screw: Ti-6Al-4V ELI (ASTM F136)

Plate: Pure Titanium (ASTM F67)

9. Cleaning and sterilization recommendations for use of non-sterile medical devices

(1) Inadequate cleaning detergent

- Strongly acidic and basic (sulfuric acid, nitric acid, hydrochloric acid, etc.) cleaning solutions are inappropriate. Washing the products for a long time in high temperature is also improper.

(2) Caution

- Use of abrasive products or equipment (sandpaper, metallic brush, etc.) is prohibited. After the cleaning, the performance and operation status of the instrument, and the presence of foreign substances in the Implant product should be verified. To do this, the cleaning facilities and cleaning methods used at each hospital should be validated

(3) Dry

- Surgical instruments and products are thoroughly dried before sterilization.

Validated mechanical alkaline cleaning and thermal disinfecting

Phase	Step	T [°C/°F]	t [min]	Water quality	Remark
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: Dr. Weigert neodischer MediClean forte

(4) Sterilization

- The product is packaged in non-sterile condition and delivered.
- Thus, the product must be sterilized before use.
- Remove all packaging material before sterilization.
- Use a sterilization and storage tray for sterilization and pre-operative storage. Steam autoclaving is recommended as a sterilization method.

* Validated sterilization process

- Steam sterilization using fractional vacuum process

- Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665

- Sterilization using fractional vacuum process at 132 °C/holding time 4min(Wrapped)

* Medical devices used with this device must be declared, certified and approved.

- Other sterilization methods may be used, but the effectiveness of sterilization should be verified before use. Validation according to the sterilization method should be performed at the hospital where sterilization is performed and sterilization time and the recommended sterilization temperature should be regularly checked.

- If paper filter is used, a new filter paper should be used for each sterilization. If there is any water left in the sterilization container or product after sterilization, it should be sterilized again.

10. Storage

Store at room temperature

11. Notice

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of Member State in which the user or patient is established.

12. Disposal

Hospitals should follow hospital medical waste disposal protocol.

Contaminated units should be decontaminated before they are discarded.

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OSTEONIC Co., Ltd