

GC-MS QP and **TQ** Series

Smart Database

**Smart Aroma Database**

1.8

**Instruction Manual**

For GCMSsolution

1.8

Read this manual thoroughly before you use the product.  
Keep this manual for future reference.

This page is intentionally left blank.

# Introduction

**Read this Instruction Manual thoroughly before using the product.**

Thank you for purchasing this product. This manual describes the installation, operation, usage cautions and accessories for this product. Read this manual thoroughly before using the product and operate the product in accordance with the instructions in this manual.

Also, keep this manual for future reference.

## Notice

- All rights are reserved, including those to reproduce this publication or parts thereof in any form without permission in writing from Shimadzu Corporation.
- Information in this manual is subject to change without notice and does not represent a commitment on the part of the vendor.
- Any errors or omissions which may have occurred in this manual despite the utmost care taken in its production will be corrected as soon as possible, although not necessarily immediately after detection.
- Shimadzu Corporation is not responsible for errors or injuries resulting from following the instructions in this document.
- Replacement parts for this product will be available for a period of seven (7) years after the product is discontinued. If the problem persists, or the symptoms are not covered in the troubleshooting section, contact your Shimadzu representative. Note, however, that the availability of parts not manufactured by Shimadzu shall be determined by the relevant manufacturers.
- Microsoft, Windows and Excel are registered trademarks of Microsoft Corporation in the United States and/or other countries. Other company names and product names mentioned in this manual are trademarks or registered trademarks of their respective companies. The TM and ® symbols are omitted in this manual.
- Microsoft® Windows® 10 Operating System is referred to as “Windows 10”.
- Microsoft® Office Excel® 2016 is referred to as “Excel”.

© 2024 Shimadzu Corporation. All rights reserved.

# Software License Agreement

PLEASE READ THIS AGREEMENT CAREFULLY BEFORE USING THE SOFTWARE.

SHIMADZU Corporation ("SMZ") is willing to license the SMZ software provided herein, together with accompanying documentation (collectively "SOFTWARE") to you only upon the condition that you accept all of the terms and condition contained in this License Agreement. By using the SOFTWARE, you agree to be bound by the terms of this Agreement. If you do not agree to all these terms of this Agreement, promptly return the unused SOFTWARE to the party (either SMZ or its reseller) from whom you acquired it to receive a refund of the amount you paid.

1. LICENSE.

SMZ grants you a non-exclusive and nontransferable license to use the SOFTWARE subject to the following terms and conditions.

2. LIMITATION OF USE.

Except as specifically authorized by SMZ, you may NOT:

- a. Use the SOFTWARE, or permit the SOFTWARE to be used, on more than one computer at any one time.
- b. Copy the SOFTWARE except one (1) archival copy of the SOFTWARE.
- c. Modify, reverse engineer, decompile, disassemble, or create derivative works based upon the SOFTWARE.
- d. Transfer, rent, lease or grant any rights in the SOFTWARE in any form to anyone else.

3. TITLE AND OWNERSHIP.

This license is not for sale and it may not be assigned or sublicensed to anyone else. Title and all associated intellectual property rights to the SOFTWARE shall remain in SMZ and/or its licensor.

4. UPGRADES.

You are entitled to receive all official software upgrades for the SOFTWARE that SMZ will release as deemed necessary by SMZ. An upgrade means certain supplemental program modules and/or information for bug fixing and/or updates to the defects and/or failures of the SOFTWARE that are acknowledged or confirmed by SMZ. An upgrade excludes hardware, network, consulting services, third party products, operation and general computer system maintenance. All supplemental program module for upgrades and enhancements furnished to you shall be deemed to be part of the SOFTWARE and subject to the terms and conditions set forth in this Agreement.

5. LIMITED WARRANTY.

SMZ warrants that for a period of one (1) year from the date of purchase, as evidenced by a copy of the receipt, the media on which SOFTWARE is furnished will be free of defects in materials and workmanship under normal use.

Except for the foregoing, SOFTWARE is provided AS IS. Your exclusive remedy and the entire liability of SMZ and its suppliers under this limited warranty will be, at SMZ's option, repair, replacement, or refund of the Software if reported (or, upon request, returned) to the party supplying the SOFTWARE to you. In no event does SMZ warrant that the Software is error free or that you will be able to operate the SOFTWARE without problems or interruptions.

EXCEPT AS SPECIFIED IN THIS WARRANTY, ALL EXPRESS OR IMPLIED CONDITIONS, REPRESENTATIONS, AND WARRANTIES INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, NONINFRINGEMENT, ARISING FROM A COURSE OF DEALING, USAGE, OR TRADE PRACTICE, ARE HEREBY EXCLUDED TO THE EXTENT ALLOWED BY APPLICABLE LAW.

6. LIMITATION OF LIABILITY.

IN NO EVENT WILL SMZ BE LIABLE FOR ANY LOST REVENUE, PROFIT OR DATA, OR FOR SPECIAL, INDIRECT, CONSEQUENTIAL, INCIDENTAL OR PUNITIVE DAMAGES, HOWEVER CAUSED REGARDLESS OF THE THEORY OF LIABILITY, ARISING OUT OF OR RELATED TO THE USE OF OR INABILITY TO USE SOFTWARE, EVEN IF SMZ HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGE. IN NO EVENT WILL SMZ'S LIABILITY TO YOU, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE), OR OTHERWISE, EXCEED THE AMOUNT YOU PAID FOR SOFTWARE.

7. TERMINATION.

This License is effective until terminated. You may terminate this License at any time by destroying all copies of SOFTWARE including any documentation. This License will terminate immediately without notice from SMZ if you fail to comply with any provision of this License. Upon termination, you must destroy all copies of SOFTWARE.

8. GENERAL

- a. This Agreement is the entire agreement. If any provision of this agreement is held invalid, the remainder of this agreement shall continue in full force and effect.
- b. This Agreement shall be construed and governed in accordance with the laws of Japan, excluding its conflict of law rules.
- c. The exclusive jurisdiction for any disputes arising out of or in connection with this Agreement shall be Kyoto District Court of Japan.
- d. The invalidity or unenforceability of any provision of this Agreement shall not affect the validity or enforceability of any other provision.

## Indications Used in This Manual

Throughout the text in this manual, warnings and other information essential when using this unit, such as cautionary or prohibited items, appear classified as per the following:

Indication	Meaning
 <b>WARNING</b>	Indicates a potentially hazardous situation which, if not avoided, could result in serious injury or possibly death.
 <b>CAUTION</b>	Indicates a potentially hazardous situation which, if not avoided, may result in minor to moderate injury or equipment damage.
 <b>NOTE</b>	Emphasizes additional information that is provided to ensure the proper use of this product.
 <b>Hint</b>	Indicates information provided to improve product performance.
 <b>Reference</b>	Indicates the location of related reference information.

### **WARNING**

When the customer uses the CD-ROM

This is a CD-ROM disk. Do not play this on an audio CD player, as the high volume may damage your hearing or the audio speakers.

# Product Warranty

Shimadzu provides the following warranty for this product.

1. **Warranty Period** Please contact your Shimadzu representative for information about the period of this warranty.
2. **Description:** If a product/part failure occurs for reasons attributable to Shimadzu during the warranty period, Shimadzu will repair or replace the product/part free of charge. However, in the case of products which are usually available on the market only for a short time, such as personal computers and their peripherals/parts, Shimadzu may not be able to provide identical replacement products.
3. **Limitation of Liability:**
  - 1) In no event will Shimadzu be liable for any lost revenue, profit or data, or for special, indirect, consequential, incidental or punitive damages, however caused regardless of the theory of liability, arising out of or related to the use of or inability to use the product, even if Shimadzu has been advised of the possibility of such damage.
  - 2) In no event will Shimadzu's liability to you, whether in contract, tort (including negligence), or otherwise, exceed the amount you paid for the product.
4. **Exceptions** Failures caused by the following are excluded from the warranty, even if they occur during the warranty period.
  - 1) Improper product handling
  - 2) Repairs or modifications performed by parties other than Shimadzu or Shimadzu designated companies
  - 3) Product use in combination with hardware or software other than that designated by Shimadzu
  - 4) Computer viruses leading to device failures and damage to data and software, including the product's basic software
  - 5) Power failures, including power outages and sudden voltage drops, leading to device failures and damage to data and software, including the product's basic software
  - 6) Turning OFF the product without following the proper shutdown procedure leading to device failures and damage to data and software, including the product's basic software
  - 7) Reasons unrelated to the product itself
  - 8) Product use in harsh environments, such as those subject to high temperatures or humidity levels, corrosive gases, or strong vibrations
  - 9) Fires, earthquakes, or any other act of nature, contamination by radioactive or hazardous substances, or any other force majeure event, including wars, riots, and crimes
  - 10) Product movement or transportation after installation
  - 11) Consumables and equivalent items  
Recording media such as are considered consumable items.

- If there is a document such as a warranty provided with the product, or there is a separate contract agreed upon that includes warranty conditions, the provisions of those documents shall apply.
- Warranty periods for products with special specifications and systems are provided separately.

# Contents

<b>1</b>	<b>Smart Aroma Database</b> .....	<b>1</b>
1.1	Overview.....	1
1.2	Operating Environment .....	2
1.3	Supported Models .....	2
1.4	Product Contents.....	3
1.4.1	Product Contents .....	3
1.4.2	Method File Details .....	4
1.4.3	Contents of Smart Database .....	6
1.4.4	Library Contents.....	6
1.4.5	Instruction Manual Included.....	7
<b>2</b>	<b>Operating Procedure</b> .....	<b>8</b>
2.1	Analysis Process Flow Using a Scan Method File .....	8
2.1.1	Creating a Folder for Analysis.....	8
2.1.2	Setting Instrument Parameters.....	9
2.1.3	Measuring and Identifying n-Alkanes .....	11
2.1.4	Adjusting Retention Times and Analyzing Samples .....	13
2.1.5	Checking Measurement Results.....	16
2.2	Library Operations .....	22
2.2.1	Library Search.....	22
2.2.2	Information Registered in the Library .....	23
2.3	Using a Smart Database to Create Method Files.....	23
2.3.1	Preparing for Analysis.....	24
2.3.2	Measuring and Identifying n-Alkanes .....	25
2.3.3	Specifying the Smart Database .....	25
2.3.4	Creating Method Files.....	28
2.3.4.1	Setting Target Compounds .....	28
2.3.4.2	Selecting Ions .....	29
2.3.4.3	Creating an MS Table .....	30
2.3.5	Editing Databases.....	36
<b>3</b>	<b>Precautions</b> .....	<b>38</b>
3.1	Changing Analytical Conditions or Using a Column that was Not Specified.....	38
3.2	Functionality for Protecting the Detector .....	38
<b>Appendix 1</b>	<b>Analytical Conditions</b> .....	<b>39</b>
<b>Appendix 2</b>	<b>List of Customer-Supplied Items</b> .....	<b>41</b>
<b>Appendix 3</b>	<b>Calculating Semi-Quantitative Results for Registered Compounds (Scan)</b> .....	<b>42</b>
Appendix 3.1	Measuring and Analyzing Standard Samples .....	43
Appendix 3.2	Creating Method Files .....	44

## Contents

---

Appendix 3.3	Analyzing Samples.....	47
Appendix 3.4	Checking Measurement Results.....	50
<b>Appendix 4</b>	<b>Using the Database with AOC-6000 Series Autosampler .....</b>	<b>51</b>
Appendix 4.1	Instrument Configuration Settings.....	51
Appendix 4.2	Analyzing n-Alkanes.....	51
Appendix 4.3	Analyzing Samples.....	52
Appendix 4.4	Examples of Method Parameter Settings .....	52
<b>Appendix 5</b>	<b>Using the Database with TD-30 Series Autosampler .....</b>	<b>62</b>
Appendix 5.1	Instrument Configuration Settings.....	62
Appendix 5.2	Analyzing n-Alkanes.....	62
Appendix 5.3	Analyzing Samples.....	63
Appendix 5.4	Examples of Method Parameter Settings .....	64
<b>Appendix 6</b>	<b>Using the Database with HS-20 Series Autosampler .....</b>	<b>66</b>
Appendix 6.1	Instrument Configuration Settings.....	66
Appendix 6.2	Analyzing n-Alkanes.....	66
Appendix 6.3	Analyzing Samples.....	67
Appendix 6.4	Examples of Method Parameter Settings .....	69
<b>Appendix 7</b>	<b>Using the Database with OPTIC-4 Series Inlet.....</b>	<b>71</b>
Appendix 7.1	Instrument Configuration Settings.....	71
Appendix 7.2	Analyzing n-Alkanes.....	71
Appendix 7.3	Analyzing Samples.....	72
Appendix 7.4	Examples of Method Parameter Settings .....	72
<b>Appendix 8</b>	<b>Connecting an Odor Sniffer .....</b>	<b>75</b>
Appendix 8.1	Flow Restrictor Settings for Connecting a Sniffer Unit.....	75
Appendix 8.2	Methods for when a Sniffer Unit is Connected.....	76
<b>Appendix 9</b>	<b>Adding Target Compounds.....</b>	<b>77</b>
Appendix 9.1	How to Add a Compound to a Method File (Scan).....	77
Appendix 9.2	How to Add a Compound to a Database File (Scan/SIM).....	79
Appendix 9.3	How to Add a Compound to a Database File (MRM).....	80
<b>Appendix 10</b>	<b>Retention Time Adjustment Using a Standard Sample.....</b>	<b>81</b>
<b>Appendix 11</b>	<b>Outputting Files for Statistical Data Analysis Software.....</b>	<b>85</b>
<b>Appendix 12</b>	<b>Installing the Database.....</b>	<b>89</b>

# 1 Smart Aroma Database

## 1.1 Overview

Optimizing columns and analytical conditions and gathering GC/MS information for compounds involved in analyzing the odors/aromas of foods, beverages, chemical products, or other samples is an extremely tedious and time-consuming process. Using non-targeted analysis to comprehensively analyze samples, without specifying specific compounds to analyze, can cause problems with tedious data analysis and inaccurate compound identification results.

1.8 **The Smart Aroma Database offers method files and database files with a wealth of information registered about columns, analytical conditions, and compounds optimized for a wide variety of odor/aroma components.**

1.8 Method files for scan mode analysis include retention index values determined from n-alkanes, so that the **automatic adjustment of retention time (AART) function can be used with a standard sample to easily adjust retention times for registered compounds.** That means target compounds can be efficiently identified based on the registered retention time, mass chromatogram, and mass spectrum information.

As a “Smart Database,” **retention times adjusted using the AART retention time adjustment function can be used to create methods for high-sensitivity targeted SIM or MRM analysis.**

Furthermore, database files are saved in Excel format to enable editing or adding compounds.



### NOTE

This instruction manual has been written assuming the reader already has a basic knowledge about using GC/MS, GC-MS/MS, and GCMSsolution products. In many instances, explanations include GCMSsolution-specific nomenclature and terminology. If using a GC/MS or GC-MS/MS system or the GCMSsolution software for the first time or if such terminology is not clear, then also refer to the corresponding GC/MS or GC-MS/MS instruction manual, the GCMS Operation Guide, and the help menu in GCMSsolution.

### 1.2 Operating Environment

Use the database in the following operating environment.

Operating System	Microsoft Windows 11 Professional 64-bit version Microsoft Windows 10 Professional 64-bit version Microsoft Windows 7 Professional SP1 32 or 64-bit version
Excel	Microsoft Excel2021 (32 or 64-bit version) Microsoft Excel2019 (32 or 64-bit version) Microsoft Excel2016 (32-bit version)
Workstation Software	GCMSsolution Ver. 4.53 SP1 or later

### 1.3 Supported Models

This product can be used with the following systems.

- GCMS-QP series: GCMS-QP2020 series or GCMS-QP2010 SE
- GCMS-TQ series: GCMS-TQ series

This product can be used with the following columns.

- SH-I-5Sil MS (30m x 0.25mm I.D. df=0.25  $\mu$ m) P/N: 221-75954-30  
Contact: Shimadzu Corporation  
This slightly polar column is suitable for targeting mainly compounds with relatively low polarity, such as terpenes.
- SH-PolarWax (60m x 0.25mm I.D. df=0.25  $\mu$ m) P/N: 227-36247-02  
Contact: Shimadzu Corporation  
This polar column is suitable for targeting mainly compounds with high polarity, such as organic acids. Due to the GC conditions used that prioritize separation, this column is ideal for comprehensive analysis.
- InertCap PureWAX (30m x 0.25mm I.D. df=0.25  $\mu$ m) P/N: 1010-68142  
Contact: GL Sciences Inc.  
This polar column is suitable for targeting mainly compounds with high polarity. It is the same column as used in GC/MS Off-Flavor Analyzer, so it can be used for both systems.



#### NOTE

In order to use the retention index values registered in the database, use a column with the specified liquid phase type, length, internal diameter, and film thickness that is made by the specified manufacturer.



### NOTE

Due to differences between column production lots of Wax column, retention index values can vary. If that occurs, use a standard sample to adjust the retention times based on the retention times predicted from the registered retention index values. (Refer to Appendix 10 Retention Time Adjustment Using a Standard Sample.)

In addition to AOC-20 and AOC-30 series autosamplers typically used for liquid injection, the database also supports the following pretreatment systems. For more information about the procedures for respective pretreatment systems, refer to Appendices 4 to 7.

- AOC-6000 series
- HS-20 series
- TD-30 series
- OPTIC-4

## 1.4 Product Contents

### 1.4.1 Product Contents

This product includes a CD-ROM. If it is missing, contact a Shimadzu representative.

P/N	Name	Quantity
225-50791-91	CD-ROM	1

The CD-ROM includes an installer program and a PDF version of the instruction manual.



### NOTE

Though the product does not include a column, be sure to use the specified type.

## 1. Smart Aroma Database

---

### 1.4.2 Method File Details

Executing the installer program (setup.exe) included with the product installs the following method files.

During installation, select the type of instrument (QP or TQ) so that files for the given instrument type are installed. (In this case, "C:" refers to the drive where GCMSsolution is installed.)

Files for the applicable columns are installed in the following folders for respective columns.

For SH-I-5Sil MS:

"C:\GCMSsolution\SmartDatabase\Aroma\SH-I-5Sil\_MS" folder

For SH-PolarWax:

"C:\GCMSsolution\SmartDatabase\Aroma\SH-PolarWax" folder

For InertCap PureWax:

"C:\GCMSsolution\SmartDatabase\Aroma\InertCap\_Pure-Wax" folder

#### QP Methods

Column	File Name	Description
SH-I-5Sil MS	Aroma_QP_SH-I-5SilMS_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_QP_SH-I-5SilMS_AART.qgm	Method file for retention time adjustment
	Aroma_QP_SH-I-5SilMS_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_QP_SH-I-5SilMS_Template.qgm	Template for creating SIM method file from Smart Database
SH-PolarWax	Aroma_QP_SH-PolarWax_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_QP_SH-PolarWax_AART.qgm	Method file for retention time adjustment
	Aroma_QP_SH-PolarWax_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_QP_SH-PolarWax_Template.qgm	Template for creating SIM method file from Smart Database
InertCap PureWAX	Aroma_QP_IC-WAX_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_QP_IC-WAX_AART.qgm	Method file for retention time adjustment
	Aroma_QP_IC-WAX_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_QP_IC-WAX_Template.qgm	Template for creating SIM method file from Smart Database

## 1. Smart Aroma Database

### TQ Methods

Column	File Name	Description
SH-I-5SiI MS	Aroma_TQ_SH-I-5SiI MS_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_TQ_SH-I-5SiI MS_AART.qgm	Method file for retention time adjustment
	Aroma_TQ_SH-I-5SiI MS_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_TQ_SH-I-5SiI MS_Template.qgm	Template for creating SIM or MRM method file from Smart Database
SH-PolarWax	Aroma_TQ_SH-PolarWax_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_TQ_SH-PolarWax_AART.qgm	Method file for retention time adjustment
	Aroma_TQ_SH-PolarWax_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_TQ_SH-PolarWax_Template.qgm	Template for creating SIM or MRM method file from Smart Database
InertCap PureWAX	Aroma_TQ_IC-WAX_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_TQ_IC-WAX_AART.qgm	Method file for retention time adjustment
	Aroma_TQ_IC-WAX_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_TQ_IC-WAX_Template.qgm	Template for creating SIM or MRM method file from Smart Database

 **Hint**

For systems with sample flow split for detection by a sniffer unit and MS, use the method files in the folders indicated below. For more details, refer to Appendix 8 Connecting to an Odor Sniffer.

“C:\GCMSsolution\SmartDatabase\Aroma\<<column name>\Sniff” folder

## 1. Smart Aroma Database

---

### 1.4.3 Contents of Smart Database

Executing the installer program (setup.exe) included with the product installs the following Smart Database files.

During installation, select the type of instrument (QP or TQ) so that files for the given instrument type are installed.

File location: "C:\GCMSsolution\SmartDatabase"

#### QP Methods

Column	File Name	Description
SH-I-5Sil MS	SmartDatabase_Aroma_QP_SH-I-5SilMS.xlsm	Smart Database file for SH-I-5Sil MS
SH-PolarWax	SmartDatabase_Aroma_QP_SH-PolarWax.xlsm	Smart Database file for SH-PolarWax
InertCap PureWAX	SmartDatabase_Aroma_QP_IC-WAX.xlsm	Smart Database file for InertCap PureWAX

#### TQ Methods

Column	File Name	Description
SH-I-5Sil MS	SmartDatabase_Aroma_TQ_SH-I-5SilMS.xlsm	Smart Database file for SH-I-5Sil MS
SH-PolarWax	SmartDatabase_Aroma_TQ_SH-PolarWax.xlsm	Smart Database file for SH-PolarWax
InertCap PureWAX	SmartDatabase_Aroma_TQ_IC-WAX.xlsm	Smart Database file for InertCap PureWAX

### 1.4.4 Library Contents

Executing the installer program (setup.exe) included with the product installs the following library files.

File location: "C:\GCMSsolution\library"

Column	File Name	Description
SH-I-5Sil MS	AROMA_5MS	Library file for SH-I-5Sil MS
SH-PolarWax	AROMA_POLARWAX	Library file for SH-PolarWax
InertCap PureWAX	AROMA_ICWAX	Library file for InertCap PureWAX

## 1. Smart Aroma Database

---

Each library comprises the following seven files. The following example is for the “AROMA.LIB\_5MS” library.

File Name	Description
AROMA_5MS.LIB	Mass spectral library file
AROMA_5MS.FOM	Includes composition formula information for registered components.
AROMA_5MS.SPC	Includes spectral information for registered components.
AROMA_5MS.NAM	Includes name information for registered components.
AROMA_5MS.COM	Includes comment information.
AROMA_5MS.FLG	Includes information about whether or not to target the component in searches.
AROMA_5MS.ADD	Includes links to structural information for registered components

- In Shimadzu GC/MS software, all libraries are indicated in terms of the “\*.LIB” file. Operating the software does not require an awareness of the other six files. However, note that the library will not function properly unless all seven files are present in the same library folder.
- The file structure for these libraries is configured as a private library. The file structure is different than for public libraries.
- These libraries include information about structural formulas. Structural formula files are included in the “Aroma\_Structure” folder. The structural formula will not function properly unless these folders are present in the following folder.  
File location: “C:\GCMSsolution\library”

### 1.4.5 Instruction Manual Included

Executing the installer program (setup.exe) included with the product installs the following PDF file of the instruction manual.


The instruction manual file location: C:\GCMSsolution\Manual\SmartDatabase\Aroma

File Name	Description
Smart Aroma Database_Manual.pdf	Smart Aroma Database Instruction Manual

## 2 Operating Procedure


### 2.1 Analysis Process Flow Using a Scan Method File

Follow the steps below to analyze samples using a scan method file included in the database. The procedure is basically the same for both QP and TQ type systems, but some of the steps could vary depending on the pretreatment system used. For more information about analysis steps using each type of pretreatment system, refer to the corresponding chapter in appendices.



#### 2.1.1 Creating a Folder for Analysis


- 1) Create a folder for saving acquired data files.  
(Example: "C:\GCMSsolution\Data\Project1")
- 2) Copy the method file that is stored in the "C:\GCMSsolution\SmartDatabase\Aroma\<column name>" folder and paste it into the folder that was created in step 1.



## 2. Operating Procedure

### 2.1.2 Setting Instrument Parameters

- 1) Specify system configuration settings as follows.  
In the [Modules Used for Analysis] field, specify the pretreatment unit to use for analysis line 1.



#### NOTE


Instrument configuration settings differ depending on the pretreatment unit. Refer to the instruction manual for the pretreatment unit to be used.

- 2) Select the column based on the method file used for analysis and install it in the instrument. Refer to 1.4.2 Method File Details regarding the relationship between methods and columns.
- 3) In the [System Configuration] window, double-click the [Column] icon for the unit to be used for analysis. Enter information for the column to be used in column property settings.



### NOTE

If the confirmation message "[1317] The hardware configuration for this method is different from the current instrument configuration. The measurement condition in the method file is according to the current instrument configuration." is displayed when a method file is loaded in the [GCMS Real Time Analysis]-[Acquisition] window, click [OK].



### 2.1.3 Measuring and Identifying n-Alkanes

Measure an n-alkane standard mixture solution so that retention times can be adjusted for each compound registered in the method file. Adjust the registered retention times for measurement target compounds based on the n-alkane identification results and retention index values for measurement target compounds (AART function).

### Hint

The following n-alkane standard mixture sample is recommended.

C7 - C30 Saturated Alkanes :

SIGMA-ALDRICH Cat#: 49451-U

(1000 µg/mL each component in hexane)

### CAUTION

Dilute standard samples under a fume hood or other exhaust system.

## 2. Operating Procedure

### <Preparing n-Alkane Standard Mixture Samples>

Prepare the standard mixture sample by diluting C7 - C30 Saturated Alkanes solution (1000 µg/mL) with hexane to a concentration of 50 µg/mL.

#### Preparation Example:

Use a microsyringe to place 500 µL of C7 - C30 Saturated Alkanes solution (1000 µg/mL) in a 10 mL volumetric flask and then fill the flask to volume with hexane. Keep the prepared solution stored in a refrigerator or freezer until ready for use.

- 1) Analyze the n-alkane solution (50 µg/mL) using the method file for adjusting retention times.



### NOTE

When AOC-30/20 series, AOC-6000 series or OPTIC-4 is used, measure n-alkane is measured by liquid injection (1 µL injection). For the TD-30 series, add n-alkane to the TD tube, and for the HS-20 series, add n-alkane to the HS vial for analysis. Refer to "Appendix 3" to "Appendix 7" for the settings of each auto sampler.

- 2) Identify the n-alkanes.



The screenshot shows the GCMS Postrun Analysis software interface. The main window displays chromatograms (TIC & MIC) and a peak table. The peak table is highlighted with a red box and labeled '4'. The table contains the following data:

Peak	Name	Ret. Time	m/z	Area	Height	SI
11	C8	2.986	129.20	625029	489593	91
12	C10	4.960	143.20	617955	419796	97
13	C11	4.854	156.20	617747	284657	84
14	C12	4.451	170.20	641763	207930	84
15	C13	3.957	184.20	743646	288318	95
16	C14	3.560	198.20	672271	230205	86
17	C15	3.200	212.20	598349	247358	97
18	C16	3.115	226.20	665691	231744	97
19	C17	3.024	240.20	529101	219290	97
20	C18	3.087	254.20	524610	200020	97
21	C19	3.027	268.20	467965	192978	95
22	C20	3.143	282.20	442754	182177	97
23	C21	2.911	296.20	420207	164827	99
24	C22	2.925	310.20	383342	149681	99
25	C23	2.978	324.20	330481	133787	98
26	C24	2.916	338.20	308337	122260	98
27	C25	2.914	352.20	282296	111699	98
28	C26	2.965	366.20	258288	98998	98
29	C27	2.978	380.20	236296	83097	96
30	C28	2.978	394.20	216699	63933	95
31	C29	2.922	408.20	194201	52216	95

Callout 3 points to the 'Peak' column in the table. Callout 4 points to the 'SI' column. Callout 5 points to the 'Compound Table' in the left sidebar.

## 2. Operating Procedure

---


1. Start the [GCMS Postrun Analysis] program and then click the [Create Compound Table] icon on the [Postrun] assistant bar.
2. Open Data Explorer and double-click the icon for the n-alkane measurement data.
3. Click the [Results] tab below the compound table on the left.
4. Identify each component by checking the mass chromatogram (MC) in the quantitation view. If the components are not identified correctly, use manual peak integration.
5. Save the data file by overwriting.

### 2.1.4 Adjusting Retention Times and Analyzing Samples

- 1) Click the [AART] icon on the [Compound] assistant bar.



- 2) Select the method with the retention times to be adjusted and then click [Open].



- 3) Check the n-alkane identification results and then click [Next]

## 2. Operating Procedure

Automatic Adjustment of Retention Time [AART] 1/2

This wizard modifies the retention time in the compound table or MS instrument parameters of the method file by using the identified result and retention index. If there are unnecessary IDs, please turn off the checks.


Identified Results of Data File

ID#	Proc	Name	Ret. Time	Ret. Index
1	<input checked="" type="checkbox"/>	C9	2.086	900
2	<input checked="" type="checkbox"/>	C10	2.960	1000
3	<input checked="" type="checkbox"/>	C11	4.854	1100
4	<input checked="" type="checkbox"/>	C12	7.431	1200
5	<input checked="" type="checkbox"/>	C13	9.557	1300
6	<input checked="" type="checkbox"/>	C14	11.340	1400
7	<input checked="" type="checkbox"/>	C15	12.900	1500
8	<input checked="" type="checkbox"/>	C16	14.315	1600
9	<input checked="" type="checkbox"/>	C17	15.624	1700
10	<input checked="" type="checkbox"/>	C18	16.857	1800
11	<input checked="" type="checkbox"/>	C19	18.027	1900
12	<input checked="" type="checkbox"/>	C20	19.143	2000
13	<input checked="" type="checkbox"/>	C21	20.211	2100
14	<input checked="" type="checkbox"/>	C22	21.235	2200
15	<input checked="" type="checkbox"/>	C23	22.219	2300
16	<input checked="" type="checkbox"/>	C24	23.166	2400
17	<input checked="" type="checkbox"/>	C25	24.078	2500
18	<input checked="" type="checkbox"/>	C26	24.959	2600
19	<input checked="" type="checkbox"/>	C27	25.874	2700
20	<input checked="" type="checkbox"/>	C28	26.909	2800
21	<input checked="" type="checkbox"/>	C29	28.118	2900
22	<input checked="" type="checkbox"/>	C30	29.575	3000

< Back   **Next >**   Cancel   Help

## 2. Operating Procedure

- 4) Confirm that the retention times have been adjusted and then click [Finish]. Use the AART function to adjust retention times in the method file.  
Clear the [Modify the time of MS Instrument parameters] checkbox and then click [Finish].



### NOTE

For scan measurement methods, clear the [Modify the time of MS Instrument parameters] checkbox. Select the checkbox if using this procedure to adjust retention times in a method file created from the database for MRM or SIM measurements. MRM and SIM measurement tables are also automatically adjusted.

- 5) Then use the method file with adjusted retention times to measure samples.

 **NOTE**

If using an AOC-6000 series or HS-20 series unit, method parameters for the pretreatment system must be specified in the created method file. For a description of that procedure, refer to Appendix 4 Using the Database with an AOC-6000 Series Autosampler or Appendix 6 Using the Database with an HS-20 Series Headspace Sampler.

 **NOTE**

Method files registered in the Smart Aroma Database include a default split ratio of 5. Changing the split ratio setting significantly can cause predicted retention times to deviate from actual retention times, especially for compounds with smaller retention index values. If the values vary significantly, use a standard sample to adjust the retention times based on the retention times predicted from the registered retention index values. (Refer to Appendix 10 Retention Time Adjustment Using a Standard Sample.)

### 2.1.5 Checking Measurement Results

Measurement results can be analyzed using either the [GCMS Postrun Analysis] program or [Labsolutions Insight GCMS] program.

- GCMS Postrun Analysis: Enables library searches for compounds other than the target compounds, adding compounds, and so on.
- Labsolutions Insight GCMS: Enables data analysis for multiple analytes, which is especially convenient if intending to apply multivariate analysis or other statistical data analysis.

 **NOTE**

Automatic peak detection can sometimes incorrectly identify or not identify peaks in the following cases. Therefore, use the standard spectrum or reference ion ratio registered in the method file or a library similarity search to check measurement results.

- If the concentration of a target component is extremely low
- If an ion for a target component overlaps with a matrix peak originating from the sample
- If the instrument condition is poor due to contamination of the glass insert, analytical column, or ion source
- If it has been a long time since auto-tuning was performed

## 2. Operating Procedure




### NOTE

Retention index values in the database are acquired using liquid injection. Depending on the pretreatment method used, predicted retention times for compounds with smaller retention index values can tend to deviate from actual retention times, due to desorption efficiency or other factors.


### GCMS Postrun Analysis

1) Peak integration and peak identification processes are performed automatically for compounds listed in the compound table, based on the following criteria.

- 1 Mass chromatogram peak integration and identification
  - \*Allowable retention time range:  $\pm 0.2$  minutes (SH-I-5Sil MS) or 0.3 minutes (SH-PolarWax, InertCap Pure-WAX)
- 2 Similarity index (SI) from a reverse search of registered standard spectra (minimum similarity score: 70)



2) Information about the odor characteristics of each compound can be confirmed in the [Description] column of the compound table.




## 2. Operating Procedure

### Hint

The compound information registered in the method file includes both mass spectral and reference ion ratio information. Results can be checked more easily by displaying and comparing results to that information.

#### How to Display Standard Spectra


Right-click in the mass spectrum view and then select [Display Setting...] on the right-click menu and select the [Display Standard Spectrum] checkbox.



#### How to Display Reference Ion Ratios

In the compound table, select the compound for the reference ion ratio to be displayed. Then right-click on the mass chromatogram view, click [Display Setting...] on the right-click menu, and select the [Reference Ion Ratio] checkbox.

The reference ion ratio information is not used for identifying peaks by automatic peak integration, but it can be used as a reference for checking results.



The image shows two screenshots. On the left, a right-click context menu is displayed with 'Display Setting...' highlighted. On the right, the 'Quantitation View Display Setting' dialog box is shown with the 'Reference Ion Ratio' checkbox checked. Below the dialog box is a table of ion ratios.

Type	m/z	Intensity	Set %	Act. %	Ref. Band
Target	43.10	1502	100.00	100.00	---
Ref. Ion1	41.10	677	66.00	45.07	50
Ref. Ion2	39.10	206	13.04	13.72	50
Ref. Ion3					
Ref. Ion4					
Ref. Ion5					


## 2. Operating Procedure

### NOTE


Due to effects from matrix components contained in the sample, retention times for compounds eluted near such peaks can deviate from predicted retention times in some cases.

If the retention time deviation is large enough that retention times extend outside the window, change the window width by the following steps.

1. On the compound table [Param's] tab page, click [Edit].
2. Enter the desired window time settings in the compound table [Proc. From] and [Proc. To] fields. (The default setting is  $\pm 1$  minutes from the predicted retention time.)



3. Click [View] in the compound table. That changes the window range to the specified duration.



In that case, it is recommended that retention times be adjusted using the retention times for detected peaks, as described in Appendix 10 Retention Time Adjustment Using a Standard Sample.

## 2. Operating Procedure




### Hint

For instructions on adding unregistered compounds, refer to Appendix 9 Adding Target Compounds.

### Labsolutions Insight GCMS

- 1) Peak integration and peak identification processes are performed automatically for compounds listed in the compound table, based on the following criteria.
  - 1 Mass chromatogram peak integration and identification
    - \*Allowable retention time range:  $\pm 0.2$  minutes (SH-I-5Sil MS) or 0.3 minutes (SH-PolarWax, InertCap Pure-WAX)
  - 2 Similarity index (SI) from a reverse search of registered standard spectra (minimum similarity score: 70)



## 2. Operating Procedure




### Hint

The compound information registered in the method file includes both mass spectral and reference ion ratio information. Results can be checked more easily by displaying and comparing results to that information.


### How to Check Standard Spectra

Standard spectra can be checked by displaying the [Compound Details] tab page.




### How to Display Reference Ion Ratios

Reference ion ratios can be checked in the [Compound Results] field.



## 2. Operating Procedure

Information about the odor/aroma characteristics of each compound can be confirmed in the [Cmment] column in the [Compound Results] field.



### Hint

For instructions on outputting data as a list, so that statistical analysis software can be used to calculate data analysis results, refer to Appendix 11 Outputting Files for Statistical Data Analysis Software.

## 2.2 Library Operations

### 2.2.1 Library Search

The library can be searched by using it as a private library for Shimadzu GC/MS software. For more details, refer to the instruction manual for the main GC/MS system, which includes instructions for library search operations.

\* Normally included in the "C:\GCMSsolution\Manual" folder.




### NOTE

The library search functionality in the library is only compatible with the [GCMS Postrun Analysis] program.

### 2.2.2 Information Registered in the Library

The library includes the following mass spectral information.



- (1) Serial numbers registered in the library
- (2) Chemical Abstracts Service (CAS) registration numbers
- (3) Compound names and corresponding information about the qualitative characteristics of odors/aromas. Odor/aroma characteristics are indicated in square brackets.
- (4) Molecular formula
- (5) Retention index values calculated for each column using n-alkanes. Data is acquired by the method specified in GC/MS analytical conditions.
- (6) Monoisotopic mass values calculated as integer values
- (7) Structural formula

## 2.3 Using a Smart Database to Create Method Files


Method files for high-sensitivity targeted analysis using the SIM or MRM mode (or alternatively Scan/SIM or Scan/MRM modes) can be created using a Smart Database.

Widely-targeted scan methods can even be recreated from database files to add a newly registered compound to the method or if the compound table is accidentally edited, for example.

This section describes how to create method files from database files. Method files are created according to the following process flow.


## 2. Operating Procedure

---



### 2.3.1 Preparing for Analysis

- 1) Create a folder for saving respective files (method and data files).  
(Example: "C:\GCMSsolution\Data\Project1")  
As an example, the following describes steps for using the Smart Database file "SmartDatabase\_Aroma\_TQ\_IC-PureWAX.xlsm" to create a method file.
- 2) Copy the necessary n-alkane method file that is stored in the "C:\GCMSsolution\SmartDatabase\Aroma\<column name>" folder and paste it into the folder that was created above.



- 3) Attach the analytical column. Then specify configuration settings and perform auto-tuning according to the procedure indicated in 2.1.2 Setting Instrument Parameters.



### NOTE

In order to use the retention index values registered in the database, use a column with the specified liquid phase type, length, internal diameter, and film thickness made by the specified manufacturer.


### 2.3.2 Measuring and Identifying n-Alkanes

Measure an n-alkane standard mixture solution so retention times registered in the method file can be adjusted. Adjust the registered retention times for target compounds based on the n-alkane identification results and retention index values for the target compounds.

Measure and identify the n-alkanes in the n-alkane standard mixture solution as described in 2.1.2 Measuring an n-Alkane Standard Mixture Sample.

### 2.3.3 Specifying the Smart Database

- 1) Click the  (Smart MRM/SIM) icon on the [Compound] assistant bar.
- 2) Select the Smart Database file (Excel file) based on the column used and then click [Open].




## 2. Operating Procedure

### NOTE

If the following [SECURITY WARNING] message is displayed on the Excel message bar, click [Enable Content].




3) Display the [Database] sheet.



The image shows a software interface with a "Create Method File" dialog box and a table below it. The dialog box has fields for "Ret. Index for AART", "n-alkane data file", "Template Method File", and "Divide Method into". The table below has columns for "Serial #", "Type", "Acq. Mode", "ISTD Group", "Level 1 Conc (13)", "Method No.", "Compound Name (E)", "Ret. Index 1", "Ret. Index 2", and "Ret. Index 3". A red box highlights the "Database" tab in the table's footer.

Serial #	Type	Acq. Mode	ISTD Group	Level 1 Conc (13)	Method No.	Compound Name (E)	Ret. Index 1	Ret. Index 2	Ret. Index 3
1	Target	Scan			1	2-Methylfuran	850		
2	Target	Scan			1	Ethyl formate	857		
3	Target	Scan			1	Butanal	880		
4	Target	Scan			1	Ethyl acetate	886		
5	Target	Scan			1	Acetal	893		
6	Target	Scan			1	2-Butanone	893		
7	Target	Scan			1	Diethyl sulfide	895		
8	Target	Scan			1	2-Methylbutanal	903		
9	Target	Scan			1	3-Methylbutanal	908		
10	Target	Scan			1	Methyl isobutyrate	914		
11	Target	Scan			1	3-Methyl-2-butanone	916		
12	Target	Scan			1	2-Ethylfuran	935		
13	Target	Scan			1	Ethyl propanoate	946		
14	Target	Scan			1	Ethyl isobutyrate	954		
15	Target	Scan			1	Diacetyl	960		

4) Click [Lang.] to specify the language for displaying compound names and other text.



The image shows a "Lang." dialog box with two dropdown menus. The first is labeled "Compound Name Lang." and is set to "English". The second is labeled "Language" and is also set to "English". There are "OK" and "Cancel" buttons at the bottom. A red box highlights the "Lang." button in the background interface.

[Compound Name lang.]: The given database only displays compound names in English.  
[Language]: Selects the language for displaying information in Smart Database windows.

### Hint

The above settings are recorded in the Smart Database.

## 2. Operating Procedure

---

- 5) Confirm that the correct instrument type is selected on the [Instrument Type] pull-down menu.

The screenshot shows a software dialog box titled "Create Method File". At the top left is a button labeled "Create Method File". To its right is a dropdown menu labeled "Instrument Type" which is open, showing three options: "TQ Series" (highlighted in blue), "TQ Series", and "QP Series". An orange box highlights the "Instrument Type" dropdown. To the right of the dropdown is a "Lang." label. Below the dropdown is a "Parameter" section with several fields: "Ret. Index for AART" with a dropdown set to "Ret. Index 1"; "n-alkane data file" with a text box containing "C:\GCMSsolution\Data\Project1\Alkane.qgd" and a browse button (".."); "Template Method File" with an empty text box and a browse button (".."); and "Divide Method into" with a dropdown set to "1" and an "Advanced" button. At the bottom left is an "Import" button.

### NOTE

Note that an error message will be displayed if the instrument type specified in [Instrument Type] does not match the instrument type specified in configuration settings for data or method files specified in [n-alkane data file] and [Template Method file] fields.

### 2.3.4 Creating Method Files

#### 2.3.4.1 Setting Target Compounds

- 1) Confirm that the data file opened in 2.3.2 Measuring and Identifying n-Alkanes is loaded.

The screenshot shows the 'Create Method File' dialog box. At the top, there is a tab labeled 'Create Method File', an 'Instrument Type' dropdown set to 'TQ Series', and a 'Lang.' button. Below this is a 'Parameter' section with the following fields: 'Ret. Index for AART' (dropdown set to 'Ret. Index 1'), 'n-alkane data file' (text box containing 'C:\GCMSsolution\Data\Project1\Alkane.qgd' and highlighted with an orange box), 'Template Method File' (empty text box), and 'Divide Method into' (dropdown set to '1'). There are 'Import' and 'Advanced' buttons at the bottom.

- 2) Of the method files stored in the “C:\GCMSsolution\SmartDatabase\Aroma\<column name>” folder, copy the method file to be used as the template method file to the folder that was created above and specify it in [Template Method File]. The data acquisition and analytical conditions in the method file specified here will be applied when a method is created.

The screenshot shows the 'Create Method File' dialog box. The 'Template Method File' field is highlighted with an orange box and contains the path 'C:\GCMSsolution\SmartDatabase\Aroma\TQ-IC-Wax\_Template.qgm'. The other fields are the same as in the previous screenshot: 'Ret. Index for AART' (dropdown set to 'Ret. Index 1'), 'n-alkane data file' (text box containing 'C:\GCMSsolution\Data\Project1\Alkane.qgd'), and 'Divide Method into' (dropdown set to '1').

#### Hint

To create a method file for scan mode analysis, specify the scan method file (“Aroma\_QP/TQ\_<column name>\_Scan.qgm”) in the [Template Method File] field.  
To create a method file for SIM or MRM mode analysis, specify the template method file (“Aroma\_QP/TQ\_<column name>\_Template.qgm”) in the [Template Method File] field.

## 2. Operating Procedure

- 3) In the [Type] field, specify “Target” only for compounds being measured and use the [Delete] key to leave the field blank for other compounds.

Serial#	Type	Acq. Mode	ISTD Group	Level1 Conc (IS)	Method No.	Compound Name (E)	Ret. Index 1
1	Target	Scan			1	2-Methylfuran	850
2	Target	Scan			1	Ethyl formate	857
3	Target	Scan			1	Butanal	880
4	Target	Scan			1	Ethyl acetate	886
5	Target	Scan			1	Acetal	893
6	Target	Scan			1	2-Butanone	893
7	Target	Scan			1	Diethyl sulfide	895
8	Target	Scan			1	2-Methylbutanal	903
9	Target	Scan			1	3-Methylbutanal	908
10	Target	Scan			1	Methyl isobutyrate	914
11	Target	Scan			1	3-Methyl-2-butanone	916
12	Target	Scan			1	2-Ethylfuran	935
13	Target	Scan			1	Ethyl propanoate	946
14	Target	Scan			1	Ethyl isobutyrate	954
15	Target	Scan			1	Diacetyl	960
16	Target	Scan			1	2-Pentanone	961



### Hint

If retention times are already known for target components from analyzing a standard sample, the method file can be created without using the AART function by using the known retention times.

(Refer to Appendix 10 Retention Time Adjustment Using a Standard Sample.)

- 4) Specify the measurement mode (select Scan, SIM, or MRM) in the [Acq. Mode] column of the target compound table.

Serial#	Type	Acq. Mode	ISTD Group	Level1 Conc (IS)	Method No.	Compound Name (E)	Ret. Index 1
1	Target	Scan			1	2-Methylfuran	850
2	Target	MRM			1	Ethyl formate	857
3	Target	SIM			1	Butanal	880
4	Target	Scan			1	Ethyl acetate	886
5	Target	Scan			1	Acetal	893
6	Target	Scan			1	2-Butanone	893
7	Target	Scan			1	Diethyl sulfide	895
8	Target	Scan			1	2-Methylbutanal	903
9	Target	Scan			1	3-Methylbutanal	908
10	Target	Scan			1	Methyl isobutyrate	914
11	Target	Scan			1	3-Methyl-2-butanone	916
12	Target	Scan			1	2-Ethylfuran	935
13	Target	Scan			1	Ethyl propanoate	946
14	Target	Scan			1	Ethyl isobutyrate	954
15	Target	Scan			1	Diacetyl	960
16	Target	Scan			1	2-Pentanone	961

### 2.3.4.2 Selecting Ions

In default settings, three ion  $m/z$  values are specified for the Scan and SIM modes and two MRM transitions for one compound. To change MRM transitions, check the mass and ion ratio values before selecting the transition.

The ion is used as a target ion if “T” is specified in the [Type] column or as a reference ion if “Ref” is specified. One target ion and up to five reference ions can be specified.

## 2. Operating Procedure

### MRM Transition

Ion1				Ion2				Ion3			
Type	m/z	CE	Rat i	Type	m/z	CE	Rat i	Type	m/z	CE	Rat i
T	82.10>54.00	9	100.00	Ref.1	81.10>53.00	6	76.60	Ref.2	82.10>39.00	18	70.48
T	72.10>57.00	3	100.00	Ref.1	72.10>44.10	6	29.54	Ref.2	72.10>43.00	3	23.44
T	70.0>55.0	6	100.00	Ref.1	88.0>61.0	6	14.32	Ref.2	70.0>43.0	18	14.77
T	73.10>45.10	6	100.00	Ref.1	103.10>75.10	6	26.44	Ref.2	103.10>47.10	12	22.51
T	72.10>43.10	6	100.00	Ref.1	72.10>57.10	3	60.62	Ref.2	57.10>41.10	3	1.87
T	90.10>75.10	9	100.00	Ref.1	90.10>62.00	6	91.84	Ref.2	75.00>47.00	6	77.31

### m/z for SIM or Scan

Ion1			Ion2			Ion3		
Type	m/z	Rat i	Type	m/z	Rat i	Type	m/z	Rat i
T	82.1	100.00	Ref.1	53.1	70.37	Ref.2	81.1	62.19
T	74.1	100.00	Ref.1	56.1	96.17	Ref.2	73.1	38.55
T	72.1	100.00	Ref.1	44.1	179.52	Ref.2	57.1	58.00
T	43.0	100.00	Ref.1	45.0	15.04	Ref.2	70.0	13.90
T	103.1	100.00	Ref.1	45.1	477.76	Ref.2	73.1	268.85
T	72.1	100.00	Ref.1	43.1	422.88	Ref.2	57.1	33.46
T	90.1	100.00	Ref.1	75.0	87.93	Ref.2	82.0	45.05

- 1) Check the registered ions or transitions.
- 2) To change or add an MRM transition, select the setting in the [Type] column pull-down menu.

### MRM Transition

Ion1				Ion2				Ion3			
Type	m/z	CE	Rat i	Type	m/z	CE	Rat i	Type	m/z	CE	Rat i
T	82.10>54.00	9	100.00	Ref.1	81.10>53.00	6	76.60	Ref.2	82.10>39.00	18	70.48
T	72.10>57.00	3	100.00	Ref.1	72.10>44.10	6	29.54	Ref.2	72.10>43.00	3	23.44
T	70.0>55.0	6	100.00	Ref.1	88.0>61.0	6	14.32	Ref.2	70.0>43.0	18	14.77
T	73.10>45.10	6	100.00	Ref.1	103.10>75.10	6	26.44	Ref.2	103.10>47.10	12	22.51
T	72.10>43.10	6	100.00	Ref.1	72.10>57.10	3	60.62	Ref.2	57.10>41.10	3	1.87
T	90.10>75.10	9	100.00	Ref.1	90.10>62.00	6	91.84	Ref.2	75.00>47.00	6	77.31
T	71.10>43.10	6	100.00	Ref.1	87.10>55.00	9	74.70	Ref.2	102.10>87.10	3	45.62
T	86.10>43.10	9	100.00	Ref.1	86.10>71.10	3	94.06	Ref.2	43.00>41.10	6	42.32
T	96.10>81.00	9	100.00	Ref.1	81.10>53.00	9	85.53	Ref.2	96.10>53.00	21	26.18
T	102.10>74.10	3	100.00	Ref.1	74.10>56.10	6	16.07	Ref.2	74.10>46.00	15	13.42
T	116.10>88.10	3	100.00	Ref.1	88.10>73.10	9	97.45	Ref.2	116.10>73.10	9	64.07
T	86.00>43.00	3	100.00	Ref.1	43.00>15.00	15	47.92	Ref.2	43.00>14.00	30	29.29
T	86.10>71.10	6	100.00	Ref.1	86.10>58.00	3	92.11	Ref.2	86.10>43.10	18	55.28
T	86.10>57.10	6	100.00	Ref.1	86.10>56.00	3	33.92	Ref.2	57.10>42.00	24	5.69
T	61.00>43.10	12	100.00	Ref.1	73.10>43.00	3	82.30	Ref.2	43.10>41.10	3	8.35

### 2.3.4.3 Creating an MS Table

Specify parameter settings necessary for creating the MS table.

- 1) Click [Create Method File] to display the [MS Table Parameter] window.

## 2. Operating Procedure

---

2) Specify [MS Table Parameter] settings.

### Creating an MRM or SIM Mode Method File

Configure the settings in the [MRM, SIM Parameter] area. For simultaneous Scan/SIM or Scan/MRM measurements, select [Scan Mode]-[ON] and configure [Scan Mode] settings.

### Creating a Scan Mode Method File

Configure the settings in the [Scan Parameter] area.

The screenshot shows a dialog box titled "MS Table Parameter" with a close button (X) in the top right corner. The dialog is divided into two main sections: "MRM, SIM Parameter" and "Scan Parameter".

**MRM, SIM Parameter**

- Loop Time (MRM, SIM): 0.30 sec
- Required Processing Time : R.T ±: 0.30 min

**Scan Mode**

ON  OFF

**Scan Parameter**

- Event Time of Scan: 0.30 sec
- Scan Range : Start m/z - End m/z: 35 - 400
- Aquisition Time: Start R.T. - End R.T.: 1 - 35 min

At the bottom right, there are two buttons: "OK" (highlighted with an orange border) and "Cancel".

 **Hint**

Configure the [MS Table Parameter] settings.

[Loop Time (MRM, SIM)] Setting

This indicates the sum of all event times for all compounds specified for a single group. For example, setting it to 0.3 seconds results in acquiring data every 0.3 seconds for each compound. To ensure good reproducibility, at least ten data points must be acquired per peak. However, specifying a short loop time results in lower sensitivity due to less time available for measuring each transition or ion.


**Recommended value: 0.2 to 0.5 seconds**

[Required Processing Time] Setting

This indicates the measurement time range, mainly based on the retention times for each compound. For example, setting it to 0.3 minutes specifies acquiring data for each compound within  $\pm 0.3$  minutes of the specified retention time (0.6 minute time span).

Check retention time deviation by using the method file that was created to measure a sample spiked with a standard sample. Then adjust the setting based on any retention time shift or peak broadening (such as tailing) caused by contaminants or other factors.

**Recommended value: 0.3 to 0.5 minutes**



The required processing time specified here is applied to all compounds being measured. To specify required processing times individually for specific compounds, specify the values in cells AB to AC in the [Required Proc. Time for Each Comp.] column of the [Database] sheet. If the [Required Proc. Time for Each Comp.] column is blank, then the required processing time settings specified in the MS table sub-window are used.

Required Proc. Time for Each Comp. R.T. - X min	Required Proc. Time for Each Comp. R.T. + Y min
0.50	0.50


## 2. Operating Procedure

### Hint

The scan parameters for each scan method in this product are indicated below. Use the values as a reference.


File Name (.qgm)	Event Time (sec)	Scan Range (start-end m/z values)	Analysis Time (min) (start-end R.T. values)
Aroma_QP/TQ_SH-I-5SiIMS_Scan	0.3	35 - 400	1 - 45
Aroma_QP/TQ_SH-PolarWax_Scan.qgm	0.3	35 - 400	3 - 90
Aroma_QP/TQ_IC-WAX_Scan.qgm	0.3	35 - 400	1 - 35

- 3) Clicking [OK] displays the [Progress Bar] window. A method file is created automatically with MS table and compound table configured.



### NOTE


If "Target" is specified for a compound without an ion or transition specified, the following type of message is displayed. If [OK] is clicked, an analytical method is created that does not measure compounds without an ion or transition specified.




## 2. Operating Procedure

### NOTE

If the number of compounds or number of ions/transitions in the MS table being created exceeds the maximum number of channels, the following type of message is displayed. If that occurs, reduce the number of compounds or the number of ions or transitions for each compound.




- 4) Name the file as desired and click [Save].




- 5) Check the time range for MRM and SIM measurements of target compounds. An MS table view is not created for scan mode measurements.

## 2. Operating Procedure



### Hint

The measurement time per transition (per ion) (hereafter referred to as “dwell time”) for each target compound specified in the created MRM or SIM analytical method can be confirmed on the [MSTableView] sheet.



The minimum dwell time value for each target compound can also be confirmed. The more target compounds that are specified, the shorter the dwell time and the lower the sensitivity. It is recommended that settings are specified for the number of compounds and number of transitions/ions that ensure a minimum dwell time of about 2 ms.

**20211130\_1.qgm**


Minimum Dwell time                      7.0 ms

## 2. Operating Procedure

### NOTE

If using an AOC-6000 series or HS-20 series unit, method parameters for the pretreatment system must be specified in the created method file. For a description of that procedure, refer to Appendix 4 Using the Database with an AOC-6000 Series Autosampler or Appendix 6 Using the Database with an HS-20 Series Headspace Sampler.

- 6) Name the Smart Database and save it in the folder that was created.



- 7) Measure samples using the method file created.

### 2.3.5 Editing Databases

Smart Database files are saved in Excel format, which means they can be freely revised and compounds can be freely added.

If a different column or GC conditions are used, a new method for use under the given conditions can be created by adjusting retention index values. Additional new compounds or transition information can be added by using the MRM Optimization Tool in GCMSSolution to register them.

#### Revising a Database to Support Different GC Analytical Conditions

Save registered retention index values in an existing database by following the procedure described up to 3.3.3 Registering Retention Index Values in Compound Tables in 3. Registering Additional New Compound Information in a Smart Database in the GCMS Operation Guide – Method Development. Then complete the process by saving the database with a different name.

## 2. Operating Procedure

---

### Registering New Compounds or Ions/Transitions

Register the information by following the series of steps described in 3. Registering Additional New Compound Information in a Smart Database in the GCMS Operation Guide – Method Development. After registering the information, complete the process by saving the database with a different name.


## 3 Precautions

### 3.1 Changing Analytical Conditions or Using a Column that was Not Specified

Retention index and calibration curve information registered in the database was obtained by analysis using the analytical conditions in Appendix 1 Analytical Conditions. Beware that changing the analytical conditions or using a column other than the type specified could prevent accurate analysis because it would result in different retention index values.

### 3.2 Functionality for Protecting the Detector

The detector voltage is specified in advance, but that can cause detector saturation for some samples. If the [The detector turns off automatically when the signal to it is saturated.] checkbox is selected, the filament is automatically switched OFF if detector saturation occurs. In such cases, that would prevent acquiring data normally, so clear the checkbox for that detector protective function.



## Appendix 1 Analytical Conditions

- SH-I-5Sil MS

GC

Injection Mode:	Split (split ratio 5)
Sample Injection Unit Temperature:	250 °C
Oven Temperature:	50 °C (5 min) → (10 °C/min) → 250 °C (10 min)
Carrier Gas:	Helium
Control Mode:	Constant pressure (83.5 kPa)
Purge Flow Rate:	3.0 mL/min

If a Sniffer Unit is Used

Control Mode:	Constant pressure (160 kPa)
Total Flow Rate:	17.2 mL/min
APC1 Pressure:	20 kPa
APC2 Pressure:	200 kPa

MS

Ion Source Temperature:	200 °C
Interface Temperature:	250 °C

- SH-PolarWax

GC

Injection Mode:	Split (split ratio 5)
Sample Injection Unit Temperature:	250 °C
Oven Temperature:	40 °C (5 min) → (3 °C/min) → 250 °C (15 min)
Carrier Gas:	Helium
Control Mode:	Constant linear velocity (25.5 cm/s)
Purge Flow Rate:	3.0 mL/min

If a Sniffer Unit is Used

Control Mode:	Constant pressure (230 kPa)
Total Flow Rate:	15.5 mL/min
APC1 Pressure:	20 kPa
APC2 Pressure:	200 kPa

MS

Ion Source Temperature:	200 °C
Interface Temperature:	250 °C

- InertCap PureWAX

## Appendix 1 Analytical Conditions

---

### GC

Injection Mode:	Split (split ratio 5)
Sample Injection Unit Temperature:	250 °C
Oven Temperature:	50 °C (5 min) → (10 °C/min) → 250 °C (10 min)
Carrier Gas:	Helium
Control Mode:	Constant pressure (83.5 kPa)
Purge Flow Rate:	3.0 mL/min

### If a Sniffer Unit is Used

Control Mode:	Constant pressure (160 kPa)
Total Flow Rate:	17.2 mL/min
APC1 Pressure:	20 kPa
APC2 Pressure:	200 kPa

### MS

Ion Source Temperature:	200 °C
Interface Temperature:	250 °C



## Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds (Scan)

Some compounds registered in the database include registered gradient correction factors based on the internal standard method. Therefore, by measuring the internal standard substance in advance, approximate quantitation values for detected compounds can be calculated based on the area value of the internal standard substance.



### NOTE

#### Applications and Precautions for Using Semi-Quantitative Results

Data measured from injecting 1  $\mu\text{L}$  of a mixture containing a standard sample (10  $\mu\text{g/mL}$ , or 50  $\mu\text{g/mL}$  for some) and an internal standard mixture solution (10  $\mu\text{g/mL}$ ) were used to calculate relative response factors (RF) based on the following formula.

$$\text{RF} = \frac{\text{Area}_{\text{Target}} \times \text{Conc}_{\text{ISTD}}}{\text{Conc}_{\text{Target}} \times \text{Area}_{\text{ISTD}}}$$

$\text{Area}_{\text{Target}}$  and  $\text{Conc}_{\text{Target}}$ : Area and concentration values for targets

$\text{Area}_{\text{ISTD}}$  and  $\text{Conc}_{\text{ISTD}}$ : Area and concentration values for internal standards

Of course, that means a different GCMS system or capillary column will result in a different RF value. Therefore, if the RF values from database obtained are used, the resulting values will be different than the true values. Furthermore, the values do not reflect the actual extraction efficiency during pretreatment, which means results can vary significantly depending on the pretreatment system used. However, the values can be used for comparing odor/aroma threshold values (minimum concentration detectable by humans) for detected compounds or as a reference guideline for adjusting standard samples.

Given that calculated concentration values are only approximate, they are less accurate than values determined by conventional quantitative methods based on preparing calibration curves. Therefore, to determine concentration values, calculate quantitative results based on a calibration curve prepared separately.

 **Hint**

The following internal standard samples are recommended.

- EPA 524.2 Fortification Solution Cat: 47358-U  
(2000 µg/mL 4-Bromofluorobenzene, 1,2-Dichlorobenzene-d4, Fluorobenzene in methanol)  
SIGMA-ALDRICH
  
- Acenaphthene-d10 solution Cat: 48417  
(2000 µg/mL component in methylene dichloride)  
SIGMA-ALDRICH

<Preparing Samples for internal standard mixture>

Dilute an internal standard/surrogate standard substance mixture (2000 µg/mL) and acenaphthene-d10 solution (2000 µg/mL) with hexane to a concentration of 10 µg/mL.

Preparation Example:


Use a microsyringe to place 50 µL each of the internal standard/surrogate standard substance mixture (2000 µg/mL) and acenaphthene-d10 solution (2000 µg/mL) in a 10 mL volumetric flask and then fill the flask to volume with hexane. Keep the prepared solution stored in a refrigerator or freezer until ready for use.

### Appendix 3.1 Measuring and Analyzing Standard Samples

n-alkanes are measured for adjusting retention times and internal standard solutions are measured for semi-quantitative analysis.

- 1) Use a method file for adjusting retention times to analyze the n-alkane solution (50 µg/mL). Use a method file for sensitivity correction to analyze the internal standard mixture solution (10 µg/mL).
- 2) Load the n-alkane data file and identify the n-alkanes as described in 2.1.2 Measuring an n-Alkane Standard Mixture Sample. After analyzing the data, save the results by overwriting.
- 3) Load the internal standard data to check and correct the identification results for each component. If any are misidentified, identify them correctly either by manual identification or by manual peak integration. After analyzing the data, save the results by overwriting.

## Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds



### Appendix 3.2 Creating Method Files

- 1) Load the data file from analyzing n-alkane data.
- 2) Specify the data as described in 2.3.3 Specifying the Smart Database.
- 3) Specify measurement target compounds as described in 2.3.4.1 Setting Target Compounds. For the template method file, specify the scan method file ("Aroma\_QP/TQ\_<column name>\_Scan.qgm").

#### NOTE

Semi-quantitative analysis is only supported for scan mode analysis.

RF values for registered compounds are calculated based on the area values for target (marked "T") compounds and ions specified in the database. If target ions are changed, RF values cannot be used for semi-quantitative analysis.

## Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds

- 4) Click [Advanced].

The screenshot shows the 'Create Method File' dialog box. At the top, there is a section for 'Instrument Type' with a dropdown menu set to 'TQ Series' and a 'Lang.' button. Below this is a 'Parameter' section with four rows of controls: 'Ret. Index for AART' with a dropdown menu set to 'Ret. Index 1'; 'n-alkane data file' with a text box containing 'C:\GCMSsolution\Data\Project1\Alkane.qgd' and a browse button; 'Template Method File' with a text box containing 'Data\Project1\Aroma\_TQ-IC-Wax\_Template.qgm' and a browse button; and 'Divide Method into' with a dropdown menu set to '1'. An 'Advanced' button is highlighted with an orange box at the bottom right of the parameter section. There is also an 'Import' button at the bottom left.

- 5) Select the [Use Semi-Quantitative Calibration Curve] checkbox. Select the [Adjust Calibration Curve Std Data for Adjustment:] checkbox and select the data file from the sensitivity correction sample analyzed.

The screenshot shows the 'Advanced Settings' dialog box. It has a title bar with a close button. The main area is titled 'Calibration Curve' and contains two checked checkboxes: 'Use Semi-Quantitative Calibration Curve' and 'Adjust Calibration Curve'. Below these is a label 'Std Data for Adjustment:' followed by a text box containing 'D:\GCMSsolution\Data\Project1\Correct.qgd' and a browse button. The 'OK' and 'Cancel' buttons are at the bottom.

- 6) Click [OK] to close the window.
- 7) Click [Create Method File] to display the [MS Table Parameter] window.

The screenshot shows the 'Create Method File' dialog box. The 'Create Method File' button at the top left is highlighted with an orange box. The 'Instrument Type' dropdown is set to 'TQ Series'. The 'Parameter' section has the following values: 'Ret. Index for AART' is 'Ret. Index 1'; 'n-alkane data file' is 'C:\GCMSsolution\Data\Project1\Alkane.qgd'; 'Template Method File' is 'n\Data\Project1\Aroma\_TQ-IC-Wax\_Scan.qgm'; and 'Divide Method into' is '1'. The 'Advanced' button is visible at the bottom right.

## Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds

---

- 8) Click [OK]. Clicking [OK] displays the [Progress Bar] window and automatically creates a method file.

MS Table Parameter

**MRM, SIM Parameter**

Loop Time (MRM, SIM) 0.30 sec

Required Processing Time : R.T ± 0.30 min

**Scan Mode**  ON  OFF

**Scan Parameter**

Event Time of Scan: 0.30 sec

Scan Range : Start m/z - End m/z 35 - 400

Aquisition Time: Start R.T. - End R.T. 1 - 35 min

OK Cancel

- 9) Name the file as desired and click [Save].

### NOTE

If using an AOC-6000 series or HS-20 series unit, method parameters for the pretreatment system must be specified in the created method file. For a description of that procedure, refer to Appendix 4 Using the Database with an AOC-6000 Series Autosampler or Appendix 6 Using the Database with an HS-20 Series Headspace Sampler.

## Appendix 3.3 Analyzing Samples


Analyze samples using the method file created for semi-quantitative analysis.

### Hint

Semi-quantitative values calculated from analytical results are indicated as absolute values (pg) injected in the GC unit.

If the sample quantity and dilution factors are entered for batch analysis, the software calculates sample concentrations automatically. Depending on the type of sample, use the following steps to change the concentration units specified in the method file.


- 1) Load the method file created by the [GCMS Analysis] program.
- 2) On the [Method] menu bar, click [Quantitative Parameters...].



- 3) Click the [Compound Table] tab. In the [Unit] column, enter a concentration unit suitable for the weight of the sample obtained.

Entry example: For solid samples, enter “pg/mg” (sample quantity acquired)

For liquid samples, enter “pg/μL” (sample quantity acquired)




- 4) Click [OK] to save the method file by overwriting.

- 1) Place the measurement samples in the autosampler.
- 2) Open the [Batch Processing] window and create a new batch table.

### Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds

- Right-click on the batch table and click [Table Style] on the right-click menu. The [Table Style] window is displayed.
- In the [Hide Items] field, select [Sample Amt.] and [Dil. Factor]. Then click [Add>>] to add those items to the [Display Items] field.



- Specify settings in the [Vial #], [Sample Name], and [Sample ID] columns. Specify the created file as the method file. To automatically calculate concentrations in samples based on quantitative results, specify [Sample Amt.] and [Dil. Factor] values.

Folder: C:\Users\shimadzu\Desktop\PIP\_ENG\SH-1MS

	Vial#	Sample Name	Sample ID	Sample Amt.	Dil. Factor	Sample Type	Analysis Type	Method File
1	1	Sample A	0001	2	100	0:Unknown	IT QT	ax_Scan.qgm
2	1	Sample A	0001	1	1	0:Unknown	IT QT	ax_Scan.qgm
3	2	Sample B	0002	2	100	0:Unknown	IT QT	ax_Scan.qgm
4	2	Sample B	0002	1	1	0:Unknown	IT QT	ax_Scan.qgm
5	3	Sample C	0003	2	100	0:Unknown	IT QT	ax_Scan.qgm
6	3	Sample C	0003	1	1	0:Unknown	IT QT	ax_Scan.qgm
7	4	Sample D	0004	2	100	0:Unknown	IT QT	ax_Scan.qgm
8	4	Sample D	0004	1	1	0:Unknown	IT QT	ax_Scan.qgm
9	5	Sample E	0005	2	100	0:Unknown	IT QT	ax_Scan.qgm
10	5	Sample E	0005	1	1	0:Unknown	IT QT	ax_Scan.qgm
11	6	Sample F	0006	2	100	0:Unknown	IT QT	ax_Scan.qgm
12	6	Sample F	0006	1	1	0:Unknown	IT QT	ax_Scan.qgm
13	7	Sample G	0007	2	100	0:Unknown	IT QT	ax_Scan.qgm
14	7	Sample G	0007	1	1	0:Unknown	IT QT	ax_Scan.qgm
15	8	Sample H	0008	2	100	0:Unknown	IT QT	ax_Scan.qgm
16	8	Sample H	0008	1	1	0:Unknown	IT QT	ax_Scan.qgm



### Hint

#### Examples of [Sample Amt.] and [Dil. Factor] Settings

The [Sample Amt.] value is the eventual quantitation value calculated by dividing the concentration value determined from the calibration curve by the value specified here. The [Dil. Factor] value is the eventual quantitation value calculated by multiplying the concentration value determined from the calibration curve by the value specified here.

#### [Sample Amt.]

Enter the acquired sample weight in the compound table, in terms of the [Unit] setting.

Example: Enter "2.5" for a sample weighing 2.5 mg.

#### [Dil. Factor]

Mainly entered for solvent extracted samples injected as a liquid, enter values as a ratio of the final extract solution divided by the sample injection volume.


Example: Assume a 2.5 mg solvent extracted sample is used to prepare 1000  $\mu\text{L}$  of a final extract solution, from which 1  $\mu\text{L}$  is injected into the GC/MS system.

Enter Sample Amt.: "2.5" (mg) and Dil Factor:  $1000 (\mu\text{L}) / 1 (\mu\text{L}) = "1000"$ .

- 6) Name and save the batch file.
- 7) Execute realtime batch analysis.

## Appendix 3.4 Checking Measurement Results

- 1) Check the identification results. The quantitative results are displayed for identified compounds.



## Appendix 4 Using the Database with AOC-6000 Series Autosampler

This chapter describes the operating procedure when using an AOC-6000 series sampler as a pretreatment system. Using an AOC-6000 series sampler requires that AOC-6000 method parameters consistent with the given sample injection method are specified in the methods file.



### Hint

For AOC-6000 series operating instructions, refer to the AOC-6000 Operation Guide and AOC-6000 Control Software Instruction Manual.

### Appendix 4.1 Instrument Configuration Settings

Configure instrument configuration settings as described in 3.2 Setting the Configuration in the AOC-6000 Control Software Instruction Manual.

### Appendix 4.2 Analyzing n-Alkanes

If using an AOC-6000 series system to analyze n-alkanes for retention time adjustment, use the liquid injection method.

- 1) Install a liquid injection syringe adapter (10  $\mu$ L capacity) on the AOC-6000 series system. For adapter installation instructions, refer to 4. PAL Tool in the AOC-6000 Operation Guide.
- 2) Specify method parameter settings for the n-alkane measurement method (“Aroma\_QP/TQ\_<column>\_AART.qgm”) as described in 4. Method Settings in the AOC-6000 Control Software Instruction Manual.  
Select “Liquid Injection” for the AOC-6000 method and specify liquid injection parameter settings indicated in Appendix 4.4.
- 3) Place an n-alkane standard mixture sample (50  $\mu$ g/mL) on the AOC-6000 rack.
- 4) Analyze the n-alkane sample as described in 5. Measuring Samples in the AOC-6000 Control Software Instruction Manual.  
Set the injection volume to 1  $\mu$ L.

 **Hint**

If an AOC-6000 series system is used for semi-quantitative analysis according to the procedure in Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds, then also analyze the internal standard substance using the same liquid injection method as used for the n-alkanes. However, note that semi-quantitative results do not reflect extraction efficiency during pretreatment, so they can differ significantly from actual quantitation values, depending on the sample injection method used.

## Appendix 4.3 Analyzing Samples

Revise the method parameter settings based on the sample injection method used.

- 1) Use the n-alkane measurement data to adjust the retention times in the method file as described in 2.1.4 Adjusting Retention Times and Analyzing Samples. If using a database to create the method, follow the procedure in 2.3 Using a Smart Database to Create Method Files, starting from 2.3.3 Specifying the Smart Database.
- 2) Replace the GC injection port insert liner based on the sample injection method used for sample analysis. Also install a syringe tool for the AOC-6000 system.
- 3) In the method adjusted or created in step 1), specify method parameter settings for the sample injection method used, according to the procedure described in 4. Method Settings in the AOC-6000 Control Software Instruction Manual. For information about AOC-6000 method parameter settings, refer to Appendix 4.4 Examples of Method Parameter Settings.
- 4) Place the measurement samples in the AOC-6000 rack and start analyzing them using the method file indicated in step 3).

## Appendix 4.4 Examples of Method Parameter Settings

Use the following examples as a reference for specifying method parameter settings for respective AOC-6000 sample injection methods. For more information about each parameter, refer to Appendix – in the AOC-6000 Control Software Instruction Manual.

- **Liquid Injection**

- Analysis Group

Syringe Tool:	Selected from the list of syringe tools for liquid injection (such as LS1 or LS2).
Pre Wash Cycles:	3
Sample Rinse Cycles:	1

## Appendix 4 Using the Database with AOC-6000 Series Autosampler

---

Sample Aspirate Flow Rate:	1 $\mu$ L/s
Sample Post Aspirate Delay:	2 s
Injector:	Selected from the list of injectors (such as Injector1 or Injector2).
Injection Flow Rate:	50 $\mu$ L/min
Post Wash Cycles:	3

### Setup Group

Gas Chromatograph:	Selected from the list of GC units (such as GC1).
Cooled Stack 1:	Not selected.
Cooled Stack 1 Temperature:	20 $^{\circ}$ C
Cooled Stack 2:	Not selected.
Cooled Stack 2 Temperature:	20 $^{\circ}$ C
Pre Wash Station:	Selected from the list of wash stations.
Post Wash Station:	Selected from the list of wash stations.
Wash Vial Depth:	40 mm
Waste Port Depth:	10 mm
Sample Vial Depth:	30 mm
Bottom Sense Sample Vial:	Off
Height From Bottom Sample Vial:	3 mm
Injection Mode:	Normal
Injection Signal Mode:	PlungerDown
Injector Penetration Depth:	40 mm
Pre Ahead Option:	Disable

### Advanced Group

Pre Wash Solvent Volume:	8 $\mu$ L
Pre Wash Solvent Position Step 1:	1
Pre Wash Solvent Position Step 2:	0
Pre Wash Solvent Position Step 3:	0
Pre Wash Solvent Position Step 4:	0
Pre Wash Aspirate Flow Rate:	1 $\mu$ L/s
Sample Vial Penetration Speed:	50 mm/s
Sample Rinse Volume:	3 $\mu$ L
Filling Strokes Count:	5
Filling Strokes Volume:	3 $\mu$ L
Filling Strokes Aspirate Flow Rate:	5 $\mu$ L/s
Filling Strokes Post Aspirate Delay:	0.2 s
Filling Strokes Post Dispense Delay:	0.5 s
Delay After Filling Strokes:	0.5 s
Air Gap Volume:	0 $\mu$ L
Injector Penetration Speed:	100 mm/s

Pre Injection Dwell Time:	0 s
Post Injection Dwell Time:	0.3 s
Post Wash Solvent Volume:	8 $\mu$ L
Post Wash Solvent Position Step 1:	1
Post Wash Solvent Position Step 2:	0
Post Wash Solvent Position Step 3:	0
Post Wash Solvent Position Step 4:	0
Post Wash Aspirate Flow Rate:	1 $\mu$ L/s



### NOTE

The settings above are for using a syringe with a metal plunger. If a syringe with a PTFE plunger is used, change the setting for [Filling Stroke Aspirate Flow Rate] from 5  $\mu$ L/s to 2  $\mu$ L/s, for example, as necessary to ensure samples are aspirated slowly enough to prevent bubbles inside the syringe.

#### • SPME Injection

##### Analysis Group

Fiber Tool:	Selected from the list of syringe tools for SPME.
Conditioning Port:	Selected from the pull-down menu.
Conditioning Temperature:	Set to the recommended value for the selected SPME fiber.
Pre Conditioning Time:	Set to the recommended value for the selected SPME fiber.
Incubation Temperature:	80 °C
Incubation Time:	5 min
Agitator Speed:	250 rpm
Sample Vial Depth:	22 mm
Sample Extract Time:	30 min
Injector:	Selected from the list of injectors (such as Injector1 or Injector2).
Sample Desorb Time:	2 min
Post Conditioning Time:	Set to the recommended value for the selected SPME fiber.
Analysis Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

### Setup Group

Gas Chromatograph:	Selected from the list of GC units (such as GC1).
Sync Before Extraction End Time:	0 min
Agitator:	Selected from the list of agitators.
Do Agitation:	True
Heat Agitator:	True
Wait For Readiness Agitator:	True
Internal Standard Station:	Not selected.
Internal Standard Position:	1
Injector Penetration Depth:	54 mm
Injection Signal Mode:	Before Expose

### Advanced Group

Internal Standard Adsorb Time:	0 min
Internal Standard Penetration Depth:	32 mm
Sample Vial Penetration Speed:	20 mm/s
Injector Penetration Speed:	100 mm/s
Agitator On Time:	5 s
Agitator Off Time:	2 s



### **NOTE**

Select the type of liquid phase for the SPME fiber based on the targeted components.

- **SPME Arrow Injection**

### Analysis Group

SPME Arrow Tool:	Selected from the list of syringe tools for SPME Arrow.
Conditioning Port:	Selected from the pull-down menu.
Conditioning Temperature:	Set to the recommended value for the selected SPME Arrow.
Pre Conditioning Time:	Set to the recommended value for the selected SPME Arrow.
Incubation Temperature:	80 °C
Incubation Time:	5 min
Agitator Speed:	250 rpm
Stirrer Speed:	250 rpm
Sample Vial Depth:	33 mm

## Appendix 4 Using the Database with AOC-6000 Series Autosampler

---

Sample Extract Time:	30 min
Injector:	Selected from the list of injectors (such as Injector1 or Injector2).
Sample Desorb Time:	2 min
Post Conditioning Time:	Set to the recommended value for the selected SPME Arrow.
Analysis Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

### Setup Group

Gas Chromatograph:	Selected from the list of GC units (such as GC1).
Wash Station:	1
Do Pre Wash:	False
Agitator:	Selected from the list of agitators.
Stirrer:	True
Do Agitation:	True
Heat Agitator and Stirrer:	True
Internal Standard Station:	Not selected.
Internal Standard Position:	1
Injector Penetration Depth:	54 mm
Injection Signal Mode:	Before Expose
Do Post Wash:	False

### Advanced Group

Pre Solvent Index:	1
Pre Wash Time:	0 min
Pre Wash Vial Depth:	35 mm
Internal Standard Adsorb Time:	0 min
Internal Standard Penetration Depth:	32 mm
Sample Vial Penetration Speed:	20 mm/s
Injector Penetration Speed:	100 mm/s
Post Solvent Index:	1
Post Wash Time:	0 min
Post Wash Vial Depth:	35 mm

 **NOTE**

Select the type of liquid phase for the SPME Arrow based on the targeted components.

- **Headspace Injection**

Analysis Group

Syringe Tool:	Selected from the list of syringe tools for HS units (such as HS1 or HS2).
Incubation Temperature:	60 °C
Incubation Time:	30 min
Syringe Temperature:	80 °C
Agitator Speed:	250 rpm
Pre Purge Time:	0 s
Injector:	Selected from the list of injectors (such as Injector1 or Injector2).
Injection Flow Rate:	25 mL/min
Post Purge Time:	60 s
Analysis Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

Setup Group

Gas Chromatograph:	Selected from the list of GC units (such as GC1).
Sync Before Incubation End Time:	0 min
Agitator:	Selected from the list of agitators.
Do Agitation:	True
Heat Agitator:	True
Wait For Readiness Agitator:	True
Sample Vial Depth:	15 mm
Heat Syringe:	True
Wait For Readiness Syringe:	True
Injection Signal Mode:	PlungerUp
Injector Penetration Depth:	40 mm
Continuous Purge:	Flase

Advanced Group

Enable Pre Filling:	True
Filling Strokes Count:	5
Filling Strokes Volume:	1.2 mL
Filling Strokes Aspirate Flow Rate:	6 mL/min
Delay After Filling Strokes:	30 s
Sample Aspirate Flow Rate:	6 mL/min
Sample Post Aspirate Delay:	0 s

Sample Vial Penetration Speed:	25 mm/s
Injector Penetration Speed:	25 mm/s
Pre Injection Dwell Time:	3 s
Post Injection Dwell Time:	10 s
Agitator On Time:	5 s
Agitator Off Time:	2 s

 **NOTE**

Depending on the version of the AOC-6000 control software, the “Headspace Injection” method has different items in the [Advanced Group]. It is recommended that the default values for the parameters of the [Advanced Group] are used as they are. The parameter examples are shown for AOC-6000 control software version 2.0 or later.

 **NOTE**

Set the injection volume to 500 µL for batch analysis.

 **注記**

Set the split ratio to 5 or higher. The splitless mode or low split ratios can cause poor peak shape due to pressure fluctuations at the injection port.

- **ITEX Injection**

Analysis Group

ITEX Tool:	Selected from the list of syringe tools for ITEX.
Trap Pre Cleaning Temperature:	250 °C
Trap Pre Cleaning Time:	300 s
Incubation Temperature:	80 °C
Incubation Time:	5 min
Agitator Speed:	250 rpm
Syringe Temperature:	90 °C
Trap Extraction Temperature:	40 °C
Extraction Strokes:	50
Injector:	Selected from the list of injectors (such as Injector1 or Injector2).

## Appendix 4 Using the Database with AOC-6000 Series Autosampler

---

Injector Aspirate Flow Rate:	10 $\mu\text{L/s}$
Desorb Temperature:	250 $^{\circ}\text{C}$
Desorb Flow Rate:	100 $\mu\text{L/s}$
Trap Post Cleaning Temperature:	250 $^{\circ}\text{C}$
Trap Post Cleaning Time:	0 s
Analysis Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

## Appendix 4 Using the Database with AOC-6000 Series Autosampler

---

### Setup Group

Gas Chromatograph:	Selected from the list of GC units (such as GC1).
Sync Before Incubation:	0 min
Agitator:	Selected from the list of agitators.
Heat Agitator:	True
Sample Vial Depth:	15 mm
Water Removal:	Disable
Water Removal Position:	Not selected.
Injection Signal Mode:	PlungerUp
Injector Penetration Depth:	33 mm

### Advanced Group

Extraction Volume:	1000 $\mu$ L
Extraction Aspirate Flow Rate:	100 $\mu$ L/s
Extraction Dispense Flow Rate:	100 $\mu$ L/s
Sample Prefill Rate:	40 %
Water Removal Trap Temperature:	60 $^{\circ}$ C
Water Removal Purge Time:	0 s
Injector Penetration Speed:	50 mm/s
Post Inject Delay:	0 s

 **NOTE**

Set the injection volume to 500  $\mu$ L for batch analysis.

 **NOTE**

Set the split ratio to 5 or higher. The splitless mode or low split ratios can cause poor peak shape due to pressure fluctuations at the injection port.



**NOTE**

Select the type of trap sorbent used based on the targeted components.

## Appendix 5 Using the Database with TD-30 Series Autosampler

This chapter describes the operating procedure when using an TD-30 series sampler as a pretreatment system. To use a TD-30 series thermal desorption unit, create a method file for the TD-30 unit.



### Hint

For TD-30 series operating instructions, refer to the TD-30 Series Instruction Manual.

### Appendix 5.1 Instrument Configuration Settings

Configure instrument configuration settings as described in 2.2 GC Set-up in the TD-30 Series Instruction Manual.

### Appendix 5.2 Analyzing n-Alkanes

To use a TD-30 series system to analyze n-alkanes for retention time adjustment, measure the n-alkanes by adding them to TD sample tubes.



### NOTE

Pre-condition sample tubes before analyzing n-alkanes and samples. For the tube conditioning procedure, refer to 5.6 Conditioning in the TD-30 Series Instruction Manual.

- 1) Use a microsyringe to add 1  $\mu\text{L}$  of the n-alkane standard mixture sample (50  $\mu\text{g}/\text{mL}$ ) to the sample tubes. Place the sample tubes in the TD-30 sample tray.
- 2) Create a TD-30 method file and batch file for measuring n-alkanes as described in 4.1 Setting the Analytical Conditions in the TD-30 Series Instruction Manual. Use Appendix 5.4 Examples of Method Parameter Settings as a reference for specifying TD-30 method parameter settings.



**Hint**

If the GCMSsolution TD Add-in is installed, a TD-30 batch file does not need to be created because the TD-30 unit is controlled by GCMSsolution batch processing.

- 3) Analyze the n-alkane sample as described in 5.1. Starting Pretreatment in the TD-30 Series Instruction Manual.



**Hint**

If a TD-30 series system is used for semi-quantitative analysis according to the procedure in Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds, then also analyze the internal standard substance using the same procedure as used for the n-alkanes.

## Appendix 5.3 Analyzing Samples

- 1) Use the n-alkane measurement data to adjust the retention times in the method file as described in 2.1.4 Adjusting Retention Times and Analyzing Samples. If using a database to create the method, follow the procedure in 2.3 Using a Smart Database to Create Method Files, starting from 2.3.3 Specifying the Smart Database.
- 2) Create a TD-30 method file and batch file for measuring samples as described in 4.1 Setting the Analytical Conditions in the TD-30 Series Instruction Manual. For TD-30 parameter settings, refer to Appendix 5.4 Examples of Method Parameter Settings.
- 3) Measure samples as described in 5.1 Starting Pretreatment in the TD-30 Series Instruction Manual.

## Appendix 5.4 Examples of Method Parameter Settings

Use the following examples as a reference for specifying method parameter settings for respective TD-30 sample injection methods.

- **TD-30R**

Pre purge flow:	20 mL/min
Pre purge time:	0 min
Tube desorption temperature:	280 °C
Tube desorption flow rate:	60 mL/min
Tube desorption time:	5 min
Trap cooling temperature:	-20 °C
Trap desorption temperature:	250 °C
Trap desorption time:	5 min
Joint temperature:	250 °C
Valve temperature:	250 °C
Transfer line temperature:	250 °C
GC Cycle Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

### Internal STD

Variable volume: –

### Restore

Restore: –

Dry Purge: –

### Stand-by Temp.

Tube: 40 °C

Trap: 50 °C

- **TD-30**

Tube desorption temperature:	280 °C
Tube desorption flow rate:	60 mL/min
Tube desorption time:	5 min
Trap cooling temperature:	-20 °C
Trap desorption temperature:	250 °C
Trap desorption time:	5 min
Joint temperature:	250 °C
Valve temperature:	250 °C
Transfer line temperature:	250 °C
GC Cycle Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

Stand-by temperature

Tube:	40 °C
Trap:	50 °C

## Appendix 6 Using the Database with HS-20 Series Autosampler

This chapter describes the operating procedure when using an HS-20 series sampler as a pretreatment system. If using an HS-20 series headspace sampler, specify HS-20 method parameters in the method file to be used.

 **Hint**

For HS-20 series operating instructions, refer to the HS-20 Operation Guide and HS-20 Control Software Instruction Manual.

 **NOTE**

Due to the operating principle involved in analysis using headspace methods, components that are difficult to vaporize tend to be difficult to detect.

 **NOTE**

Semi-quantitative calculations described in Appendix 3 are not available when using the HS-20 series system.

### Appendix 6.1 Instrument Configuration Settings

Configure instrument configuration settings as described in the HS-20Control Software Instruction Manual.

### Appendix 6.2 Analyzing n-Alkanes

If using an HS-20 series system to analyze n-alkanes for retention time adjustment, add n-alkanes to an empty vial.

- 1) Add 5  $\mu\text{L}$  of the n-alkane standard mixture sample (1000  $\mu\text{g}/\text{mL}$ ) to an HS-20 vial. Place the vial in the HS-20 sample tray.
- 2) Specify method parameter settings according to the procedure 2.2 Method Development in

the HS-20 Control Software Instruction Manual. Use Appendix 6.4 Examples of Method Parameter Settings as a reference for specifying HS-20 method parameter settings.

- 3) Measure the n-alkane sample as described in 2.3 Batch Processing in the HS-20 Control Software Instruction Manual.

### Appendix 6.3 Analyzing Samples


- 1) Use the n-alkane measurement data to adjust the retention times in the method file as described in 2.1.4 Adjusting Retention Times and Analyzing Samples. If using a database to create the method, follow the procedure in 2.3 Using a Smart Database to Create Method Files, starting from 2.3.3 Specifying the Smart Database.

## Appendix 6 Using the Database with HS-20 Series Autosampler


### NOTE

If using an HS-20 system to measure n-alkanes, note that it cannot detect n-alkanes with all carbon numbers. If the vial is heated to 60 °C, then carbon numbers up to C18 can be detected by the trap mode and C16 by the loop mode. If undetected alkanes are misidentified, then delete those identification results according to the following procedure.


1. On the compound table [Results] tab page, select the alkanes that were not detected.



2. Right-click to select [Delete Result] on the right-click menu.



3. Click [No].





**NOTE**

AART-predicted retention times will tend to deviate more for compounds that have a retention index that is after detected alkanes.

- 2) Specify method parameter settings for the method adjusted or created in step 1), according to the procedure described in 2.2 Method Development in the HS-20 Control Software Instruction Manual. Use Appendix 6.4 Examples of Method Parameter Settings as a reference for specifying HS-20 method parameter settings.
- 3) Place the measurement samples in the HS-20 rack and start analysis using the method file from step 3).

## Appendix 6.4 Examples of Method Parameter Settings

Use the following examples as a reference for specifying method parameter settings for respective HS-20 sample injection methods.

- **HS-20 (Loop Mode)**

Mode:	Loop
Oven Temp.:	60 °C
Sample Line Temp.:	100 °C
Transfer Line Temp.:	100 °C
Shaking Level:	Off
Multi Injection Count:	1
Pressurize Gas Pressure:	80 kPa
Equilibrating Time:	10 min
Pressurizing Time:	1 min
Pressure Equilib. Time:	0.1 min
Load Time:	1 min
Load Equilib. Time:	0.1 min
Injection Time:	3 min
Needle Flush Time:	5 min
GC Cycle Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

- **HS-20 (Trap Mode)**

Mode:	Trap
Oven Temp.:	60 °C
Sample Line Temp.:	100 °C

## Appendix 6 Using the Database with HS-20 Series Autosampler

---

Transfer Line Temp.:	100 °C
Trap Cooling Temp.:	-10 °C
Trap Desorb Temp.:	280 °C
Trap Equilib. Temp.:	25 °C
Shaking Level:	Off
Multi Injection Count:	5
Pressurize Gas Pressure:	80 kPa
Dry Purge Gas Pressure:	60 kPa
Equilibrating Time:	10 min
Pressurizing Time:	1 min
Pressure Equilib. Time:	0.1 min
Load Time:	1 min
Load Equilib. Time:	0.1 min
Dry Purge Time:	10 min
Injection Time:	3 min (n-Alkanes: 30 min)
Needle Flush Time:	5 min
GC Cycle Time:	Sum of GC program time, cooling time, and equilibration time, plus an extra margin.

## Appendix 7 Using the Database with OPTIC-4 Series Inlet

This chapter describes the operating procedure when using an OPTIC-4 series inlet unit as a pretreatment system. To use an OPTIC-4 unit, create an OPTIC-4 method file.



### Hint

For OPTIC-4 unit operating instructions, refer to the OPTIC-4 Instruction Manual.



### Hint

For instructions on operating it in combination with an AOC-6000 series optional LINEX-2 system, refer to the separate LINEX Method Instructions for AOC-6000 Control Software Instruction Manual.

### Appendix 7.1 Instrument Configuration Settings

Configure instrument configuration settings as described in 6. Basic Operations in the OPTIC-4 Instruction Manual.

### Appendix 7.2 Analyzing n-Alkanes

If using an OPTIC-4 unit to analyze n-alkanes for retention time adjustment, use the liquid injection method.

- 1) Install an OPTIC-4 insert liner for liquid injection.
- 2) Create a TD-30 method file and batch file for measuring n-alkanes as described in 7. Creating and Executing Methods in the OPTIC-4 Instruction Manual. Use Appendix 7.4 Examples of Method Parameter Settings as a reference for specifying OPTIC-4 method parameter settings.
- 3) Use the OPTIC-4 method created in step 1) and the GCMS method created in step 2) to analyze the n-alkane standard mixture sample (50 µg/mL).

 **Hint**

If an OPTIC-4 series system is used for semi-quantitative analysis according to the procedure in Appendix 3 Calculating Semi-Quantitative Results for Registered Compounds, then also analyze the internal standard substance using the same liquid injection method as used for the n-alkanes.

## Appendix 7.3 Analyzing Samples

Revise the method parameter settings based on the sample injection method used.

- 1) Use the n-alkane measurement data to adjust the retention times in the method file as described in 2.1.4 Adjusting Retention Times and Analyzing Samples. If using a database to create the method, follow the procedure in 2.3 Using a Smart Database to Create Method Files, starting from 2.3.3 Specifying the Smart Database.
- 2) Replace the OPTIC-4 port insert liner based on the sample injection method used for sample analysis.
- 3) In the method adjusted or created in step 1), specify method parameter settings for the sample injection method used, according to the OPTIC-4 Instruction Manual. For information about OPTIC-4 method parameter settings, refer to Appendix 7.4 Examples of Method Parameter Settings.
- 4) Start analyzing the measurement samples using the method file indicated in step 3).

## Appendix 7.4 Examples of Method Parameter Settings

Use the following parameter settings as a reference for specifying method parameter settings for respective OPTIC-4 sample injection methods. For more information about respective parameters, refer to 7.4 Method Parameters in the OPTIC-4 Instruction Manual.

- **Expert**

General

Equilibration Time: 5 sec  
End Time: GC program duration

Inlet Temperature

Initial Temperature: 40 °C  
Delay Time: 0 sec  
Ramp Rate 1: 50 °C/sec

## Appendix 7 Using the Database with an OPTIC-4 Series Inlet

---

Hold Temperature 1: 250 °C  
Solvent Colling Effect: No  
Cooling Valve Mode: No

### Septum Purge Flow

Septum Purge Flow: 3 mL/min

### Solvent Vending

Vent Mode: Fixed Time  
Vent Time: 0

If Using SH-I-5Sil MS or InertCap Pure-Wax Column

### Column Flow / Inlet Pressure

Carrier Control Mode: Pressure  
Zero LINEX Head Pressure: Yes  
Initial Inlet Pressure 1: 83.5 kPa  
Start Inlet Pressure 1: 83.5 kPa  
End Inlet Pressure 1: 83.5 kPa

### Split Flow

Direct Split Value Control: No  
Initial Split Flow: 5 mL/min  
Split Flow 1: 5 mL/min

If Using SH-PolarWax Column

### Column Flow / Inlet Pressure

Carrier Control Mode: Flow  
Zero LINEX Head Pressure: Yes  
Initial Column Flow 1: 1 mL/min  
Start Column Flow: 1 mL/min  
End Column Flow: 1 mL/min

### Split Flow

Direct Split Value Control: No  
Initial Split Flow: 5 mL/min  
Split Flow 1: 5 mL/min

## Appendix 7 Using the Database with an OPTIC-4 Series Inlet

---

### Cryotrap (only if cryogenics is available)

Cryotap Low Temperature:	-150 °C
Low Temperature Hold Time:	300 °C
Cryotrap High Temperature:	60 °C/sec
Cryotrap Heat Ramp Rate:	250 °C

## Appendix 8 Connecting an Odor Sniffer




### NOTE

The retention index values included in the product were acquired using only a standalone GC/MS system. Splitting flows between an odor sniffer unit and an MS unit can cause a large shift in the retention index values that varies depending on the compound.

### Appendix 8.1 Flow Restrictor Settings for Connecting a Sniffer Unit

If connecting a sniffer unit to a GC/MS system, use a flow restrictor with a 4-way valve installed and a digital flow controller (APC) to maintain a constant split ratio of components eluted from the column. The following split ratios are recommended for balancing the sensitivity of the MS unit with the flow rate of the sniffer unit.

	Internal Diameter (mm)	Length (m)	Connected to	Flow Rate Split Ratio (MS: Sniffer)
Flow Restrictor 1 (MS Flow)	0.15	1.7	MS	About 0.5:1
Flow Restrictor 2 (Sniffer Flow)	0.25	2.0	Sniffer unit	



## Appendix 8.2 Methods for when a Sniffer Unit is Connected

If a sniffer unit is connected, it applies pressure to the column outlet, which requires changing the injection port pressure setting in regular GC-MS methods. To analyze samples using a system with a sniffer unit connected, use a method file located in the folder below.


“C:\GCMSsolution\SmartDatabase\Aroma\\Sniff” folder

Column	File Name	Description
SH-I-5Sil MS	Aroma_QP/TQ_Sniff_SH-I-5SilMS_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_QP/TQ_Sniff_SH-I-5SilMS_AART.qgm	Method file for retention time correction
	Aroma_QP/TQ_Sniff_SH-I-5SilMS_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_QP/TQ_Sniff_SH-I-5SilMS_Template.qgm	Template for creating SIM method file from Smart Database
SH-PolarWax	Aroma_QP/TQ_Sniff_SH-PolarWax_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_QP/TQ_Sniff_SH-PolarWax_AART.qgm	Method file for retention time correction
	Aroma_QP/TQ_Sniff_SH-PolarWax_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_QP/TQ_Sniff_SH-PolarWax_Template.qgm	Template for creating SIM method file from Smart Database
InertCap PureWAX	Aroma_QP/TQ_Sniff_IC-WAX_Scan.qgm	Method file for analyzing aromas/odors (Scan)
	Aroma_QP/TQ_Sniff_IC-WAX_AART.qgm	Method file for retention time correction
	Aroma_QP/TQ_Sniff_IC-WAX_Correct.qgm	Method file for sensitivity correction analysis
	Aroma_QP/TQ_Sniff_IC-WAX_Template.qgm	Template for creating SIM method file from Smart Database


## Appendix 9 Adding Target Compounds

### Appendix 9.1 How to Add a Compound to a Method File (Scan)

- 1) Use a method file with adjusted retention times to measure a sample that contains the additional compound to be registered.
- 2) In the [GCMS Postrun Analysis] program, open the data file.
- 3) Click the compound table [Edit] button.



- 4) Left-click to select the last row of the compound table, where nothing is registered.
- 5) Select the peaks for the compound to register and, if necessary, subtract the background.
- 6) Left-click to select the last row of the compound table, where nothing is registered. Then right-click to click [Set Compound Information] on the right-click menu.



- 7) Select [Set name] and register the compound name.

## Appendix 9 Adding Target Compounds

Register the target and reference ion values. Generally, the ions with the highest relative intensity are registered automatically. To change a registered ion, click the [m/z] setting and then click the arrow button. That registers the mass spectrum.

Set Compound Information

Retention Time: 29.575 min  
Ret. Index: 0  
Type: Target


Target Ion:  
 TIC  MIC1  MC  
Decimal for mass: None

Compound Name  
 RT:29.575  
 Set name  
A-123456 >>

	Type	m/z	Rel Inten.
1	Target Ion	57	100.00
2	Ref. Ion	71	73.67
3	Ref. Ion	85	62.01
4	Ref. Ion	43	47.09
5	Not used	99	31.80
6	Not used	113	22.83
7	Not used	55	19.31
8	Not used	127	16.85
9	Not used	97	16.26
10	Not used	83	14.43

Buttons: Insert, Overwrite, Add To Tail, Cancel, Help

- 8) In the mass spectrum, double-click on the peaks to register the corresponding ions to be specified.



- 9) Click [Overwrite] or [Add to Tail] to register the compound.  
10) Save the method file with a new name.  
11) If necessary, overwrite the data file.


### NOTE

Compound table information is saved in the method file. If it was not saved, open the data file again and save the method file.

## Appendix 9.2 How to Add a Compound to a Database File (Scan/SIM)

Compound information registered in method files is registered in a database file in the product. That database file can be used to create method files. Therefore, newly registered compounds can be managed more conveniently by registering them in a database file.

- 1) Register the compounds using the same procedure as described in 1. How to Add a Compound to a Method File.
- 2) Open the data file and click [Quantitative Parameters] on the [Quantitative] toolbar.
- 3) Click the [Quantitative] tab and change the [Quantitative Method] setting to [External Standard].
- 4) Click the compound table [Edit] button. Then delete all compounds except the compound to be registered.



- 5) Save the method file with a different name.
- 6) Open the database file.
- 7) Click [Import] to open the method file saved in step 5). The compound is registered automatically.

**Create Method File**
Instrument Type  Lang.

**Parameter**

Ret. Index for AART

n-alkane data file  ..

Template Method File  ..

Divide Method into

## **Appendix 9.3    How to Add a Compound to a Database File (MRM)**

Register the information by following the series of steps described in 3. Registering Additional New Compound Information in a Smart Database in the GCMS Operation Guide – Method Development. After registering the information, complete the process by saving the database with a different name.

## Appendix 10 Retention Time Adjustment Using a Standard Sample

1. Use the method file with adjusted retention times to measure the standard sample.



2. Identify target compounds.

The screenshot shows the GCMS Postrun Analysis software interface. The workflow is as follows:

1. Click the GCMS Postrun Analysis icon.
2. Click the Create Compound Table icon.
3. In the File Name dialog box, select STD.qgd.
4. In the main software interface, click the Results tab. The table below shows the list of compounds and their corresponding mass chromatograms.


Ret	Name	SI	Description	Ret. Time	Ret. Index	av/c	Area	Height
1.1	1,2-Dichloroethane	50	1,2-Dichloroethane	1.155	172	74.10	2400	943
2.1	Chloroform	50	Chloroform	1.325	172	74.10	2400	943
2.2	1,1,1-Trichloroethane	50	1,1,1-Trichloroethane	1.325	172	74.10	2400	943
2.3	1,1,2-Trichloroethane	50	1,1,2-Trichloroethane	1.325	172	74.10	2400	943
2.4	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.5	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.6	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.7	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.8	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.9	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.10	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.11	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.12	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.13	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.14	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.15	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.16	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.17	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.18	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.19	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.20	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.21	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.22	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.23	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.24	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.25	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.26	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.27	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.28	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.29	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.30	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.31	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.32	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.33	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.34	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.35	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.36	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.37	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.38	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.39	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.40	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.41	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.42	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.43	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.44	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.45	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.46	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.47	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.48	1,1,1,1-Tetrachloroethane	50	1,1,1,1-Tetrachloroethane	1.325	172	74.10	2400	943
2.49	1,1,1,2-Tetrachloroethane	50	1,1,1,2-Tetrachloroethane	1.325	172	74.10	2400	943
2.50	1,1,2,2-Tetrachloroethane	50	1,1,2,2-Tetrachloroethane	1.325	172	74.10	2400	943

- 1) Start the [GCMS Postrun Analysis] program and then click the [Create Compound Table] icon on the [Postrun] assistant bar.
- 2) Open Data Explorer and click on the measurement data for the standard sample.
- 3) Click the [Results] tab below the compound table on the left.
- 4) Identify each component by checking the mass chromatogram (MC) in the quantitation view. If the components are not identified correctly, use manual peak integration.


### **Hint**

If the target compound peaks extend outside the quantitation view window, use the following procedure to specify a wider window for all compounds.

- 1) On the toolbar, click [Quantitative]-[Quantitative Parameters].





- 2) Set the [Processing Time]-[Std. Ret. Time +/- ] setting to 2 min or other suitable value and then click [OK].





- 3) That changes the window width for all compounds in the compound table based on the specified time range.

3. Use the [Wizard (Modify)] to apply the retention times from the identification results so they are reflected in the compound table.

1 

2 

3 

4 

- 1) Click the [Wizard (Modify)] icon on the [Compound] assistant bar.
- 2) Clear the indicated checkboxes.
- 3) Click [Next].
- 4) Click [Next].
- 5) Click [Finish].

 **NOTE**

About the [Modify Compound Table (Assist Reference Ion Ratio)] Window (1 of 3)  
Normally, these checkboxes are cleared in step 2. Select these checkboxes if the reference ion ratios specified in the compound table need to be reset based on the data currently loaded.

Reference Ion



Set Reference Ion ratios to existing reference ions (where no ratio was specified before)

Overwrite existing Reference Ion ratios

Decimal for mass:

None

## 4. Saving the Compound Table




- 1) Click the [View] icon located on the right side above the compound table.
- 2) Click the [Edit] tab on the left side below the compound table.  
That allows checking that the retention times are adjusted to the retention times in identification results.
- 3) Click the [Save Compound Table] icon on the [Compound] assistant bar.
- 4) Select the method file used for analysis.
- 5) Click [Save].

## Appendix 11 Outputting Files for Statistical Data Analysis Software


To use statistical data analysis software for multivariate analysis or other statistical analysis, the area and area ratio values for respective peaks in multiple sets of data must be output in list form. Such data can be output as a CSV file using the following procedure in LabSolutions Insight.

### 1. Load the data to output.


- 1) Click [File]-[Open]  on the menu bar to open the method file to be used for data analysis.
- 2) Click [File]-[Import]  on the menu bar to load the data file to be output.
- 3) Check the data analysis results.
- 4) On the [View] menu, specify the [Summary] view. Right-click the caption in the summary pane and select [Properties].



- 5) In the [Summary Results View Style] window, click the [Comparison] icon.




- 6) Use table style settings to specify the data to be output (Area, etc.) in the [Selected Columns] field. Then click the  (OK) icon.




## 2. Output an Export File.

- 1) Click [File]-[Export] on the menu bar and then select [Results to File].
- 2) Select the [Output to File] checkbox and then click the [...] (browse) button.



- 3) Specify the file name and click [Save as type]-[Comma-separated Values (\*.csv)].



- 4) Click [Save]. That saves the data as a CSV file.

## Appendix 12 Installing the Database

Follow the steps below to install the database.

The explanations given here assume Windows 10 is used.



### NOTE

GCMSsolution must already be installed before the database can be installed.


- 1) Insert the Smart Aroma Database installation disk into the CD-ROM drive. The installer program launches automatically.



### NOTE

The dialog box "Do you want to allow the following program to make changes to this computer?" may be displayed. In this case, click the [Yes] to proceed.


- 2) Click [Next>].




## Appendix 12 Installing the Database

---

- 3) Select the MS model to be used. Click [Next>].



- 4) In the [Ready to Install the Program] window, click [Install] to start installing the various database files in the specified location.



Method files and database files are normally installed in the "C:\GCMSsolution\SmartDatabase" folder.

## Appendix 12 Installing the Database

---

- 5) When installation is finished, the [InstallShield Wizard Completed] window is displayed. Click [Finish] to complete the process.

