

GCMS-TQ NX Series

GCMSsolution

Instruction Manual

Read this manual thoroughly before you use the product.
Keep this manual for future reference.

This page is intentionally left blank.

Introduction

Read this Instruction Manual thoroughly before using the product.

Thank you for purchasing Shimadzu analytical instrument workstation “GCMSsolution” (hereafter referred to as “the software” or “GCMSsolution”).

This manual describes the procedures for operating this product. Read this manual thoroughly before using the product and operate the product in accordance with the instructions in this manual.

Keep this manual for future reference.

This manual assumes that the reader is knowledgeable of basic operations of Windows. For the operation of Windows, refer to the instruction manual that comes with that product.

For information about how to use the GCMS-TQ NX series gas chromatograph mass spectrometer, its operating precautions, and so on, refer to the Instruction Manual for each instrument.

Important

- If the user or installation location changes, ensure that this Instruction Manual is transferred with the product.
- If this manual is lost or damaged, immediately contact your Shimadzu representative to request a replacement.
- To ensure safe operation, contact your Shimadzu representative for product installation, adjustment, or reinstallation (after the product is moved).

© 2018 Shimadzu Corporation All rights reserved.
Original version is approved in English.

Notice

- Information in this manual is subject to change without notice and does not represent a commitment on the part of the vendor.
- Any errors or omissions which may have occurred in this manual despite the utmost care taken in its production will be corrected as soon as possible, although not necessarily immediately after detection.
- All rights are reserved, including those to reproduce this manual or parts thereof in any form without permission in writing from Shimadzu Corporation.
- Microsoft, Windows, and Excel are registered trademarks of Microsoft Corporation in the United States and/or other countries.
Adobe, Adobe logo and Adobe Acrobat, Adobe Acrobat Reader, Acrobat Distiller are trademarks or registered trademarks of Adobe Systems Incorporated in the United States and/or other countries.
Other company names and product names mentioned in this manual are trademarks or registered trademarks of their respective companies.
The TM and ® symbols are omitted in this manual.
- Microsoft® Windows® 7 Operating System is referred to as “Windows 7”.
Microsoft® Windows® 10 Operating System is referred to as “Windows 10”.
- Microsoft® Office Excel® 2016 is referred to as “Excel”.
- GCMS-TQ and GCMSsolution are registered trademarks of Shimadzu Corporation in the United States and Europe (Office for Harmonization in the Internal Market).
- Replacement parts for this product will be available for a period of seven (7) years after the product is discontinued. Thereafter, such parts may cease to be available. Note, however, that the availability of parts not manufactured by Shimadzu shall be determined by the relevant manufacturers.

Indications Used in This Manual

This instruction manual uses the notation described below.

Notation	Meaning
 CAUTION	Indicates a potentially hazardous situation which, if not avoided, may result in minor to moderate injury or equipment damage.
 NOTE	Indicates additional information that is provided to ensure the proper use of this product.
 Reference	Indicates the location of related information.
 Hint	Indicates information provided to improve product performance.
[]	Indicates items displayed on the screen, such as buttons, menu selections, settings, windows, and icons. Example: Click [OK].

Product Warranty

Shimadzu Corporation provides the following warranty for this product.

Details

1. Period: Please contact your Shimadzu representative for information about the period of this warranty.
2. Description: If a product/part failure occurs for reasons attributable to Shimadzu during the warranty period, Shimadzu will repair or replace the product/part free of charge. However, in the case of products which are usually available on the market only for a short time, such as personal computers and their peripherals/parts, Shimadzu may not be able to provide identical replacement products.
3. Limitation of Liability
 - 1) In no event will Shimadzu be liable for any lost revenue, profit or data, or for special, indirect, consequential, incidental or punitive damages, however caused regardless of the theory of liability, arising out of or related to the use of or inability to use the product, even if Shimadzu has been advised of the possibility of such damage.
 - 2) In no event will Shimadzu's liability to you, whether in contract, tort (including negligence), or otherwise, exceed the amount you paid for the product.
4. Exceptions: Failures caused by the following are excluded from the warranty, even if they occur during the warranty period.
 - 1) Improper product handling
 - 2) Repairs or modifications performed by parties other than Shimadzu or Shimadzu designated companies
 - 3) Product use in combination with hardware or software other than that designated by Shimadzu
 - 4) Computer viruses leading to device failures and damage to data and software, including the product's basic software
 - 5) Power failures, including power outages and sudden voltage drops, leading to device failures and damage to data and software, including the product's basic software
 - 6) Turning OFF the product without following the proper shutdown procedure leading to device failures and damage to data and software, including the product's basic software
 - 7) Reasons unrelated to the product itself
 - 8) Product use in harsh environments, such as those subject to high temperature or humidity levels, corrosive gasses, or strong vibrations
 - 9) Fires, earthquakes or any other act of nature, contamination by radioactive or hazardous substances, or any other force majeure event, including wars, riots, and crimes
 - 10) Product movement or transportation after installation
 - 11) Consumable items
Note: Recording media, such as floppy disks and CD-ROMs are considered consumable items.

* If there is a document such as a warranty provided with the product, or there is a separate contact agreed upon that includes warranty conditions, the provisions of those documents shall apply.

The warranty period for products with special specifications or for system products is specified separately.

Software License Agreement

PLEASE READ THIS AGREEMENT CAREFULLY BEFORE USING THE SOFTWARE.

SHIMADZU Corporation ("SMZ") is willing to license the SMZ software provided herein, together with accompanying documentation (collectively "SOFTWARE") to you only upon the condition that you accept all of the terms and condition contained in this License Agreement. By using the SOFTWARE, you agree to be bound by the terms of this Agreement. If you do not agree to all these terms of this Agreement, promptly return the unused SOFTWARE to the party (either SMZ or its reseller) from whom you acquired it to receive a refund of the amount you paid.

1. License: SMZ grants you a non-exclusive and nontransferable license to use the SOFTWARE subject to the following terms and conditions.
2. Limitation of Use: Except as specifically authorized by SMZ, you may NOT:
 - 1) Use the SOFTWARE, or permit the SOFTWARE to be used, on more than one computer at any one time.
 - 2) Copy the SOFTWARE except one (1) archival copy of the SOFTWARE.
 - 3) Modify, reverse engineer, decompile, disassemble, or create derivative works based upon the SOFTWARE.
 - 4) Transfer, rent, lease or grant any rights in the SOFTWARE in any form to anyone else.
3. Title and Ownership: This license is not for sale and it may not be assigned or sublicensed to anyone else. Title and all associated intellectual property rights to the SOFTWARE shall remain in SMZ and/or its licensor.
4. Upgrades: You are entitled to receive all official software upgrades for the SOFTWARE that SMZ will release as deemed necessary by SMZ. An upgrade means certain supplemental program modules and/or information for bug fixing and/or updates to the defects and/or failures of the SOFTWARE that are acknowledged or confirmed by SMZ. An upgrade excludes hardware, network, consulting services, third party products, operation and general computer system maintenance. All supplemental program module for upgrades and enhancements furnished to you shall be deemed to be part of the SOFTWARE and subject to the terms and conditions set forth in this Agreement.
5. Limited Warranty: SMZ warrants that for a period of one (1) year from the date of purchase, as evidenced by a copy of the receipt, the media on which SOFTWARE is furnished will be free of defects in materials and workmanship under normal use.

Except for the foregoing, SOFTWARE is provided AS IS. Your exclusive remedy and the entire liability of SMZ and its suppliers under this limited warranty will be, at SMZ's option, repair, replacement, or refund of the Software if reported (or, upon request, returned) to the party supplying the SOFTWARE to you. In no event does SMZ warrant that the Software is error free or that you will be able to operate the SOFTWARE without problems or interruptions.

EXCEPT AS SPECIFIED IN THIS WARRANTY, ALL EXPRESS OR IMPLIED CONDITIONS, REPRESENTATIONS, AND WARRANTIES INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, NONINFRINGEMENT, ARISING FROM A xvi Introduction COURSE OF DEALING, USAGE, OR TRADE PRACTICE, ARE HEREBY EXCLUDED TO THE EXTENT ALLOWED BY APPLICABLE LAW.
6. Limitation of Liability: IN NO EVENT WILL SMZ BE LIABLE FOR ANY LOST REVENUE, PROFIT OR DATA, OR FOR SPECIAL, INDIRECT, CONSEQUENTIAL, INCIDENTAL OR PUNITIVE DAMAGES, HOWEVER CAUSED REGARDLESS OF THE THEORY OF LIABILITY, ARISING OUT OF OR RELATED TO THE USE OF OR INABILITY TO USE SOFTWARE, EVEN IF SMZ HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGE. IN NO EVENT WILL SMZ'S LIABILITY TO YOU, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE), OR OTHERWISE, EXCEED THE AMOUNT YOU PAID FOR SOFTWARE.

- 
7. Termination: This License is effective until terminated. You may terminate this License at any time by destroying all copies of SOFTWARE including any documentation. This License will terminate immediately without notice from SMZ if you fail to comply with any provision of this License. Upon termination, you must destroy all copies of SOFTWARE.
8. General:
- 1) This Agreement is the entire agreement. If any provision of this agreement is held invalid, the remainder of this agreement shall continue in full force and effect.
 - 2) This Agreement shall be construed and governed in accordance with the laws of Japan, excluding its conflict of law rules.
 - 3) The exclusive jurisdiction for any disputes arising out of or in connection with this Agreement shall be Kyoto District Court of Japan.
 - 4) The invalidity or unenforceability of any provision of this Agreement shall not affect the validity or enforceability of any other provision.



Contents

1 Description

1.1	Features.....	1
1.2	Common Operations.....	2
1.2.1	Starting Up the Program	2
1.2.2	Explanation of Window Parts	3
1.2.3	How to Use Data Explorer.....	4
1.2.4	Opening Help	5
1.2.5	Using Online Manuals	6
1.2.6	File Formats	7
1.3	Operation Flow Chart.....	8

2 Before Starting Data Acquisition

2.1	System Configuration	11
2.1.1	Opening the [System Configuration] Sub-Window	11
2.1.2	Checking GC Communication Settings	12
2.1.3	Registering Instruments	13
2.1.4	Changing CID Gas Settings	16
2.1.5	Setting Criteria for Auto-Tuning Result Judgment.....	18
2.2	Executing a System Check.....	20
2.2.1	Opening the [System Check] Sub-Window	20
2.2.2	Executing a System Check	22
2.3	Executing Auto-Tuning	23
2.3.1	Executing Auto-Tuning.....	24
2.3.2	Checking Tuning Results	26
2.3.3	Printing Tuning Reports	28

3 Data Acquisition

3.1	Flow of Data Acquisition	29
3.2	[Acquisition] Window.....	29
3.2.1	Opening the [Acquisition] Window.....	29
3.3	Setting Data Acquisition Methods	31
3.3.1	Creating New Method Files.....	31
3.3.2	AOC Instrument Parameters.....	32

3.3.3	GC Instrument Parameters	33
3.3.4	MS Instrument Parameters	36
	3.3.4.1 Creating MS Tables Based on Compound Tables	40
	3.3.4.2 Creating FASST Measurement Methods	43
3.3.5	Saving Method Files	45
3.4	Starting Single Run.....	46
3.4.1	Registering Samples	46
3.4.2	Preparations	47
3.4.3	Starting Single Run	47
3.4.4	Stopping Single Run	48
3.4.5	Extending the Data Acquisition Time	48
3.4.6	Snapshot	49
3.5	Monitoring Instrument Status.....	50
3.5.1	Instrument Monitor	50
3.6	Ecology Mode.....	51
3.6.1	Setting the Ecology Mode	51
3.6.2	Entering the Ecology Mode	53
3.6.3	Canceling the Ecology Mode	54
3.6.4	Setting the PC or Display to the Power Save Mode.....	54

4 Continuous Data Acquisition

4.1	Opening the [Batch Table] Window	55
4.2	Creating Batch Tables	55
4.2.1	Batch Table Wizard	55
4.2.2	Editing Batch Tables	59
4.3	Starting Continuous Data Acquisition	60
4.3.1	Starting Continuous Data Acquisition.....	60
4.3.2	Pausing Continuous Data Acquisition	61
4.3.3	Stopping Continuous Data Acquisition	61
4.3.4	Partially Executing Continuous Data Acquisition	62
4.4	Using Batch Queue to Perform Additional Data Acquisition.....	62
4.4.1	Creating Additional Batch Files	62
4.4.2	Adding Batch Files	63
4.5	Displaying Remaining Time of Continuous Data Acquisition.....	65
4.6	Entering the Ecology Mode after Continuous Data Acquisition	66
4.7	Errors During Batch Analysis.....	67

5 Postrun Analysis

5.1	[Data Analysis] Window	71
5.1.1	Opening the [Data Analysis] Window	71
5.2	Displaying Chromatograms and Spectra	72
5.2.1	Chromatogram View.....	74
5.2.2	Spectrum View	75
5.2.3	Quantitative View	77
5.2.4	Displaying the Spectrum of a Specific Retention Time	78
5.2.5	Displaying the Chromatogram of a Specific m/z	78
5.2.6	Displaying the Standard Spectrum in [Spectrum View].....	79
5.2.7	Comparing the Quantitative Chromatogram with Other Measurement Data....	80
5.2.8	Comparing Chromatograms with the Intensity Axis of the Quantitative Chromatogram Fixed.....	80
5.3	Qualitative Processing	82
5.3.1	Flow of Qualitative Processing	82
5.3.2	Performing Qualitative Peak Integration Automatically	82
5.3.3	Checking the Peak Integration Results	85
5.3.4	Performing Qualitative Peak Integration Manually	85
5.3.5	Performing Average/Subtraction Processing on Spectra.....	89
5.3.6	Searching Spectra from the Library.....	94
5.3.7	Printing Graph Images	97
5.4	Quantitative Processing	99
5.4.1	Flow of Quantitative Processing.....	99
5.4.2	Creating Compound Tables	99
5.4.2.1	Creating Compound Tables Using the Wizard	99
5.4.2.2	Creating Compound Tables from MS Tables in the Method.....	105
5.4.2.3	Copying from Spreadsheet Format Files to Create Compound Tables.....	108
5.4.3	Creating Calibration Curves	111
5.4.4	Performing Quantitative Processing on Unknown Samples.....	115
5.4.5	Checking the Quantitative Calculation Results	116
5.4.6	Performing Quantitative Peak Integration Manually	117
5.4.7	Performing Identification Manually	118
5.4.8	Detecting Unidentified Peaks	118
5.4.9	Calculating S/N	120
5.5	Saving New Parameters to Method Files	122
5.6	Printing the Content of Data Files	123
5.7	Performing Continuous Postrun Analysis	124

6 Quant Browser

6.1	[Quant Browser] Window	127
6.1.1	Opening the [Quant Browser] Window	127
6.2	Checking the Quantitative Results	128
6.2.1	Opening Method Files	128
6.2.2	Opening Data Files to Display Quantitative Results	129
6.2.3	Opening Batch Files to Display Quantitative Results	130
6.2.4	Editing [Quantitative Result View]	132
6.2.5	Fixing the Intensity Axis of [Chromatogram View] to Check Quantitative Results	133
6.3	Performing Postrun Analysis Collectively on Multiple Data Files	135
6.3.1	Editing Compound Tables	135
6.3.2	Setting the Peak Integration Parameters	135
6.3.3	Performing Peak Integration Collectively on Multiple Data	136
6.4	Correcting Calibration Curves	136
6.5	Printing Summary Reports	137
6.6	Printing Chromatogram Images	137

7 Data Comparison

7.1	Opening the [Data Comparison] Window	139
7.2	Displaying Chromatograms and Spectra	140
7.2.1	Displaying the Spectrum of a Specific Retention Time	141
7.2.2	Comparing Spectra with the Intensity Axis of the Spectra Fixed	141
7.2.3	Displaying the Chromatogram of a Specific m/z	142
7.2.4	Moving Chromatograms	143
7.2.5	Expanding/Reducing Chromatograms	144
7.2.6	Performing Calculations on Chromatograms and Spectra	144
7.3	Printing Graph Images	146

8 Reports

8.1	Opening the [Report] Window.....	147
8.2	Creating Report Formats	148
8.2.1	Manually Setting Report Formats.....	148
8.2.2	Using a Template to Set a Report Format.....	152
8.3	Printing Reports	154
8.3.1	Printing Data Processing Results in the [Data Analysis] Window	154
8.3.2	Printing Reports at Batch Processing.....	156

Appendix A Maintenance

A.1	Installing the Software	157
A.1.1	Preparations for Installation	157
A.1.2	Installing GCMSsolution	157
A.2	Uninstalling the Software	163
A.3	Connecting to the Instrument	164
A.3.1	Connect the PC to the MS Unit	164
A.3.2	Connect the PC to the GC Unit	164
A.3.3	Installation of the Driver for MS or GC	169
A.4	Easy sTop	172
A.5	Resetting the Number of Uses/Operation Time	175
A.6	Changing the Replacement Guidelines for the Septa and Glass Inserts	176
A.7	Column Maintenance	178
A.8	Errors and Remedies	179
A.8.1	The Instrument Is Not Properly Recognized	179
A.8.2	Messages Displayed at Startup	179
A.8.3	Windows Does Not Start Up	179

Appendix B Before PDF Output

B.1	Before PDF Output	181
B.2	Installing the Postscript Printer Driver	181
B.3	Setting PDF Output	187

Appendix C Exporting Files

C.1	Exporting from the [Data Explorer] Sub-Window	189
C.2	Exporting Data to AIA Files from the [Data Analysis] Window	190
C.3	Exporting Data by Batch Processing	191
C.4	Exporting from the [Library Editor] Window	192
C.5	Exporting Batch Files from the [Batch Table] Window	194

Appendix D Specifications

D.1	Workstation	197
------------	-------------------	-----

Index	199
--------------------	------------

1

Description

1.1.6.1

The GCMSsolution software for gas chromatograph mass spectrometers runs on Windows.

This software controls the GCMS-TQ NX series and peripherals such as the autosampler and acquires data. It also analyzes acquired data and prints out analysis results.

1

1.1 Features

This software has many functions for acquiring and analyzing data by using the GCMS-TQ NX series gas chromatograph mass spectrometer. The following introduces the main features of this software.

1. Automatic Adjustment of Retention Time Function (AART)

Even if the retention time of a compound has fluctuated, the software automatically corrects for this simply by performing data acquisition on the reference sample (n-alkane).

2. Automatic Scan/SIM Parameter Creation Function (COAST)

With the COAST function, method files for MRM, simultaneous scan/MRM measurements and product ion scan can be created automatically from the Compound Table.

3. Compound Identification Support Function

This function helps identify compounds quickly and accurately by displaying chromatograms from multiple sets of data, such as data from samples with standard substances added, on the same window for easy comparison.

4. Quantitation Browser with Multianalyte Quantitation Function

This browser allows easily loading large quantities of analytical results and performing quantitative calculations for all the data simultaneously or adjusting or confirming chromatograms for identified components. It also allows printing area values, concentration values, and other quantitation results, which can improve the efficiency of quality control processes.

5. Similarity Search Function Based on Retention Indices

In addition to searches based on mass spectra similarity, searches can be further refined based on retention indices. This ensures components can be identified even among compounds with similar mass spectra.

6. Highly Flexible Report Customization

Chromatograms, spectra, quantitation results, or a variety of other report items can be freely laid out in any location or size. This allows creating a report format that is tailored to what you want.

7. GLP/GMP Support

The software is fully equipped with security, audit trail, hardware/software verification and QA/QC functions for GLP/GMP compliance.

8 Ecology Mode

Implementation of the ecology mode means that power and carrier gas consumption can be reduced when the instrument is standing by for data acquisition.

1.2 Common Operations

This software comprises the following programs:

Program Name	Contents
[GCMS Realtime Analysis] program	Performs instrument start/stop, system configuration and data acquisition.
[GCMS Analysis Editor] program	Creates and edits method files or batch files during data acquisition.
[GCMS Postrun Analysis] program	Performs data processing such as qualitative processing, quantitative processing, and printing of reports.
[GCMS Browser] program	Performs data processing on multiple data, such as qualitative processing, quantitative processing, and printing of reports.

1.2.1 Starting Up the Program

- 1 Double-click a GCMSsolution program icon  (e.g. (GCMS Real Time Analysis)) on the desktop.

The [Login] window opens.

- 2 Select the user ID from the [User ID] list, and enter the password.




 **NOTE**

When logging in for the first time, select [Admin] in the [User ID] list and leave [Password] blank.

- 3 Click [OK].
[GCMS Real Time Analysis] program starts up.

1.2.2 Explanation of Window Parts



No.	Name	Explanation
①	Title Bar	Displays information such as name of the currently running program, window name, currently loaded file name, and logged in user name.
②	(Close) button	Clicking this button exits the program.
③	Menu Bar	Displays the menus that are enabled for the currently used window and the rights of the currently logged in user.
④	Tool Bar	Displays the icons of frequently used menu items and the icons for operating analytical instruments.
⑤	Window	Different windows such as [Acquisition] and [Batch Table] are displayed in this frame of the program. The currently displayed window is switched by clicking icons on the assistant bar or tabs displayed under each window.
⑥	Output window	Displays an operation history of data acquisition and error messages that occur.
⑦	Data Explorer	Displays files in the currently selected folder according to the type of file selected by the tabs underneath. The content of files is displayed by dragging-and-dropping the file in the [Data Explorer] sub-window onto the window.
⑧	Assistant Bar	Icons corresponding to commands are arranged according to the general order of operations. Normally, click an icon on this bar to select a respective operation.

1.2.3 How to Use Data Explorer


GCMSsolution displays a list of files stored to specified projects (folders) in the [Data Explorer] sub-window.

Files displayed in the [Data Explorer] sub-window can be moved or copied to other folders, or deleted. Tabs are also provided for each file type so that files of the type according to the clicked tab can be extracted and displayed.



NOTE

When the [Data Explorer] sub-window is not displayed, click [Data Explorer] on the [View] menu.



No.	Name	Explanation
①	 (Create New Project) button	Creates a new project at the same level of the folder hierarchy as the currently referenced project (folder). Files can be also copied from the current project when creating new projects.
②	 (Select Project) button	Changes the currently referenced project.
③	Project in	Displays the path of the currently referenced project. The paths of up to ten recently referenced projects are stored in the combo box list. The reference project displayed in Data Explorer can be selected and changed from this list.
④	(File View)	Displays a list of files contained in the currently referenced project. Double-clicking a file in this list starts up the specific program related to that file to load the file. Files can be loaded by dragging-and-dropping onto specific windows.
⑤	(Preview)	<p>Previews the TIC of a file selected in the file view.</p> <p> NOTE When the preview is not displayed, select [Data Preview] on the menu displayed by right-clicking the file view.</p>
⑥	(File Type Tab)	Extracts and displays files of the same type as the selected tab from the files contained in the currently referenced project.

1.2.4 Opening Help


Display and refer to Help when you are not sure of software operations or the terms displayed in a window.

Help	Click [Help] in the sub-window. The Help for the currently open sub-window is displayed.
 (Help)	Click  (Help) on the toolbar. [GCMSsolution Help] opens.
[Help] menu	Click [Contents] on the [Help] menu. [GCMSsolution Help] opens.
[F1] key	Press the [F1] key on your keyboard. The Help for the currently open sub-window is displayed.

1

Keyword Search


- 1 Open Help.
- 2 Click the [Index] tab.



NOTE

Keywords can be searched from the entire text of the Help file on the [Index] tab page.

3 Select the keyword in [Title] in the sub-window, and click [Display].



The content of the selected topic is displayed.

NOTE


When there are multiple matching keywords in multiple topics, the [Topic Found] sub-window opens. Select the matching keyword in the list of titles in this sub-window, and click the [Display] button.

1.2.5 Using Online Manuals

When this software is installed, this Instruction Manual, System User's Guide and other online manuals are also installed in PDF format.

1 Double-click the desired instruction manual in the Manual folder in the folder where this software is installed.

Adobe Acrobat Reader starts up and the instruction manual opens.



No.	Explanation
1	You can jump to the desired page by clicking the hierarchically structured bookmarks (table of contents).

No.	Explanation
②	You can search for terms that you want to look up.
③	You can jump to the page of the related item by clicking the references or the terms in blue.

1.2.6 File Formats

This software handles files in the following formats.

File type	Icon	Extension	File contents
Data file		.qgd	In addition to the raw data acquired (e.g., chromatograms and spectra), the following information is saved. <ul style="list-style-type: none"> • Calculation results such as area values and concentrations • Status information such as the oven temperature and error status at the time data is acquired • Contents of method files used in analysis (including configuration settings used for analysis) • Contents of report format file (when reports are output) • Contents of batch files (when batch processing is performed) • Contents of tuning file used in analysis
Method file		.qgm	Analysis conditions, peak integration parameters, compound tables, etc. are saved. Because the configuration settings are saved when the method is edited, the configuration settings are checked when the method file is loaded to ensure that they agree with the current settings. Created calibration curves are also saved in the method file.  NOTE If a method file created in GCMSsolution Ver. 2 is opened in GCMSsolution Ver. 4, SIM events are loaded as Q3SIM events and scan events are loaded as Q3 scan events.
Report format file		.qgr	The report format information used to output a report, such as layout information and detailed settings, is saved. Once a report format file has been created, it can be used repeatedly to output reports of the same format.
Batch file		.qgb	Batch tables used to perform automatic sequential processing are saved. The same files can be used in both the [GCMS Real Time Analysis] program and the [GCMS Postrun Analysis] program.
Tuning file		.qgt	The conditions used to perform instrument adjustment (tuning) and the tuning results are saved.  NOTE Tuning files created in GCMSsolution Ver. 2 cannot be opened in GCMSsolution Ver. 4 or used for data acquisition.
Library file		.lib *1	These files are used to register the compound information and spectral data used to perform similarity searches. The libraries consist of public libraries (e.g., NIST and Wiley) and private libraries.
Browsing file	None	.qqq	These files store information such as compound information displayed in [Quantitative Results View] and the names of method and data files loaded in the [Quant Browser] window.
Layout file	None	.lyt	These files store information such as data file names and display layouts loaded in [Data Browser].


*1 The library is actually composed of multiple files. The constitutive files are somewhat different for commercial and private libraries. When copying library files, ensure that all associated files are copied.
Commercial library: .lib, .c2s, .fom, .nam, .ncv, .spc, .str
Private library: .lib, .com, .flg, .fom, .nam, .spc, .add

Reference

For details, refer to the instruction manual for the library.

1.3 Operation Flow Chart

This section describes a simple operation procedure in GCMSsolution in the form of a flow chart.



Quantitative Analysis



1 Create method file

- Set instrument parameters. [P.31](#)
- Single run using a standard sample [P.46](#)



2 Create the Compound Table

- Create the Compound Table using standard sample data
- Create the Compound Table for the MRM/SIM acquisition mode. [P.108](#)
 - Create the Compound Table for the scan acquisition mode. [P.99](#)
- Create the Compound Table using the Compound Table Wizard. [P.99](#)



3 Continuous data acquisition

- Perform continuous data acquisition on standard samples and unknown samples of each concentration for calibration curves using the method file used for creating the Compound Table.
- Create the Batch Table. [P.55](#)
 - Start continuous data acquisition. [P.60](#)



4 Quantitation Browser

- Check/correct calibration curves. [P.136](#)
- Quantitative calculation of multiple unknown concentration samples [P.135](#)



5 Reports

- Print graph image. [P.137](#)



4 Data analysis

- Check/correct calibration curves. [P.111](#)
- Quantitative calculation of unknown concentration samples. [P.115](#)



5 Reports

- Print data reports. [P.123](#)

This page is intentionally left blank.

2

Before Starting Data Acquisition

This chapter describes preparations for using the instrument before actually starting data acquisition.

2


2.1 System Configuration

Register the instrument configuration of the GC to be used for data acquisition together with the autosampler, injection unit, column, detector, and other units. Registered details are saved to a system configuration file unique to each system, and need not be set each time that GCMSsolution is started up. The system configuration file must be set again when the column has been replaced or the instrument configuration has been changed.

2.1.1 Opening the [System Configuration] Sub-Window

- 1 Start up the [GCMS Real Time Analysis] program, and click the  (System Configuration) icon on the [Real Time] assistant bar.

The [System Configuration] sub-window opens.




No.	Name	Explanation
1	Available Modules	Available instruments are registered. The model name of autosamplers, injection units, detectors, additional heaters, and additional flows are displayed by clicking  on the left.
2	Modules Used for Analysis	Displays registered instrument configurations for each system.
3		Moves the instrument selected in the [Available Modules] list to the [Modules Used for Analysis] list, and adds the instrument as the instrument to be used for data acquisition.
4		Moves the instrument selected in the [Modules Used for Analysis] list to the [Available Modules] list, and deletes the instrument from the instrument to be used for data acquisition.
5	Properties	A detailed setting sub-window of the instrument selected at [Modules Used for Analysis] opens.

No.	Name	Explanation
6	Print	Prints the information of the set instrument configuration and its detailed settings.
7	Set	Applies the information of the set instrument configuration and its detailed settings, and closes the [System Configuration] sub-window.
8	Audit Trail	<p>Clicking this button opens the [Audit Trail] sub-window. The [Audit Trail] sub-window displays the audit trail results until the [Set] button was previously clicked.</p> <p>To enable the audit trail of the system configuration, click [Audit Trail Settings For Config File] on the [Tools] menu in the [GCMS Real Time Analysis] program, and select [Input reason when saving Config] in the sub-window that is displayed.</p> <p> NOTE Once a history is set to be left behind, this setting cannot be canceled.</p>

2.1.2 Checking GC Communication Settings

1 Double-click [Instrument 1] at [Modules Used for Analysis].




The [Instrument] sub-window opens.

2 Check the PC-GC communication settings.

When [USB] is selected, enter the serial number of GC in [Serial Number].
When [Ethernet] is selected, enter the IP address of GC in [IP Address].
(The figure below shows the setting example for USB.)

Reference


See "A.3.2 Connect the PC to the GC Unit".



2


2.1.3 Registering Instruments

1 Select the instruments to be used for data acquisition at [Available Modules], and click



The instruments are registered to [Modules Used for Analysis].

2 Double-click the instrument whose settings are to be changed at [Modules Used for Analysis].




The sub-window for setting the details of the instrument opens.

3 Set each parameter, and click [OK].

Almost all parameters can be used at their default settings.

[Column] Tab Page

Register several types of columns, and select which one is to be used. Register the information of the column to be used, and be sure to select that column.



[MS] Tab Page

Make sure that the type of MS detector ion source to be used is set correctly.

When using an MS/MS acquisition mode, such as MRM, Product Ion Scan, Precursor Ion Scan, or Neutral Loss Scan, select the [Use CID Gas] checkbox.

When changing the CID gas operating pressure, set the pressure in the [CID Gas (Initial value of tuning)] field.

When using the instrument as a single GCMS without using an MS/MS acquisition mode, deselect the [Use CID Gas] checkbox. See ["2.1.4 Changing CID Gas Settings"](#).


To automatically judge auto-tuning results, select [Check Auto Tuning Result].

 Reference

See ["2.1.5 Setting Criteria for Auto-Tuning Result Judgment"](#).

**NOTE**

This setting is available only when the ion source of the MS detector is EI.


**NOTE**

When the CID gas operating pressure or CID gas use/non-use setting has been changed, be sure to execute auto-tuning before performing data acquisition.

 Reference

See ["2.3.1 Executing Auto-Tuning"](#).

4

Repeat steps 2 and 3 to set other instruments.

5


When registration is over, click [Set] to close the [System Configuration] sub-window.

2.1.4 Changing CID Gas Settings


When setting whether or not to use CID gas in data acquisition, follow the procedure below. When CID gas is set not to be used, data acquisition can be performed in the Q3 scan or Q3 SIM acquisition mode with the instrument as a single GCMS.

1 In the [System Configuration] sub-window, set so that CID gas is not used.

- 1 Start up the [GCMS Real Time Analysis] program, and click the  (System Configuration) icon on the [Real Time] assistant bar. The [System Configuration] sub-window opens.
- 2 Double-click [MS] at [Modules Used for Analysis].



- 3 When not using CID gas, deselect the [Use CID Gas] checkbox.



- 4 Click [OK] to close the [System Configuration] sub-window. CID gas control is turned OFF, and CID gas no longer flows.

2 Create a new tuning file, and execute auto-tuning.



NOTE

When the CID gas operating pressure or CID gas use/non-use setting has been changed, be sure to execute auto-tuning before performing data acquisition.



Reference

See "[2.3.1 Executing Auto-Tuning](#)".

3

Start data acquisition.

Specify the tuning file obtained by performing auto-tuning by the new CID gas setting, and execute data acquisition.


2




NOTE

When CID gas is set not to be used, the following message is sometimes displayed when preparing for data acquisition by single run or when continuous data acquisition is started.

- When data acquisition is started using a method file set with events other than Q3 scan or Q3 SIM scan to MS parameters



- When data acquisition is started using a tuning file that has not been tuned after use/non-use of CID gas has been changed in the system configuration



Reference

Clicking [Help] in the error sub-window displays the error message Help. Remove the cause of the error by referring to causes and remedies.


2.1.5 Setting Criteria for Auto-Tuning Result Judgment

When the ion source of the MS detector is EI, auto-tuning results can be judged automatically. Configure the settings using the procedure below.


- 1 Start up the [GCMS Real Time Analysis] program, and click the  (System Configuration) icon on the [Real Time] assistant bar.

The [System Configuration] sub-window opens.

- 2 Double-click [MS] at [Modules Used for Analysis].



- 3 To automatically judge auto-tuning results, select [Check Auto Tuning Result].



4

1.1.6.2

Click [Test Values of System Check/Auto Tuning] and open the [Auto Tuning] tab page.
Set the criteria for judging auto-tuning results.

2

No.	Name	Explanation	Default
①	FWHM	Specifies the upper limit of the shift of the FWHM of spectrum peaks at m/z 69, 219, and 502 from the set value.	0.1
②	Detector Gain	Specifies the upper limit of the detector voltage set by auto-tuning. NOTE Although the detector voltage on the actual instrument is a negative value, the software handles the parameter as a positive value.	2
③	Vacuum Leak (Intensity Ratio of 69/28)	Specifies the lower limit of the intensity ratio between m/z 69 (PFTBA) and m/z 28 (Nitrogen). NOTE <ul style="list-style-type: none"> The reference value of the criteria differs depending on the pretreatment unit. Set the criteria in accordance with the pretreatment unit by following the values below. <ul style="list-style-type: none"> Pyrolyzer, OPTIC-4: intensity ratio 1 Other pretreatment units: intensity ratio 2 When using the advanced flow technology system (such as detector splitting system or backflush), the intensity for m/z 28 may be larger than when using the standard system. Therefore, it is not recommended to select [Vacuum Leak (Intensity Ratio of 69/28)] for check. 	2
④	Maximum Mass Shift	Specifies the upper limit of the mass axis shift.	0.1
⑤	Relative Intensity Ratio for High m/z (502)	Specifies the lower limit of the relative intensity for m/z 69 in the high mass range (m/z 502).	2

5

Click [OK] to close the [System Configuration] sub-window.


Refer to "2.3 Executing Auto-Tuning" to execute auto-tuning.

2.2 Executing a System Check


The system check judges whether or not there are problems in the instrument status. It is recommended that a system check is executed before the series of data acquisition operations are begun.

2.2.1 Opening the [System Check] Sub-Window

- 1 Start up the [GCMS Real Time Analysis] program, and click [System Check] on the [Instrument] menu.



The [System Check] sub-window opens.




No.	Name	Explanation
①	Check	Select the items to be checked. The operation time/number of uses of consumables and instrument status are checked.
	Maintenance	Checks the operation time/number of uses of GC and MS consumables and maintenance parts. If the result of the check is "Fail," replace the consumables and reset operation status.
	GC Check	Checks the status of the GC system. Fail: If repair or adjustment is required, contact your local Shimadzu representative. Replace any parts that need replacing.  NOTE GC system check criteria can be set by clicking [Rough Standard for Exchange] in each of the instrument property sub-windows of the [System Configuration] sub-window. These values are merely for reference and are not manufacturers' guaranteed values.
	MS Check	Checks setting values relating to MS data acquisition precision. Fail: If repair or adjustment is required, contact your local Shimadzu representative. Replace any parts that need replacing or perform auto-tuning if required.  NOTE MS system check criteria can be set by clicking [Test Value of System Check] in each of the instrument property sub-windows of the [System Configuration] sub-window. These values are merely for reference and are not manufacturers' guaranteed values.
②	Run	Executes the system check.
③	View Result	Displays the check results in the [System Check Result] sub-window.
④	Advanced	Sets GC and MS detailed check items.  NOTE This item can be set only when a user with [Modify system check settings] rights is logged in. Some items can be executed only when EI is selected as the ion source. Also, some items can be executed only when the MS has finished starting up. Items that cannot be executed are grayed out.
⑤	Run Auto Tuning As Needed / Run Auto-tuning	When [Run Auto Tuning As Needed] is selected, auto-tuning is executed when one of the [Signal Intensity Ratio], [Maximum Mass Shift] and [FWHM] MS check items is judged as "Fail." When [Run Auto-tuning] is selected, auto-tuning is executed regardless of the system check results.  NOTE <ul style="list-style-type: none"> Auto-tuning conditions can be set by clicking the  (Auto Tuning Condition) icon on the [Tuning] assistant bar in the [Tuning] window. Execution of auto-tuning using defaults is recommended. When auto-tuning is executed, previously made calibration curves sometimes can no longer be used.
⑥	Report Out	Automatically prints a result report after the system check is executed.

No.	Name	Explanation
7	(Progress bar)	Indicates the progress of the system check.

2.2.2 Executing a System Check


1 Set the system check items, and click [Run].



The system check is started. During execution, the time elapsed is incremented and the progress bar advances.

When the check ends without any problems, "Pass" is displayed above the progress bar. If a problem is found in the result, "Fail" is displayed.

2 To check the system check results, click [View Result].



The [System Check Result] sub-window is displayed.

Item	Judgment	Ratio	Actual/Test
<Maintenance GC>			
LCD Back Light	Pass	[+] 16%	7451/46380 Hour
Fan Motor	Pass	[+] 12%	7223/61320 Hour
Injection Unit Septum			
SPL1	Fail	[*****] 607%	607/100 times
SPL2	Pass	[] 0%	0/100 times
Injection Unit Insert			
SPL1	Fail	[*****] 121%	607/500 times
SPL2	Pass	[] 0%	0/500 times
CRG coolant Used Time			
<Maintenance MS>			
Filament #1		[****] 44%	436/1000 Hour
Filament #2		[] 7%	68/1000 Hour
Ion Source		[***] 31%	467/1500 Hour
Detector		[] 9%	511/6000 Hour
Rotary Pump 1		[+] 10%	1548/15000 Hour
Rotary Pump 1 Oil		[*****] 52%	1552/3000 Hour
<GC Check>			
Temp. Sensor in over 300C			
Column Oven	Pass	[] 0%	35/26280 Hour

2

Reference

When the system check result is [Fail], click [Help] and refer to GC check and MS check topics.

NOTE

- The system check results are saved to a file named in accordance with the following rule at \Systemcheck\Log in the folder where GCMSSolution is installed:
System No._YYYYMMDDHHMM.qgc
To check the results of system checks executed so far, click [Load] in the [System Check Result] sub-window, and select the file with the date of the result to be checked.
- The [Reset Consumables] sub-window is displayed by clicking the [Reset] button that is displayed when [Advanced] is clicked.
When a consumable is replaced, reset the operation time or number of uses.

2.3 Executing Auto-Tuning

There are two ways to tune the MS control parameters, auto-tuning (automatic adjustment) and tuning by [Peak Monitor] (manual adjustment).

Auto-tuning automatically adjusts the MS control parameters in accordance with the auto-tuning conditions so that desirable spectrum peaks can be acquired.

Peak monitor directly changes and tunes the individual control parameters of the instrument while observing the mass peaks.

This is used, for example, to check for any vacuum leaks.

This section describes auto-tuning.

Reference

When the ion source of the MS detector is EI, auto-tuning results can be judged automatically.


See "2.1.5 Setting Criteria for Auto-Tuning Result Judgment".

2.3.1 Executing Auto-Tuning


- 1 Click the  (Tuning) icon on the [Real Time] assistant bar in the [GCMS Real Time Analysis] program.

The [Tuning] window opens.

- 2 To create a new tuning file, click [New Tuning File] on the [File] menu.



When the ion source of the MS detector is EI, the [Select Tuning Mode] sub-window opens. Select the tuning mode according to the concentration of the component to be measured, and click [OK]. Measurement can be performed by the appropriate dynamic range since the tuning file is created by the emission current matched to the selected mode.




Set 20 uA, 60 uA or 150 uA as the emission current when [High Conc.], [Normal] or [High Sens.] is selected, respectively.

NOTE

When using an existing tuning file to perform tuning, click [Open Tuning File] on the [File] menu to display the [Open Tuning File] sub-window. Then, select the desired tuning file, click the [Open] button, and proceed to step 5.

- 3 Click the  (Auto Tuning Condition) icon on the [Tuning] assistant bar. The [Tuning Information] sub-window opens.



4 Set the items to auto-tune and their auto-tuning conditions, and click [OK].



NOTE


Execution of auto-tuning using default settings is recommended.

5


Click the (Start Auto Tuning) icon on the [Tuning] assistant bar.

The [Auto Tuning] sub-window is displayed, and auto-tuning is started. In the [Auto Tuning] sub-window, you can check the tuning status.


By default, auto-tuning will be performed with CID gas turned ON and then OFF when the [Use CID Gas] checkbox is selected.



When the [Print Result after Auto Tuning] checkbox on the [Tuning] menu is selected, a tuning report will be printed after auto-tuning ends.




- 6** Click [Save Tuning File As] on the [File] menu, and enter the tuning file name to save the file.



2.3.2 Checking Tuning Results

- 1** Click the  (Tuning Result View) icon on the [Tuning] assistant bar. The tuning results are displayed in the [Tuning] window.



No.	Name	Explanation
①	Ion Source	The ion source currently set in the system configuration is displayed. When the ion source is NCI, buttons for switching display of the tuning results for each of the NCI, SEI and SCI ionization modes are displayed.
②	Tuning Date	Displays the date and time that tuning was performed.
③	Acquisition Mode	Select the acquisition mode for displaying the tuning results.
④	CID Gas	Displays the tuning results for the CID gas status (ON/OFF) selected here.
⑤	Q1 Resolution Q3 Resolution	Select the resolution of Q1 or Q3 at SIM measurement. Selectable items change according to the acquisition mode.  Reference For details, refer to Help.


No.	Name	Explanation
⑥	[MS] tab page [GC] tab page	The [MS] tab page displays the MS monitor values that indicate instrument status during tuning and the MS control parameters saved to the tuning file. The [GS] tab page displays the GC instrument status (column information, temperature, pressure, flowrate) during tuning.
⑦	[Peak Profile] tab page [Spectrum] tab page	The [Peak Profile] tab page displays the peak profile and spectrum for a specific m/z of the standard sample (PFTBA) that was acquired when tuning was performed. The [Spectrum] tab page displays the spectrum of the standard sample (PFTBA) that was acquired when tuning was performed.
⑧	Judgment result display	Displays the result of auto-tuning judgment at the Q3 scan with CID gas set OFF when auto-tuning results are set to be judged in the System Configuration.

**NOTE**

When the [Tuning] window is opened, the previously saved tuning file or the tuning file used in the previous data acquisition is automatically loaded. Unless otherwise specified, the tuning file loaded here is used in system checks or in data acquisition.

2.3.3 Printing Tuning Reports

- 1 Click the  icon on the toolbar. The tuning results are output to the printer.



Reference

When tuning results are not set to be judged in the System Configuration, check the auto-tuning results using the procedure below.

- 1 Check that the FWHM (full width at half maximum) values are in the range 0.5 to 0.7.
- 2 Check that the detector voltage does not exceed 2 kV.
- 3 Check that the relative intensity ratio for m/z 502 is at least 2 %.
- 4 Check that the peak intensity for m/z 69 is at least twice that for m/z 28.

NOTE

If any irregularities are discovered above, possible causes could include a vacuum leak, poor column connections, or contaminated ion source. Remedy the problem by referring to the maintenance Help (MS Navigator) that is displayed by clicking [Maintenance] on the [Help] menu in the [GCMS Real Time Analysis] program.

3

Data Acquisition

This chapter describes how to use the instrument to set the data acquisition parameters and perform a single run.

3.1 Flow of Data Acquisition

1. Create a method.

Set the operating conditions for the GC, MS, autosampler, and other units. If a method file is already created, check its content.

Reference

See ["3.3 Setting Data Acquisition Methods"](#).

2. Register the sample.

Register the sample name, sample ID, data file name, and other information.

Reference

See ["3.4.1 Registering Samples"](#).

3. Inject sample and perform data acquisition.

Ready the instrument for data acquisition, inject the sample manually or automatically, and start data acquisition.

See ["3.4.2 Preparations"](#) and ["3.4.3 Starting Single Run"](#).

4. Complete acquisition

Acquired data is saved and data acquisition is complete.

Reference

For details on data processing from here on, see Chapters ["5 Postrun Analysis"](#) and ["8 Reports"](#)

Operations from data acquisition through to data processing can be also performed automatically.

See ["4 Continuous Data Acquisition"](#).

3.2 [Acquisition] Window


The [Acquisition] window has views for displaying sample information, chromatograms and spectra, and [Instrument Parameters View] for displaying and setting the parameters of each instrument.

3.2.1 Opening the [Acquisition] Window

1

Start up the [GCMS Real Time Analysis] program, and click the  (Data Acquisition) icon on the [Real Time] assistant bar.

The [Acquisition] window opens.



No.	Name	Explanation
①	Sample information	Displays the sample name, sample ID, and data comment. Items displayed at the sample information view are set on the [Setting 3] tab page in the [Display Settings] sub-window that is displayed by clicking [Display Settings] on the [View] menu.
②	Chromatogram View	Displays the chromatogram acquired on the instrument in real time. Oven temperature, carrier gas flow rate, and other information can also be displayed overlaying chromatograms. When performing acquisition on two lines, the chromatogram of each detector is displayed.
③	Spectrum View	Displays the MS spectrum acquired on the instrument in real time. Displays the spectrum obtained by measurement using the data acquisition conditions of event 1.
④	Instrument Method View	Sets the data acquisition conditions (method) of the GC, MS, autosampler, and other units. When 2 lines are set in the system configuration, set the instrument parameters for each line. Switch the line by using the [Select Line] button.
⑤	Instrument Monitor	Displays each of the parameters of the analytical instrument in real time. Currently, this is used for judging whether or not the instrument can start data acquisition immediately. When 2 lines are set in the system configuration, set the instrument parameters for each line.

**NOTE**

Set the display proportion of [Chromatogram View] and [Spectrum View] and the status curve of the instruments to display in [Chromatogram View] on the [Display Settings] sub-window that appears by selecting [Display Settings] on the [View] menu.

3.3 Setting Data Acquisition Methods

Set the parameters to be used for data acquisition in [Instrument Parameters View]. Parameters that are set are saved to a method file.

This section describes how to set the parameters in [Instrument Parameters View].


3.3.1 Creating New Method Files

1 Click [New Method File] on the [File] menu.




Hint

The [Setting Options] sub-window opens by clicking [Option] on the [Tools] menu. Selecting [Prompt on File New] in the [Method File] area on the [File New] tab page allows you to use a template file when creating a method file.




GCMSsolution includes method template files suitable for frequently used acquisition modes, as shown below. Setting items on MS Tables that are not used in some acquisition modes are hidden in method template files.



 **NOTE**

- If [Instrument Parameters View] is not displayed, click the  (Method Detail) icon on the [Acquisition] assistant bar.
- To open an existing method file, drag-and-drop the method file onto the [Acquisition] window from the [Method] tab page in the [Data Explorer] sub-window.




3.3.2 AOC Instrument Parameters

Set the injection conditions when using the autosampler to inject samples.

This sub-window is displayed when an autosampler is registered at [Modules Used for Analysis] in the system configuration.

1

Click the [Sampler] tab, and set the autosampler parameters.


 **NOTE**

Set the parameters in more detail in the [Advanced] sub-window that is displayed by clicking [Advanced]. Normally, data acquisition can be performed at the default settings.


 **Reference**

For details about each of the parameters, refer to Help.

3.3.3 GC Instrument Parameters


Set the gas chromatograph (GC) operating conditions.

1 Click the [GC] tab.



2 Set the parameters.

■ [Injection Port] window



- 1 Click the injection unit from the unit tree used.
- 2 Set the temperature of the injection unit taking the boiling point of the target component into consideration.
- 3 Select the injection mode.



NOTE

- Select split mode when there is a large volume of target component to inject. (Amount of target component: 10 ng to 100 ng min.)
- Select splitless mode when there is a small volume of target component to inject. (Amount of target component: less than 10 ng)

- 4 Set the carrier gas.



Reference

For details about each of the parameters, refer to Help.

- 5 Create flow rate and pressure programs so that the target component separates from other peaks.

■ [Column] window


1 Click the column unit from the unit tree used.

2 Set the initial temperature of the column oven.

3 Create a column oven temperature program so that the target component separates from other peaks.

■ [General] window

When varying the parameters for controlling the GC by time during measurement, set a prerun program or a time program.



- 1 Click the general unit from the unit tree used.
- 2 At [Time], enter the retention time during which processing is executed, and at [Device], select the unit on which processing is executed.
At [Event], select the processing event.
Set values for processing at [Set Value], if required.



NOTE

- A "prerun program" is a time program for controlling each unit before data acquisition.
- A "time program" is a time program for controlling each unit during data acquisition.

3.3.4 MS Instrument Parameters


Set the mass spectrometer (MS) operating conditions.

There are also other ways of setting MS tables. Refer to the following sections.


Reference

- ["3.3.4.1 Creating MS Tables Based on Compound Tables"](#)
- ["3.3.4.2 Creating FASST Measurement Methods"](#)
- GCMS Operation Guide - Basic Operation Guide
- GCMS Operation Guide - Method Development Guide

1 Click the [MS] tab.



2 Set the MS parameters.




No.	Name	Explanation
①	Instrument parameters	<p>Sets the ion source temperature, interface temperature, detector voltage, and other parameters.</p> <p> Reference</p> <p>For details about each of the parameters, refer to Help.</p>
②	MS Table	<p>Set detailed data acquisition settings for each event. The group No. and event No. of events are indicated joined by a hyphen "-" on the left edge of the table. If the start time and end time of adjacent events are each set to the same value, the events will automatically belong to the same group.</p> <p>Select the acquisition mode from the following eight modes: MRM, Product Ion Scan, Precursor Ion Scan, Neutral Loss Scan, Q3 Scan, Q3 SIM, Q1 Scan or Q1 SIM</p> <p> Reference</p> <p>Available parameters differ according to the acquisition mode. For details, refer to Help.</p> <p> NOTE</p> <p>Up to 128 groups can be registered to one method file. Up to 64 events can be registered to one group. The maximum total number of events for all groups is 2048 in the case of SIM type events and 128 in the case of scan type events.</p>

NOTE

With this instrument, data acquisition settings can be successively switched during a single data acquisition. Each data acquisition setting is called an "event". When multiple events have been set to the same measurement time range (group), operation switches to the data acquisition settings of the next event after the event time set to each event has elapsed. When measurement of the last event ends, operation returns to the data acquisition settings of the first event.

In the example below, for ten minutes starting from one minute, data acquisition is performed in the order event 1 -> event 2 -> event 3 -> event 1 -> event 2 -> event 3 and so forth. The cycle from event 1 to event 3 is called Loop. For 20 minutes starting from ten minutes, data acquisition is performed in the order event 1 -> event 2 -> event 1 -> event 2 and so forth.




MS Table Setting Example (MRM Acquisition Mode)

Set the parameters for measuring by SIM for both Q1 and Q3.

Compound Name	Start Time (min)	End Time (min)	Acq. Mode	Event Time (min)	Ch1 m/z	Ch2 m/z	Ch3 m/z	Q1 Resolution	Q3 Resolution
1-1 A	3.00	10.00	MRM	0.003	273.10>100	5.00		Low	Low
1-2 B	3.00	10.00	MRM	0.300	455.20>164	5.00		Low	Low
1-3 C	3.00	10.00	MRM	0.300	308.10>162	10.00		Low	Low

- 1 Hide unwanted columns in the MS Table, and set the number of ions (i.e. number of channels) to be measured.

The [Table Style] sub-window opens by right-clicking on the MS Table, and clicking [Table Style] on the displayed menu.



In this sub-window, hide unwanted columns and set the number of channels to display in the MS Table at [Number of Ch].

- 2 At the row of the event to be set, set the compound name of the component to be measured. For events in the MRM acquisition mode, set the data acquisition settings of one component per one event.
- 3 Set the time that data acquisition will start after the sample is injected, and set the time that data acquisition will end.
- 4 Click the [Acq. Mode] cell, and select [MRM] from the list that is displayed.
- 5 At the [Event Time] cell, set the time that spectral data will be acquired at each event.
- 6 Set the m/z of the precursor ion and the m/z of the product ion that is acquired by destruction of the precursor ion on each channel.

**NOTE**

Enter in the format "precursor ion m/z > product ion m/z ". Up to 16 ions can be set.

- 7 Enter the **collision energy of the channel**. 1.1.5.1

3 Repeat steps 1 through 6 to set the other components.

Acquire Data without Using CID Gas(Q3Scan)

Even when the option of using CID gas is selected for the system configuration, selecting the [Acquire Data without Using CID Gas(Q3Scan)] checkbox enables data acquisition in Q3 Scan and Q3 SIM modes with CID gas turned OFF. Without effects of CID gas, better spectrum peaks can be acquired.

	Compound Name	Start Time (min)	End Time (min)	Acq. Mode	Event Time(sec)	Scan Speed	Start m/z	End m/z	Precursor m/z	Product m/z
1-1		3.00	10.00	Q3 Scan	0.300	1666	45.00	500.00		
		0.00	0.00	MRM	0.000					

**NOTE**

This function is valid only in Q3 Scan or Q3 SIM acquisition mode.

However, when an MS/MS acquisition mode (MRM, Product Ion Scan, Precursor Ion Scan, or Neutral Loss Scan) is included in the MS table, CID gas will not be turned OFF even with the [Acquire Data without Using CID Gas(Q3Scan)] checkbox selected.


**NOTE**

This parameter can be selected only when the option of using CID gas is selected for the system configuration.


MS Program

When varying the parameters for controlling the MS by time during measurement, set the MS program.

1 Select the [Use MS Program] checkbox, and click [Set].



2 Make an MS program, and click [OK].



1 At [Time (min)], enter the retention time during which processing is executed, and at [Command], select processing content.

Set values for processing at [Value], if required.

2 Click [OK].



NOTE

- The [GCMS Data File Open] sub-window opens by clicking the [Load Chromatogram] button. When the data file is selected and [Open] is clicked, the TIC of that data file is displayed.
- When a chromatogram is clicked with the [Time (min)] cell selected, the retention time for the clicked part of the chromatogram is set to the cell.

3.3.4.1 Creating MS Tables Based on Compound Tables

This section describes how to copy data acquisition settings, which are set to spreadsheet format files (e.g. Excel) that are preset with data acquisition settings, to Compound Tables, and use COAST to create MRM/SIM events for MS parameters based on the Compound Table.

COAST is used to create MS Tables for MRM or SIM measurement.

The MS Table for measuring one compound per event can be created.

1 Copy the data acquisition settings of the spreadsheet format file to the Compound Table.


- 1 Select and copy the content in the spreadsheet format file.
Make sure that the columns in the compound table and in the spreadsheet format file are arranged in the same order.

ID#	Name	Type	ISTD Group	m/z	Ret. Time	Ret. Index	Unit	Ref. Ions
4	Octafluoronaphthalene	Target	1	272>241	4.78			272>222-272>172
5	Dichlorvos	Target	1	185>109	5.91			185>93-185>63
6	Fenobucarb	Target	1	150>121	9.665			150>103-150>77
7	Hexachlorobenzene	Target	1	284>249	10.805			284>214-284>212
8	Simazine	Target	1	201>173	11.095			201>186-201>158
9	Propyzamide	Target	1	173>145	11.555			173>109-173>74
10	Diazinon	Target	1	304>179	11.61			304>162-304>195
11	Chlorothalonil	Target	1	266>231	11.71			266>170-266>168
12	Iprobenfos	Target	1	204>91	12.12			204>122-204>171
13	Fenitrothion	Target	1	277>260	13.165			277>109-277>125
14	Thiobencarb	Target	1	257>100	13.5			257>72-257>224
15	Stearic acid, methyl ester	Target	1	298>101	14.885			298>185-298>143
16	Isoprothiolane	Target	1	290>204	15.265			290>118-290>162
17	Isoxathion	Target	1	313>177	15.77			313>130-313>159
18	Chlornitrofen	Target	1	317>287	16.68			317>236-317>224
19	EPN	Target	1	157>110	17.625			157>77-157>133


- 2 Click [Quantitative Parameters] on the [Method] menu in the [Acquisition] window, and click the [Compound Table] tab in the [Quantitative Parameters] sub-window that is displayed.

ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit
1		Target	1	TIC	0.000	0	

- 3 Click at the start of the row in the Compound Table to copy, and select [Paste] on the right-click menu.



Rows are automatically added, and the content is pasted.



- 4 Click [OK] to close the [Quantitative Parameters] sub-window.
The screen returns to the [Acquisition] window.

2 Use COAST to create the MS Table based on the Compound Table.

- 1 Click [Creation of Automatic MS Table [COAST]] on the [Method] menu.
The [Creation of Automatic MS Table [COAST]] sub-window is displayed.
The retention time set to the Compound Table and the MRM/SIM Table with components divided into measurement groups based on the processing time of peak integration parameters are displayed at the bottom of the sub-window.
- 2 Adjust the end time so that there are 20 channels or less for each group.

Creation of Automatic MRM or SIM Table [COAST]

Here, you can create MRM, SIM, or Product Ion Scan table automatically by using a compound table. Please select data acquisition mode, and press "Update Table" button.

MRM/SIM Mode : All Compounds of the table are measured by MRM or SIM mode.
 FASST Mode : Only checked compounds in the table are measured by Scan and MRM or Scan and SIM mode synchronously. Not checked compounds are measured by Scan mode.
 Product Ion Scan Mode : All Compounds of the table are measured by Product Ion Scan mode.

Select Data Acquisition Mode: MRM/SIM FASST Product Ion Scan Requested Loop Time: 0.300 sec

ID#	Name	Ret. Time	m/z	Ref. Ions	FASST
1	Dichlorvos	5.910	185.00>109.00	2	<input type="checkbox"/>
2	Fenobucarb	9.665	150.00>121.00	2	<input type="checkbox"/>
3	HCB	11.095	284.00>249.00	2	<input type="checkbox"/>
4	Simazine	11.585	201.00>173.00	2	<input type="checkbox"/>
5	Propyzamide	11.610	173.00>145.00	2	<input type="checkbox"/>
6	Diazinon	11.710	304.00>179.00	2	<input type="checkbox"/>
7	Chlorothalonil	12.120	266.00>231.00	2	<input type="checkbox"/>
8	Iprobenfos	13.165	204.00>91.00	2	<input type="checkbox"/>
9	Fenitrothion	13.500	277.00>260.00	2	<input type="checkbox"/>
10	Thiobencarb	15.265	257.00>100.00	2	<input type="checkbox"/>
11	Stearic acid	15.770	298.00>101.00	2	<input type="checkbox"/>
12	Isoprothiolan	16.680	290.00>204.00	2	<input type="checkbox"/>
13	Isoxathion	17.825	313.00>177.00	2	<input type="checkbox"/>

End Time (min)	m/z-1	m/z-2	m/z-3	m/z-4	m/z-5	m/z-6	m/z-7	m/z-8	m/z-9	m/z-10
1	7.79	317>267	157>110	185>109	185>93	185>63	0	0	0	0
2	10.38	150>121	150>103	150>77	0	0	0	0	0	0
3	12.64	284>249	201>186	201>158	201>173	173>109	173>74	173>145	304>162	304>195
4	14.38	204>91	277>109	277>125	277>260	257>72	257>224	0	0	0
5	17.25	257>100	290>118	290>162	298>101	313>130	313>159	290>204	317>236	317>224
6	18.32	313>177	157>77	157>139	0	0	0	0	0	0

Start Time: 0.00 min Load Chromatogram Update Table

OK Cancel Help



NOTE

When there is a group division where peaks are closely concentrated, set the end time so that there is a gap between peaks.

3 Click [OK].

An MS table is created with one compound assigned to one event.

Compound Name	Start Time (min)	End Time (min)	Acq. Mode	Event Time (min)	Ch1 m/z	Ch1 CE	Ch2 m/z	Ch2 CE	Ch3 m/z	Ch3 CE	Q1 Description
1-1 Dichlorvos	3.00	7.79	MRM	0.300	185.00>109	15.00	185.00>93.0	15.00	185.00>63.0	15.00	Low
2-1 Fenobucarb	7.79	10.38	MRM	0.300	150.00>121	15.00	150.00>103	15.00	150.00>77.0	15.00	Low
3-1 Simazine	10.38	12.64	MRM	0.300	201.00>173	15.00	201.00>186	15.00	201.00>158	15.00	Low
3-2 Propyzamide	10.38	12.64	MRM	0.300	173.00>145	15.00	173.00>109	15.00	173.00>74.0	15.00	Low
3-3 Diazinon	10.38	12.64	MRM	0.300	304.00>179	15.00	304.00>162	15.00	304.00>195	15.00	Low
3-4 Chlorothalonil	10.38	12.64	MRM	0.300	266.00>231	15.00	266.00>170	15.00	266.00>168	15.00	Low
3-5 Iprobenfos	10.38	12.64	MRM	0.300	204.00>91.0	15.00	204.00>122	15.00	204.00>171	15.00	Low
4-1 Fenrothion	12.64	14.38	MRM	0.300	277.00>260	15.00	277.00>109	15.00	277.00>125	15.00	Low
4-2 Thioencarb	12.64	14.38	MRM	0.300	257.00>100	15.00	257.00>72.0	15.00	257.00>224	15.00	Low
5-1 Isothiolane	14.38	17.25	MRM	0.300	290.00>204	15.00	290.00>118	15.00	290.00>162	15.00	Low
5-2 Isoxathion	14.38	17.25	MRM	0.300	313.00>177	15.00	313.00>130	15.00	313.00>159	15.00	Low
5-3 CNP	14.38	17.25	MRM	0.300	317.00>287	15.00	317.00>236	15.00	317.00>224	15.00	Low
6-1 EPN	17.25	19.32	MRM	0.300	157.00>110	15.00	157.00>77.0	15.00	157.00>139	15.00	Low

At this time, the event numbers of the Compound Table are also set.

ID#	Ret. Index	Unit	Ref. Ions	Conc. 1	Event	STD Spec
1	0		272.00>222.00-2	1	1	
2	0		185.00>93.00-18	1	1	
3	0		150.00>103.00-1	1	1	
4	0		284.00>214.00-2	1	1	
5	0		201.00>186.00-2	1	2	
6	0		173.00>109.00-1	1	3	
7	0		304.00>162.00-3	1	4	
8	0		266.00>170.00-2	1	5	
9	0		204.00>122.00-2	1	6	
10	0		277.00>109.00-2	1	1	
11	0		257.00>72.00-25	1	2	
12	0		298.00>185.00-2	1	1	
13	0		290.00>118.00-2	1	2	
14	0		313.00>130.00-3	1	3	
15	0		317.00>287.00-2	1	4	

3.3.4.2 Creating FASST Measurement Methods

FASST (Fast Automated SCAN/SIM Type) measurement refers to switching between the scan acquisition mode and MRM (or SIM) acquisition mode in a single measurement group.

The purpose of FASST measurement is broadly classified into two types.

Scan measurement is suited to qualitative analysis of components since mass spectra over a broad range can be obtained. This is also used for quantitative analysis since components can be identified more reliably than by MRM/SIM measurement. Yet, the difference in sensitivity required for each component increases as the number of target components to be batch-analyzed increases, with the result that sensitivity is sometimes insufficient by scan measurement. In cases such as this, high-sensitivity data acquisition is possible by performing FASST measurement on components that require higher sensitivity. Also, by scan measurement, elution of unknown components can be checked from the mass spectra of masses that have not been measured by MRM/SIM measurement.

This section describes how to create a method file for FASST measurement from a Compound Table when a method file for MRM/SIM measurement is already available.

1 2

In the [Acquisition] window, open the method file for MRM/SIM measurement.

Click [Creation of Automatic MS Table [COAST]] on the [Method] menu.

The [Creation of Automatic MS Table [COAST]] sub-window is displayed.

The retention time set to the Compound Table and the MRM/SIM Table with components divided into measurement groups based on the processing time of peak integration parameters are displayed at the bottom of the sub-window.

3 Select [FASST] at [Select Data Acquisition Mode], select the [FASST] column checkboxes of the compounds to undergo FASST measurement. Click [Update Table].

Here, you can create MRM, SIM, or Product Ion Scan table automatically by using a compound table. Please select data acquisition mode, and press "Update Table" button.

MRM/SIM Mode : All Compounds of the table are measured by MRM or SIM mode.
 FASST Mode : Only checked compounds in the table are measured by Scan and MRM or Scan and SIM mode synchronously. Not checked compounds are measured by Scan mode.
 Product Ion Scan Mode : All Compounds of the table are measured by Product Ion Scan mode.

Select Data Acquisition Mode: MRM/SIM FASST Product Ion Scan Requested Loop Time: 0.300 sec

ID#	Name	Ret. Time	m/z	Ref. Ions	FASST
1	Dichlorvos	5.910	185.00>109.00	2	<input checked="" type="checkbox"/>
2	Fenobucarb	9.665	150.00>121.00	2	<input checked="" type="checkbox"/>
3	HCB	11.095	284.00>249.00	2	<input checked="" type="checkbox"/>
4	Simazine	11.555	201.00>173.00	2	<input checked="" type="checkbox"/>
5	Propyzamide	11.610	173.00>145.00	2	<input checked="" type="checkbox"/>
6	Diazinon	11.710	304.00>179.00	2	<input checked="" type="checkbox"/>
7	Chlorothalonil	12.120	266.00>231.00	2	<input checked="" type="checkbox"/>
8	Iprobenfos	13.165	204.00>91.00	2	<input checked="" type="checkbox"/>
9	Fenitrothion	13.500	277.00>260.00	2	<input checked="" type="checkbox"/>
10	Thiobencarb	15.265	257.00>100.00	2	<input checked="" type="checkbox"/>
11	Stearic acid	15.770	298.00>101.00	2	<input checked="" type="checkbox"/>
12	Isoprothiolan	16.680	290.00>204.00	2	<input checked="" type="checkbox"/>
13	Isoxathion	17.825	313.00>177.00	2	<input checked="" type="checkbox"/>

Start Time: 0.00 min

MRM/SIM Table

End Time (min)	m/z-1	m/z-2	m/z-3	m/z-4	m/z-5	m/z-6	m/z-7	m/z-8	m/z-9	m/z-10
1	7.79	317>287	157>110	185>109	185>93	185>63	0	0	0	0
2	10.38	150>121	150>103	150>77	0	0	0	0	0	0
3	12.64	284>249	201>186	201>158	201>173	173>109	173>74	173>145	304>162	304>195
4	14.38	204>91	277>109	277>125	277>260	257>72	257>224	0	0	0
5	17.25	257>100	290>118	290>162	290>101	313>130	313>159	290>204	317>236	317>224
6	18.32	313>177	157>77	157>139	0	0	0	0	0	0

You can change the end time and m/z of the current row by clicking on the graph.

4 If necessary, change the end time setting of the group.

5 Click [OK] to complete COAST.

An MS Table set with FASST measurement events is created.

GCMS-TQ Series

Ion Source Temp.: 230 °C
 Interface Temp.: 25 °C
 Solvent Cut Time: 2 min
 Detector Voltage: Relative to the Tuning Result Absolute
 Do Not Use CID Gas in Q3 Scan or Q3 SIM Mode
 Use MS Program: Set...
 GC Program Time: 19.67 min


Compound Name	Start Time (min)	End Time (min)	Acq. Mode	Event Time(sec)	Ch1 m/z	Ch1 CE	Ch2 m/z	Ch2 CE	Ch3 m/z	Ch3 CE	Q1 Resolution
1-1	3.00	7.79	Q3 Scan	0.300							
1-2	Dichlorvos	3.00	MRM	0.300	185.00>109	15.00	185.00>93.0	15.00	185.00>63.0	15.00	Low
1-3	3.00	7.79	Q3 Scan	0.300							
2-1	Fenobucarb	7.79	MRM	0.300	150.00>121	15.00	150.00>103	15.00	150.00>77.0	15.00	Low
2-2	7.79	10.38	Q3 Scan	0.300							
3-1	Simazine	10.38	MRM	0.300	201.00>173	15.00	201.00>186	15.00	201.00>158	15.00	Low
3-2	Propyzamide	10.38	MRM	0.300	173.00>145	15.00	173.00>109	15.00	173.00>74.0	15.00	Low
3-3	Diazinon	10.38	MRM	0.300	304.00>179	15.00	304.00>162	15.00	304.00>195	15.00	Low

NOTE

- Quantitation is performed to the compounds whose [FASST] column checkboxes are not selected.
- The measurement mass range and event time for scan measurement events are automatically set. Confirm the [Start m/z], [End m/z], and [Scan Speed] in the MS Table, and if necessary, change the settings for the measurement mass range and event time. If [Start m/z], [End m/z], and [Scan Speed] are not displayed in the MS Table, select [Table Style] on the right-click menu on the MS Table, and add them to the display items.


3.3.5 Saving Method Files

- 1 Click [Save Method File As] on the [File] menu.



The [Save Method File As] sub-window opens.

- 2 Select the folder to save the method file in, enter a file name, and click [Save]. Specified parameters are saved as a method file.



3

3.4 Starting Single Run

There are two ways of acquiring data, by single run (data acquisition executed only once) or by continuous data acquisition (data acquisition performed continuously from multiple samples).


This section describes how to perform single run acquisition.

3.4.1 Registering Samples

Before executing single run, set the name or ID of the sample to be used for data acquisition, data file name to save acquired data to, vial No., report settings, format file name, and other settings.

Click [Advanced] to configure the report settings and format file name.

- 1** Click the  (Sample Login) icon on the [Acquisition] assistant bar.



- 1** Enter the data file name.

 **NOTE**

- When only the file name is entered, the sample is saved to the same folder as the method file open in the [Acquisition] window.
- To prepare for data acquisition when a data file of the same name already exists, a confirmation sub-window is displayed to prompt whether to overwrite the file or automatically rename the file.

- 2** Enter the position of the sample vial.

 **NOTE**

When "0" or "-1" is set, data acquisition is performed without samples being injected from the autosampler.

- 3** Enter the sample injection volume.


3.4.2 Preparations

Download the [Instrument Parameters View] parameters to the instrument to ready it for data acquisition.

- 1 Click the  (Download) icon on the [Acquisition] assistant bar.

 **NOTE**

- When the  (Download) icon is clicked without the parameters set at [Instrument Parameters View] saved to a method file, a message box for selecting whether or not to save the currently set method to the method file will be displayed.



To save the method file, click [Yes]. If the method file is opened and edited, it will be overwritten. Otherwise, the [Save Method File As] sub-window opens. Select the folder, enter the name of the method file, and save the file.

To ready the system without saving the method file, click [No]. At this time, too, the set parameters are downloaded to the instrument.

To cancel readying of the instrument, click [Cancel].

- When a sampler is connected to the instrument, single run is automatically started if the GC status after it is readied for data acquisition is [Ready].

3.4.3 Starting Single Run

Start single run data acquisition.

- 1 When a sampler (AOC or HS) is not connected to the instrument, click the  (Start) icon on the [Acquisition] assistant bar.

During data acquisition, the chromatogram and spectrum are displayed, and the status of [Instrument Monitor] is [Run].

When data acquisition ends, the status of [Instrument Monitor] returns to [Ready].

 **NOTE**

- Data acquisition ends when the GC program time or the MS measurement time set to the method, whichever is longer, elapses.
- Data acquisition can be also started by pressing the [START] button on the GC unit. The same applies when a sampler is connected.

 **NOTE**

While the [GCMS Real Time Analysis] program is in operation, do not insert or remove devices (e.g. flash memory, printers) connected on the USB interface or turn the power ON/OFF. Doing so might impede communications with the instrument.

3.4.4 Stopping Single Run

Stop single run data acquisition.

- 1 Click the  (Stop) icon on the [Acquisition] assistant bar.




When the  (Stop) icon is clicked to stop measurement midway, the MS program stops at that point, but the GC program does not stop. This is to prevent sample from remaining in the column. If it is OK to stop GC operation, press the [STOP] button on the GC unit.

3.4.5 Extending the Data Acquisition Time


The time may be extended during data acquisition.

This section describes how to extend the data acquisition time.

- 1 Click [Extend] on the [Acquisition] menu.



- 2 Enter the extended time for both GC and MS in minute increments.




NOTE

The data acquisition time can be extended only when one event is set in one group to MS parameters and when one of the Q3 Scan, Q1 Scan, Product Ion Scan, Precursor Ion Scan, or Neutral Loss Scan acquisition modes is selected.

3.4.6 Snapshot

Data obtained so far since start of measurement can be displayed and analyzed during data acquisition.

- 1 Click the  (Snapshot) icon on the [Acquisition] assistant bar.
Data obtained so far since start of data acquisition are displayed in the [Data Analysis] window.



- 2 Analyze the data in the [Data Analysis] window.
- 3 To load acquired data from a continuing data acquisition after snapshot has been executed, click  (Update) on the toolbar in the [Data Analysis] window.
The latest data are displayed in the [Data Analysis] window.


3

3.5 Monitoring Instrument Status

This section describes how to monitor the current status of the instrument.

3.5.1 Instrument Monitor

[Instrument Monitor] displays the status of instruments and the measured values of various parameters. In the [Monitor Settings] sub-window, parameter settings can be changed without changing the instrument methods in the method file.



No.	Name	Explanation
①	Instrument status	The color changes according to the current status of the GC, MS and HS. Reference For details, refer to Help.
②	Line1/Line2	Switches the line to be monitored. This item is displayed only when 2 lines are set in the system configuration.
③	Units	Monitors the information of units on the selected line.
④	Vacuum	Displays the vacuum level as a graph. The triangular mark moves according to the measured values. [H.Vac.] is not displayed when the ion gauge is set to [Ion Gauge] in the system configuration.
⑤	GC Consumables MS Consumables	The color of these items changes according to the degree of consumption. Reference For details, refer to Help.
⑥	Detail	Clicking this item opens the [Monitor Settings] sub-window. Here, units' measured values and settings for other parameters, and the current status and estimated replacement time of consumables can be checked. Reference For details, refer to Help.

No.	Name	Explanation
7	Ecology Mode	Units can be turned OFF and carrier gas flow rate can be reduced when the instrument is standing by for data acquisition. The original status is automatically returned to after the ecology mode is canceled.

3.6 Ecology Mode

1.6.1

The ecology mode can be used to reduce power consumption or consumed carrier gas and CID gas when the instrument is standing by for data acquisition. When the instrument enters the ecology mode, the shutdown method is downloaded to the instrument. When the ecology mode is canceled, the instrument is returned to the mode before it entered the ecology mode. There are two ways of entering the ecology mode, manual and automatic. This section describes manual operations.


3

Reference


For details on how to automatically enter the ecology mode, see ["4.6 Entering the Ecology Mode after Continuous Data Acquisition"](#).

3.6.1 Setting the Ecology Mode

- 1 Click [Ecology Settings] on the [Tools] menu in the [Acquisition] window.



- 2 Set the content of the ecology mode, and click [OK].




- 1 To download the shutdown method file to the instrument, select the [Shutdown Method File] checkbox, and specify the name of the method file.

NOTE


To edit the parameters of the shutdown method file, click the [Edit] button.

This starts up the [GCMS Analysis Editor] program, and the method file is loaded to the [Acquisition] window. After parameters have been edited, save the method file.


- 2 Click the button of the line for the parameters in the shutdown method file to download to the instrument.
- 3 Select the parameters to download to the instrument, and click [OK].



- 4 To control the unit when the instrument enters the ecology mode, click the button for the line corresponding to the unit.



- 5 Select the checkboxes of the commands to download to the instrument, and click [OK].



- 6 In the [Ecology Settings] sub-window, and click [OK].

3

3.6.2 Entering the Ecology Mode 1.6.1

- 1 Click the  (Ecology Mode) button at [Instrument Monitor] in the [GCMS Real Time Analysis] program.




NOTE

This button cannot be clicked while data acquisition, tuning and system check are in progress or while the instrument is starting up or is shutting down.

- 2 A message confirming whether or not to enter the ecology mode is displayed. Click [Yes].

The ecology settings are downloaded to the instrument, and the instrument enters the ecology mode.



NOTE

While the instrument is in the ecology mode, the [Ecology Mode] sub-window is displayed, and other sub-window operations in the [GCMS Real Time Analysis] program are not possible until the ecology mode is canceled.

3.6.3 Canceling the Ecology Mode

1

In the [Ecology Mode] sub-window, click [Cancel].

When the ecology mode is canceled, the instrument returns to the mode before it entered the ecology mode.



3.6.4 Setting the PC or Display to the Power Save Mode

To set the ecology mode when the PC or display is not in use, set the mode in the Windows power options.

 **NOTE**

- The instrument will not automatically switch to the power save mode while the instrument is running. The instrument will switch to the power save mode after the current process is finished and a specified amount of time elapses.
- The power save mode can be activated manually, depending on power option settings. However, do not manually switch to the power save mode while the instrument is running. Doing so will cause an error communicating with the instrument.
- This feature is not compatible with pretreatment units other than HS-20 or AOC-20.

4

Continuous Data Acquisition

Continuous data acquisition refers to data acquisition performed continuously on multiple samples. To perform continuous data acquisition, first start by making a Batch Table.

4.1 Opening the [Batch Table] Window

- 1 Start up the [GCMS Real Time Analysis] program, and click the  (Batch Processing) icon on the [Acquisition] assistant bar.
The [Batch Table] window opens.

4

4.2 Creating Batch Tables


To acquire data continuously from multiple samples, set the sample information, vial #, method file, and data file name to a Batch Table.

This section describes how to create a Batch Table.

4.2.1 Batch Table Wizard

Batch Tables can be made easily by using the Batch Table Wizard.

- 1 Click the  (Wizard) icon on the [Batch] assistant bar.
The [Batch Table Wizard] sub-window opens. Set the following parameters.




- 1 Click [New] at [Batch Table].
To add a row to the currently displayed batch file, click [Append].


- 2 Set the data acquisition pattern for the standard sample and unknown sample.
This description is for an example of how to create a Batch Table used to perform data acquisition on both a standard sample and unknown sample.
- 3 Set [Method File] to be used for continuous data acquisition.
- 4 Set whether to perform quantitative processing or qualitative processing on the data.
- 5 When settings are finished, click [Next].
The [Batch Table Wizard Standard Sample (1)] sub-window opens.

2 Set the standard sample information.

This sub-window is displayed only when the sample type is set to [Standard & Unknown] or [Standard Only] in the previous sub-window.




- 1 Enter [Vial #] where injection of the standard sample will begin.
- 2 Set [# of Calib. Levels] and [Average Count] that injection is performed from the same vial.
- 3 Set the injection volume.
- 4 Enter [Sample Name] and [Sample ID] of the standard sample.
When [Auto-increment] is selected, [Sample Name] and [Sample ID] will automatically be assigned to a continuous number.
- 5 Click [Next].
The [Batch Table Wizard Standard Sample (2)] sub-window opens.



- 6 Enter [Data File Name].
To automatically generate a data file name, select [Auto-increment].

**NOTE**

Set the format of the data file name on the [Data Filename] tab page in the [Settings] sub-window that is displayed by clicking the  (Settings) icon on the [Batch] assistant bar.



Batch Table settings matching [Selected Items] are joined by a hyphen "-", and are appended a continuous number in accordance with the number assignment rule selected at [Auto-increment Format].


- 7 To print a report, select [Report Out], and set the report format file name.

**Reference**

To create a new report format, refer to ["8.2 Creating Report Formats"](#).

- 8 Click [Next].
The [Batch Table Wizard Unknown Sample (1)] window opens.

3 Set the unknown sample information in the same way as the standard sample.




1 Enter [Sample Name] and [Sample ID] of the unknown sample.

When [Auto-increment] is selected, [Sample Name] and [Sample ID] will automatically be assigned to a continuous number.

2 Click [Next].

The [Batch Table Wizard Unknown Sample (2)] window opens.



3 Enter [Data File Name].

To automatically generate a data file name, select [Auto-increment].

4 To print a report, select [Report Out], and set the report format file name.

Reference

To create a new report format, refer to ["8.2 Creating Report Formats"](#).

5 When settings are finished, click [Finish].

The [Batch Table Wizard] sub-window closes, and the Batch Table is created.

NOTE

If there are values other than those in the Wizard that require setting, directly edit the Batch Table.

4.2.2 Editing Batch Tables


In daily data acquisition, partially editing and executing an existing batch file is handy. This section describes an example of how to batch-edit specified rows based on the vial #, sample name, data file name, and other information of pre-edited rows. This section describes two functions, [Fill Series] and [Fill Down], that are often used for directly editing Batch Tables.

Appending Settings with a Continuous Number

The [Vial#], [Sample Name], [Sample ID], and [Data File] items in the Batch Table are assigned a continuous number.

This section describes how to assign a continuous number to [Vial#].

- 1 Drag-and-drop the batch file from the [Batch] tab page in the [Data Explorer] sub-window onto the [Batch Table] window.



- 2 Edit the vial No. of the first row to set with a continuous number.

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	Report
1	2	OFN-100pg-01	0:Unknown	IT QT	TEST_Ultra.qgm	Repeat_000.qgd	1	<input type="checkbox"/> Print	
2	1	OFN-100pg-02	2-Control	IT QT ILM	TEST_Ultra.qgm	Repeat_001.qgd	1	<input type="checkbox"/> Print	
3	1	OFN-100pg-03	2-Control	IT QT ILM	TEST_Ultra.qgm	Repeat_002.qgd	1	<input type="checkbox"/> Print	
4	1	OFN-100pg-04	2-Control	IT QT ILM	TEST_Ultra.qgm	Repeat_003.qgd	1	<input type="checkbox"/> Print	
5	1	OFN-100pg-05	2-Control	IT QT ILM	TEST_Ultra.qgm	Repeat_004.qgd	1	<input type="checkbox"/> Print	
6	1	OFN-100pg-06	2-Control	IT QT ILM	TEST_Ultra.qgm	Repeat_005.qgd	1	<input type="checkbox"/> Print	

- 3 Right-click the edited cell, and click [Fill Series].

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	Report
1	2			IT QT	TEST_Ultra.qgm	Repeat_000.qgd	1	<input type="checkbox"/> Print	
2				IT QT ILM	TEST_Ultra.qgm	Repeat_001.qgd	1	<input type="checkbox"/> Print	
3				IT QT ILM	TEST_Ultra.qgm	Repeat_002.qgd	1	<input type="checkbox"/> Print	
4				IT QT ILM	TEST_Ultra.qgm	Repeat_003.qgd	1	<input type="checkbox"/> Print	
5				IT QT ILM	TEST_Ultra.qgm	Repeat_004.qgd	1	<input type="checkbox"/> Print	
6				IT QT ILM	TEST_Ultra.qgm	Repeat_005.qgd	1	<input type="checkbox"/> Print	

The vial numbers, and the numbers appended to sample names and data file names will be set with a continuous number.

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	Report Fi
1	2	OFN-100pg-01	0:Unknown	IT QT	TEST_Ultra.qgm	Repeat_000.qgd	1	<input type="checkbox"/> Print	
2	3	OFN-100pg-02	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_001.qgd	1	<input type="checkbox"/> Print	
3	4	OFN-100pg-03	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_002.qgd	1	<input type="checkbox"/> Print	
4	5	OFN-100pg-04	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_003.qgd	1	<input type="checkbox"/> Print	
5	6	OFN-100pg-05	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_004.qgd	1	<input type="checkbox"/> Print	
6	7	OFN-100pg-06	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_005.qgd	1	<input type="checkbox"/> Print	

**NOTE**

[Fill Series] can be also executed for the selected range of cells.

Copying Settings

Settings can be copied to each individual column of a Batch Table.

This section describes how to copy and set [Sample Name].

1 Enter [Sample Name] of the first row.

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	Report Fi
1	2	STD	0:Unknown	IT QT	TEST_Ultra.qgm	Repeat_000.qgd	1	<input type="checkbox"/> Print	
2	3	OFN-100pg-02	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_001.qgd	1	<input type="checkbox"/> Print	
3	4	OFN-100pg-03	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_002.qgd	1	<input type="checkbox"/> Print	
4	5	OFN-100pg-04	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_003.qgd	1	<input type="checkbox"/> Print	
5	6	OFN-100pg-05	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_004.qgd	1	<input type="checkbox"/> Print	
6	7	OFN-100pg-06	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_005.qgd	1	<input type="checkbox"/> Print	

2 Right-click the cell at the first row of the [Sample Name] column to copy, and click [Fill Down].

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	Report
1	2	STD	Fill Series		Ultra.qgm	Repeat_000.qgd	1	<input type="checkbox"/> Print	
2	3	OFN-100pg-02			Ultra.qgm	Repeat_001.qgd	1	<input type="checkbox"/> Print	
3	4	OFN-100pg-03	Fill Down		Ultra.qgm	Repeat_002.qgd	1	<input type="checkbox"/> Print	
4	5	OFN-100pg-04			Ultra.qgm	Repeat_003.qgd	1	<input type="checkbox"/> Print	
5	6	OFN-100pg-05	Cut	Ctrl+X	Ultra.qgm	Repeat_004.qgd	1	<input type="checkbox"/> Print	
6	7	OFN-100pg-06	Conv	Ctrl+C	Ultra.qgm	Repeat_005.qgd	1	<input type="checkbox"/> Print	

The content of the first row will be copied to all cells.

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	Report Fi
1	2	STD	0:Unknown	IT QT	TEST_Ultra.qgm	Repeat_000.qgd	1	<input type="checkbox"/> Print	
2	3	STD	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_001.qgd	1	<input type="checkbox"/> Print	
3	4	STD	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_002.qgd	1	<input type="checkbox"/> Print	
4	5	STD	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_003.qgd	1	<input type="checkbox"/> Print	
5	6	STD	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_004.qgd	1	<input type="checkbox"/> Print	
6	7	STD	2:Control	IT QT ILM	TEST_Ultra.qgm	Repeat_005.qgd	1	<input type="checkbox"/> Print	

**NOTE**

[Fill Down] can be also executed for the selected range of cells.

4.3 Starting Continuous Data Acquisition

This section describes how to start and stop continuous data acquisition, pause continuous data acquisition to edit the Batch Table, and partially execute continuous data acquisition.

4.3.1 Starting Continuous Data Acquisition

1

Click the  (Start) icon on the [Batch] assistant bar. This starts continuous data acquisition according to the Batch Table.

4.3.2 Pausing Continuous Data Acquisition

Continuous data acquisition can be paused so that rows of the Batch Table that have not undergone data acquisition can be edited. This section describes how to delete a row of the Batch Table that has not undergone data acquisition.

- 1 During continuous data acquisition, click the  (Pause/Restart) icon on the [Batch] assistant bar.

Data acquisition of the currently executing row is continued as it is.

- 2 Select the row on the Batch Table to be deleted, right-click the selected row, and click [Delete Row].

Folder: C:\GCMSsolution\Sample\Sample2010\QP2010 Ultra

	Vial#	Sample Name	Sample Type	Analysis Ty	Method File	Data File	Inj. Vo	Report Out	F
1	6	OFN-100pg-01	0:Unknown		EI_TEST_Ultra.qcm	EN_SN_TEST_001.qcd	1	<input type="checkbox"/> Print	
2	6	OFN-100pg-02	2-Control	IT QT ILM	EI_TEST_Ultra.qcm	EN_SN_TEST_002.qcd	1	<input type="checkbox"/> Print	
3	6	OFN-100pg-03	2-Control	IT QT ILM	EI_TEST_Ultra.qcm	EN_SN_TEST_003.qcd	1	<input type="checkbox"/> Print	
4	6	OFN-100pg-04	2-Control	IT QT ILM	EI_TEST_Ultra.qcm	EN_SN_TEST_004.qcd	1	<input type="checkbox"/> Print	
5	6	OFN-100pg-05	2-Control			EN_SN_TEST_005.qcd	1	<input type="checkbox"/> Print	
6	6	OFN-100pg-06	2-Control			EN_SN_TEST_006.qcd	1	<input type="checkbox"/> Print	

Context menu options: Fill Series, Fill Down, Cut (Ctrl+X), Copy (Ctrl+C), Paste (Ctrl+V), Clear, Select All, Copy Row, Add Row, Insert Row, Paste Row, Delete Row.

The selected row is deleted.

- 3 Click the  (Save) icon on the toolbar.
The batch file is overwritten and saved.

- 4 To resume continuous data acquisition, click the  (Pause/Restart) icon again on the [Batch] assistant bar.

Data acquisition is resumed.

NOTE


To cancel data acquisition without resuming it, click the  (Stop) icon.

4.3.3 Stopping Continuous Data Acquisition

- 1 Click the  (Stop) icon on the [Batch] assistant bar.
Continuous data acquisition is stopped.

4.3.4 Partially Executing Continuous Data Acquisition

- 1 Select the row No. of the Batch Table to perform partial execution.
- 2 Click the  (Start) icon on the [Batch] assistant bar.
The [Select Batch Execution Range] sub-window opens.
- 3 Check [Start], and click [Selected Row(s)].



Batch processing of the selected rows is partially executed.



NOTE

Partial execution is not available when TurboMatrixHS is connected.


4.4 Using Batch Queue to Perform Additional Data Acquisition

This software enables data to be acquired consecutively using different batch files. Batch files used for continuous data acquisition are registered to a memory area called "batch queue."

4.4.1 Creating Additional Batch Files

- 1 Start up the [GCMS Analysis Editor] program, and click the  (Batch Processing) icon on the [Real Time] assistant bar.
The [Batch Table] window opens.

2 Create the batch file to add to the batch queue, and save the file under a new name.



4

3 Exit the [GCMS Analysis Editor] program.



NOTE


- Data acquisition is not started if there are duplicate data file names or a non-existent method file has been specified.
- The batch queue function operates only when the [GCMS Analysis Editor] program has been exited.

4.4.2 Adding Batch Files

1


Start up the [GCMS Real Time Analysis] program.

During data acquisition, two windows, the [Acquisition] and [Batch Table] windows are displayed simultaneously.




- 2 Click the [Batch Table] window.**
The content of the toolbar, menu bar and assistant bar is changed.

- 3 Select [Batch Queue] on the [Batch] menu.**



The [Batch Queue] sub-window opens.

- 4 Click [Add], and select the batch file to add.**



The [Batch Queue] sub-window displays a list of registered batch files, and data acquisition is continuously performed on the files in order from the top of the queue. The order of registered batch files can be changed by the [Move Up] / [Move Down] buttons, and registered batch files can be deleted from the list by the [Remove] button.


- 5 After editing the batch file list, click [OK].**
The batch file is added to the batch queue.

4.5 Displaying Remaining Time of Continuous Data Acquisition

Acquisition 1.1.6.2


The remaining time of the current continuous data acquisition can be displayed after the acquisition is started. When batch files are set to the batch queue, the remaining time for completing all batch files is displayed.

- 1 Open the [Batch Table] window, and click [Option] on the [Tools] menu.



The [Setting Options] sub-window opens.

- 2 Open the [Batch Processing] tab page.




3 Select [Display Remaining Time to Batch Processing End.].

4 Click [OK] to close the [Setting Options] sub-window.

5 Click the  (Start) icon on the [Batch] assistant bar.

The remaining time for completing the acquisition is displayed in the instrument monitor.



NOTE

The remaining time display can be configured only when all items below meet in the System Configuration.

- Number of analytical lines: 1
- Autosampler: AOC-20i or AOC-20i+s
- Injection port: SPL
- Detector: MS


4.6 Entering the Ecology Mode after Continuous Data Acquisition 1.6.1

The ecology mode is a function for reducing the amount of electrical power and carrier gas consumed during standby for data acquisition. The instrument can be switched to the ecology mode after all continuous data acquisition operations end.


1 Open the batch file to be used for continuous data acquisition.

2 Select [Enters Ecology Mode when ending Realtime Batch] on the [Batch] menu.

A checkmark is added to the selected menu item.



- 3** Click the  (Start) icon on the [Batch] assistant bar. A message for confirming the ecology mode is displayed. Click [Yes].



The instrument switches to the ecology mode after all continuous data acquisition operations including batch queue end.

**NOTE**

Though the setting to enter the ecology mode after batch processing ends can be cancelled by repeating Step 2, it is recommended that this setting remain enabled at all times.

**NOTE**

In the ecology mode, the [Ecology Mode] sub-window is displayed. To perform operations in other sub-windows in the [GCMS Real Time Analysis] program, cancel the ecology mode.

To cancel the ecology mode, click the [Cancel] button in the [Ecology Mode] sub-window.

After the ecology mode is canceled, the instrument returns to the mode before it switched to the ecology mode.



4.7 Errors During Batch Analysis

This section describes main instrument-related errors that occur during batch analysis, and whether batch analysis is continued or canceled according to each error.

Error Code	Error Message	Real Time Batch
0x0D47	The specified vial does not exist.	Batch analysis is continued.
0x0D55	The specified vial does not exist.	Same as above
0x0D95	Over 25,000 hours of usage. Be sure to perform maintenance on turbo molecular pump.	Same as above
0x0D9F	Filament in ion gauge consumed.	Same as above
0x0DA2	Vacuum level deterioration. Filament turned off.	Same as above
0x0DA3	Filament turned on again.	Same as above
0x0DAD	Filament 1 consumed.	Same as above

Error Code	Error Message	Real Time Batch
0x0DAE	Filament 2 consumed.	Same as above
0x0DB6	Ion gauge cannot be controlled normally.	Same as above
0x0DA7	Vacuum level deterioration.	Same as above
0x0D98	Filament 1 consumed. Filament switched.	Batch analysis is canceled.
0x0D99	Filament 2 consumed. Filament switched.	Same as above
0x0DA0	Detector saturated.	Same as above(*1)
0x0D9C	Trap/Total current is insufficient.	Same as above
0x0D9D	R.F. power source error.	Same as above
0x0D9E	Detector high-voltage power supply error.	Same as above
0x0DB1	Total current is insufficient.	Same as above
0x0DB0	Trap current is insufficient.	Same as above
0x0DB2	Conversion dynode power source error.	Same as above
----	(AOC-20i/s-related error)(*2)	Same as above
0x0D80 - 0x0D91	(Instrument communications-related error)(*3)	Same as above
0x0D92	Vacuum system is not ready.	Same as above
0x0F00	Connection to instrument failed.	Same as above
0x0F01	Process stopped due to active instruments.	Same as above
0x0F02	Received error status.	Same as above
0x0F03	Received unexpected status.	Same as above
0x0F04	Failed to initialize instruments.	Same as above
0x0F06	Communication Hardware is not connected.	Same as above
0x4300	CID gas flow rate may be insufficient or unstable.	Same as above
0x4301	CID gas flow rate may be insufficient or unstable.	Same as above
0x4302	CID gas flow rate may be insufficient or unstable.	Same as above
0x4303	The connection with the GC unit failed.	Same as above
0x4304	Detected carrier gas leak. Cooling down column oven, and then turning off flow controller.	Same as above
0x4305	Detected carrier gas leak. Cooling down column oven, and then turning off flow controller.	Same as above
0x4306	Vacuum level doesn't come to be enough. Vacuum system will shut down automatically.	Same as above
0x4307	An abnormality occurred in the vacuum system. The vacuum system will stop automatically.	Same as above
0x4308	Turbo molecular pump 1 isn't ready.	Same as above
0x4309	Ion gauge can not be controlled normally.	Same as above
0x430A	Filament in ion gauge consumed.	Same as above

Error Code	Error Message	Real Time Batch
0x430B	Ion gauge can not be controlled normally.	Same as above
0x430C	Ion source temperature can not be controlled normally.	Same as above
0x430D	Ion source temperature sensor error.	Same as above
0x430E	Ion source temperature can not be controlled normally.	Same as above
0x430F	Ion source temperature can not be controlled normally.	Same as above
0x4310	Conversion dynode power source error.	Same as above
0x4311	Detector high-voltage power supply error.	Same as above
0x4312	Detector saturated.	Same as above
0x4313	An error occurred in the high-frequency power supply (Q1 Main rod).	Same as above
0x4314	An error occurred in the high-frequency power supply (Collision Cell).	Same as above
0x4315	An error occurred in the high-frequency power supply (Q3 Main rod)	Same as above
0x4316	The temperature may be unstable (Q1 RF capacitor).	Same as above
0x4317	The temperature may be unstable (Q3 RF capacitor).	Same as above
0x4318	The temperature may be unstable (Q1 RF capacitor).	Same as above
0x4319	The temperature sensor is disconnected (Q1 Main rod).	Same as above
0x431A	The temperature may be unstable (Q3 RF capacitor).	Same as above
0x431B	The temperature sensor is disconnected (Q3 Main rod).	Same as above
0x431C	Temperature of DI probe can not be controlled normally.	Same as above
0x431D	DI probe is not set or something wrong.	Same as above
0x431E	Temperature of DI probe can not be controlled normally.	Same as above
0x431F	Temperature of DI probe can not be controlled normally.	Same as above
0x4320	Trap/Total current is insufficient.	Same as above
0x4321	Filament 1 consumed. Filament switched.	Same as above
0x4322	Filament 1 consumed. Data acquisition stopped.	Same as above
0x4323	Filament 2 consumed. Filament switched.	Same as above
0x4324	Filament 2 consumed. Data acquisition stopped.	Same as above
0x4325	Ion gauge can not be controlled normally. Data acquisition will stop.	Same as above

Error Code	Error Message	Real Time Batch
0x4326	Ion gauge can not be controlled normally. Data acquisition will stop.	Same as above
0x4327	Q1 RF Heater has stopped.	Batch analysis is continued.
0x4328	Q3 RF Heater has stopped.	Same as above
0x4329	Q1 RF Fan has stopped.	Same as above
0x432A	Q3 RF Fan has stopped.	Same as above
0x432B	Main Fan has stopped.	Same as above
0x432D	Detected carrier gas leak. Stop ion source heater.	Batch analysis is canceled.
0x433D	Filament 1 consumed.	Same as above
0x433E	Filament 2 consumed.	Same as above
0x433F	Vacuum level. Deterioration. Heater in MS turned off.	Same as above
0x4340	Vacuum level. Deterioration. Heater in MS turned off.	Same as above
0x4341	Vacuum level. Deterioration. High voltage and filament are immediately switched off.	Same as above
0x4342	Vacuum level. Deterioration. High voltage and filament are immediately switched off.	Same as above
0x4343	Vacuum level. Deterioration. Ion gauge is immediately switched off.	Same as above
0x4347	Vacuum level. Deterioration. High voltage and filament are immediately switched off.	Same as above

(*1)When the [The batch analysis does not stop even when the signal to the detector is saturated.] checkbox is selected on the [Batch Processing] tab of the [Setting Options] dialog box, this error is not counted as an error.

(*2) Refer to Troubleshooting section in GC-2030 Instruction Manual (PDF file) for details.

(*3) Refer to Error Message Help for details.



NOTE

There are other errors in addition to the above that cancel batch analysis or allow it to continue. Refer to GCMS Help for more information.

5

Postrun Analysis

This chapter describes how to display the results of acquired data and how to set data processing parameters during postrun analysis. In the [Data Analysis] window, one data file is analyzed.

Reference

Multiple data files are analyzed in the [Quant Browser] window. For details, see "6 Quant Browser".

5.1 [Data Analysis] Window

The [Data Analysis] window has the following views.

- [Chromatogram View] which comprises a full chromatogram that displays the entire chromatogram after data acquisition, and a zoomed chromatogram that displays an enlarged chromatogram
- [Spectrum View] for displaying the spectrum of the retention time (or retention time range) specified in [Chromatogram View]
- [Quantitative View] for displaying the quantitative chromatogram of the compound selected in [Compound Table View], the calibration curves, and reference ion information
- [Compound Table View] for displaying the settings or identification results of the Compound Table in the data file.


5

5.1.1 Opening the [Data Analysis] Window

1 Start up the [GCMS Postrun Analysis] program.

2 Click the  (Qualitative) icon or  (Quantitative) icon on the [Postrun] assistant bar.

The [Data Analysis] window opens.



No.	Name	Explanation
①	Chromatogram View	The upper level displays the TIC or MIC of all measurement times after data acquisition. The lower level displays the enlarged TIC, MIC and MC. The sample information is displayed on the right side.
②	Spectrum View	Displays the spectrum of the retention time (or retention time range) specified in [Chromatogram View]. When data is acquired with multiple events, spectra for all events registered to the group to which a specified retention time belongs are displayed in a stack.
③	Quantitative View	Displays the quantitative chromatogram of the compound selected in [Compound Table View], the calibration curves, and reference ion information.
④	Compound Table View	Displays the settings or identification results of the Compound Table in the data file, and the settings or grouping results of the Grouping Table.

5.2 Displaying Chromatograms and Spectra

This section describes how to display chromatograms and spectra in the [Data Analysis] window.

- 1 Click the  (Qualitative) icon on the [Postrun] assistant bar.

2 Drag-and-drop the data file from the [Data] tab page in the [Data Explorer] sub-window onto the [Data Analysis] window.

The screenshot shows the GCMS Postrun Analysis software. On the left, the 'Data Explorer - Data' window lists several data files. An arrow indicates that the file 'pept_14_003.qgd' is being dragged from this list to the 'Data Analysis - Untitled' window. The 'Data Analysis' window displays a Total Ion Chromatogram (TIC & MIC) with a prominent peak at 3.436 minutes. Below the chromatogram is a mass spectrum plot showing relative intensity versus m/z. A table on the right side of the window lists peak data including ID#, Name, Type, ISTD G, m/z, Ret. Time, Ret. Index, and Unit.


The content of the data files is displayed in the [Data Analysis] window.

The screenshot shows the GCMS Postrun Analysis software. On the left, the 'Data Explorer - Data' window lists several data files. The 'Data Analysis - scan_pept.QGD' window on the right displays a Total Ion Chromatogram (TIC & MIC) with multiple peaks. Below the chromatogram is a mass spectrum plot showing relative intensity versus m/z. A table on the right side of the window lists peak data including ID#, Name, Type, ISTD G, m/z, Ret. Time, Ret. Index, and Unit.

ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit
1	Dichlorine	Target	0	181.00	6.550	0	ug/L
2	Fenbutacab	Target	0	150.00	9.667	0	ug/L
3	Simazine	Target	0	201.00	10.992	0	ug/L
4	Phosphamide	Target	0	173.00	11.433	0	ug/L
5	Diazinon	Target	0	304.00	11.692	0	ug/L
6	TPN	Target	0	266.00	11.692	0	ug/L
7	Imidacloprid	Target	0	244.00	12.017	0	ug/L
8	Fenitrothion	Target	0	277.00	13.075	0	ug/L
9	Thiobencab	Target	0	100.00	13.367	0	ug/L
10	Isochlorogenic	Target	0	290.00	15.133	0	ug/L
11	Isochlorogenic	Target	0	313.00	15.592	0	ug/L
12	OMP	Target	0	317.00	16.467	0	ug/L
13	EPN	Target	0	157.00	17.642	0	ug/L

5.2.1 Chromatogram View

[Chromatogram View] displays chromatograms and the instrument status curves.



No.	Name	Explanation
①	Chromatogram information	Displays the range indicating TIC and MIC. Displays the range of the group number of the zoomed chromatogram.
②	Chromatogram label	Displays the event No., m/z and magnification of the currently displayed chromatogram.  NOTE When a reference data file is open, the reference data file No. and m/z of the chromatogram are also displayed.
③	Graph comment	Displays the retention time, scan number and intensity of the mouse position, the column oven temperature when sample is injected from the GC, and the DI temperature when sample is injected from the DI.
④	Peak	Moves the vertical line to the peak top of the adjacent peak, and displays the spectrum of that retention time in [Spectrum View] when the peaks of the chromatogram are detected.
⑤	Scan	Displays the spectrum of the retention time of the immediately preceding or following scan in [Spectrum View].
⑥	Group	Displays the chromatogram of the immediately preceding or following group in the zoomed chromatogram.
⑦	Full chromatogram	Displays the TIC or MIC of all measurement times after data acquisition. Dragging the pointer across the chromatogram displays an enlarged view of the range indicated by dragging.  NOTE When a reference data file is open, the chromatogram of the reference data file is also displayed.
⑧	Zoomed chromatogram	Displays the chromatogram according to the Fragment Table. Dragging the pointer across the chromatogram displays an enlarged view of the range indicated by dragging.  NOTE When a reference data file is open, the chromatogram of the reference data file is also displayed.

No.	Name	Explanation
9	Sample information	Displays the sample information. "There is an error log" is displayed when an error is logged during data acquisition. When a reference data file is open, the reference data file No. and file name are displayed.
10	(vertical line)	Displays a vertical line at the retention time by which spectra are extracted.
11	Peak top comment	Displays the information of the detected peak. Items displayed in the [Chromatogram View Properties] sub-window can be selected. Some of the items in reference data files cannot be displayed as peak top comments. Reference For details, refer to Help.

NOTE

Set the view ratio between the full chromatogram and zoomed chromatogram, and full chromatogram and sample information display/hide selection in the [Chromatogram View Display Setting] sub-window.


To open the [Chromatogram View Display Setting] sub-window, click [Display Settings] - [Chromatogram View] on the [View] menu.

5

5.2.2 Spectrum View

[Spectrum View] displays the spectrum of the retention time (or retention time range) specified in [Chromatogram View].

In the case of data files acquired by measuring with multiple events, spectra for all events registered to the group to which a specified retention time belongs are displayed in a stack.




No.	Name	Explanation
1	Spectrum information	Displays the event No. of the spectrum, acquisition mode, mass according to acquisition mode, retention time, scan No., and retention index. NOTE Mass according to acquisition mode refers to the mass of the precursor ion when the spectrum in the product ion scan mode is displayed, the mass of the product ion in the precursor ion scan mode, and the "Losses of" mass in the neutral loss scan mode.

No.	Name	Explanation
②	Graph comment	Displays the mass, peak intensity and relative intensity at the mouse position.
③	Base peak	Displays m/z and intensity values for base peaks.
④	Regist	This is displayed on the upper right of a spectrum for each event. The currently displayed spectra for the corresponding events are registered to the spectrum process table.
⑤	Spectrum	<p>Displays the spectrum of the retention time at the vertical line of [Chromatogram View] according to the spectrum format parameters.</p> <p> NOTE</p> <p>In the case of events in the MRM and product ion scan acquisition modes, the horizontal axis of the displayed spectrum becomes the scale of the Q3.</p> <p>In the case of events in the precursor ion scan and neutral loss scan acquisition modes, the horizontal axis of the displayed spectrum becomes the scale of the Q1.</p>
⑥	Active mark	This is displayed on the left side of the active spectrum. A similarity search can be performed on the spectrum, the library can be registered to the spectrum, and the MS Table can be displayed.

 **NOTE**

Whether or not a spectrum is displayed for each event when measurement has been performed by multiple events can be set in the [Spectrum Graph Display Setting] sub-window. Also, the number of spectra displayed in a single sub-window and display/hide of the standard spectrum of the Compound Table are set in this sub-window.



By default, [The number of display spectrum] is set to "1". When analyzing data scanned using multiple events (data acquired by Product Ion Scan, Precursor Ion Scan, or Neutral Loss Scan), set a high value, "10" for example, to display all events in the [Spectrum View].


To open the [Spectrum Graph Display Setting] sub-window, click [Display Setting] - [Spectrum View] on the [View] menu.

 **Reference**

For details on how to display the standard spectrum, see ["5.2.6 Displaying the Standard Spectrum in \[Spectrum View\]"](#).

5.2.3 Quantitative View

[Quantitative View] displays the quantitative chromatogram, calibration curves, and reference ion information of the compound selected in [Compound Table View].



5

No.	Name	Explanation
①	Quantitative chromatogram	Displays the target ion of the compound selected in [Compound Table View], and the chromatogram of the reference ion. NOTE When a reference data file is open, its quantitative chromatogram is displayed in a stack under the measurement data file.
②	Chromatogram label	Displays the <i>m/z</i> of the currently displayed quantitative chromatogram according to the ID chromatogram Fragment Table.
③	Peak top comment	Displays the information of the detected peak. Items to be displayed can be selected in the [Chromatogram View Properties] sub-window.
④	Calibration curve	Displays the calibration curve of the compound selected in [Compound Table View]. NOTE When a reference data file is open, sample information is displayed in a stack under the measurement data file instead of the calibration curve.
⑤	Reference ion information	Displays the <i>m/z</i> , area/height/mass, setting% and measurement% of the intensity ratio, and allowable width of the target ion and reference ion. This information can be edited by clicking Edit in [Compound Table View].




NOTE

Set the quantitative chromatogram and calibration curve, and display/hide selection of the reference ion information in the [Quantitation View Display Setting] sub-window.


To open the [Quantitation View Display Setting] sub-window, click [Display Settings] - [Quantitation View] on the [View] menu.

5.2.4 Displaying the Spectrum of a Specific Retention Time

- 1 In [Chromatogram View], place the mouse pointer over the retention time to be extracted as the spectrum, and double-click it.




The spectrum at that retention time position is displayed in [Spectrum View].




5.2.5 Displaying the Chromatogram of a Specific m/z

- 1 Double-click the m/z of the chromatogram to display in [Spectrum View].



That m/z is registered to the Fragment Table, and the chromatogram is displayed in the zoomed chromatogram.




5.2.6 Displaying the Standard Spectrum in [Spectrum View]

The standard spectrum currently registered to the Compound Table can be displayed at the lowermost level of [Spectrum View].

1 Click [Spectrum View] - [Display Setting] on the [View] menu.


The [Spectrum Graph Display Setting] sub-window opens.



2 Select [Display Standard Spectrum], and click [OK].

The standard spectrum is displayed under the measured spectrum.

The range of the mass axis is aligned with the measurement spectrum of event 1 and the standard spectrum, whichever is larger.



3 Drag on the spectrum to enlarge it.

Both the measured spectrum and standard spectrum are displayed zoomed.

NOTE

To initialize zoom, right-click the spectrum, and select [Initialize Zoom] on the displayed menu. Both the measured spectrum and standard spectrum are initialized.

5.2.7 Comparing the Quantitative Chromatogram with Other Measurement Data

The chromatograms of up to three reference data files can be displayed in [Chromatogram View] and [Quantitative View], and the chromatograms of multiple data, such as standard spiked samples, can be compared on the same sub-window.

Opening Reference Data Files

- 1 Click [Open Reference Data File] on the [File] menu.
- 2 In the [Open Data File] sub-window, select the reference data file, and click [Open].
The chromatogram is displayed in [Chromatogram View] and [Quantitative View].

 **NOTE**

- GCMS data files and data files measured with GC-2010 detectors can be selected. Data files measured with GC-2030 detectors cannot be selected.
- Analysis processing cannot be performed on the chromatograms of reference data files.

Closing Reference Data Files


- 1 On the [Close Reference Data File] sub-menu of the [File] menu, click the file to close.
 **NOTE**
Even if the reference data file is closed, the data file open in the [Data Analysis] window does not close.

5.2.8 Comparing Chromatograms with the Intensity Axis of the Quantitative Chromatogram Fixed

When the reference data file is open, the quantitative chromatograms of the measurement data file and reference data file are displayed in a stack in [Quantitative View].

The intensity axis of quantitative chromatograms displayed in [Quantitative View] can be fixed to the identified peaks in a specified data file. By fixing the intensity axis, whether or not the chromatogram contains the target component and whether or not the peak of the component is at the detection limit or below can be checked at a glance.

- 1 In [Quantitative View], right-click the chromatogram to fix the intensity axis, and click [Fix the Intensity Axis to this Data.].
Select the data file used as the detection limit such as data file for a calibration curve level 1.



The intensity axis of the quantitative chromatogram is fixed to the intensity of the selected data file.


5

2 Check presence of the target component.

- 1 In [Compound Table View], click the row of compound.

ID#	Name	Tune	ISTD G	m/z	Ret Time	Ret Index	Unit	Ref Ions	Conc 1
1	Dichlorvos	Target	0	185.00	6.250	0	ug/L	109.00	100
2	Fenobucarb	Target	0	150.00	3.667	0	ug/L	121.00	100
3	Simazine	Target	0	201.00	10.992	0	ug/L	186.00	100

- 2 In [Quantitative View], check whether or not the peak is at the detection limit or below.
In the figure below, it can be seen that the peak is at the detection limit or above since the chromatogram is displayed exceeding the maximum value of the intensity axis of the graph.



5.3 Qualitative Processing

This section describes how to set peak integration parameters, perform qualitative peak integration on chromatograms to detect peaks, search spectra from the library, and identify components.



NOTE

Qualitative processing (e.g. peak detection or registration of spectra to the Spectrum Process Table) can only be performed on event 1.

5.3.1 Flow of Qualitative Processing

1. Peak detection

Detect peaks to extract the spectrum for searching from the library. Peaks are detected manually and automatically.



Reference

See ["5.3.2 Performing Qualitative Peak Integration Automatically"](#) and ["5.3.4 Performing Qualitative Peak Integration Manually"](#).

2. Spectrum search from library

Compare spectra registered to the library file with the spectra of detected peaks to display library spectra as a list in order of highest similarity. Components are identified by comparing the spectra of detected peaks with the list.



Reference


See ["5.3.6 Searching Spectra from the Library"](#).

5.3.2 Performing Qualitative Peak Integration Automatically


This section describes how to set peak integration parameters to automatically detect peaks.

1


Click **[Peak Integrate]** on the **[Qualitative]** menu, and select the *m/z* to perform peak integration with on the sub-menu.



2 In the [Qualitative Parameters] sub-window, set the peak integration parameters, and click [OK].



The chromatogram of the specified m/z undergoes peak integration and the peak is detected. The baseline, peak detection mark, and peak top comment are displayed on the detected peak.



5


Integration Time Program

If the peak integration parameters must be changed according to the retention time to detect peaks, set an integration time program.

1
2

Click  at the peak integration parameters.

Create an integration time program, and click [OK].



- At [Time], enter the retention time at which to execute the command.
- At [Command], select the processing content.
- Set values for processing at [Value], if required.

Reference

For details about each of the parameters, refer to Help.


NOTE

- The [Load Chrom] sub-window opens by clicking the [Load Chrom] button. When a data file is selected and [Open] is clicked, TIC is displayed at the top of the sub-window. When a chromatogram is clicked with the [Time] cell selected, the retention time for the clicked part of the chromatogram is set to the cell.
- Clicking the [Simulate] button executes a simulation of peak integration on the chromatogram using the content set in the peak integration program.

5.3.3 Checking the Peak Integration Results

Detected peaks can be checked in the Peak Table. The Peak Table is divided into tabs for each m/z .

- 1 Click the  (Qualitative Table) icon on the [Qualitative] assistant bar.
- 2 Select the tab of the desired m/z to check the result.



Peak#	Ret Time	Start Tm	End Tm	m/z	Area	Area%	Height	Height%	A/H	Mark	Name
1	6.249	6.217	6.300	TIC	5186839	6.64	4605017	8.50	1.13		
2	6.450	6.417	6.492	TIC	422599	0.54	385796	0.71	1.10		
3	9.666	9.625	9.733	TIC	9888348	12.80	7628646	14.08	1.31		
4	10.992	10.950	11.058	TIC	6578994	8.43	4751729	8.77	1.39		
5	11.434	11.392	11.492	TIC	5286172	6.75	3517725	7.23	1.34		
6	11.893	11.492	11.658	TIC	10204642	13.07	6973621	12.88	1.46	V	
7	11.892	11.850	11.967	TIC	6906138	8.85	4753300	8.78	1.45		
8	12.012	11.975	12.092	TIC	8406725	10.77	5101016	9.42	1.85		
9	13.079	13.042	13.100	TIC	2709332	3.47	2022774	3.74	1.34		
10	13.270	13.217	13.342	TIC	9630070	12.41	5965101	11.02	1.62		
11	15.132	15.092	15.208	TIC	6425176	8.23	4032305	7.45	1.59		
12	15.930	15.542	15.642	TIC	2408346	3.09	1524894	2.82	1.59		
13	16.464	16.417	16.517	TIC	2100303	2.69	1338162	2.47	1.57		
14	17.642	17.683	17.883	TIC	1764218	2.28	1152385	2.13	1.53		

A compound name can be displayed as a peak top comment on the chromatogram by entering the name in the [Name] cell.


5.3.4 Performing Qualitative Peak Integration Manually

This section describes how to perform peak integration manually when peaks cannot be detected even by changing the peak integration parameters.

Detection of Peaks

- 1 In the zoomed chromatogram, click the chromatogram label of the m/z whose peak is to be detected.

The corresponding chromatogram is displayed by a bold line for a preset duration.



2 Move the mouse cursor over the chromatogram, and drag from the start point to the end point using the right mouse button with the [Shift] key held down.

During the drag operation, the mouse pointer changes to a vertical line.

The base line is drawn joining the start and end points, and the peaks are detected.




NOTE

- This operation can be also executed by clicking [Manual Peak Integrate] on the [Qualitative] menu, and selecting the m/z to perform peak integration with on the sub-menu.
- By dragging from the start point to the end point using the right mouse button with the [Ctrl] key held down, peaks are detected with the baseline kept horizontal.


Deletion of Peaks

1 In the zoomed chromatogram, click the chromatogram label of the m/z whose peak is to be deleted.

The corresponding chromatogram is displayed by a bold line for a preset duration.



2 Double-click between the peak start time to end time with the [Shift] key held down.



The peak is deleted.



5

Movement of Detection Point


1 In the zoomed chromatogram, move the mouse pointer to near the peak detection point (start point or end point).

The mouse pointer changes to a triangle.



2 Click the detection point and drag-and-drop onto the time where the detection point is to be moved.

During the drag operation, the mouse pointer changes to a vertical line, and the baseline correction line (red) moves.



The detection point moves to the dropped retention time, and peak integration is performed.

NOTE


The baseline can be easily changed when the peaks are tailing or leading.

Do not move the detection point beyond the adjacent peak.


Splitting of Peaks

1 In the zoomed chromatogram, click the chromatogram label of the m/z whose peak is to be split.


The corresponding chromatogram is displayed by a bold line for a preset duration.



- 2** Click [Split Peak] on the [Qualitative] menu.
The mouse pointer changes to a cross.



- 3** Click at the retention time where the peak is to be split.
The peak is split at that retention time and is detected as respective peaks.



- 4** Click [Split Peak] on the [Qualitative] menu again to end peak splitting.

5

5.3.5 Performing Average/Subtraction Processing on Spectra

When a peak is selected and its spectrum is displayed in [Chromatogram View], average/subtraction processing is performed on the spectrum according to the spectrum processing parameters, and the resulting spectrum is displayed.

The same operation can be also performed when making a Spectrum Process Table from a Peak Table or when displaying the spectrum of a compound selected in [Compound Table View].

This section describes how to manually perform average/subtraction processing on a spectrum to obtain a clearer spectrum.

A processed spectrum can be easily reproduced later if it is registered to the Spectrum Process Table beforehand.


Averaging of Spectra

- 1** Click  (Average Spectrum) on the toolbar.

- 2** In [Chromatogram View], drag from the start to end time to be averaged.
During the drag operation, the mouse pointer changes to a vertical line.



The average spectrum within the dragged time range is displayed in [Spectrum View].




NOTE


Averaging near the top of peaks on the chromatogram is more effective.

Subtraction of Spectra

- 1** Click  (Subtract Spectrum) on the toolbar.
- 2** In [Chromatogram View], double-click at the time to perform subtraction.




The spectrum obtained by subtracting the spectrum at the specified time in background is displayed in [Spectrum View].



**NOTE**

In the case of the following peaks, set the part marked by arrows as the background processing position.



Averaging and Subtraction of Spectra


1

Click  (Average Subtract Spectrum) on the toolbar.


2

In [Chromatogram View], drag the time range to perform subtraction on.

5




The spectrum obtained by subtracting the spectrum in the dragged time range in background is displayed in [Spectrum View].

**NOTE**

- Subtracting near the start or end of the peaks on the chromatogram is more effective.
- To further perform subtraction on already subtracted peaks, either click  (Subtract Spectrum) or double-click at the time to perform subtraction with the [Shift] key held down, or click  (Average_Subtract Spectrum) and drag the time range to perform subtraction in with the [Shift] key held down.

Display of Spectrum Calculation Results


The top of [Spectrum View] displays which time range was calculated.



Display Text String	Explanation
[10.000]	Displays a spectrum with retention time 10.000 min.
[10.000 -> 10.100]	Averages from retention time 10.000 min through 10.100 min.
[10.000] - [10.500]	Subtracts the spectrum at 10.500 min from the spectrum with retention time 10.000 min.
[10.000 -> 10.100] - [10.500]	Averages from retention time 10.000 min through 10.100 min, and subtracts the spectrum at 10.500 min from that range.
[10.000 -> 10.100] - [10.500 -> 11.000]	Averages from retention time 10.000 min through 10.100 min, averages from retention time 10.500 min through 11.000 min from that, and subtracts the result of averaging as background.
[10.000 -> 10.100] - [10.500] - [10.040]	Averages from retention time 10.000 min through 10.100 min, subtracts the spectrum at 10.500 min from that range, and further subtracts the spectrum at 10.040 min.

Registering Spectra in [Spectrum View] to the Spectrum Process Table


1 Display the spectrum for the event to register, and click the  icon.



The spectrum of [Spectrum View] is added to the end of the Spectrum Process Table.

Checking the Spectrum Process Table

- 1 Click the  (Qualitative Table) icon on the [Qualitative] assistant bar.
- 2 Select the [Spectrum Process] tab.



	Spectrum			Background			Search	Report	Event	Name
	Ret. Time	Start Tm	End Tm	Ret. Time	StartRT	EndRT				
1	...	6.433	6.450	6.442	6.408	6.475		✓	1	
2	...	6.633	6.650	6.642	6.608	6.675		✓	1	
3	...	6.933	6.950	6.942	6.908	6.975		✓	1	
4	...	8.562	8.578	8.570	8.537	8.603		✓	1	
5	...	9.895	9.912	9.903	9.870	9.953		✓	1	
6	...	11.228	11.245	11.237	11.203	11.278		✓	1	
7	...	11.670	11.687	11.678	11.645	11.720		✓	1	
8	...	11.828	11.845	11.837	11.803	11.862		✓	1	
9	...	12.145	12.162	12.153	12.120	12.195		✓	1	
10	...	12.262	12.278	12.270	12.195	12.312		✓	1	
11	...	12.445	12.462	12.453	12.312	12.503		✓	1	
12	...	13.358	13.375	13.367	13.333	13.417		✓	1	
13	...	13.525	13.542	13.533	13.492	13.575		✓	1	


5.3.6 Searching Spectra from the Library

You can compare the spectra of peaks detected in a chromatogram with the spectra in library files that have a high similarity to identify components, or filter library spectra using a specific search index as the search key.


This section describes how to perform a similarity search for spectra in [Spectrum View] and in the Spectrum Process Table.

Setting Similarity Search Parameters

- 1 Click [Qualitative Parameters] on the [Qualitative] menu.



- 2 On the [Similarity Search] tab page, set the search parameters, and click [OK].



- 1 Up to five library files can be set. The minimum similarity of the library spectrum to search can be set to each library file.

Reference

For details about each of the parameters, refer to Help.

NOTE

When a library file name is not set in the similarity search parameters, the similarity search is performed using the default library. The default library is set in the sub-window that opens by clicking [Set Default Library] on the [Tools] menu in the [Data Analysis] window.

Performing a Similarity Search on the Spectra in [Spectrum View]

- 1** Click the  (Similarity Search) icon on the [Qualitative] assistant bar. The [Similarity Search Results] sub-window is displayed.



Right double-clicking on mass spectra also executes a similarity search.

The screenshot shows the 'Similarity Search Results' window. At the top is a table with columns: Hit, Similarity, Register, Compound Name, Mol Wt, Formula, and Library. The first row is highlighted. Below the table is a button 'Copy marked compound name to the Spectrum Process Table' and a 'Spectrum Process Table #1' field. Below this are three mass spectra plots. The top plot is the target spectrum. The middle plot is labeled '1 : 255 : Propyzamide \$\$ Benzamide, 3,5-dichloro-N-(1,1-dimethyl-2-propynyl)- \$\$ Kerb \$\$ Kerb 50w \$\$ N-(1,1-Dimethylpropynyl)-3,5-dichlorobenzan' and includes a chemical structure. The bottom plot is labeled '2 : 255 : Propyzamide \$\$ Benzamide, 3,5-dichloro-N-(1,1-dimethyl-2-propynyl)- \$\$ Kerb \$\$ Kerb 50w \$\$ N-(1,1-Dimethylpropynyl)-3,5-dichlorobenzan' and also includes a chemical structure. Numbered callouts 1-10 point to: 1. Similarity Search icon; 2. Similarity column; 3. Register column; 4. Compound Name column; 5. Copy button; 6. Spectrum Process Table field; 7. Target mass spectrum plot; 8. Reference mass spectrum plot; 9. Reference mass spectrum plot; 10. A vertical scrollbar on the left side of the spectra area.


No.	Name	Explanation
1	Search results list	Displays the results of the similarity search in order of highest similarity.
2	Similarity	Displays how similarity among spectra with 100 used as the match reference.
3	Register	Select the target compound checkbox when registering a compound name to the Spectrum Process Table.

No.	Name	Explanation
4	Copy marked compound name to the Spectrum Process Table	Copies compound names whose [Register] checkbox is selected to the Spectrum Process Table.  NOTE This button is displayed only when a search has been made from the Spectrum Process Table.
5	(Spectrum Processing Table # spin buttons)	Switches the row numbers of the Spectrum Process Table. The target spectrum switches to the spectrum of the row after switching, and the hit list and search results are also switched according to the content of that spectrum.  NOTE This button is displayed only when a search has been made from the Spectrum Process Table.
6	Register Target Spectrum to Spectrum Process Table	Adds the target spectrum to the end of the Spectrum Process Table.  NOTE This button is displayed only when a search has been made from [Spectrum View].
7	Target spectrum	Displays the target spectrum.
8	Library spectrum	Displays the spectrum of the compound selected at the search results list.
9	Comparison/Information/Subtraction	The display content is switched by clicking [Compare]/[Information] / [Subtracted] on the [View] menu. [Compare] displays the library spectrum. [Information] displays the compound information for the library spectrum. [Subtracted] displays the spectrum obtained by subtracting the library spectrum from the target spectrum.
10	(Library spectrum spin buttons)	These are used to switch the library spectrum.

To print the results of the similarity search, click [Search Results] on the [Report] menu in the [Similarity Search Results] sub-window.

Performing a Similarity Search on Spectra in the Spectrum Process Table

- 1 Click the  (Qualitative Table) icon on the [Qualitative] assistant bar. The [Qualitative Table] sub-window is displayed.



	Spectrum			Background			Search	Report	Event	Name
	Ret. Time	Start Tm	End Tm	Ret. Time	StartRT	EndRT				
1	...	6.433	6.450	6.442	6.408	6.475		<input checked="" type="checkbox"/>	1	
2	...	6.633	6.650	6.642	6.608	6.675		<input checked="" type="checkbox"/>	1	
3	...	6.933	6.950	6.942	6.908	6.975		<input checked="" type="checkbox"/>	1	
4	...	9.562	9.578	9.570	9.537	9.603		<input checked="" type="checkbox"/>	1	
5	...	9.895	9.912	9.903	9.870	9.953		<input checked="" type="checkbox"/>	1	
6	...	11.228	11.245	11.237	11.203	11.278		<input checked="" type="checkbox"/>	1	
7	...	11.670	11.687	11.678	11.645	11.720		<input checked="" type="checkbox"/>	1	
8	...	11.828	11.845	11.837	11.803	11.862		<input checked="" type="checkbox"/>	1	
9	...	12.145	12.162	12.153	12.120	12.195		<input checked="" type="checkbox"/>	1	
10	...	12.262	12.278	12.270	12.195	12.312		<input checked="" type="checkbox"/>	1	
11	...	12.445	12.462	12.453	12.312	12.503		<input checked="" type="checkbox"/>	1	
12	...	13.358	13.375	13.367	13.333	13.417		<input checked="" type="checkbox"/>	1	
13	...	13.525	13.542	13.533	13.492	13.575		<input checked="" type="checkbox"/>	1	

- 2** On the [Spectrum Process] tab page, click [Search All Table] on the [Similarity Search] menu.



The similarity search is executed, and "Done" is displayed at the [Search] row.

	Spectrum			Background			Search	Report	Event	Name
	Ret. Time	Start Tm	End Tm	Ret. Time	StartRT	EndRT				
1	---	6.433	6.450	6.442	6.408	6.475	Done	✓	1	
2	---	6.633	6.650	6.642	6.608	6.675	Done	✓	1	
3	---	6.933	6.950	6.942	6.908	6.975	Done	✓	1	
4	---	8.562	8.578	8.570	8.537	8.603	Done	✓	1	
5	---	9.895	9.912	9.903	9.870	9.953	Done	✓	1	
6	---	11.228	11.245	11.237	11.203	11.278	Done	✓	1	
7	---	11.670	11.687	11.678	11.645	11.720	Done	✓	1	
8	---	11.828	11.845	11.837	11.803	11.862	Done	✓	1	
9	---	12.145	12.162	12.153	12.120	12.195	Done	✓	1	
10	---	12.262	12.278	12.270	12.195	12.312	Done	✓	1	
11	---	12.445	12.462	12.453	12.312	12.503	Done	✓	1	
12	---	13.358	13.375	13.367	13.333	13.417	Done	✓	1	
13	---	13.525	13.542	13.533	13.492	13.575	Done	✓	1	

NOTE

- Correct search results are not displayed when a similarity search is performed on spectra obtained by MS/MS measurement.
- If the [Copy Compound Name of Hit #1 to Spectrum Process Table] is selected, the compound name of the compound having the highest similarity is copied to the [Name] cell after the similarity search is performed.

- 3** Double-click the similarity-searched row.

The [Similarity Search Results] sub-window is displayed.

NOTE

When a similarity search is performed from the Spectrum Process Table, the search results are also saved when the data file is saved.

5.3.7 Printing Graph Images


This section describes how to print chromatograms and spectra that have been analyzed in the [Data Analysis] window.

- 1** Display or analyze chromatograms or spectra in the [Data Analysis] window.

- 2** Select [Print Image] on the [File] menu, and click [Preview].

The [Report] window opens.

3 Check the chromatogram and spectra in the preview sub-window, and click [Print].



NOTE

- To edit a report format, select [Print Image] on the [File] menu, and click [Edit Format]. The [Report] window opens, and the format can be edited. If an edited format is saved, it can be used when next printing graph images.
- To return an edited format to the default format, select [Print Image] on the [File] menu, and select [Initialize Format]. This returns the format to its default state when GCMSsolution was installed.

5.4 Quantitative Processing

This section describes how to set quantitative parameters to perform quantitative processing.

5.4.1 Flow of Quantitative Processing

1. Creation of Compound Tables

Create the Compound Table required for setting peak integration parameters, making calibration curves and performing quantitative processing.

There are several ways of creating a Compound Table.

Reference

See the following sections:

- In case of data obtained by measurement by event 1 in Q3 Scan acquisition mode, etc., see ["5.4.2.1 Creating Compound Tables Using the Wizard"](#).
- In case of data acquired in MRM or SIM acquisition mode, see ["5.4.2.2 Creating Compound Tables from MS Tables in the Method"](#).
- See ["5.4.2.3 Copying from Spreadsheet Format Files to Create Compound Tables"](#).

2. Creation of calibration curves

Make calibration curves based on the Compound Table created at 1.

Reference

See ["5.4.3 Creating Calibration Curves"](#).

3. Data acquisition of unknown samples

Acquire data of unknown samples whose concentrations are to be calculated.

Reference

See ["3 Data Acquisition"](#).

4. Execution of quantitative calculation

Perform quantitative peak integration on the data file of the unknown samples to identify compounds and calculate their concentrations.

Reference

See ["5.4.4 Performing Quantitative Processing on Unknown Samples"](#), ["5.4.6 Performing Quantitative Peak Integration Manually"](#) and ["5.4.7 Performing Identification Manually"](#).

5. Confirmation of quantitative results

Check the identification results and quantitative calculation results.

Reference

See ["5.4.5 Checking the Quantitative Calculation Results"](#).


5.4.2 Creating Compound Tables

5.4.2.1 Creating Compound Tables Using the Wizard

This section describes how to create Compound Tables from data measured by Event 1 in Q3 scan mode, etc., using the Compound Table Creation Wizard.


The event numbers for Compound Tables created with the Compound Table Creation Wizard will all be "1".

- 1 Load the data file to the [Data Analysis] window.
- 2 Click the  (Create Compound Table) icon on the [Postrun] assistant bar.
- 3 Click the  (Wizard(New)) icon on the [Compound] assistant bar.
- 4 Select "1" or "2", and click [Next].




No.	Name	Explanation
①	Use current Spectrum Process Table	Creates a Compound Table using the spectra already registered to the Spectrum Process Table. Before starting up the Compound Table Wizard, register the required spectra to the Spectrum Process Table.
②	Integration of TIC	Performs peak integration on TIC, and creates a Compound Table using the spectra of detected peaks. Clicking the [Integration Param] button opens the sub-window for setting the qualitative peak integration parameters. Clicking the [Spectrum Format] button opens the sub-window for setting the spectrum process parameters.  Reference For a detailed description of the parameters, refer to Help.

- 5 Check the number of spectra registered to [Spectrum Process Table], and click [Next].




6 Select the [Proc] column checkboxes of the spectra to be processed, and click [Next].



No.	Explanation
①	Select the checkboxes for the retention times of the spectra to be used for creating the Compound Table. The peaks near the retention times registered here are targeted for calculation.
②	The spectra of the rows selected in the table at 1 can be checked.

5

7 Set quantitative method, calibration curve and other quantitative parameters, and click [Next].




No.	Explanation
①	Sets the quantitative calculation method and type of calibration curve.
②	Sets the concentration unit or format of the concentration value.
③	By setting [Grouping], quantitative calculation of groups can be performed in addition to quantitative calculation for each compound.

Reference

For details about the parameters, refer to Help.

8 Set the concentration and measurement ion, and click [Next].



No.	Name	Explanation
①	Standard	Sets the concentration of the standard sample. When the concentration differs according to the component, correct each concentration after the Wizard ends.
②	# of Reference Ions	Enter the number of reference ions to be used for peak identification.
③	Decimal for mass	Determines the number of digits past the decimal point for the target ion and reference ion. Selecting [1 Decimal] improves sensitivity.
④	Internal Standard	Sets the concentration of the internal standard when the internal standard method is selected as the quantitative calculation method.

9 Set the type, compound name, target ion, and reference ion of each component, and click [Next] when the information pertaining to all compounds has been entered.

Compound Table Wizard 6/7

Edit all fields, as necessary. To change the type, place the cursor in the type column and select a new type from the drop down list in the field.

	Type	m/z	Rel. Inten.
1	Target Ion	109	100.00
2	Ref. Ion	185	23.71
3	Ref. Ion	79	22.16
4	Not used	145	8.83
5	Not used	187	7.92
6	Not used	220	5.11
7	Not used	222	3.40
8	Not used	76	3.32
9	Not used	93	3.22
10	Not used	110	3.12

1: ID# field
2: Type dropdown menu
3: Table
4: Compound Name field

- 1 Switch the ID of the compound to set.
- 2 Select [Target] from the [Type] list.




NOTE

Select [I.S.] for the internal standard substance when the internal standard method is selected as the quantitative calculation method. In addition, the same number must be set to the [ISTD Group] cells of the internal standard substance and the compounds that use that substance for the calculation in the Compound Table.

ID #	ISTD G	m/z
1	1	109.00
2	1	121.00
3	1	201.00
4	1	173.00
5	1	179.00
6	1	266.00
7	1	91.00
8	2	125.00
9	2	100.00
10	2	118.00
11	2	105.00

3 Change the type and mass.

- To change the type, click the cell of the type to change, and select one of [Target Ion], [Ref. Ion] and [Not used].
- To change the mass, click the arrow button that is displayed by clicking the cell of the mass to change. The mass spectrum is displayed. Drag near the spectrum peak to enlarge it. The mass can be set by double-clicking the spectrum peak. The mass spectrum also is switched by changing the compound ID.




- 4 The name registered to the Spectrum Process Table is automatically entered. If the entry field is blank, enter the compound name.

**NOTE**

Clicking the button on the right if the similarity search has already been executed on the Spectrum Process Table displays a list of names of library registered substances found in the search.

10

In the confirmation sub-window, click [Finish].



The Compound Table will be created. If necessary, correct the content of the Compound Table.


ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit	Ref. Ions	Conc. 1	E
1	RT-6.250	Target	0	109.00	6.250	0		185.00-79.00	1	1
2	RT-9.667	Target	0	121.00	9.667	0		150.00-91.00	1	1
3	RT-10.992	Target	0	201.00	10.992	0		186.00-68.00	1	1
4	RT-11.433	Target	0	173.00	11.433	0		175.00-145.00	1	1
5	RT-11.592	Target	0	179.00	11.592	0		137.00-152.00	1	1
6	RT-11.892	Target	0	266.00	11.892	0		264.00-268.00	1	1
7	RT-12.008	Target	0	91.00	12.008	0		204.00-123.00	1	1
8	RT-13.075	Target	0	125.00	13.075	0		109.00-277.00	1	1
9	RT-13.267	Target	0	100.00	13.267	0		72.00-125.00	1	1
10	RT-15.133	Target	0	118.00	15.133	0		162.00-189.00	1	1
11	RT-15.592	Target	0	105.00	15.592	0		177.00-77.00	1	1
12	RT-16.467	Target	0	317.00	16.467	0		319.00-236.00	1	1
13	RT-17.642	Target	0	157.00	17.642	0		169.00-141.00	1	1

**NOTE**

To correct the content, click Edit at the top right of the Compound Table to set to the edit mode. When editing is finished, click View to set to the display mode.

11 Click the  (Save Compound Table) icon on the [Compound] assistant bar.
The method file used when this data was acquired is selected.

12 Click [Save].




5

5.4.2.2 Creating Compound Tables from MS Tables in the Method

This section describes how to create Compound Tables from data acquired in the MRM or SIM acquisition mode.

Complete the Compound Table by setting the retention time of the Compound Table while observing the chromatogram peaks, and setting the reference ion ratio while observing the spectrum.

- 1** Load the data file to the [Data Analysis] window.
- 2** Click [Update compound table m/z by MRM/SIM event] on the [Compound Table] menu.
A confirmation message is displayed.
- 3** Click [OK].



The Compound Table enters the edit mode.

If the Compound Table contains the same compounds as the MS Table's MRM or SIM events, update the m/z value of the Compound Table with the m/z setting values of the MS Table.

If there are no identical compounds, add them as new compounds to the Compound Table.

NOTE


Update the m/z value set in MS Table Ch1 as the quantitation ion, and update the m/z values set in Ch2 to Ch6 as the reference ions.

4 Set [Ret. Time] in the Compound Table.

- 1 Click the [Ret. Time] cell.

ID#	Name	Type	ISTD G	m/z	Ret Time	Ret. Index	Unit
1	Dichlorvos	Target	0	185.00:93.00	0.000	0	
2	Fenobucarb	Target	0	150.00:125.00	0.000	0	
3	Simazine	Target	0	201.00:173.00	0.000	0	
4	Propyzamide	Target	0	173.00:150.00	0.000	0	
5	Diazinon	Target	0	304.00:179.00	0.000	0	
6	TPN	Target	0	266.00:250.00	0.000	0	
7	Iprobenfos	Target	0	204.00:150.00	0.000	0	
8	Fenitrothion	Target	0	277.00:260.00	0.000	0	
9	Thiobencarb	Target	0	100.00:80.00	0.000	0	
10	Isoprothiolan	Target	0	290.00:260.00	0.000	0	
11	Isoxathion	Target	0	313.00:300.00	0.000	0	
12	CNP	Target	0	317.00:290.00	0.000	0	
13	EPN	Target	0	157.00:110.00	0.000	0	

- 2 Click the peak top of the chromatogram.



The retention time of the clicked position is set to the [Ret. Time] cell.

ID#	Name	Type	ISTD G	m/z	Ret Time	Ret. Index	Unit
1	Dichlorvos	Target	0	185.00:93.00	5.398	0	
2	Fenobucarb	Target	0	150.00:125.00	0.000	0	
3	Simazine	Target	0	201.00:173.00	0.000	0	
4	Propyzamide	Target	0	173.00:150.00	0.000	0	
5	Diazinon	Target	0	304.00:179.00	0.000	0	
6	TPN	Target	0	266.00:250.00	0.000	0	
7	Iprobenfos	Target	0	204.00:150.00	0.000	0	
8	Fenitrothion	Target	0	277.00:260.00	0.000	0	
9	Thiobencarb	Target	0	100.00:80.00	0.000	0	
10	Isoprothiolan	Target	0	290.00:260.00	0.000	0	
11	Isoxathion	Target	0	313.00:300.00	0.000	0	
12	CNP	Target	0	317.00:290.00	0.000	0	
13	EPN	Target	0	157.00:110.00	0.000	0	

5 Set the reference ion ratio.

- 1 Click  icon at the right end of the [Ref. Ions] cell.

ID#	Ret. Time	Ret. Index	Unit	Ref. Ions
1	5.398	0		185.00:85.00
2	0.000	0		
3	0.000	0		201.00:138.00
4	0.000	0		173.00:140.00
5	0.000	0		304.00:152.00
6	0.000	0		266.00:220.00
7	0.000	0		204.00:100.00
8	0.000	0		277.00:109.00
9	0.000	0		100.00:68.00
10	0.000	0		290.00:240.00
11	0.000	0		313.00:260.00
12	0.000	0		317.00:240.00
13	0.000	0		157.00:141.00

The [Reference Ion Ratio] sub-window is displayed.

- 2 While observing the spectrum in [Spectrum View], set the relative intensity ratio of the reference ion when the quantitation ion is 100.

Reference Ion Ratio - ID#: 1 - Target m/z: 185.00>93.00

Ref. #	m/z	Ratio	Allowance
1	185.00>85.00	30.00	Default
2	TIC	0.00	Default
3	TIC	0.00	Default
4	TIC	0.00	Default
5	TIC	0.00	Default

6 Repeat steps 4 and 5 to set the retention times and reference ion ratios of all compounds.

7 For all compounds, set the standard sample concentrations for each level.

ID#	Ret. Index	Unit	Ref Ions	Conc. 1	Conc. 2	Conc. 3	S
1	0		180.00>85.00	10	500	1000	
2	0		130.00>110.00	10	500	1000	
3	0		180.00>138.00	10	500	1000	
4	0		150.00>140.00	10	500	1000	
5	0		280.00>152.00	10	500	1000	
6	0		240.00>220.00	10	500	1000	
7	0		180.00>100.00	10	500	1000	
8	0		250.00>109.00	10	500	1000	
9	0		90.00>68.00	10	500	1000	
10	0		270.00>240.00	10	500	1000	
11	0		290.00>260.00	10	500	1000	
12	0		285.00>240.00	10	500	1000	
13	0		130.00>141.00				

5

Reference

Set the number of concentration levels using the quantitative parameters. See "Changing the Quantitative Parameters" in ["5.4.3 Creating Calibration Curves"](#) for the procedures for revising the quantitative parameters.


8 Click the  icon in the Compound Table.

The Compound Table returns to the display mode.

9 Click the  (Save Compound Table) icon on the [Compound] assistant bar.

Select the method file used when this data was acquired.

10 Click [Save].




5.4.2.3 Copying from Spreadsheet Format Files to Create Compound Tables

This section describes how to copy information set to spreadsheet format files that are preset with data acquisition settings to a Compound Table.


- 1 Load the data file to the [Data Analysis] window.
- 2 Click the  (Quantitative) icon on the [Postrun] assistant bar. The [Data Analysis] window is displayed in the quantitative processing layout.
- 3 Click the  icon in [Compound Table View]. The Compound Table enters the edit mode.
- 4 Copy items from the spreadsheet format file to the Compound Table.

- 1 Select and copy the content in the spreadsheet format file. Make sure that the columns in the compound table and in the spreadsheet format file are arranged in the same order.



ID#	Name	Type	ISTD (Group)	m/z	Ret. Time	Ret. Index	Unit	Ref. Ions
4	Octafluoronaphthalene	Target	1	272>241	4.78			272>222-272>172
5	Dichlorvos	Target	1	185>109	5.91			185>93-185>63
6	Fenobucarb	Target	1	150>121	9.665			150>103-150>77
7	Hexachlorobenzene	Target	1	284>249	10.805			284>214-284>212
8	Simazine	Target	1	201>173	11.095			201>186-201>158
9	Propyzamide	Target	1	173>145	11.555			173>109-173>74
10	Diazinon	Target	1	304>179	11.61			304>162-304>195
11	Chlorothalonil	Target	1	266>231	11.71			266>170-266>188
12	Iprobenfos	Target	1	204>91	12.12			204>122-204>171
13	Fenitrothion	Target	1	277>260	13.165			277>109-277>125
14	Thiobencarb	Target	1	257>100	13.5			257>72-257>224
15	Stearic acid, methyl ester	Target	1	298>101	14.885			298>185-298>143
16	Isoprothiolane	Target	1	290>204	15.265			290>118-290>162
17	Isoxathion	Target	1	313>177	15.77			313>130-313>159
18	Chlornitrofen	Target	1	317>287	16.68			317>236-317>224
19	EPN	Target	1	157>110	17.825			157>77-157>139

- 2 Click the [Name] cell in the first row of the Compound Table, and paste the content.



ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit
1	Dichlorvos	Target	1	185.00>93.00	6.439	0	
2	Fenobucarb	Target	1	150.00>125.00	9.906	0	
3	Simazine	Target	1	201.00>173.00	11.236	0	
4	Propyzamide	Target	1	173.00>150.00	11.681	0	
5	Diazinon	Target	1	304.00>179.00	11.832	0	
6	TPN	Target	1	266.00>250.00	12.153	0	
7	Iprobenfos	Target	1	204.00>150.00	12.265	0	
8	Fenitrothion	Target	1	277.00>260.00	13.326	0	
9	Thiobencarb	Target	1	100.00>80.00	13.532	0	
10	Isoprothiolan	Target	1	290.00>260.00	15.389	0	
11	Isoxathion	Target	1	313.00>300.00	15.840	0	
12	CNP	Target	1	317.00>290.00	16.727	0	
13	EPN	Target	1	157.00>110.00	17.903	0	

Rows are automatically added and the content is pasted.

- 5 Set the reference ion ratio.

- 1 Click the [Ref. Ions] cell.

ID#	Ret. Time	Ret. Index	Unit	Ref. Ions	Conc. 1	Event	STD Sp
1	6.439	0		185.00>85.00	1	1	
2	9.906	0			1	1	
3	11.236	0		201.00>138.00	1	1	
4	11.681	0		173.00>140.00	1	1	
5	11.832	0		304.00>152.00	1	1	
6	12.153	0		266.00>220.00	1	1	
7	12.265	0		204.00>100.00	1	1	
8	13.326	0		277.00>109.00	1	1	
9	13.532	0		100.00>68.00	1	1	
10	15.389	0		290.00>240.00	1	1	
11	15.840	0		313.00>260.00	1	1	
12	16.727	0		317.00>240.00	1	1	
13	17.903	0		157.00>141.00	1	1	

The [Reference Ion Ratio] sub-window is displayed.

- Set the relative intensity ratio and allowable width for the reference ion when the quantitation ion is 100.

Ref. #	m/z	Ratio	Allowance
1	185.00>85.00	30.00	Default
2	TIC	0.00	Default
3	TIC	0.00	Default
4	TIC	0.00	Default
5	TIC	0.00	Default

- Repeat steps 1 and 2 to set the reference ion ratios and allowable widths of all compounds.



Hint

To collectively set the mass, ion ratio, and allowable width for the reference ion, copy a text string in the following format to the [Ref. Ions] cell in the Compound Table.

Mass of reference ion 1-Mass of reference ion 2\r\nIon ratio of reference ion 1-Ion ratio of reference ion 2\r\nAllowable width of reference ion 1-Allowable width of reference ion 2

Specify the settings values for reference ions 1 to 5, using a "-" as a separator. The "\r\n" is an item (mass, ion ratio, allowable width) delimiter.

Example:

Copying "185>85-127>80\r\n60-40\r\n10-20" to the [Ref. Ions] cell will configure the following reference ion settings.

"185>85-127>80\r\n60-40\r\n10-20"

ID#	Ret. Time	Ret. Index	Unit	Ref. Ions	Conc. 1	Event	STD Sp
1	6.439	0			1	1	
2	9.906	0			1	1	
3	11.236	0			1	1	
4	11.681	0			1	1	
5	11.832	0			1	1	
6	12.153	0			1	1	
7	12.265	0			1	1	
8	13.326	0			1	1	
9	13.532	0			1	1	
10	15.389	0			1	1	
11	15.840	0			1	1	
12	16.727	0			1	1	
13	17.903	0			1	1	

Ref. #	m/z	Ratio	Allowance
1	185.00>85.00	60.00	10
2	127.00>80.00	40.00	20
3	TIC	0.00	Default
4	TIC	0.00	Default
5	TIC	0.00	Default

6 For all compounds, set the standard sample concentrations for each level.

ID#	Ret. Index	Unit	Ref. Ions	Conc. 1	Conc. 2	Conc. 3	S
1	0		160.00>85.00	10	500	1000	
2	0		130.00>110.00	10	500	1000	
3	0		180.00>138.00	10	500	1000	
4	0		150.00>140.00	10	500	1000	
5	0		280.00>152.00	10	500	1000	
6	0		240.00>220.00	10	500	1000	
7	0		180.00>100.00	10	500	1000	
8	0		250.00>109.00	10	500	1000	
9	0		90.00>68.00	10	500	1000	
10	0		270.00>240.00	10	500	1000	
11	0		290.00>260.00	10	500	1000	
12	0		285.00>240.00	10	500	1000	
13	0		130.00>141.00	10	500	1000	



NOTE

Set the number of concentration levels using the quantitative parameters. See "Changing the Quantitative Parameters" in ["5.4.3 Creating Calibration Curves"](#) for the procedures for revising the quantitative parameters.


7 Click the icon in the Compound Table.

The Compound Table returns to the display mode.

8 Click the (Save Compound Table) icon on the [Compound] assistant bar.

Select the method file used when this data was acquired.

9 Click [Save].




5.4.3 Creating Calibration Curves

This section describes how to create a calibration curve using the method file used for creating the Compound Table.

Opening the [Calibration Curve] Window

- 1 Start up the [GCMS Postrun Analysis] program, and click the  (Calibration Curve) icon on the [Postrun] assistant bar.

The [Calibration Curve] window opens.



No.	Name	Explanation
1	Calibration Curve View	Displays the calibration curve graph and calibration curve information of each component (each group in the case of group quantitation). It also displays a table of the level of each calibration point, preset concentration values and area values (height values).
2	Compound Table View	Displays the settings or quantitation results of the Compound Table (Grouping Table) of the method.
3	Chromatogram View	Displays the quantitative chromatogram of the compound selected in [Compound Table View] for the data file selected in [Data File Tree].
4	Data File Tree	Displays the data files of the standard samples that comprise each calibration point at each level. Calibration points can be manually added by dragging-and-dropping the data files of the standard samples onto the specified level.

Creating Calibration Curves

- 1 Load a method file that is set with a Compound Table.
- 2 Click the [Data] tab page in the [Data Explorer] sub-window, and drag-and-drop the data file for the standard sample onto the data file tree.



NOTE

When there are two or more levels of calibration points, drag-and-drop each of the data file corresponding to each level.

The screenshot shows the GCMSSolution Analyst interface. The 'Data Explorer - Data' window is open, displaying a file tree for 'C:\GCMSSolution\Sample'. A file named 'ipest_std_001.qgd' is selected and being dragged towards the 'Level 1' area of the calibration curve plot. The plot is currently blank, with axes for 'Level' and 'Conc.'. A table on the right side of the window lists various compounds and their properties.

ID#	Name	Type	ISTD G	m/z
1	Diazepam	Target	0	150.0
2	Fenofibrate	Target	0	150.0
3	Simvastatin	Target	0	201.0
4	Propofolone	Target	0	173.0
5	Diazepam	Target	0	304.0
6	PHN	Target	0	266.0
7	Isofenflone	Target	0	234.0
8	Fenofibrate	Target	0	277.0
9	Fluticasone	Target	0	180.0
10	Isofenflone	Target	0	290.0
11	Isofenflone	Target	0	313.0
12	ENP	Target	0	377.0
13	EPN	Target	0	157.0

3

- Click the  (Peak Integration for All Data) icon on the [Calibration] assistant bar. The calibration curve is produced and displayed.

The screenshot shows the GCMSSolution Analyst interface with the calibration curve completed. The 'Data Explorer - Method' window is open, displaying a file tree for 'C:\GCMSSolution\Sample'. A file named 'ipest_std_001.qgd' is selected. The calibration curve plot shows a linear relationship between 'Level' and 'Conc.'. The plot includes a linear regression line and the following equation: $Y = aX^2 + bX + c$, $a = 1.332225e-002$, $b = 400.2147$, $c = 1028.9899$, $R^2 = 1.0$, $R = 0.999999994$. The plot also shows the 'Area' and 'Conc.' for each level. A table on the right side of the window lists various compounds and their properties.

Level	Conc.	Area1
1	100.00	40.025
2	500.00	208.112
3	1,000.00	433.227


Type	m/z	Intern.	Act%	Set%	Disp.
Target	150.00	25412	100.00	100.00	
Int1	109.00	159764	451.16	421.00	

- 4** Click the  (Save) icon on the toolbar, and save the method file.

Changing the Quantitative Parameters

If necessary, change the quantitative parameters.

- 1** Click the  (Quantitative Parameters) icon on the [Calibration] assistant bar. The [Quantitative Parameters] sub-window is displayed.
- 2** Set the parameters on the [Peak Integration] tab page so that the intended compound peaks can be detected.
- 3** On the [Quantitative] tab page, set the calibration curve.



5

No.	Name	Explanation
1	Quantitative Method	<ul style="list-style-type: none"> External Standard: Quantitation is performed using a calibration curve obtained from the absolute quantity (concentration) and the area or height value of the target compound in a standard sample. Internal Standard: An internal standard is added to the sample, the sample is analyzed, and quantitation is performed using the relationship between the relative sensitivity and the quantitative ratio with respect to the internal standard compound.
2	# of Calib. Levels	Input the number of concentration levels of the calibration curve.

No.	Name	Explanation
③	Curve Fit Type	<p>Specifies how to plot the calibration curve.</p> <ul style="list-style-type: none"> • Linear: Determines the calibration curve as a straight line from the obtained values. • Point to point: Points are connected by a broken line. No formula is displayed for point to point calibration curves. • Quadratic: Curve is fit to each point using a quadratic equation. This requires at least three points on the calibration curve. For two points or less, the curve is calculated as linear. • Mean RF: First, it determines straight lines passing through the origin and each point. Then it finds the simple average of the slopes for each line. Consequently, the resulting calibration curve always passes through the origin.
④	Zero	Select either [Not Forced] or [Force Through]. Normally, select [Not Forced].
⑤	Weighted Regression	<p>A typical least squares method of plotting calibration curves could result in a quantitation error that is larger the lower the concentration at the calibration point. In general, when the calibration curve has a large dynamic range (maximum concentration is at least 50 times higher than the minimum quantitation limit), formulas are weighted to reduce the weight of higher concentration points of the calibration curve. Typically, formulas are optimized by checking the correlation coefficient and contribution ratio.</p> <ul style="list-style-type: none"> • [1/C²]: Formulas are weighted by the inverse of the concentration value squared. • [1/C]: Formulas are weighted by the inverse of the concentration value. • [1/A²]: Formulas are weighted by the inverse of the area value squared (or height value when a height is specified for the data used). • [1/A]: Formulas are weighted by the inverse of the area

4 After the settings are changed, click [OK].

The calibration curve will be revised in accordance with the new parameter settings.


Printing Calibration Curves

The procedure for printing calibration curves that are produced is as follows.

1 Select [Print Image] on the [File] menu, and click [Preview].

The [Report] window opens.

2 Check the calibration curve in the preview sub-window, and click [Print].




5

5.4.4 Performing Quantitative Processing on Unknown Samples

This section describes how to perform quantitative calculation on unknown samples using the method file used for creating the calibration curve.

- 1 In the [Data Analysis] window, open the data file to undergo quantitative processing.
- 2 Click the  (Quantitative) icon on the [Postrun] assistant bar.
- 3 Click the  (Load Method) icon on the [Quantitative] assistant bar. The [Load Method] sub-window opens.
- 4 Select the method file that has already been used to create a calibration curve, and click [Open].



The calibration curve and the Compound Table are loaded to the data file.

- 5** Click the  (Peak Integration) icon on the [Quantitative] assistant bar. Peak integration, identification processing and quantitative calculation are executed.

5.4.5 Checking the Quantitative Calculation Results

Check the quantitative calculation results on the [Results] tab page in [Compound Table View].

- 1** Click the [Results] tab in [Compound Table View].

ID#	Name	Conc	Ret. Time	Type	m/z	Area
1	Dichlorvos	Ratio of reference ion does not match.				
2	Fenobucarb	Ratio of reference ion does not match.				
3	Simazine	0.00000	10.990	Target	201.00	957657
4	Propyzamide	0.00000	11.434	Target	173.00	946982
5	Diazinon	0.00000	11.590	Target	304.00	524013
6	EPN	No peak is found in Window/Band range.				
7	Iprobenfos	0.00000	12.011	Target	204.00	1479624
8	Fenitrothion	0.00000	13.072	Target	277.00	311707
9	Thobencarb	0.00000	13.271	Target	100.00	4086151
10	Isoxathion	0.00000	15.129	Target	290.00	413007
11	Isoxathion	0.00000	15.595	Target	313.00	106667
12	CNP	Ratio of reference ion does not match.				
13	EPN	0.00000	17.642	Target	157.00	393079

The Identified Results Table is displayed.



NOTE

The following text strings are displayed according to cause for compounds that failed to be identified.

Display Text String	Explanation
No peak is detected.	The peak has not been detected by quantitative peak integration.
No Peak is found in Window/Band range.	The peak has not been detected within the retention time width of the identification parameters.
Ratio of reference ion does not match.	The difference between the setting value of the reference ion ratio and measured value exceeds the allowable width, preventing the peak from being identified.
Under the minimum similarity index.	The peak has not been identified due to the calculated similarity value being the setting value or less when mass pattern matching is set in the identification parameters.

Display Text String	Explanation
No peak is identified.	The automatic identification results have been deleted manually.

**NOTE**

When a row on the Compound Table is selected, the quantitative chromatogram of its compound, or calibration curve and reference ion information are displayed in [Quantitation View].


5.4.6 Performing Quantitative Peak Integration Manually

This section describes how to perform peak integration manually when peaks cannot be identified or detected even by setting the peak integration parameters.

1 Move the mouse pointer over the chromatogram, and drag from the start point to the end point using the right mouse button with the [Shift] key held down.

During the drag operation, the mouse pointer changes to a vertical line.


The base line is drawn joining the start and end points, and a peak is detected and identified.

**NOTE**


- The above operation can be also executed by clicking [Manual Peak Integrate] on the [Quantitative] menu.
- By dragging from the start point to the end point using the right mouse button with the [Ctrl] key held down, peaks are detected with the baseline kept horizontal.
- When [Peak Picked] is selected at [Ref. Ions Based On] in the identification parameters, peak integration is performed on the chromatograms of the quantitation ion and reference ion. When [Spectrum] is selected, peak integration is performed on only the chromatogram of the quantitation ion.
To perform peak integration on only the chromatogram of a specific m/z , click the chromatogram label of that m/z and perform peak integration by following the description in step 1.

5.4.7 Performing Identification Manually

- 1 Right-click in the range between the start and end times of the currently detected peak with the [Shift] and [Ctrl] keys held down.



A peak is identified.




NOTE

- The above operation can be also executed by clicking [Manual Identification] on the [Quantitative] menu.
- When the outside of the range between the start and end times of the peak is right-clicked, the peak of the quantitation ion having the retention time closest to the right-clicked position is identified.

5.4.8 Detecting Unidentified Peaks

This section describes how to detect peaks outside the target compound registered to the Compound Table and perform a similarity search on its spectrum. This feature is used for analyzing contamination and unforeseen components.

1 Select [TIC] or [MIC] on the [Pickup Unidentified Peaks] sub-menu on the [Quantitative] menu.




Unidentified peaks in the selected chromatogram are detected and registered to the Identified Results Table.

ID#	Name	Conc	Ret. Time	Type	m/z	Area
1	Dichlorvos					
2	Fenobucarb					
3	Simazine	0.00000	10.990	Target	201.00	857657
4	Propyzamide	0.00000	11.434	Target	173.00	946982
5	Diazinon	0.00000	11.590	Target	304.00	524013
6	TPN					
7	Iprobenfos	0.00000	12.011	Target	204.00	1479624
8	Fenitrothion	0.00000	13.072	Target	277.00	311707
9	Thiobencarb	0.00000	13.271	Target	100.00	4086151
10	Isoprothiolan	0.00000	15.129	Target	290.00	413007
11	Isosathion	0.00000	15.585	Target	313.00	106667
12	CNP					
13	EPN	0.00000	17.642	Target	157.00	393079
14	Tentatively Id	0.00000	9.666	Unknown	TIC	9988348
15	Tentatively Id	0.00000	11.892	Unknown	TIC	6906138
16	Tentatively Id	0.00000	6.248	Unknown	TIC	5186829
17	Tentatively Id	0.00000	16.464	Unknown	TIC	2100203
18	Tentatively Id	0.00000	6.450	Unknown	TIC	422599

2 Select [Similarity Search for Compound Table] on the [Quantitative] menu.

- 3** Select [Unknown] on the [Compound Table Search] tab page in the [Quantitative Parameters] sub-window, and click [OK].



The similarity search is executed, and [Finish] is displayed at the [Search] column.



NOTE

Double-clicking the row of a compound displays the result of the similarity search in the [Similarity Search Results] sub-window.


5.4.9 Calculating S/N

This section describes how to set parameters for calculating noise and S/N.


1
2

Click the  (Quantitative Parameters) icon on the [Quantitative] assistant bar.

Click [Noise Calculation] on the [Peak Integration] tab page.



3 Set the calculation method and noise range, and click [OK].



No.	Explanation
1	Select one of [ASTM], [RMS] or [Peak to Peak].
2	Selecting [Specified Range] calculates noise from the ranges set to [From] and [To]. If [Calculation Method] is set to [RMS] or [Peak to Peak], selecting [Auto] automatically searches for the range having the smallest noise in the time width set at [Noise Range]. When [Calculation Method] is set to [ASTM], noise is calculated by the range obtained by subtracting the peak detected time from the processing time.

Reference

For details about the parameters, refer to Help.

4 Click the (Peak Integration) icon on the [Quantitative] assistant bar.

Peak integration, identification and quantitative calculation are performed on all IDs.

5 Click the [Results] tab in [Compound Table View].


The S/N calculation results are displayed in the Identified Results Table.

ID#	SI	S/N	Noise From	Noise To	Noise Calc
1	84	3495.94	5.910	6.010	PtoP
2	94	7042.18	9.177	9.277	PtoP
3	78	7109.76	10.713	10.813	PtoP
4	48	1391.01	11.630	11.730	PtoP
5	32	2664.68	11.250	11.350	PtoP
6	72	4358.33	12.270	12.370	PtoP
7	72	4474.33	11.554	11.654	PtoP
8	41	1725.23	13.220	13.320	PtoP
9	81	10287.63	12.767	12.867	PtoP
10	51	2777.95	14.687	14.787	PtoP
11	82	405.44	15.377	15.477	PtoP
12	82	1421.89	16.050	16.150	PtoP
13	72	1696.59	17.220	17.320	PtoP


5.5 Saving New Parameters to Method Files

Save new settings to a method file when data processing parameters (e.g. qualitative parameters or quantitative parameters) have been changed.

- 1 Click **[Save Method As]** on the **[File]** menu.
The **[Save Method As]** sub-window opens.



- 2 Enter the file name, and click **[Save]**.
The **[Select Method Parameters]** sub-window opens.
- 3 Set the method parameters in the **[Select Method Parameters]** sub-window, and click **[OK]**.



- 1 By **[Current Settings]**, the currently displayed latest parameter settings are saved. By **[Acquired Settings]**, the parameter settings that were used at data acquisition are saved. In this example, select **[Current Settings]**.
- 2 Select the checkboxes of the parameters to save.

5.6 Printing the Content of Data Files

This section describes how to print the content of data files in the [Data Analysis] window.

- 1 Load the data file to the [Data Analysis] window, and click [Report] on the [File] menu.**
The [Data Report] window is displayed. When a report format is specified during data acquisition or batch processing, that report format is set.

 **NOTE**

The [Data Report] window is displayed also when the  (Report) icon on the [Qualitative] assistant bar or on the [Quantitative] assistant bar is clicked.

 **Reference**

For details on how to paste new report items, see "[8.2 Creating Report Formats](#)".

- 2 Place the report items, and click the  (Preview) icon.**

 **NOTE**


Report format files can be also dragged-and-dropped from the [Report] tab page in the [Data Explorer] sub-window onto the [Data Report] window.

- 3 Check the report in the preview sub-window, and click [Print].**

5.7 Performing Continuous Postrun Analysis


This section describes how to create a Batch Table from a data file and perform postrun analysis.

- 1 Click the  (Batch Processing) icon on the [Postrun] assistant bar.
- 2 Click the  (Select Data Files) icon on the [Batch] assistant bar.
The [Select Data File] sub-window is displayed.
- 3 Select the data file to register to the Batch Table, and click .



- 4** At [Selected Data File], select the data file, change the order of the data files by the [Up] and [Down] buttons, and click [OK].

The Batch Table is created.



- 5** In the [Data Analysis] sub-window that opens by clicking the [Analysis Type] cell, set the analysis details, and click [OK].



Reference

For details, refer to Help.

- 6** Click the (Start) icon on the [Batch] assistant bar. Continuous postrun analysis is executed.

NOTE

A batch file made when performing continuous data acquisition can be loaded to perform continuous postrun analysis. At this time, settings relating to data acquisition are ignored.

This page is intentionally left blank.

6

Quant Browser

By using the [Quant Browser] window, the quantitative results (concentration, area, height, etc.) of multiple data files can be displayed for each compound set to the Compound Table in the method. Quantitative results in the list undergo statistical calculation (average, %RSD, maximum, minimum, standard deviation) as necessary before being displayed. These results can be printed out as a summary report.

This chapter describes how to check and edit the quantitative results of multiple data files, and perform post-run analysis collectively on multiple data files.


6.1 [Quant Browser] Window

The [Quant Browser] window has four views: [Quantitative Result View] for displaying the quantitative calculation results of each data file, [Compound Table View] for displaying the Compound Table and the quantitative results of each data file, [Chromatogram View] for displaying chromatograms and sample information, and [Calibration Curve View] for displaying the calibration curves of compounds.

6.1.1 Opening the [Quant Browser] Window

- 1 Start up the [GCMS Browser] program, and click the  (Quant Browser) icon on the [Browser] assistant bar.

The [Quant Browser] window opens.



The screenshot displays the GCMS Browser interface with the [Quant Browser] window open. The window is divided into four main panes:

- 1 Assistant Bar:** Located on the left side, it contains various icons for navigation and analysis, including 'Top', 'Modify Calibration Curve', 'Quantitative Parameters', 'Peak Integration for All Data', 'Peak Table', 'Summary Report', and 'Guide'.
- 2 Quantitative Result View:** A table showing quantitative data for multiple samples. The table has columns for Data File Name, Sample Type, Level, Ret. Time, and Area. The data is summarized as follows:

Data File Name	Sample Type	Level	Ret. Time	Area
PEST_STD_000.QGD	Standard	1	8.607	97.541
PEST_STD_000.QGD	Standard	2	8.607	516.485
PEST_STD_000.QGD	Standard	3	8.607	1,108.323
pest_unknown_001.qgd	Unknown	---	8.607	300,898
Average			8.607	504.331
%RSD			0.00000	86.568166
Maximum			8.607	1,108.323
Minimum			8.607	97.541
Std. Dev.			0.000000	436,596.3863
- 3 Compound Table View:** A table listing compounds with columns for DB Name, Type, ISTD G, m/z, Ret. Time, and Ret. The data is summarized as follows:

DB Name	Type	ISTD G	m/z	Ret. Time	Ret.
Dichlorvos	Target	0	185.00	6.250	
Fenobacarb	Target	0	150.00	9.867	
Simazine	Target	0	201.00	10.952	
Phoxonyl	Target	0	173.00	11.423	
Diazinon	Target	0	304.00	11.652	
TPH	Target	0	266.00	11.862	
Imidacloprid	Target	0	204.00	12.017	
Fenitrothion	Target	0	277.00	13.075	
Thiofencarb	Target	0	100.00	13.257	
Inoprotolol	Target	0	250.00	15.153	
Isoaathion	Target	0	313.00	15.552	
CHP	Target	0	317.00	16.467	
EPN	Target	0	157.00	17.642	
- 4 Chromatogram View:** A chromatogram plot showing intensity versus time. The x-axis is labeled 'Time' and ranges from 9.25 to 10.00. The y-axis is labeled 'Intensity' and ranges from 0.0 to 7.5 (x100,000). A single prominent peak is visible at approximately 9.50 minutes. Below the plot is a table with columns for Type, m/z, Inten., Act%, Set%, and Disp. The data is summarized as follows:

Type	m/z	Inten.	Act%	Set%	Disp.
Target	150.00	221411	100.00	100.00	
Ref Ion	121.00	679396	306.82	297.00	

No.	Name	Explanation
①	Quantitative Result View	Displays the quantitative calculation results (e.g. concentration, retention time, area) of multiple data files, and also displays their statistical values (average, %RSD, maximum, minimum, standard deviation). Clicking  changes the ID# and switches the displayed compound.
②	Compound Table View	Displays the Compound Table and Group Table in the currently open method file. The identification results and grouping results of the data file selected in [Quantitative Result View] are displayed.
③	Calibration Curve View	Displays the calibration curve of the ID selected in [Compound Table View] in the currently open method file.
④	Chromatogram View	Displays the chromatogram of the ID selected in [Compound Table View] in the data file selected in [Quantitative Result View].

**NOTE**

- The [Quant Browser] window comprises views that are separated by several splitters. These views can be resized by dragging these splitters with the mouse.
- In the [Quant Browser] window, you can collectively check the quantitative calculation results of up to 1024 data files using the same method file.
- When [Save Browsing File As] on the [Layout] menu in the [Quant Browser] window is selected, the name of the currently loaded method file and data file, file sort order, layout information, and other details can be saved as a browsing file (file extension *.qgq).
- When a file is being edited in a different window, it becomes [Read Only] and cannot be edited. To edit the file, close the file in the other window, and open the file again. When there are read only files, neither the method file can be saved nor parameters can be edited.

6.2 Checking the Quantitative Results

In the [Quant Browser] window, multiple data can be easily displayed.

This section describes how to display data in each of the views in the [Quant Browser] window.

6.2.1 Opening Method Files


This section describes how to open a method file and display the content of the data file (sample type is [Standard]) that comprises the calibration curve contained in that method file in each view.

1


Click the [Method] tab in the [Data Explorer] sub-window.

2

Drag-and-drop the method file from the [Data Explorer] sub-window onto the [Quant Browser] window.




The content of the standard sample data file and method file are displayed in the [Quant Browser] window.




6.2.2 Opening Data Files to Display Quantitative Results

This section describes how to open a data file of an unknown sample and display the content of the data file in each view.

- 1 Click the [Data] tab in the [Data Explorer] sub-window.
- 2 Drag-and-drop the data file of an unknown sample from the [Data Explorer] sub-window onto the [Quant Browser] window.



The contents of the data file and method file currently set to the data file are displayed in the [Quant Browser] window.




NOTE

When the data file is dragged-and-dropped from the [Data Explorer] sub-window onto the [Quant Browser] window, multiple files can be selected simultaneously. At this time, if the method file is still not opened, the content of the method file set to the top data file is displayed. When the loaded method file holds calibration curve information, the content of the data file of the standard sample that was used when making that calibration curve is also displayed.


6.2.3 Opening Batch Files to Display Quantitative Results

This section describes how to open batch files to display the contents of the method file and data file currently set to batch files in each view.

- 1 Click the [Batch] tab in the [Data Explorer] sub-window.
- 2 Drag-and-drop the batch file from the [Data Explorer] sub-window onto the [Quant Browser] window.



Rows whose sample type is [Standard] in the batch file are searched for, and where the first [Standard] row is found, the method file used in that row and all of the data files using that method file are loaded.




NOTE

- If no row whose sample type is [Standard] is found in the batch file, a search is performed from the 1st row of the batch file for an [Unknown] sample row, and where the first [Unknown] sample row is found, the method file used in that row and all of the data files that use that method file are loaded.
- To display the content of the method file and data file using a batch file, the data file and method file must be saved beforehand to the same folder as the batch file.
- Method files and data files can be edited in the [Quant Browser] window. Batch files cannot be changed.

6.2.4 Editing [Quantitative Result View]

Switching the Displayed Data Type

The data type can be selectively displayed by one of the [All], [Standard], [Unknown], [Control], and [Spiked]

sample types by clicking the  buttons on the toolbar.

NOTE

[Standard] refers to the data that comprises the calibration curve in the method file. Even if the sample type of the data file is standard sample, the sample type is handled as [Unknown] when the data is not the data comprising the calibration curve of the currently open method file.

Switching Compounds for Displaying the Quantitative Results List

Quantitative results are displayed in [Quantitative Result View] for each compound currently displayed in the Compound Table. The compound is switched by the  (ID#) spin buttons at the top of [Quantitative Result View] or by selecting the row of the target compound in [Compound Table View].

NOTE


- The error between the concentration value found by quantitative calculation performed on data with sample type [Standard] and the set concentration value at the same level in the Compound Table can be displayed.
- The calibration curve of the target compound is displayed in [Calibration Curve View].

Displaying Statistical Calculation Results

When the [Statistic] cell checkboxes in the Quantitative Results Table are selected, statistical calculation (average, %RSD, maximum, minimum, standard deviation) is performed between data with selected cells, and the statistical calculation result items are added to the lowest row of the table. Statistical calculation can be performed on cells other than data file name, sample name, sample ID, sample type, set concentration, level No., mark, and statistical calculation.

NOTE

- [Quantitative Result View] is sorted by file name by clicking [Quantitative Result] - [Sort by Filename] on the [View] menu.
- Items displayed in the [Quantitative Result View] can be selected in the [Table Style] sub-window by clicking [Quantitative Result] - [Table Style] the [View] menu.



6.2.5 Fixing the Intensity Axis of [Chromatogram View] to Check Quantitative Results

The intensity axis of chromatograms displayed in [Chromatogram View] can be fixed to the identified peaks in a specified data file. By fixing the intensity axis, whether or not the chromatogram contains the target component and whether or not the peak of the component is at the detection limit or below can be checked at a glance.

- 1 In [Quantitative Result View], select the data file whose intensity axis is to be fixed.
Select the data file used as the detection limit such as data file for a calibration curve level 1.


- 2 Right-click the selected row, and click [Fix the Intensity Axis to This Data.].

The intensity axis of the quantitative chromatogram is fixed to the intensity of the selected data file.

6



To fix the intensity axis on a data file that has not been loaded on the quantitation browser, click [Intensity Range Setting] - [Chromatogram] on the [View] menu. In the [ID Chromatogram Axis Setting] sub-window, select [Fix to identified peak], and specify [Data File].



3 Check presence of the target component.

- 1 Select the compound ID#1 by the spin buttons in [Quantitative Result View].


Quantitative Result View ID# 1 Dichlorvos					
	Data Filename	Sample Type	Level	Ret. Time	Area
1	PEST_STD_003.QGD	Standard	1	6.211	40,060
2	PEST_STD_002.QGD	Standard	2	6.207	208,667

The intensity axis of [Chromatogram View] is already fixed to the intensity of compound ID#1 of the data file that was used for fixing the intensity axis.

- 2 Click on the toolbar.
The quantitative result for the unknown sample is displayed.
- 3 In [Quantitation Result View], select the data file of the unknown sample.

Quantitative Result View ID# 1 Dichlorvos					
	Data Filename	Sample Type	Level	Ret. Time	Area
1	pest_unknown_001.qgd	Unknown	6.237	101
	Average			
	%RSD			
	Maximum			6.297	101
	Minimum			6.237	101
	Std. Dev.			

- 4 In [Chromatogram View], check whether or not the peak is at the detection limit or below.
In the figure below, it can be seen that the peak is at the detection limit or above since the chromatogram is displayed exceeding the maximum value of the intensity axis of the graph.



- 5 Repeat steps 3 and 4, and check the peaks of data files of other unknown samples.
- 6 Select the compound ID#2 by the spin buttons in [Quantitation Result View].
At this time, the intensity axis of [Chromatogram View] is already fixed to the intensity of compound ID#2 of the data file that was used for fixing the intensity axis.

- Repeat steps 3 to 5, and check the peaks of data files of unknown samples for compound ID#3 onwards.

6.3 Performing Postrun Analysis Collectively on Multiple Data Files

This section describes how to edit [Compound Table View] and perform postrun analysis collectively on multiple data files.

6.3.1 Editing Compound Tables

- Click  (Edit) in [Compound Table View].
- Click the [Param's] tab, select the compound to be changed, and change the setting of [Conc. 1].

In this example, change the [Conc. 1] setting from "100" to "50."

ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit	Ref. Ions	Conc. 1	Conc. 2	Conc. 3	Event
1	Dichlorvos	Target	0	185.00	6.250	0	ug/L	109.00	50	500	1000	1
2	Fenobucarb	Target	0	150.00	9.667	0	ug/L	121.00	50	500	1000	1
3	Simazine	Target	0	201.00	10.992	0	ug/L	186.00	50	500	1000	1
4	Propyzamide	Target	0	173.00	11.433	0	ug/L	255.00	50	500	1000	1
5	Diazinon	Target	0	304.00	11.592	0	ug/L	179.00	50	500	1000	1
6	TPN	Target	0	266.00	11.892	0	ug/L	264.00	50	500	1000	1
7	Iprobenfos	Target	0	204.00	12.017	0	ug/L	91.00	50	500	1000	1
8	Fenitrothion	Target	0	277.00	13.075	0	ug/L	260.00	50	500	1000	1
9	Thiobencarb	Target	0	100.00	13.267	0	ug/L	257.00	50	500	1000	1
10	Isoprothiolan	Target	0	290.00	15.133	0	ug/L	162.00	50	500	1000	1
11	Isoxathion	Target	0	313.00	15.592	0	ug/L	177.00	50	500	1000	1
12	CNP	Target	0	317.00	16.467	0	ug/L	319.00	50	500	1000	1
13	EPN	Target	0	157.00	17.642	0	ug/L	169.00	50	500	1000	1

- Click  (View) in [Compound Table View].
The calibration curve is recreated. At the same time, the quantitative results are re-calculated, and [Quantitative Result View] is updated.

NOTE

- To cancel editing in [Compound Table View], right-click [Compound Table View], and click [Cancel Edit] on the displayed menu.
- Select the items that are displayed in the Compound Table in the [Table Style] sub-window. The [Table Style] sub-window is displayed by right-clicking a Compound Table, and clicking [Table Style] on the displayed menu.

6.3.2 Setting the Peak Integration Parameters

- Click the  (Quantitative Parameters) icon on the [Quantitative] assistant bar.
- Click the [Peak Integration] tab in the [Quantitative Parameters] sub-window, and set each parameter as required.

Reference

For details about the parameters, refer to Help.

6.3.3 Performing Peak Integration Collectively on Multiple Data

- 1 Click the  (Peak Integration for All Data) icon on the [Quant Browser] assistant bar. Peak integration is performed on all compound IDs of all data files, and quantitative calculation follows after that.

 **NOTE**


- To perform analysis on only data files selected in [Quantitative Result View], select [Peak Integration] - [Peak Integrate for all IDs] on the [Process] menu.
- To perform analysis on only the data files currently selected in [Quantitative Result View], select [Peak Integration] - [Peak Integrate by ID] on the [Process] menu.

6.4 Correcting Calibration Curves

- 1 Click the  (Modify Calibration Curve) icon on the assistant bar. The [Calibration Curve] window opens for correcting the calibration curve of the method file currently loaded on the [Quant Browser] window.
- 2 Change the parameters or perform manual peak integration or manual identification to correct the calibration curve.
- 3 Check the calibration curve, and click  (Save) on the toolbar. The calibration curve is saved to the currently open method file.
- 4 Return to the [Quant Browser] window. When the [Quant Browser] window is displayed, quantitative calculation is performed again based on the corrected calibration curve.

6.5 Printing Summary Reports

- 1 Click the  (Summary Report) icon on the assistant bar. The [Report] window opens with the default report format loaded.



NOTE

To change the content of the summary report, double-click the summary report item, and edit the items in the [GCMS Summary(Compound) Properties] sub-window. To save changes made to

summary report items, click the  (Save) button on the toolbar. Changes are saved to the file "MS Quantitative Browser Report.qgr," and this file is used when next printing a summary report.

- 2 Click the  (Preview) icon on the [Report] assistant bar.
- 3 Check the report in the preview sub-window, and click [Print].

6.6 Printing Chromatogram Images

This section describes how to print the chromatogram displayed in the [Quant Browser] window.

- 1 Display the chromatogram in the [Quant Browser] window.
- 2 Select [Print Current Data] on the [File] menu, and click [Preview].
Check the chromatogram in the preview sub-window, and click [Print].

7


Data Comparison

7.1 Opening the [Data Comparison] Window

In the [Data Comparison] window, displays the chromatograms and spectra of the multiple data files, and compares/performs them.

- 1 Start up the [GCMS Postrun Analysis] program, and click the  (Data Comparison) icon on the assistant bar.

The [Data Comparison] window opens.



No.	Item	Explanation
1	Data File Information	Displays the information of the currently displayed data file. To close a data file, right-click on the data file information and click [Close] on the displayed menu.
2	Chromatogram View	The chromatogram of the currently open data file is displayed as a "full chromatogram" or a "zoomed chromatogram." Three view methods are available, Overlay, Stack and Base Shift. [Chromatogram View] displays the MC (MIC and TIC also) of the specified mass extracted from data. The chromatogram of the desired channel is displayed in the case of GC data files.
3	[m/z] Box	Enter the mass of the chromatogram to display in this box. The [m/z] box has a history function. A list of previously set masses is displayed in the pull down list, and the mass can be selected from this list.
4	Event	Set the event No. of the chromatograms and spectra of all data to display. Either enter the event No. directly or select the event No. by clicking the spin buttons.
5	Compound Name	Displays the compound name corresponding to the event specified at [Event] and the group to which the retention time specified on the chromatogram belongs for the instrument method of the data file currently selected by the [Data] button.
6	Move/Parallel Move Buttons	By clicking these buttons, the chromatogram of the selected data can be moved and expanded/reduced as desired in both the time and intensity axis directions.

No.	Item	Explanation
7	Data Selection Buttons	Select the data No. to set to a selected state.
8	Spectrum View	[Spectrum View] displays the spectrum of the retention time specified on the chromatogram. When multiple data sets are loaded, the spectra are displayed in multiple levels. A similarity search can be performed from the menu that is displayed by right-clicking on the spectrum to be processed.


NOTE

- Chromatograms can be displayed shifted by an equal interval by right-clicking on the chromatogram and selecting [Base Shift] on the displayed menu.
- Chromatograms can be pasted to other applications as image files by right-clicking on the chromatogram and selecting [Copy] on the displayed menu.
- To close all data, click [Close] on the [File] menu.

7.2 Displaying Chromatograms and Spectra

This section describes how to display chromatograms and spectra on the [Data Comparison] window.

- 1 From the [Data] tab page in the [Data Explorer] sub-window, drag-and-drop the multiple data files onto the [Data Comparison] window.




The content of the data files is displayed in the [Data Comparison] window.

NOTE

- Data files measured with GC-2030 detectors cannot be selected. Data files measured with GC-2010 detectors can be selected.

7.2.1 Displaying the Spectrum of a Specific Retention Time

- 1 In [Chromatogram View], double-click with the mouse pointer over the retention time to extract the spectrum.



The spectrum of event 1 at that retention time position is displayed in [Spectrum View].



NOTE

To display the chromatogram or spectrum of other events, switch to the desired chromatogram or spectrum using the [Event] spin buttons.

7


7.2.2 Comparing Spectra with the Intensity Axis of the Spectra Fixed

The intensity axis of all spectra displayed in [Spectrum View] can be fixed to the spectrum having the greatest intensity. By fixing the intensity axis in this way, the spectra of data files acquired by altered conditions can be compared to find the optimum data acquisition settings (CE or product ion).


- 1 In [Chromatogram View], double click at the retention time position of the compound whose spectrum is to be checked.

The spectrum of the specified retention time is displayed.

- 2 Check the compound name in [Chromatogram View].




3 Right-click on [Spectrum View], and select [Fix Intensity Axis] on the menu that is displayed.



The intensity axis of the spectrum is fixed to match the spectrum having the greatest intensity.

4 Compare the spectra to find the optimum CE or precursor ion.



5 Select the compound of the next event using the [Event] spin buttons at [Event].


The chromatogram and spectrum of the data file are displayed.

The intensity axis of [Spectrum View] is fixed to the spectrum having the greatest intensity.

6 Find the optimum CE or precursor ion.

7.2.3 Displaying the Chromatogram of a Specific m/z

1 At [m/z], enter the m/z for which the chromatogram is to be displayed.



The chromatogram of event 1 is displayed.




NOTE

To display the chromatogram of other events, switch to the desired chromatogram using the [Event] spin buttons.

7.2.4 Moving Chromatograms

This section describes how to move chromatograms in parallel with the time axis direction.

1 Select the No. of the data to move by the data selection button.



The following operations are performed on the chromatogram having this data No.


2 Click (Move Left/Right).




NOTE

To move parallel in the intensity axis direction, click the  button.

3 Move the mouse over the chromatogram, and drag from the start point to the end point in the move direction.



The chromatogram moves to that point.




7.2.5 Expanding/Reducing Chromatograms

The chromatogram of selected data can be expanded or reduced as desired.

1 Set the desired chromatogram to a selected state.

2 Click the  (Base Point) button, and click the expanding/reducing base point on the chromatogram to specify that point.


The base point is set by a cross line. This base point will be the origin from where expanding/reducing will be performed, and only that point will be processed as fixed and will not move. If the base point button is already in a clicked state, click the base point button again to deselect that point, and click the button again to set the base point.



3 Click  (Zoom Left/Right) or .

4 Drag the chromatogram to the desired expanding/reducing position.

The chromatogram is expanded or reduced to that point.



NOTE


- Parallel move processing can be used in combination with expanding/reducing processing.
- To cancel a parallel move to return data to its original state, right-click on the chromatogram, and click [Initialize Process] on the displayed menu.

7.2.6 Performing Calculations on Chromatograms and Spectra

Subtraction can be performed on the chromatogram and spectrum of the currently displayed first and second data files.

1 Double-click the icon of the data file to load, and load the two data files to perform calculations on.


2 Click the (Operation Mode) button on the toolbar to switch to the calculation mode.



The comparison chromatogram is displayed under the zoomed chromatogram in [Chromatogram View], and the comparison spectrum is displayed at the lowermost level of [Spectrum View]. Even if three or more data files are loaded, the calculation results of data files 1 and 2 will be displayed at all times.

3 Click the subtraction operation ($1-2$ or $2-1$) to execute.


By $1-2$, the result of subtracting chromatogram 2 from chromatogram 1 is displayed at the lower level. By $2-1$, the result of subtracting chromatogram 1 from chromatogram 2 is displayed at the lower level.



7.3 Printing Graph Images

The currently displayed chromatogram and spectrum can be printed.

- 1 Click **[Print Image]-[Preview]** on the **[File]** menu.
The **[Report]** window opens with the default report format in a set state.
- 2 Check the report in the preview sub-window, and click **[Print]**.



8

Reports

This chapter describes how to use the report function to print chromatograms and quantitative calculation results.

With the report function, you combine report items, such as sample information, chromatograms, and quantitative results, to create reports in various formats.


8.1 Opening the [Report] Window

- 1 Start up the [GCMS Postrun Analysis] program, and click the  (Report Format) icon on the [Postrun] assistant bar.

The [Report] window opens.



The [Report] window can also be opened from the  (Report Format) icon on the [Browser] assistant bar of the [GCMS Browser] program.




No.	Name	Explanation
①	Toolbar (report creation tools)	Various buttons corresponding to the currently displayed window are displayed.
②	Toolbar (report item tools)	Icons corresponding to the items to be placed are displayed on this toolbar. Items, such as data file chromatograms and quantitative result information, can be pasted to reports using the item buttons on the toolbar.

No.	Name	Explanation
③	Layout View	This is the report paper screen. Click an item button on the toolbar, and drag on screen to place the desired item to print.

**NOTE**

The buttons on the toolbar can be turned into a floating palette by either double-clicking an area enclosing the respective buttons or by dragging the buttons to the outside of the toolbar.



Use this feature, for example, to make the toolbar as narrow as possible. To return the buttons to inside the toolbar, either double-click the title bar of the palette or drag the palette to inside the toolbar.

8.2 Creating Report Formats

In [Report] window, either set new report formats or output reports using pre-created templates. Results (e.g. spectrum information) to output must be processed in advance and saved to a data file.

8.2.1 Manually Setting Report Formats


This section describes an example of how to make a new report format in the [Report] window using chromatogram items.

- 1 Click the  (Report Format) icon on the [Postrun] assistant bar or the [Browser] assistant bar.

The [Report] window opens.


- 2 Select [Page Setup] on the [File] menu.

The [Page Setup] sub-window opens.




In this sub-window, set the paper size to match the printer paper source, the orientation (portrait or landscape), and the margins of the report.

- 3** Click  (Chromatogram) on the toolbar, and drag it onto the layout view. Drag from one corner to the diagonally opposite corner.



- 4** Correct the display positions or display items on each tab page of the [GCMS Chromatogram Properties] sub-window, and click [OK].




The sub-window closes, and the newly set items are placed on the layout.


NOTE

- To re-open the [GCMS Chromatogram Properties] sub-window, double-click inside the item frame.
- When making layouts that straddle multiple pages, click the  (Insert) button on the toolbar to insert a page after the currently displayed page.


To display the header and footer, select [Header/Footer] on the [View] menu on the menu bar, and set these in the [Header/Footer] sub-window.



- 5** To check the print content of report format files, click the [Data] tab in the [Data Explorer] sub-window, and drag-and-drop the target data file onto the format window.




The chromatogram information of the data file is displayed in the [Report] window.



NOTE


- To delete an item, right-click on the item to delete, and click [Delete] on the displayed menu.
- When the data file is not displayed in the [Data Explorer] sub-window, click  (Select Project), and specify the folder containing the data file.

- 6 Click the  (Preview) icon on the [Report] assistant bar.
- 7 Check the report in the preview sub-window, and click [Close].
- 8 Click [Save As] on the [File] menu.
The [Save Report Format File As] sub-window opens.
- 9 Select the folder to save the report format file in, enter a file name, and click [Save].
The newly set report format is saved as a report format file.




8.2.2 Using a Template to Set a Report Format

- 1 Click [Option] on the [Tools] menu.
The [Setting Options] sub-window opens.
- 2 Select [Prompt on File New] at [Report Format File] on the [File New] tab page, and click [OK].




3 Click the  (New) button on the toolbar.
The [File New] sub-window opens.

4 Select [Use Template], then select the template file from the list of templates, and click [OK].



The content of the template file is loaded to the [Report] window.



- 5** Edit the report format and save the template file under a new name.

8.3 Printing Reports


8.3.1 Printing Data Processing Results in the [Data Analysis] Window

This section describes how to print reports from data files on the [Data Analysis] window.

- 1** Start up the [GCMS Postrun Analysis] program, and load the data file to the [Data Analysis] window.
The same report is output in either of the qualitative and quantitative sub-windows.
- 2** Click the  (Report) icon on the [Qualitative] assistant bar or the [Quantitative] assistant bar.
The [Data Report] sub-window opens.

3 Click the  (Report) tab in the [Data Explorer] sub-window.


4 Drag-and-drop the target report format file onto the [Data Report] window.



8.3.2 Printing Reports at Batch Processing

This section describes how to print reports during batch processing.

- 1 Open the [Batch Table] window, and load the batch file.
- 2 Select the [Report Output] cell of the row to be printed, and set the name of the report format file to print in the [Report File] cell.



- 3 Click the  (Start) icon on the [Batch] assistant bar. The report is automatically output when processing ends.

A.1 Installing the Software

This chapter describes the procedure for installing the software.

A.1.1 Preparations for Installation

Check the following points before installing the software.

- Installation disk
The installation disk is provided on a DVD-ROM.
When the software installation program is executed, the files are decompressed automatically and copied to the hard disk on the PC.
- The software runs on Windows. Make sure that the OS running on the PC is appropriate.
- Before re-installing GCMSsolution, be sure to uninstall it first.

Reference

See "[A.2 Uninstalling the Software](#)".

- Use LabSolutions to analyze the data measured with GC detectors.

A.1.2 Installing GCMSsolution

1 Turn ON the PC, start up Windows and insert the GCMSsolution installation disk into the DVD-ROM drive.

The [Autorun] window opens.



The [Do you want the following program to allow changes to this computer?] dialog box may be displayed. In this case, click the [Yes] button to proceed.

NOTE

If the above screen does not open automatically, click the [Start] menu on the taskbar. Enter "E:\autorun.exe", and press the [Enter] key. (In this example, "E:" indicates the DVD-ROM drive.)




2 Click [Install] in the [Autorun] window.

After a time the [Welcome to the InstallShield Wizard for GCMSsolution] window opens.




NOTE

- It takes about one minute before the [Welcome to the InstallShield Wizard for GCMSsolution] window is displayed.
- When .NET Framework 4.6.1 or later is not installed, the installation window of .NET Framework opens.




Click [Install].

When the message "The installation of GCMSsolution requires a reboot" is displayed, click [Yes].



After the PC is rebooted, repeat the installation procedure from step 1.

- 3** Click **[Next]** in the **[Welcome to the InstallShield Wizard for GCMSsolution]** window.
The **[User Authentication Mode Select]** window opens.



- 4** Select the user authentication mode in the **[User Authentication]** window.
Select **[User Database in GCMSsolution]** if CLASS-Agent is not to be used. When using CLASS-Agent, first read the following notes and select the required mode.


 **NOTE**

- When the optional CLASS-Agent software is used, user information can be managed using the Shimadzu user authentication tool. Use of the Shimadzu user authentication tool allows user information to be shared with other software (e.g. CLASS-Agent or GCsolution) that utilize the user authentication tool.
- When using the Shimadzu user authentication tool to manage user information, be sure to install Shimadzu user authentication tool Ver. 1.09 or later. If the software is run on a 64-bit version of Windows 7 or Windows 10, be sure to install Shimadzu user authentication tool Ver. 1.10 or later.
- If the Shimadzu user authentication tool is set to be used when the system is already being run using GCMSsolution, the user information will be switched to the information that is being managed by the Shimadzu user authentication tool. For this reason, users who have been using GCMSsolution and who have not been registered to the Shimadzu user authentication tool must be re-registered. Register the required users by referring to "2.4 Registering (Changing/Deleting) Users" in the GCMSsolution Administration Manual.
Also, password and other user information set by the Shimadzu user authentication tool are enabled when a user ID or user name for the same name is set to GCMSsolution and the Shimadzu user authentication tool.


 **Reference**

For details on installation of the Shimadzu user authentication tool and system configuration, refer to the Shimadzu user authentication tool Instruction Manual.


- 5** Click **[Next]** in the **[User Authentication]** window.
The **[Customer Information]** window opens.



- 6** In the **[Customer Information]** window, enter the name of the user at **[User Name]** and the department that the user belongs to at **[Organization]**, and click **[Next]**.
The **[Destination Folder]** window opens. GCMSSolution is now installed to the folder that is displayed as the GCMSSolution installation directory.



- 7** To change the folder to install GCMSSolution to, click **[Change]** in the **[Change Current Destination Folder]** window.
The **[Change Current Destination Folder]** window opens.



Select the folder to install GCMSSolution to.

 **NOTE**


- Unless you have good reason to do otherwise, leave the currently displayed folder as the installation directory.
- When multiple Windows users use GCMSsolution, the following operation may be impossible if a new data folder is created in a location other than the Data folder (or its sub-folders) at the installation destination.
 - Creating sub-folders
 - Saving files
 - File reference by other users

In such case, change the folder security settings using the appropriate Windows function.

Click [OK]. The display returns to the [Destination Folder] window.


8 Click [Next] in the [Destination Folder] window.

The [Ready to Install the Program] window opens.



9 Click [Install] in the [Ready to Install the Program] window.


The [Installing GCMSsolution] window opens, and GCMSsolution is installed.



A

NOTE

To cancel installation midway, click [Cancel].
A window for confirming whether or not to cancel installation opens.



To continue with installation, click [No]. To cancel installation, click [Yes].

10 Install the driver.

NOTE


For Windows 10, the driver installation window does not open. Proceed to step 11.

- 1 The [Would you like to install this device software?] window opens. Click [Install].




The USB driver for MS is installed.

- 2 If the following message is displayed, click [Install].




11 When the display returns to the GCMSSolution installer window, click [Finish].




A.2 Uninstalling the Software

- 1 **Select [Control Panel] on the [Start] menu on the taskbar.**
For Windows 10, press the [X] key with the [Windows] key held down, and click [Control Panel] from the displayed list.


- 2 **Click [Uninstall a Program].**



- 3 **At [Name] select [GCMSsolution], and click [Uninstall] on the right-click menu.**



- 4 **In the [Are you sure you want to uninstall GCMSsolution?] confirmation window, click [Yes].**



NOTE

To re-install GCMSsolution after uninstalling it, first restart the PC, and then install GCMSsolution.

A.3 Connecting to the Instrument

This section describes how to connect the respective MS unit and GC unit to the computer, as well as how to install the MS or GC driver.

A.3.1 Connect the PC to the MS Unit

The MS unit connects to the computer via a USB cable.

Connect the USB cable from a USB port on the computer to the USB port on the back panel of the MS unit. Make sure that the MS is turned ON.



NOTE

Always use a USB cable specified by Shimadzu.

To reduce the risk of the cable becoming disconnected, connect to a USB port on the back of the computer.


A.3.2 Connect the PC to the GC Unit

The GC unit connects to the computer via a USB cable or a Ethernet cable.

■ Connecting via a Ethernet cable


- 1** Connect the Ethernet cable to the connector on the back panel of the GC unit.
- 2** Turn ON the GC power.
The GC unit starts.
- 3** Press  (Home) - [Function] - [Configuration] - [Transmission Settings] on the GC unit.
The [Transmission Settings] window opens.

- 4** Press the [Ethernet] tab, set [IP Address], [Subnet Mask], [Default Gateway], and press [Apply].




- 5** Turn ON the GC power.
The GC unit starts.

- 6** Start up the [GCMS Real Time Analysis] program, and click the [System Configuration] icon on the [Real Time] assistant bar.
The [System Configuration] sub-window opens.



A

- 7** Double-click [Instrument 1] at [Modules Used for Analysis].
The [Instrument] sub-window opens.



The screenshot shows a dialog box titled "Instrument" with the following fields and options:

- Instrument Name: Instrument1
- Instrument Type: GCMS-TQ Series
- Communication: USB
- GC Type: GC-2030
- GC Communication: Ethernet
- IP Address: 0 0 0 0
- Description: (empty text area)

Buttons at the bottom: OK, Cancel, Help.

- 8** Select [GC-2030] at [GC Type].

- 9** Select [Ethernet] at [GC Communication].

- 10** Enter the IP address of GC in [IP Address].

Enter the IP address set in the [Transmission Settings] window on the GC unit in step 4.

- 11** Click [OK].

Close the [Instrument] sub-window.

 **Reference**

For details on the system configurations, see ["2.1 System Configuration"](#).

■ Connecting via a USB cable

- 1 **Connect the USB cable from a USB port on the computer to the USB port on the back panel of the GC unit.**



NOTE

Always use a USB cable specified by Shimadzu.

To reduce the risk of the cable becoming disconnected, connect to a USB port on the back of the computer.


- 2 **Turn ON the GC power.**
The GC unit starts.

- 3 **Press  (Home) - [Function] - [Configuration] - [Transmission Settings] on the GC unit.**
The [Transmission Settings] window opens.

- 4 **Press the [USB] tab and then press [Apply].**

**A**

5 Follow steps 5 to 7 in "■ Connecting via a Ethernet cable"



Instrument

Instrument Name : Instrument1

Instrument Type : GCMS-TQ Series

Communication : USB

GC Type : GC-2030

GC Communication : USB

Serial No : C1234567890

Description :

OK Cancel Help

6 Select [USB] at [GC Communication].

7 Enter the serial number of GC in [Serial Number].

Enter the serial number displayed in the [Transmission Settings] window on the GC unit in step 4.

8 Click [OK].

Close the [Instrument] sub-window.

Reference

For details on the system configurations, see ["2.1 System Configuration"](#).

A.3.3 Installation of the Driver for MS or GC

In Windows 10 or Windows 7, the driver is installed automatically when the MS or GC unit and USB cable are connected. If the driver is not installed automatically, install the driver according to the following procedure.




NOTE

GCMSSolution must be installed on the PC before installing the driver for MS.

1 Connect the MS or GC to the PC using USB cable while Windows is running.


2 From the Windows Start Menu, select [Control Panel]. Click [System and Security].

For Windows 10, press the [X] key with the [Windows] key held down, select [Control Panel] from the displayed list, and click [System and Security].



3 Click [Device Manager].

For Windows 10, click [System] and click [Device Manager] from the list displayed on the left side of the sub-window.



[Device Manager] window opens.


A

4 Double-click the [GCMS-TQ Series] or [GC-2030] icon in [Other devices].




NOTE

Although the explanation below uses the figures for MS driver installation, perform a similar operation for the GC driver as well.




The [GCMS-TQ Series Properties] sub-window opens.

5 Click the [Update Driver] on the [Driver] tab.



The [Update Driver Software] sub-window opens.

6 Click the [Browse my computer for device software].



- 7** At [Search for driver software in this location], set "C:\GCMSsolution\Program" for the MS driver, and "C:\GCMSsolution\Common" for the GC driver, and click [Next]. ("C:" refers to the drive on which GCMSsolution is installed.)




Installation of device driver starts.

If the [Would you like to install this device software?] window opens. Click [Install].



- 8** Click [Close].




A

A.4 Easy sTop

Use of Easy sTop enables septa or glass inserts to be replaced without venting the vacuum system. The time required for instrument re-stabilization is significantly reduced, and re-tuning is usually not necessary.

Easy sTop sets the temperature of the injection port, column oven and interface to 70 °C or less. It may take up to 30 minutes until replacement of the glass inserts or septa is possible.



1 Double-click one of the consumables icons at the instrument monitor.



The [Consumable] tab page in the [Monitor Settings] sub-window opens.

2 Click [Easy sTop].


The [Easy sTop] sub-window is displayed, and the temperature of the injection unit, column oven and interface drops. When each of the temperatures falls to 70 °C or below, the [Push the replace button] status is displayed in the [Easy sTop] sub-window.



A

3 Click [Replace], and then start replacing the or glass insert of the injection unit.

For details on the replacement procedure, refer to the septum replacement procedure and insert replacement procedure in the [MS Navigator] sub-window.




NOTE

Clicking [Replace] stops supply of the carrier gas.

Leaving supply of carrier gas stopped for a long time sometimes lowers the column performance. So, complete replacement as soon as possible.

4 After replacement, click [Complete] in the [Easy sTop] sub-window.

The temperature of the injection unit, column oven and interface returns to the original set point.



5 Reset the number of uses of [Septum] or [Glass Insert].

For details on how to reset the number of uses, refer to step 3 onwards of "A.5 Resetting the Number of Uses/Operation Time".


A.5 Resetting the Number of Uses/Operation Time

GCMSsolution's instrument monitor counts the number of uses and operation time of each of the consumables.

When consumables have been replaced or the **ion source** has been cleaned, for example, reset the numbers of uses and operation time by following the procedure below.


1 Double-click one of the consumables icons at the instrument monitor.

The [Consumable] tab page in the [Monitor Settings] sub-window opens.




2 Click [Reset Consumables].

The [Reset Consumables] sub-window opens.



A

- 3** Either click [Clear] for the consumable that has just undergone maintenance, or enter "0" in the operation time cell and click [OK].



The display returns to the previous sub-window.

- 4** Click [OK] in the [Monitor Settings] sub-window.
The sub-window closes.


A.6 Changing the Replacement Guidelines for the Septa and Glass Inserts

The frequency that the septum must be replaced changes according to the needle thickness of the syringe in use. Recommended syringes can be used about 100 times, and recommended gas-tight syringes can be used about 30 times.

Also, the frequency that glass inserts must be replaced changes according to the sample. Set replacement guidelines according to the sample.

- 1** Click the  (System Configuration) icon on the [Real Time] assistant bar.
The [System Configuration] sub-window opens.


2 Double-click the [SPL1] icon at [Modules Used for Analysis].



The system configuration sub-window of the injection unit opens.

3 Enter [Septum Used Counts] and [Insert Used Counts] settings.

To return settings to their defaults, click [Initialize].




4 Click [OK].


The display returns to the [System Configuration] sub-window.

A

5 Click [Set].



The replacement guidelines for the septa and glass inserts are changed.



A.7 Column Maintenance

The retention time of measurement data sometimes differs from previous data, for example, when the column has been cut to change its length or when a column from a different lot is used. In cases such as this, the MS instrument parameters or retention times in the Compound Table must be adjusted so that the target compound can be reliably measured and identified.

In such cases, use the AART (automatic adjustment of retention time) function to batch-adjust all retention times based on retention indices assigned to target compounds.

1. Perform single run on the reference compound (mainly, n- alkane).
2. Identify the reference compound.
3. Execute the AART function, and adjust the retention times of the method file.
4. Measure the standard sample of the target compound, and adjust the retention time.

Reference

For a more detailed description of operations, refer to the AART Instruction Manual in the \Manual folder of the GCMSsolution installation folder.

A.8 Errors and Remedies

A.8.1 The Instrument Is Not Properly Recognized

When the software is started up, the instrument sometimes cannot be recognized properly and the following error messages are sometimes displayed:

[0d80] MS Communication Hardware is not connected.

[0d00] GC Communication Hardware is not connected.

If this happens, click [Help] and check connections and settings according to the procedure described in the error message that is displayed.

A.8.2 Messages Displayed at Startup

This section describes messages that are displayed when the [Real Time Analysis] program is started up.

[1317] The hardware configuration for this method is different from the current instrument configuration. The measurement condition in the method file is modified according to the current instrument configuration.

This message is displayed when the method file used for the previous data acquisition and the method file set to the Batch Table differ from those for the current instrument configuration. To display the method file or execute realtime batch, the system configuration information currently saved in the method files must be matched to the current instrument configuration.

To match the system configuration information currently saved in the method files to the current instrument configuration, click [OK].

Clicking [Cancel] cancels loading of method files and start of realtime batch.

When the system configuration information has been matched to the current instrument configuration, be sure to check the various parameters since defaults are entered to units in a different configuration.

A.8.3 Windows Does Not Start Up

If an error message is displayed when Windows does not start up normally after the PC is turned ON, check the instruction manual provided with the PC.

A

This page is intentionally left blank.

Appendix **B**

Before PDF Output

The GCMSsolution provides the capability to output a report in the PDF (Portable Document Format) when it is used with Adobe Acrobat, a tool for creation, management, and editing of PDF documents.

B.1 Before PDF Output

To use this function, the following product must be set up beforehand:

Adobe Acrobat (Version 9.2 or later)

For details on how to install this product, perform "Standard Setup" by following the on-screen instructions that are automatically launched when the Acrobat product CD-ROM is inserted into the CD-ROM drive.



NOTE

Adobe Acrobat is not included in this software package. It must be purchased separately.

Adobe Acrobat Reader that is provided on the GCMSsolution installation disk only allows PDF electronic documents to be viewed and does not provide a function for creating PDF files.


B.2 Installing the Postscript Printer Driver

1 Click [Devices and Printers] on the Windows [Start] menu.

For Windows 10, click the Start button on the taskbar, and click [Settings]. Click [Devices] in the displayed sub-window, and then click [Devices and printers] in [Related settings] of the sub-window that is displayed.

2 Click [Add a printer].


The [What type of printer do you want to install?] window opens.



B


3 Click [Add a local printer].

For Windows 10, this windows is not displayed. Click [The printer that I want isn't listed] in the [Add a device] sub-window. Then, select [Add a local printer or network printer with manual settings] in the sub-window that is displayed.




The [Choose a printer port] window opens.

4 Select [Use an existing port] and [FILE: Print to File] for the port.




The [Install the printer driver] window opens.

5 Click [Have Disk].



The [Install From Disk] window opens.


6 Click [Browse].




The [Locate File] window opens.

7 Select "PDF_Setup.inf" from the "\\PDFOutput\English" folder in the GCMS installation disk, and click [Open].


For Windows 10 or Windows 7 (64 bit version), select "PDF_Setup.inf" from the "\\PDFOutput\English\x64" folder in the GCMS installation disk.



8 In the [Install From Disk] window, check that the path to "PDF_Setup.inf" is correctly selected, and click [OK].




9 In the [Install the printer driver] window, select "LabSolutions PDF Printer," and click [Next].




10

In the [Type a printer name] window, check the printer name selected in Step 9 is correctly entered, and click [Next].


If the printer needs to be renamed, rename it at this step.



When the [User Account Control] window opens, enter the password of the administrator account, and click [Yes].




When the [Would you like to install this device software?] window is displayed, click [Install].



This starts installation of the driver.

11


Click **[Finish]**.



12

When installation is completed, load the data file to the **[Data Analysis]** window in the **[GCMS Postrun Analysis]** program, and click **[Save Data File As]** on the **[File]** menu.

The **[Save Data File As]** sub-window opens.



Check whether or not PDF output is set to enabled. If "PDF File (*.pdf)" is selected and **[Save]** is clicked, the data report will be saved in PDF format.


NOTE

If a data file not set with a data report format is converted and saved as a PDF file, a PDF containing a white sheet without any data will be generated.


B.3 Setting PDF Output

This section describes how to set output to PDF file when outputting a report from GCMSsolution.

1 Click [Option] on the [Tools] menu.



2 On the [PDF Export] tab page, select the [Redirect Report to PDF file] checkbox, and select [Normal Speed Mode (Acrobat Distiller)].



Reference

For details about parameters, refer to Help.

3 Click [OK].

This page is intentionally left blank.

With GCMSsolution, chromatograms, spectra, data analysis results, and other information can be exported to other file formats.

Window/Sub-Window	File Type	Convertible File Format (Extender)
[Data Analysis] window	Data files	ASCII files (*.txt)
		AIA files (*.cdf)
[Data Explorer] window	Data files	ASCII files (*.txt)
		AIA files (*.cdf)
		JCAMP files (*.jdx)(Note 1)
		mzData files (*.mzData)
		mzXML files (*.mzXML)
[Library Editor] window	Library files	Text files (*.txt)
		JCAMP files (*.jdx)
[Batch table] window	Batch files	ASCII files (*.txt)

(Note1) GC/MS/MS data files cannot be converted to this format.

C.1 Exporting from the [Data Explorer] Sub-Window

In the [Data Explorer] sub-window, data files can be converted to files of various formats.



NOTE

GC/MS/MS data files can be converted to ASCII format, AIA format, mzData format, and mzXML format.

1

Select the data files (multiple files can be selected) in the [Data Explorer] sub-window, and select the file format on the [File Convert] sub-menu on the right-click menu.


The [To *** File] sub-window is displayed.



NOTE


"""" stands for the file format selected in the sub-menu.

2 Select the folder as the save directory, and click [OK].




NOTE

To cancel file conversion, click [Terminate] in the progress sub-window.




When converting to ASCII format, select the items to export in the [Output Items] sub-window that opens by clicking [Output Items].



C.2 Exporting Data to AIA Files from the [Data Analysis] Window

Data files opened in the [Data Analysis] window can be converted to AIA files.

- 1 Load data to the [Data Analysis] window.
- 2 Click [Save Data File As] on the [File] menu.
- 3 Select [AIA (Andi) File (*.cdf)] at [Save as type] in the [Save Data File As] sub-window.



- 4 Specify the AIA file name at [File name], and click [Save].


C.3 Exporting Data by Batch Processing

Analysis results can be exported to a file each time that processing of one row of the Batch Table ends or when all batch processing ends.

- 1 Open the [Batch Table] window, and load the batch file.
- 2 Click the  (Settings) icon on the [Batch] assistant bar. The [Settings] sub-window is displayed.

C

3 On the [ASCII Convert] tab page, select the [Output ASCII File(s)] checkbox.



NOTE

When converting to AIA format, specify the appropriate settings on the [File Convert] tab page.

4 Set the parameters, and click [OK].

Reference

For details about parameters, refer to Help.

5 Click the (Start) icon on the [Batch] assistant bar.

Batch processing is started, and data is converted to the ASCII file format after it is analyzed.

C.4 Exporting from the [Library Editor] Window

Compound information can be exported from the [Library Editor] window.

Library spectra can be converted to JCAMP format only in the case of private libraries.

1 Start up the [GCMS Postrun Analysis] program, and click the (Library Editor) icon on the [Postrun] assistant bar.

The [Library Editor] window is displayed.

- 2 Click [Open Library] on the [File] menu.
- 3 In the [Open Library File] sub-window, specify the library file, and click [Open].
- 4 Set the criteria in the Library Search Parameters Table.

The screenshot shows the 'Library Editor' window with the following data:

Index	Parameter	Upper/Lower	# of Hit
1	Serial Number	1-50	107886
2	No Setting	No need to set	0
3	No Setting	No need to set	0
4	No Setting	No need to set	0
5	No Setting	No need to set	0
6	No Setting	No need to set	0
7	No Setting	No need to set	0


Hit#	Compd Name	Mol Wt	Formula
1	Hydrogen \$\$ o-Hydrogen \$\$ p-Hydrogen \$\$ M	2	H2
2	Methane \$\$ Marsh gas \$\$ Methyl hydride \$\$	16	CH4
3	Ammonia \$\$ Ammonia gas \$\$ Nitro-SI \$\$ Sprk	17	HN
4	Water \$\$ Distilled water \$\$ Ice \$\$ Water vapo	18	H2O
5	Acetylene \$\$ Ethyne \$\$ Ethane \$\$ Naclylen \$	26	C2H2
6	Hydrogen cyanide \$\$ Hydrocyanic acid \$\$ Ae	27	CHN
7	Diborane \$\$ Diborane(6) \$\$ Boron hydride (B2	28	B2H6
8	Carbon monoxide \$\$ Carbon oxide (CO) \$\$ Ca	28	CO
9	Ethene \$\$ Acetylene \$\$ Bicarbonellted hydrog	28	C2H4
10	Nitrogen \$\$ Nitrogen gas \$\$ Nitrogen-14 \$\$ N	28	N2

Mass Spectrum Plot: 1 : 2 : Hydrogen \$\$ o-Hydrogen \$\$ p-Hydrogen \$\$ Molecular hydrogen \$\$ H2 \$\$ UN 1049 \$\$ UN 1966 \$\$

Y-axis: Inten (x10,000) from 0.00 to 1.00. X-axis: m/z from 1.0 to 2.0. A peak is visible at m/z 2.0.

Summary: CAS#: 1333-74-0 Mol Wt: 2 Serial#: 1
 Compd Name: Hydrogen \$\$ o-Hydrogen \$\$ p-Hydrogen \$\$ Molecular hydrogen \$\$ H2 \$\$ UN 1049 \$\$ UN 1966 \$\$
 Compd Form: H2 Class Flag: No Class Flags.


- 5 Click [Start] on the [Index Search] menu.
The compounds in the library are filtered using the search parameters.



- 6 Click [Export Search Results] on the [Index Search] menu.
The [Export Index Search results] sub-window opens.



7 Specify the text file name at [File name], and click [Save].



The compound information is converted to a text file.



NOTE


To convert spectra to JCAMP format, select [JCAMP File(*.jdx)] at [Save as type].

C.5 Exporting Batch Files from the [Batch Table] Window

The content of Batch Tables can be converted to a text file.


1 Load the batch file in the [Batch Table] window.

2 Click [Save Batch File As] on the [File] menu.



The [Save Batch File As] sub-window opens.

3 Select [Text File (*.txt)] at [Save as type].



4 Specify the text file name at [File name], and click [Save]. The Batch Table is converted to a text file.

This page is intentionally left blank.

Appendix **D**

Specifications

D.1 Workstation

PC	CPU	Intel® Core™ i5-7500 (3.4 GHz) equivalent or greater
	Internal memory	8 GB or more
	Hard Disk	500 GB or more
	OS	Windows 10 Pro 64-bit version
Display	21.5" wide color LCD, number of pixels 1920 x 1080 dpi	
Printer	Laser printer, A4 or larger	

This page is intentionally left blank.



Index

A

Acquisition Mode	36, 75
AIA file	189
ASCII file	189
Assistant Bar	3
Audit Trail	12
Auto Tuning	17, 23
Auto Tuning Condition	24
Autosampler	32
Average Spectrum	89
Averate Subtract Spectrum	91

C

Calculating S/N	120
Calibration Curve	111, 136
Calibration curve	77
Calibration Curve View	111, 128
Chromatogram	74, 78
Chromatogram View	30, 72, 74, 111, 128, 139
CID gas	16
CLASS-Agent	159
COAST	1, 41, 43
Column	14, 178
Comound Table View	128
Compound Table	40, 99, 105
Compound Table Creation Wizard	99
Compound Table View	72, 111
Consumable	50, 175
Continuous data acquisition	55
Creation of Automatic MRM or SIM Table	43

D

Data Acquisition	29
Data Analysis	71
Data Comparison	139
Data Explorer	3, 4
Data file	7
Data File Name	57
Data File Tree	111
Download	47

E

Easy sTop	172
Ecology Mode	51, 53, 66
Errors During Batch Analysis	67
Event	36, 43, 44, 75, 76, 139, 142

F

FASST	43
Fix the Intensity Axis	80, 141

G

Gas Chromatograph	33
GCMS Analysis Editor	2
GCMS Analysis Editor	62
GCMS Browser	2
GCMS Postrun Analysis	71
GCMS Real Time Analysis	29
GCMS Realtime Analysis	2
GCMSPostrun Analysis	2
Glass Insert	176
Group	36

I

Instrument Method View	30
Instrument Parameter	36
Instrument Parameters View	31
Instrument status	50

K

Keyword	5
---------------	---

L

Login	2
-------------	---

M

Manual Identification	118
Manual Peak Ingetrate	117

Manual Peak Integrate.....	86
Mass Spectrometer.....	36
Movement of Detection Point	87
MS detector.....	15
MS Program.....	39
MS Table.....	37, 40

O

Output window.....	3
--------------------	---

P

Pause/Restart	61
PDF Output	181, 187
Printing Graph Image.....	97, 137, 146

Q

Qualitative	72
-------------------	----

R

Reference Data File.....	80
Reference ion information	77

S

Sample information.....	30, 75
Sample Login	46
Select Tuning Mode.....	24
Septum	176
Shimadzu user authentication tool.....	159
Shutdown Method File	51
Similarity Search Parameter.....	94
Single Run	46
Snapshot.....	49
Specification.....	197
Spectrum.....	75
spectrum	78
Spectrum Process Table.....	93, 96
Spectrum View	30, 72, 75, 140
Start.....	47, 60
Stop.....	48, 61
Subtract Spectrum	90
Summary Report	137
System Check	20
System Check Result	21
System Check Results.....	22
System Configuration	11

T

Tuning File	17
Tuning file	7
Tuning Report.....	28

U

Uninstalling the Software	163
---------------------------------	-----

V

Vacuum Level.....	50
-------------------	----