SERION ELISA *classic*

Yersinia IgA/IgG/IgM

Instructions - English

(Version V 138.15)

virion\serion

	<u>Updates</u>
Current version No.:	V 138.15
Previous version:	V 138.14
Update in section:	8.2; 9.1

SERION ELISA *classic* Yersinia IgA/IgG/IgM

CONTENTS

- 1 INTENDED USE
- 2 DIAGNOSTIC RELEVANCE
- 3 TEST PRINCIPLE SERION ELISA *classic*
- 4 KIT COMPONENTS
- 5 MATERIAL REQUIRED BUT NOT SUPPLIED
- 6 STORAGE AND STABILITY
- 7 TEST PROCEDURE SERION ELISA *classic*
 - 7.1 General information
 - 7.2 Sample Preparation and Storage
 - 7.3 Preparation of Kit Reagents
 - 7.4 Overview - Test Procedure
 - 7.5 Manual Test Procedure
 - 7.6 Automated Test Procedure
 - 7.7 SERION ELISA *controls* (external Positive Control / Accuracy Control)
- 8 TEST EVALUATION
 - 8.1 SERION ELISA *classic*
 - 8.2 Automated Evaluation / Software
 - 8.3 Borderline Ranges
 - 8.4 Limits of Quantification
 - 8.5 Qualitative Evaluation with SERION ELISA *classic*
 - 8.6 Criteria of Validity
 - 8.7 Interpretation of Results
 - 8.8 Reference Range of Healthy Individuals
- 9 PERFORMANCE CHARACTERISTICS
 - 9.1 Sensitivity and Specificity
 - 9.2 Reproducibility
 - 9.3 Cross-reactivities
 - 9.4 Interfering substances
- 10 SAFETY MEASURES
 - 10.1 Statements of Warning
 - 10.2 Disposal
 - 10.3 Limitation of the test
- 11 LITERATUR / REFERENCES / RÉFÉRENCES / BIBLIOGRAFIA / BIBLIOGRAFÍA / REFERÊNCIAS / ΑΝΑΦΟΡΕΣ / ΟΔΚΑΖΥ / BIBLIOGRAFIA / ΛΙΤΕΡΑΤΥΡΑ / REFERENCER / REFERENSER / REFERENCIE / REFERENCE / REFERANSER / REFERENCIÁK

SERION ELISA *classic* Yersinia IgA/IgG/IgM

Enzyme-immunoassay for determination of human antibodies

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1 INTENDED USE

The SERION ELISA *classic* Yersinia IgA, IgG and IgM tests are quantitative and qualitative immunoassays for the detection of human antibodies in serum or plasma directed against plasmid encoded virulence markers of *Yersinia*. The SERION ELISA *classic* Yersinia are recommended for the differentiation of acute and recent infections in case of gastrointestinal diseases and reactive arthritis.

2 DIAGNOSTIC RELEVANCE

Y. enterocolitica and *Y. pseudotuberculosis* are important human pathogens which are responsible for a range of intestinal syndromes e.g. enteritis, diarrhoeae, pseudoappendicitis and, in immune compromised patients, haemopathies and sepsis. Post infection complications, such as reactive arthritis, *erythema nodosum* and other rheumatic diseases, may occur. These complications are often apparent in individuals carrying the HLA-B27 gene.

Yersinia are ubiquitous, but primarily distributed in temperate and subtropical areas. They are found in the intestinal tracts and faeces of warm blooded animals (wild and domesticated). Transmission is primarily through contaminated foodstuffs.

The virulence of the bacterium is dependent upon plasmid-encoded proteins with a molecular weight of 23 – 51 kD. These virulence-associated proteins, referred to as YOPs (Yersinia Outer Proteins) have been characterized and are accepted as important diagnostic markers. The Yersinia outer proteins are genus-specific so antibodies against all relevant serotypes are detected.

The dynamics of the various antibody classes are important for the differential diagnosis of rheumatoid disease such as rheumatoid arthritis.

While a *Yersinia* infection results in production of IgG antibodies which may remain elevated for several years, the YOP-specific IgA and IgM responses are more fleeting and in a normal infection fall off within a few months. By contrast, in the case of post-infection complications the IgA titer remains elevated (possibly for years), often at a very high level, while the IgM titer drops off within a few months.

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The test systems are used for the detection of antibodies against all *Yersinia* serovars pathogenic for humans.

3 TEST PRINCIPLE SERION ELISA *classic*

The ELISA (Enzyme Linked Immunosorbent Assay) is an immunoassay, which is particularly suited to the determination of antibodies in the field of infectious serology. The reaction is based on the specific interaction of antibodies with their corresponding antigen. The test strips of the SERION ELISA *classic* microtiter plate are coated with specific antigens of the pathogen of interest. If antibodies in the patient's serum sample are present, they bind to the fixed antigen. A secondary antibody, which has been conjugated with the enzyme alkaline phosphatase, detects and binds to the immune complex. The colourless substrate p-nitrophenylphosphate is then converted into the coloured product p-nitrophenol. The signal intensity of this reaction product is proportional to the concentration of the analyte in the sample and is measured photometrically.

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4 KIT COMPONENTS

Test Components	Pieces / Volume
Break apart microtiter test strips each with eight antigen coated single wells. (altogether 96) [MTP] , 1 frame. The coating material is inactivated.	12 pieces
Standard serum (ready-to-use) [STD] , Human serum in protein containing phosphate buffer; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface antigen) and anti-HCV Ab; preservative: < 0.1 % sodium azide; colouring: Amaranth O	2 x 2 ml
Negative control serum (ready-to-use) [NEG] , Human serum in protein containing phosphate buffer; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface antigen) and anti-HCV Ab; preservative: < 0.1 % sodium azide; colouring: Lissamin Green V	2 ml
Anti-human IgA, IgG or IgM conjugate (ready-to-use) [APC] , Anti-human IgA, IgG or IgM polyclonal antibody, conjugated to alkaline phosphatase, stabilised with protein stabilisation solution; preservative: < 0.1 % methylisothiazolone, < 0.1 % bromnitrodioxane	13 ml
Washing solution concentrate (sufficient for 1000 ml) [WASH] , Sodium chloride solution with Tween 20 and 30 mM Tris/HCl, pH 7.4; preservative: < 0.1 % sodium azide	33.3 ml
Dilution buffer (ready-to-use) [DILB] , Protein containing phosphate buffer with Tween 20; preservative: < 0.1 % sodium azide; colouring: 0.01 g/l Bromphenol blue	2 x 50 ml
Stopping solution (ready-to-use) [STOP] , < 0.1 N sodium hydroxide, 40 mM EDTA	15 ml
Substrate (ready-to-use) [pNPP] , Para-nitrophenylphosphate in solvent free buffer; preservative: < 0.1 % sodium azide	13 ml
Quality control certificate with standard curve and evaluation table [INFO] , (quantification of antibodies in IU/ml or U/ml)	2 pages

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5 MATERIAL REQUIRED BUT NOT SUPPLIED

- Common laboratory equipment
- For the IgM detection: SERION RF Absorbent (Order No. Z200 (20 ml))
- Photometer for microtiter plates with filter, wavelength 405 nm, recommended reference wavelength 620 nm - 690 nm (e.g. 650 nm)
- Microtiter plate washer
- Incubator 37 °C
- Moist chamber
- Distilled water
- Click-Clips (Order Nr. VT120.1)
- Optional: SERION ELISA *control*

6 STORAGE AND STABILITY

Reagent	Storage	Stability
Microtiter strips (coated with antigen)	unopened after opening at 2 – 8 °C in closed aluminum bag with desiccant	see expiry date 6 months
Control sera / Standard sera	unopened after opening at 2 – 8 °C	see expiry date 6 months
Conjugate	unopened after opening at 2 – 8 °C	see expiry date 6 months
Dilution buffer	unopened after opening at 2 – 8 °C	see expiry date 6 months
Washing solution	unopened / after opening at 2 – 8 °C working dilution at 2 – 8 °C working dilution at room temperature	see expiry date 2 weeks 1 week
Substrate	unopened after opening at 2 – 8 °C	see expiry date 6 months
Stopping solution	unopened after opening at 2 – 8 °C	see expiry date 6 months

7 TEST PROCEDURE SERION ELISA *classic*

7.1 General information

Optimum results can only be achieved if the instructions are strictly followed. Only use SERION ELISA *classic* reagents when using SERION ELISA *classic* immunoassays. The components must not be exchanged for reagents of other manufacturers. Standard and control sera of SERION ELISA *classic* immunoassays are defined exclusively for the test kit to be used and must not be used in other lots. Washing solution, substrate and stop solution can be used for all SERION ELISA *classic* immunoassays irrespective of the lot and the test.

Each SERION ELISA *classic* test contains a ready-to-use sample dilution buffer. In some cases the use of special dilution buffers is necessary to guarantee consistent quality and reliable results. The dilution buffers can be used irrespective of the lots.

There are three different conjugate concentrations for each immunoglobulin class (IgA, IgG, IgM), indicated on the label as + (low), ++ (medium) and +++ (high). Conjugates with the same concentration and of the same immunoglobulin class are interchangeable and can be used for other SERION ELISA *classic* immunoassays irrespective of the lot and the test. Dilution or alteration of the reagents may result in a loss of sensitivity. Use aseptic techniques when removing aliquots from the reagent tubes to avoid contamination.

Reproducibility of test results is dependent on thorough mixing of the reagents. Agitate the flasks containing control sera before use and also all samples after dilution (e.g. by using a vortex mixer).

Be sure to pipette carefully and comply with the given incubation times and temperatures. Significant time differences between pipetting the first and last well of the microtiter plate when dispensing samples and control sera, conjugate or substrate can result in different pre-incubation times, which may influence the precision and reproducibility of the results. Avoid exposure of reagents to strong light during storage and incubation.

Adequate washing avoids test unspecificities. Therefore, the washing procedure should be carried out carefully. All of the flat bottom wells should be filled with equal volumes of washing buffer. At the end of the procedure ensure that the wells are free of all washing buffer in order to avoid uncontrolled dilution effects. Avoid foaming!

Reagents must be tightly closed after use to avoid evaporation and contamination. Take care not to mix-up the caps of the bottles and/or vials.

The SERION ELISA *classic* immunoassay is only valid if the lot-specific validation criteria on the quality control certificate are fulfilled.

7.2 Sample Preparation and Storage

Lipaemic, hemolytic or icteric samples (serum or plasma) should only be tested with caution. Obviously contaminated samples should not be tested. Serum or plasma (EDTA, citrate, heparin) collected according to standard laboratory methods are suitable samples. Samples must not be thermally inactivated.

7.2.1 Dilution of Samples

Before running the test, patient samples (V_1) must be diluted in dilution buffer (V_2) as follows:

SERION ELISA *classic* Yersinia IgA

$V_1 + V_2 = 1 + 100$	add	10 μ l	patient's sample
	each to	1000 μ l	dilution buffer

SERION ELISA *classic* Yersinia IgG

$V_1 + V_2 = 1 + 19$	add	10 μ l	patient's sample
	each to	190 μ l	dilution buffer

After dilution and before pipetting into the microtiter plate the samples must be mixed thoroughly to prepare a homogenous solution.

SERION ELISA *classic* Yersinia IgM

Interference with rheumatoid factors

Rheumatoid factors are autoantibodies mainly of the IgM class, which preferably bind to IgG immune complexes. The presence of non-specific IgM antibodies (rheumatoid factors) can lead to false-positive results in the IgM assay. Furthermore, the possibility exists, that weak-binding pathogen-specific IgM antibodies may be displaced by stronger-binding IgG antibodies leading to a false-negative IgM result. Therefore it is necessary to pretreat samples with rheumatoid factor-absorbent prior to IgM detection (SERION RF Absorbent, Order No.: Z200 (20 ml/100 tests)). Rf-absorption is performed by incubation of the patient's sample in Rf-dilution buffer for 15 minutes at room temperature or over night at 4 °C. The test procedure is described in a separate instruction manual.

Before running the test, rheumatoid factor-absorbent (V_1) must be diluted 1+4 in dilution buffer (V_2).

$V_1 + V_2 = V_3 (1 + 4)$	add	200 μ l	Rf-absorbent
	each to	800 μ l	dilution buffer

Patient's samples (V_4) must be diluted in this Rf-dilution buffer (V_3):

$V_4 + V_3 = 1 + 100$	add	10 μ l	patient's sample
	each to	1000 μ l	Rf-dilution buffer

7.2.2 Sample Storage

The patient's samples should not be stored for more than 7 days at 2 – 8 °C. Extended storage is possible at ≤ -20 °C. Avoid repeated freezing and thawing of samples. Diluted samples can be stored at 2 – 8 °C for one week.

7.3 Preparation of Kit Reagents

Bring all reagents to room temperature before testing.

7.3.1 Microtiter Test Strips

The microtiter test strips labeled with abbreviations for pathogen and immunoglobulin class are packed with a desiccant in an aluminum bag. To open the aluminum bag of the microtiter plate please cut off the top of the marked side only, in order to guarantee proper resealing. Take unrequired cavities out of the frame and put them back into the aluminum bag. Close bag carefully to ensure airtight conditions. Do not use the strips if the aluminum bag is damaged or if the bag with remaining strips and desiccant was not properly resealed.

7.3.2 Negative control Sera / Standard Sera (ready-to-use)

Negative control and standard sera are ready-to-use and must not be diluted any further. For each test run - independent of the number of microtiter test strips to be used – negative control and standard sera must be included. Standard sera should be set up in duplicate. Do not treat negative control and standard sera with Rf-absorbent.

7.3.3 Anti-human IgA, IgG or IgM AP-Conjugate (ready-to-use)

The required conjugate concentration (+, ++, +++) is indicated on the quality control certificate. Please refer also to the specification on the label. Avoid contamination.

7.3.4 Washing Solution (Concentrate)

Dilute washing buffer concentrate (V_1) 1:30 with aqua dest. to a final volume of V_2 . Bottles used for the working dilution should be cleaned regularly. Discard cloudy solutions.

Example:

Buffer concentrate (V_1)	Final volume (V_2)
33.3 ml	1000 ml
1.0 ml	30 ml

7.3.5 Dilution Buffer for Samples (ready-to-use)

Discard cloudy solutions.

7.3.6 Substrate (ready-to-use)

Substrate in unopened bottle may have a slightly yellow coloring, which does not reduce the quality of the product! Avoid contamination.

7.3.7 Stopping Solution (ready-to-use)


7.4 Overview - Test Procedure

SERION ELISA *classic* Yersinia IgA/IgG/IgM quantitative

In case of IgM detection absorption of rheumatoid factor, see No. 7.2.1;
Incubation 15 minutes at room temperature or over night at 4°C

sample dilution¹
(patient's samples)
IgG: 1 + 19
IgA/IgM: 1 + 100

Pipette diluted samples and ready-to-use negative control /
standard sera into the microtest wells (100 µl)

¹Special dilution buffers for the following SERION ELISA *classic* tests:
Borrelia burgdorferi IgG, IgM, Dengue Virus superior IgM and Epstein-Barr Virus EA IgG.

²For manual use:
tap plate at the end of the wash procedure on paper towel.

7.5 Manual Test Procedure

1. Place the required number of **cavities in the frame** and prepare a protocol sheet.
2. Add each **100 µl of diluted sample or ready-to-use negative control/standard sera** into the appropriate wells of microtiter test strips. Spare one well for substrate blank, e.g.:

Well	Quantitative ELISA
A1	substrate blank
B1	negative control
C1	standard serum
D1	standard serum
E1	patient 1 ...
F1	patient 2 ...

3. **Sample incubation** for 60 minutes (+/- 5 min.) at 37 °C (+/- 1°C) in moist chamber
4. After incubation **wash** all wells with washing solution (by automated washer or manually):
 - aspirate or shake out the incubation solution
 - fill each well with 300 µl washing solution
 - aspirate or shake out the washing buffer
 - repeat the washing procedure 3 times (altogether 4 times!)
 - dry by tapping the microtiter plate on a paper towel
5. **Addition of conjugate**
Add 100 µl of the ready-to-use IgA/IgG/IgM conjugate to the appropriate wells (except substrate blank)
6. **Conjugate incubation** for 30 minutes (+/- 1 min.) at 37 °C (+/- 1 °C) in moist chamber.
7. After incubation **wash** all wells with washing solution (see above).
8. **Addition of substrate**
Add 100 µl of ready-to-use substrate solution to each well (including well for substrate blank!)
9. **Substrate incubation** for 30 minutes (+/- 1 min.) at 37 °C (+/- 1 °C) in moist chamber. Ensure dark incubation.
10. **Stopping of the reaction**
Add 100 µl stopping solution to each well, shake microtiter plate gently to mix.
11. **Read extinction**
Read optical density (OD) within 60 minutes at 405 nm against substrate blank, reference wave length between 620 nm and 690 nm (e.g. 650 nm).

7.6 Automated Test Procedure

SERION ELISA are validated for use with Immunomat (using the following consumables: VT124, VT111, VT112) and suited for processing on similar analyzers. For processing on the Immunomat the current software version including reagent check has to be used. The automated processing is performed analogous to manual use. Please note, that under special working-conditions internal laboratory adaptations of the substrate incubation times may be necessary.

7.7 SERION ELISA controls (external Positive Control / Accuracy Control)

For the periodic verification of the test method, in order to fulfil the requirements of laboratory internal quality management systems, we recommend using SERION ELISA *controls* to determine precision and accuracy of SERION ELISA *classic* test runs. SERION ELISA *controls* are separately available and the usage is described in specific instruction manuals. SERION ELISA *controls* are not available in all countries and the customer should consult the local distributor.

8 TEST EVALUATION

8.1 SERION ELISA *classic* Yersinia IgA/IgG/IgM

The mathematical curve fitting for antibody quantification with SERION ELISA *classic* immunoassays is based on the 4-parameter logistic (4 PL) function.

$$Activity (U / ml) = e^{C - \frac{1}{B} \ln\left(\frac{D-A}{OD(Patient)*F-A}\right)}$$

The 4 parameters A, B, C, and D are representative for the exact shape of the standard curve:

Parameter A:	Lower asymptote (OD)
Parameter B:	Slope of the curve
Parameter C:	Inflection point
Parameter D:	Upper asymptote (OD)

Institut Virion\Serion GmbH establishes a lot-specific 4 PL standard curve for each SERION ELISA *classic* immunoassay in multiple test runs under optimal test conditions. The four parameters are indicated on the quality control certificate of each individual SERION ELISA *classic* test.

For the adaption of the test level to the given 4 PL standard curve the correction factor F is calculated by dividing the standard reference OD value indicated on the quality control certificate with the measured, and consequently test run-specific, standard OD value.

$$F = \frac{STD \text{ reference } OD \text{ value}}{\text{measured } STD \text{ OD value}}$$

By multiplying the OD values obtained from patient samples with the correction factor F, the level of each individual test run is adjusted to the given 4 PL standard curve. Thereby,

interassay deviations are compensated for and antibody activities can be directly evaluated from the 4 PL standard curve.

After subtraction of the substrate blank from all measured OD values and calculation of the mean OD value of the standard serum (STD), tested in duplicate, the evaluation of antibody activities from the optical measurement signals (OD) of patient samples can be performed with the 4PL function presented above.

8.2 Automated Evaluation / Software

Institut Virion\Serion GmbH recommends the use of the SERION *easyANALYZE* software for the automated evaluation of optical measurement signals.

8.3 Borderline Ranges

The borderline ranges of the SERION ELISA *classic* Yersinia IgA/IgG/IgM tests are specified on the quality control certificates and indicate the range of borderline test results. Values below this range indicate a negative test result; values above the borderline range are interpreted positive.

8.4 Limits of Quantification

The limits of quantification are specified on the quality control certificate of the SERION ELISA *classic* Yersinia IgA/IgG/IgM. The linearity of dilution within this range has been demonstrated in comprehensive evaluation studies. In case a patient sample shows a test result above the upper limit of quantification, the sample may be tested at a higher dilution. The resulting antibody activity must then be multiplied by the additional dilution factor.

8.5 Qualitative Evaluation with SERION ELISA *classic* Yersinia IgA/IgG/IgM

For the SERION ELISA *classic* test evaluation a lot-specific quality control certificate with standard curve and an evaluation table is included in the test kit so that the obtained OD values may be assigned to the corresponding antibody activities. The substrate blank must be subtracted from all OD values prior to evaluation. Mean OD value of the standard serum (STD), tested in duplicate, has to be used.

Method 1:

In the first line of the table, several ranges of OD values for the standard serum are depicted covering the whole standard validity range. According to the measured mean OD value of the standard serum, the corresponding column can be chosen. This column contains the information of upper and lower cut-off OD values to allow evaluation of the patient sample. OD values below the lower cut-off are evaluated negative and values above the upper cut-off are evaluated positive. Implementation of the correction factor F is not necessary in the context of the evaluation table.

Method 2:

Qualitative Evaluation

To fix the cut-off ranges multiply the mean value of the measured standard OD with the numerical data of the quality control certificate (see special case formulas), e.g.:

OD = 0.502 x MW(STD) with upper cut-off

OD = 0.352 x MW(STD) with lower cut-off

If the measured mean absorbance value of the standard serum is 0.64 OD, the range of the cut-off is in between 0.225-0.321 OD.

Calculation example:

Standard serum mean OD = 0.64

Upper cut-off: OD = 0.502 x 0.64 = 0.321

Lower cut-off: OD = 0.352 x 0.64 = 0.225

8.6 Criteria of Validity

The substrate blank must be < 0.25 OD.

The negative control must be negative.

The mean OD value (after subtraction of the substrate blank!) of the standard serum must be within the validity range, which is given on the lot specific quality control certificate.

The variation of OD values of the standard serum must not be higher than 20 %.

If these criteria are not met, the test is not valid and must be repeated.

8.7 Interpretation of Results

A positive test result confirms the presence of specific antibodies. A negative result indicates that no clinically relevant antibodies against the pathogen are present in the patient's sample, but does not exclude the possibility of an acute infection. In case of a borderline result a reliable evaluation is not possible. A definitive diagnosis can only be achieved by testing paired serum samples, taken at one to two weeks intervals, in parallel.

IgG, IgA and IgM antibodies are detectable following infection with *Yersinia* and IgG antibodies may persist for years. The published literature states that the seroprevalence in the normal population in Europe is approx. 30 to 45 %. In an uncomplicated acute, primary infection levels of IgA and IgM antibodies fall within a few weeks or months to undetectable levels. Secondary immunological complications such as arthritis and *erythema nodosum* are characterized by a persistent elevated IgA titer while IgM falls as usual.

8.8 Reference Range of Healthy Individuals

Testing of random blood donor sera, collected in the region of southern Germany, with SERION ELISA *classic* Yersinia resulted in the following distribution:

SERION ELISA <i>classic</i>	Number	negative	borderline	positive
Yersinia IgA	105	91 (86.7 %)	7 (6.7 %)	7 (6.7 %)
Yersinia IgG	105	63 (60.0 %)	9 (8.6 %)	33 (31.4 %)
Yersinia IgM	105	103 (98.1 %)	2 (1.9 %)	-

9 PERFORMANCE CHARACTERISTICS

9.1 Sensitivity and Specificity

The performance characteristics of the SERION ELISA *classic* Yersinia IgA and IgG tests were evaluated by the analysis of more than 150 serum samples from patients with suspected yersiniosis, healthy blood donors and samples from interlaboratory tests. The performance characteristics of the SERION ELISA *classic* IgM was evaluated using 22 serum samples from patients with suspected infection. An ELISA and an Immunoblot of a leading European manufacturer was used as reference test. The immunoassays of leading European manufacturers were used as reference tests. Sera classified as borderline were not included in the calculation of sensitivity and specificity.

	Sensitivity	Specificity
SERION ELISA <i>classic</i> Yersinia IgA	94.8 %	89.7 %
SERION ELISA <i>classic</i> Yersinia IgG	97.7 %	97.8 %
SERION ELISA <i>classic</i> Yersinia IgM	> 99 %	90.0 %

9.2 Reproducibility

SERION ELISA *classic* Yersinia IgA

Sample	Mean Value (OD)	Intraassay CV (%)	Mean Value (OD)	Interassay CV (%)
Serum 1	0.400	3.1	0.491	4.1
Serum 2	0.866	2.7	0.990	3.0
Serum 3	2.289	3.0	2.697	2.9

SERION ELISA classic Yersinia IgG

Sample	Mean Value (OD)	Intraassay CV (%)	Mean Value (OD)	Interassay CV (%)
Serum 1	0.755	3.9	0.634	7.1
Serum 2	1.030	4.1	0.940	6.9
Serum 3	1.518	3.7	1.480	5.0

SERION ELISA classic Yersinia IgM

Sample	Mean Value (OD)	Intraassay CV (%)	Mean Value (OD)	Interassay CV (%)
Serum 1	0.272	2.9	0.296	2.8
Serum 2	0.626	1.2	0.605	2.8
Serum 3	1.045	2.6	1.028	2.6

9.3 Cross-reactivities

SERION ELISA classic Yersinia IgA

To determine detection of cross-reactive antibodies directed against different parameters sera were analyzed with SERION ELISA *classic* Yersinia IgA and a commercially available anti-Yersinia IgA ELISA. Positive sera (10 sera each) for *Helicobacter pylori* IgA and *Chlamydia pneumoniae* IgA have been tested as well as sera positive for rheumatoid factor (RF) and anti-nuclear antibodies (ANA). Within this internal evaluation potential cross-reactivities with six RF positive serum samples have been observed. The majority of reactivities have been confirmed by positive or borderline results in the reference assay. Other cross-reactivities cannot be ruled out in general.

SERION ELISA classic Yersinia IgG

To determine detection of cross-reactive antibodies directed against different parameters sera were analyzed with SERION ELISA *classic* Yersinia IgG and a commercially available anti-Yersinia IgG ELISA. Positive sera (10 sera each) for *Helicobacter pylori* IgG, *Francisella tularensis* IgG, *Borrelia burgdorferi* IgG and *Chlamydia pneumoniae* IgG have been tested as well as sera positive for rheumatoid factor (RF) and anti-nuclear antibodies (ANA). Within this internal evaluation potential cross-reactivities with one *Helicobacter pylori* IgG, one *Chlamydia pneumoniae* IgG, one ANA and seven RF positive serum samples have been observed. All reactivities have been confirmed by positive or borderline results in the reference assay..Other cross-reactivities cannot be ruled out in general.

SERION ELISA *classic* Yersinia IgM

To determine detection of cross-reactive antibodies directed against different parameters sera were analyzed with SERION ELISA *classic* Yersinia IgM and a commercially available anti-Yersinia IgM ELISA. Positive sera (10 sera each) for Helicobacter pylori IgM, Francisella tularensis IgM, Borrelia burgdorferi IgM and Chlamydia pneumoniae IgM have been tested as well as sera positive for rheumatoid factor (RF) and anti-nuclear antibodies (ANA). Within this internal evaluation potential cross-reactivities with one Francisella tularensis IgM and one Chlamydia pneumoniae IgM positive serum samples have been observed. Other cross-reactivities cannot be ruled out in general.

9.4 Interfering substances

SERION ELISA *classic* Yersinia IgA/IgG/IgM

To determine the influence of interfering substances, sera with different reactivities were analyzed with SERION ELISA *classic* Yersinia IgA/IgG/IgM. No interferences have been detected for sera with concentrations up to 2.00 g/L hemoglobin, 11.50 g/L lipemia/triglyceride or 0.201 g/L bilirubin (conjugated and unconjugated).

10 SAFETY MEASURES

10.1 Statements of Warning

The SERION ELISA *classic* is designed for use by qualified personnel who are familiar with good laboratory practice. All kit reagents and human specimens should be handled carefully, using established good laboratory practice.

- This kit contains human blood components. Although all control- and cut-off sera have been tested and found negative for anti-HIV-ab, HBs-Ag (*Hepatitis B-Virus-surface Antigen*) and anti-HCV-ab, they should be considered potentially infectious.
- Do not pipette by mouth.
- Do not smoke, eat or drink in areas in which specimens or kit reagents are handled.
- Wear disposable gloves, laboratory coat and safety glasses while handling kit reagents or specimens. Wash hands thoroughly afterwards.
- Patient's material and other potentially infectious material should be decontaminated after the test run.
- Reagents should be stored safely and be inaccessible to unauthorized access e.g. children.

10.2 Disposal

Please observe the relevant statutory requirements!

10.3 Limitation of the test

Please note that diagnosis should never be solely based on serological data. Rather, serological results have to be interpreted in the context of the clinical picture and other diagnostic findings.

11 LITERATUR / REFERENCES / RÉFÉRENCES / BIBLIOGRAFIA / BIBLIOGRAFÍA / REFERÊNCIAS / ΑΝΑΦΟΡΕΣ / ΟΔΚΑΖΥ / BIBLIOGRAFIA / ЛИТЕРАТУРА / REFERENCER / REFERENSER / REFERENCIE / REFERENCE / REFERANSER / REFERENCIÁK

- [1] Baert, F. (1994) Yersiniosis: The Clinical Spectrum. *Acta Clinica Belgica* 49, 76-85.
- [2] Bottone, E.J. (1997) *Yersinia enterocolitica*: The Charisma Continues. *Clin. Microbiol. Rev.* 10, 257-76.
- [3] Cremer, J. *et al.* (1993) Immunoblotting of *Yersinia* plasmid-encoded released proteins: A tool for serodiagnosis. *Electrophoresis* 14, 952-9.
- [4] Cremer, J. *et al.* (1995) Plasmidkodierte sekretorische Proteine von *Yersinia enterocolitica* als Antigene in einem IgG- und IgA-spezifischen immundiagnostischen Elisa. *Klin. Lab.* 41, 55-9.
- [5] Gransfors, K. (1980) Persistence of IgM, IgG, and IgA Antibodies to *Yersinia* in *Yersinia* Arthritis. *J. Infect. Dis.* 141, 424-9.
- [6] Gransfors, K. *et al.* (1988) IgM, IgG, and IgA Antibodies in *Yersinia* Infection. *J. Infect. Dis.* 157, 601-2.
- [7] Hammer, M. *et al.* (1995) Postenteritische reaktive Arthritiden und Spondarthritiden. *Deutsches Ärzteblatt* 92, 2738-49.
- [8] Lahesmaa-Rantala, R. *et al.* (1989) Avidity of antibodies against released proteins of *Yersinia* spp: comparison of patients with or without reactive arthritis. *Ann. Rheum. Dis.* 48, 1003-6.
- [9] Mäki-Ikola, O. *et al.* (1994) *Yersinia*-specific antibodies in serum and synovial fluid in patients with *Yersinia* triggered reactive arthritis: *Ann. Rheum. Dis.* 53; 535-9.
- [10] Mäki-Ikola, O. (1997) High frequency of *Yersinia* antibodies in healthy populations in Finland and Germany. *Rheumatol. Int.* 16, 227-9.
- [11] Ruckdeschel, K., Deuretzbacher, A., and Haase, R. (2008) Crosstalk of signalling processes of innate immunity with *Yersinia* Yop effector functions. *Immunobiol.* 213, 261-9.
- [12] Skurnik, M. (2007) My Life with *Yersinia*. *Adv. Exp. Med. Biol.* 603, 44-73.
- [13] Tomaso, H., Mooseder, G., Al Dahouk, S., Bartling, C., Scholz, H.C., Strauss, R. Treu, T.M. and Neubauer, H. (2006) Seroprevalence of Anti-*Yersinia* Antibodies in Healthy Austrians. *Eur. J. Epidemiol.* 21, 77-81.

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

Dear Sir or Madam,

IVDR certificate

In February 2025 we have informed you about the achievement of the IVDR certification for our class B and C products of the following product groups:

- SERION ELISA *classic* (ESRxxx)
- SERION ELISA *antigen* (ESR200)
- SERION ELISA *control* (BCxxx)
- SERION *immunoCONTROL* (ICxxx)
- SERION ELISA AI *control* (CLSxxx)
- SERION ELISA Avidity Control (BRxxxAVID)
- SERION ELISA Avidity Reagent (BxxxAVID)
- SERION RF Absorbent (Z200)

It is certified that requirements according to REGULATION (EU) 2017/746 on In Vitro Diagnostic Medical Devices Annex IX Chapter I, Section 2 and 3 and Chapter III are fulfilled for the covered products of class B as well as class C products.

With this mailing we would like to inform you that the current IVDR certificates are being transferred to another Notified Body, TÜV Süd (number 0123).

The IVDR certificate will be issued by TÜV Süd with the corresponding number 0123. Please note that the content of the certificate remains unchanged and therefore the certificate transfer does not affect the certification status of our products. A change will apply for the labelling of the devices as the number of the new Notified Body will be implemented for labels, Instructions for Use (IFU) and Declarations of Conformity (DoC).

virion\serion

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19, 97076 Würzburg, Germany

Phone +49 931 3045 0 Fax +49 931 3045 100 Mail diagnostics@virion-serion.de Web www.virion-serion.de

The effective date for certificate transfer is **June 30, 2025**.

The certificate is shared on our homepage here: [Quality and Certificates | Virion\Serion \(virion-serion.de\)](https://www.virion-serion.de/Quality-and-Certificates-Virion-Serion)

From July 1, 2025, newly produced product lots will be labelled with the new number 0123.

Updated DoCs, IFUs and labels will be made available upon request after July 1, 2025. Please let us know which information is required from your side.

The change of labelling will be implemented successively for the individual products and depends on the production schedule. Within a transition period of 18 months all products will carry the new labelling.

The CE with number of Notified Body will be removed from component labels of device components as this is not mandatory for components.

EN ISO 13485 certificate

The EN ISO 13485 certificate has been transferred to TÜV Süd (number 0123) and validity date is from 2025-05-30. The certificate number is Q5 131304 0001 Rev. 00. The scope has been extended regarding the instruments for infection diagnostics.

The certificate is shared on our homepage here: [Quality and Certificates | Virion\Serion \(virion-serion.de\)](https://www.virion-serion.de/Quality-and-Certificates-Virion-Serion)

If you have any further questions or comments, don't hesitate to contact us.

We thank you for your understanding and cooperation.

With best regards,

Institut Virion\Serion GmbH

EU Quality Management System Certificate

Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices, Annex IX Chapter I

Certificate No. V13 131304 0003 Rev. 00

Manufacturer: **Institut Virion\Serion GmbH**
Friedrich-Bergius-Ring 19
97076 Würzburg
GERMANY

SRN Manufacturer - DE-MF-000007582

The quality management system has been evaluated in accordance with Regulation (EU) 2017/746, Annex IX Chapter I with a positive result.

Details on devices covered by the quality management system are described on the following page(s). The report referenced below summarises the results of the assessment and includes reference to relevant CS, harmonised standards and test reports.

The certified quality management system is subject to periodical surveillance.

If class A devices in sterile conditions are covered by this certificate, the audit was limited to the aspects relating to establishing, securing, and maintaining sterile conditions.

If class B or C excluding self-/near-patient-testing, or class C companion diagnostics devices are covered by this certificate, the quality management system assessment was accompanied by the assessment of technical documentation for devices selected on a representative basis. The periodical surveillance includes further assessment of the technical documentation on the basis of representative samples.

If class D devices, class B or C self-/near-patient testing, or class C companion diagnostics devices are covered by this certificate, an EU Technical Documentation Assessment Certificate in accordance with Annex IX Chapter II is required before placing them on the market.

All applicable requirements of the Testing, Certification, Validation and Verification Regulations TÜV SÜD Group have to be complied with.

For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:V13 131304 0003 Rev. 00

Report No.: 713371585
Valid from: 2025-07-01
Valid until: 2029-05-15

Marta Carnielli
Head of Certification IVD

Issue date: 2025-07-01

EU Quality Management System Certificate

Regulation (EU) 2017/746 on in Vitro Diagnostic Medical Devices, Annex IX Chapter I

Certificate No. V13 131304 0003 Rev. 00

Classification: Class B
Device Group: IVR 0503 - Infectious agent detection: Presence of, or exposure to an infectious agent
Intended Purpose: Devices intended to be used to detect the presence of, or exposure to an infectious agent

Classification: Class C
Device Group: W0105 + IVP 3007 - Infectious diseases
Intended Purpose: IVD Reagents for Infectious diseases

The validity of this certificate depends on conditions and/or is limited to the following: -none-

Rev.	Dated	Report	Description
00	2025-07-01	713371585	Initial issuance

Yersinia IgA Serum 1:100

YerA5.0.2

No password protection.

Written by Virion\Serion GmbH.

Date created: 2007 m. geguže 15 d., antradienis 10:26:02.

Last edited: 2024 m. liepa 8 d., pirmadienis 10:50:11.

Plate ID Check: *

Generic 96 well plate type.

Maximum tolerable ambient temperature is 28,0°C.

Combination group: 0

	1	2	3	4	5	6	7	8	9	10	11	12
A	SLW1	T5	T13	T21	T29	T37	T45	T53	T61	T69	T77	T85
B	NC1	T6	T14	T22	T30	T38	T46	T54	T62	T70	T78	T86
C	STD1	T7	T15	T23	T31	T39	T47	T55	T63	T71	T79	T87
D	STD1	T8	T16	T24	T32	T40	T48	T56	T64	T72	T80	T88
E	T1	T9	T17	T25	T33	T41	T49	T57	T65	T73	T81	T89
F	T2	T10	T18	T26	T34	T42	T50	T58	T66	T74	T82	T90
G	T3	T11	T19	T27	T35	T43	T51	T59	T67	T75	T83	T91
H	T4	T12	T20	T28	T36	T44	T52	T60	T68	T76	T84	T92

Pipette

Predilution

Pipette into wells E1-H12.

Action on error: Raise alarm and stop

Aspirate 210ul from DILB (speed 2).

Aspirate 10ul from patient sample (speed 1).

Dispense 210ul into dilution plate ("Dilution 1:101") (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Dilb

Pipette into wells E1-H12.

Action on error: Raise alarm and stop

Aspirate 800ul from DILB (speed 2).

Dispense 800ul into dilution plate ("Dilution 1:101") (speed 5).

Incubate

Incubate for 1 (+15,-0) minutes at room temperature.

Incubate in ambient light.

--- Scheduling break ---

Pipette

NC

Pipette into wells B1.

Action on error: Raise alarm and stop

Aspirate 110ul from Yersinia A NC (speed 3).

Dispense 100ul into microplate (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Std

Pipette into wells C1-D1.

Action on error: Raise alarm and stop

Aspirate 110ul from Yersinia A Std (speed 3).

Dispense 100ul into microplate (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Sample 1:101

Pipette into wells E1-H12.

Action on error: Raise alarm and stop

Aspirate 110ul from dilution plate ("Dilution 1:101") (speed 3).

Dispense 100ul into microplate (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Incubate

Incubate for 60 (+5,-5) minutes at 37,0+/-1,0°C.

Wash

Perform 4 x 300ul wash cycles using Wash buffer on strips 1-12.

Dispense rate 75.

With aspirate sweep.

Partial plate mode: as quick as possible.

Perform 2 aspirate cycle(s) on strips 1-12.
With aspirate sweep.

Partial plate mode: as quick as possible.

Dispense

Add 8x100ul (+50ul oversoak) of APC (Conjugate) IgA (3/3) into wells B1-H12.

Incubate

Incubate for 30 (+1,-1) minutes at 37,0+/-1,0°C.

Wash

Perform 4 x 300ul wash cycles using Wash buffer on strips 1-12.

Dispense rate 75.

With aspirate sweep.

Partial plate mode: as quick as possible.

Perform 2 aspirate cycle(s) on strips 1-12.

With aspirate sweep.

Partial plate mode: as quick as possible.

Dispense

Add 8x100ul (+50ul oversoak) of pNPP (Substrate) (2/3) into wells A1-H12.

Incubate

Incubate for 25 (+1,-1) minutes at 37,0+/-1,0°C.

Dispense

Add 8x100ul (+50ul oversoak) of STOP (2/3) into wells A1-H12.

Read

Reader settings

Shake for 5 seconds.

Using wavelengths: 405nm,620nm.

Average inputs. Use arithmetic mean.

Report results to 3 decimal places.

Auto scientific mode for:

Average blank mode.

Blank calculation on final results only.

Validation criteria

SLW<0.25

0.5*RefValue<=Std<=1.7*RefValue

0.8<=STD1.1/Std1.2<=1.25

NC<lcutoff*Std

Full V.C. report.

Do not display results if V.C. fails.

Qualitative Settings

If 'Sample<lcutoff*Std' Then Result:='neg'

If 'Sample>ucutoff*Std' Then Result:='pos'

Default result := ?

- pos Indicates a positive result
- ? Indicates an equivocal result
- neg Indicates a negative result

No histogram required

Average inputs. Use arithmetic mean.

Report results to 3 decimal places.

Auto scientific mode for:

Quantitative

STD1 = 0 U/mL

Average standard replicates.

Virion data model.

Show predicted standard values.

Display a concentration graph.

X axis settings: label: 'auto.', minimum=5, maximum=100, units=auto., show gridlines, logarithmic

Y axis settings: label: 'auto.', minimum=0, maximum=3, units=1, show gridlines

Dilutions:

T1 = 1	T2 = 1	T3 = 1	T4 = 1
T5 = 1	T6 = 1	T7 = 1	T8 = 1

T9 = 1	T10 = 1	T11 = 1	T12 = 1
T13 = 1	T14 = 1	T15 = 1	T16 = 1
T17 = 1	T18 = 1	T19 = 1	T20 = 1
T21 = 1	T22 = 1	T23 = 1	T24 = 1
T25 = 1	T26 = 1	T27 = 1	T28 = 1
T29 = 1	T30 = 1	T31 = 1	T32 = 1
T33 = 1	T34 = 1	T35 = 1	T36 = 1
T37 = 1	T38 = 1	T39 = 1	T40 = 1
T41 = 1	T42 = 1	T43 = 1	T44 = 1
T45 = 1	T46 = 1	T47 = 1	T48 = 1
T49 = 1	T50 = 1	T51 = 1	T52 = 1
T53 = 1	T54 = 1	T55 = 1	T56 = 1
T57 = 1	T58 = 1	T59 = 1	T60 = 1
T61 = 1	T62 = 1	T63 = 1	T64 = 1
T65 = 1	T66 = 1	T67 = 1	T68 = 1
T69 = 1	T70 = 1	T71 = 1	T72 = 1
T73 = 1	T74 = 1	T75 = 1	T76 = 1
T77 = 1	T78 = 1	T79 = 1	T80 = 1
T81 = 1	T82 = 1	T83 = 1	T84 = 1
T85 = 1	T86 = 1	T87 = 1	T88 = 1
T89 = 1	T90 = 1	T91 = 1	T92 = 1
NC1 = 1			

Result annotation:
 IF(x(y)<5,Result:='<5',,)
 IF(x(y)>100,Result:='>100',,)
 IF(y<=A*(Std/RefValue),Result:='<5',,)
 IF(y>=D*(Std/RefValue),Result:='>100',,)
 IF(y<0,0,Result:='<5',,)
 IF(y>3,0,Result:='>100',,)
 Average inputs. Use arithmetic mean.
 Report results to 1 decimal places.
 Auto scientific mode for:
 Header:

Report Settings

Header:	&I	&D	&Q
Body:		Laboratory Details Assay Header Lot Specific Values Validation criteria Quantitative Results 1 Combined Report	
Footer:	&V &F		&S &P

Combined Report
 Table of results :- Patient ID, Well Location, Reader value, Quant. 1 value, Qual. value, Flag
 Table order: STD, NC, EC, T, PC, CO

ASTM E 1394 Quant. 1 mean

End of protocol

Yersinia IgG Serum 1:20

YerG5.0.2

No password protection.

Written by Virion\Serion GmbH.

Date created: 2007 m. geguže 15 d., antradienis 10:26:02.

Last edited: 2024 m. liepa 8 d., pirmadienis 10:52:02.

Plate ID Check: *

Generic 96 well plate type.

Maximum tolerable ambient temperature is 28,0°C.

Combination group: 0

	1	2	3	4	5	6	7	8	9	10	11	12
A	SLW1	T5	T13	T21	T29	T37	T45	T53	T61	T69	T77	T85
B	NC1	T6	T14	T22	T30	T38	T46	T54	T62	T70	T78	T86
C	STD1	T7	T15	T23	T31	T39	T47	T55	T63	T71	T79	T87
D	STD1	T8	T16	T24	T32	T40	T48	T56	T64	T72	T80	T88
E	T1	T9	T17	T25	T33	T41	T49	T57	T65	T73	T81	T89
F	T2	T10	T18	T26	T34	T42	T50	T58	T66	T74	T82	T90
G	T3	T11	T19	T27	T35	T43	T51	T59	T67	T75	T83	T91
H	T4	T12	T20	T28	T36	T44	T52	T60	T68	T76	T84	T92

Pipette

Predilution

Pipette into wells E1-H12.

Action on error: Raise alarm and stop

Aspirate 210ul from DILB (speed 2).

Aspirate 40ul from patient sample (speed 1).

Dispense 240ul into dilution plate ('Dilution 1:16') (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Dilb

Pipette into wells E1-H12.

Action on error: Raise alarm and stop

Aspirate 400ul from DILB (speed 2).

Dispense 400ul into dilution plate ('Dilution 1:16') (speed 5).

Incubate

Incubate for 1 (+15,-0) minutes at room temperature.

Incubate in ambient light.

--- Scheduling break ---

Pipette

NC

Pipette into wells B1.

Action on error: Raise alarm and stop

Aspirate 110ul from Yersinia G NC (speed 3).

Dispense 100ul into microplate (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Std

Pipette into wells C1-D1.

Action on error: Raise alarm and stop

Aspirate 110ul from Yersinia G Std (speed 3).

Dispense 100ul into microplate (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Sample 1:20

Pipette into wells E1-H12.

Action on error: Raise alarm and stop

Aspirate 30ul from DILB (speed 2).

Aspirate 80ul from dilution plate ('Dilution 1:16') (speed 3).

Dispense 100ul into microplate (speed 3).

Passive wash: waste: 0,0 seconds, wash: 0,0 seconds, volume: 0ul.

Incubate

Incubate for 60 (+5,-5) minutes at 37,0+/-1,0°C.

Wash

Perform 4 x 300ul wash cycles using Wash buffer on strips 1-12.

Dispense rate 75.

With aspirate sweep.
 Partial plate mode: as quick as possible.

Perform 2 aspirate cycle(s) on strips 1-12.
 With aspirate sweep.
 Partial plate mode: as quick as possible.

Dispense
 Add 8x100ul (+50ul oversoak) of APC (Conjugate) IgG (3/3) into wells B1-H12.

Incubate
 Incubate for 30 (+1,-1) minutes at 37,0+/-1,0°C.

Wash
 Perform 4 x 300ul wash cycles using Wash buffer on strips 1-12.
 Dispense rate 75.
 With aspirate sweep.
 Partial plate mode: as quick as possible.

Perform 2 aspirate cycle(s) on strips 1-12.
 With aspirate sweep.
 Partial plate mode: as quick as possible.

Dispense
 Add 8x100ul (+50ul oversoak) of pNPP (Substrate) (2/3) into wells A1-H12.

Incubate
 Incubate for 25 (+1,-1) minutes at 37,0+/-1,0°C.

Dispense
 Add 8x100ul (+50ul oversoak) of STOP (2/3) into wells A1-H12.

Read
 Reader settings

Shake for 5 seconds.
 Using wavelengths: 405nm,620nm.
 Average inputs. Use arithmetic mean.
 Report results to 3 decimal places.
 Auto scientific mode for:

Average blank mode.
 Blank calculation on final results only.

Validation criteria
 SLW<0.25
 $0.5 * \text{RefValue} \leq \text{Std} \leq 1.7 * \text{RefValue}$
 $0.8 \leq \text{STD}1.1 / \text{Std}1.2 \leq 1.25$
 NC<lcutoff*Std
 Full V.C. report.
 Do not display results if V.C. fails.

Qualitative Settings
 If 'Sample<lcutoff*Std' Then Result:='neg'
 If 'Sample>ucutoff*Std' Then Result:='pos'
 Default result := ?

pos Indicates a positive result
 ? Indicates an equivocal result
neg Indicates a negative result

No histogram required
 Average inputs. Use arithmetic mean.
 Report results to 3 decimal places.
 Auto scientific mode for:

Quantitative
 STD1 = 0 U/mL
 Average standard replicates.
 Virion data model.
 Show predicted standard values.
 Display a concentration graph.
 X axis settings: label: 'auto.', minimum=5, maximum=500, units=auto., show gridlines, logarithmic
 Y axis settings: label: 'auto.', minimum=0, maximum=3, units=1, show gridlines

Dilutions:
 T1 = 1 T2 = 1 T3 = 1 T4 = 1

T5 = 1	T6 = 1	T7 = 1	T8 = 1
T9 = 1	T10 = 1	T11 = 1	T12 = 1
T13 = 1	T14 = 1	T15 = 1	T16 = 1
T17 = 1	T18 = 1	T19 = 1	T20 = 1
T21 = 1	T22 = 1	T23 = 1	T24 = 1
T25 = 1	T26 = 1	T27 = 1	T28 = 1
T29 = 1	T30 = 1	T31 = 1	T32 = 1
T33 = 1	T34 = 1	T35 = 1	T36 = 1
T37 = 1	T38 = 1	T39 = 1	T40 = 1
T41 = 1	T42 = 1	T43 = 1	T44 = 1
T45 = 1	T46 = 1	T47 = 1	T48 = 1
T49 = 1	T50 = 1	T51 = 1	T52 = 1
T53 = 1	T54 = 1	T55 = 1	T56 = 1
T57 = 1	T58 = 1	T59 = 1	T60 = 1
T61 = 1	T62 = 1	T63 = 1	T64 = 1
T65 = 1	T66 = 1	T67 = 1	T68 = 1
T69 = 1	T70 = 1	T71 = 1	T72 = 1
T73 = 1	T74 = 1	T75 = 1	T76 = 1
T77 = 1	T78 = 1	T79 = 1	T80 = 1
T81 = 1	T82 = 1	T83 = 1	T84 = 1
T85 = 1	T86 = 1	T87 = 1	T88 = 1
T89 = 1	T90 = 1	T91 = 1	T92 = 1
NC1 = 1			

Result annotation:
 IF(x(y)<5,Result:='<5',,)
 IF(x(y)>500,Result:='>500',,)
 IF(y<=A*(Std/RefValue),Result:='<5',,)
 IF(y>=D*(Std/RefValue),Result:='>500',,)
 IF(y<0,0,Result:='<5',,)
 IF(y>3,0,Result:='>500',,)

Average inputs. Use arithmetic mean.
 Report results to 1 decimal places.
 Auto scientific mode for:
 Header:

Report Settings

Header:	&I	&D	&Q
Body:		Laboratory Details Assay Header Lot Specific Values Validation criteria Quantitative Results 1 Combined Report	
Footer:	&V &F		&S &P

Combined Report
 Table of results :- Patient ID, Well Location, Reader value, Quant. 1 value, Qual. value, Flag
 Table order: STD, NC, EC, T, PC, CO

ASTM E 1394 Quant. 1 mean

End of protocol