Thank you for purchasing a SleepSense® sensor.

Description

Snoring sensors respond to snoring and other sounds in the audio range picked up through the skin, and convert them to a small analog voltage that provides a clear, reliable indication of the presence of these sounds. Piezo or dynamic sensors do not require a battery or interface. Electret microphones are powered by the system. They are easily activated by plugging the output connector into the appropriate input on the head box.

- Dynamic snore sensors offer a wide frequency range, and sensitivity to other respiratory sounds besides snoring.
- Piezo sensors are robust, and relatively insensitive to ambient noise.
- Electret microphones offer a wide frequency response and high sensitivity.

Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

Applying the Sensor

- The sensor is attached to the patient's neck.

To position the sensor, place the sensor with the opening side or the protruding side (for the Piezo sensor) against the skin. The sensor is intended for facial skin to hold the sensor in place. Ensure that vibrations are clearly transmitted to the sensor during snoring.

- Ask the patient to simulate snoring, and check that clear, strong signals are being transmitted before leaving the patient and starting the study.

Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the complete sensor assembly is thoroughly dry before reusing it.

2.5 p.d. Jutiklį galima valyti plastikui tinkamu antikoroziniu dezinfekantu

Technical Specifications:

Sensor:	Dynamic / Piezo	Electret
Cable:	White insulation wire	
Sensor design:	Round flexible TPE disc	
Recommended filter settings:	High pass: 50 HZ / Low pass: 500 HZ or as recommended in the system's instructions	
Signal Output:	Approximately 1 mV for loud snoring	Approximately 50 mV for loud snoring
Operating conditions:	5°C (40°F) - 40°C (104°F)	
Storage temperature:	-20°C (-4°F) - 60°C (140°F)	
Operating and storage humidity:	5% - 95% (Non-condensing)	

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