HI-TORQUE INFILTRAC™

Guide Wire

CE
2797

en English:	4
de Deutsch:	6
fr Français :	8
es Español:	10
it Italiano:	12
pt Português:	14
sv Svenska:	16
nl Nederlands:	18
da Dansk:	20
el Ελληνικά:	22
pl Polski:	24
hu Magyar:	26
cs Český:	28
tr Türkçe:	30
sk Slovensky:	32
bg Български:	34
ro Română:	36
ru Русский:	38
fi Suomi:	40
no Norsk:	42
et Eesti keel:	44
lv Latviešu:	46
lt Lietuvių:	48
sl Slovenščina:	50
hr Hrvatski:	52
sr Srpski:	54
uk Українська:	56
mk Македонски:	58

PPL2132783 (2021-05-12)

en: Graphical Symbols for Medical Device Labeling; **de:** Grafische Symbole zur Beschriftung von Medizinprodukten; **fr:** Symboles graphiques pour l'étiquetage des dispositifs médicaux ; **es:** Símbolos gráficos utilizados en el etiquetado de productos sanitarios; **it:** Simboli grafici per l'etichettatura di dispositivi medicali; **pt:** Símbolos Gráficos para Rotulagem de Dispositivos Médicos; **sv:** Grafiska symboler för märkning av medicinsk utrustning; **nl:** Grafische symbolen voor de etikettering van medische hulpmiddelen; **da:** Grafiske symboler beregnet til mærkning af medicinske produkter; **el:** Γραφικά σύμβολα για τη σήμανση ιατρικών συσκευών; **pl:** Symbole graficzne do oznaczania sprzętu medycznego; **hu:** Grafikus szimbólumok orvosi eszközök címkézéséhez; **cs:** Grafické symboly k označení zdravotnického prostředku; **tr:** Tıbbi Cihaz Etiketleri için Grafik Semboller; **sk:** Grafické symboly na označovanie zdravotníckej pomôcky; **bg:** Графични символи за етикетиране на медицински уреди; **ro:** Simboluri grafice pentru etichetarea dispozitivelor medicale; **ru:** Графические символы на этикетках медицинских изделий; **fi:** Lääketieteellisten laitteiden tuotetaroissa esiintyvät symbolit; **no:** Grafiske symboler for merking av medisinsk utstyr; **et:** Meditsiiniseadme märgistuse graafilised sümbolid; **lv:** Grafiskie simboli medicīnisko ierīču apzīmēšanai; **It:** Medicinos prietaisų etikečių grafiniai simboliai; **sl:** Grafični simboli za označevanje medicinskega pripomočka; **hr:** Grafički simboli za označavanje medicinskog uređaja; **sr:** Grafički simboli za označavanje medicinskog sredstva; **uk:** Графічні позначки на етикетках медичних виробів; **mk:** Графички симболи за означување медицински уреди

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HI-TORQUE INFILTRAC™ Guide Wire

CAUTIONS

This medical device should be used only by physicians trained in angiography and percutaneous transluminal coronary angioplasty (PTCA), and / or percutaneous transluminal angioplasty (PTA).

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

Refer to the instructions supplied with any interventional devices to be used in conjunction with the HI-TORQUE™ Guide Wire, for their intended uses, contraindications, and potential complications.

A Summary of the Safety and Clinical Performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the Basic UDI-DI 8717648CGW0001UZ. This is the SSCP location after the launch of the European Database on Medical Devices/EUDAMED.

This IFU as per European Regulation (EU) 2017/745 can be found on the Abbott eLabeling website at vascular.eifu.abbott.

DESCRIPTION

This Instructions for Use pertains to the HI-TORQUE INFILTRAC™ and the HI-TORQUE INFILTRAC PLUS Guide Wires.

The HI-TORQUE™ Family of Guide Wires are medical devices that are indicated for use in PTCA and PTA to treat patients with coronary or peripheral artery disease as specified in the indications section of these instructions. These devices are steerable guide wires available in several lengths and diameters. The distal tip is shapeable or, as an option, a preshaped "J" tip is available. Refer to the product label for product specifications (e.g., wire length, diameter, and length of tip radiopacity).

Coronary guide wires facilitate PTCA and PTA, which offer clinical benefits for treatment of coronary artery disease (CAD) and multiple types of diseases in peripheral vessels.

HI-TORQUE Extendible Guide Wires: Some HI-TORQUE Guide Wires have a modified proximal end that permits the attachment of the DOC™ Guide Wire Extension. Refer to the product label for Guide Wire Extension system compatibility. Joining the guide wire extension to the guide wire facilitates the exchange of one interventional device for another, while maintaining guide wire position in the anatomy. After the interventional device exchange has been completed, the extension can be detached and the guide wire can be used in its original capacity.

.014 HI-TORQUE Guide Wires with Proximal Markers: Brachial and femoral markers located on the proximal segment of the 0.014" (0.36 mm) guide wire aid in gauging guide wire position relative to the guiding catheter tip when using bare wire technique. These guide wires are compatible with guiding catheters that are at least 90 cm (brachial) or 100 cm (femoral) long.

HI-TORQUE Guide Wires with Distal Radiopaque Markers: Some HI-TORQUE Guide Wires have distal radiopaque markers. Refer to the product drawing on the label to determine the presence and location of the markers. Markers are represented by the following designation: ■. A series of spaced markers are provided proximal to the coil as references in determining lesion length.

HI-TORQUE Guide Wires with Hydrophilic Coating: Refer to the product label for presence of a hydrophilic coating.

HOW SUPPLIED

Sterile – Sterilized with electron-beam radiation. Refer to the product label for the sterilization method. Non-pyrogenic. Do not use if the package is open or damaged.

This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Contents – One (1) guide wire per pouch.

Storage – Store in a dry, dark, cool place.

INTENDED PURPOSE

Guide wires are intended to facilitate navigation across vascular anatomies and structures in order to deliver diagnostic and therapeutic devices.

INDICATIONS FOR USE

This guide wire is indicated to facilitate the delivery of catheter-based interventional devices during PTCA and PTA. This guide wire may be used with compatible stent devices during therapeutic procedures. The guide wire may be used to reach and cross a target lesion, provide a pathway within the vessel structure, facilitate the substitution of one diagnostic or interventional device for another, and to distinguish the vasculature. This guide wire may also be used to cross or assist in crossing *de novo* chronic total coronary occlusions (CTO).

CONTRAINDICATIONS

Not indicated for use in the cerebral vasculature.

WARNINGS

This device is not designed for use with atherectomy devices. The safety and effectiveness of the use of this HI-TORQUE™ Guide Wire with atherectomy devices are unestablished.

The lifetime of the device is the labeled shelf life. This device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and / or reuse. The safety and effectiveness of this device have not been established after being reprocessed for multiple uses.

This device has a hydrophilic coating at the distal end of the device for a length of 25.6 cm, and a silicone based hydrophobic coating at the proximal end for a length of 133 cm (on 190 cm wires) and 250 cm (on 300 cm wires), and a polytetrafluoroethylene (PTFE) hydrophobic coating at the proximal end for a length of 156.5 cm (on 190 cm wires) and 266.5 cm (on 300 cm wires).

Please refer to section PREPARATION FOR USE for further information on how to prepare and use this device to ensure it performs as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

The safety and effectiveness of the device have not been established, or is unknown, in vascular regions other than those specifically indicated.

When the guide wire is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation.

Persons with known history of allergies to any of the components of this device listed below may suffer an allergic reaction to this guide wire. Prior to its use on the patient, the patient should be counseled on the materials contained in the device and a thorough history of allergies must be discussed. This device contains stainless steel, platinum-nickel alloy, tin-silver alloy, gold-tin alloy, MICROGLIDE™ silicone coating, PTFE coating and polyvinylpyrrolidone (PVP) coating.

Carefully observe the instructions under "Do Not" and "Do" below. Failure to do so may result in vessel trauma, guide wire damage, guide wire tip separation, or stent damage. If resistance is observed at any time, determine the cause under fluoroscopy and take remedial action as needed. Use the most suitable guide wire for the lesion being treated.

Do Not:

- Push, auger, withdraw, or torque a guide wire that meets resistance.
- Torque a guide wire if the tip becomes entrapped within the vasculature.
- Allow the guide wire tip to remain in a prolapsed condition.
- Deploy a stent such that it will entrap the wire between the vessel wall and the stent.

Do:

- Advance or withdraw the guide wire slowly.
- Use the radiopaque marker of the interventional device to confirm position.
- Examine the tip movement under fluoroscopy before manipulating, moving, or torquing the guide wire.
- Observe the wire under fluoroscopy for tip buckling, which is a sign of resistance.
- Maintain continuous flush while removing and reinserting the guide wire to prevent air from entering the catheter system.
- Perform exchanges slowly to prevent air entry and / or trauma.
- When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and that the tip is parallel to the vessel wall.
- Use extreme caution when moving a guide wire through a non-endothelialized stent, or through stent struts, into a bifurcated vessel. Use of this technique involves additional patient risks, including the risk that the wire may become caught on the stent strut.
- Consider that if a secondary wire is placed in a bifurcation branch, this wire may need to be retracted prior to stent deployment, because there is additional risk that the secondary wire may become entrapped between the vessel wall and the stent.

For the HI-TORQUE INFILTRAC™ Guide Wire family only: the HI-TORQUE INFILTRAC Guide Wire family of guide wires have distal ends of varying stiffness. Operate these guide wires carefully so as to not injure the blood vessel, observing the information in these instructions. The higher torque performance, stiffer distal ends, and / or higher advancement force may present a higher risk of perforation or injury than a guide wire with a more pliable distal end. Therefore, use the guide wire with the least stiff distal end that will treat the lesion, and use extreme care to minimize the risk of perforation or other damage to blood vessels.

LIMITATIONS

These HI-TORQUE™ Guide Wires do not have a specific limitation in duration of use but are designed for vascular procedures that vary in duration depending on lesion complexity.

Percutaneous transluminal coronary angioplasty (PTCA) should only be performed at centers where emergency coronary artery bypass graft (CABG) surgery is available.

PRECAUTIONS

Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and / or inaccurate torque response.

Confirm the compatibility of the guide wire diameter with the interventional device before actual use.

Please refer to the indication on the label and the Instructions for Use to confirm the appropriate vasculature that this guide wire may be used in. Failure to abide by the above recommendation may result in size mismatch of blood vessel and guide wire, which can result in vessel injury, such as, but not limited to, perforation, dissection, rupture, and avulsion.

The safety and effectiveness of this HI-TORQUE™ Guide Wire have not been established in the pediatric population.

Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system, because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit guide wire movement. It is recommended that the user determine the source of resistance, exercise caution when removing the device and / or other components as a unit and exchange the device for a new one to complete the procedure.

Never attach the torque device to the modified portion of the proximal end of the extendible guide wire; otherwise, guide wire damage may occur, preventing the ability to attach the DOC™ Guide Wire Extension.

Avoid abrasion of the hydrophilic coating. Do not withdraw or manipulate the hydrophilic-coated wire through a metal cannula or sharp-edged object. Manipulation, advancement, and / or withdrawal through a metal device may result in destruction and / or separation of the outer coating, which may cause coating material to remain in the vasculature. This in turn may lead to unintended adverse events requiring additional intervention.

When wet, a hydrophilic coating increases the lubricity of the guide wire surface.

Do not soak the device for longer than four hours when the device is not in use. Avoid pre-soaking devices for longer than instructed, as this may impact the coating performance.

The coating swells when exposed to aqueous media, but does not have any impact on device use.

The integrity and performance of the device coating can be negatively impacted by preparation with incompatible media or solvents. Please take note of the following important recommendations:

- Avoid wiping the device with dry gauze as this may damage the device coating.
- Avoid excessive wiping of the coated device.
- Avoid using alcohol, antiseptic solutions, or other solvents (e.g., use of medications and rotoflush) to pre-treat the device because this may cause unpredictable changes in the coating which could negatively affect the safety and performance of the guide wire.

Attempting to alter the shape of devices by bending, twisting, or similar methods beyond instructed methods may compromise the coating integrity, and that damage to the coating may not always be noticeable to the naked eye.

One or more components of this device may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0.

Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

After use, this product may be a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

ADVERSE EVENTS

Potential adverse events associated with use of this device may include the following, but are not limited to:

- Allergic reaction or hypersensitivity to latex, contrast agent, anesthesia, device materials, and drug reactions to anticoagulation, or antiplatelet drugs
- Perforation
- Dissection
- Bleeding
- Stenosis / restenosis
- Occlusion
- Arrhythmia
- Myocardial infarction
- Embolism
- AV fistula
- Aneurysm / pseudoaneurysm
- Infection
- Vasoconstriction
- Vasospasm
- Hypotension
- Hypertension
- Death

REPORTING OF ADVERSE EVENTS

The user and / or patient should report any serious incident that has occurred in relation to this HI-TORQUE™ Guide Wire to the manufacturer and the competent authority of the European Member State in which the user and / or patient is established.

PREPARATION FOR USE

Contraindications, warnings, limitations and intended uses of interventional devices compatible with HI-TORQUE™ Guide Wires are described in the instructions supplied with the respective devices.

Prior to the interventional procedure, carefully examine all equipment to be used, including the interventional device, for defects. Do not use if the device is damaged or defective.

1. Prepare the interventional device according to the manufacturer's instructions. Be sure to flush the guide wire lumen before introducing the guide wire.
2. Remove the guide wire from the dispenser by pushing the exposed section of the guide wire into the dispenser until the guide wire tip and a portion of the core exit the end of the hoop. Then grasp the core of the wire to remove it totally from the dispenser. Avoid damaging the fragile guide wire tip. Do not grasp the tip of the wire while removing it from the dispenser.
3. If indicated, the guide wire tip may be carefully shaped using standard tip-shaping practices. Do not use a shaping instrument with a sharp edge.

HI-TORQUE Guide Wires with Hydrophilic Coating

1. Before removing the guide wire from the wire dispenser, inject normal saline into the hub end of the dispenser to thoroughly wet the complete surface of the guide wire.
2. Carefully remove the guide wire from the dispenser, as suggested above in **PREPARATION FOR USE**, step 2. If the guide wire cannot be removed easily from the dispenser, inject more normal saline and attempt to remove the guide wire again.
3. Do not reinsert the guide wire into the dispenser once it has been removed.
4. If the surface of the hydrophilic-coated wire becomes dry, wetting the surface with normal saline will renew the hydrophilic effect. Be sure to thoroughly rewet the guide wire before reintroducing it into an interventional device.
5. After the guide wire is withdrawn from the body, it should be wiped clean with saline-soaked gauze and kept wet.

DIRECTIONS FOR USE

Over-the-Wire Type Systems (Preload Technique)

1. Carefully insert the guide wire through the guide wire lumen hub of the interventional device.
2. Advance the guide wire until the tip of the guide wire is just proximal to the tip of the interventional device.
3. If using a guiding catheter, engage the guiding catheter and insert the interventional device / guide wire assembly through the hemostatic valve. Advance the system through the guiding catheter until it is just proximal to the tip of the guiding catheter.
4. Tighten the hemostatic valve to create a seal around the interventional device. Ensure intentional guide wire movement is still permitted.
5. Attach the torque device to the guide wire, if desired.
6. Under fluoroscopy, advance the guide wire out of the interventional device while securing the interventional device in place. Use the torque device to steer the guide wire across the lesion.
7. Secure the guide wire in place while tracking the interventional device over it and into the lesion.
8. If a different tip configuration or guide wire is indicated, carefully remove the guide wire while observing guide wire movement under fluoroscopy.
9. Reshape the guide wire tip according to standard practice or prepare the next guide wire to be used.
10. Reinsert the guide wire following steps 1 through 7 of this section.

Rail Type Systems (Bare Wire Technique)

1. Engage the guiding catheter and then insert a guide wire introducer through the hemostatic valve attached to the guiding catheter.
2. Carefully insert the distal tip of the guide wire through the introducer and into the guiding catheter.
3. If a metal guide wire introducer was used, be sure to remove it before withdrawing or further manipulating the wire.
4. Advance the guide wire to the appropriate proximal marker. When the proximal marker is aligned with the knurled knob of the hemostatic valve, the guide wire tip is just proximal to the guiding catheter tip.

Note: Use the most distal marker as a distance gauge when using a 90 cm brachial guiding catheter, and the most proximal marker as the distance gauge when using a 100 cm femoral guiding catheter.

5. Attach the torque device.
6. Under fluoroscopy, advance the guide wire out of the guiding catheter and into the selected vessel. Use the torque device to steer the guide wire across the lesion.

7. If a different tip configuration or guide wire is indicated, the guide wire may be removed as follows:

- a. Open the hemostatic valve and the flush line on the coronary manifold. Slowly withdraw the guide wire while observing guide wire movement under fluoroscopy.
- b. Close the hemostatic valve and coronary manifold flush line.

8. Reshape the guide wire tip according to standard practice or prepare the next guide wire.

9. Reinsert the guide wire following steps 2 through 6 of this section.

10. Remove the torque device and the guide wire introducer from the guide wire.

11. Secure the guide wire while tracking the interventional device over it and into the lesion.

INSTRUCTIONS FOR INTERVENTIONAL DEVICE EXCHANGE PROCEDURE

Over-the-Wire Type Systems (Preload Technique)

1. If using a HI-TORQUE™ Exchange Guide Wire, proceed to step 3.
2. If using a HI-TORQUE Extendible Guide Wire, extend the guide wire using the DOC™ Guide Wire Extension. Refer to the instructions included with the extension for specific information on the use of the HI-TORQUE Guide Wire as an exchange-length guide wire.
3. Maintain guide wire position while withdrawing the interventional device over the exchange-length guide wire.
4. Prepare the other interventional device per manufacturer's instructions.
5. Load the device onto the guide wire and advance it over the HI-TORQUE Exchange Guide Wire and across the lesion.
6. Proceed according to standard medical practices.

Rail Type Systems (Bare Wire Technique)

1. Maintain guide wire position while withdrawing the device over the guide wire.
2. Prepare the next interventional device per manufacturer's instructions.
3. Load the device onto the guide wire and advance it over the HI-TORQUE Guide Wire and across the lesion.
4. Proceed according to standard medical practices.

Reference Abbott website for patent markings: www.abbott.com/patents

™ Indicates a trademark of the Abbott group of companies.

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© 2021 Abbott. All Rights Reserved.