ENGLISH

IMPORTANT INFORMATION Please Read Before Use

CODMAN® Disposable Perforator

 Only

STERILE	R
---------	---

FAILURE TO READ AND FOLLOW THESE INSTRUCTIONS MAY RESULT IN SERIOUS PATIENT INJURY OR DEATH.

Recommendation

A second sterile perforator or alternative method must always be available as backup in the event sterile technique is compromised. Any resterilization of this device can cause malfunction resulting in serious patient injury or death.

Indications

The CODMAN® Disposable Perforator is for use in perforating the cranium. When properly used, it is designed to automatically disengage once perforation is accomplished and when pressure is removed from the drill point.

Description

The CODMAN Disposable Perforator is a single use device. It is a disposable perforator with a Hudson end and is available in three color-coded sizes:

14mm (blue)- catalog no. 26-1221
11mm (green)- catalog no. 26-1222
9mm (yellow)- catalog no. 26-1223

WARNINGS

It is extremely important that this section be read carefully before attempting to use the perforator.

- The surgeon must always be aware of the possibility that the perforator may not disengage.
- Exercise extreme caution when the perforator nears the point of complete perforation. The underlying bone may not be sound and can fracture if excessive pressure is used, causing the perforator to penetrate the brain.
- The perforator can nick or tear dura or brain at the time of disengagement if conditions such as adherent dura, high intracranial pressure, or other underlying abnormalities exist in the area of perforation.
- The perforator can nick or tear dura or brain if it is not held perpendicular to the skull at the point of perforation.

Integra LifeSciences Production Corporation
11 Cabot Boulevard
Mansfield, MA 02048 USA

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CE
2797

LCN 205434-001 Rev D 10/19 1295483-2
© 2019 Integra LifeSciences Corporation. All rights reserved.

5. Excess deviation from the perpendicular can cause malfunction that could result in patient injury or death.

6. The CODMAN Disposable Perforator is designed for perforating normal, healthy cranial bone. Exercise extreme caution when performing skull perforations on infants, children, the elderly, and diseased bone because of the varying consistencies and thicknesses of bone.

7. Integra does not recommend use of the Disposable Perforator in a previously drilled burr hole or in surrounding areas. The composition of the restoration differs from healthy bone, thereby increasing the risk of penetration into the brain.

8. Exercise extreme caution when making perforations in thin bone (<1.75mm) such as the Temporal Bone.

9. If the perforator is stopped or prematurely disengages in the middle of a burr hole, there is the possibility that it may not reengage when attempting to drill again in the same burr hole.

10. The CODMAN Disposable Perforator (Figure 1) is to be used for one procedure only. Resterilization by any method could cause a malfunction that could result in patient injury or death. Do not use if perforator has been resterilized by any means.

11. Do not use the perforator if the label is damaged or missing from the outer plastic sleeve.

PRECAUTIONS

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.

Sterility

This product is for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Do not reuse. Discard after one procedure. A component of the device is specifically designed to fail if resterilized to prohibit reuse, which could cause impairment of the device's function.

Integra will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the individual package is not opened or damaged, the product is sterile.

Testing Procedure

PERFORM THE FOLLOWING TEST PROCEDURE BEFORE EACH PERFORATION. FAILURE TO CONDUCT THIS PROCEDURE COULD RESULT IN SERIOUS PATIENT INJURY OR DEATH.

CAUTION: Place gauze between the perforator point and thumb when performing the test to eliminate the possibility of puncturing or tearing the glove.

1. Visually examine the perforator and its packaging for any obvious defects, damage, or distortion. Read and inspect the label on the perforator. Do not use perforator if brown stripe is visible on perforator label as this indicates the perforator has been resterilized by ethylene oxide contrary to the labeling.

2. Do not use perforator if deformation of the colored sleeve is present as this indicates that the perforator has been resterilized using steam sterilization.

3. With adequate pressure applied on the perforator point, engagement should occur as the Hudson End is rotated. (See Figure 2.)

If engagement does not occur, do not use perforator.

4. When engagement occurs, thumb pressure on the perforator point should produce a smooth, positive spring action. (See Figure 3.) If smooth, positive spring action does not occur, do not use perforator.

5. The Hudson End should rotate smoothly within the perforator body when the unit is in the disengaged position. (See Figure 4.) If the Hudson End does not rotate smoothly, do not use the perforator.

General Operating Procedure

- Insert perforator securely into driver.
- Press the perforator on the skull.
- Manually rotate the Hudson Snap-Lock on the driver until the perforator is in the engaged position.
- With the driver and attached perforator positioned perpendicular to the skull, apply power and maintain constant but controlled pressure until disengagement occurs.
- Repeat steps 1, 2, and 3 of the Testing Procedure and steps 2, 3, and 4 of the General Operating Procedure for each perforation.

PRODUCT INFORMATION DISCLOSURE

INTEGRA LIFESCIENCES CORPORATION ("INTEGRA") HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA WARRANTS THAT THESE PRODUCTS SHALL CONFORM TO THE PRODUCT LIMITED WARRANTY AS PROVIDED IN THE PRODUCT LABELING OR APPLICABLE PRODUCT CATALOG. THIS WARRANTY IS EXCLUSIVE, AND INTEGRA DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THESE PRODUCTS. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Codman is a registered trademark of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.

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