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3 Levels / 1 Box

INTENDED USE

i-Smart Electrolyte Quality Control is an assayed quality control material used for monitoring the performance of i-Smart 30 analyzers that measure electrolytes (Na⁺, K⁺, Cl⁻ and Ca²⁺) by potentiometry using ion-selective electrodes.

PRODUCT DESCRIPTION

i-Smart Electrolyte Quality Control is provided in three levels for performance monitoring within the clinical range. It is a buffered aqueous solution containing electrolytes (Na⁺, K⁺, Cl⁻ and Ca²⁺) and preservatives. The solution of each level is packaged in sealed plastic bottles, each containing 10mL of solution. This control contains no human or biological materials.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 2 to 8 °C. Following the first opening of the bottle, all analytes will be stable for ten (10) weeks when stored tightly capped at 2 to 8 °C.

This product is shipped under ambient temperature. The product should not be exposed to very extreme temperatures (below -23 °C or greater than 45 °C) for more than fourteen (14) consecutive days.

DIRECTIONS FOR USE

Open the QC bottle and gently squeeze out the liquid into a sample cup (at least 60 µL). Immediately introduce the liquid to the analyzer following the instructions in the QC sampling section of i-Smart analyzers' operator's manuals. After use, discard all materials in a safe and acceptable manner that is in compliance with all country, state, and local requirements.

EXPECTED VALUES

Expected values for this lot of product are listed on the reverse side of this sheet.

The expected values are provided as a guide in evaluating the performance of the i-Smart 30 analyzers and it is recommended that each laboratory establish its own means and acceptance ranges.

Incorrect sampling, storage, or other mishandling may cause the readings to deviate from the expected values.

LIMITATIONS

1. This product should not be used past the expiration date.
2. This product is not for use as a calibration standard.

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Medical Technology Promedt Consulting GmbH
Altenhofstrasse 80, 66386 St. Ingbert, Germany

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i-SENS, Inc.
[Manufacturing facility] i-SENS, Inc.