AQT90 FLEX

TnI Test Kit

Contents

en	Tnl Test Kit.....	2
bg	Набор тестове за Tnl.....	8
cs	Souprava testů Tnl.....	15
da	Testkit til Tnl.....	21
de	Tnl-Testkit.....	27
el	Κιτ εξέτασης Tnl.....	33
es	Kit de Test de Tnl.....	40
et	Tnl mõõtekomplekt.....	47
fi	Tnl-testipakkaus.....	53
fr	Kit de tests Tnl.....	59
hr	Testni komplet Tnl.....	65
hu	Tnl tesztkészlet.....	71
it	Kit per il test Tnl.....	77
lt	Tnl tyrimo rinkinys.....	84
lv	Tnl testēšanas komplekts.....	90
nl	Tnl-testkit.....	96
no	Tnl-testkit.....	102
pl	Zestaw testowy do oznaczania Tnl.....	108
pt	Kit de Testes de Tnl.....	115
ro	Kit de testare pentru Tnl.....	121
ru	Набор тестов на Tnl.....	128
sk	Súprava na testovanie Tnl.....	135
sv	Tnl-testkit.....	141
tr	Tnl Test Kiti.....	147
zh	Tnl 检测试剂盒.....	153
References		159

en

TnI Test Kit

[REF] 942-903

Intended use

For *in vitro* diagnostic use.

The TnI Test is an *in vitro* diagnostic assay for the quantitative determination of cardiac TnI in EDTA or lithium-heparin whole-blood or plasma specimens on the AQT90 FLEX analyzer in point of care and laboratory settings. It is indicated for use as an aid in the diagnosis of myocardial infarction and in the risk stratification of patients with acute coronary syndromes with respect to their relative risk of mortality.

The TnI CAL Cartridges are *in vitro* diagnostic reagents intended for the calibration adjustment of the TnI Test on the AQT90 FLEX analyzer by establishing points of reference to estimate TnI values.

NOTE: The cartridges are ready for use.

Test

TnI is a protein of the muscular thin filaments and part of the troponin-tropomyosin complex in striated muscle. Troponin complex consists of three components: troponin C (TnC), troponin T (TnT) and troponin I (TnI). This ternary troponin complex is a calcium-sensitive molecular apparatus that regulates the interaction of actin and myosin [1,2].

Both TnI and TnT exist in skeletal and cardiac isoforms. The fact that the cardiac isoforms are extremely specific to the heart tissue is the basis for their increased use as specific markers of cardiac muscle damage. The main part of plasma TnI is found in the form of a complex, mainly with TnC. Only a small part of circulating TnI is found in a non-complexed form [3,4].

Plasma concentrations of TnI in healthy individuals are extremely low or undetectable, which constitutes the basis for the usefulness of TnI in the detection of myocardial infarction (MI). TnI is released from necrotic cardiac myocytes into the bloodstream within hours (~4-8 hours) after the onset of chest pain. The peak TnI concentration is generally reached in 12-48 hours. Plasma TnI levels remain elevated several days after severe MI [5,6,7,8].

The exceptional cardiac specificity and high clinical sensitivity of the troponins shown in clinical studies have markedly increased their role in the diagnostics of MI. According to the Third universal definition of myocardial infarction published in October 2012 cardiac biomarkers like troponin or CKMB are central in the diagnosis of myocardial infarction, together with ECG and the patient history [9].

Calibration adjustment

The assay has been factory-calibrated using a 9-point calibration. The factory-defined calibration data is provided in the form of a barcode. This data is used as

a lot-specific reference curve, which is adjusted to the correct analyzer-specific signal level by using a CAL Cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and analyzer-specific calibration curve for the assay. When a CAL Cartridge with a new lot number is inserted, the analyzer requests the factory-defined calibration data so that it can do the calibration adjustment. If the user does not insert a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after a calibration adjustment is done with a CAL Cartridge with the same lot number.

See the AQT90 FLEX *instructions for use* document for details.

Principles of operation

The AQT90 FLEX analyzer is a fully automatic continuous-access analyzer utilizing immunoassay technology and time-resolved fluorometric detection. See the AQT90 FLEX *instructions for use* document for details.

The AQT90 FLEX analyzer is operated via a touch screen. The analyzer aspirates samples from capped sample tubes that are inserted manually.

The analyzer does all assay steps automatically and reports quantitative results. See the AQT90 FLEX *instructions for use* document for details.

Contents

- 10 TnI Test Cartridges for 160 tests
- One TnI CAL Cartridge for one calibration adjustment
- A factory-defined calibration data sheet with barcode

Components

- Each Test Cartridge contains 16 analyte-specific assay cups
- The CAL Cartridge contains 8 analyte-specific background cups and 8 cups with added antigen

Storage

Store unopened at 2-8 °C (36-46 °F).

Shelf life

If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the expiry date shown on the pouch label.

On-board stability

The value (in days) is printed on the cartridge pouches alongside the following symbol:

Chemical composition

Reactive ingredients:

- Mouse monoclonal anti-TnI capturing and tracing antibodies: approx. 400 ng/cup
- Native TnI (only in the 8 CAL Cartridge cups with added antigen): approx. 0.03 ng/ cup

Other ingredients:

- Bovine serum albumin
- Bovine γ -globulin
- Mouse IgG (blocker substances for heterophilic interferences)

Main buffer component: TRIS (Tris hydroxymethylaminomethane)

Toxic ingredient: Sodium azide <0.1 %

Data in cartridge barcodes

- Checksum
- Barcode type
- Parameter code
- Lot no.
- Cartridge ID
- Expiry date
- Barcode version
- On-board stability

NOTE: The barcodes are lot-specific.

About using the cartridges

WARNING – Risk of incorrect results

Cartridges are sensitive to humidity and must therefore not be taken out of their sealed pouches until just before use. Only use cartridges that have been kept according to storage specifications and are not damaged.

WARNING – Risk of incorrect results

Do not use a cartridge that has already been used on another analyzer, otherwise you risk obtaining incorrect results.

WARNING – Risk of infection

This product contains materials from human source and should be handled and disposed of as potentially infectious. The plasma from human donors used in the manufacture of this product has been tested and found non-reactive to HbsAg, HIV-1, HIV-2 and HCV.

Product calibrator traceability

The calibration of the TnI assay is traceable to NIST SRM 2921.

Sample collection and preparation

Sample tubes approved for use with the analyzer

See the *AQT90 FLEX instructions for use* document for details.

Anticoagulants

Use either EDTA or lithium heparin.

Sample size

The minimum sample volume is 2 mL (independent of the number of tests requested).

Sample type

Use whole-blood or plasma samples.

Sample collection

Collect blood samples by venipuncture.

WARNING – Risk of incorrect results

For sample collection, do not use tubes that contain a gel.

Sample handling and storage

Samples that cannot be analyzed immediately after collection must be handled and stored correctly before they are analyzed. The table provides an overview.

Sample type	Handling and storage	Analyze within...
Whole blood	1. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	2 hours ^a
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	2 hours ^a
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Refrigerate at 2-8 °C (36-46 °F).	24 hours ^a
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Freeze at -18 °C (<-1 °F).	>24 hours

^a Sample stability.

NOTE: Never freeze whole-blood samples.

To thaw frozen plasma samples

1. Let the sample thaw at room temperature.

2. Gently invert the sample tube 5 times to mix the sample.
3. Centrifuge the thawed sample to remove possible micro clots.

NOTE: Do not freeze and thaw a sample more than once.

Maximum sample age

This is the period of time after a sample is collected within which whole-blood samples kept at an ambient temperature of 18-25 °C (65-77 °F) must be analyzed. See 'Sample handling and storage' for details.

Sample analysis

Samples must be analyzed immediately after they are collected. When this is not possible, samples must be handled and stored correctly and analyzed within the time period given in 'Sample handling and storage'.

WARNING – Risk of falsely low results

Analyze samples within the time periods given in 'Sample handling and storage', otherwise you risk falsely low TnI results.

See the *AQT90 FLEX instructions for use* document for detailed instructions.

NOTE: Keep the analyzer connected to the mains power supply so that it is always ready for use.

Results

Results are shown on the screen and can be printed as they become available. Results are saved in a data log. If the analyzer detects an error condition, a message is attached to the analysis and no result is given.

WARNING – Risk of incorrect interpretation of results

Current international guidelines should always be followed when interpreting TnI results.

WARNING – Risk of incorrect results

Always interpret with caution any assay results that are inconsistent with other available clinical data and a patient's medical history. If samples contain heterophilic antibodies, immunoassays could give falsely elevated or depressed results. Blockers have been added to this assay to minimize potential interference from heterophilic antibodies.

Calibration adjustment

See the *AQT90 FLEX instructions for use* document for details.

Results are saved in a data log and can be printed. If a calibration adjustment fails, a message is shown.

NOTE: After every calibration adjustment, liquid quality control (LQC) measurements must be done.

Quality control

Liquid quality control (LQC) measurements must be done:

- after a calibration adjustment
- in accordance with local, national and/or international regulations or accreditation requirements

See the *AQT90 FLEX instructions for use* document for details.

LQC material

See the *AQT90 FLEX Ordering information* document for details.

Performance characteristics

Analytical specificity

The cross-reactivity was studied between skeletal troponin T, skeletal troponin I, cardiac troponin C and cardiac troponin T at concentrations of 1,000 µg/L (ng/mL) to determine the analytical specificity of the TnI assay. All observed cross-reactivities were <0.01 %.

Analytical sensitivity and reportable range

Analytical sensitivity	TnI µg/L (ng/mL)
Limit of blank (LoB)	<0.005
Limit of detection (LoD)	<0.010

The reportable range of the assay is 0.010-25 µg/L (ng/mL).

Imprecision

Five levels of frozen native or antigen-spiked lithium-heparin plasma pools were analyzed at 4 point-of-care (POC) test sites, on one analyzer. Three Test Kit lots were used per site. Each level was analyzed in duplicate in 2 separate runs for a minimum of 20 test days to determine the within-run, within-lab and between-lot imprecision on plasma.

Plasma pool	TnI ^a µg/L (ng/mL)	Within-run		Within-lab		N ^b
		SD ^c µg/L (ng/mL)	CV ^d (%)	SD µg/L (ng/mL)	CV (%)	
Sites 1, 2, 3 & 4	0.021	0.0025	11.9	0.0027	12.9	370
	0.035	0.0023	6.7	0.0027	7.7	368
	0.37	0.008	2.0	0.016	4.4	378
	9.2	0.14	1.5	0.28	3.1	380
	24	0.4	1.5	0.7	3.1	378

Plasma pool	Tnl ^a µg/L (ng/mL)	Within-run		Within-lab		N ^b
		SD ^c µg/L (ng/mL)	CV ^d (%)	SD µg/L (ng/mL)	CV (%)	
Site 1	0.021	0.0020	9.5	0.0022	10.4	94
	0.035	0.0030	8.5	0.0030	8.5	90
	0.39	0.008	2.1	0.018	4.6	96
	9.3	0.11	1.2	0.26	2.8	96
	24	0.4	1.7	0.7	2.9	96
Site 2	0.020	0.0015	7.1	0.0018	8.7	92
	0.035	0.0015	4.4	0.0024	6.8	92
	0.35	0.006	1.8	0.014	4.1	92
	9.1	0.14	1.5	0.28	3.1	92
	24	0.3	1.2	0.8	3.4	92
Site 3	0.022	0.0039	17.8	0.0042	18.9	98
	0.034	0.0019	5.5	0.0022	6.3	98
	0.37	0.008	2.0	0.017	4.6	98
	9.1	0.14	1.5	0.28	3.1	100
	23	0.3	1.5	0.8	3.5	98
Site 4	0.020	0.0017	8.4	0.0022	10.8	86
	0.035	0.0027	7.8	0.0034	9.9	88
	0.36	0.008	2.2	0.016	4.3	92
	9.2	0.17	1.9	0.30	3.3	92
	24	0.4	1.6	0.6	2.7	92

^a Mean concentration; ^b Number of measurements; ^c Standard deviation; ^d Coefficient of variation

Plasma pool	Tnl ^a µg/L (ng/mL)	Between-lot		N ^b
		SD ^c µg/L (ng/mL)	CV ^d (%)	
Sites 1, 2, 3 & 4	0.021	0.0009	4.2	370
	0.035	0.0011	3.2	368
	0.37	0.008	2.1	378
	9.2	0.24	2.6	380
	24	0.6	2.5	378

^a Mean concentration; ^b Number of measurements; ^c Standard deviation; ^d Coefficient of variation

Five levels of fresh native or antigen-spiked lithium-heparin whole blood samples were analyzed at 3 POC test sites on one analyzer. One Test Kit lot was used per site. Each level was analyzed in 5 replicates in 5 separate runs within a 2-hour period to determine the within-run and within-lab imprecision on whole blood.

Whole blood	Tnl ^a µg/L (ng/mL)	Within-run		Within-lab		N ^b
		SD ^c µg/L (ng/mL)	CV ^d (%)	SD µg/L (ng/mL)	CV (%)	
Site 1	0.023	0.0025	10.7	0.0025	10.7	25
	0.030	0.0031	10.2	0.0036	11.8	34
	0.51	0.014	2.8	0.017	2.8	25
	3.4	0.10	2.8	0.10	2.8	25
	16	0.9	5.4	0.9	5.4	25
Site 2	0.025	0.0018	7.1	0.0018	7.1	22
	0.041	0.0026	6.4	0.0026	6.4	25
	0.46	0.012	2.6	0.013	2.7	25
	4.9	0.10	2.0	0.10	2.0	25
	18	0.4	2.4	0.5	2.7	25
Site 3	0.024	0.0021	8.6	0.0022	9.4	25
	0.040	0.0023	5.8	0.0028	7.1	25
	0.43	0.011	2.5	0.012	2.7	25
	4.8	0.10	2.1	0.10	2.1	25
	18	0.4	2.1	0.4	2.2	25

^a Mean concentration; ^b Number of measurements; ^c Standard deviation; ^d Coefficient of variation

Reference values

Whole-blood (lithium-heparin and EDTA) samples and plasma (lithium-heparin and EDTA) samples were obtained from 231 apparently healthy individuals (106 women and 125 men) and analyzed using the TnI assay. The 99th percentile was determined to be 0.023 µg/L (ng/mL).

NOTE: This value should only be used as an example. Each laboratory should establish its own reference ranges.

Clinical sensitivity and specificity

Clinical sensitivity and specificity of the TnI assay was evaluated analyzing a clinically characterized panel consisting of lithium-heparin plasma specimens of 390 patients. Of this patient population, 79 patients were according to World Health Organisation (WHO) MONICA criteria diagnosed as definite AMI patients and 311 were diagnosed as possible or non-AMI patients.

The specimens included in the analysis were collected at the following time points:

- On admission to the hospital (0-2 hours)
- 6-9 hours after admission

The biochemical marker assay used was Elecsys Troponin T (3rd generation) and the cut-off for AMI was 0.10 µg/L.

Cut-off (99 th percentile)	Hours after admission			
	0-2 hours		6-9 hours	
	Sens. ^a	Spec. ^b	Sens.	Spec.
0.023 µg/L (ng/mL)	65 %	91 %	91 %	88 %

^a Sensitivity; ^b Specificity

Hook effect

No hook effect was found when TnI concentrations up to 1,470 µg/L (ng/mL) were measured.

Carry over

Carry over from a sample with high TnI value (~40 µg/L (ng/mL)) to an adjacent negative sample was determined to be <500 ppm (<0.020 µg/L (ng/mL)).

Interference tests

Blockers have been added to this assay to minimize potential interference from heterophilic antibodies, for example HAMA (human anti-mouse antibodies) and rheumatoid factor (RF). Tests showed that the blockers have the desired effect.

Hemolytic, lipemic and icteric samples were shown not to interfere with the assay.

The interference was tested by using plasma samples with approximately 0.04 µg/L (ng/mL) and 0.8 µg/L (ng/mL) of TnI spiked with the interfering substance at

concentrations about 3 times the upper therapeutic range. Blood component interference was tested in the same manner, but in plasma as well as whole blood samples.

The following concentrations of potentially interfering substances were found to have no notable effect on the TnI assay (interference ≤10 %).

Drugs and solvents	Tested concentration mg/L
Abciximab	90
Acetaminophen	200
Acetylcysteine	1,700
Acetylsalicylic acid	700
Allopurinol	40
Amiodarone	12.3
Ampicillin	60
Ascorbic acid	60
Atenolol	10
Caffeine	60
Captopril	5
Cinnarizine	1
Clopidogrel	200
Cocaine	1
Cyclosporine	6
Diclofenac	50
Digoxin	0.006
Diltiazem	6
Dipyridamole	0.005
Dopamine	1
Enalaprilate	0.3
Erythromycin	60
Ethanol	4,000
Furosemide	60
Heparin, low molecular weight	3,000 IU/L
Hydralazine	2.4
Ibuprofen	500

Drugs and solvents	Tested concentration mg/L
Isosorbide dinitrate	0.15
Methyldopa	17
Metronidazole	120
Milrinone	0.63
Nicotine (+/–)	1
Nicotinic acid	35
Nifedipine	0.40
Nitrofurantoin	4
Nitroglycerin	0.02
Oxazepam	5
Oxytetracycline	12
Phenobarbital	80
Phenylbutazone	600
Phenytoin	50
Probenecid	600
Procainamide	24
Propranolol	2
Quinidine	12
Simvastatin	0.024
Sulfamethoxazole	400
Tetracycline	15
Theophylline	40
Trimethoprim	40
Verapamil	2
Warfarin	10

Substances in blood	Tested concentration
Bilirubin unconjugated	270 mg/L
Bilirubin conjugated	270 mg/L
Biotin	3 µg/L
Cholesterol	5 g/L
Creatinine	50 mg/L
Fibrinogen	10 g/L

Substances in blood	Tested concentration
Glucose	55 mmol/L
Hemoglobin	10 g/L
Hemolysis index	50 mg/dL
Human albumin	100 g/L
Human IgG	25 g/L ^a
Triglycerides (Intralipid)	15 g/L
Urea	2 mg/L

^a Human IgG concentrations above 25 g/L may cause negative interference.

Method comparison

The Tnl assay (y) was compared to the commercially available cardiac troponin I assay Tnl-Ultra for the ADVIA Centaur immunoassay system (x). Lithium-heparin plasma samples in the range 0.010-25 µg/L (as measured on the AQT90 FLEX analyzer) were used.

The linear regression line and correlation coefficient were found to be:

$$y = 0.148x - 0.012; R^2 = 0.984 (n = 220)$$

Patents

Radiometer products may be covered by one or more patents and patent applications, see www.radiometer.com/en/legal/patents.

Explanation of graphical symbols

See the *AQT90 FLEX instructions for use* document.

References

See 'References' on the last page of this insert.

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References

1. Filatov VL, Katrukha AG, Bulargina TV, Gusev NB. Troponin: Structure, properties, and mechanism of functioning. *Biochemistry (Moscow)* 1999; 64(9): 969-85.
2. Farah CS, Reinach FC. The troponin complex and regulation of muscle contraction. *FASEB J* 1995; 9: 755-67.
3. Bodor GS, Porterfield D, Voss EM, Smith S, Apple FS. Cardiac troponin I is not expressed in fetal and healthy or diseased adult human skeletal tissue. *Clin Chem* 1995; 41: 1710-15.
4. Katrukha AG, Bereznikova AV, Esakova TV et al. Troponin I is released in bloodstream of patients with acute myocardial infarction not in free form but as complex. *Clin Chem* 1997; 43: 1379-85.
5. Wu AHB. Cardiac markers. In: Wild D, ed. *The immunoassay handbook*. London: Nature Publishing Group, 2001: 623-34.
6. Collinson PO, Boa FG, Gaze DC. Measurement of cardiac troponins. *Ann Clin Biochem* 2001; 38(5): 423-49.
7. Apple FS, Falahati A, Paulsen PR, Miller EA, Sharkey SW. Improved detection of minor ischemic myocardial injury with measurement of serum cardiac troponin I. *Clin Chem* 1997; 43: 2047-51.
8. Manuel Del Rey JM, Madrid AH, Valino JM et al. Cardiac troponin I and minor cardiac damage: biochemical markers in a clinical model of myocardial lesions. *Clin Chem* 1998; 44: 2270-76.
9. Thygesen K, Alpert JS, Jaffe AS, Simoons ML, Chaitman BR, White HD, the Writing Group on behalf of the Joint ESC/ACCF/AHA/WHF Task Force for the Universal Definition of Myocardial Infarction. Third universal definition of myocardial infarction. *Eur Heart J* 2012; 33(20): 2551-2567.

Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Denmark
www.radiometer.com

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.

© 2016 Radiometer Medical ApS. All rights reserved.
Code number: 990-872
Version: 201608X