
AQT90 FLEX

PCT (8-test) Test Kit

Contents

en	PCT (8-test) Test Kit.....	2
bg	Набор тестове PCT (8 теста).....	7
cs	Souprava testů pro PCT (8 testů).....	13
da	Testkit til PCT (8-test).....	19
de	PCT-(8-Test)-Testkit.....	24
el	Κιτ εξέτασης PCT (8 εξετάσεις).....	30
es	Kit de Test de PCT (8 tests).....	36
et	PCT mõõtekomplekt (8 mõõtmist).....	42
fi	PCT-testaussarja (8 testiä).....	48
fr	Kit de tests PCT (test-8).....	53
hr	Testni komplet PCT (8 testova).....	59
hu	PCT tesztkészlet (8 teszt).....	65
it	Kit per il test PCT (8 test).....	70
lt	PCT (8 tyrimų) tyrimų rinkinys.....	76
lv	PCT (8 testu) testēšanas komplekts.....	81
nl	PCT-(8 tests)-testkit.....	86
no	PCT-(8-test)-testkit.....	92
pl	Zestaw testowy do oznaczania PCT (8 testów).....	98
pt	Kit de testes (8 testes) de PCT.....	104
ro	Kit de testare pentru PCT (8 teste).....	110
ru	Набор тестов на PCT (8 тестов).....	116
sk	Súprava testov PCT (8 testov).....	122
sv	PCT-testkit (8-test).....	127
tr	PCT (8 test) Test Kiti.....	133
		139

en

REF 942-970

Intended use

For *in vitro* diagnostic use.

The PCT (8-test) Test Kit is an *in vitro* diagnostic assay for the quantitative determination of PCT in EDTA or lithium-heparin whole-blood or plasma specimens on the AQT90 FLEX analyzer in point-of-care and laboratory settings. It is indicated for use as an aid in the diagnosis of sepsis.

The PCT CAL Cartridges are *in vitro* diagnostic reagents intended for the calibration adjustment of the procalcitonin (PCT) Test on the AQT90 FLEX analyzer by establishing points of reference to estimate PCT values.

NOTE: The cartridges are ready for use.

Test

Procalcitonin (PCT) is a 116-amino acid polypeptide prohormone of calcitonin. PCT is currently considered by many researchers to be the most useful biomarker for severe systemic inflammation, infection and sepsis although not totally specific for sepsis [1,2]. It is a biomarker for differentiating patients with sepsis from those with non-infectious causes of systemic inflammations. PCT concentrations increase earlier and PCT has a shorter half-life than CRP, which is another frequently used biomarker for sepsis. Due to the favorable kinetics, PCT has the potential advantage of earlier disease diagnosis and better monitoring of disease progression. Early diagnosis and treatment are critical for survival rates.

Interpretation of PCT levels for patients in intensive care units (ICU) according to the literature [1,2,3]:

PCT level (ng/mL)	Interpretation
<0.05	Healthy
<0.50	Local infection, sepsis not likely
0.50-2.0	Systemic infection suggestive of sepsis
>2.0-10	High risk for sepsis
>10	Septic shock

Limitations

Evaluation of the PCT Test was done with EDTA and lithium heparin whole-blood and plasma samples in plastic tubes. Use of other anticoagulants or serum samples is not recommended. PCT Test results should be evaluated together with all available clinical and laboratory data. If results are not consistent with clinical observations, additional information may be required for diagnosis.

NOTE: In neonates or in patients with special medical conditions such as polytrauma, major surgery and severe burns, PCT levels may be elevated independently of bacterial infection. In these particular cases, levels usually decrease rapidly [4,5].

Calibration adjustment

The assay has been factory-calibrated using a 11-point calibration. The factory-defined calibration data is provided in the form of a barcode. This data is used as a lot-specific reference curve, which is adjusted to the correct analyzer-specific signal level by using a CAL Cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and analyzer-specific calibration curve for the assay. When a CAL Cartridge with a new lot number is inserted, the analyzer requests the factory-defined calibration data so that it can do the calibration adjustment. If the user does not insert a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after a calibration adjustment is done with a CAL Cartridge with the same lot number.

See the *AQT90 FLEX instructions for use* document for details.

Principles of operation

The AQT90 FLEX analyzer is a fully automatic continuous-access analyzer utilizing immunoassay technology and time-resolved fluorometric detection. See the *AQT90 FLEX instructions for use* document for details.

The AQT90 FLEX analyzer is operated via a touch screen. The analyzer aspirates samples from capped sample tubes that are inserted manually.

The analyzer does all assay steps automatically and reports quantitative results. See the *AQT90 FLEX instructions for use* document for details.

Contents

- 10 PCT Test Cartridges for 80 tests
- One PCT CAL Cartridge for one calibration adjustment
- One *PCT Factory Defined Calibration Data* sheet with barcode

Components

- Each Test Cartridge contains 8 analyte-specific assay cups
- The CAL Cartridge contains 8 analyte-specific background cups and 8 cups with added antigen

Storage

Store unopened at 2-8 °C (36-46 °F).

Shelf life

If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the expiry date shown on the pouch label.

On-board stability

The value (in days) is printed on the cartridge pouches alongside the following symbol:

Chemical composition

Reactive ingredients:

- Mouse monoclonal anti-PCT capturing antibodies: approx. 400 ng/cup
- Mouse monoclonal anti-PCT tracing antibodies: approx. 270 ng/cup
- Recombinant PCT (only in the 8 cups with added antigen): approx. 0.13 ng/cup

Other ingredients:

- Bovine serum albumin
- Bovine γ -globulin
- Mouse IgG (blocker substances for heterophilic interferences)

Main buffer component: TRIS (Tris hydroxymethylaminomethane)

Toxic ingredient: Sodium azide <0.1 %

Data in cartridge barcodes

- Checksum
- Barcode type
- Parameter code
- Lot no.
- Cartridge ID
- Expiry date
- Barcode version
- On-board stability

NOTE: The barcodes are lot-specific.

About using the cartridges

WARNING – Risk of incorrect results

Cartridges are sensitive to humidity and must therefore not be taken out of their sealed pouches until just before use. Only use cartridges that have been kept according to storage specifications and are not damaged.

WARNING – Risk of incorrect results

Do not use a cartridge that has already been used on another analyzer, otherwise you risk obtaining incorrect results.

Product calibrator traceability

The calibration of the PCT assay is traceable to in-house reference calibrators, which have been value-assigned to correlate with the B·R·A·H·M·S PCT sensitive KRYPTOR assay.

Sample collection and preparation

Sample tubes approved for use with the analyzer

See the *AQT90 FLEX instructions for use* document for details.

Anticoagulants

Use either EDTA or lithium heparin.

Sample size

The minimum sample volume is 2 mL (independent of the number of tests requested).

Sample type

Use whole-blood or plasma samples.

Sample collection

Collect blood samples by venipuncture.

WARNING – Risk of incorrect results

For sample collection, do not use tubes that contain a gel.

Sample handling and storage

Samples that cannot be analyzed immediately after collection must be handled and stored correctly before they are analyzed. The table provides an overview.

Sample type	To handle and store the sample	Analyze within this time period
Whole blood	1. Keep at an ambient storage temperature 18-25 °C (65-77 °F).	3 hours
Plasma	1. Separate the plasma immediately after the sample is collected. 2. Refrigerate at 2-8 °C (36-46 °F).	24 hours

Maximum sample age

This is the period of time after a sample is collected within which whole-blood samples kept at an ambient temperature of 18-25 °C (65-77 °F) must be analyzed.

See 'Sample handling and storage' for details.

Sample analysis

Samples must be analyzed immediately after they are collected. When this is not possible, samples must be handled and stored

correctly and analyzed within the time period given in 'Sample handling and storage'.

See the *AQT90 FLEX instructions for use* document for detailed instructions.

NOTE: Keep the analyzer connected to the mains power supply so that it is always ready for use.

Results

Results are shown on the screen and can be printed as they become available. Results are saved in a data log. If the analyzer detects an error condition, a message is attached to the analysis and no result is given.

 **WARNING – Risk of incorrect interpretation of results**

Current international guidelines should always be followed when interpreting PCT results.

 **WARNING – Risk of incorrect results**

Always interpret with caution any assay results that are inconsistent with other available clinical data and a patient's medical history. If samples contain heterophilic antibodies, immunoassays could give falsely elevated or depressed results. Blockers have been added to this assay to minimize potential interference from heterophilic antibodies.

Calibration adjustment

See the *AQT90 FLEX instructions for use* document for details. Results are saved in a data log and can be printed. If a calibration adjustment fails, a message is shown.

NOTE: After every calibration adjustment, liquid quality control (LQC) measurements must be done.

Quality control

- Liquid quality control (LQC) measurements must be done:
- after a calibration adjustment
 - in accordance with local, national and/or international regulations or accreditation requirements

See the *AQT90 FLEX instructions for use* document for details.

LQC material

See *Ordering information* in the *AQT90 FLEX instructions for use* document for details.

Performance characteristics

Analytical specificity

The cross reactivity was studied for calcitonin, katacalcin, α -human CGRP and β -human CGRP at concentrations of 200 ng/mL (μ g/L) to determine the analytical specificity of the PCT assay. All observed cross reactivities were ≤ 0.05 %.

Analytical sensitivity and reportable range

The limit of blank (LoB) determination was based on 180 replicates of a blank sample. Three analyzers and 3 Test Kit lots were used.

The limit of detection (LoD) was based on a precision profile with 40 replicates of each of 9-11 low-level samples. Six analyzers and 3 Test Kit lots were used.

The limit of quantitation (LoQ) was based on a precision profile with 40 replicates of each of 16 low-level samples. Six analyzers and 3 Test Kit lots were used. The LoQ was based on a $CV_{\text{within-lab}}$ criterion of less than 20 %.

Analytical sensitivity	PCT (ng/mL (μ g/L))	
	Plasma	Whole blood
LoB	0.033	0.033
LoD	0.059	0.072
LoQ	0.082	0.12
Concentration with $CV_{\text{within-lab}} = 10$ %	0.13	0.19

Sample type	Reportable range
Plasma	0.082-100 ng/mL (μ g/L)
Whole blood	0.12-100 ng/mL (μ g/L)

Reference ranges

Lithium-heparin whole blood and plasma samples obtained from 252 apparently healthy individuals (128 women and 124 men, aged 18 years or older) were analyzed with the PCT assay.

The 95th percentile was determined to be <0.12 ng/mL (μ g/L) for whole blood and 0.087 ng/mL (μ g/L) for plasma for both men and women.

NOTE: These values should only be used as examples. Each laboratory should establish its own reference ranges.

Imprecision

Within-run, within-lab and between-lot imprecision on plasma were determined by analyzing 5 levels of frozen native or antigen-spiked lithium-heparin plasma pools at 2 test sites. Four analyzers were used across sites and 3 Test Kit lots were used at each site. Each level was analyzed in duplicate in 2 separate runs per day for a minimum of 20 test days.

Plasma pool	Mean concentration of PCT (ng/mL (µg/L))	Within-run		Within-lab		N*
		SD (ng/mL (µg/L))	CV (%)	SD (ng/mL (µg/L))	CV (%)	
Both sites combined	0.092	0.0110	11.9	0.0115	12.5	323
	0.30	0.010	3.5	0.011	3.8	324
	0.52	0.015	3.0	0.017	3.3	324
	2.0	0.04	1.8	0.04	2.0	323
	67	1.0	1.5	1.3	1.9	324
Site 1	0.094	0.0093	9.9	0.0102	10.9	164
	0.30	0.011	3.6	0.012	3.9	164
	0.52	0.015	2.9	0.018	3.5	164
	2.0	0.04	1.9	0.04	2.2	163
	68	1.1	1.7	1.5	2.1	164
Site 2	0.089	0.0124	13.9	0.0127	14.2	159
	0.29	0.010	3.5	0.011	3.7	160
	0.52	0.016	3.0	0.016	3.1	160
	2.0	0.03	1.7	0.04	1.8	160
	66	0.9	1.4	1.1	1.7	160

*N = number of measurements

Plasma pool	Mean concentration of PCT (ng/mL (µg/L))	Between lots		No. of measurements
		SD (ng/mL (µg/L))	CV (%)	
Both sites combined	0.092	0.0045	4.9	323
	0.30	0.006	2.0	324
	0.52	0.013	2.5	324
	2.0	0.03	1.3	323
	67	1.3	1.9	324

Within-run and within-lab imprecision on whole-blood were determined by analyzing 5 levels of antigen-spiked lithium-heparin whole blood samples at 2 test sites. Four analyzers and 3 Test Kit lots were used across sites. Each level was analyzed in 5 replicates in 5 separate runs within a 3-hour period.

Whole blood	Mean concentration of PCT (ng/mL (µg/L))	Within-run		Within-lab		N*
		SD (ng/mL (µg/L))	CV (%)	SD (ng/mL (µg/L))	CV (%)	
Site 1	0.32	0.019	5.8	0.0188	5.8	50
	0.51	0.020	3.9	0.020	4.0	50
	1.9	0.07	3.4	0.07	3.6	50
	7.5	0.14	1.9	0.14	1.9	50
	55	1.3	2.4	1.3	2.4	50
Site 2	0.29	0.015	5.2	0.0148	5.2	50
	0.52	0.022	4.2	0.022	4.2	50
	1.9	0.04	2.1	0.04	2.3	50
	4.1	0.10	2.5	0.10	2.6	50
	41	0.6	1.5	0.6	1.6	50

*N = number of measurements

Hook effect

No hook effect was found when PCT concentrations up to ~1600 ng/mL (µg/L) were measured.

Carry over

Carry over from a sample with high PCT value (~1600 ng/mL (µg/L)) to an adjacent negative sample was determined to be <0.25 ng/mL (µg/L).

Interfering substances

Blockers have been added to this assay to minimize potential interference from heterophilic antibodies, for example HAMA (human anti-mouse antibodies) and rheumatoid factor (RF). Assessment indicated that the blockers have the desired effect.

Potential interfering substances (at concentrations listed in the table) were found to have no notable effect on the PCT assay (interference ≤10 %). Interference was tested with plasma samples containing approximately 0.5 ng/mL (µg/L) and 2 ng/mL (µg/L) of PCT spiked with the substance of interest at concentrations approximately 3 times the upper limit of the therapeutic range.

Drugs and solvents	Tested concentration
Acetaminophen	200 mg/L
Acetylsalicylic acid	652 mg/L

Drugs and solvents	Tested concentration
Ampicillin	53 mg/L
Beclomethasone	0.64 mg/L
Budesonide	533 mg/L
Cefotaxim	180 µg/mL
Ciclesonide	0.64 mg/L
Cyclosporine	6 mg/L
Diclofenac	50 mg/L
Dobutamine	1 µg/mL
Dopamine	0.9 mg/L
Erythromycin	60 mg/L
Flunisolide	0.464 mg/L
Fluticasone	0.88 mg/L
Furosemide	60 mg/L
Heparin, low molecular weight	3000 IU/L
Ibuprofen	500 mg/L
Imipenem	600 µg/mL
Metronidazole	120 mg/L
Mometasone	0.88 mg/L
Nitrofurantoin	4 mg/L
Noradrenaline	2 µg/mL
Oxytetracycline	12 mg/L
Tetracycline	15 mg/L
Triamcinolone	0.048 mg/L
Vancomycin	100 mg/L

The following substance was found to have a notable effect on the PCT assay at the concentrations shown in the table (interference >10 %).

Interfering substance	Interfering concentrations		Therapeutic concentration (2 h after an oral dose) [6]
	PCT 0.5 ng/mL (µg/L)	PCT 2.0 ng/mL (µg/L)	
Trimethoprim	5.9 mg/L	10.6 mg/L	3-5 mg/L

Endogenous blood components (at concentrations listed in the table) were found to have no notable effect on the PCT assay (interference ≤10 %). Interference was tested with plasma and whole blood samples containing approx. 0.5 ng/mL (µg/L) and 2 ng/mL (µg/L) of PCT spiked with the substance of interest at the concentration shown in the table. The concentrations represent increases to the native level.

Substances in blood	Tested concentration
Bilirubin unconjugated	200 mg/L
Bilirubin conjugated	330 mg/L
Biotin	3 µg/L
Cholesterol	5 g/L
Fibrinogen	10 g/L
Hemoglobin	2 g/L
Human albumin	60 g/L
Triglycerides	33 g/L

No sign of interference for hemolysis up to 300 mg/dL. For high degrees of hemolysis, that is >300 mg/dL, a maximum interference of –11 % was observed at PCT level 0.5 ng/mL.

Method comparison

The PCT assay (y) for the AQT90 FLEX analyzer was compared to a B·R·A·H·M·S PCT sensitive KRYPTOR assay for the B·R·A·H·M·S KRYPTOR System (x).

EDTA and lithium-heparin whole-blood and plasma samples were analyzed on the AQT90 FLEX analyzer and plasma samples on the B·R·A·H·M·S KRYPTOR System. Samples in the range 0.085-96 ng/mL (µg/L) for plasma and 0.12-87 ng/mL (µg/L) for whole blood were used (as measured on the AQT90 FLEX analyzer).

The Passing-Bablok regression equation for whole blood (versus plasma on the B·R·A·H·M·S KRYPTOR System) was found to be:

$$y = 1.05x + 0.007 \quad (n = 199)$$

The Passing-Bablok regression equation for plasma (versus plasma on the B·R·A·H·M·S KRYPTOR System) was found to be:

$$y = 1.02x - 0.007 \quad (n = 211)$$

Patents

Radiometer products may be covered by one or more patents and patent applications, see www.radiometer.com/en/legal/patents.

Explanation of graphical symbols

See the *AQT90 FLEX instructions for use* document.

References

See 'References' on the last page of this insert.

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1. Becker KL, Snider R, Nylen ES. Procalcitonin assay in systemic inflammation, infection, and sepsis: Clinical utility and limitations. *Crit Care Med* 2008 Mar; 36(3): 941-52.
2. McGee KA, Bauman NA. Procalcitonin, Clinical utility in Diagnosing Sepsis. *Clin Lab News* 2009 July; 35(7): 1-8.
3. O'Grady NP, Barie PS, Bartlett JG, Bleck T, Carroll K, Kalil AC *et al*. Guidelines for evaluation of new fever in critically ill adult patients: 2008 update from the American College of Critical Care Medicine and the Infectious Diseases Society of America. *Crit Care Med* 2008 Apr; 36(4): 1330-49.
4. Chiesa C, Panero A, Rossi N, Stegagno M, De Giusti M, Osborn JF, Pacifico L. Reliability of procalcitonin concentrations for the diagnosis of sepsis in critically ill neonates. *Clin Infect Dis* 1998 Mar; 26(3): 664-72.
5. Uzzan B, Cohen R, Nicolas P, Cucherat M, Perret GY. Procalcitonin as a diagnostic test for sepsis in critically ill adults and after surgery or trauma: a systematic review and meta-analysis. *Crit Care Med* 2006; 34(7): 1996-2003.
6. Stevens RC, Laizure SC, Sanders PL, Stein DS. Multiple-Dose Pharmacokinetics of 12 Milligrams of Trimethoprim and 60 Milligrams of Sulfamethoxazole per Kilogram of Body Weight per Day in Healthy Volunteers. *Antimicrob Agents Chemother* 1993 Mar; 37(3): 448-52.

Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Denmark
www.radiometer.com

Code number: 996-391
Version: 201504C