DIAPATH

Constant Innovation in Anatomic Pathology

INSTRUCTION FOR USE

C030X

EN

IT

FR

ES

DE

Mayer's Haematoxylin

Staining solution

EN

Manufacturer

Diapath S.p.A.

Via Savoldini, 71 - 24057 Martinengo (BG) - Italy

Tel. +39 0363 986411

Fax +39 0363 948000

www.diapath.com

info@diapath.com

EMDN code: W01030708

Code	Format	UDI code
C0301	125 ML	08033895006903
C0302	500 ML	08033895006910
C0303	1 LT	08033895003544
C0304	2 LT	08033895006934
C0305	2.5 LT	08033895003568
C0306	5 LT	08033895006958

Intended use

Disposable, non-sterile in vitro diagnostic device for clinical use for nuclear staining.

The device, for diagnostic use or research, is intended to be used by qualified personnel in laboratories who prepare cytological and histological samples.

In terms of staining, the result obtained from the combination of the different components can be intended for inspection under a microscope and interpretation of the result by experienced and qualified personnel.

The device for manual or automatic use associated with an automatic colouring machine.

The device is intended to be used on preparations previously prepared according to the methods set forth by the standard sample processing technique used by the laboratory.

The device is disposable, not recyclable.

Description of the reagents and any limitations of use

Haematoxylin	CAS No. 517-28-2	EC No. 208-237-3	C.I. 75290
Potassium aluminium sulfate	CAS No. 7784-24-9	EC No. 215-691-6	
Deionised water	CAS No. 7732-18-5	EC No. 231-791-2	
Sodium iodate	CAS No. 7681-55-2	EC No. 231-672-5	
Chloral hydrate	CAS No. 302-17-0	EC No. 206-117-5	
Citric acid	CAS No. 77-92-9	EC No. 201-069-1	

Device	Stability when the product is closed	Stability after opening	Ready to use	Storage temperature
C0301	2 years	2 years	YES	15 – 30°C
C0302	2 years	2 years	YES	15 – 30°C
C0303	2 years	2 years	YES	15 – 30°C
C0304	2 years	2 years	YES	15 – 30°C
C0305	2 years	2 years	YES	15 – 30°C
C0306	2 years	2 years	YES	15 – 30°C

Diapath S.p.A.

Via Pietro Savoldini, 71 - 24057 Martinengo (BG) Italy

T +39 0363 986411 | F +39 0363 948000

info@diapath.com | diapath.com

Rev.01 14/06/2024

Information concerning interfering substances that may affect the performance of the device

Not applicable

List of materials provided and of specific materials that are necessary but not provided

The materials required for the analysis, but not provided, are:

- Distilled water
- Running water
- Reagents for tissue dehydration (xylene and substitutes, graded alcohols)
- Eosin
- Differentiating solutions and toning solutions, as needed
- Mounting medium
- Glass coverslip or replacement (at user's discretion)
- Positive control slides

Additional equipment, if provided for in the protocol of use:

- Automatic colouring machine

Information on the combined use of the device

The device is intended for use within a Haematoxylin and eosin (H&E) staining protocol, in combination with specific reagents, previously listed in the section on materials that are necessary but not supplied.

Principle of the Test

Haematoxylin, as an alkaline stain with positive charge, stains the negatively charged (basophilic) cell components blue-purple, as they are nucleic acids contained in the nucleus.

Instructions for use

Examples of staining protocols for histological sections

Reagent	Time
Xylene	03:00
Xylene	03:00
Xylene	03:00
Ethanol 99	02:00
Ethanol 99	02:00
Ethanol 95	02:00
Running water	01:00
Mayer's Haematoxylin	05:00 – 10:00
Running water	05:00
Ethanol 70	01:00
Alcoholic eosin 0.5%	03:00
Ethanol 95	01:00
Ethanol 99	02:00
Ethanol 99	02:00
Xylene	03:00
Xylene	03:00
Xylene	03:00

Reagent	Time
Xylene	03:00
Xylene	03:00
Ethanol 99	02:00
Ethanol 99	02:00
Ethanol 95	02:00
Running water	01:00
Mayer's Haematoxylin	05:00 – 10:00
Running water	05:00
Aqueous Eosin	01:00 - 04:00
Running water	00:30 – 01:00
Ethanol 95	01:00
Ethanol 99	02:00
Ethanol 99	02:00
Xylene	03:00
Xylene	03:00
Xylene	03:00

It is good practice to follow all the manufacturer's instructions in order to obtain the desired result, as an incorrect procedure may hinder the quality of the staining.

All operational procedures, prior to staining, i.e. fixation, dehydration, clearing, paraffin impregnation and finally cutting, are supervised by the operator and must be validated according to normal routine operations within the reference laboratory.

It is neither self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel, since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified and trained personnel.

refer to the Bibliography section

Interpretation of results

Nuclei:	Blue-Purple
Cytoplasm, connective tissue, intercellular substances:	Pink-Red

Quality control

In order to ensure the correct performance of haematoxylin, before using it routinely, it is essential to carry out a quality control, including tissues fixed and treated according to the laboratory's operating procedures.

Analytical performance characteristics

Not applicable

Conditions for sample collection, handling and preparation

Preparation of the histological and cytological sample follows the activities expected in the normal workflow of a pathology laboratory.

Do not handle until all safety precautions have been read and understood.

Storage

Store the device at a temperature between +15°C and +30°C.

Do not expose the device to environmental conditions other than those prescribed, as failure to comply with storage conditions would lead to a change in the chemical and physical characteristics of the device, resulting in underperforming staining.

If properly stored and in unopened packaging, the product is stable as stated in the section *Description of the reagents and any limitations of use*. Do not use after the expiry date.

If the product is not stored as instructed, its yield may vary and it will have to be tested by the operator.

Stability when handling the device

The use of the device is under the absolute control of the skilled and qualified operator.

The environmental conditions in the laboratory must be controlled:

- Temperature: +15°C /+30°C
- Relative humidity ≤80%

Instructions for disposal / Safe final disposal of the device

Specific instructions on how to safely dispose of the device after use are to be found in the local regulations of the country in which the device is used.

It is good practice to dispose of the device as special waste in dedicated containers and in accordance with internal disposal procedures.

For information on the content of the reagents, please refer to the safety data sheet. If the product has expired and/or is unused, dispose of it in accordance with local waste regulations, based on the hazard classification on the label, after checking for any contamination. In some cases, analytical evaluation may be required so as to correctly classify the waste and determine the hazard characteristics to be attributed to it

Warnings, precautions and restrictions when using the device

The device is to be used only for its intended use.

For proper handling, it is recommended to:

- Wear personal protective equipment such as a lab coat, gloves and goggles. Protect your eyes and face
- The operator may come into contact with an irritating device
- Work with the device under a fume hood
- Avoid inhaling the fumes of vapours while handling the device
- In order to obtain the expected result, it is advisable to follow the indications given in the protocol of use
- Keep the work area clean in order to avoid compromising the diagnostic results

Device malfunction or degradation

The device retains its properties throughout its estimated useful life if the storage requirements are guaranteed and complied with.

The device does not generate risks from any debris due to wear and tear, degradation products and processing residues.

In case of malfunction, isolate the device, replace it with a new one and dispose of the non-functioning device. Should this occur, please contact the manufacturer by sending an e-mail to claims@diapath.com

Error	Possible cause	Possible corrective action
Loss of staining specificity	Failure to observe the instructions on correct preservation of the device after opening or unsuitable storage conditions	Consult the storage instructions and temperatures beforehand for suitable maintenance of the device, also provided on the label
Unsuitable handling during use	Underperforming device	Staining protocol provided on the IFU given to personnel

Warnings from the user, with the need to report any serious incidents that have occurred involving the device

In case of a serious accident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Legend of symbols used in labelling and in the instructions for use

Pictogram	Pictogram name	Explanatory notes
	Manufacturer	This indicates the manufacturer of the device
	Expiry date	This indicates the date after which the device must not be used
	Batch number	This indicates the manufacturer's batch required for batch identification
	Catalogue code	This indicates the code assigned by the manufacturer for identification
	Distributor	This indicates the company that distributes the medical device
	Non-sterile	This indicates that the device has not undergone a sterilisation process
	Keep away from light	This indicates that the device needs to be protected from light sources.
	Temperature limit	This indicates the temperature limit that the medical device can be safely exposed to
	Do not re-use	This indicates a device that is to be used once only, disposable
	Consult instructions for use	This indicates that the user needs to consult the instructions
	Attention	This indicates that the user needs to be careful when operating the device or that the situation requires greater awareness in order to avoid unwanted consequences
	In vitro diagnostic medical device	This indicates a device that is intended to be used as an in vitro diagnostic medical device
	CE Marking	The CE mark certifies that the product complies with the applicable EU directives

	Sufficient content for <n> tests	This indicates the number of tests that can be performed with the device
	Date of manufacture	This indicates the date of manufacture of the medical device
	Unique Device Identification	This indicates a carrier containing the fixed and variable information related to the UDI

Substance/mixture and hazard identification

Substance or mixture	Hazard pictograms	Indication of the hazard	Cautionary advice
MAYER'S HAEMATOXYLIN		The product is not considered hazardous according to EC Regulation 1272/2008 (CLP).	

References

- "Laboratory methods in histotechnology. Stain methods of the Arm Forces Institute of Pathology". Prophet, Mills, Arrinton, Sobin. Editor McGRAW Hil. Washington D.C. 1968
- "Tecnica di anatomia patologica". M.Melis, F. Carpino, U. Di Tondo. Edi Ermes. 1989
- "Theory and practice of Histological Techniques". by Bancroft, J.D. and Gamble, M., Publisher Churchill Livingstone. New York. 2002

Table of revisions

Rev.	Date	Purpose of the revision
1st issue	03/03/2023	First issue
01	14/06/2024	Update UDI-DI code and staining protocols in the "instructions for use" section

(*) If the request is due to a change, please provide a detailed description of the reason for revising the document and the changes that have been made since the last revision; if available, please enclose recalls, non-conformities, NC AC AP statistics, etc.

You can find this document in the reserved area of the www.diapath.com website. If you want to download the document and keep up to date on the latest available version, you will need to log in after registering.

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