QUALIRIS QC FM-FDP

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
QUALIRIS QC FM-FDP

Control plasmas for external quality assessment of fibrin monomers and fibrin/fibrinogen degradation products assays

- Kit containing 6 x 1-ml Vials

(REF 01188)

May 2018

English

1/ INTENDED USE

The **Qualiris QC FM-FDP** kit provides unassayed plasmas intended for external quality assessment of fibrin monomers and fibrin/fibrinogen degradation products assays used with the Qualiris by Stago programme.

2/ KIT REAGENTS

Kit of 6 vials of lyophilized human plasmas numbered from 1 to 6.

WARNING - POTENTIAL BIOHAZARDOUS MATERIAL

These reagents contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.

3/ PRECAUTIONS

Store at 2-8 °C. For *in vitro* use only. These reagents are to be used by certified medical laboratory personnel only. The disposal of waste materials must be carried out according to current local regulations. Exercise great care in the handling of these reagents and of patient samples.

4/ REAGENTS PREPARATION

Reconstitute each vial with exactly 1 ml of distilled water. Allow the solution to stand at room temperature (18-25 °C) for 30 minutes. Then, swirl the vial gently before use.

5/ REAGENTS STABILITY

The reagents in intact vials are stable until the expiration date indicated on the box label, when stored at 2-8 °C.

Once reconstituted, reagents remain stable for 4 hours at 20 ± 5 °C.

Do not freeze.

6/ REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Reagents for performing tests mentioned in Section 1.
- Common clinical laboratory equipment and materials (distilled water...).

7/ DIRECTIONS FOR USE

The **Qualiris QC FM-FDP** plasmas are to be used in the same manner as patients' plasmas.

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Significant changes are indicated by dotted lines in the margin.

DIAGNOSTICA STAGO S.A.S.
3 allée Thérèse
92600 Asnières sur Seine (France)
+33 (0)1 46 88 20 20
webmaster@stago.com

Information and/or pictures contained in this document are protected by copyrights and other intellectual property rights. © 2018, Diagnostica Stago, all rights reserved. Diagnostica Stago's logos and products names are registered trademarks.