30 years already!

Over thirty years ago, SATELEC pioneered the piezo-ceramic ultrasonic scaler. As a result of collaborating with dental professionals worldwide, we are proud to offer today the most advanced ultrasonic scaling devices and a variety of more than 70 inserts encompassing all ultrasonic applications (including endodontic and apical surgery procedures). All our scalers work on the piezo principle. Piezo means smooth, quiet and more comfortable vibration. By utilizing only one handpiece, the powerful and compact **Suprasson® P5 Booster** removes deposits with increased patient comfort and reduced operator fatigue - all in less time.

Daily hygiene procedures are often strenuous and time consuming. The piezo scaling technique is like using a hand instrument, but without the effort!

Incorporating the routine use of ultrasonics is the backbone of periodontal treatment. Hygiene maintenance is deliverable in less time, maximizing production and reducing the number of patient appointments.

Technical specifications

51.1.

Supply voltage:	110 V- or 220/230 V- 50/60 Hz
Equipment classification:	Class II - BF type
Operating conditions:	Intermittent duty 5min / 10min
Max. power consumption:	40 VA
Vibration frequency:	27 to 32 kHz
Water pressure:	1 to 5 bar (15 to 73 psi)
Overall dimensions:	Height 60 mm
.	Width 155 mm
.	Depth 150 mm
Weight:	1.2 kg

References

Autoclavable handpiece:	F12200
Universal wrench:	F00406
Easy Torque M:	F01400
Easy Torque C:	F01405



Manufactured in compliance with regulations in force 93/42/EEC.

Scaling tips

These four supra- and sub-gingival scaling tips are ideal for treating commonly encountered cases (pockets <3-4mm). All tips come with the new Spray Control System (Satelec exclusive).

The tip must be oriented tangentially to the surface being treated and must be moved in a back-and-forth sweeping motion, over the whole surface, without excessive lateral pressure. Removing thick deposits may require using the Boost mode for short periods. (19)

POWER SETTINGS

		P-Max Newtron Prophy Max Newtron	P5 Newtron	SP Newtron	Suprasson P5 Booster
51.2.	No. 1 tip	8 to 10	14 to 15	7 to 8	7 to 9
	No. 2 tip	8 to 10	14 to 15	7 to 8	7 to 10
	No. 10X tip	6 to 8	12 to 14	6 to 7	7 to 8
51.4.	No. 10Z tip	6 to 8	12 to 14	6 to 7	7 to 8

51.2.

No. 1 tip: "Universal". Recommended for treating simple cases and gross supra-gingival scaling. Tangential orientation with respect to the tooth surface. To-and-fro sweeping motion to "push" and "pull" calculus without damaging the enamel.

No. 2 tip: "Voluminous calculus". Recommended for removing voluminous supra-gingival deposits. Apply flat end to surface of teeth.


No. 10X tip: "Interproximal". Recommended for the treatment of interproximal spaces and for supra-gingival scaling. Anatomic shape allows quick, efficient hand movements.

51.4.

No. 10Z tip: "Sub-gingival". New design for the 10Z for sub-gingival scaling, this instrument is recommended for scaling medium pockets (<4mm). The new version is more powerful as its amplitude has been increased by 75% and the working edges are longer.

Color Coding System - Newtron Technology
Identification of power modes by color coding.

51.2.
51.4.



Scaling tips

The tip must be oriented tangentially to the surface being treated, parallel to the main axis of the tooth, and must be moved in a back-and-forth sweeping motion without excessive lateral pressure. Begin the treatment at the center of the dental face. Push or pull the calculus towards the mesial and distal margins.

The Twiny® is a patented tip with remarkable versatility. It greatly reduces scaling time and removes pathogen toxins without damaging the cementum. The Twiny is recommended for both supra- and sub-gingival scaling. (14)

Thick deposits may require using the Boost mode for short periods.

POWER SETTINGS

		P-Max Newtron Prophy Max Newtron	P5 Newtron	SP Newtron	Suprasson P5 Booster
51.3.	No. 1-S tip	8 to 10	14 to 15	7 to 8	7 to 10
	No. 3 tip	8 to 10	14 to 15	7 to 8	7 to 8
	No. 10P tip	8 to 10	14 to 15	7 to 8	7 to 8
	Twiny tip	8 to 10	14 to 15	7 to 8	7 to 10

51.3.

No. 1-5 tip: "Slim". This tip was developed for supra- and sub-gingival scaling. It has improved load resistance, and superior amplitude and power compared to the original tip. The more effective lateral edges make it particularly suitable for scaling the interproximal spaces.

No. 3 tip: "Stains". For removing discolorations and stains (tobacco, tea, coffee, etc.). The tip is used by placing the rounded end in contact with the surface to be treated.

No. 10P tip: "Shallow pockets". Fine tip designed for scaling shallow pockets (<2-3mm). It provides greater irrigation than the other tips.

Twiny tip: "Double spray". Very powerful "surface" tip. The concave, symmetrical surfaces amplify the cavitation effect. It can be used to remove large deposits, polish treated surfaces, and irrigate pockets thanks to its double spray (pockets <3mm).

Color Coding System - Newtron Technology
Identification of power modes by color coding.

51.3.

1-10

Perio mini-tips

These mini-tips make it possible to perform a full-mouth periodontal debridement treatment. All the tips have a universal curette shape. The distal surface of the last segment is at a 90° angle to the shaft. The end of the mini-tip has a quarter-loop to prevent blunt trauma. The double guiding edges remove calculus and deposits with no risk of damage to connective fibers. (8, 9, 10)



H4 tips are designed for sub-gingival scaling, non-surgical treatment of deep, narrow pockets (>4-6mm), and open treatments (flap surgery). Low pressure combined with a very low amplitude provide tactile sensitivity unequalled by any other ultrasonic tip.

The mini-tips operate in the opposite direction to manual curettes (ie. using a push stroke without scraping). They should be moved from the sulcular opening towards the bottom of the pocket without excessive lateral pressure (0.3 to 0.5N). The calculus is removed by methodically crisscrossing all root surfaces, while the biofilm is mechanically disrupted and eliminated by irrigation. Residual deposits are removed from pockets and cementum surfaces are decontaminated by the bacteriostatic effect of the micro-cavitation. (5, 8, 9)

POWER SETTINGS

		P-Max Newtron Prophy Max Newtron	P5 Newtron	SP Newtron	Suprasson P5 Booster
51.5.	H3 tip	1 to 5	1 to 2	1 to 2	1 to 4
	H4L tip	1 to 5	1 to 2	1 to 2	1 to 4
	H4R tip	1 to 5	1 to 2	1 to 2	1 to 4

51.5.

H3 tip: "Anterior teeth". Ideal tip to begin with, for the treatment of anterior teeth. The guiding edge is used parallel to the pocket.

H4L tip: "Premolars and molars". First tip in the sequence. Left angled. Recommended for the treatments of all surfaces and furcations:

- **Mandibular:** The lingual and mesial surfaces of sector 3, then the buccal and distal surfaces of sector 4. Pivoted at 43.
- **Maxillary:** The buccal and distal surfaces of sector 1, then the palatal and mesial surfaces of sector 2. Pivoted at 13.

H4R tip: "Premolars and molars". Second tip in the whole-mouth sequence. Right angled. Follows the H4L tip in the treatment sequence. Recommended for the treatment of all surfaces and furcations:


- **Mandibular:** Buccal and distal surfaces of sector 3, then lingual and mesial surfaces of sector 4. Pivoted at 43.
- **Maxillary:** The buccal and distal surfaces of sector 1, then the palatal and mesial surfaces of sector 2. Pivoted at 13.

Tip equivalents

- **H3:** Anterior teeth and premolars. Replaces manual curettes nos. 1-2, 3-4, and 5-6.
- **H4L and H4R:** Premolars and molars. Replace manual curettes nos. 7-8, 9-10, 11-12, and 13-14.

51.5.

Color Coding System - Newtron Technology
Identification of power modes by color coding.

Tips

Mastering ultrasound does not stop at devising generators and control software, it also means designing the appropriate accessories that will give them optimum vibration power and guarantee their reliability.


Our instruments are made from the best alloys, exclusive to **Satelec**, using patented procedures.

We offer you a range of different coatings (diamond, thermal treatments) or materials (carbon) so you can choose the most suitable and efficient instrument for each material, giving you better performance.

Satelec's expertise in tuning all its ultrasonic instruments is second to none. Our transducers are probably the best in the world and certainly the only ultrasonic motors capable of driving all these tips at maximum performance, whatever their type or composition.

Advantages

- **Corrosion resistant:** the steel alloy we have chosen does not corrode and is resistant to decontamination products as well as to irrigation solutions. **51.2.**
51.3.
51.4.
51.5.
- **Tensile strength:** the manufacturing procedures are designed, and the tensile properties adjusted, to suit the treatments for which they are intended.
- **Sterilization:** the tips are suitable for autoclave sterilization, to comply with ISO-11134, as are the handpieces and accessories (wrenches, storage kits, etc.): **51.2.**
51.3.
51.4.
51.5.
Temperature: 134°C;
Pressure: 2 bars (29 P.S.I.);
Sterilization time: 18 minutes.

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Medical devices and electromedical devices for dental art, electrosurgery and bone and bone substitute surgery equipments, with associated sterile and non-sterile accessories

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LNE/G-MED certifies that, on the basis of the results contained in the file referenced P159650-1, the quality system - for design, manufacturing, and final inspection - of medical devices listed here aboved complies with the requirements of the Directive 93/42/EEC, annex II excluding section 4.

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On behalf of the Certification Director
Cécile VAUGELADE
G-MED Certification Division Manager

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LNE/G-MED
0459

 On behalf of the Certification Director
Cécile VAUGELADE
 G-MED Certification Division Manager

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Periodontal debridement and root planning

The micro tips H3 and H4 with curette have shape double guiding edges to remove calculus and deposits with no risk of damage to connective fibers. They are designed for deep and narrow pockets debridement.

The diamond-coated micro-tips H1 and H2 are used in the second step for root planning, treatment of furcations and abscesses and to remove granulation tissues. They must be inserted delicately, without pressure (approximately 0.3N), to avoid over-instrumentation.

The micro-tips operate in the opposite direction to manual curettes (ie. using a push stroke without scraping), from the sulcular opening towards the bottom of the pocket without excessive lateral pressure. Residual deposits are removed from pockets and cementum surfaces are decontaminated by the bacteriostatic effect of micro-cavitation.

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


Scaling tips from SATELEC are made of an exclusive alloy that respects the enamel.

To be effective, the tips must be oriented tangentially to the surface being treated and must be moved in a back-and-forth sweeping motion without excessive lateral pressure, for more patient comfort.

Scaling tips are 20% to 50% faster to remove dental plaque compared with manual instruments and cause less damage to tooth surfaces.

Precise control of power and irrigation is a guarantee of better efficiency.

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